

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17

DRAFT ICCVAM RECOMMENDATIONS
***In Vitro* Acute Toxicity Test Methods**

March 17, 2006

17
18
19
20
21
22
23
24
25
26
27
28
29
30

[This Page Intentionally Left Blank]

30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45

TABLE OF CONTENTS

Page Number

Section 1.0 Draft ICCVAM Recommendations for *In Vitro* Acute Toxicity Test Methods..... 5

Appendix A Draft Performance Standards for *In Vitro* Acute Toxicity Methods..... A-3

Appendix B Draft Recommended Test Method Protocols B-1

Appendix B1 Test Method Protocol for the BALB/c 3T3 Neutral Red Uptake (NRU) Cytotoxicity Test B1-1

Appendix B2 Test Method Protocol for the Normal Human Epidermal Keratinocyte (NHK) Neutral Red Uptake (NRU) Cytotoxicity Test B2-1

45

46

47

48

49

50

51

52

53

54

55

56

57

58

[This Page Intentionally Left Blank]

59 **1.0 Draft ICCVAM Recommendations for *In Vitro* Acute Toxicity Test Methods**

60

61 **1.1 Draft Recommended Test Method Uses**

62

63 This independent validation study evaluated two *in vitro* neutral red uptake (NRU) basal
64 cytotoxicity assays: the BALB/c 3T3 (clone A31) mouse fibroblast NRU (hereafter referred
65 to as the 3T3 NRU) test method and the normal human keratinocyte NRU (hereafter referred
66 to as NHK NRU) test method. The objective of the study was to determine their ability to
67 estimate rodent acute oral toxicity LD₅₀ values to be used in a weight-of-evidence approach
68 to set the starting dose for *in vivo* acute oral toxicity tests. Based on the results of this
69 validation study, ICCVAM proposes the following draft recommendations:

- 70 1. The 3T3 and NHK NRU test methods are not sufficiently accurate to predict the acute
71 oral toxicity of substances for the purposes of hazard classification (see Section 6 of
72 the *In Vitro* Acute Toxicity Test Methods Background Review Document).
- 73 2. For the purposes of acute oral toxicity testing, the 3T3 and NHK NRU test methods
74 may be used in a weight-of-evidence approach to determine the starting dose for the
75 current acute oral *in vivo* toxicity protocols (i.e., the Up-and-Down Procedure [UDP]
76 and Acute Toxic Class [ATC])
- 77 3. Consistent with the U.S. Government Principles on the Use of Animals in Research,
78 Testing, and Education (National Research Council 1996), and the U.S. Public Health
79 Service Policy on Humane Care and Use of Laboratory Animals (PHS 2002)¹, *in vitro*
80 basal cytotoxicity test methods as part of a weight-of-evidence approach to estimate
81 the starting dose for acute oral *in vivo* toxicity test methods should be considered and
82 used where appropriate before testing is conducted using animals. For some types of
83 substances, this approach will reduce the number of animals needed. In some testing
84 situations, the approach may also reduce the numbers of animals that die or need to be
85 humanely killed.

¹ National Research Council. 1996. Guide for the Care and Use of Laboratory Animals. Washington, DC: National Academy Press.

PHS. 2002. Public Health Service (PHS) Policy on Humane Care and Use of Laboratory Animals.

- 86 4. Substances with specific toxic mechanisms that are not expected to be active in 3T3
87 or NHK cells (e.g., those that are neurotoxic, cardiotoxic, interfere with energy
88 utilization, or alkylate proteins and other macromolecules) will likely be
89 underpredicted by these *in vitro* basal cytotoxicity test methods. Therefore, until such
90 time as a more predictive testing approach is developed, the results from basal
91 cytotoxicity testing with such substances may not be appropriate.
- 92 5. The regression formula used to determine starting doses should be the revised
93 Registry of Cytotoxicity (RC) regression line [with IC₅₀ values in µg/mL and LD₅₀
94 values in mg/kg] developed with the RC chemicals using rat LD₅₀ data only and
95 excluding chemicals with mechanisms of action that are not expected to be active in
96 *in vitro* basal cytotoxicity test methods.
- 97 6. The performance of other *in vitro* basal cytotoxicity test methods that are based on
98 similar scientific principles and that measure or predict the same biological response
99 (i.e., basal cytotoxicity and the rat acute oral LD₅₀ value, respectively) should be
100 demonstrated to meet or exceed the accuracy and reliability of the 3T3 and NHK
101 NRU test methods.
- 102 7. Compared to the NHK NRU test method, the 3T3 NRU test method appears to be less
103 labor intensive and less expensive to conduct; therefore, the 3T3 NRU cytotoxicity
104 test method is recommended for general use.

105

106

107 1.2 Draft Recommended Future Studies

108

- 109 1. Additional data should be collected using the 3T3 and/or the NHK NRU test methods
110 to evaluate their usefulness for predicting the *in vivo* acute oral toxicity of chemical
111 mixtures.
- 112 2. Additional high quality comparative *in vitro* basal cytotoxicity data should be
113 collected in tandem with *in vivo* rat acute oral toxicity test results to supplement the
114 high quality validation database started by this study. Periodic evaluations of the
115 expanded database should be conducted to further characterize the usefulness and

- 116 limitations of using *in vitro* cytotoxicity data as part of a weight-of-evidence approach
117 to estimate starting doses.
- 118 3. Additional efforts should be conducted to identify additional *in vitro* tests and other
119 methods necessary to achieve accurate acute oral hazard classification; specifically,
120 studies should be conducted to investigate the potential use of *in vitro* cell-based test
121 methods that incorporate mechanisms of action and evaluations of ADME
122 (absorption, distribution, metabolism, excretion) to provide improved estimates of
123 acute toxicity hazard categories.
- 124 4. The *in vivo* database of reference substances used in this validation study should be
125 used to evaluate the utility of other nonanimal approaches to estimate starting doses
126 for acute oral systemic toxicity tests (e.g., widely available software that uses
127 quantitative structure-activity relationships [QSAR]).
- 128 5. Standardized procedures to collect information pertinent to an understanding of the
129 mechanisms of lethality should be included in future *in vivo* rat acute oral toxicity
130 studies. Such information will likely be necessary to support the further development
131 of predictive mechanism-based *in vitro* methods.
- 132 6. An expanded list of reference substances with estimated rat LD₅₀ values substantiated
133 by high quality *in vivo* data should be developed for use in future *in vitro* test method
134 development and validation studies.

135

136 **Appendix A** provides *Draft Performance Standards for In Vitro Acute Toxicity Methods* that
137 are based on ICCVAM guidelines (ICCVAM 2003²). **Appendix B** provides two draft
138 recommended *in vitro* NRU basal cytotoxicity protocols that are revised versions of the
139 Phase III protocols used in the validation study. **Appendix B-1** is the *Test Method Protocol*
140 *for the BALB/c 3T3 Neutral Red Uptake (NRU) Cytotoxicity Test*. **Appendix B-2** is the *Test*
141 *Method Protocol for the Normal Human Epidermal Keratinocyte (NHK) Neutral Red Uptake*
142 *(NRU) Cytotoxicity Test*.

² ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised, and Alternative Test Methods. NIH Publication No. 03-4508. National Institute for Environmental Health Sciences, Research Triangle Park, NC. Available: <http://iccvam.niehs.nih.gov/>. [accessed 2 June 2005].

143

144

145

146

147

148

149

150

151

Appendix A

152

153

154

Draft Performance Standards for *In Vitro*

155

Acute Toxicity Methods

156

156

157

158

159

160

161

162

163

164

165

166

167

168

169

170

[This Page Intentionally Left Blank]

171	TABLE OF CONTENTS		
172			Page Number
173	LIST OF TABLES		A-5
174	LIST OF ACRONYMS AND ABBREVIATIONS		A-6
175	INTERAGENCY COORDINATING COMMITTEE ON THE VALIDATION OF		
176	ALTERNATIVE METHODS (ICCVAM) - DESIGNATED AGENCY		
177	REPRESENTATIVES		A-8
178	INTERAGENCY COORDINATING COMMITTEE ON THE VALIDATION OF		
179	ALTERNATIVE METHODS		
180	(ICCVAM) ACUTE TOXICITY WORKING GROUP (ATWG)		A-9
181	NATIONAL TOXICOLOGY PROGRAM (NTP) INTERAGENCY CENTER FOR		
182	THE EVALUATION OF ALTERNATIVE TOXICOLOGICAL METHODS		
183	(NICEATM)		A-10
184	PREFACE		A-11
185	EXECUTIVE SUMMARY		A-13
186	1.0 PURPOSE AND BACKGROUND OF PERFORMANCE STANDARDS		A-15
187	1.1 Introduction		A-15
188	1.2 Elements of ICCVAM Performance Standards		A-16
189	1.3 ICCVAM Process for the Development of Performance Standards		A-17
190	1.4 ICCVAM Development of Recommended Performance Standards		
191	for <i>In Vitro</i> Acute Toxicity Test Methods		A-18
192	1.4.1 Current Regulatory Testing Requirements for Acute Systemic		
193	Toxicity		A-18
194	1.4.2 Test Methods for Assessing Acute Systemic Toxicity		A-19
195	1.4.3 Intended Regulatory Uses for <i>In Vitro</i> Cytotoxicity Test		
196	Methods.....		A-20
197	1.4.4 Similarities and Differences in the Endpoints of <i>In Vitro</i>		
198	Cytotoxicity Test Methods and <i>In Vivo</i> Acute Oral Toxicity Test		
199	Methods.....		A-21
200	2.0 <i>IN VITRO</i> ACUTE TOXICITY TEST METHODS		A-22
201	2.1 Background		A-22
202	2.2 Principles of <i>In Vitro</i> Basal Cytotoxicity Assays to Predict Starting		
203	Doses for Acute Oral Toxicity Tests		A-23
204	2.3 Essential Test Method Components for <i>In Vitro</i> Basal Cytotoxicity		
205	Assays to Predict Starting Doses for Acute Oral Toxicity (Lethality)		
206	Tests		A-24

207	2.3.1	<i>In Vitro</i> Cell Culture Conditions	A-25
208	2.3.2	Application of the Test Substances	A-25
209	2.3.3	Control Substances	A-27
210	2.3.4	Viability Measurements	A-28
211	2.3.5	Interpretation of Results.....	A-29
212	2.3.6	Test Report.....	A-30
213	2.4	Reference Substances for <i>In Vitro</i> Basal Cytotoxicity Assays to Predict	
214		Starting Doses for Acute Oral Toxicity Tests.....	A-32
215	2.5	Accuracy and Reliability.....	A-35
216	2.5.1	Accuracy	A-35
217	2.5.2	Reliability.....	A-37
218	3.0	REFERENCES	A-41
219			

219
220
221
222
223
224
225
226
227
228
229
230
231
232
233
234
235
236
237

LIST OF TABLES

Table 1-1	Summary of Current U.S. Legislation for Using Acute Systemic Toxicity Data for Product Labeling	A-19
Table 1-2	Regulatory Classification Systems for Acute Oral Toxicity	A-20
Table 2-1	Recommended Reference Substances for Evaluation of <i>In Vitro</i> Basal Cytotoxicity Methods for Predicting Starting Dose for Acute Systemic Toxicity Tests	A-34
Table 2-2	Linear Regression Analyses of <i>In Vitro</i> and <i>In Vivo</i> Results	A-36
Table 2-3	Prediction of GHS Toxicity Category by RC Rat-Only Weight Regression Excluding Chemicals with Specific Mechanisms of Toxicity.....	A-38
Table 2-4	Summary of CV Results for the 3T3 and NHK NRU Test Methods.....	A-40

237
238
239

LIST OF ACRONYMS AND ABBREVIATIONS

240	ANOVA	Analysis of Variance
241	ATC	Acute Toxicity Class
242	BRD	Background review document
243	CASRN	Chemical Abstracts Service Registry Number
244	CPSC	U.S. Consumer Product Safety Commission
245	CS	Calf serum
246	CV	Coefficient of Variation
247	°C	Degrees Celsius
248	DMSO	Dimethyl sulfoxide
249	DOT	U.S. Department of Transportation
250	DPBS	Dulbecco's Phosphate Buffered Saline
251	ECVAM	European Center for the Validation of Alternative Methods
252	EPA	U.S. Environmental Protection Agency
253	ETOH	Ethanol
254	FDA	U.S. Food and Drug Administration
255	FDP	Fixed Dose Procedure
256	FL	Fluorescein leakage
257	FR	Federal Register
258	GHS	Globally Harmonized System of Classification and Labelling of
259		Chemicals (UN 2005).
260	IC ₅₀	Test substance concentration producing 50% inhibition of the endpoint
261		measured
262	ICCVAM	Interagency Coordinating Committee on the Validation of Alternative
263		Methods
264	LD ₅₀	Lethal dose that produces lethality in 50% of test animals
265	LDH	Lactate dehydrogenase

266	MTT	[3-(4,5-dimethylthiazol-2yl)2,5-diphenyl tetrazolium bromide]
267	NCS	Newborn calf serum
268	NHK	Normal human epidermal keratinocytes
269	NICEATM	National Toxicology Program Center for the Evaluation of Alternative
270		Toxicological Methods
271	NIEHS	National Institute of Environmental Health Sciences
272	NIOSH	U.S. National Institute for Occupational Safety and Health
273	NR	Neutral red
274	NRR	Neutral red release
275	NRU	Neutral red uptake
276	NTP	U.S. National Toxicology Program
277	OD	Optical density
278	OECD	Organisation for Economic Cooperation and Development
279	OSHA	U.S. Occupational Safety & Hazards Administration
280	RC	Registry of Cytotoxicity
281	RTECS	Registry of Toxic Effects for Chemical Substances
282	3T3	BALB/c mouse fibroblasts, clone A31
283	UDP	Up and Down Procedure
284	UN	United Nations
285	XTT	[Sodium 3,3,-[(Phenylamino)carbonyl]-3,4-Tetrazolium-Bis(4-
286		methoxy-6-nitro)benzenesulfonic acid hydrate]
287	ZEBET	German Center for Documentation and Evaluation of Alternative
288		Methods to Animal Experiments
289		
290		

290
291
292
293

**INTERAGENCY COORDINATING COMMITTEE ON THE
VALIDATION OF ALTERNATIVE METHODS (ICCVAM)
Designated Agency Representatives**

Agency for Toxic Substances and Disease Registry

- Moiz Mumtaz, Ph.D.

Consumer Product Safety Commission

- Marilyn L. Wind, Ph.D. (Vice-Chair)
- * Kailash C. Gupta, D.V.M., Ph.D.
- * Patricia Bittner, M.S.
- * Kristina Hatelid, Ph.D.

Department of Agriculture

- Jodie Kulpa-Eddy, D.V.M.
- ◇ Elizabeth Goldentyer, D.V.M.

Department of Defense

- Robert E. Foster, Ph.D.
- ◇ Patty Decot
- * Harry Salem, Ph.D.
- * John M. Frazier, Ph.D.

Department of Energy

- Marvin Stodolsky, Ph.D.

Department of the Interior

- Barnett A. Rattner, Ph.D.
- ◇ Sarah Gerould, Ph.D.

Department of Transportation

- George Cushmac, Ph.D.
- ◇ Steve Hwang, Ph.D.

Environmental Protection Agency**Office of Science Coordination and Policy**

- Karen Hamernik, Ph.D.

Office of Research and Development

- ◇ Julian Preston, Ph.D.
- * Suzanne McMaster, Ph.D.

OECD Test Guidelines Program

- * Jerry Smrcek, Ph.D.

Office of Pesticides Programs

- * Amy Rispin, Ph.D.
- * Deborah McCall

Food and Drug Administration

- Leonard M. Schechtman, Ph.D. (Chair)
- Principal Agency Representative
- ◇ Alternate Principal Agency Representative
- * Other Designated Agency Representatives

Office of Science and Health Coordination

- ◇ Richard Canady, Ph.D., D.A.B.T.

Center for Drug Evaluation and Research

- * Abigail C. Jacobs, Ph.D.

Center for Devices and Radiological Health

- * Raju Kammula., D.V.M., Ph.D., D.A.B.T.
- * Melvin E. Stratmeyer, Ph.D.

Center for Biologics Evaluation and Research

- * Richard McFarland, Ph.D., M.D.
- * Ying Huang, Ph.D.

Center for Food Safety and Nutrition

- * David G. Hattan, Ph.D.
- * Robert L. Bronaugh, Ph.D.

Center for Veterinary Medicine

- * Devaraya Jagannath, Ph.D.
- * M. Cecilia Aguila, D.V.M.

National Center for Toxicological Research

- * William T. Allaben, Ph.D.

Office of Regulatory Affairs

- * Lawrence A. D'Hoostelaere, Ph.D.

National Cancer Institute

- Alan Poland, M.D.
- ◇ Vacant

National Institute of Environmental Health Sciences

- William S. Stokes, D.V.M., D.A.C.L.A.M.
- ◇ John R. Bucher, Ph.D., D.A.B.T.
- * Rajendra S. Chhabra, Ph.D., D.A.B.T.
- * Jerrold J. Heindel, Ph.D.

National Institute for Occupational Safety and Health

- Paul Nicolaysen, V.M.D.
- ◇ K. Murali Rao, MD, Ph.D.

National Institutes of Health

- Margaret D. Snyder, Ph.D.

National Library of Medicine

- Vera Hudson, M.S.
- ◇ Jeanne Goshorn, M.S.

Occupational Safety and Health Administration

- Surender Ahir, Ph.D.

294

294	INTERAGENCY COORDINATING COMMITTEE ON THE	
295	VALIDATION OF ALTERNATIVE METHODS (ICCVAM) ACUTE	
296	TOXICITY WORKING GROUP (ATWG)	
297		
298		
299		
300	Consumer Product Safety Commission	325
301	(CPSC)	326 Food and Drug Administration (FDA)
302		327
303	Patricia Brundage, Ph.D.	328 Leonard Schechtman, Ph.D. (ICCVAM Chair)
304	Kailash Gupta, Ph.D.	329 Kenneth Hastings, Ph.D.
305	Cassandra Prioleau, Ph.D.	330 Abigale Jacobs, Ph.D.
306	Marilyn Wind, Ph.D. (ICCVAM Co-	331 Suzanne Morris, Ph.D.
307	Chair)	332 David Morse, Ph.D.
308		333 Thomas Umbreit, Ph.D.
309	Department of Energy (DOE)	334
310		335 National Institute for Occupational
311	Po-Yung Lu, Ph.D.	336 Safety & Health (NIOSH)
312		337
313	Environmental Protection Agency	338 Stephen Reynolds, Ph.D.
314	(EPA)	339
315		340 National Institute of Environmental
316	Karen Hamernik, Ph.D.	341 Health Sciences (NIEHS)
317	Masih Hashim, Ph.D.	342
318	Marianne Lewis, Ph.D.	343 Rajendra Chhabra, Ph.D.
319	Elizabeth Margosches, Ph.D.	344 William Stokes, D.V.M., D.A.C.L.A.M.
320	Debbie McCall, Ph.D.	345 (Executive Director, ICCVAM)
321	John Redden, Ph.D.	346 Raymond Tice, Ph.D., (Deputy Director,
322	Amy Rispin, Ph.D.	347 NICEATM)
323		348
324		
349		

349 **NATIONAL TOXICOLOGY PROGRAM (NTP) INTERAGENCY**
350 **CENTER FOR THE EVALUATION OF ALTERNATIVE**
351 **TOXICOLOGICAL METHODS (NICEATM)**
352
353

354 **Bradley Blackard, M.S.P.H.**
355 ILS, Inc.
356
357 **Sue Brenzel**
358 ILS, Inc.
359
360 **Thomas Burns, M.S.**
361 ILS, Inc.
362
363 **Jeffery Charles, Ph.D., M.B.A.,**
364 **D.A.B.T.**
365 ILS, Inc.
366
367 **Linda Litchfield**
368 ILS, Inc.
369
370 **Deborah McCarley**
371 NIEHS
372
373
374
375
376
377
378
379
380
381
382
383
384
385
386
387
388

374 **Michael Paris**
375 ILS, Inc.
376
377 **William Stokes, D.V.M., Diplomate,**
378 **ACLAM (Director)**
379 NIEHS
380
381 **Judy Strickland, Ph.D., D.A.B.T.**
382 ILS, Inc.
383
384 **Raymond Tice, Ph.D. (Deputy Director)**
385 NIEHS
386
387

389

PREFACE

390

391 The National Toxicology Program (NTP) Interagency Center for the Evaluation of
392 Alternative Toxicological Methods (NICEATM) collaborated with the European Centre for
393 the Validation of Alternative Methods (ECVAM), a component of the European
394 Commission's Joint Research Centre, to further characterize the usefulness of *in vitro* basal
395 cytotoxicity test methods as predictors of starting doses for acute oral systemic toxicity (i.e.,
396 lethality) assays. NICEATM and ECVAM designed a multi-laboratory validation study to
397 evaluate the performance of two standardized *in vitro* basal cytotoxicity test methods using
398 72 reference substances with the ZEBET approach of using the Registry of Cytotoxicity (RC)
399 regression model to estimate starting doses for acute oral systemic toxicity tests. Based on
400 the procedures described in the *Guidance Document on Using In Vitro Data to Estimate In*
401 *Vivo Starting Doses for Acute Toxicity* (hereafter referred to as *Guidance Document*
402 (ICCVAM 2001), the validation study used two mammalian cell types (i.e., BALB/c 3T3
403 mouse fibroblasts [3T3] and primary normal human epidermal keratinocytes [NHK]) with a
404 neutral red uptake (NRU) cell viability endpoint.

405

406 NICEATM developed draft performance standards that could be used to evaluate the
407 acceptability of test methods that are based on similar scientific principles and that measure
408 or predict the same biological or toxic effect as the *in vitro* basal cytotoxicity test methods.
409 The Interagency Coordinating Committee on the Validation of Alternative Methods
410 (ICCVAM) and its Acute Toxicity Working Group (ATWG) reviewed the Background
411 Review Document (BRD) for the NICEATM/ECVAM validation study and the draft
412 performance standards. After commenting and recommending revisions, ICCVAM and the
413 ATWG recommended that these standards be presented to and reviewed by an Expert Peer
414 Panel convened for evaluation of the validation study. ICCVAM, in collaboration with
415 NICEATM, subsequently proposed and sought public comment on draft performance
416 standards for these types of test methods. Following consideration of public and advisory
417 committee comments, ICCVAM will finalize recommended performance standards for *in*
418 *vitro* basal cytotoxicity test methods.

419

420 Performance standards provide the basis by which a validated and accepted proprietary (i.e.,
421 copyrighted, trademarked, registered) or non-proprietary test method has been determined to
422 have sufficient accuracy and reliability for a specific testing purpose. In addition,
423 performance standards should assist other test developers in the validation of test methods
424 that are similar in structure and function, and facilitate acceptance of test methods that adhere
425 to the applicable performance standards.

426

427 This document is available online at <http://iccvam.niehs.nih.gov>; printed copies are available
428 on request from the NICEATM (NIEHS, P.O. Box 12233, MD EC-17, Research Triangle
429 Park, NC 27709; telephone: 919-541-3398, fax: 919-541-0947, e-mail:
430 iccvam@niehs.nih.gov).

431

432 We gratefully acknowledge the significant contributions of the ICCVAM agency
433 representatives and members of the ICCVAM Acute Toxicity Working Group (ATWG) in
434 the preparation of this document, and the NICEATM staff that assisted throughout the
435 process. We also appreciate the constructive suggestions from interested stakeholders in
436 response to the *Federal Register* notice.

437

438 Leonard M. Schechtman, Ph.D.

439 Chair, ICCVAM

440

441 William S. Stokes, D.V.M., Diplomate, A.C.L.A.M.

442 Director, NICEATM

443 Executive Director, ICCVAM

444

444

EXECUTIVE SUMMARY

445

446 The purpose of performance standards is to communicate the basis by which validated new
447 proprietary (e.g., copyrighted, trademarked, registered) and nonproprietary test methods have
448 been determined to have sufficient accuracy and reliability for specific testing purposes.

449 Performance standards can then be used to evaluate the accuracy and reliability of other test
450 methods that are based on similar scientific principles and that measure or predict the same
451 biological or toxic effect. The three elements of performance standards are essential test
452 method components (i.e., structural, functional, and procedural elements of a validated test
453 method that a proposed, mechanistically and functionally similar test method should adhere
454 to), a minimum list of reference chemicals for assessing the accuracy and reliability of the
455 proposed test method, and the accuracy and reliability values that should be achieved by the
456 proposed test method using the minimum list of reference chemicals.

457

458 The Interagency Coordinating Committee on the Validation of Alternative Methods
459 (ICCVAM) has reviewed and evaluated information presented in the Background Review
460 Document (BRD) for a validation study that evaluated the performance of two standardized
461 *in vitro* basal cytotoxicity test methods for predicting starting doses for acute oral systemic
462 toxicity tests. The study, a collaboration of the National Toxicology Program (NTP)
463 Interagency Center for the Evaluation of Alternative Toxicological Methods (NICEATM)
464 and the European Centre for the Validation of Alternative Methods (ECVAM), tested 72
465 reference substances using *in vitro* neutral red uptake (NRU) test methods with two
466 mammalian cell types (i.e., BALB/c 3T3 mouse fibroblasts [3T3] and primary normal human
467 epidermal keratinocytes [NHK]). The 3T3 and NHK NRU test methods are not sufficiently
468 accurate to predict the acute oral toxicity of substances for the purposes of hazard
469 classification. However, for the purposes of acute oral toxicity testing, the 3T3 and NHK
470 NRU test methods may be used in a weight-of-evidence approach to determine the starting
471 dose for the current acute oral *in vivo* acute systemic toxicity protocols for the Up-and-Down
472 Procedure (UDP; OECD 2001a; EPA 2002a) and Acute Toxic Class (ATC; OECD 2001b)
473 method. The performance of other *in vitro* basal cytotoxicity test methods that are based on
474 similar scientific principles and that measure or predict the same biological response (i.e.,

475 basal cytotoxicity and the rat acute oral LD₅₀ [median lethal dose], respectively) should be
476 demonstrated to meet or exceed the accuracy and reliability of the 3T3 and NHK NRU test
477 methods.

478

479 ICCVAM, the ICCVAM Acute Toxicity Working Group (ATWG), and NICEATM have
480 drafted proposed performance standards for the 3T3 and NHK NRU test methods. This
481 document describes those proposed performance standards. The extent to which proposed *in*
482 *vitro* basal cytotoxicity test methods should demonstrate comparable performance to these
483 two *in vitro* NRU cytotoxicity test methods should be considered on a case-by-case basis.
484 While it would be desirable for proposed test methods to have reliability and accuracy values
485 at least as good as that of these two NRU test methods, some flexibility might be acceptable
486 to the extent that it would not compromise the ultimate protection of human and animal
487 health.

488

489 To demonstrate technical proficiency with the validated test method, ICCVAM recommends
490 that the user evaluate his/her ability to calculate IC₅₀ values (i.e., the test concentration that
491 produces 50% inhibition of the endpoint measure) for at least 12 of the 25 reference
492 substances (minimum of two unclassified chemicals and two from each from the five GHS
493 hazard categories) listed in **Table 2-1**. The resulting IC₅₀ values should be within 2.5
494 standard deviations of the reported RC IC₅₀ values in the table and a linear regression
495 calculated with the IC₅₀ values from the proposed test method and LD₅₀ values in the table
496 should not be different from a linear regression calculated using the IC₅₀ and LD₅₀ values
497 from the table. Intralaboratory Coefficient of Variation (CV) for the IC₅₀ of the reference test
498 substances should not exceed 129% for reference substance test chemicals and the mean for
499 the substances tested should not exceed 30%. The proposed *in vitro* basal cytotoxicity test
500 method should provide correct predictions of GHS acute oral toxicity category using the
501 recommended IC₅₀-LD₅₀ regression for at least 38% of the substances tested.

502

503

504 **1.0 PURPOSE AND BACKGROUND OF PERFORMANCE STANDARDS**

505

506 **1.1 Introduction**

507

508 Prior to the acceptance of a new test method for regulatory testing applications, validation
509 studies are conducted to assess its reliability (i.e., the extent of intra- and inter-laboratory
510 reproducibility) and its relevance (i.e., the ability of the test method to correctly predict or
511 measure the biological effect of interest) (ICCVAM 1997, 2003; OECD 1996, 2002). The
512 purpose of performance standards is to communicate the basis by which new proprietary (i.e.,
513 copyrighted, trademarked, registered) and nonproprietary test methods have been determined
514 to have sufficient relevance and reliability for specific testing purposes. These performance
515 standards, based on test methods accepted by regulatory agencies, can be used to evaluate the
516 reliability and relevance of other test methods that are based on similar scientific principles
517 and measure or predict the same biological or toxic effect. Two *in vitro* basal cytotoxicity
518 test methods, the BALB/c 3T3 mouse fibroblast (3T3) neutral red uptake (NRU) assay and
519 the normal human keratinocyte (NHK) NRU assay, underwent a validation process to
520 evaluate the correlation between *in vitro* cytotoxicity and acute lethality and the feasibility of
521 using *in vitro* NRU assays to predict starting doses for *in vivo* acute oral systemic toxicity
522 assays.

523

524 This section describes the three elements of performance standards identified by ICCVAM
525 (ICCVAM 2003) and the ICCVAM process used to develop performance standards during a
526 test method evaluation. These test method performance standards are proposed as initial
527 standards that can be used to evaluate future *in vitro* basal cytotoxicity test methods. If other
528 *in vitro* basal cytotoxicity test methods are adequately validated and demonstrate
529 significantly improved performance, then the test method performance standards may be
530 revised accordingly.

531

532

532 1.2 Elements of ICCVAM Performance Standards

533

534 Performance standards are standards based on a validated test method that provide a basis for
535 evaluating the comparability of a proposed test method that is mechanistically and
536 functionally similar (ICCVAM 2003). The three elements of performance standards are:

- 537 • **Essential test method components:** These consist of essential structural,
538 functional, and procedural elements of a validated test method that should be
539 included in the protocol of a proposed, mechanistically and functionally similar
540 test method. Essential test method components include unique characteristics of
541 the test method, critical procedural details, and quality control measures.
542 Adherence to essential test method components will help to assure that a
543 proposed test method is structurally and functionally similar to the
544 corresponding validated test method.
- 545 • **A minimum list of reference substances:** Reference substances are used to
546 assess the accuracy and reliability of a proposed, mechanistically and
547 functionally similar test method. These substances are a representative subset
548 of those used to demonstrate the reliability and the accuracy of the validated test
549 method. To the extent possible, this subset of substances should:
 - 550 – be representative of the range of responses that the validated test method is
551 capable of measuring or predicting
 - 552 – have produced consistent results in the validated test method and in the *in*
553 *vivo* reference test method and/or the species of interest
 - 554 – have well-defined chemical structures
 - 555 – be readily available
 - 556 – not be associated with excessive hazard or prohibitive disposal costs
 - 557 – have performance characteristics (accuracy, sensitivity, specificity, false
558 negative and false positive rates) of the validated test method for the subset
559 of reference chemicals should approximate the performance values obtained
560 during the validation process for all appropriate substances

561

562 These reference substances are the minimum number that should be used to
563 evaluate the performance of a proposed, mechanistically and functionally similar
564 test method. Reference substances should not be used to develop the prediction
565 model for the proposed test method. However, if reference substances have been
566 used for this purpose, then they should be replaced with other substances of the
567 same chemical class and biological activity for the endpoint of interest and for
568 which adequate reference data are available. Similarly, if any of the recommended
569 reference substances are unavailable, other substances for which adequate reference
570 data are available could be substituted. Again, to the extent possible, the substituted
571 substance(s) should be of the same chemical class as the original reference
572 substance(s). If desired, additional substances representing other chemical or
573 product classes and for which adequate reference data are available can be used to
574 more comprehensively evaluate the accuracy of the proposed test method.
575 However, none of these additional substances should have been used to develop the
576 proposed test method.

- 577 • **Accuracy and reliability values:** These are the accuracy and reliability
578 characteristics that the proposed test method should be comparable to when
579 evaluated using the minimum list of reference chemicals.

580

581 **1.3 ICCVAM Process for the Development of Performance Standards**

582

583 The process followed by ICCVAM for developing performance standards for new test
584 methods is as follows:

- 585 • NICEATM and the appropriate ICCVAM working group develop proposed
586 performance standards for consideration during the ICCVAM evaluation
587 process. If performance standards are proposed by a test method sponsor, they
588 will be considered by ICCVAM at this stage. Generally, the proposed
589 performance standards are based on the information and data provided in the
590 test method submission or on other available applicable data.
- 591 • The ICCVAM/NICEATM Peer Review Panel evaluates the proposed
592 performance standards for completeness and appropriateness during its

593 evaluation of the validation status of the proposed test method. The proposed
594 performance standards, as well as the test method submission, are made
595 available to the public for comment prior to and during the Peer Review Panel
596 meeting.

597 • The appropriate ICCVAM working group, with the assistance of NICEATM,
598 prepares the final performance standards for ICCVAM approval, taking into
599 consideration the recommendations of the Peer Review Panel and public
600 comments.

601

602 Performance standards recommended by ICCVAM are incorporated into ICCVAM test
603 method evaluation reports, which are then provided to U.S. Federal agencies and made
604 available to the public. Regulatory authorities can then reference the performance standards
605 in the ICCVAM report when they communicate their acceptance of a new test method. In
606 addition, performance standards adopted by U.S. Federal regulatory authorities can be
607 provided in guidelines issued for new test methods. Availability of ICCVAM test method
608 evaluation reports are announced routinely in the *Federal Register*, NTP Newsletters, and by
609 e-mail to ICCVAM/NICEATM listserv groups.

610

611 **1.4 ICCVAM Development of Recommended Performance Standards for** 612 ***In Vitro* Acute Toxicity Test Methods**

613

614 1.4.1 Current Regulatory Testing Requirements for Acute Systemic Toxicity

615 The major regulatory requirement for acute systemic toxicity testing is for the hazard
616 classification and labeling of products, which is intended to protect handlers and consumers
617 from toxic hazards. The LD₅₀ results (i.e., median lethal dose) from acute systemic toxicity
618 tests are used to place substances in various toxicity categories that, in turn, invoke the
619 associated hazard phrases to be used on product labels. **Table 1-1** shows the current U.S.
620 legislation requiring the use of acute systemic toxicity testing for product labeling and the
621 substances regulated. **Table 1-2** shows the statutory protocol requirements and classification
622 systems used by each U.S. regulatory agency. Also included is an international guideline for
623 labeling, the Harmonized Integrated Classification System for Human Health and

624 Environmental Hazards of Chemical Substances and Mixtures (OECD 2001c), which
 625 provides guidance to regulatory agencies on the use of the Globally Harmonized System of
 626 Classification and Labelling of Chemicals (GHS; UN 2005) as a method for an
 627 internationally comprehensible system for hazard communication.

628 **Table 1-1 Summary of Current U.S. Legislation for Using Acute Systemic Toxicity**
 629 **Data for Product Labeling**

Legislation (Year of Initial Enactment)	U.S. Regulatory Agency	Substance
Federal Insecticide, Fungicide and Rodenticide Act (1947)	EPA	Pesticides
Federal Hazardous Substances Act (1964)	CPSC	Household products
Occupational Safety and Health Act (1970)	OSHA	Occupational materials
Federal Hazardous Material Transportation Act (1975)	DOT	Transported substances

630 Abbreviations: EPA = U.S. Environmental Protection Agency; CPSC = U.S. Consumer Product Safety
 631 Commission; OSHA = U.S. Occupational Safety and Health Administration; DOT = U.S. Department of
 632 Transportation. [Note: The U.S. Food and Drug Administration (FDA) does not require data for acute lethality
 633 testing, and in fact, discourages the use of animals for such testing (FDA 1993).]

634
 635

636 1.4.2 Test Methods for Assessing Acute Systemic Toxicity

637 The current internationally recognized test methods for acute systemic toxicity testing are the
 638 UDP (OECD 2001a; EPA 2002a), the ATC method (OECD 2001b), and the Fixed Dose
 639 Procedure (FDP; OECD 2001d). Information on signs of acute toxicity and target organs can
 640 be obtained using any of the three test methods. All three methods are sequential tests in
 641 which the outcome of testing one or more animals at the first dose is used to determine the
 642 second dose that should be tested. The FDP differs from the UDP and ATC in that it
 643 involves testing more animals per dose and the primary endpoint of interest is evident
 644 toxicity³ rather than lethality. The ATC method provides a range for the LD₅₀ for
 645 classification purposes. The UDP generally provides a point estimate of the LD₅₀ with a
 646 confidence interval (EPA 2002a).

647

³ *Evident toxicity* is a general term describing clear signs of toxicity following administration of test substance, such that an increase to the next highest fixed dose would result in the development of severe toxic signs and probably mortality (ICCVAM 2000).

648 **Table 1-2 Regulatory Classification Systems for Acute Oral Toxicity**

Regulatory Agency (Authorizing Act)	Animals	Endpoint	Classification
EPA (Federal Insecticide, Fungicide and Rodenticide Act)	Use current EPA or OECD protocol	Death ¹	I – LD ₅₀ ≤ 50 mg/kg II – 50 < LD ₅₀ ≤ 500 mg/kg III – 500 < LD ₅₀ ≤ 5000 mg/kg IV – LD ₅₀ > 5000 mg/kg
CPSC (Federal Hazardous Substances Act)	White rats, 200-300 g	Death ¹ within 14 days for ≥ half of a group of ≥ 10 animals	Highly toxic – LD ₅₀ ≤ 50 mg/kg Toxic – 50 mg/kg < LD ₅₀ < 5 g/kg
OSHA (Occupational Safety and Health Act)	Albino rats, 200-300 g	Death ¹ , duration not specified.	Highly toxic - LD ₅₀ ≤ 50 mg/kg Toxic – 50 < LD ₅₀ < 500 mg/kg
DOT (Federal Hazardous Material Transportation Act)	Male and female young adult albino rats	Death ¹ within 14 days of half the animals tested. Number of animals tested must be sufficient for statistically valid results.	<i>Packing Group I</i> LD ₅₀ ≤ 5 mg/kg <i>Packing Group II</i> 5 < LD ₅₀ ≤ 50 mg/kg <i>Packing Group III</i> LD ₅₀ < 500 mg/kg (liquid) LD ₅₀ < 200 mg/kg (solid)
OECD Guidance for Use of GHS (2001a)	Protocol not specified	Protocol not specified	I - LD ₅₀ ≤ 5 mg/kg II - 5 < LD ₅₀ ≤ 50 mg/kg III - 50 < LD ₅₀ ≤ 300 mg/kg IV - 300 < LD ₅₀ ≤ 2000 mg/kg V - 2000 < LD ₅₀ ≤ 5000 mg/kg Unclassified - LD ₅₀ > 5000 mg/kg

649 ¹Guidance Document on the Recognition, Assessment and Use of Clinical Signs as Humane Endpoints for
650 Experimental Animals Used in Safety Evaluation calls for humane killing of moribund animals (OECD 2000).
651 Moribund animals that are humanely euthanized are accepted as deaths.
652 Abbreviations: EPA = U.S. Environmental Protection Agency; CPSC = U.S. Consumer Product Safety
653 Commission; OECD = Organisation for Economic Co-operation and Development; OSHA = U.S. Occupational
654 Safety and Health Administration; DOT = U.S. Department of Transportation; GHS = Globally Harmonized
655 System of Classification and Labelling of Chemicals (UN 2005)

656

657 **1.4.3 Intended Regulatory Uses for *In Vitro* Cytotoxicity Test Methods**

658 *In vitro* cytotoxicity test methods are not recommended for the replacement of acute oral
659 toxicity tests in animals. Rather, such test methods are intended to serve as adjuncts for *in*
660 *vivo* acute systemic toxicity test methods. To select a starting dose for a test substance, the
661 current test guidelines for acute oral systemic toxicity recommend using information on
662 structurally-related substances and the results of any other toxicity tests, including *in vitro*
663 cytotoxicity results (OECD 2001a, b, d; EPA 2002a). *In vitro* basal cytotoxicity test methods
664 are intended to be used as part of the weight-of-evidence approach to select starting doses for
665 the UDP and ATC assays in order to reduce and refine the use of animals for *in vivo* acute

666 toxicity testing. Since the estimation of the true LD₅₀ is irrelevant to setting doses for
667 measuring evident toxicity, the use of *in vitro* basal cytotoxicity test methods for setting
668 starting doses for the FDP was not considered in the NICEATM/ECVAM validation study.

669

670 1.4.4 Similarities and Differences in the Endpoints of *In Vitro* Cytotoxicity Test Methods 671 and *In Vivo* Acute Oral Toxicity Test Methods

672 The endpoint measured in the *in vitro* NRU cytotoxicity test methods is cell death (neutral
673 red [NR] is taken up only by live cells) and the major endpoint of interest is the concentration
674 at 50% inhibition of NRU (i.e., the IC₅₀). The endpoint measured in acute systemic toxicity
675 assays is usually animal death. Cell death and animal death may be similar since animals are
676 comprised of organ systems consisting of tissues, which are comprised of cells. All cells,
677 regardless of whether they are in animals or *in vitro* cell cultures, have similar cellular
678 mechanisms of energy production and utilization and maintenance of cell membrane
679 integrity. Animal death and death of cells in culture due to toxicity are similar in that both
680 involve some type of cellular injury. For the animal, the cellular injury produces tissue and
681 organ injury to the most sensitive target organ, which may then cause the death of the whole
682 organism. Cell death in a culture system involves the death of a single cell type. Cell death
683 and animal death may be produced by the same mechanisms, such as disruption of membrane
684 structure or function, inhibition of mitochondrial function, disturbance of protein turnover,
685 disruption of energy production, etc. (Gennari et al. 2004).

686

687 Animal and cell culture systems are different with respect to how a substance or toxin is
688 delivered to the cell and how it is distributed, metabolized, and excreted. After oral
689 administration, animals must absorb the toxin from the gastrointestinal tract, which involves
690 the passage of membranes. The toxin may or may not be heavily bound to serum proteins;
691 this would reduce the availability of the toxin to the target organ. The toxin may then be
692 metabolized during and/or after distribution to the target organs and then the toxin or its
693 metabolites are excreted. In a cell culture system, the only membranes that must be passed
694 are those of the target cell and cellular organelles. No absorption and distribution by other
695 cellular systems is required. Cell culture systems may or may not include serum proteins,
696 which could reduce the availability of toxin to act as its target site. Excretion from the cell

697 culture milieu cannot occur since cell culture systems have no excretory system. The
698 cultured cells are exposed to substances for the entire duration of exposure in the test system.

699

700 Animal and cell culture systems may also be different with respect to the target on which a
701 toxin acts. If a toxin acts in a specialized organ system in a whole animal, it may not produce
702 a toxic effect by the same mechanism in cultured cells that are derived from tissue different
703 from the target organ. For example, a neurotoxin that acts by a neuroreceptor-mediated
704 pathway in animals, would be expected to produce toxicity by a different mechanism in 3T3
705 or NHK cells, which are derived from fibroblasts, and skin cells, respectively. Even if a
706 neurotoxin were applied to neuronal cells in culture, the cultured cells may not respond in the
707 same way as neuronal cells in a whole animal. Cultured cells may not retain the same
708 functionality as cells *in vivo*.

709

710 **2.0 IN VITRO ACUTE TOXICITY TEST METHODS**

711

712 **2.1 Background**

713

714 Pre-validation and validation studies have been completed to evaluate the ability of the 3T3
715 and NHK NRU test methods to be used in a weight of evidence approach (and not as stand-
716 alone tests) to estimate acute rat oral LD₅₀ values, which are to be used as the basis for
717 selecting the starting dose for acute systemic toxicity (i.e., lethality) studies. This section
718 briefly describes the principles of *in vitro* basal cytotoxicity test methods followed by the
719 recommended performance standards that would be used to evaluate test methods that are
720 functionally and mechanistically similar to the 3T3 and NHK NRU test methods. The
721 performance standards consist of 1) essential test method components, 2) reference
722 substances, and 3) the comparable accuracy and reliability that should be achieved.

723 Theoretically, any *in vitro* test capable of determining basal cytotoxicity could be used to
724 estimate the starting dose for acute systemic toxicity testing using the UDP (EPA 2002c,
725 OECD 2001a) or the ATC method (OECD 2001b).

726

727

728

729 **2.2 Principles of *In Vitro* Basal Cytotoxicity Assays to Predict Starting**
730 **Doses for Acute Oral Toxicity Tests**

731

732 The basic steps of an *in vitro* basal cytotoxicity assay are as follows:

- 733 • the test substance is dissolved in an appropriate solvent and applied as a solution
734 to a sub-confluent cell culture monolayer or to exponentially growing cultures
735 of non-adherent cells
- 736 • the test substance is incubated with the cells for a specified period of time
- 737 • the test substance is removed and an endpoint indicative of cell viability or
738 cytotoxicity is measured
- 739 • the IC₅₀ value is calculated (i.e., the concentration at which cell viability or
740 growth is inhibited by 50% compared to control values)

741

742 Many different *in vitro* basal cytotoxicity methods can be used to estimate the acute rat oral
743 LD₅₀ value (the calculated value of the oral dose that produces lethality in 50% of test
744 animals) and, thus, to predict the starting dose for an acute lethality assay. *In vitro* basal
745 cytotoxicity data determined in various primary cells, as well as in various permanent non-
746 differentiated finite or transformed cell lines, generally exhibit the same dose response
747 cytotoxicity relationship in response to the same xenobiotic, regardless of the type of toxic
748 endpoints investigated. The following endpoints are sufficiently characteristic of basal
749 cytotoxicity (Spielmann et al. 1999; Halle 1998):

750

- 751 • Inhibition of cell proliferation: cell number, cell protein, DNA content, DNA
752 synthesis, colony formation
- 753 • Cell viability - metabolic markers: metabolic inhibition test, mitochondrial
754 reduction of tetrazolium salts
- 755 • Decreased cell viability - membrane markers: NRU into cell lysosomes, Trypan
756 Blue exclusion, cell detachment for monolayer cultures
- 757 • Differentiation markers: functional differentiation within cell clusters,
758 morphological differentiation within cell clusters, intracellular morphology

759

760 Markers of the release of intracellular components, such as the enzyme lactate
761 dehydrogenase (i.e., LDH release test), or of dye introduced into the cells previous to
762 chemical exposure, as occurs, for example, in the fluorescein leakage [FL] test or the Neutral
763 Red Release [NRR] test, are not considered to be characteristic for basal cytotoxicity because
764 they specifically detect damage of the outer cell membrane and generally are associated with
765 short-term chemical exposure (ICCVAM 2001). A chemical that specifically damages only
766 cell membranes, however, will be detected correctly in one of the tests for basal cytotoxicity
767 listed above.

768

769 Investigators using an *in vitro* basal cytotoxicity system for prediction of the *in vivo* starting
770 dose for acute oral toxicity/lethality studies must be able to demonstrate that the assay is
771 valid for its intended use. This includes demonstrating that any modification to the existing
772 validated reference test method does not adversely affect its performance characteristics. *In*
773 *vitro* systems may be used to test solids, liquids, and emulsions of any chemical or product
774 class. The liquids can be aqueous or nonaqueous; solids can be soluble or insoluble in water.
775 The samples may be pure chemicals, dilutions, formulations, or waste. Test substances must
776 be soluble in cell culture medium, dimethyl sulfoxide (DMSO), or ethanol (ETOH). The test
777 method endpoint (i.e., percent of control values) will be used to generate an IC₅₀ value in
778 µg/mL and this value will be used in the regression formula developed to estimate the LD₅₀
779 value in mg/kg.

780

781 **2.3 Essential Test Method Components for In Vitro Basal Cytotoxicity** 782 **Assays to Predict Starting Doses for Acute Oral Toxicity (Lethality)** 783 **Tests**

784

785 Essential test method components consist of essential structural, functional, and procedural
786 elements of a validated test method that should be included in the protocol of a
787 mechanistically and functionally similar proposed test method. These components include
788 unique characteristics of the test method, critical procedural details, and quality control

789 measures. Adherence to these components will help assure that a proposed test method is
790 based on the same concepts as the corresponding validated test method.

791

792 The following is a description of the essential test method components for *in vitro* basal
793 cytotoxicity assays to predict starting doses for acute oral toxicity/lethality tests.

794

795 2.3.1 In Vitro Cell Culture Conditions

796 • Use a mammalian cell line (or primary cells) that divides rapidly with doubling
797 times of less than 30 hours under standard culture conditions, preferably with
798 calf serum [CS], newborn calf serum [NCS]), or serum-free medium (ICCVAM
799 2001).

800 • Propagate cells in sterile tissue culture vessels (e.g., flasks) and then subculture
801 cells to other sterile tissue culture vessels (e.g., 96 well-plates) for use in testing.
802 Initial cell seeding should be done at a density that allows rapid growth
803 throughout the exposure period. However, cell density should not reach
804 confluency by the end of the test exposure period.

805 • Maintain appropriate cell culture growth conditions throughout the testing
806 period (e.g., 37°C ± 1°C, 90% ± 10% humidity, 5.0% ± 1% CO₂/air). The cell
807 cultures should be free of contamination with bacteria, mycoplasma, or fungi.

808

809 Cell culture media should be prequalified by the testing laboratory via a standardized
810 protocol before initiating the test to guarantee that the media provide cells with appropriate
811 nutrients to meet the growth criteria needed for the test method.

812

813 2.3.2 Application of the Test Substances

814 *Test Substance Preparation*

815 • Test substance solutions should be prepared in cell culture medium within one
816 hour (unless otherwise known stability conditions of the substance require
817 different parameters) before application to the cell cultures.

818 • Standard protocol methods for solubility procedures can include mixing the test
819 substance by vortexing, sonication, warming, and stirring. Test substances

820 should be fully solubilized (i.e., no visual observation of test substance in the
821 dosing solution) before application.

822 • One inherent limitation to *in vitro* cytotoxicity is the testing of volatile
823 chemicals since the material may evaporate before application to the cells or
824 may not remain in the test vessel when incubated. If volatility is predicted or
825 identified for a test substance (e.g., by detection of cross-contamination of the
826 high concentrations of test substance in culture with lower concentrations or
827 controls in the test vessel), measures can be employed to test moderately
828 volatile chemicals (e.g., cover the test plate with a CO₂ permeable plastic film
829 cover/sealer).

830

831 *Cytotoxicity Test*

- 832 • Each cytotoxicity test should contain a range of test substance concentrations
833 such that the IC₅₀ value can be determined with at least one cytotoxic point
834 between 0 – 50% viability and at least one cytotoxic point between 50 – 100%
835 viability.
- 836 • A minimum of three adequate data points should be collected for each test
837 substance concentration. Note: The NICEATM/ECVAM validation study
838 required the testing of six replicates for each test substance concentration with
839 at least four successful replicates].
- 840 • Blanks (i.e., culture vessels without cells) should be available for assessing
841 background interference when measuring the endpoint.
- 842 • Cell monolayers in tissue culture vessels should be adequately covered (e.g., a
843 minimum of 100 µL of test substance solution per well in a 96-well test plate).
- 844 • The chemical exposure period should be at least the duration of one cell cycle,
845 i.e., 24 – 72 hours (Riddell et al. 1986). [Note: The NICEATM/ECVAM
846 validation study required an exposure period of 48 hours for 3T3 and NHK
847 cells; the cell cycle duration for these cells is generally 17-19 and 10-22 hours
848 in the log phase, respectively.]
- 849 • At the end of the exposure period, most endpoints require washing the test
850 substance from the cells with an appropriate buffering solution (e.g., Dulbecco's
851 Phosphate Buffered Saline [DPBS]) before applying the endpoint material (e.g.,

852 neutral red dye). Washing cells to remove the test substance is the default
853 recommendation unless it is known that washing would interfere with
854 measurement of the endpoint.

855

856 2.3.3 Control Substances

857 *Vehicle Controls:* The vehicle controls provide the reference for 100% cell growth in the test
858 vessel and, thus, the vehicle (or solvent) must be compatible with the cell culture system (i.e.,
859 not cause cytotoxicity or reduce cell growth through other mechanisms) and should not alter
860 the properties of the test substance. The vehicle controls should contain the solvent at the
861 concentration applied to the cells. For example, DMSO and ethanol at a final concentration
862 $\leq 0.5\%$ [v/v] were demonstrated to be compatible with cell growth for 3T3 and NHK cells in
863 the NICEATM/ECVAM validation study. If the compatibility of the solvent with the cell
864 culture system is unknown, cultures with and without the solvent must be included in each
865 experiment.

866

867 *Positive Controls:* The purpose of a positive control chemical is to demonstrate that the cell
868 culture system is responding with adequate sensitivity to a cytotoxic agent for which the
869 magnitude of the cytotoxic response is well characterized. The positive control chemical
870 should be tested concurrently with (and independent of) the test substance. The positive
871 control should be well characterized for its cytotoxicity potential and each test should
872 generate a response that is comparable to the historic IC_{50} range generated by the laboratory.
873 A laboratory should perform a minimum of 10 cytotoxicity tests using the positive control
874 over a number of days to develop a minimum historical database of IC_{50} data. Typically, for
875 biologically based test methods, suggested acceptable ranges for the positive control response
876 are within two to three standard deviations of the historical mean response, but developers of
877 proprietary test methods may establish tighter ranges. Sodium lauryl sulfate (SLS) is an
878 effective positive control substance for use in *in vitro* basal cytotoxicity test methods. [Note:
879 The NICEATM/ECVAM validation study used SLS as the positive control and required 2.5
880 standard deviations of the historical mean response as the acceptable range.]

881

882 *Benchmark Controls:* Benchmark controls may be useful to demonstrate that the test method
883 is functioning properly for detecting the cytotoxic potential of chemicals of a specific
884 chemical class or a specific range of responses, or for evaluating the relative cytotoxic
885 potential of a cytotoxic test substance. Appropriate benchmark controls should have the
886 following properties:

- 887 • consistent and reliable source(s) for the chemical
- 888 • structural and functional similarity to the class of the substance being tested
- 889 • known physical/chemical characteristics
- 890 • supporting data on known effects in animal models
- 891 • known potency in the range of response (including moderate response)

892

893 2.3.4 Viability Measurements

- 894 • Only standardized, quantitative methods should be used to measure cell
895 viability. The protocol should be compatible with laboratory apparatus such as
896 spectrophotometers that allow a quick and precise measurement of the endpoint.
- 897 • Non-specific dye binding must not interfere with the viability measurement. A
898 measurement endpoint that is well established and that has good interlaboratory
899 reproducibility should be used (ICCVAM 2001).
- 900 • A detailed concentration-response experiment should be conducted using a
901 progression factor that yields graded effects between no effect and total
902 cytotoxicity. Any desired toxicity measure can be derived from a well-designed
903 concentration-response experiment.
- 904 • Preference should be given to endpoints that determine either cell proliferation
905 or cell viability (e.g., NRU, MTT [3-(4,5-dimethylthiazol-2-yl)2,5-diphenyl
906 tetrazolium bromide], XTT [Sodium 3,3-[(Phenylamino)carbonyl]-3,4-
907 Tetrazolium-Bis(4-methoxy-6-nitro)benzenesulfonic acid hydrate]) (ICCVAM
908 2001).
- 909 • Simple endpoints such as total protein content are not recommended, as they
910 may under-predict the toxicity of certain test chemicals by including protein
911 from dead cells.

- 912 • A lack of information and a low level of accuracy characterize experiments that
913 seek only to identify the highest tolerated dose or the lowest cytotoxic dose.

914

915 Colorimetric endpoints (e.g., NRU) should have the optical density (OD)
916 spectrascopically-measured at the appropriate wavelength (e.g., 540 nm for NRU) and
917 OD values for blanks should be subtracted from the vehicle control and test substance
918 ODs.

919

920 2.3.5 Interpretation of Results

921 *IC₅₀ Determination*

922 The endpoint values obtained for each test sample can be used to calculate the percentage of
923 cell viability or growth relative to the negative (vehicle) control, which is arbitrarily set at
924 100%. The cell viability criteria used to determine an IC₅₀ value must be clearly defined and
925 documented, and be shown to be appropriate. In general, such criteria are established during
926 test optimization, tested during a prevalidation phase, and confirmed in a validation study.

927

928 *Regression Formula*

929 The regression formula for *in vitro* acute toxicity test methods in the NICEATM/ECVAM
930 validation study was determined from the chemical IC₅₀ values (µg/mL) and rat LD₅₀ values
931 (mg/kg) from the Registry of Cytotoxicity (RC) database. The RC contains acute oral LD₅₀
932 values for rats and mice obtained from the Registry of Toxic Effects for Chemical Substances
933 (RTECS[®]) and IC₅₀ values from *in vitro* cytotoxicity assays using multiple cell lines and
934 cytotoxicity endpoints for chemicals with known molecular weights (Halle 1998). Of the
935 347 chemicals in the RC, 282 chemicals have rat LD₅₀ values. Those chemicals without rat
936 data were excluded. Chemicals with mechanisms of action not expected to be active in
937 culture cells used for basal cytotoxicity test methods were excluded from the database prior
938 to deriving the regression line. The regression formula based upon rat data and excluding
939 chemicals with specific mechanisms other than cytotoxicity was developed from the data for
940 232 chemicals.

941

$$942 \quad \log(\text{LD}_{50} \text{ in mg/kg}) = 0.357 \times \log(\text{IC}_{50} \text{ in } (\mu\text{g/mL})) + 2.194$$

943

944 Before using a candidate *in vitro* basal cytotoxicity test to predict starting doses, the
945 correlation between the *in vitro* test and the *in vivo* test must be established quantitatively.
946 This can be accomplished by testing in the new assay a minimum of 12 reference substances
947 from **Table 2-1** that cover all six hazard classification categories (i.e., the entire range of
948 cytotoxicity) and that produce the same regression formula as the total database. After
949 testing, the IC₅₀ data are used to calculate a regression formula (least square method) for the
950 selected reference substances using the corresponding LD₅₀ values provided in **Table 2-1**.
951 The resulting regression is compared against the original regression. If the regressions are
952 not statistically significantly different (at $p < 0.05$; e.g., comparison of slope and intercept),
953 then the test is considered suitable to generate IC₅₀ data to use with the recommended
954 regression formula for estimating starting doses for acute oral systemic toxicity/lethality
955 tests.

956

957 2.3.6 Test Report

958 The test report should include the following information, if relevant to the conduct of the
959 study:

960 *Test Substances and Control Substances*

- 961 • Chemical name(s) such as Chemical Abstracts Service Registry Number
962 (CASRN) and molecular weight (if known), followed by other names, if known
- 963 • Formulation (if available) of the test substance if the material is a mixture
- 964 • Purity and composition of the substance or preparation (in percentage[s] by
965 weight)
- 966 • Physicochemical properties such as physical state, volatility, pH, stability,
967 chemical class, water solubility relevant to the conduct of the study
- 968 • Treatment of the test/control substances prior to testing, if applicable (e.g.,
969 vortexing, sonication, warming; solvent used)
- 970 • Stability, if known

971 *Justification of the In Vitro Test Method and Protocol Used*

972 *Test Method Integrity*

- 973 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
974 test method over time
- 975 • If the test method employs proprietary components, documentation on the
976 procedure used to ensure their integrity from “lot-to-lot” and over time
- 977 • The procedures that the user may employ to verify the integrity of the
978 proprietary components

979 *Criteria for an Acceptable Test*

- 980 • Acceptable concurrent positive control ranges based on historical data
- 981 • Acceptable negative and solvent/vehicle control data

982 *Test Conditions*

- 983 • Cell system used
- 984 • Calibration information for measuring device used for measuring cell viability
985 (e.g., spectrophotometer)
- 986 • Details of test procedure used
- 987 • Test doses used
- 988 • Description of any modifications of the test procedure
- 989 • Reference to historical data of the model
- 990 • Description of evaluation criteria used

991 *Results*

- 992 • Tabulation of data from individual test samples (e.g., OD values and calculated
993 percentage cell viability data for the test substance and the positive, negative,
994 and benchmark controls, reported in tabular form, including data from replicate
995 repeat experiments as appropriate, and means and \pm the standard deviation for
996 each trial)
- 997 • Calculated IC₅₀ value
- 998 • Calculated starting dose (i.e., LD₅₀ value) using IC₅₀ value in regression formula
- 999 • Regression formula (prediction model) used

1000 *Description of Other Effects Observed*

1001 *Discussion of the Results*

1002 *Conclusion*

1003

1004

1005 **2.4 Reference Substances for In Vitro Basal Cytotoxicity Assays to Predict**
1006 **Starting Doses for Acute Oral Toxicity Tests**

1007

1008 Reference substances are used to assess the accuracy and reliability of a proposed,
1009 mechanistically and functionally similar test method and are a representative subset of those
1010 used to demonstrate the reliability and the accuracy of the validated test method. These
1011 substances are:

- 1012 • representative of the range of responses that the validated test method is capable
1013 of measuring or predicting
- 1014 • have produced consistent results in the validated test method
- 1015 • reflect the accuracy of the validated test method
- 1016 • have well-defined chemical structures
- 1017 • are readily available
- 1018 • are not associated with excessive hazard or prohibitive disposal costs

1019

1020 To demonstrate technical proficiency with the validated test method, the user should evaluate
1021 his/her ability calculating IC₅₀ values of at least 12 of the reference substances (minimum of
1022 two unclassified chemicals and two from each from the five GHS hazard categories) listed in
1023 **Table 2-1**. The **resulting** IC₅₀ values should be within 2.5 standard deviations of the
1024 reported RC IC₅₀ values in the table.

1025

1026 This subset of 25 reference substances was chosen from the 72 reference substances used in
1027 the NICEATM/ECVAM *in vitro* acute toxicity test methods validation study and are
1028 common to the RC database. Criteria used for selecting the substances in this subset are:

- 1029 • they must be in common with the RC database and have rat LD₅₀ data
- 1030 • they cannot be the same substances as those that were removed from the revised
1031 RC regression formula because of mechanism of action
- 1032 • they cannot be volatile or insoluble (determined by the NICEATM/ECVAM
1033 study)

1034

1035 The low number of substances in the $5 < LD_{50} \leq 50$ mg/kg and $300 < LD_{50} \leq 2000$ mg/kg
1036 GHS categories is due to the exclusionary criteria described above. Reference substances
1037 that exhibited solubility difficulties (with 3T3 medium) or were volatile are included as a
1038 secondary subset and are recommended for investigational purposes only.

1039

1040 The substances in this list represent the following types of chemical classes: alcohols;
1041 amides; boron compounds; cadmium compounds; carboxylic acids; copper compounds;
1042 cyclic hydrocarbons; heterocyclics; mercury compounds; organometalics; potassium
1043 compounds; sodium compounds; and sulfur compounds. Chemicals that exhibited solubility
1044 difficulties (with 3T3 medium) or were volatile are included as a secondary subset and are to
1045 be used for investigational purposes only.

1046

1047

1047 **Table 2-1 Recommended Reference Substances for Evaluation of In Vitro**
 1048 **Basal Cytotoxicity Methods for Predicting Starting Dose for Acute**
 1049 **Systemic Toxicity Tests**

Chemical	CASRN ^a	LD ₅₀ ^b (mg/kg)	IC ₅₀ ^c (µg/mL)
LD₅₀ ≤ 5 mg/kg			
Mercury II chloride	7487-94-7	1	4.07
Triethylenemelamine	51-18-3	1	0.159
Cycloheximide	66-81-9	2	0.166
Busulfan	55-98-1	2	11.3
Phenylthiourea	103-85-5	3	82.2
5 < LD₅₀ ≤ 50 mg/kg			
Triphenyltin hydroxide	76-87-9	44	0.0180
Sodium bichromate	10588-01-9	50	0.244
50 < LD₅₀ ≤ 300 mg/kg			
Hexachlorophene	70-30-4	61	3.21
Cadmium II chloride	10108-64-2	88	1.17
Sodium oxalate	62-76-0	155	59.0
Sodium fluoride	7681-49-4	180	77.7
Diquat dibromide	85-00-7	231	55.1
Cupric sulfate pentahydrate	7758-99-8	300	82.4
300 < LD₅₀ ≤ 2000 mg/kg			
Acetylsalicylic acid	50-78-2	1000	409.0
Propranolol HCl ³	350-60-90	470	35.5
2000 < LD₅₀ ≤ 5000 mg/kg			
Acetaminophen	103-90-2	2404	409.7
Potassium chloride	7447-40-7	2602	6113
Dimethylformamide	68-12-2	2800	8335
Sodium chloride	7647-14-5	2998	4436
Chloramphenicol	56-75-7	3393	255.3
Lactic acid	50-21-5	3730	5946
Trichloroacetic acid	76-03-9	4999	1338
LD₅₀ > 5000 mg/kg			
Ethylene glycol	107-21-1	8567	34454
Dibutyl phthalate	84-74-2	11998	211.6
Glycerol	56-81-5	12691	57477

Chemical	CASRN ^a	LD ₅₀ ^b (mg/kg)	IC ₅₀ ^c (µg/mL)
Secondary Subset			
Precipitating Chemicals ¹			
LD₅₀ ≤ 5 mg/kg			
Arsenic trioxide	1327-53-3	20	1
Parathion	56-38-2	2	27.1
300 < LD₅₀ ≤ 2000 mg/kg			
Giberrellic Acid	77-06-5	6305	797
Volatile Chemicals ²			
300 < LD₅₀ ≤ 2000 mg/kg			
Phenol	108-95-2	414	283.3
LD₅₀ > 5000 mg/kg			
Ethanol	64-17-5	14008	17464
2-Propanol	67-63-0	5843	10038

1050
1051
1052
1053
1054
1055
1056
1057
1058
1059
1060

^aChemical Abstracts Service Registry Number

^bThe calculated value of the oral dose that produces lethality in 50% of test animals (rats and mice).

^cReference substance concentration (geometric mean) producing 50% inhibition of the endpoint measured (i.e., cell viability).

¹Reference substances expected to precipitate at cytotoxic concentrations.

²Reference substances expected to contaminate neighboring wells at high concentrations.

³Propranol HCl data were not used in developing the regression formula because the compound lacked rat data (RC LD₅₀ value = 466 mg/kg [from mouse]). NICEATM found a rat LD₅₀ value of 470 mg/kg in the literature and decided that this reference substance would be an acceptable representative of the 300 < LD₅₀ ≤ 2000 mg/kg category.

1061 2.5 Accuracy and Reliability

1062

1063 The third element of the performance standards is the determination of accuracy (also known
1064 as relevance) and reliability values. The proposed test method, functionally and
1065 mechanistically similar to the 3T3 NRU test method, will use selected reference substances
1066 to assess accuracy and reliability.

1067

1068 2.5.1 Accuracy

1069 When evaluated using the minimum list of recommended reference substances (**Table 2-1**),
1070 the proposed test method should have performance characteristics that are comparable to the
1071 performance of the validated 3T3 NRU test method. Accuracy is defined as the closeness of
1072 agreement between a test method result and an accepted reference value (ICCVAM 2003).
1073 Substances, ranging in toxicity activity from strong to weak, and representing relevant
1074 chemical classes across all GHS hazard categories are included so that the performance of the

1075 proposed test method can be determined and compared to that of the validated reference test
1076 method.

1077

1078 Although the *in vitro* basal cytotoxicity test methods are not intended as replacements for
1079 acute systemic toxicity assays, the accuracy of these assays to predict LD₅₀ values in the
1080 correct GHS acute oral toxicity hazard classification category (UN 2005) was evaluated.
1081 This accuracy evaluation characterizes the extent that additional test methods will be
1082 necessary to achieve accurate *in vitro* predictions of acute oral toxicity hazards for regulatory
1083 classification and labeling purposes.

1084

1085 The log IC₅₀ data from the 3T3 and NHK NRU test methods and the corresponding reference
1086 log LD₅₀ values (mg/kg) values were used to calculate linear regressions for each validation
1087 study laboratory. The slopes for all regressions were statistically significant ($p < 0.0001$),
1088 which indicated a significant relationship between IC₅₀ and LD₅₀ values. The higher adjusted
1089 R² values for the 3T3 regressions, compared with the NHK regressions indicated that the 3T3
1090 IC₅₀ data provided a better fit to the LD₅₀ data than the corresponding NHK data (see **Table**
1091 **2-2**).

1092

1093 **Table 2-2 Linear Regression Analyses of *In Vitro* and *In Vivo* Results¹**

Laboratory	Weight Unit Regressions ²		
	Slope	Intercept	Adjusted R ²
3T3 NRU Test Method			
Lab 1	0.509	1.552	0.420
Lab 2	0.453	1.513	0.307
Lab 3	0.515	1.542	0.421
Combined ³	0.516	1.498	0.409
NHK NRU Test Method			
Lab 1	0.425	1.679	0.319
Lab 2	0.375	1.798	0.276
Lab 3	0.424	1.704	0.318
Combined ³	0.424	1.720	0.322

1094

1095

1096

1097

1098

1099

¹Data for 70 reference substances in the 3T3 assay and 71 reference substances in the NHK assay.

²log IC₅₀ in µg/mL; log LD₅₀ in mg/kg.

³Single regression for the test method using the geometric mean of the laboratory-specific IC₅₀ values for each reference substance.

1100 **Table 2-3** shows that the accuracy of the 3T3 and NHK NRU test methods for predicting
1101 GHS acute oral toxicity categories was similar. The accuracy for the 3T3 regression was
1102 46% while the accuracy for the NHK regression was 38%.

1103

1104 2.5.2 Reliability

1105 Reliability is the degree to which a test method can be performed reproducibly within and
1106 among laboratories over time (ICCVAM 2003) and was assessed in the NICEATM/ECVAM
1107 validation study by determining both intra- and inter-laboratory reproducibility.

1108 Intralaboratory reproducibility is the agreement of results produced when qualified people
1109 within the same laboratory perform the test method using the same test protocol at different
1110 times. Interlaboratory reproducibility is the agreement of results from different qualified
1111 laboratories using the same protocol and reference substances. Interlaboratory
1112 reproducibility indicates the extent to which a test method can be transferred successfully
1113 among laboratories. Intra- and inter-laboratory reproducibility of the 3T3 and NHK NRU
1114 test methods were determined using Analysis of Variance (ANOVA) and Coefficient of
1115 Variation (CV) analysis. Interlaboratory reproducibility of the 3T3 and NHK NRU test
1116 methods was also assessed by comparing the laboratory-specific *in vitro-in vivo* regressions
1117 (for each test method) to one another.

1118 **Table 2-3 Prediction of GHS Toxicity Category¹ by RC Rat-Only Weight Regression Excluding Chemicals with**
 1119 **Specific Mechanisms of Toxicity**

Reference Rodent LD ₅₀ ²	3T3 NRU-Predicted Toxicity Category						Total	Accuracy	Toxicity Overpredicted	Toxicity Underpredicted
	< 5	5 – 50	50 – 300	300-2000	2000-5000	> 5000				
< 5	0	0	2	2	0	0	4 ³	0%	100%	0%
5 – 50	0	1	4	2	0	0	7 ⁴	14%	86%	0%
50 – 300	0	0	4	1	0	0	5 ⁵	80%	20%	0%
300 – 2000	0	1	1	7	0	0	9 ⁶	78%	0%	22%
2000 – 5000	0	0	0	3	6	0	9 ⁷	67%	0%	33%
> 5000	0	0	0	5	4	3	12 ^{8,9}	25%	0%	75%
Total	0	2	11	20	10	3	46	46%	24%	30%
Predictivity	0%	50%	36%	35%	60%	100%				
Category Underpredicted	0%	50%	9%	40%	40%	0%				
Category Overpredicted	0%	0%	55%	25%	0%	0%				
Reference Rodent LD ₅₀ ²	NHK NRU-Predicted Toxicity Category						Total	Accuracy	Toxicity Overpredicted	Toxicity Underpredicted
	< 5	5 – 50	50 – 300	300 – 2000	2000 – 5000	> 5000				
< 5	0	0	2	2	0	0	4 ³	0%	100%	0%
5 – 50	0	1	4	2	0	0	7 ⁴	14%	86%	0%
50 – 300	0	1	3	1	0	0	5 ⁵	60%	20%	20%
300 – 2000	0	1	0	8	0	0	9 ⁶	89%	0%	11%
2000 – 5000	0	0	0	5	4	0	9 ⁷	44%	0%	56%
> 5000	0	0	0	4	7	2	13 ⁹	15%	0%	85%
Total	0	3	9	22	11	2	47	38%	23%	38%
Predictivity	0%	33%	33%	36%	36%	100%				
Category Underpredicted	0%	67%	0%	41%	64%	0%				
Category Overpredicted	0%	33%	67%	23%	0%	0%				

1120 ¹GHS-Globally Harmonized System of Classification and Labelling of Chemicals with LD₅₀ in mg/kg (UN 2005). The RC rat-only weight regression
 1121 excluding chemicals with specific mechanisms of toxicity is $\log LD_{50} \text{ (mg/kg)} = \log IC_{50} \text{ (}\mu\text{g/mL)} \times 0.357 + 2.194$.
 1122
 1123

1124		
1125	< 5:	$LD_{50} \leq 5 \text{ mg/kg}$
1126	5 – 50:	$5 < LD_{50} \leq 50 \text{ mg/kg}$
1127	50 – 300:	$50 < LD_{50} \leq 300 \text{ mg/kg}$
1128	300 – 2000:	$300 < LD_{50} \leq 2000 \text{ mg/kg}$
1129	2000 – 5000:	$2000 < LD_{50} \leq 5000 \text{ mg/kg}$
1130	> 5000:	$LD_{50} > 5000 \text{ mg/kg}$
1131	Numbers in table represent number of reference substances.	
1132	² In vivo reference LD ₅₀ values	
1133		

Goodness of fit F-tests indicated that the laboratory-specific regressions for both test methods were not significantly different from one another ($p = 0.796$ for comparison of the 3T3 regressions and $p = 0.985$ for comparison of the NHK regressions).

The reliability of the proposed test method for the reference substances should be comparable to or better than that of the validated 3T3 NRU test method. **Table 2-4** illustrates intra- and inter-laboratory CV results obtained from the NICEATM/ECVAM *in vitro* acute toxicity test methods validation study.

Table 2-4 Summary of CV Results for the 3T3 and NHK NRU Test Methods

	3T3 NRU Test Method	N	NHK NRU Test Method	N
Intralaboratory CV				
Range	1-122%	202	1-129%	208
Mean CV (and Range)	26%	202	26%	208
Lab 1	23% (2-95%)	68	23% (2-76%)	69
Lab 2	33% (1-98%)	66	42%(1-129%)	69
Lab 3	21% (1-122%)	68	14% (1-38%)	70
Interlaboratory CV				
Range	2-135%	68	1-99%	69
Mean	46%	68	28%	69

CV - coefficient of variation; N - number of values

3.0 REFERENCES

EPA. 2002. Health Effects Test Guidelines OPPTS 870.1100 Acute Oral Toxicity. EPA 712–C–02–190. Washington, DC: U.S. Environmental Protection Agency.

FDA. 1993. Toxicological Principles for the Safety Assessment of Direct Food Additives and Color Additives Used in Food. Center for Food Safety and Applied Nutrition. Washington, DC: U.S. Food and Drug Administration.

Gennari A, van den Berghe C, Casati S, Castell J, Clemedson C, Coecke S, et al. 2004. Strategies to Replace In Vivo Acute Systemic Toxicity Testing. The Report and Recommendations of ECVAM Workshop 50. *Altern Lab Anim* 32:437-459.

Halle W. 1998. Toxizitätsprüfungen in Zellkulturen für eine Vorhersage der akuten Toxizität (LD₅₀) zur Einsparung von Tierversuchen. *Life Sciences/Lebenswissenschaften*, Volume 1, 94 pp., Jülich: Forschungszentrum Jülich. English translation: Halle W. 2003. The Registry of Cytotoxicity: Toxicity testing in cell cultures to predict acute toxicity (LD50) and to reduce testing in animals. *Altern Lab Anim* 31:89-198.

ICCVAM. 1997. Validation and Regulatory Acceptance of Toxicological Test Methods. A Report of the *ad hoc* Interagency Coordinating Committee on the Validation of Alternative Methods. NIH Publication No. 97-3981. Research Triangle Park, NC. National Institute of Environmental Health Sciences. Available: <http://iccvam.niehs.nih.gov/docs/guidelines/validate.pdf> [accessed 14 March 2006].

ICCVAM. 2001. Guidance Document On Using In Vitro Data To Estimate In Vivo Starting Doses For Acute Toxicity. NIH Publication No. 01-4500. National Institute for Environmental Health Sciences, Research Triangle Park, NC. Available: <http://iccvam.niehs.nih.gov/>. [accessed 14 March 2006].

ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised, and Alternative Test Methods. NIH Publication No. 03-4508. National Institute for Environmental Health Sciences, Research Triangle Park, NC. Available: <http://iccvam.niehs.nih.gov/>. [accessed 14 March 2006].

OECD. 1996. Final Report of the OECD Workshop on Harmonisation of Validation and Acceptance Criteria for Alternative Toxicological Methods. (Solna Report) Paris: Organisation for Economic Co-operation and Development. Available: <http://www.oecd.org> [accessed 14 March 2006].

OECD. 2000. Guidance Document on the Recognition, Assessment and Use of Clinical Signs as Humane Endpoints for Experimental Animals Used in Safety Evaluation. Environmental Health and Safety Monograph Series on Testing and Assessment No. 19. Paris France: OECD. Available: http://www.oecd.org/document/22/0,2340,en_2649_34377_1916054_1_1_1_1,00.html

[accessed 14 March 2006]

OECD. 2001a. Guideline for Testing of Chemicals, 425, Acute Oral Toxicity – Up-and-Down Procedure. Paris France: OECD. Available: <http://www.oecd.org> [accessed 14 March 2006].

OECD. 2001b. Guideline For Testing of Chemicals, 423, Acute Oral Toxicity – Acute Toxic Class Method. Paris France: OECD.

OECD. 2001c. OECD Series on Testing and Assessment, 33, Harmonized Integrated Classification System for Human Health and Environmental Hazards of Chemical Substances and Mixtures. ENV/JM/MONO(2001)6. Paris France: OECD. Available: <http://www.oecd.org> [accessed 12 September 2004].

OECD. 2001d. Guideline for Testing of Chemicals, 420, Acute Oral Toxicity – Fixed Dose Method. Paris France: OECD. Available: <http://www.oecd.org> [accessed 12 September 2004].

OECD. 2002. Report of the Stockholm Conference on Validation and Regulatory Acceptance of New and Updated Methods in Hazard Assessment. Paris: Organisation for Economic Co-operation and Development. Available: <http://www.oecd.org> [accessed 14 March 2006].

Riddell RJ, Panacer DS, Wilde SM, Clothier RH, Balls M. 1986. The importance of exposure period and cell type in *in vitro* cytotoxicity tests. *Altern Lab Anim* 14:86-92.

Spielmann H, Genschow E, Liebsch M, Halle W. 1999. Determination of the starting dose for acute oral toxicity (LD₅₀) testing in the up and down procedure (UDP) from cytotoxicity data. *Altern Lab Anim* 27:957-966.

UN. 2005. Globally Harmonized System of Classification and Labelling of Chemicals (GHS), First Revised Edition. [ST/SG/AC.10/30/Rev.1]. United Nations, New York and Geneva. Available: http://www.unece.org/trans/danger/publi/ghs/ghs_rev01/01files_e.html [accessed 14 March 2006].

Appendix B

Draft Recommended Test Method Protocols

B-1	Test Method Protocol for the BALB/c 3T3 Neutral Red Uptake (NRU) Cytotoxicity Test	B-1-3
B-2	Test Method Protocol for the Normal Human Epidermal Keratinocyte (NHK) Neutral Red Uptake (NRU) Cytotoxicity Test.....	B-2-3

[This Page Intentionally Left Blank]

Appendix B-1

Test Method Protocol for the BALB/c 3T3 Neutral Red Uptake (NRU) Cytotoxicity Test

This draft recommended 3T3 NRU test method protocol is substantially the same as the protocol used in Phase III of the NICEATM/ECVAM validation study. Revisions were made based on recommendations from NICEATM and the study directors involved in the study. The changes are as follows:

- Explanations and directions for the use of the revised Hill function for determining IC50 values are included in the protocol.
- The range for relative humidity values for the cell culture incubators was changed from 90 % \pm 5 % humidity to 90 % \pm 10 % humidity.
- An additional step was added to the test substance solubility protocol to allow testing of higher concentrations of test material.
- The spreadsheet templates used in the NICEATM/ECVAM validation study are incorporated into this protocol as an annex (ANNEX I).
- The stand-alone solubility protocol is incorporated into this protocol as an annex (ANNEX II).

[This Page Intentionally Left Blank]

Test Method Protocol

The BALB/c 3T3 Neutral Red Uptake (NRU) Cytotoxicity Test A Test for Basal Cytotoxicity

I. PURPOSE

This test method is used to evaluate the cytotoxicity of test substances using the BALB/c 3T3 Neutral Red Uptake (NRU) *in vitro* cytotoxicity test. The data generated from the *in vitro* cytotoxicity assays are used to evaluate the effectiveness of the assay to predict the starting doses for rodent acute oral systemic toxicity assays. This test method protocol outlines the procedures for performing the basal cytotoxicity test and is the result of the *in vitro* validation study organized by NICEATM and the European Centre for the Validation of Alternative Methods (ECVAM).

If changes or modifications are made to this protocol, the testing laboratory should prove that the results are comparable to those obtained when using the original protocol.

II. TEST SYSTEM

The NRU cytotoxicity assay procedure is a cell survival/viability chemosensitivity assay based on the ability of viable cells to incorporate and bind neutral red (NR), a supravital dye. NR is a weak cationic dye that readily penetrates cell membranes by non-ionic diffusion and accumulates intracellularly in lysosomes. Alterations of the cell surface or the sensitive lysosomal membrane lead to lysosomal fragility and other changes that gradually become irreversible. Such changes brought about by the action of xenobiotics result in a decreased uptake and binding of NR. It is thus possible to distinguish between viable, damaged, or dead cells, which is the basis of this assay.

Healthy mammalian cells, when maintained in culture, continuously divide and multiply over time. A toxic chemical, regardless of site or mechanism of action, will interfere with this process and result in a reduction of the growth rate as reflected by cell number. Cytotoxicity is expressed as a concentration dependent reduction of the uptake of the NR after chemical exposure thus providing a sensitive, integrated signal of both cell integrity and growth inhibition.

III. KEY PERSONNEL

A. Laboratory

- 1) Study Director
- 2) Laboratory Technician(s)

46 **B. Testing Facility**

- 47 1) Scientific Advisor
 48 2) Quality Assurance Director
 49 3) Safety Manager
 50 4) Facility Management
 51

52 **IV. DEFINITIONS**

- 53
 54 3 **Hill function:** a four parameter logistic mathematical model relating the
 55 concentration of test substance to the response being measured in a sigmoidal
 56 shape.
 57

$$58 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log EC_{50} - X)\text{HillSlope}}}$$

59 where Y= response, X is the logarithm of dose (or concentration), Bottom is
 60 the minimum response, Top is the maximum response, logEC50 is
 61 logarithm of X at the response midway between Top and Bottom, and
 62 HillSlope describes the steepness of the curve. When Top = 100 and
 63 Bottom = 0, the EC₅₀ is the concentration at 50% viability (i.e., the IC₅₀).
 64

- 65 4 **Documentation:** all methods and procedures will be noted in a Study Workbook;
 66 logs will be maintained for general laboratory procedures and equipment (e.g.,
 67 media preparation, test substance preparation, incubator function); all optical
 68 density data obtained from the spectrophotometer plate reader will be saved in
 69 electronic and paper formats; all calculations of IC_x values and other derived data
 70 will be in electronic and paper format; all data will be archived.
 71
 72 5 **IC₅₀:** test substance concentration producing 50% inhibition of the endpoint
 73 measured (i.e., cell viability).
 74

75 **V. IDENTIFICATION OF CONTROL SUBSTANCES**

76

77 **A. Positive Control (PC)**

78 Sodium Lauryl Sulfate (SLS)

79

80 **B. Vehicle Control (VC)**

81 Assay medium (DMEM containing 5% NBCS, 4 mM L-Glutamine, 100 IU/mL
 82 Penicillin, 100 µg/mL Streptomycin)
 83

83

84 **C. Solvent Control**

85 VC control with solvent (i.e., assay medium, dimethyl sulfoxide [DMSO], or
 86 ethanol [ETOH])

87 (DMSO is the preferred solvent for substances that are not water [i.e., assay
88 medium] soluble.)
89

90 VI. PROCEDURES

91 A. Materials

92 1. Cell Line

93 BALB/c 3T3 cells, clone A31
94 (e.g., CCL-163, American Type Culture Collection [ATCC], Manassas, VA,
95 USA)
96

97 2. Technical Equipment

98
99 [Note: Suggested brand names/vendors are listed in parentheses. Equivalents
100 may be used.]
101

102 4 Incubator: 37°C ± 1°C, 90 % ± 10 % humidity, 5.0 % ± 1 % CO₂/air

103 5 Laminar flow clean bench/cabinet (standard: "biological hazard")

104 6 Waterbath: 37°C ± 1°C

105 7 Inverse phase contrast microscope

106 8 Sterile glass tubes with caps (e.g., 5 mL)

107 9 Centrifuge

108 10 Laboratory balance

109 11 96-well plate spectrophotometer (i.e., plate reader) equipped with 540 nm ±
110 10 nm filter

111 12 Shaker for microtiter plates

112 13 Cell counter or hemocytometer

113 14 Pipetting aid

114 15 Pipettes, pipettors (multi-channel and single channel; multichannel repeater
115 pipette), dilution block

116 16 Cryotubes

117 17 Tissue culture flasks (e.g., 75 - 80 cm², 25 cm²)

118 18 96-well flat bottom tissue culture microtiter plates (e.g., Nunc # 167 008;
119 Falcon tissue culture-treated)

120 19 pH paper (wide and narrow range)

121 20 Multichannel reagent reservoir

122 21 Waterbath sonicator

123 22 Magnetic stirrer

124 23 Antistatic bar ionizer/antistatic gun (optional for neutralizing static on 96-
125 well plates)

126 24 Dry heat block (optional)

127 25 Adhesive film plate sealers (e.g., Excel Scientific SealPlate™, Cat # STR-
128 SEAL-PLT or equivalent)

129 26 Vortex mixer

130 27 Filters/filtration devices
131

132 [Note: Tissue culture flasks and microtiter plates should be prescreened to
133 ensure that they adequately support the growth of 3T3 cells. Multi-channel
134 repeater pipettes may be used for plating cells in the 96-well plates, dispensing
135 plate rinse solutions, NR medium, and desorb solution. Do not use the repeater
136 pipette for dispensing test substances to the cells.]
137

138 3. Chemicals, Media, and Sera

- 139 a) Dulbecco's Modification of Eagle's Medium (DMEM) without L-
140 Glutamine; should have high glucose [4.5gm/l] (e.g., ICN-Flow Cat. No. 12-
141 332-54)
142 b) L-Glutamine 200 mM (e.g., ICN-Flow # 16-801-49)
143 c) New Born Calf Serum (NCS) (e.g., Biochrom # SO 125)
144 d) 0.05 % Trypsin/0.02 % EDTA solution (e.g., SIGMA T 3924, ICN-Flow, #
145 16891-49)
146 e) Phosphate buffered saline (PBS) without Ca^{2+} and Mg^{2+} (for trypsinization)
147 f) Hanks' Balanced Salt Solution (HBSS) without Ca^{2+} and Mg^{2+} (CMF-
148 HBSS)
149 g) Dulbecco's Phosphate Buffered Saline (D-PBS) [formulation containing
150 calcium and magnesium cations; glucose optional] (for rinsing)
151 h) Penicillin/streptomycin solution (e.g. ICN-Flow # 16-700-49)
152 i) Neutral Red (NR) Dye – tissue culture-grade; liquid form (e.g., SIGMA N
153 2889); powder form (e.g., SIGMA N 4638)
154 j) DMSO, U.S.P. analytical grade (Store under nitrogen @ -20°C)
155 k) ETOH, U.S.P. analytical grade (100 %, non-denatured for test substance
156 preparation; 95 % can be used for the desorb solution)
157 l) Glacial acetic acid, analytical grade
158 m) Distilled H_2O or any purified water suitable for cell culture and NR desorb
159 solution (sterile)
160 n) Sterile/non-sterile paper towels (for blotting 96-well plates)
161

162 [Note: Due to lot variability of NCS, first check a lot for growth stimulating
163 properties with 3T3 cells (approximately 20-24 hours doubling time) and then
164 reserve a sufficient amount of NCS.]
165

166 B. Preparations of Media and Solutions

167 [Note: All solutions (except NR stock solution, NR medium and NR desorb),
168 glassware, pipettes, etc., shall be sterile and all procedures should be carried out
169 under aseptic conditions and in the sterile environment of a laminar flow cabinet
170 (biological hazard standard). All methods and procedures will be adequately
171 documented.]
172

173 1. Media

174 DMEM (buffered with sodium bicarbonate) supplemented with (final
175 concentrations in DMEM are quoted):
176

- 177 a) Freeze Medium; contains 2X concentration of NCS and DMSO of final
 178 freezing solution
 179 40 % NCS
 180 20 % DMSO
 181
 182 b) Routine Culture Medium
 183 10 % NCS
 184 4 mM Glutamine
 185
 186 c) Chemical Dilution Medium
 187 4 mM Glutamine
 188 200 IU/mL Penicillin
 189 200 µg/mL Streptomycin
 190
 191 d) NR Dilution Medium
 192 5 % NCS
 193 4 mM Glutamine
 194 100 IU/mL Penicillin
 195 100 µg/mL Streptomycin
 196

197 [Note: The Chemical Dilution Medium with test substance will dilute the serum
 198 concentration of the Routine Culture Medium in the test plate to 5 %. Serum
 199 proteins may mask the toxicity of the test substance, but serum cannot be totally
 200 excluded because cell growth is markedly reduced in its absence.]
 201

202 Completed media formulations should be kept at approximately 2-8°C and
 203 stored for no longer than two weeks.
 204

205 2. Neutral Red (NR) Stock Solution

206 The liquid tissue culture-grade stock NR Solution is the first choice (e.g.,
 207 SIGMA #N2889, 3.3 mg/mL). Store liquid tissue culture-grade NR Stock
 208 Solution at the storage conditions and shelf-life period recommended by the
 209 manufacturer.
 210

211 A stock solution can be made with powder NR dye and water (e.g., 0.25 g NR
 212 Dye powder in 100 mL H₂O) if the liquid stock form is not available. The stock
 213 should be stored in the dark at room temperature for up to two months.
 214

215 3. Neutral Red (NR) Medium

216 EXAMPLE:

217 0.758 mL (3.3 mg NR dye/mL sol.) NR Stock Solution
 218 99.242 mL NR Dilution Medium (pre-warmed to
 219 37°C)
 220

221 The final concentration of the NR Medium is **25 µg NR dye/mL** and aliquots will be
222 prepared on the day of application.

223

224 [Note: The NR Medium shall be filtered (e.g., Millipore filtering, 0.2 – 0.45 µm pore
225 size) to reduce NR crystals. Aliquots of the NR Medium should be maintained
226 at 37°C (e.g., in a waterbath) before adding to the cells and used within 60
227 minutes of preparation but also used within 15 minutes after removing from
228 37°C storage. Examine the solution for crystals.]

229

230

4. ETOH/Acetic Acid Solution (NR Desorb)

231 1 % Glacial acetic acid solution

232 50 % ETOH

233 49 % H₂O

234

235

C. Methods

236

1. Cell Maintenance and Culture Procedures

237 BALB/c 3T3 cells are routinely grown as a monolayer in tissue culture grade
238 flasks (e.g., 75 - 80 cm²) at 37°C ± 1°C, 90 % ± 10 % humidity, and 5.0 % ± 1
239 % CO₂/air. The cells should be examined on a daily (i.e., on workdays) basis
240 under a phase contrast microscope, and any changes in morphology or their
241 adhesive properties noted in a Study Workbook.

242

243

2. Receipt of Cryopreserved BALB/c 3T3 Cells

244 Upon receipt of cryopreserved BALB/c 3T3 cells, the vial(s) of cells shall be
245 stored in a liquid nitrogen freezer until needed.

246

247

3. Thawing Cells

248 Thaw cells by putting ampules into a waterbath at 37°C ± 1°C. Leave for as brief
249 a time as possible.

250

251 3 Resuspend the cells in pre-warmed Routine Culture Medium and
252 transfer into pre-warmed Routine Culture Medium in a tissue-culture
253 flask.

254 4 Incubate at 37°C ± 1°C, 90 % ± 10 % humidity, and 5.0 % ± 1 %
255 CO₂/air.

256 5 When the cells have attached to the bottom of the flask (within 4 to
257 24 hours), decant the supernatant and replace with fresh pre-warmed
258 (37°C) medium. Culture as described above.

259 6 Passage at least two times before using the cells in a cytotoxicity test.

260

261 A fresh batch of frozen cells from the stock lot of cells should be thawed out and cultured
262 approximately every two months. This period resembles a sequence of about 18

263 passages.

264

265

4. Routine Culture of BALB/C 3T3 Cells

266

When cells exceed 50 % confluence (but less than 80 % confluent) they should be removed from the flask by trypsinization:

267

268

269

3 Decant medium, briefly rinse cultures with 5 mL PBS or Hanks' BSS (without Ca^{2+} , Mg^{2+}) per 25 cm^2 flask (15 mL per 75 cm^2 flask). Wash cells by gentle agitation to remove any remaining serum that might inhibit the action of the trypsin.

270

271

272

273

4 Discard the washing solution. Repeat the rinsing procedure and discard the washing solution.

274

275

5 Add 1-2 mL trypsin-EDTA solution per 25 cm^2 to the monolayer for a few seconds (e.g., 15-30 seconds).

276

277

6 Remove excess trypsin-EDTA solution and incubate the cells at room temperature.

278

279

7 After 2-3 minutes, lightly tap the flask to detach the cells into a single cell suspension.

280

281

282

5. Cell Counting

283

After detaching the cells, add 0.1-0.2 mL of pre-warmed (37°C) Routine Culture Medium/ cm^2 to the flask (e.g., 2.5 mL for a 25 cm^2 flask). Disperse the monolayer by gentle trituration to obtain a single cell suspension for exact counting. Count a sample of the cell suspension obtained using a hemocytometer or cell counter (e.g., Coulter counter).

284

285

286

287

288

289

6. Subculture of Cells

290

After determination of cell number, the culture can be sub-cultured into other flasks or seeded into 96-well microtiter plates (see **Section VI.E.1** for 96-well test plate configuration). BALB/c 3T3 cells are routinely passaged at suggested cell densities as listed in the table (approximate doubling time is 20-24 hours). Laboratories must determine and adjust the final density to achieve appropriate growth.

291

292

293

294

295

296

297

Table 1. Cell Density Guidelines for Subculturing

298

Days in Culture	Seeding Density (cells/ cm^2)	Total Cells per 25 cm^2 flask	Total Cells per 75 cm^2 flask
2	16800	4.2×10^5	1.26×10^6
3	8400	2.1×10^5	6.3×10^5
4	4200	1.05×10^5	3.15×10^5

299

300

[Note: It is important that cells have overcome the lag growth phase when they are used for the test.]

301

302

303 **7. Freezing Cells** (procedure required only if current stock of cells is depleted)

304 Stocks of BALB/c 3T3 cells can be stored in sterile, freezing tubes in a
305 liquid nitrogen freezer. DMSO is used as a cryoprotective agent.

306

307 4 Centrifuge trypsinized cells at approximately 200 x g.

308 5 Suspend the cells in cold Routine Culture Medium (half the final
309 freezing volume) so a final concentration of $1-5 \times 10^6$ cells/mL can be
310 attained.

311 6 Slowly add cold Freeze Medium to the cells so that the solvent will
312 equilibrate across the cell membranes. Bring the cell suspension to the
313 final freezing volume. The final cell suspension will be 10 % DMSO.
314 Aliquot the cell suspension into freezing tubes and fill to 1.8 mL.

315 7 Place the tubes into an insulated container (e.g., styrofoam trays) and
316 place in a freezer (-70 to -80°C) for 24 hours (~freezing rate of
317 1°C/minutes). The laboratory needs to ensure that the freezing protocol
318 is applicable to the 3T3 cells and that the cells are viable when removed
319 from cryopreservation.

320 8 Place the frozen tubes into liquid nitrogen for storage.

321

322

8. Preparation of Cells for Assays

323 • Cultured cells that will be used in seeding the 96-well plates should be fed
324 fresh medium the day before subculturing to the plates. On the day of plate
325 seeding, prepare a cell suspension of $2.0 - 3.0 \times 10^4$ cells/mL in Routine
326 Culture Medium. Using a multi-channel pipette, dispense 100 μ l Routine
327 Culture Medium only into the peripheral wells (blanks) of a 96-well tissue
328 culture microtiter plate (See **Section VI.E.1**). In the remaining wells,
329 dispense 100 μ l of a cell suspension of $2.0 - 3.0 \times 10^4$ cells/mL (= $2.0 -$
330 3.0×10^3 cells/well). The seeding density should be noted to ensure that the
331 cells in the control wells are not overgrown after three days (i.e., 24 hour
332 incubation in step **b** and 48 hour exposure to test substances). Prepare one
333 plate per substance to be tested.

334 • Incubate cells for 24 hours \pm 2 hours ($37^\circ\text{C} \pm 1^\circ\text{C}$, 90 % \pm 10 % humidity,
335 5.0 % \pm 1 % CO_2 /air) so that cells form a less than half (< 50%) confluent
336 monolayer. This incubation period assures cell recovery and adherence and
337 progression to exponential growth phase.

338 • Examine each plate under a phase contrast microscope to assure that cell
339 growth is relatively even across the microtiter plate. This check is
340 performed to identify experimental and systemic cell seeding errors. Record
341 observations in the Study Workbook.

342

343

9. Determination of Doubling Time

344 a) Establish cells in culture and trypsinize cells as per **Section VI.C.4** for
345 subculture. Resuspend cells in NR Dilution Medium (5 % NBCS/NCS).
346 Seed cells at 4200 cells/cm².

- 347 b) Seed five sets of cell culture vessels in triplicate for each cell type (e.g., 15
348 tissue culture dishes [60mm x 15mm]). Use appropriate volume of culture
349 medium for the culture vessels. Note number of cells placed into each
350 culture dish. Place dishes into the incubators ($37^{\circ}\text{C} \pm 1^{\circ}\text{C}$, $90\% \pm 10\%$
351 humidity, $5.0\% \pm 1\%$ CO_2/air).
- 352 c) After 4 - 6 hours (use the same initial measurement time for each subsequent
353 doubling time experiment), remove three culture dishes and trypsinize cells.
354 Count cells using a cell counter or hemocytometer. Study Director may
355 determine cell viability by dye exclusion (e.g., Trypan Blue; Nigrosin). Use
356 appropriate size exclusion limits if using a Coulter counter. Determine the
357 total number of cells and document. Repeat sampling at 24 hours, 48 hours,
358 72 hours, and 96 hours post inoculation. Change culture medium at 72
359 hours or sooner in remaining dishes if indicated by pH drop.
- 360 d) Plot cell concentration (per mL of medium) on a log scale against time on a
361 linear scale. Determine lag time and population doubling time. Additional
362 dishes and time are needed if the entire growth curve is to be determined
363 (lag phase, log phase, plateau phase).
- 364

365 **D. Preparation of Test Substances**

366 [Note: Preparation under red or yellow light is recommended to preserve substances
367 that degrade upon exposure to light.]

368

369 Test substance solubility should be determined by following the procedures
370 outlined in **ANNEX II** of this protocol.

371

372 **1. Test Substances in Solution**

- 373 a) Allow test substances to equilibrate to room temperature before dissolving
374 and diluting.
- 375 b) Prepare test substance immediately prior to use rather than preparing in bulk
376 for use in subsequent tests. Ideally, the solutions must not be cloudy nor
377 have noticeable precipitate. Each stock dilution should have at least 1-2 mL
378 total volume to ensure adequate solution for the test wells in a single 96-well
379 plate. The Study Director may store an aliquot (e.g., 1 mL) of the highest
380 2X stock solution (e.g., low solubility substances) in a freezer (e.g., -70°C)
381 for use in future substance analyses.
- 382 c) For substances dissolved in DMSO or ETOH, the final DMSO or ETOH
383 concentration for application to the cells must be 0.5 % (v/v) in the vehicle
384 controls and in all of the eight test concentrations.
- 385 d) The stock solution for each test substance should be prepared at the highest
386 concentration found to be soluble in the solubility test conducted per
387 **ANNEX II: Test Method Procedure - Solubility Determination of Test**
388 **Substances**. Thus, the highest test concentration applied to the cells in each
389 range finding experiment is:
- 390 4 0.5 times the highest concentration found to be soluble in the solubility
391 test, if the substance was soluble in Chemical Dilution Medium, or

- 392 5 1/200 the highest concentration found to be soluble in the solubility test
393 if the substance was soluble in ETOH or DMSO.
394 e) The seven lower concentrations in the range finding experiment are
395 prepared by successive dilutions that decrease by one log unit each. The
396 following example illustrates the preparation of test substance in solvent and
397 the dilution of dissolved test substance in Chemical Dilution Medium before
398 application to 3T3 cells.
399

400 Example: Preparation of Test Substance in Solvent Using a Log Dilution

401 Scheme

402 If DMSO is determined to be the preferred solvent at Tier 3 of the solubility test
403 (i.e., 200,000 µg/mL), dissolve the substance in DMSO at 200,000 µg/mL for
404 the chemical stock solution.

- 405
- 406 1) Label eight tubes 1 – 8. Add 0.9 mL solvent (e.g., DMSO) to tubes 2 -- 8.
 - 407 2) Prepare stock solution of 200,000 µg test substance/mL solvent in tube # 1.
 - 408 3) Add 0.1 mL of 200,000 µg/mL dilution from tube #1 to tube #2 to make a
409 1:10 dilution in solvent (i.e., 20,000 µg/mL).
 - 410 4) Add 0.1 mL of 20,000 µg/mL dilution from tube #2 to tube #3 to make
411 another 1:10 dilution (i.e., 1:100 dilution from stock solution) in solvent
412 (i.e., 2,000 µg/mL)
 - 413 5) Continuing making serial 1:10 dilutions in the prepared solvent tubes.
 - 414 6) Since each concentration is 200 fold greater than the concentration to be
415 tested, make a 1:100 dilution by diluting 1 part dissolved test substance in
416 each tube with 99 parts of Chemical Dilution Medium (e.g., 0.1 mL test
417 substance in DMSO + 9.9 mL Chemical Dilution Medium) to derive the
418 eight 2X concentrations for application to 3T3 cells. Each 2X test substance
419 concentration will then contain 1 % v/v solvent. The 3T3 cells will have
420 0.05 mL Routine Culture Medium in the wells prior to application of the test
421 substance. By adding 0.05 mL of the appropriate 2X test substance
422 concentration to the appropriate wells, the test substance will be diluted
423 appropriately (e.g., highest concentration in well will be 1,000 µg/mL) in a
424 total of 0.1 mL and the solvent concentration in the wells will be 0.5% v/v.
 - 425 7) A test substance prepared in Chemical Dilution Medium, DMSO, or ETOH
426 may precipitate upon transfer into the Routine Culture Medium. The 2X
427 dosing solutions should be evaluated for precipitates and the results
428 recorded in the Study Workbook. It is permissible to test all of the dosing
429 solutions in the dose range finding assay and main experiments. However,
430 doses containing test substance precipitates should be avoided because it
431 creates doubt about the concentration of test substance exposed to the cells.
432

433 Document all test substance preparations in the Study Workbook.
434

435 **2. pH of Test Substance Solutions**

436 Prior to or immediately after application of the test substance to the 96-well
437 plate, measure the pH of the highest 2X dosing concentration of the test
438 substance (i.e., C1 in the test plate, see **Figure 1**) in culture medium. Use pH
439 paper (e.g., pH 0 - 14 to estimate and pH 5 - 10 to determine more precise
440 value; or Study Director's discretion) for measurements. The pH paper should
441 be in contact with the solution for approximately one minute. Document the pH
442 and note the color of the 2X concentration medium (i.e., in the Microsoft Excel
443 ® template; see **ANNEX 1** for an example template). Medium color for all
444 dosing dilutions should be noted in the Study Workbook. Do not adjust the pH.
445


446 **3. Concentrations of Test Substance**

447 a) Range Finder Experiment

- 448 • Test eight concentrations of the test substance by diluting the stock
449 solution using log dilutions (e.g., 1:10, 1:100, 1:1000, etc.).
450 • If a range finder experiment does not generate enough cytotoxicity, then
451 higher doses should be attempted. If cytotoxicity is limited by
452 solubility, then more stringent solubility procedures to increase the stock
453 concentration (to the maximum concentration specified in **Section**
454 **VI.D.3.b.**) should be employed.
455 • Place the test substance concentration into an incubator (37°C ± 1°C, 90
456 % ± 10 % humidity, 5.0 % ± 1 % CO₂/air) and stir or rock for up to 3
457 hours, if necessary, to facilitate dissolution. For stocks prepared in
458 medium, vessel caps should be loose to allow for CO₂ exchange.
459 Proceed with dosing solution preparation and dosing.
460 • If a range finding test produces a biphasic curve, then the doses selected
461 for the subsequent main experiments should cover the most toxic dose-
462 response range (see **Example 1** – the most toxic range is 0.001 – 0.1
463 µg/mL) that reduces viability to 50 %.

464
465
466

Example 1 – Biphasic Curve



467
468

469

b) Main Experiment (Definitive Assay)

470

- Depending on the slope of the concentration-response curve estimated from the range finder, the dilution/progression factor in the concentration series of the main experiment should be smaller (e.g., dilution factor of $\sqrt[6]{10} = 1.47$). Cover the relevant concentration range around the IC_{50} (> 0 % and < 100 % effect) preferably with several points of a graded effect, but with a minimum of two points, one on each side of the estimated IC_{50} value, avoiding too many non-cytotoxic and/or 100 %-cytotoxic concentrations.
- Determine which test substance concentration is closest to the IC_{50} value. Use that value as a central concentration and adjust dilutions higher and lower in equal steps for the definitive assay.

471

472

473

474

475

476

477

478

479

480

481

482

Maximum Doses to be Tested in the Main Experiments

483

If minimal or no cytotoxicity was measured in the dose range finding assay, a maximum dose for the main experiments will be established as follows:

484

485

486

For test substances prepared in Chemical Dilution Medium

487

- The highest test substance concentration that may be applied to the cells in the main experiments will be either 100 mg/mL, or the maximum soluble dose.
- Test substance will be weighed into a glass tube and the weight will be documented. A volume of Chemical Dilution Medium will be added to the vessel so that the concentration is 200,000 µg/mL (200 mg/mL).
- The solution is mixed using the mechanical procedures that produced solubility when performing the solubility test (See ANNEX II). If complete solubility is achieved in medium, then 7 additional serial stock dosing solutions may be prepared from the 200 mg/mL 2X stock. If the test substance is insoluble in medium at 200 mg/ml, proceed by adding medium, in small incremental amounts, to attempt to dissolve the substance by using the sequence of mechanical procedures specified in ANNEX II.

488

489

490

491

492

493

494

495

496

497

498

499

500

- 501 • More stringent solubility procedures may be employed if needed based
 502 on results from the range finder experiment (**Section VI.D.3.a.**). The
 503 highest soluble stock solution will be used to prepare the 7 additional
 504 serial stock dosing solutions.

505

506 For test substances prepared in either DMSO or ETOH

- 507 • The highest test substance concentration that may be applied to the cells
 508 in the main experiments will be ≤ 2.5 mg/mL or less, depending upon
 509 the maximum solubility in solvent.
- 510 • Weigh the test substance into a glass tube and document the weight.
 511 Add the appropriate solvent (determined from the original solubility test)
 512 to the vessel so that the concentration is 500,000 $\mu\text{g/mL}$ (500 mg/mL).
- 513 • Mix the solution using the sequence of mechanical procedures specified
 514 in **ANNEX II**. If complete solubility is achieved in the solvent, then
 515 seven additional serial stock dosing solutions may be prepared from the
 516 500 mg/mL 200X stock. If the test substance is insoluble in solvent at
 517 500 mg/ml, proceed by adding solvent, in small incremental amounts, to
 518 attempt to dissolve the substance by again using the sequence of mixing
 519 procedures. The highest soluble stock solution will be used to prepare
 520 the seven additional serial stock dosing solutions.

521

522 If precipitates are observed in the 2X dilutions, continue with the
 523 experiment and make the appropriate observations and documentation.

524

525

c) Test Substance Dilutions

526 The dosing factor of 3.16 ($=\sqrt[2]{10}$) divides a log into two equidistant steps,
 527 2.15 ($=\sqrt[3]{10}$) into three steps, 1.47 ($=\sqrt[6]{10}$) into six steps, 1.78 ($=\sqrt[4]{10}$) into
 528 four steps, and 1.21 ($=\sqrt[12]{10}$) into 12 steps.

529

530

EXAMPLE:

531

10						31.6						100
10				21.5				46.4				100
10		14.7		21.5		31.6		46.4		68.1		100
10	12.1	14.7	17.8	21.5	26.1	31.6	38.3	46.4	56.2	68.1	82.5	100

532

533

An example of decimal geometric concentration series for factor 1.47:
*Dilute 1 volume of the highest concentration by adding 0.47 volumes of
 534 diluent. After equilibration, dilute 1 volume of this solution by adding 0.47
 535 volumes of diluent...(etc.).*

536

537

538

538
539
540
541
542
543
544
545

E. Test Procedure

1. 96-Well Plate Configuration

The 3T3 NRU assay for test substances will use the 96-well plate configuration as shown in **Figure 1**.

Figure 1. 96-Well Plate Configuration for Positive Control (PC) and Test Substance Assays

	1	2	3	4	5	6	7	8	9	10	11	12
A	VCb	VCb	C _{1b}	C _{2b}	C _{3b}	C _{4b}	C _{5b}	C _{6b}	C _{7b}	C _{8b}	VCb	VCb
B	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
C	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
D	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
E	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
F	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
G	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
H	VCb	VCb	C _{1b}	C _{2b}	C _{3b}	C _{4b}	C _{5b}	C _{6b}	C _{7b}	C _{8b}	VCb	VCb

546
547
548
549
550
551
552
553

VC1 and VC2 = VEHICLE CONTROL
C₁ – C₈ = Test Substances or PC (SLS) at eight concentrations
(C₁ = highest, C₈ = lowest)
b = BLANKS (Test substance or PC, but contain **no** cells)
VCb = VEHICLE CONTROL BLANK (contain **no** cells)

2. Application of Test Substance

554
555
556
557
558
559
560
561
562
563
564
565
566
567

6.0 Two optional methods for rapidly applying the 2X dosing solutions onto the 96-well plates may be utilized.

- 1) Add each of the 2X dosing solutions into labeled, sterile reservoirs (e.g., Corning/Costar model 4870 sterile polystyrene 50 mL reagent reservoirs; or Corning/Transtar model 4878 disposable reservoir liners, 8-channel; or other multichannel reservoirs).
- 2) Use a “dummy” plate (i.e., an empty sterile 96-well plate) prepared to hold the dosing solutions immediately prior to treatment of the test plate (with cells). The test substance and control dosing solutions should be dispensed into the dummy plate in the same pattern/order as will be applied to the plate containing cells. More volume than needed for the test plate (i.e., greater than 50 µl/well) should be in the wells of the dummy plate.

568
569
570
571
572
573
574
575
576
577
578
579
580
581
582
583
584
585
586
587
588
589
590
591
592
593
594
595
596
597
598
599
600
601
602
603
604
605
606
607
608

At the time of treatment initiation, a multi-channel micropipettor is used to transfer the 2X dosing solutions, from the reservoirs or dummy plate, to the appropriate wells on the treatment plate (as described in step c. below). These methods will ensure that the dosing solutions can be transferred rapidly to the appropriate wells of the test plate to initiate treatment times and to minimize the range of treatment initiation times across a large number of treatment plates, and to prevent “out of order” dosing. Do not use a multichannel repeater pipette for dispensing test substance to the plates.

- 7.0 After 24 hours \pm 2 hours incubation of the cells, remove Routine Culture Medium from the cells by careful inversion of the plate (i.e., “dump”) over an appropriate receptacle. Gently blot the plate on a sterile paper towel so that the monolayer is minimally disrupted. Do not use automatic plate washers for this procedure nor vacuum aspiration.
- 8.0 Immediately add 50 μ L of fresh pre-warmed Routine Culture Medium to all of the wells, including the blanks. Fifty microliters (50 μ L) of dosing solution will be rapidly transferred from the 8-channel reservoir (or dummy plate) to the appropriate wells of the test plate using a single delivery multi-channel pipettor. For example, the VC may be transferred first (into columns 1, 2, 11, and 12), followed by the test substance dosing solutions from lowest to highest dose, so that the same pipette tips on the multi-channel pipettor can be used for the whole plate. [The Vehicle Control blank (VCb) wells (column 1, column 12, wells A2, A11, H2, H11) will receive the Vehicle Control dosing solutions (which should include any solvents used)]. Blanks for wells A3 – A10 and H3 – H10 shall receive the appropriate test substance solutions for each concentration (e.g., wells A3 and H3 receive C₁ solution).
- 9.0 Incubate cells for 48 hours \pm 0.5 hours (37°C \pm 1°C, 90 % \pm 10 % humidity, and 5.0 % \pm 1 % CO₂/air).
- 10.0 **Positive Control:** For each set of test substance plates used in an assay, prepare a separate plate of positive control concentrations. If multiple sets of test substance plates are set up, then clearly designate the positive control plates for each set; each set will be an individual entity. The Study Director will decide how many test substance plates will be run with a positive control plate. This plate will follow the same schedule and procedures as used for the test substance plates (including appropriate test substance concentrations in the appropriate wells and meeting test acceptance criteria – see **sections VI.E.1, E.2, and E.5**).

3. Microscopic Evaluation

609
610
611
612
613

After at least 46 hours of treatment, examine each plate under a phase contrast microscope to identify systematic cell seeding errors and growth characteristics of control and treated cells. Record any changes in morphology of the cells due to the cytotoxic effects of the test substance, but do not use these records for any quantitative measure of cytotoxicity. Undesirable growth characteristics of

614 control cells may indicate experimental error and may be cause for rejection of
 615 the assay. Use the following Visual Observations Codes in the description of
 616 cell culture conditions. Numerical scoring of the cells (see **Section VI.E.3**)
 617 should be determined and documented in the Study Workbook and in the
 618 appropriate section of the Microsoft Excel® template.

619
 620
 621

Visual Observations Codes

Note Code	Note Text
1	Normal Cell Morphology
2	Low Level of Cell Toxicity
3	Moderate Level of Cell Toxicity
4	High level of Cell Toxicity
1P	Normal Cell Morphology with Precipitate
2P	Low Level of Cell Toxicity with Precipitate
3P	Moderate Level of Cell Toxicity with Precipitate
4P	High level of Cell Toxicity with Precipitate
5P	Unable to View Cells Due to Precipitate

622
 623

4. Measurement of NRU

- 624 a) Carefully remove (i.e., “dump”) the medium with test substance and rinse
 625 the cells very carefully with 250 µL pre-warmed D-PBS. Remove the
 626 rinsing solution by dumping and remove excess by gently blotting on paper
 627 towels. Add 250 µL NR medium (to all wells including the blanks) and
 628 incubate (37°C ± 1°C, 90 % ± 10 % humidity, and 5.0 % ± 1 % CO₂/air) for
 629 3 hours ± 0.1 hour. Observe the cells briefly during the NR incubation (e.g.,
 630 between 2 and 3 hours – Study Director’s discretion) for NR crystal
 631 formation. Record observations in the Study Workbook. Study Director
 632 can decide to reject the experiment if excessive NR crystallization has
 633 occurred.
- 634 b) After incubation, remove the NR medium, and carefully rinse cells with 250
 635 µl pre-warmed D-PBS.
- 636 c) Decant and blot D-PBS from the plate.
- 637 d) Add exactly 100 µl NR Desorb (ETOH/acetic acid) solution to all wells,
 638 including blanks.
- 639 e) Shake microtiter plate rapidly on a microtiter plate shaker for 20 – 45
 640 minutes to extract NR from the cells and form a homogeneous solution.
 641 Plates should be protected from light by using a cover during shaking.
- 642 f) Plates should be still for at least five minutes after removal from the plate
 643 shaker (or orbital mixer). If any bubbles are observed, assure that they have
 644 been ruptured prior to reading the plate. Measure the absorption (within 60
 645 minutes of adding NR Desorb solution) of the resulting colored solution at
 646 540 nm ± 10 nm in a microtiter plate reader (spectrophotometer), using the
 647 blanks as a reference.

648 [Note: A mean $OD_{540 \pm 10nm}$ of 0.031 - 0.065 for the VC blanks is a target
649 range of ODs but not a test acceptance criterion (range = mean OD \pm 2.5
650 standard deviations; mean = 0.048; SD = 0.007; N = 233).] Save raw data in
651 the Microsoft Excel® template.
652

653 5. Quality Check of 3T3 NRU Assay

654 a) Test Acceptance Criteria

655 All acceptance criteria (i.e., criteria 1, 2, and 3) must be met for a test to be
656 acceptable.
657

658 1) The PC (SLS) IC_{50} must be within \pm two and a half (2.5) standard
659 deviations (SD) of the historical mean established by the Test Facility,
660 and must meet criteria 2 and 3, and must have an r^2 (coefficient of
661 determination) value calculated for the Hill model fit (i.e., from
662 PRISM® software) \geq 0.85. NICEATM/ECVAM study generated the
663 following PC data:

- 664 • IC_{50} mean = 41.5 μ g/mL; SD = 4.8 (n = 233)
- 665 • range for IC_{50} mean \pm 2.5 SD = 29.5 μ g/mL – 53.5 μ g/mL

666 2) The left and right mean of the VCs do not differ by more than 15% from
667 the mean of all VCs.

- 668 • At least one calculated cytotoxicity value $>$ 0 % and \leq 50 % viability
669 and at least one calculated cytotoxicity value $>$ 50 % and $<$ 100 %
670 viability must be present.
671

672 *Exception:* If a test has only one point between 0 and 100 % **and** the
673 smallest dilution factor (i.e., 1.21) was used **and** all other test acceptance
674 criteria were met, then the test will be considered acceptable.
675

676 **Stopping Rule for Insoluble Substances:** If the most rigorous solubility
677 procedures have been performed and the assay cannot achieve adequate
678 toxicity to meet the test acceptance criteria after three definitive trials, then
679 the Study Director may end all testing for that particular substance.
680

681 [Note: A corrected mean $OD_{540 \pm 10nm}$ of 0.183 - 0.769 for the VCs is a target
682 range of ODs but not a test acceptance criterion (range = mean OD \pm 2.5
683 standard deviations; mean = 0.476; SD = 0.117; N = 233).]
684

685 b) Checks for Systematic Cell Seeding Errors

686 To check for systematic cell seeding errors, untreated VCs are placed both at
687 the left side (row 2) and the right side (row 11 for the test plates) of the 96-
688 well plate. Aberrations in the cell monolayer for the VCs may reflect a
689 volatile and toxic test substance present in the assay. If volatility is
690 suspected, then proceed to **Section VI.E.6**. Checks for cell seeding errors
691 may also be performed by examining each plate under a phase contrast
692 microscope to assure that cell quantity is consistent.

693
694
695
696
697
698
699
700
701
702
703
704
705
706
707
708
709
710
711
712
713
714
715
716
717
718
719
720
721
722
723
724
725
726
727
728
729
730
731
732
733
734
735
736
737
738

6. Testing Volatile Substances

Although this test method is not suitable for highly volatile substances, mildly volatile substances may be tested with some success. Volatile test substances may generate vapors from the treatment medium during the test substance treatment incubation period. These vapors may become resorbed into the treatment medium in adjacent wells, such that culture wells nearest the highest doses may become contaminated by exposure. If the test substance is particularly toxic at the doses tested, the cross contamination may be evident as a significant reduction in viability in the VC cultures (i.e., VC1) adjacent to the highest test substance doses.

If potential test substance volatility is suspected (e.g., for low density liquids) or if the initial range finder test (non-sealed plate) results show evidence of toxic effects in the control cultures (i.e., > 15 % difference in viability between VC1 [column 2] and VC2 [column 11]), then seal the subsequent test plates using the following procedure.

- a) Plates and substances will be prepared as usual according to **Sections VI.D and VI.E.**
- b) Immediately after the 96-well culture plate has been treated with the suspected volatile substance (**Section VI.E.2.b**), apply the adhesive plate sealer (e.g., using a hand, microplate roller, etc.) directly over the culture wells. Assure that the sealer adheres to each culture well (well tops should be dry). Place the 96-well plate cover over the sealed plate and incubate the plate under specified conditions (**Section VI.E.2.b**). [Note: Do not jam the plate lid over the film to avoid deforming the sealer and causing the sealer to detach from culture wells. Loose fit of the plate lid is acceptable.]
- c) At the end of the treatment period, the plate sealer should be carefully removed to avoid spillage. Continue with the NRU assay as per **Section VI.E.4.**

F. Data Analysis

- The Study Director will use good biological/scientific judgment for determining “unusable” wells that will be excluded from the data analysis and provide explanations for the removal of any data from the analysis.
- A calculation of cell viability expressed as NRU is made for each concentration of the test substance by using the mean NRU of the six replicate values (minimum of four acceptable replicate well) per test concentration (blanks will be subtracted). This value is compared with the mean NRU of all VC values. Relative cell viability is then expressed as percent of untreated VC. If achievable, the eight concentrations of each substance tested will span the range of no effect up to total inhibition of cell viability.
- Data from the microtiter plate reader should be transferred to a spreadsheet template (e.g., Microsoft Excel®) that will automatically determine cell

- 739 viability, calculate IC₅₀ values by linear interpolation, and perform statistical
 740 analyses (including statistical identification of outliers) (see ANNEX 1 for an
 741 example spreadsheet template).
- 742 • A Hill function analysis should be performed using statistical software (e.g.,
 743 GraphPad PRISM® 3.0) and a template to calculate IC₂₀, IC₅₀, and IC₈₀ values
 744 (and the associated confidence limits) for each test substance.
 - 745 • Dose-responses for which the toxicity plateaus as concentration increases do not
 746 fit the Hill function well when Bottom =0. To obtain a better model fit,
 747 unconstrain the Bottom parameter so that the model calculates the Bottom
 748 value. However, when Bottom ≠ 0, the EC₅₀ reported by the Hill function ≠
 749 50% viability since the Hill function defines EC₅₀ as the point midway between
 750 Top and Bottom. To obtain the appropriate IC₅₀ when Bottom ≠ 0, use the
 751 following rearranged Hill function:

$$X = \log EC_{50} - \frac{\log\left(\frac{Top - Bottom}{Y - Bottom} - 1\right)}{HillSlope}$$

- 753
- 754 • X is the logarithm of concentration at 50% response, logEC₅₀ is logarithm of
 755 concentration at the response midway between Top and Bottom, Top is the
 756 maximum response, Bottom is the minimum response, Y = 50 (i.e., 50%
 757 response), and HillSlope describes the steepness of the curve.

758

759 [Note: IC₅₀ values are used in a regression formula to predict the LD₅₀ value of a
 760 test substance as an estimate of the starting dose for an acute oral toxicity test.]

761 VII. REFERENCES

- 762
- 763 Hackenberg, U. and H. Bartling. 1959. Messen und Rechnen im pharmakologischen
 764 Laboratorium mit einem speziellen Zahlensystem (WL24-System). Arch. Exp. Pathol.
 765 Pharmakol. 235: 437-463.
- 766
- 767 Spielmann, H., S. Gerner, S. Kalweit, R. Moog, T. Wirmserberger, K. Krauser, R.
 768 Kreiling, H. Kreuzer, N.P. Luepke, H.G. Miltenburger, N. Müller, P. Murmann, W.
 769 Pape, B. Siegmund, J. Spengler, W. Steiling, and F.J. Wiebel. 1991. Interlaboratory
 770 assessment of alternatives to the Draize eye irritation test in Germany. Toxicol. In
 771 Vitro 5: 539-542.
- 772
- 773 *Test Method Protocol for Solubility Determination. In Vitro Cytotoxicity Validation*
 774 *Study. Phase III. August 29, 2003. Prepared by The National Toxicology Program*
 775 *(NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods*
 776 *(NICEATM). (<http://iccvam.niehs.nih.gov/methods/invitro.htm>)*
 777
- 778
- 779 U. S. Environmental Protection Agency. 1996. Product Properties Test Guidelines.
 780 OPPTS 803.7840. Water Solubility: Column Elution Method; Shake Flask Method.
 781 EPA712-C-96-041, Prevention, Pesticides and Toxic Substances, Washington DC.

782
783
784
785
786
787
788

Test Method Protocol for the 3T3 Neutral Red Uptake Cytotoxicity Test. A Test for Basal Cytotoxicity for an In Vitro Validation Study. November 4, 2004. Prepared by the National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods (NICEATM). (<http://iccvam.niehs.nih.gov/methods/invitro.htm>)

788
789
790

ANNEX 1
Microsoft EXCEL® Example Spreadsheet Template

Test Facility :	A				Study Number.:	A1							
Chemical Code :	SLS				96-Well Plate ID :	A11							
2nd Chem. Code*:	11				Experiment ID :	XX							
96-WELL PLATE MAP													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	
B	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
C	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
D	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
E	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
F	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
G	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank	
H	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	
RAW ABSORBANCE DATA (OD₅₅₀)													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	0.044	0.044	0.045	0.045	0.045	0.046	0.051	0.057	0.057	0.043	0.041	0.044	
B	0.042	0.456	0.043	0.043	0.130	0.300	0.395	0.414	0.418	0.402	0.401	0.042	
C	0.043	0.407	0.042	0.041	0.130	0.294	0.383	0.382	0.413	0.375	0.385	0.044	
D	0.043	0.438	0.042	0.043	0.147	0.337	0.409	0.404	0.438	0.436	0.391	0.047	
E	0.044	0.448	0.041	0.045	0.132	0.321	0.429	0.414	0.416	0.420	0.441	0.042	
F	0.045	0.411	0.040	0.042	0.127	0.375	0.397	0.402	0.422	0.447	0.403	0.043	
G	0.041	0.405	0.043	0.040	0.124	0.361	0.444	0.442	0.425	0.448	0.405	0.044	
H	0.041	0.041	0.048	0.042	0.042	0.044	0.042	0.042	0.040	0.044	0.041	0.041	
Max	0.045	0.456	0.043	0.045	0.147	0.375	0.444	0.442	0.438	0.448	0.441	0.047	
Min	0.041	0.405	0.040	0.040	0.124	0.294	0.383	0.382	0.413	0.375	0.385	0.041	
Next Max	0.044	0.448	0.042	0.043	0.132	0.361	0.429	0.414	0.425	0.447	0.405	0.044	
Next Min	0.042	0.407	0.041	0.041	0.127	0.300	0.395	0.402	0.416	0.402	0.391	0.042	
Rmax	-0.250	-0.157	-0.333	-0.400	-0.652	-0.173	-0.246	-0.467	-0.520	-0.014	-0.643	-0.500	
Rmin	0.250	0.039	0.333	0.200	0.130	0.074	0.197	0.333	0.120	0.370	0.107	0.167	
CORRECTED ABSORBANCE (Sample OD₅₅₀ - Mean Blank OD₅₅₀)													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	0.001	0.001	-0.002	0.002	0.002	0.001	0.005	0.008	0.009	-0.001	-0.002	0.001	
B	-0.001	0.413	-0.004	-0.001	0.087	0.255	0.349	0.365	0.370	0.359	0.358	-0.001	
C	0.000	0.364	-0.005	-0.003	0.087	0.249	0.337	0.333	0.365	0.332	0.342	0.001	
D	0.000	0.395	-0.005	-0.001	0.104	0.292	0.363	0.355	0.390	0.393	0.348	0.004	
E	0.001	0.405	-0.006	0.002	0.089	0.276	0.383	0.365	0.368	0.377	0.398	-0.001	
F	0.002	0.368	-0.007	-0.001	0.084	0.330	0.351	0.353	0.374	0.404	0.360	0.000	
G	-0.002	0.362	-0.004	-0.004	0.081	0.316	0.398	0.393	0.377	0.405	0.362	0.001	
H	-0.002	-0.002	0.002	-0.001	-0.001	-0.001	-0.005	-0.008	-0.009	0.001	-0.002	-0.002	
Mean Blank =	0.043		0.047	0.044	0.044	0.045	0.047	0.050	0.049	0.044			
RELATIVE VIABILITY (% OF VEHICLE CONTROL)													
	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B		110.7%	-0.9%	-0.1%	23.2%	68.4%	93.4%	97.7%	99.0%	96.1%	96.0%		
C		97.6%	-1.2%	-0.7%	23.2%	66.7%	90.2%	89.1%	97.7%	88.9%	91.7%		
D		105.9%	-1.2%	-0.1%	27.7%	78.3%	97.2%	95.0%	104.4%	105.2%	93.3%		
E		108.6%	-1.5%	0.4%	23.7%	74.0%	102.5%	97.7%	98.5%	100.9%	106.7%		
F		98.7%	-1.7%	-0.4%	22.4%	88.5%	94.0%	94.5%	100.1%	108.2%	96.5%		
G		97.1%	-0.9%	-0.9%	21.6%	84.7%	106.5%	105.2%	100.9%	108.4%	97.1%		
H													

791

802
803

TEST CHEMICAL							
Test Facility : A				Study Number : A1			
Chemical Code : SLS				96-Well Plate ID : A11			
2 nd Chem. Code* : 11				Experiment ID : XX			
* Testing Facility Accession Code, if applicable							
PREPARATION OF TEST CHEMICAL							
Solvent:		Medium		Dilution factor:		1.4	
Solvent Conc. (% v/v) in dosing solutions :				N/A		Highest Stock Conc.: 20,000 µg/mL	
Aids used to dissolve : <input type="checkbox"/> Vortexing <input type="checkbox"/> sonication <input type="checkbox"/> heating to 37C							
pH (highest medium stock or 2X dosing solution) : 8.0							
Medium Clarity/Color (highest 2X dosing solution): clear red						If ppt, note lowest conc.:	
Concentration Series (µg/mL)							
C1	C2	C3	C4	C5	C6	C7	C8
100	71.4	51.0	36.4	26.0	18.6	13.3	9.49
Positive Control (SLS) 100 - 9.49 µg/mL							
CELL LINE/TYPE							
Name: BALB/c 3T3		Supplier: ATCC		Lot No. not provided			
Passage No.: 69		Passage No. in Assay: 75		Proliferating/frozen 24-May-02			
CELL CULTURE CONDITIONS							
Medium: DMEM		Supplier:		Lot No.:			
Serum: NCS		Supplier:		Lot No.:			
Serum Conc.:		Growth Medium: 10%		Treatment Medium: 0%			
TEST ACCEPTANCE CRITERIA							
No. of values >50% and <100%: 3		No. of values >0% and <50%: 1		Accept? YES			
VC: % Difference between Col 2 and mean VC.:				-3%		Accept? YES	
PC: Hill Function R ² Value of SLS: 0.99						Accept? YES	
PC: IC ₅₀ of SLS: 43.2 µg/mL						Accept? YES	
TIMELINE							
Cell Seeding Date		Dose Application Date		OD ₅₅₀ Determination Date			
TEST RESULTS							
VC: Mean Corrected OD ₅₅₀ : 0.373				Hill Function R ² Value: 0.9869			
log IC ₂₀ : 1.551E+00 µg/mL		log IC ₅₀ : 1.635E+00 µg/mL		log IC ₈₀ : 1.718E+00 µg/mL			
IC ₂₀ : 3.56E+01 µg/mL		IC ₅₀ : 4.32E+01 µg/mL		IC ₈₀ : 5.22E+01 µg/mL			
Test Chemical F.W. : 288.4							
IC ₂₀ : 0.12331183 mM		IC ₅₀ : 0.1496252 mM		IC ₈₀ : 0.18113599 mM			

804
805
806
807
808

ANNEX II

TEST METHOD PROCEDURE
Solubility Determination of Test Substances

PROPOSAL

This procedure was designed to identify the solvent that would provide the highest soluble concentration of a test substance so there would be uniform availability of the substance to cells used for *in vitro* basal cytotoxicity testing. The solubility exercises can be performed in a routine and repeatable manner and provide guidelines to effectively prepare test substances for toxicity testing in the NRU test methods.

TEST SYSTEM

The solubility test procedure is based on attempting to dissolve substances in various solvents with increasingly rigorous mechanical techniques. The solvents to be used, in the order of preference, are cell culture medium, DMSO, and ETOH. Determination of whether a test substance has dissolved is based entirely on visual observation for the purposes of this protocol. A test substance has dissolved if the solution is clear and shows no signs of cloudiness or precipitation.

PROCEDURES

Preparation of the 3T3 medium will follow all procedures in the 3T3 NRU protocol.

Materials – see Section VI.A

Preparations of Media and Solutions – see Section VI.B

All solutions glassware, pipettes, etc., should be sterile and all procedures should be carried out under aseptic conditions and in the sterile environment of a laminar flow cabinet (biological hazard standard). All methods and procedures should be adequately documented.

Determination of Solubility

- Solubility should be determined in a step-wise procedure that involves attempting to dissolve a test substance at a relatively high concentration with the sequence of mechanical procedures specified in **Mechanical Procedures**. **Table 1** and **Figures 1** and **2** illustrate the step-wise procedures. The hierarchy of preference of solvent for dissolving test substances is medium, DMSO, and then ETOH. If the substance does not dissolve in the solvent, the volume of solvent is increased so as to decrease the test substance concentration by a factor of 10, and then the sequence of mechanical procedures are repeated in an attempt to solubilize the substance at the lower concentrations.
- For testing solubility in medium, the starting concentration is 200,000 µg/ml (i.e., 200 mg/mL) in Tier 1, but for DMSO and ETOH the starting concentration is 200,000 µg/ml (i.e., 200 mg/mL) in Tier 3.

853

854 *Methods*

- 855 A. Tier 1 begins with testing 200 mg/mL in Chemical Dilution Medium (see **Table 1**).
- 856 1. Weigh approximately 100 mg (100,000 µg) of the test substance into a glass tube.
- 857 Document the test substance weight.
- 858 2. Add approximately 0.5 mL of medium into the tube so that the concentration is
- 859 200,000 µg/ml (200 mg/mL).
- 860 3. Mix the solution as specified in **Mechanical Procedures**. If complete solubility is
- 861 achieved, then additional solubility procedures are not needed.
- 862
- 863 B. If the test substance is insoluble in Tier 1 at 200 mg/mL, then proceed to Tier 2.
- 864 1. Weigh approximately 10 mg (10,000 µg) of the test substance into a glass tube.
- 865 Document the substance weight.
- 866 2. Add approximately 0.5 mL of medium into the tube so that the concentration is
- 867 20,000 µg/ml (20 mg/mL).
- 868 3. Mix the solution as specified in **Mechanical Procedures**. If complete solubility is
- 869 achieved, then additional solubility procedures are not needed.
- 870
- 871 C. If the test substance is insoluble in Chemical Dilution Medium, proceed to Tier 3.
- 872 1. Add enough medium, approximately 4.5 mL, to attempt to dissolve the substance at 2
- 873 mg/mL by using the sequence of mixing procedures. If the test substance dissolves in
- 874 medium at 2 mg/mL, no further procedures are necessary.
- 875 2. If the test substance does NOT dissolve in medium, weigh out approximately 100 mg
- 876 test substance in a second glass tube and add enough DMSO to make the total volume
- 877 approximately 0.5 mL (for 200 mg/mL) and mix the solution as specified in
- 878 **Mechanical Procedures**.
- 879 3. If the test substance does not dissolve in DMSO, weigh out approximately 100 mg
- 880 test substance in another glass tube and add enough ETOH to make the total volume
- 881 approximately 0.5 mL (for 200 mg/mL) and mix the solution as specified in
- 882 **Mechanical Procedures**.
- 883 4. If the substance is soluble in either solvent, no additional solubility procedures are
- 884 needed.
- 885
- 886 D. If the substance is NOT soluble in Chemical Dilution Medium, DMSO, or ETOH at Tier
- 887 3, then continue to Tier 4 in **Table 1**.
- 888 1. Add enough solvent to increase the volume of the three (or four) Tier 2 solutions by
- 889 10 and attempt to solubilize again using the sequence of mixing procedures. If the
- 890 test substance dissolves, no additional solubility procedures are necessary.
- 891 2. If the test substance does NOT dissolve, continue with Tier 5 and, if necessary, Tier 6
- 892 using DMSO and ETOH.
- 893 3. Tier 5 begins by diluting the Tier 4 samples with DMSO or ETOH to bring the total
- 894 volume to 50 mL. The mixing procedures are again followed to attempt to solubilize
- 895 the substance.

896 4. Tier 6 is performed, if necessary, by weighing out another two samples of test
897 substance at ~10 mg each and adding ~50 mL DMSO or ETOH for a 200 µg/mL
898 solution, and following the mixing procedures.
899

900 Example

- 901 • If complete solubility is not achieved at 20,000 µg/mL in Chemical Dilution Medium at
902 Tier 2 using the mixing procedures, then the procedure continues to Tier 3 by diluting the
903 solution to 5 mL with medium and mixing again.
- 904 • If the substance is not soluble in Chemical Dilution Medium, two samples of ~ 100 mg
905 test substance are weighed to attempt to solubilize in DMSO and ETOH at 200,000
906 µg/mL (i.e., 200 mg/mL). Solutions are mixed following the sequence of procedures
907 prescribed in **Mechanical Procedures** in an attempt to dissolve.
- 908 • If solubility is not achieved at Tier 3, then the solutions prepared in Tier 3 are diluted by
909 10 so as to test 200 µg/mL in media, and 20,000 µg/mL in DMSO and ETOH. This
910 advances the procedure to Tier 4. Solutions are again mixed in an attempt to dissolve.
- 911 • If solubility is not achieved in Tier 4, the procedure continues to Tier 5, and to Tier 6 if
912 necessary (see **Figures 1 and 2** and **Table 1**).
913

914 **MECHANICAL PROCEDURES**

915 A. The following hierarchy of mixing procedures will be followed to dissolve the test
916 substance:

- 917
- 918 1. Add test substance to solvent as in Tier 1 of **Table 1**. (Test substance and solvent
919 should be at room temperature.)
- 920 2. Gently mix at room temperature. Vortex the tube (1 –2 minutes).
- 921 3. If test substance has not dissolved, use waterbath sonication for up to 5 minutes.
- 922 4. If test substance is not dissolved after sonication, then warm solution to 37°C for 5 -
923 60 minutes. This can be performed by warming tubes in a 37°C waterbath or in a
924 CO₂ incubator at 37°C. The solution may be stirred during warming (stirring in a
925 CO₂ incubator will help maintain proper pH).
- 926 5. Proceed to Tier 2 (and Tiers 3-6, if necessary of **Table 1** and repeat procedures 2-4).
927

928 B. The preference of solvent for dissolving test substances is Chemical Dilution Medium,
929 DMSO, and then ETOH. Thus, if all solvents for a particular tier are tested
930 simultaneously and a test substance dissolves in more than one solvent, then the choice of
931 solvent follows this hierarchy. For example, if, at any tier, a substance were soluble in
932 Chemical Dilution Medium and DMSO, the choice of solvent would be medium. If the
933 substance were insoluble in medium, but soluble in DMSO and ETOH, the choice of
934 solvent would be DMSO.
935

936

943

FIGURE 1 SOLUBILITY STEP-WISE (TIERED) PROCEDURE**TIER 1**

STEP 1:	200 mg/mL test substance (TS) in 0.5 mL Chemical Dilution Medium <ul style="list-style-type: none"> • if TS soluble in medium, then <u>STOP</u>. • if TS insoluble in medium, then go to STEP 2.
---------	---

TIER 2

STEP 2:	20 mg/mL TS in 0.5 mL Chemical Dilution Medium <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 3.
---------	---

TIER 3

STEP 3:	200 mg/mL TS in DMSO <ol style="list-style-type: none"> 4 if TS soluble, then <u>STOP</u>. 5 if TS insoluble, test at 200 mg/mL in ETOH. 6 if TS soluble, then <u>STOP</u>. 7 If TS insoluble, go to STEP 4.
---------	--

TIER 4

STEP 4:	0.2 mg/mL TS in medium (one or both) – increase volume from STEP 2 by 10 (i.e., to 50 mL) <ul style="list-style-type: none"> • if TS soluble in both media, then <u>STOP</u>. • if TS insoluble in one medium, test at 20 mg/mL in DMSO – increase volume from STEP 3 by 10 (i.e., to 5 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 20 mg/mL in ETOH – increase volume from STEP 3 by 10 (i.e., to 5 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 5.
---------	--

TIER 5

STEP 5:	2 mg/mL TS in DMSO – increase volume from STEP 4 by 10 (i.e., to 50 mL) <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 2 mg/mL in ETOH – increase volume from STEP 4 by 10 (i.e., to 50 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 6.
---------	---

TIER 6

STEP 6:	0.2 mg/mL TS in 50 mL DMSO <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 0.2 mg/mL in 50 mL ETOH <ul style="list-style-type: none"> • <u>STOP</u>
---------	--

944

Figure 2 Solubility Flow Chart

Tier	1		2		3		4		5		6
Concentration in Medium	Start Here 200 mg/mL	Incomplete solubility →	20 mg/mL	Incomplete solubility →	2 mg/mL		0.20 mg/mL				
Concentration in DMSO					200 mg/mL		20 mg/mL		2 mg/mL		0.2 mg/mL
Concentration in ETOH					200 mg/mL	Incomplete solubility →	20 mg/mL	Incomplete solubility →	2 mg/mL	Incomplete solubility →	0.2 mg/mL End
Concentration on Cells	100 mg/mL		10 mg/mL		1 mg/mL		0.1 mg/mL		0.01 mg/mL		0.001 mg/mL

Notes: 3T3 Medium - Dulbecco's Modification of Eagle's Medium, with supplements, for 3T3 mouse fibroblasts

[This Page Intentionally Left Blank]

42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65
66
67
68
69
70
71

Appendix B-2

Test Method Protocol for the Normal Human Epidermal Keratinocyte (NHK) Neutral Red Uptake (NRU) Cytotoxicity Test

This draft recommended NHK NRU test method protocol is substantially the same as the protocol used in Phase III of the NICEATM/ECVAM validation study. Revisions were made based on recommendations from NICEATM and the study directors involved in the study. The changes are as follows:

- Explanations and directions for the use of the revised Hill function for determining IC50 values are included in the protocol.
- The range for relative humidity values for the cell culture incubators was changed from 90 % ± 5 % humidity to 90 % ± 10 % humidity.
- An additional step was added to the test substance solubility protocol to allow testing of higher concentrations of test material.
- The spreadsheet templates used in the NICEATM/ECVAM validation study are incorporated into this protocol as an annex (ANNEX I).
- The stand-alone solubility protocol is incorporated into this protocol as an annex (ANNEX II).
- The stand-alone NHK media prequalification protocol is incorporated into this protocol as an annex (ANNEX III).

71

72

73

74

75

76

77

78

79

80

81

82

83

[This Page Intentionally Left Blank]

84

TEST METHOD PROTOCOL

THE NORMAL HUMAN KERATINOCYTE (NHK) NEUTRAL RED UPTAKE (NRU) CYTOTOXICITY TEST A Test for Basal Cytotoxicity

I. PURPOSE

This test method is used to evaluate the cytotoxicity of test substances using the Normal Human Keratinocyte (NHK) Neutral Red Uptake (NRU) *in vitro* cytotoxicity test. The data generated from the *in vitro* cytotoxicity assays are used to evaluate the effectiveness of the cytotoxicity assay to predict the starting doses for rodent acute oral systemic toxicity assays. This test method protocol outlines the procedures for performing the basal cytotoxicity test and is the result of the *in vitro* validation study organized by NICEATM and the European Centre for the Validation of Alternative Methods (ECVAM).

If changes or modifications are made to this protocol, the testing laboratory should prove that the results are comparable to those obtained when using the original protocol.

II. TEST SYSTEM

The NRU cytotoxicity assay procedure is a cell survival/viability chemosensitivity assay based on the ability of viable cells to incorporate and bind neutral red (NR), a supravital dye. NR is a weak cationic dye that readily penetrates cell membranes by non-ionic diffusion and accumulates intracellularly in lysosomes. Alterations of the cell surface or the sensitive lysosomal membrane lead to lysosomal fragility and other changes that gradually become irreversible. Such changes brought about by the action of xenobiotics result in a decreased uptake and binding of NR. It is thus possible to distinguish between viable, damaged, or dead cells, which is the basis of this assay.

Healthy mammalian cells, when maintained in culture, continuously divide and multiply over time. A toxic substance, regardless of site or mechanism of action, will interfere with this process and result in a reduction of the growth rate as reflected by cell number. Cytotoxicity is expressed as a concentration dependent reduction of the uptake of the NR after substance exposure thus providing a sensitive, integrated signal of both cell integrity and growth inhibition.

III. KEY PERSONNEL

A. Laboratory

- 1) Study Director
- 2) Laboratory Technician(s)

- 130 **B. Testing Facility**
 131 1) Scientific Advisor
 132 2) Quality Assurance Director
 133 3) Safety Manager
 134 4) Facility Management
 135

136 IV. DEFINITIONS

- 137
 138 **A. Hill function:** a four parameter logistic mathematical model relating the
 139 concentration of test substance to the response being measured in a sigmoidal
 140 shape.
 141

$$142 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log EC_{50} - X)\text{HillSlope}}}$$

143
 144 where Y= response, X is the logarithm of dose (or concentration), Bottom is
 145 the minimum response, Top is the maximum response, logEC50 is logarithm
 146 of X at the response midway between Top and Bottom, and HillSlope
 147 describes the steepness of the curve. When Top = 100 and Bottom = 0, the
 148 EC₅₀ is the concentration at 50% viability (i.e., the IC₅₀).
 149

- 150 **B. Documentation:** all methods and procedures will be noted in a Study Workbook; logs
 151 will be maintained for general laboratory procedures and equipment (e.g., media
 152 preparation, test substance preparation, incubator function); all optical density data
 153 obtained from the spectrophotometer plate reader will be saved in electronic and
 154 paper formats; all calculations of IC_x values and other derived data will be in
 155 electronic and paper format; all data will be archived.
 156
 157 **C. IC₅₀:** test substance concentration producing 50% inhibition of the endpoint measured
 158 (i.e., cell viability).
 159

160 V. IDENTIFICATION OF CONTROL SUBSTANCES

- 161
 162 **A. Positive Control (PC)**
 163 Sodium Lauryl Sulfate (SLS)
 164
 165 **B. Vehicle Control (VC)**
 166 Keratinocyte assay medium
 167
 168 **C. Solvent Control**
 169 VC control with solvent (i.e., keratinocyte assay medium, dimethyl sulfoxide
 170 [DMSO], or ethanol [ETOH])
 171 (DMSO is the preferred solvent for substances that are not water [i.e., assay medium]
 172 soluble.)

173

174 **VI. PROCEDURES**

175

176 **A. Materials**

177

1. Cell Line

178

Normal Human Epidermal Keratinocytes (NHK)

179

Non-transformed cells; from cryopreserved primary or secondary cells (e.g.,

180

Clonetics #CC-2507 or equivalent - Cambrex [Cambrex Bio Science, 8830

181

Biggs Ford Road, Walkersville, MD). Cells will be Clonetics NHK cells.

182

183

2. Technical Equipment

184

[Note: Suggested brand names/vendors are listed in parentheses. Equivalentents may
be used.]

185

186

8 Incubator: 37°C ± 1°C, 90 % ± 10 % humidity, 5.0 % ± 1 % CO₂/air

187

9 Laminar flow clean bench (standard: "biological hazard")

188

10 Waterbath: 37°C ± 1°C

189

11 Inverse phase contrast microscope

190

12 Sterile glass tubes with caps (e.g., 5 mL)

191

13 Centrifuge

192

14 Laboratory balance

193

15 96-well plate spectrophotometer (i.e., plate reader) equipped with 540 nm ± 10
nm filter

194

16 Shaker for microtiter plates

195

17 Cell counter or hemocytometer

196

18 Pipetting aid

197

19 Pipettes, pipettors (multi-channel and single channel; multichannel repeater
pipette), dilution block

198

20 Cryotubes

199

21 Tissue culture flasks (75 - 80 cm², 25 cm²)

200

22 96-well flat bottom tissue culture microtiter plates (e.g., Nunc # 167 008;

201

Corning/COSTAR tissue culture-treated)

202

23 pH paper (wide and narrow range)

203

q) Multichannel reagent reservoir

204

r) Waterbath sonicator

205

28 Magnetic stirrer

206

29 Antistatic bar ionizer/antistatic gun (optional for neutralizing static on 96-well
plates)

207

30 Dry heat block (optional)

208

31 Adhesive film plate sealers (e.g., Excel Scientific SealPlate™, Cat # STR-
SEAL-PLT or equivalent)

209

32 Vortex mixer

210

33 Filters/filtration devices

211

212

213

214

215

216

217 [Note: Tissue culture flasks and microtiter plates should be prescreened to ensure
218 that they adequately support the growth of NHK. Multi-channel repeater pipettes
219 may be used for plating cells in the 96-well plates, dispensing plate rinse
220 solutions, NR medium, and desorb solution. Do not use the repeater pipette for
221 dispensing test substances to the cells.]
222

223

3. Chemicals, Media, and Sera

- 224 9 Keratinocyte Basal Medium without Ca^{++} (e.g., KBM®, Clonetics CC-3104)
225 that is completed by adding supplements (e.g., KBM® SingleQuots®,
226 Clonetics CC-4131) to achieve the proper concentrations of epidermal growth
227 factor, insulin, hydrocortisone, antimicrobial agents, bovine pituitary extract, and
228 calcium (e.g., Clonetics Calcium SingleQuots®, 300 mM CaCl_2 , Clonetics CC-
229 4202).
230 10 HEPES Buffered Saline Solution (HEPES-BSS) (e.g., Clonetics # CC-5022)
231 11 0.025 % Trypsin/EDTA solution (e.g., Clonetics # CC-5012)
232 12 Trypsin Neutralizing Solution (TNS) (e.g., Clonetics # CC-5002)
233 13 Phosphate Buffered Saline (PBS)
234 14 Dulbecco's Phosphate Buffered Saline (D-PBS) [formulation containing calcium
235 and magnesium cations; glucose optional] (for rinsing)
236 15 Neutral Red (NR) Dye – tissue culture-grade; liquid form (e.g., SIGMA N
237 2889); powder form (e.g., SIGMA N 4638)
238 16 DMSO, U.S.P analytical grade (Store under nitrogen @ -20°C)
239 17 ETOH, U.S.P. analytical grade (100 %, non-denatured for test substance
240 preparation; 95 % can be used for the desorb solution)
241 18 Glacial acetic acid, analytical grade
242 19 Hanks' Balanced Salt Solution without Ca^{2+} or Mg^{2+} (CMF-HBSS) (e.g.,
243 Invitrogen # 14170)
244 20 Distilled H_2O or any purified water suitable for cell culture and NR desorb
245 solution (sterile)
246 21 Sterile/non-sterile paper towels (for blotting 96-well plates)
247

248

B. Preparations of Media and Solutions

249 [Note: All solutions (except NR stock solution, NR medium and NR desorb),
250 glassware, pipettes, etc., shall be sterile and all procedures should be carried out
251 under aseptic conditions and in the sterile environment of a laminar flow cabinet
252 (biological hazard standard). All methods and procedures will be adequately
253 documented.]
254

255

1. Media

256 [Note: This protocol is based on the use of Clonetics KBM® medium and
257 supplements. Other media may be acceptable if proper cell growth conditions can
258 be maintained as per this protocol. Prequalify candidate media by using the
259 keratinocyte medium prequalification in ANNEX III.]
260

- 261 a) Routine Culture Medium/Treatment Medium
 262
 263 KBM® (Clonetics CC-3104) supplemented with KBM® SingleQuots®
 264 (Clonetics CC-4131) and Clonetics Calcium SingleQuots® (CC-4202) to make
 265 500 mL medium. Final concentrations of supplements in medium are:
 266
 267 0.0001 ng/mL Human recombinant epidermal growth factor
 268 5 µg/mL Insulin
 269 0.5 µg/mL Hydrocortisone
 270 30 µg/mL Gentamicin
 271 15 ng/mL Amphotericin B
 272 0.10 mM Calcium
 273 30 µg/mL Bovine pituitary extract

274

275 *Complete media should be kept at 2-8°C and stored for no longer than two*
 276 *weeks.*

277

278 NOTE: KBM® SingleQuots® contain the following stock concentrations and
 279 volumes:

280

281	0.1 ng/mL	hEGF	0.5 mL
282	5.0 mg/mL	Insulin	0.5 mL
283	0.5 mg/mL	Hydrocortisone	0.5 mL
284	30 mg/mL	Gentamicin, 15 µg/mL Amphotericin-B	0.5 mL
285	7.5 mg/mL	Bovine Pituitary Extract (BPE)	2.0 mL

286

287 Clonetics Calcium SingleQuots® are 2 mL of 300 mM calcium.

288

289 165 µl of solution per 500 mL calcium-free medium equals 0.10 mM calcium
 290 in the medium.

291

292

2. Neutral Red (NR) Stock Solution

293 The liquid tissue culture-grade stock NR Solution will be the first choice for
 294 performing the assay (e.g., SIGMA #N2889, 3.3 mg/mL). Store liquid tissue
 295 culture-grade NR Stock Solution at the storage conditions and shelf-life period
 296 recommended by the manufacturer.

297

298 A stock solution can be made with powder NR dye and water (e.g., 0.33 g NR
 299 Dye powder in 100 mL H₂O) if the liquid stock form is not available. The stock
 300 should be stored in the dark at room temperature for up to two months.

301

302

3. Neutral Red (NR) Medium

303

EXAMPLE:

304

1.0 mL (3.3 mg NR dye/mL) NR Stock Solution

305 99 ml 99.0 mL Routine Culture Medium (pre-warmed to 37°C)

306

307 The final concentration of the NR Medium is **33 µg NR dye/mL** and aliquots will
308 be prepared on the day of application.

309

310 [Note: The NR Medium shall be filtered (e.g., Millipore filtering, 0.2 – 0.45 µm
311 pore size) used to reduce NR crystals. Aliquots of the NR Medium should be
312 maintained at 37°C (e.g., in a waterbath) before adding to the cells and used
313 within 60 minutes of preparation but also used within 15 minutes after removing
314 from 37°C storage. Examine the solution for crystals.]

315

316 4. ETOH/Acetic Acid Solution (NR Desorb)

317 1 % Glacial acetic acid solution

318 50 % ETOH

319 49 % H₂O

320

321 C. Methods

322

1. Cell Maintenance and Culture Procedures

323 NHK cells are routinely grown as a monolayer in tissue culture grade flasks (e.g.,
324 25 cm²) at 37°C ± 1°C, 90 % ± 10 % humidity, and 5.0 % ± 1 % CO₂/air. The
325 cells should be examined on a daily (i.e., on workdays) basis under a phase
326 contrast microscope, and any changes in morphology or their adhesive properties
327 must be noted in a Study Workbook.

328

329 2. Receipt of Cryopreserved Keratinocytes

330 Upon receipt of cryopreserved keratinocytes, the vial(s) of cells shall be stored in
331 a liquid nitrogen freezer until needed.

332

333 3. Thawing Cells and Establishing Cell Cultures

334 4 Thaw cells by putting ampules into a water bath at 37°C for as brief a time as
335 possible. Do not thaw cells at room temperature or by hand. Seed the thawed
336 cells into culture flasks as quickly as possible and with minimal handling.

337 5 Slowly (taking approximately 1-2 minutes) add 9 mL of pre-warmed Routine
338 Culture Medium to the cells suspended in the cryoprotective solution and
339 transfer cells into flasks containing pre-warmed Routine Culture Medium (See
340 **Table 1**).

341 6 Incubate the cultures at 37°C ± 1°C, 90 % ± 10 % humidity, 5.0 % ± 1 %
342 CO₂/air until the cells attach to the flask (within 4 to 24 hours), at which time the
343 Routine Culture Medium should be removed and replaced with fresh Routine
344 Culture Medium.

345 7 Unless otherwise specified, the cells should be incubated at $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$, $90\% \pm$
346 10% humidity, $5.0\% \pm 1\%$ CO_2 /air and fed every 2-3 days until they exceed
347 50% confluence (but less than 80% confluent).
348
349

349
350**Table 1. Guidelines for Establishing Cell Cultures**

Cells/25 cm ² flask (in approximately 5 mL) 1 flask each cell concentration	6.25 x 10 ⁴ (2500/cm ²)	1.25 x 10 ⁵ (5000/cm ²)	2.25 x 10 ⁵ (9000/cm ²)
Approximate Time to Subculture	96+ hours	72 - 96 hours	48 - 72 hours
Cells to 96-Well Plates	6 – 8 plates	6 – 8 plates	6 – 8 plates

351
352
353
354

Cell growth guidelines – actual growth of individual cell lots may vary.

4. Subculture of NHK Cells to 96-Well Plates

355
356
357
358
359
360
361
362
363
364
365
366
367
368
369
370
371
372
373
374
375
376
377
378
379
380
381
382
383
384
385
386

[Note: It is important that cells have overcome the lag growth phase when they are used for the test. Keratinocytes will be passaged only into the 96-well plates and will not be subcultured into flasks for use in later assays]

- 8 When the keratinocyte culture in a 25 cm² flask exceeds 50 % confluence (but less than 80 % confluent), remove the medium and rinse the culture twice with 5 mL HEPES-BSS. The first rinse may be left on the cells for up to 5 minutes and the second rinse should remain on the cells for approximately 5 minutes. Discard the washing solutions.
- 9 Add 2 mL trypsin/EDTA solution to each flask and remove after 15 to 30 seconds. Incubate the flask at room temperature for 3 to 7 minutes. When more than 50 % of the cells become dislodged, rap the flask sharply against the palm of the hand.
- 10 When most of the cells have become detached from the surface, rinse the flask with 5 mL of room temperature TNS. If more than one flask is subcultured, the same 5 mL of TNS may be used to rinse a total of up to two flasks.
- 11 Then rinse the flask with 5 mL CMF-HBSS and transfer the cell suspension to a centrifuge tube.
- 12 Pellet the cells by centrifugation for 5 minutes at approximately 220 x g. Remove the supernatant by aspiration.
- 13 Resuspend the keratinocyte pellet by gentle trituration (to have single cells) in Routine Culture Medium. It is important to obtain a single cell suspension for exact counting. Count a sample of the cell suspension using a hemocytometer or cell counter.
- 14 Prepare a cell suspension $1.6 - 2.0 \times 10^4$ cells/mL in Routine Culture Medium. Using a multi-channel pipette, dispense 125 μ l Routine Culture Medium only into the peripheral wells (blanks) of a 96-well tissue culture microtiter plate. In the remaining wells, dispense 125 μ l of the cell suspension ($2 \times 10^3 - 2.5 \times 10^3$ cells/well). Prepare one plate per substance to be tested (see **Figure 1, Section VI.E.1**).
- 15 Incubate cells ($37^\circ\text{C} \pm 1^\circ\text{C}$, $90\% \pm 10\%$ humidity, and $5\% \pm 1\%$ CO₂/air) so that cells form a 20+ % monolayer (~48-72 hours). This incubation period

387 assures cell recovery and adherence and progression to exponential growth
388 phase.
389 16 Examine each plate under a phase contrast microscope to assure that cell
390 growth is relatively even across the microtiter plate. This check is performed
391 to identify experimental and systemic cell seeding errors. Record
392 observations in the Study Workbook.
393
394

5. Determination of Doubling Time

- 395 a) Establish cells in culture and trypsinize cells as per **Section VI.C.4** for
396 subculture. Resuspend cells in appropriate culture medium. Use **Table 1** to
397 determine seeding densities.
398 b) Seed five sets of cell culture vessels in triplicate for each cell type (e.g., 15
399 tissue culture dishes [60mm x 15mm]). Use appropriate volume of culture
400 medium for the culture vessels. Note number of cells placed into each culture
401 dish. Place dishes into the incubators ($37^{\circ}\text{C} \pm 1^{\circ}\text{C}$, $90\% \pm 10\%$ humidity, 5.0
402 $\% \pm 1\%$ CO_2/air).
403 c) After 4-6 hours (use the same initial measurement time for each subsequent
404 doubling time experiment), remove three culture dishes and trypsinize cells.
405 Count cells using a cell counter or hemocytometer. Cell viability may be
406 determined by dye exclusion (e.g., Trypan Blue; Nigrosin). Determine the
407 total number of cells and document. Repeat sampling at 24 hours, 48 hours,
408 72 hours, and 96 hours post inoculation. Change culture medium at 72 hours
409 or sooner in remaining dishes if indicated by pH drop.
410 d) Plot cell concentration (per mL of medium) on a log scale against time on a
411 linear scale. Determine lag time and population doubling time. The doubling
412 time will be in the log (exponential) phase of the growth curve. Additional
413 dishes and time are needed if the entire growth curve is to be determined (lag
414 phase, log phase, plateau phase).
415

D. Preparation of Test Substances

416 [Note: Preparation under red or yellow light is recommended to preserve substances
417 that degrade upon exposure to light.]
418

419
420 Test substance solubility should be determined by following the procedures outlined
421 in **ANNEX II** of this protocol.
422

1. Test Substance in Solution

- 423
424 • Allow test substances to equilibrate to room temperature before dissolving and
425 diluting.
426 • Prepare test substance immediately prior to use and not in bulk for use in
427 subsequent tests. Ideally, the solutions must not be cloudy nor have
428 noticeable precipitate. Each stock dilution should have at least 1-2 mL total
429 volume to ensure adequate solution for the test wells in a single 96-well plate.
430 The Study Director may store an aliquot (e.g., 1 mL) of the highest 2X stock
431 solution (e.g., low solubility substances) in a freezer (e.g., -70°C) for use in
432 future substance analyses.

- 433
- 434
- 435
- 436
- 437
- 438
- 439
- 440
- 441
- 442
- 443
- 444
- 445
- 446
- 447
- 448
- 449
- 450
- For substances dissolved in DMSO or ETOH, the final DMSO or ETOH concentration for application to the cells must be 0.5 % (v/v) in the vehicle controls and in all of the eight test concentrations.
 - The stock solution for each test substance should be prepared at the highest concentration found to be soluble in the solubility test (*ANNEX II: Test Method Procedure - Solubility Determination of Test Substances*). Thus, the highest test concentration applied to the cells in each range finding experiment is:
 - 0.5 times the highest concentration found to be soluble in the solubility test, if the substance was soluble in medium, or
 - 1/200 the highest concentration found to be soluble in the solubility test if the substance was soluble in ETOH or DMSO.
 - The seven lower concentrations in the range finding experiment are prepared by successive dilutions that decrease by one log unit each. The following example illustrates the preparation of test substance in solvent and the dilution of dissolved test substance in medium before application to NHK cells.

451 *Example: Preparation of Test Substance in Solvent Using a Log Dilution Scheme*
452 If DMSO is determined to be the preferred solvent at Tier 3 of the solubility test
453 (i.e., 200,000 µg/mL), dissolve the substance in DMSO at 200,000 µg/mL for the
454 substance stock solution.

- 455
- 456 8) Label eight tubes 1 – 8. Add 0.9 mL solvent (e.g., DMSO) to tubes 2 -- 8.
- 457 9) Prepare stock solution of 200,000 µg test substance/mL solvent in tube # 1.
- 458 10) Add 0.1 mL of 200,000 µg/mL dilution from tube #1 to tube #2 to make a
459 1:10 dilution in solvent (i.e., 20,000 µg/mL).
- 460 11) Add 0.1 mL of 20,000 µg/mL dilution from tube #2 to tube #3 to make
461 another 1:10 dilution (i.e., 1:100 dilution from stock solution) in solvent (i.e.,
462 2,000 µg/mL)
- 463 12) Continuing making serial 1:10 dilutions in the prepared solvent tubes.
- 464 13) Since each concentration is 200 fold greater than the concentration to be
465 tested, make a 1:100 dilution by diluting 1 part dissolved substance in each
466 tube with 99 parts of culture medium (e.g., 0.1 mL of test substance in DMSO
467 + 9.9 mL culture medium) to derive the eight 2X concentrations for
468 application to NHK cells. Each 2X test substance concentration will then
469 contain 1 % v/v solvent. The NHK cells will have 0.125 mL of culture
470 medium in the wells prior to application of the test substance. By adding
471 0.125 mL of the appropriate 2X test substance concentration to the appropriate
472 wells, the test substance will be diluted appropriately (e.g., highest
473 concentration in well will be 1,000 µg/mL) in a total of 0.250 mL and the
474 solvent concentration in the wells will be 0.5% v/v.
- 475 14) A test substance prepared in DMSO or ETOH may precipitate upon transfer
476 into the Routine Culture Medium. The 2X dosing solutions should be
477 evaluated for precipitates and the results recorded in the Study Workbook. It

478 will be permissible to test all of the dosing solutions in the dose range finding
479 assay and main experiments. However, doses containing test substance
480 precipitates should be avoided because it creates doubt about the
481 concentration of test substance exposed to the cells.

482

483

Document all test substance preparations in the Study Workbook.

484

485

2. pH of Test Substance Solutions

486

487

488

489

490

491

492

493

494

495

Prior to or immediately after application of the test substance to the 96-well plate, measure the pH of the highest 2X dosing concentration of the test substance (i.e., C1 in the test plate, see Figure 1) in culture medium. Use pH paper (e.g., pH 0 – 14 to estimate and pH 5 – 10 to determine more precise value; or Study Director's discretion). The pH paper should be in contact with the solution for approximately one minute. Document the pH and note the color of the 2X concentration medium (i.e., in the Microsoft Excel® template; see ANNEX 1 for an example template). Medium color for all dosing dilutions should be noted in the Study Workbook. Do not adjust the pH.

496

3. Concentrations of Test Substance

497

a) Range Finder Experiment

498

499

500

501

502

503

504

505

506

507

508

509

510

511

512

513


514

515

- Test eight concentrations of the test substance by diluting the stock solution using log dilutions (e.g., 1:10, 1:100, 1:1000, etc.).
- If a range finder experiment does not generate enough cytotoxicity, then higher doses should be attempted. If cytotoxicity is limited by solubility, then more stringent solubility procedures to increase the stock concentration (to the maximum concentration specified in **Section VI.D.3.b.**) should be employed.
- Place the highest test substance concentration into an incubator (37°C ± 1°C, 90 % ± 10 % humidity, 5.0 % ± 1 % CO₂/air) and stir or rock for up to 3 hours, if necessary, to facilitate dissolution. For stocks prepared in medium, vessel caps should be loose to allow for CO₂ exchange. Proceed with dosing solution preparation and dosing.
- If a range finding test produces a biphasic curve, then the doses selected for the subsequent main experiments should cover the most toxic dose-response range (see Example 1 – the most toxic range is 0.001 – 0.1 µg/mL) that reduces viability to 50%.

515
516

Example 1 – Biphasic Curve



517
518
519
520
521
522
523
524
525
526
527
528
529
530
531
532

b) Main Experiment (Definitive Assay)

- Depending on the slope of the concentration-response curve estimated from the range finder, the dilution/progression factor in the concentration series of the main experiment should be smaller (e.g., dilution factor of $\sqrt[6]{10} = 1.47$). Cover the relevant concentration range around the IC_{50} (> 0 % and < 100 % effect) preferably with several points of a graded effect, but with a minimum of two points, one on each side of the estimated IC_{50} value, avoiding too many non-cytotoxic and/or 100 %-cytotoxic concentrations.
- Determine which test substance concentration is closest to the IC_{50} value. Use that value as a central concentration and adjust dilutions higher and lower in equal steps for the definitive assay.

Maximum Doses to be Tested in the Main Experiments

If minimal or no cytotoxicity was measured in the dose range finding assay, a maximum dose for the main experiments will be established as follows:

533
534
535
536
537
538
539
540
541
542
543
544
545
546
547
548

For test substances prepared in Routine Culture Medium

- The highest test substance concentration that may be applied to the cells in the main experiments will be either 100 mg/mL, or the maximum soluble dose.
- Test substance will be weighed into a glass tube and the weight will be documented. A volume of Routine Culture Medium will be added to the vessel so that the concentration is 200,000 µg/mL (200 mg/mL).
- The solution is mixed using the mechanical procedures specified in ANNEX II of this protocol. If complete solubility is achieved in medium, then seven additional serial stock dosing solutions may be prepared from the 200 mg/mL 2X stock. If the test substance is insoluble in medium at 200 mg/mL, proceed by adding medium, in small incremental amounts, to

549 attempt to dissolve the substance by using the sequence of mixing
 550 procedures specified in ANNEX II.
 551 • More stringent solubility procedures may be employed if needed based on
 552 results from the range finder experiment (**Section VI.D.3.a.**). The highest
 553 soluble stock solution will be used to prepare the seven additional serial
 554 stock dosing solutions.

555
 556 For test substances prepared in either DMSO or ETOH

557 11.0 The highest test substance concentration that may be applied to the
 558 cells in the main experiments will be ≤ 2.5 mg/mL, depending upon the
 559 maximum solubility in solvent.

560 12.0 Test substance will be weighed into a glass tube and the weight will be
 561 documented. A volume of the appropriate solvent (determined from the
 562 original solubility test) will be added to the vessel so that the concentration
 563 is 500,000 $\mu\text{g/mL}$ (500 mg/mL).

564 13.0 The solution is mixed as specified ANNEX II. If complete solubility
 565 is achieved in the solvent, then seven additional serial stock dosing
 566 solutions may be prepared from the 500 mg/mL 200X stock. If the test
 567 substance is insoluble in solvent at 500 mg/ml, proceed by adding solvent,
 568 in small incremental amounts, to attempt to dissolve the substance by
 569 using the sequence of mixing procedures. The highest soluble stock
 570 solution will be used to prepare the seven additional serial stock dosing
 571 solutions.

572
 573 If precipitates are observed in the 2X dilutions, continue with the experiment,
 574 make the appropriate observations and documentation, and report data to the
 575 SMT.

576
 577 c) Test Substance Dilutions

578 The dosing factor of 3.16 ($=\sqrt[2]{10}$) divides a log into two equidistant steps,
 579 2.15 ($=\sqrt[3]{10}$) into three steps, 1.47 ($=\sqrt[6]{10}$) into six steps, 1.78 ($=\sqrt[4]{10}$) into
 580 four steps, and 1.21 ($=\sqrt[12]{10}$) into 12 steps.

581
 582 EXAMPLE:

10						31.6						100
10				21.5				46.4				100
10		14.7		21.5		31.6		46.4		68.1		100
10	12.1	14.7	17.8	21.5	26.1	31.6	38.3	46.4	56.2	68.1	82.5	100

584
 585 An example of decimal geometric concentration series for factor 1.47: *Dilute*
 586 *1 volume of the highest concentration by adding 0.47 volumes of diluent.*
 587 *After equilibration, dilute 1 volume of this solution by adding 0.47 volumes of*
 588 *diluent...(etc.).*
 589
 590

590
591
592
593
594
595
596
597

E. Test Procedure

1. 96-Well Plate Configuration

The NHK NRU assay for test substances will use the 96-well plate configuration shown in **Figure 1**.

Figure 1. 96-Well Plate Configuration for Positive Control (PC) and Test Substance Assays

	1	2	3	4	5	6	7	8	9	10	11	12
A	VCb	VCb	C ₁ b	C ₂ b	C ₃ b	C ₄ b	C ₅ b	C ₆ b	C ₇ b	C ₈ b	VCb	VCb
B	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
C	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
D	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
E	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
F	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
G	VCb	VC1	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	VC2	VCb
H	VCb	VCb	C ₁ b	C ₂ b	C ₃ b	C ₄ b	C ₅ b	C ₆ b	C ₇ b	C ₈ b	VCb	VCb

598
599
600
601
602
603
604

VC1 and VC2 = VEHICLE CONTROL
 C₁ – C₈ = Test Substances or PC (SLS) at eight concentrations
 (C₁ = highest, C₈ = lowest)
 b = BLANKS (Test substance or PC, but contain **no** cells)
 VCb = VEHICLE CONTROL BLANK (contain **no** cells)

605
606
607
608
609
610
611
612
613
614
615
616
617

2. Application of Test Substance

- Two optional methods for rapidly applying the 2X dosing solutions onto the 96-well plates may be utilized.
 - Add each of the 2X dosing solutions into labeled, sterile reservoirs (e.g., Corning/Costar model 4870 sterile polystyrene 50 mL reagent reservoirs or Corning/Transtar model 4878 disposable reservoir liners, 8-channel; or other multichannel reservoirs).
 - Use a “dummy” plate (i.e., an empty sterile 96-well plate) prepared to hold the dosing solutions immediately prior to treatment of the test plate (with cells). The test substance and control dosing solutions should be dispensed into the dummy plate in the same pattern/order as will be applied to the plate containing cells. More volume than needed for the test

618 plate (i.e., greater than 125 μl /well) should be in the wells of the dummy
619 plate.

620

621 At the time of treatment initiation, a multi-channel micropipettor is used to
622 transfer the 2X dosing solutions, from the reservoirs or dummy plate, to the
623 appropriate wells on the treatment plate (as described in step c. below). These
624 methods will ensure that the dosing solutions can be transferred rapidly to the
625 appropriate wells of the test plate to initiate treatment times and to minimize
626 the range of treatment initiation times across a large number of treatment
627 plates, and to prevent “out of order” dosing. Do not use a multichannel
628 repeater pipette for dispensing test substance to the plates.

629

630 2. After 48 - 72 hours (i.e., after cells attain 20+ % confluency [see **Section**
631 **VI.C.4(h)**]) incubation of the cells, add 125 μl of the appropriate
632 concentration of test substance, the PC, or the VC (see **Figure 1** for the
633 plate configuration) directly to the test wells. Do not remove Routine
634 Culture Medium for re-feeding the cells. The dosing solutions will be
635 rapidly transferred from the 8-channel reservoir (or dummy plate) to the
636 test plate using a single delivery multi-channel pipettor. For example, the
637 VC may be transferred first (into columns 1, 2, 11, and 12), followed by
638 the test substance dosing solutions from lowest to highest dose, so that the
639 same pipette tips on the multi-channel pipettor can be used for the whole
640 plate. [The Vehicle Control blank (VCb) wells (column 1, column 12,
641 wells A2, A11, H2, H11) will receive the Vehicle Control dosing solutions
642 (which should include any solvents used)]. Blanks for wells A3 – A10
643 and H3 – H10 shall receive the appropriate test substance solution for each
644 concentration (e.g., wells A3 and H3 receive C₁ solution.) Incubate cells
645 for 48 hours \pm 0.5 hours (37°C \pm 1°C, 90 % \pm 10 % humidity, and 5.0 % \pm
646 1 % CO₂/air).

647 3. **Positive Control:** For each set of test substance plates used in an assay,
648 prepare a separate plate of positive control concentrations. If multiple sets
649 of test substance plates are set up, then clearly designate the positive
650 control plates for each set; each set will be an individual entity. The Study
651 Director will decide how many test substance plates will be run with a
652 positive control plate. This plate will follow the same schedule and
653 procedures as used for the test substance plates (including appropriate
654 substance concentrations in the appropriate wells and meeting test
655 acceptance criteria see **Sections VI.E.1, E.2, and E.5**).

656

657

3. Microscopic Evaluation

658 After at least 46 hours treatment, examine each plate under a phase contrast
659 microscope to identify systematic cell seeding errors and growth characteristics of
660 control and treated cells. Record any changes in morphology of the cells due to
661 the cytotoxic effects of the test substance, but do not use these records for any
662 quantitative measure of cytotoxicity. Undesirable growth characteristics of
663 control cells may indicate experimental error and may be cause for rejection of

664 the assay. Use the following Visual Observations Codes in the description of cell
 665 culture conditions. Numerical scoring of the cells (see **Section VI.E.3**) should be
 666 determined and documented in the Study Workbook and in the appropriate
 667 section of the Microsoft Excel® template.

668
 669
 670

Visual Observations Codes

Note Code	Note Text
1	Normal Cell Morphology
2	Low Level of Cell Toxicity
3	Moderate Level of Cell Toxicity
4	High level of Cell Toxicity
1P	Normal Cell Morphology with Precipitate
2P	Low Level of Cell Toxicity with Precipitate
3P	Moderate Level of Cell Toxicity with Precipitate
4P	High level of Cell Toxicity with Precipitate
5P	Unable to View Cells Due to Precipitate

671
 672

4. Measurement of NRU

- 673 a) Carefully remove (i.e., “dump”) the Routine Culture Medium (with test
 674 substance) and rinse the cells very carefully with 250 µL pre-warmed D-PBS.
 675 Remove the rinsing solution by dumping and remove excess by gently
 676 blotting on paper towels. Add 250 µL NR medium (to all wells including the
 677 blanks) and incubate (37°C ± 1°C, 90 % ± 10 % humidity, and 5.0 % ± 1 %
 678 CO₂/air) for 3 hours ± 0.1 hour. Observe the cells briefly during the NR
 679 incubation (e.g., between 2 and 3 hours – Study Director’s discretion) for NR
 680 crystal formation. Record observations in the Study Workbook. Study
 681 Director can decide to reject the experiment if excessive NR crystallization
 682 has occurred.
- 683 b) After incubation, remove the NR medium, and carefully rinse cells with 250
 684 µL pre-warmed D-PBS.
- 685 c) Decant and blot D-PBS from the plate. (Optionally: centrifuge the reversed
 686 plate.)
- 687 d) Add exactly 100 µL NR Desorb (ETOH/acetic acid) solution to all wells,
 688 including blanks.
- 689 e) Shake microtiter plate rapidly on a microtiter plate shaker for 20 – 45 minutes
 690 to extract NR from the cells and form a homogeneous solution. Plates should
 691 be protected from light by using a cover during shaking.
- 692 f) Plates should be still for at least five minutes after removal from the plate
 693 shaker (or orbital mixer). If any bubbles are observed, assure that they have
 694 been ruptured prior to reading the plate. Measure the absorption (within 60
 695 minutes of adding NR Desorb solution) of the resulting colored solution at
 696 540 nm ± 10 nm in a microtiter plate reader (spectrophotometer), using the
 697 blanks as a reference.

698 [Note: A mean $OD_{540 \pm 10nm}$ of 0.043 - 0.059 for the VC blanks is a target
699 range of ODs but not a test acceptance criterion (range = mean OD \pm 2.5
700 standard deviations; mean = 0.054; SD = 0.003; N = 114)]. Save raw data in
701 the Microsoft Excel® template.

702
703

5. Quality Check of Assay

704
705
706
707

a) Test Acceptance Criteria

All acceptance criteria (i.e., criteria 1, 2, and 3) must be met for a test to be acceptable.

708
709
710
711
712
713
714
715
716
717
718
719

1) The PC (SLS) IC_{50} must be within two and a half (2.5) standard deviations (SD) of the historical mean established by the Test Facility and must meet criteria 2 and 3, and must have an r^2 (coefficient of determination) value calculated for the Hill model fit (i.e., from PRISM® software) \geq 0.85.

NICEATM/ECVAM study generated the following PC data:

- IC_{50} mean = 3.11 μ g/mL; SD = 0.72 (n = 114)
- range for IC_{50} mean \pm 2.5 SD = 1.31 μ g/mL – 4.91 μ g/mL

2) The left and the right mean of the VCs do not differ by more than 15 % from the mean of all VCs.

3) At least one calculated cytotoxicity value $>$ 0 % and \leq 50 % viability and at least one calculated cytotoxicity value $>$ 50 % and $<$ 100 % viability must be present.

720
721
722
723

Exception: If a test has only one point between 0 and 100 % **and** the smallest dilution factor (i.e., 1.21) was used **and** all other test acceptance criteria were met, then the test will be considered acceptable.

724
725
726
727
728
729

Stopping Rule for Insoluble Substances: If the most rigorous solubility procedures have been performed and the assay cannot achieve adequate toxicity to meet the test acceptance criteria after three definitive trials, then the Study Director may end all testing for that particular substance.

730
731
732
733

[Note: A corrected mean $OD_{540 \pm 10nm}$ of 0.205 - 1.645 for the VCs is a target range of ODs but will not be a test acceptance criterion (range = mean OD \pm 2.5 standard deviations; mean = 0.685; SD = 0.175; N = 114).]

734
735
736
737
738
739
740
741
742

b) Checks for Systematic Cell Seeding Errors

To check for systematic cell seeding errors, untreated VCs are placed both at the left side (row 2) and the right side (row 11 for the test plates) of the 96-well plate. Aberrations in the cell monolayer for the VCs may reflect a volatile and toxic test substance present in the assay. If volatility is suspected, then proceed to **Section VI.E.6**. Checks for cell seeding errors may also be performed by examining each plate under a phase contrast microscope to assure that cell quantity is consistent.

743

743 **7. Testing Volatile Substances**

744 Although this test method is not suitable for highly volatile substances, mildly
745 volatile substances may be tested with some success. Volatile test substances
746 may generate vapors from the treatment media during the test substance treatment
747 incubation period. These vapors may become resorbed into the treatment medium
748 in adjacent wells, such that culture wells nearest the highest doses may become
749 contaminated by exposure. If the test substance is particularly toxic at the doses
750 tested, the cross contamination may be evident as a significant reduction in
751 viability in the VC cultures (i.e., VC1) adjacent to the highest test substance
752 doses.

753
754 If potential test substance volatility is suspected (e.g., for low density liquids) or if
755 the initial range finder test (non-sealed plate) results show evidence of toxic
756 effects in the control cultures (i.e., > 15 % difference in viability between VC1
757 [column 2] and VC2 [column 11]), then seal the subsequent test plates by the
758 following procedure.

759
760 a) Plate Sealer Method

761 d) Plates and substances will be prepared as usual according to **Sections**
762 **VI.D and VI.E.**

763 e) Immediately after the 96-well culture plate has been treated with the
764 suspected volatile substance (**Section VI.E.2.b**), apply the adhesive plate
765 sealer (e.g., using a hand, microplate roller, etc.) directly over the culture
766 wells. Assure that the sealer adheres to each culture well (well tops should
767 be dry). Place the 96-well plate cover over the sealed plate and incubate
768 the plate under specified conditions (**Section VI.E.2.b**). [Note: Do not
769 jam the plate lid over the film to avoid deforming the sealer and causing
770 the sealer to detach from culture wells. Loose fit of the plate lid is
771 acceptable.]

772 f) At the end of the treatment period, the plate sealer should be carefully
773 removed to avoid spillage. Continue with the NRU assay as per **Section**
774 **VI.E.4.**

775

776 **F. Data Analysis**

- 777 • The Study Director will use good biological/scientific judgment for determining
778 “unusable” wells that will be excluded from the data analysis and provide
779 explanations for the removal of any data from the analysis.
- 780 • A calculation of cell viability expressed as NRU is made for each concentration of
781 the test substance by using the mean NRU of the six replicate values (minimum of
782 four acceptable replicates wells) per test concentration. This value is compared
783 with the mean NRU of all VC values. Relative cell viability is then expressed as
784 percent of untreated VC. If achievable, the eight concentrations of each substance
785 tested will span the range of no effect up to total inhibition of cell viability.
- 786 • Data from the microtiter plate reader should be transferred to a spreadsheet
787 template (e.g., Microsoft Excel®) that will automatically determine cell viability,
788 calculate IC₅₀ values by linear interpolation, and perform statistical analyses

789 (including statistical identification of outliers). The template should also calculate
 790 the concentrations associated with 20 %, 50 %, and 80 % viability using the Hill
 791 slope and EC₅₀ (i.e., IC₅₀) from the Hill function analysis (see ANNEX 1 for an
 792 example spreadsheet template).
 793 • A Hill function analysis should be performed using statistical software (e.g.,
 794 GraphPad PRISM® 3.0) and a template to calculate IC₂₀, IC₅₀, and IC₈₀ values
 795 (and the associated confidence limits) for each test substance.
 796 • Dose-responses for which the toxicity plateaus as concentration increases do not
 797 fit the Hill function well when Bottom = 0. To obtain a better model fit,
 798 unconstrain the Bottom parameter so that the model calculates the Bottom value.
 799 However, when Bottom ≠ 0, the EC₅₀ reported by the Hill function ≠ 50%
 800 viability since the Hill function defines EC₅₀ as the point midway between Top
 801 and Bottom. To obtain the appropriate IC₅₀ when Bottom ≠ 0, use the following
 802 rearranged Hill function:
 803

$$X = \log EC_{50} - \frac{\log\left(\frac{Top - Bottom}{Y - Bottom} - 1\right)}{HillSlope}$$

804
 805 • X is the logarithm of concentration at 50% response, logEC₅₀ is logarithm of
 806 concentration at the response midway between Top and Bottom, Top is the
 807 maximum response, Bottom is the minimum response, Y = 50 (i.e., 50%
 808 response), and HillSlope describes the steepness of the curve.
 809

810 [Note: IC₅₀ values are used in a regression formula to predict the LD₅₀ value of a
 811 test substance as an estimate of the starting dose for an acute oral toxicity test.]
 812
 813

814 VII. REFERENCES

815
 816 Clonetics Normal Human Keratinocyte Systems Instructions for Use, AA-1000-4-
 817 Rev.03/00. (<http://www.clonetics.com>).
 818
 819 Hackenberg, U. and H. Bartling. 1959. Messen und Rechnen im pharmakologischen
 820 Laboratorium mit einem speziellen Zahlensystem (WL24-System). Arch. Exp. Pathol.
 821 Pharmakol. 235: 437-463.
 822
 823 Triglia, D., P.T. Wegener, J. Harbell, K. Wallace, D. Matheson, and C. Shopsis. 1989.
 824 Interlaboratory validation study of the keratinocyte neutral red bioassay from
 825 Clonetics Corporation. In *Alternative Methods in Toxicology*, Volume 7. A.M.
 826 Goldberg, ed., pp. 357-365. Mary Ann Liebert, Inc., New York.
 827
 828 *Test Method Protocol for Solubility Determination. In Vitro Cytotoxicity Validation*
 829 *Study. Phase III. August 29, 2003. Prepared by The National Toxicology Program*
 830 *(NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods*
 831 *(NICEATM).* (<http://iccvam.niehs.nih.gov/methods/invitro.htm>)

832 U. S. Environmental Protection Agency. 1996. Product Properties Test Guidelines.
833 OPPTS 803.7840. Water Solubility: Column Elution Method; Shake Flask Method.
834 EPA712-C-96-041, Prevention, Pesticides and Toxic Substances, Washington DC.

835

836 Test Method Protocol for the NHK Neutral Red Uptake Cytotoxicity Test. A
837 Test for Basal Cytotoxicity for an In Vitro Validation Study. November 4,
838 2004. Prepared by the National Toxicology Program (NTP) Interagency
839 Center for the Evaluation of Alternative Toxicological Methods (NICEATM).
840 (<http://iccvam.niehs.nih.gov/methods/invitro.htm>)

841

842 Personal communication, Technical Services, Cambrex Bio Science,
843 Walkersville, MD

844

845

846

847

848

849

849
850

ANNEX 1
Microsoft EXCEL® Example Spreadsheet Template

Test Facility : A						Study Number.: A1						
Chemical Code : SLS						96-Well Plate ID : A11						
2nd Chem. Code*: 11						Experiment ID : XX						
96-WELL PLATE MAP												
	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank
B	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
C	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
D	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
E	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
F	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
G	Blank	VC1	C1	C2	C3	C4	C5	C6	C7	C8	VC2	Blank
H	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank
RAW ABSORBANCE DATA (OD₅₅₀)												
	1	2	3	4	5	6	7	8	9	10	11	12
A	0.044	0.044	0.045	0.045	0.045	0.046	0.051	0.057	0.057	0.043	0.041	0.044
B	0.042	0.456	0.043	0.043	0.130	0.300	0.395	0.414	0.418	0.402	0.401	0.042
C	0.043	0.407	0.042	0.041	0.130	0.294	0.383	0.382	0.413	0.375	0.385	0.044
D	0.043	0.438	0.042	0.043	0.147	0.337	0.409	0.404	0.438	0.436	0.391	0.047
E	0.044	0.448	0.041	0.045	0.132	0.321	0.429	0.414	0.416	0.420	0.441	0.042
F	0.045	0.411	0.040	0.042	0.127	0.375	0.397	0.402	0.422	0.447	0.403	0.043
G	0.041	0.405	0.043	0.040	0.124	0.361	0.444	0.442	0.425	0.448	0.405	0.044
H	0.041	0.041	0.048	0.042	0.042	0.044	0.042	0.042	0.040	0.044	0.041	0.041
Max	0.045	0.456	0.043	0.045	0.147	0.375	0.444	0.442	0.438	0.448	0.441	0.047
Min	0.041	0.405	0.040	0.040	0.124	0.294	0.383	0.382	0.413	0.375	0.385	0.041
Next Max	0.044	0.448	0.042	0.043	0.132	0.361	0.429	0.414	0.425	0.447	0.405	0.044
Next Min	0.042	0.407	0.041	0.041	0.127	0.300	0.395	0.402	0.416	0.402	0.391	0.042
Rmax	-0.250	-0.157	-0.333	-0.400	-0.652	-0.173	-0.246	-0.467	-0.520	-0.014	-0.643	-0.500
Rmin	0.250	0.039	0.333	0.200	0.130	0.074	0.197	0.333	0.120	0.370	0.107	0.167
CORRECTED ABSORBANCE (Sample OD₅₅₀ - Mean Blank OD₅₅₀)												
	1	2	3	4	5	6	7	8	9	10	11	12
A	0.001	0.001	-0.002	0.002	0.002	0.001	0.005	0.008	0.009	-0.001	-0.002	0.001
B	-0.001	0.413	-0.004	-0.001	0.087	0.255	0.349	0.365	0.370	0.359	0.358	-0.001
C	0.000	0.364	-0.005	-0.003	0.087	0.249	0.337	0.333	0.365	0.332	0.342	0.001
D	0.000	0.395	-0.005	-0.001	0.104	0.292	0.363	0.355	0.390	0.393	0.348	0.004
E	0.001	0.405	-0.006	0.002	0.089	0.276	0.383	0.365	0.368	0.377	0.398	-0.001
F	0.002	0.368	-0.007	-0.001	0.084	0.330	0.351	0.353	0.374	0.404	0.360	0.000
G	-0.002	0.362	-0.004	-0.004	0.081	0.316	0.398	0.393	0.377	0.405	0.362	0.001
H	-0.002	-0.002	0.002	-0.001	-0.001	-0.001	-0.005	-0.008	-0.009	0.001	-0.002	-0.002
Mean Blank =	0.043		0.047	0.044	0.044	0.045	0.047	0.050	0.049	0.044		
RELATIVE VIABILITY (% OF VEHICLE CONTROL)												
	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		110.7%	-0.9%	-0.1%	23.2%	68.4%	93.4%	97.7%	99.0%	96.1%	96.0%	
C		97.6%	-1.2%	-0.7%	23.2%	66.7%	90.2%	89.1%	97.7%	88.9%	91.7%	
D		105.9%	-1.2%	-0.1%	27.7%	78.3%	97.2%	95.0%	104.4%	105.2%	93.3%	
E		108.6%	-1.5%	0.4%	23.7%	74.0%	102.5%	97.7%	98.5%	100.9%	106.7%	
F		98.7%	-1.7%	-0.4%	22.4%	88.5%	94.0%	94.5%	100.1%	108.2%	96.5%	
G		97.1%	-0.9%	-0.9%	21.6%	84.7%	106.5%	105.2%	100.9%	108.4%	97.1%	
H												

851
852

858

859

TEST CHEMICAL															
Test Facility : A			Study Number : A1												
Chemical Code : SLS			96-Well Plate ID : A11												
2 nd Chem. Code* : 11			Experiment ID : XX												
* Testing Facility Accession Code, if applicable															
PREPARATION OF TEST CHEMICAL															
Solvent:			Medium			Dilution factor:				1.4					
Solvent Conc. (% v/v) in dosing solutions :			N/A			Highest Stock Conc.:				20,000 µg/mL					
Aids used to dissolve : <input type="checkbox"/> Vortexing <input type="checkbox"/> sonication <input type="checkbox"/> heating to 37C															
pH (highest medium stock or 2X dosing solution) : 8.0															
Medium Clarity/Color (highest 2X dosing solution): clear red						If ppt, note lowest conc.:									
Concentration Series (µg/mL)															
C1		C2		C3		C4		C5		C6		C7		C8	
100		71.4		51.0		36.4		26.0		18.6		13.3		9.49	
Positive Control (SLS) 100 - 9.49 µg/mL															
CELL LINE/TYPE															
Name: BALB/c 3T3			Supplier: ATCC			Lot No. not provided									
Passage No.: 69			Passage No. in Assay: 75			Proliferating/frozen 24-May-02									
CELL CULTURE CONDITIONS															
Medium: DMEM			Supplier:			Lot No.:									
Serum: NCS			Supplier:			Lot No.:									
Serum Conc.:			Growth Medium: 10%			Treatment Medium: 0%									
TEST ACCEPTANCE CRITERIA															
No. of values >50% and <100%:			3			No. of values >0% and <50%:			1		Accept?	YES			
VC: % Difference between Col 2 and mean VC.:						-3%				Accept?	YES				
PC: Hill Function R ² Value of SLS:						0.99				Accept?	YES				
PC: IC ₅₀ of SLS:						43.2 µg/mL				Accept?	YES				
TIMELINE															
Cell Seeding Date			Dose Application Date			OD ₅₅₀ Determination Date									
TEST RESULTS															
VC: Mean Corrected OD ₅₅₀ :						0.373				Hill Function R ² Value:		0.9869			
log IC ₂₀ :		1.551E+00 µg/mL		log IC ₅₀ :		1.635E+00 µg/mL		log IC ₈₀ :		1.718E+00 µg/mL					
IC ₂₀ :		3.56E+01 µg/mL		IC ₅₀ :		4.32E+01 µg/mL		IC ₈₀ :		5.22E+01 µg/mL					
Test Chemical F.W. : 288.4															
IC ₂₀ : 0.12331183 mM			IC ₅₀ : 0.1496252 mM			IC ₈₀ : 0.18113599 mM									

859
860
861
862
863
864
865
866
867
868
869
870
871
872
873
874
875
876
877
878
879
880
881
882
883
884
885
886
887
888
889
890
891
892
893
894
895
896
897
898
899
900
901
902
903

ANNEX II

TEST METHOD PROCEDURE Solubility Determination of Test Substances

PROPOSAL

This procedure was designed to identify the solvent that would provide the highest soluble concentration of a test substance so there would be uniform availability of the substance to cells used for *in vitro* basal cytotoxicity testing. The solubility exercises can be performed in a routine and repeatable manner and provide guidelines to effectively prepare test substances for toxicity testing in the NRU test methods.

TEST SYSTEM

The solubility test procedure is based on attempting to dissolve substances in various solvents with increasingly rigorous mechanical techniques. The solvents to be used, in the order of preference, are cell culture medium, DMSO, and ETOH. Determination of whether a substance has dissolved is based entirely on visual observation for the purposes of this protocol. A substance has dissolved if the solution is clear and shows no signs of cloudiness or precipitation.

PROCEDURES

Preparation of the keratinocyte medium and supplements will follow all procedures in the NHK NRU protocol.

Materials – see Section VI.A

Preparations of Media and Solutions – see Section VI.B

All solutions glassware, pipettes, etc., should be sterile and all procedures should be carried out under aseptic conditions and in the sterile environment of a laminar flow cabinet (biological hazard standard). All methods and procedures should be adequately documented.

Determination of Solubility

- Solubility should be determined in a step-wise procedure that involves attempting to dissolve a test substance at a relatively high concentration with the sequence of mechanical procedures specified in **Mechanical Procedures**. **Table 1** and **Figures 1** and **2** illustrate the step-wise procedures. The hierarchy of preference of solvent for dissolving test substances is medium, DMSO, and then ETOH. If the substance does not dissolve in the solvent, the volume of solvent is increased so as to decrease the test substance concentration by a factor of 10, and then the sequence of mechanical procedures are repeated in an attempt to solubilize the substance at the lower concentrations.
- For testing solubility in medium, the starting concentration is 200,000 µg/ml (i.e., 200 mg/mL) in Tier 1, but for DMSO and ETOH the starting concentration is 200,000 µg/ml (i.e., 200 mg/mL) in Tier 3.

904

905 *Methods*906 A. Tier 1 begins with testing 200 mg/mL in Routine Culture Medium (see **Table 1**).

907 1. Weigh approximately 100 mg (100,000 µg) of the test substance into a glass tube.

908 Document the substance weight.

909 2. Add approximately 0.5 mL of medium into the tube so that the concentration is

910 200,000 µg/ml (200 mg/mL).

911 3. Mix the solution as specified in **Mechanical Procedures**. If complete solubility is

912 achieved, then additional solubility procedures are not needed.

913

914 B. If the test substance is insoluble in Tier 1 at 200 mg/mL, then proceed to Tier 2.

915 1. Weigh approximately 10 mg (10,000 µg) of the test substance into a glass tube.

916 Document the substance weight.

917 2. Add approximately 0.5 mL of medium into the tube so that the concentration is

918 20,000 µg/ml (20 mg/mL).

919 3. Mix the solution as specified in **Mechanical Procedures**. If complete solubility is

920 achieved, then additional solubility procedures are not needed.

921

922 C. If the test substance is insoluble in Routine Culture Medium, proceed to Tier 3.

923 1. Add enough medium, approximately 4.5 mL, to attempt to dissolve the substance at 2

924 mg/mL by using the sequence of mixing procedures. If the test substance dissolves in

925 medium at 2 mg/mL, no further procedures are necessary.

926 2. If the test substance does NOT dissolve in medium, weigh out approximately 100 mg

927 test substance in a second glass tube and add enough DMSO to make the total volume

928 approximately 0.5 mL (for 200 mg/mL) and mix the solution as specified in

929 **Mechanical Procedures**.

930 3. If the test substance does not dissolve in DMSO, weigh out approximately 100 mg

931 test substance in another glass tube and add enough ETOH to make the total volume

932 approximately 0.5 mL (for 200 mg/mL) and mix the solution as specified in

933 **Mechanical Procedures**.

934 4. If the substance is soluble in either solvent, no additional solubility procedures are

935 needed.

936

937 D. If the substance is NOT soluble in Routine Culture Medium, DMSO, or ETOH at Tier 3,

938 then continue to Tier 4 in **Table 1**.

939 1. Add enough solvent to increase the volume of the three (or four) Tier 2 solutions by

940 10 and attempt to solubilize again using the sequence of mixing procedures. If the

941 test substance dissolves, no additional solubility procedures are necessary.

942 2. If the test substance does NOT dissolve, continue with Tier 5 and, if necessary, Tier 6

943 using DMSO and ETOH.

944 3. Tier 5 begins by diluting the Tier 4 samples with DMSO or ETOH to bring the total

945 volume to 50 mL. The mixing procedures are again followed to attempt to solubilize

946 the substance.

947 4. Tier 6 is performed, if necessary, by weighing out another two samples of test
948 substance at ~10 mg each and adding ~50 mL DMSO or ETOH for a 200 µg/mL
949 solution, and following the mixing procedures.

950
951

Example

- 952 • If complete solubility is not achieved at 20,000 µg/mL in Routine Culture Medium at
953 Tier 2 using the mixing procedures, then the procedure continues to Tier 3 by diluting the
954 solution to 5 mL with medium and mixing again.
- 955 • If the substance is not soluble in Routine Culture Medium, two samples of ~ 100 mg test
956 substance are weighed to attempt to solubilize in DMSO and ETOH at 200,000 µg/mL
957 (i.e., 200 mg/mL). Solutions are mixed following the sequence of procedures prescribed
958 in **Mechanical Procedures** in an attempt to dissolve.
- 959 • If solubility is not achieved at Tier 3, then the solutions prepared in Tier 3 are diluted by
960 10 so as to test 200 µg/mL in media, and 20,000 µg/mL in DMSO and ETOH. This
961 advances the procedure to Tier 4. Solutions are again mixed in an attempt to dissolve.
- 962 • If solubility is not achieved in Tier 4, the procedure continues to Tier 5, and to Tier 6 if
963 necessary (see **Figures 1 and 2** and **Table 1**).

964

MECHANICAL PROCEDURES

- 965
- 966 A. The following hierarchy of mixing procedures will be followed to dissolve the test
967 substance:
- 968 1. Add test substance to solvent as in Tier 1 of **Table 1**. (Test substance and solvent
969 should be at room temperature.)
 - 970 2. Gently mix at room temperature. Vortex the tube (1 –2 minutes).
 - 971 3. If test substance has not dissolved, use waterbath sonication for up to 5 minutes.
 - 972 4. If test substance is not dissolved after sonication, then warm solution to 37°C for 5 -
973 60 minutes. This can be performed by warming tubes in a 37°C waterbath or in a
974 CO₂ incubator at 37°C. The solution may be stirred during warming (stirring in a
975 CO₂ incubator will help maintain proper pH).
 - 976 5. Proceed to Tier 2 (and Tiers 3-6, if necessary of **Table 1** and repeat procedures 2-4).
- 977
- 978 B. The preference of solvent for dissolving test substances is Routine Culture Medium,
979 DMSO, and then ETOH. Thus, if all solvents for a particular tier are tested
980 simultaneously and a test substance dissolves in more than one solvent, then the choice of
981 solvent follows this hierarchy. For example, if, at any tier, a substance were soluble in
982 Routine Culture Medium and DMSO, the choice of solvent would be medium. If the
983 substance were insoluble in medium, but soluble in DMSO and ETOH, the choice of
984 solvent would be DMSO for both assays.

985

986

987

987 **Table 1 Determination of Solubility in Routine Culture Medium, DMSO, or**
 988 **ETOH**
 989

Tier	1	2	3	4	5	6
Total Volume Routine Culture Medium	0.5 mL	0.5 mL	5 mL	50 mL		
Concentration of Test Substance Tier 1: Add ~ 100 mg to a tube. Add enough medium to equal Tier 1 volume. If insoluble, go to Tier 2. Tier 2: Add ~10 mg to another tube. Add enough medium to equal the first volume. Dilute to subsequent volumes if necessary.	200,000 µg/mL (200 mg/mL)	20,000 µg/mL (20 mg/mL)	2,000 µg/mL (2 mg/mL)	200 µg/mL (0.20 mg/mL)		
Total Volume DMSO/			0.5 mL	5 mL	50 mL	
Concentration of Test Substance (Add ~100 mg to a large tube. Add enough DMSO or ETOH to equal the first volume. Dilute with subsequent volumes if necessary.)			200,000 µg/mL (200 mg/mL)	20,000 µg/mL (20 mg/mL)	2,000 µg/mL (2 mg/mL)	
Total Volume DMSO/ETOH						50 mL
Concentration of Test Substance (Add ~10 mg to a large tube. Add enough DMSO or ETOH to equal 50 mL.)						200 µg/mL (0.2 mg/mL)
EQUIVALENT CONCENTRATION ON CELLS	100,000 µg/mL (100 mg/mL)	10,000 µg/mL (10 mg/mL)	1000 µg/mL (1 mg/mL)	100 µg/mL (0.1 mg/mL)	10 µg/mL (0.01 mg/mL)	1 µg/mL (0.001 mg/mL)

990
 991 [NOTE: The amounts of test substance weighed and Routine Culture Medium added may be
 992 modified from the amounts given above, provided that the targeted concentrations specified
 993 for each tier are tested.]
 994

994

FIGURE 1 SOLUBILITY STEP-WISE (TIERED) PROCEDURE**TIER 1**

STEP 1:	200 mg/mL test substance (TS) in 0.5 mL Routine Culture Medium <ul style="list-style-type: none"> • if TS soluble in medium, then <u>STOP</u>. • if TS insoluble in medium, then go to STEP 2.
---------	---

TIER 2

STEP 2:	20 mg/mL TS in 0.5 mL Routine Culture Medium <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 3.
---------	---

TIER 3

STEP 3:	200 mg/mL TS in DMSO <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 200 mg/mL in ETOH. <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • If TS insoluble, go to STEP 4.
---------	--

TIER 4

STEP 4:	0.2 mg/mL TS in medium (one or both) – increase volume from STEP 2 by 10 (i.e., to 50 mL) <ul style="list-style-type: none"> • if TS soluble in both media, then <u>STOP</u>. • if TS insoluble in one medium, test at 20 mg/mL in DMSO – increase volume from STEP 3 by 10 (i.e., to 5 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 20 mg/mL in ETOH – increase volume from STEP 3 by 10 (i.e., to 5 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 5.
---------	--

TIER 5

STEP 5:	2 mg/mL TS in DMSO – increase volume from STEP 4 by 10 (i.e., to 50 mL) <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 2 mg/mL in ETOH – increase volume from STEP 4 by 10 (i.e., to 50 mL). <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, then go to STEP 6.
---------	---

TIER 6

STEP 6:	0.2 mg/mL TS in 50 mL DMSO <ul style="list-style-type: none"> • if TS soluble, then <u>STOP</u>. • if TS insoluble, test at 0.2 mg/mL in 50 mL ETOH <ul style="list-style-type: none"> • <u>STOP</u>
---------	--

995

Figure 2 Solubility Flow Chart

Tier	1		2		3		4		5		6
Concentration in Medium	Start Here 200 mg/mL	Incomplete solubility →	20 mg/mL	Incomplete solubility →	2 mg/mL		0.20 mg/mL				
Concentration in DMSO					200 mg/mL		20 mg/mL		2 mg/mL		0.2 mg/mL
Concentration in ETOH					200 mg/mL	Incomplete solubility →	20 mg/mL	Incomplete solubility →	2 mg/mL	Incomplete solubility →	0.2 mg/mL End
Concentration on Cells	100 mg/mL		10 mg/mL		1 mg/mL		0.1 mg/mL		0.01 mg/mL		0.001 mg/mL

Flowchart details: The chart shows a 12-column grid. The top row (Tier) is shaded. The second row (Concentration in Medium) starts at 200 mg/mL in Tier 1. Arrows labeled 'Incomplete solubility' point from Tier 1 to 2 (20 mg/mL), and from Tier 2 to 3 (2 mg/mL). From Tier 3, an arrow points to Tier 4 (0.20 mg/mL). The third row (Concentration in DMSO) has values: 200 mg/mL (Tier 3), 20 mg/mL (Tier 4), 2 mg/mL (Tier 5), and 0.2 mg/mL (Tier 6). Downward arrows labeled 'Incomplete solubility' are present in the DMSO row for Tiers 3, 4, 5, and 6. The fourth row (Concentration in ETOH) has values: 200 mg/mL (Tier 3), 20 mg/mL (Tier 4), 2 mg/mL (Tier 5), and 0.2 mg/mL (Tier 6, labeled 'End'). Horizontal arrows labeled 'Incomplete solubility' point from Tier 3 to 4, 4 to 5, and 5 to 6. Downward arrows labeled 'Incomplete solubility' are present in the ETOH row for Tiers 3, 4, and 5. The bottom row (Concentration on Cells) is shaded and has values: 100 mg/mL (Tier 1), 10 mg/mL (Tier 2), 1 mg/mL (Tier 3), 0.1 mg/mL (Tier 4), 0.01 mg/mL (Tier 5), and 0.001 mg/mL (Tier 6).

Notes: NHK medium - Keratinocyte Growth Medium (e.g., KGM® from Cambrex) for normal human keratinocytes.

ANNEX III

TEST METHOD PROCEDURE

Prequalification of Normal Human Epidermal Keratinocyte Growth Medium

This annex provides the guidelines and testing requirements for prequalifying manufacturer lots of Keratinocyte Basal Medium and the medium supplements for use with the Test Method Protocol for the NHK Neutral Red Uptake (NRU) Cytotoxicity Test. The medium and supplements should be tested so as to demonstrate their ability to perform adequately in the recommended assay.

The Testing Facility should request the quality control (QC) test data from the manufacturer for each potential lot of medium and supplements. Based upon the QC test data, purchase and test the one or two most current lots of medium and supplements that appear to have the potential to support NHK cultures according to the requirements of the aforementioned protocol.

TEST SYSTEM

The NHK NRU test is performed to analyze NHK growth characteristics and the *in vitro* toxicity of sodium lauryl sulfate (SLS), as measured by the IC_{50} , with each NHK medium/supplements being tested.

Every combination of medium/supplements expected to be used should be tested. Potential medium testing/supplement combinations are:

- One lot of medium/one lot of supplements: Test the lot of medium using the lot of supplements.
- Two or more lots of medium/one lot of supplements: Test each lot of medium using the one lot of supplements.
- One lot of medium/two or more lots of supplements: Test the lot of medium using each lot of supplements.

NHK cultures should be established using each medium/supplement combination to be tested, and should be subcultured on three different days into 96-well plates (1 plate per day) for three subsequent SLS cytotoxicity tests using each test medium/supplement combination along with a control medium (if available) for which performance has been previously established.

PROCEDURES

Prequalification of the keratinocyte medium and supplements will follow all procedures in the NHK NRU protocol.

Materials – see Section VI.A

Preparations of Media and Solutions – see Section VI.B

Methods – see Section VI.C

- NHK cultures should be established with cryopreserved cells seeded into individual tissue culture 25 cm² flasks using a proven medium/supplement combination (i.e., the “control” medium) and each test medium/supplement combination.
- Suspend freshly thawed cells initially into 9 mL of control medium and then add the cell suspension to 25 cm² culture flasks containing pre-warmed control or test medium. Cell seeding densities (1 flask/density/medium) of 1 x 10⁴, 5 x 10³, and 2.5 x 10³ are recommended.
- The cells should be subcultured on three different days into 96-well plates for three subsequent NRU tests (three test plates total [one plate per day] for each medium/supplement combination and each control).

Flask	Culture: 1 Test Plate and 1 Control Plate	Application of SLS
10 ⁴ cells/mL)	Day A	Day X
10 ³ cells/mL)	Day B	Day Y
2.5 x 10 ³ cells/mL)	Day C	Day Z

- Subculturing the cells and application of the SLS will follow procedures in the protocol in reference to appropriate cell confluency. Cell numbers should be recorded for each flask prior to subculturing to the 96-well plates.

[Note: Use of a control medium assumes that the Testing Facility has recent experience with a medium/supplement combination proven to support adequate NHK growth and provide adequate sensitivity to SLS. It is not absolutely necessary to use a control medium.]

Doubling Time – see Section VI.C.5

A doubling time experiment may be considered as an additional quality assurance check.

Preparation of SLS – see Sections VI.D.1.a, b, and d

Preparation of SLS concentrations/dilutions should follow the main experiment (definitive assay) procedures specifically for testing compounds in Routine Culture Medium as outlined in **Section VI.D.3.b**. The concentrations/dilutions should be the same or similar to those used previously with control medium/supplements. SLS concentration ranges used by three laboratories in the NICEATM/ECVAM validation study were 20.0 µg/mL – 1.4 µg/mL and 10.0 µg/mL – 0.6 µg/mL.

Test Procedure – see Sections VI.E.1, E.2, and E.4

The C₁ test concentration will be the highest SLS concentration and C₈ the lowest concentration. Cells cultured in control medium and in each test medium/supplement combination should be tested in parallel for their sensitivity to SLS (see **Methods**). Each of the three test plates of the new medium/supplement combinations is considered a replicate test plate.

Microscopic Evaluation – see Section VI.E.3

Changes in morphology of the cells due to cytotoxic effects of the SLS (prior to measurement of NRU) should be recorded as per procedures outlined in **Section VI.E.3**. In addition to the general microscopic evaluation of the cell cultures, the Study Director should make the following specific observations:

General culture observations

- rate of proliferation (e.g., rapid, fair, slow)
- percent confluence (e.g., daily estimate);
- number of mitotic figures (e.g., average per field);
- contamination (present/not present)

Cell morphology observations

- overall appearance (e.g., good, fair, poor)
- colony formation (e.g., tight/defined, fair, loose/migrating)
- distribution (e.g., even/uneven)
- abnormal cells (e.g., enlarged, vacuolated, necrotic, spotted, blebby - [average per field])

Data Analysis and Test Evaluation see Sections VI.E.5 and VI.F

Test Acceptance Criteria in section VI.E.5 will be used to determine acceptability of a test plate. Other criteria that should be considered by the Study Director includes the following:

- Mean corrected OD₅₄₀₋₅₅₀ of the VCs [Note: The target range for corrected mean OD_{540 ± 10nm} = 0.248 - 1.123 for the VCs, but it is not a test acceptance criterion (range = mean OD ± 2.5 standard deviations; mean = 0.685; SD = 0.175; N = 114).]
- Cell morphology and confluence of the VCs at the end of the 48 hour treatment
- Doubling time

The Study Director should utilize all observed growth characteristics and test results in addition to comparison of results to the media manufacturer's QC data to determine whether the medium/supplements combinations perform adequately. The Testing Facility should request that the manufacturer reserve a portion of an acceptable lot based on estimates of media need.

[This Page Intentionally Left Blank]