X10A SAFETY CHECKLIST				N.	BROOKHAVEN NATIONAL LABORATORY ATIONAL SYNCHROTRON LIGHT SOURCE
Rev: B	Effective:5/24/0	7	Page 1 of	1 N	fumber: LS-SCL-0021
Reviewed by:	Thomas McDonald	Reviewed by	γ:	·	Approved by: J. Aloi

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

## X10A CAN ONLY OPERATE IN THE FOCUSED MONOCHROMATIC BEAM CONFIGURATION

	1. Valid Padlock Index and log, no pad	llocks open					
	2. Bremsstrahlung shield (BS #1) in place and banded as per photo						
	3. Exclusion Zone (EZ #1) wire mesh in place.						
	4. Exclusion Zone (EZ #2) wire mesh i	in place					
	5. Bremsstrahlung shield (BS #2) in pla	ace and banded as per photo					
	6. Exclusion Zone (EZ #3) wire mesh i	1 1					
	7. Bremsstrahlung shield (BS #3) in pla	1					
	8. Exclusion Zone (EZ #4) wire mesh i	1 1					
	9. Slit view port covered with lead	r					
	10. Mirror shields (3) in place, scatter s	shielding as per photo					
	11. Exclusion Zone (EZ #5) wire mesh in						
	12. Bremsstrahlung shield (BS #4) in p	L. Company of the com					
	13. Exclusion Zone (EZ #6) wire mesh						
	14. Exclusion Zone (EZ #7) wire mesh	1					
	15. Bremsstrahlung shield (BS #5) in p	1					
	16. Exclusion Zone (EZ #8) wire mesh						
	17. Exclusion Zone (EZ #9) wire mesh	±					
	18. Slit view port covered with lead	an piuce					
	19 Monochromator leaded glass winder	ows (3) shielded as per photo					
	20. Exclusion Zone (EZ #10) wire mes						
	21. Bremsstrahlung shield (BS #6) in p	•					
		The property of the property o					
OPO	CO / Beamline Staff:	Date:					
	USER SYSTEM CHECKS						
	<ul> <li>A. Beamline Isolation Valve O</li> </ul>	pen					
	O B. Mirror Isolation valve open						
	O C. Monochromator isolation va	llve open					
	Beamline Rep.:	Date:					