

## Post-BioAnalyzer DNase Treatment and Cleanup of RNA

**\*\*Note:** RNeasy Mini Protocol for RNA Cleanup should be used on samples that the NMG requests **DNase treatment** of samples after QC of RNA

Qiagen RNeasy Mini Kit – Cat # 74903 or 74904

RNase-Free DNase Set – Qiagen Cat# 79254 (50 digestions)

- o a maximum of 100µg RNA can be used in the RNA cleanup protocol
- o add β-ME to Buffer RLT before use (10µl/ml RLT)
- o prepare DNase I stock solution before using for the first time. Dissolve the solid DNase I (1500 Kunitz units) in 550 µl RNase-free water provided. Mix gently by inverting the tube. **Do not vortex.**
- o for long-term storage of DNase I, divide it into single-use aliquots, and store at –20°C for up to 9 months. Do not refreeze the aliquots after thawing.
- o perform all steps of the protocol at RT. During the procedure, work quickly.
- o all centrifugation steps are performed at 20-25°C, ≥10,000 rpm.

1. Adjust sample to a volume of 100µl with **RNase-free water**.
2. Add **350µl Buffer RLT** and mix thoroughly.
3. Add **250µl 96-100% EtOH** to the diluted RNA, and mix thoroughly by pipetting.
4. Apply the **sample (700µl) to an RNeasy mini column** placed in a 2ml collection tube (supplied). Centrifuge **15 s**, ≥10,000 rpm. Discard the flow-through and collection tube.
5. Transfer column to a new 2 ml collection tube. Pipet **350 µl Buffer RW1** (Buffer **BR3**, if blood RNA) into column, and centrifuge for **15 s** at ≥10,000 rpm. Discard the flow-through.
6. Add 10 µl DNase I stock solution to 70 µl Buffer RDD for each digestion (example: for 10 samples, add 100 µl DNase I stock solution to 700 µl Buffer RDD). Mix by **gently flicking** the tube, and centrifuge briefly to collect residual liquid from the sides of the tube. **Do not vortex.**
7. Pipet **80µl DNase I** incubation mix directly onto the spin-column membrane, and place on the benchtop (20-30°C) for **15 min**.
8. Pipet **350µl Buffer RW1** (or **BR3**, if blood RNA) into column. Centrifuge for **15 s**, ≥10,000 rpm. Discard the flow-through.
9. Pipet **500µl Buffer RPE** into column. Centrifuge **15 s**, ≥10,000 rpm. Discard the flow-through.
10. Add another **500µl Buffer RPE** onto the column. Centrifuge **2 min**, ≥10,000 rpm.
11. Transfer the RNeasy column to a new 2 ml collection tube. Centrifuge **1 min**, full speed.
12. Transfer the RNeasy column to a new 1.5 ml collection tube (supplied). To **elute**, pipet **30 - 50µl RNase-free water** directly onto the RNeasy silica-gel membrane. Centrifuge **1 min**, ≥10,000 rpm.
13. If the expected RNA yield is >30µg, repeat the elution step (step 8) with a 2<sup>nd</sup> volume of RNase-free water. Elute into the same collection tube.
  - To obtain a higher total RNA concentration, use the 1<sup>st</sup> eluate in this 2<sup>nd</sup> elution step (pipet the 1<sup>st</sup> eluate from the bottom of the collection tube and pipet it directly onto the RNeasy silica-gel membrane). Centrifuge 1 min, ≥10,000 rpm.
14. Measure the concentration and 260/280 on the NanoDrop and record.