

## QCP7

### AUTOMATIC PIPET QUALITY CONTROL

#### 1.0 PURPOSE

To describe the procedure for operational accuracy check and quality control of laboratory pipettes.

#### 2.0 RESPONSIBILITIES

The Laboratory Manager or designee is responsible for assuring that this procedure is implemented on a quarterly basis.

Laboratory staff is responsible for following this procedure.

#### 3.0 PROCEDURE

##### 3.1 Equipment

3.1.1 Micro pipettes: adjustable and nonadjustable  
All varieties 20 - 1,000  $\mu\text{L}$

3.1.2 Macro pipettes: adjustable  
All varieties 1 -10 mL

3.1.3 Analytical balance  
Mettler AG 204 analytical balance or equivalent

3.1.4 Small and medium size weighing pans

3.1.5 Pipette tips, 4 sizes

3.1.6 Reagent grade water

##### 3.2 Nonadjustable exact micro pipette

3.2.1 Place small weighing pan on balance and tare.

3.2.2 Pipette reagent grade water into pan.

3.2.3 Record weight.

3.2.4 Repeat procedure for each pipette.

### 3.3 Adjustable exact micro pipette

3.3.1 Place small weighing pan on balance and tare.

3.3.2 Select pipette setting as appropriate to evaluate full adjustment

Example: P200: 50  $\mu$ L, 100  $\mu$ L, 150  $\mu$ L, 200  $\mu$ L

P1000: 200  $\mu$ L, 250  $\mu$ L, 500  $\mu$ L, 1000  $\mu$ L

3.3.3 Pipette reagent grade water into pan.

3.3.4 Record weight.

3.3.5 Repeat procedure one time for each pipette setting.

3.3.6 If weight deviates >3 percent, remove from service.

### 3.4 Adjustable gross macro pipette

3.4.1 Place medium weighing pan on balance and tare.

3.4.2 Select pipette setting.

**Note: 5 mL: 1 mL, 3 mL, 5 mL**

**10 mL: 5 mL, 8 mL, 10 mL**

3.4.3 Pipette reagent grade water into pan.

**Note: Weights are not recorded.**

3.4.4 Repeat procedure one time for each pipette setting.

3.4.5 If more than 10 percent deviation, affix a tag with an explanation of this problem, and remove from service for cleaning, part replacement and/or further check.

3.4.6 Repeat check procedure after cleaning.