

ALKALINE DUSTS

7401

NaOH, KOH, LiOH, (NaOH) and basic salts	MW : 40.00 (NaOH); 56.11 (KOH) 23.95 (LiOH)	CAS: 1310-73-2 1310-58-3 1310-65-2	RTECS: WB490000 TT2100000 (KOH) OJ6307070 (LiOH)
---	---	--	--

METHOD: 7401, Issue 2

EVALUATION: FULL

Issue 1: 15 February 1984

Issue 2: 15 August 1994

OSHA : 2 mg/m³ (NaOH)
NIOSH: C 2 mg/m³/15 min (NaOH); Group I Pesticide
ACGIH: C 2 mg/m³ (NaOH)

PROPERTIES: basic, hygroscopic, caustic solids and aerosols; VP not significant

SYNONYMS: alkali; caustic soda; lye; sodium hydroxide; potassium hydroxide

SAMPLING	MEASUREMENT
<p>SAMPLER: FILTER (1-μm PTFE membrane)</p> <p>FLOW RATE: 1 to 4 L/min</p> <p>VOL-MIN: 70 L @ 2 mg/m³ -MAX: 1000 L</p> <p>SHIPMENT: routine</p> <p>SAMPLE STABILITY: at least 7 days @ 25 °C [1,2]</p> <p>BLANKS: 2 to 10 field blanks per set</p>	<p>TECHNIQUE: ACID-BASE TITRATION</p> <p>ANALYTE: OH⁻ (alkalinity)</p> <p>EXTRACTION: 5.00 mL 0.01 N HCl, 15 min under nitrogen with stirring</p> <p>TITRATION: 0.01 N NaOH under nitrogen, endpoint by pH electrode</p> <p>CALIBRATION: 0.01 N NaOH standardized with 0.01 N HCl</p> <p>RANGE: 0.14 to 1.9 mg (as NaOH) per sample [1]</p> <p>ESTIMATED LOD: 0.03 mg per sample (as NaOH) [1] (7 x 10⁻⁴ moles of alkalinity)</p> <p>PRECISION (\hat{S}_r): 0.033 @ 0.38 to 1.5 mg NaOH per sample [1]</p>
ACCURACY	
<p>RANGE STUDIED: 0.76 to 3.9 mg/m³ [1] (360-L samples)</p> <p>BIAS: 5.6%</p> <p>OVERALL PRECISION ($\hat{S}_{r,T}$): 0.062 [1]</p> <p>ACCURACY: \pm 16.2%</p>	

APPLICABILITY: The working range is 0.4 to 5.4 mg/m³ for a 360-L air sample. The method measures total alkalinity of alkali hydroxides, carbonates, borates, silicates, phosphates, and other basic salts, expressed as equivalents of NaOH.

INTERFERENCES: Carbon dioxide in the air may react with alkali on the filter to produce carbonates but does not interfere when titrated. The carbonates will produce the equivalent amount of strong alkali that was consumed on the filter [1]. Acid aerosols may neutralize the sample, if present, producing a negative interference.

OTHER METHODS: This revises Methods S381 [2] and P&CAM 241 [3].

REAGENTS:

1. Sodium carbonate, primary standard grade.
2. Hydrochloric acid stock solution, 0.1 N. Standardize with sodium carbonate primary standard.
3. Dilute hydrochloric acid, 0.01 N. Dilute 10.0 mL 0.1 N stock HCl to 100 mL in a volumetric flask with distilled water.
4. Water, distilled, CO₂-free. Boil and cool under N₂ or bubble nitrogen through distilled water for 30 min. Store with an Ascarite trap.
5. Nitrogen, compressed.
6. Sodium hydroxide, 50% w/v.* Dissolve 50 g NaOH in CO₂-free distilled water and dilute to 100 mL.
7. Stock sodium hydroxide, 0.1 N. Dilute 8 mL 50% NaOH to 1.0 L with CO₂-free distilled water. Store under Ascarite or other CO₂-absorbing trap.
8. Working sodium hydroxide solution, 0.01 N. Dilute 10 mL stock (0.1 N NaOH) to 100 mL with CO₂-free distilled water.
9. Standard buffer solutions, pH 4 and 7.

* See Special Precautions

EQUIPMENT:

1. Sampler: 37-mm diameter PTFE membrane filter (Millipore, Fluoropore or equivalent), 1.0- μ m pore size, supported by a cellulose backup pad in a cassette filter holder.
2. Personal sampling pump, 1 to 4 L/min, with flexible connecting tubing.
3. pH meter with pH electrode and recorder.
4. Titration vessel, 150 to 200 mL beaker, flask or jar with cover containing openings for the pH electrode and N₂ inlet and outlet.
5. Stirrer, magnetic, and stir bar.
6. Glass rod, ca. 5-mm diameter and 10 cm long to hold filter under liquid surface in titration vessel.
7. Pipets, 5- and 10-mL.
8. Volumetric flasks, 100-mL and 1-L.
9. Burets, 50-mL, readable to 0.1 mL.
10. Tweezers.

SPECIAL PRECAUTIONS: NaOH solutions are corrosive to tissue [4]. Handle with care.

SAMPLING:

1. Calibrate each personal sampling pump with a representative sampler in line.
2. Sample at an accurately known flow rate between 1 and 4 L/min for a sample size of 70 to 1000 L. Do not exceed a filter loading of ca. 2 mg total dust.

SAMPLE PREPARATION:

3. Transfer the sample filter to a titration vessel with tweezers. Place the filter face down in the titration vessel.
4. Place the end of a glass rod in the center of the filter to maintain the filter below the liquid surface during the analysis.
5. Cover the titration vessel, add 5.00 mL 0.01 N HCl, start the magnetic stirrer and N₂ purge (ca. 0.1 L/min).
6. Allow to stand 15 min (with stirring).

CALIBRATION AND QUALITY CONTROL:

7. Calibrate the pH meter with pH 4 and pH 7 buffer solutions.
8. Standardize aliquots of the 0.1 N HCl stock solution with sodium carbonate in triplicate [3].
 - a. Dry 3 to 5 g primary standard grade Na₂CO₃ at 250 °C for 4 h. Cool in a desiccator.
 - b. Weigh ca. 2.5 g Na₂CO₃ to the nearest mg. Dissolve and dilute to exactly 1 L with CO₂-free distilled water. The concentration is ca. 0.05 N Na₂CO₃.
 - c. Place 5.00 mL 0.05 N Na₂CO₃ solution into a titration vessel and titrate potentiometrically to a pH of 5.

- d. Remove electrodes, rinse them into the titration vessel, and bubble N₂ gas through contents of the titration vessel for 3 to 5 min to remove dissolved CO₂.
- e. Proceed with the titration to the inflection point.
- f. Calculate the normality of the stock HCl solution:

$$N_{\text{HCl}} = \frac{(\text{g Na}_2\text{CO}_3 \text{ weighed})(\text{mL Na}_2\text{CO}_3 \text{ solution used in titration})}{(52.99)(\text{mL HCl used})}$$

9. Standardize the working (ca. 0.01 N) NaOH solution against the standardized HCl solution by following steps 8.c. and 8.e, substituting the standardized HCl stock solution for the Na₂CO₃ solution and the 0.01 N NaOH solution for the 0.1 N HCl solution. Calculate the normality of the NaOH titrating solution.

$$N_{\text{NaOH}} = (N_{\text{HCl}})(\text{mL HCl used})/(\text{mL NaOH used})$$

10. Prepare at least three spiked media blank samples to check analytical recoveries at levels expected on the field samples.

MEASUREMENT:

11. Back-titrate the excess HCl in the samples, blanks and spiked blank solutions with the standardized NaOH solution while maintaining the N₂ purge.
12. Observe the pH meter. Calculate the endpoint (mL of 0.01 N NaOH used).

CALCULATIONS:

13. Using the normality (N) and volumes of NaOH in the titration of the sample (V_{NaOH-s}) and average blank filter (V_{NaOH-b}), and the volume of air sampled, V(L), calculate the concentration, C (mg/m³), of alkalinity (as NaOH with equivalent weight = 40.0):

$$C = \frac{(V_{\text{NaOH-b}} - V_{\text{NaOH-s}}) \cdot N \cdot 4 \times 10^4}{V}, \text{ mg/m}^3.$$

EVALUATION OF METHOD:

Method S381 [1] was issued on July 8, 1977, and was validated using generated atmospheres of NaOH over the range 0.76 to 3.9 mg/m³ using a 360-L sample [1,6]. Overall precision, \hat{S}_{RT} , was 0.062 with an average recovery of 105%, representing an insignificant bias. The NaOH concentration was independently verified by analyzing additional samples for sodium by atomic absorption spectrophotometry. An overall average collection efficiency of greater than 99% was found by sampling 1 to 360 L in an atmosphere of 7.50 mg/m³ (NaOH) using backup Fluoropore filters and also by collecting two 615-L samples using backup impingers filled with water.

REFERENCES:

- [1] Backup Data Report No. S381, Sodium Hydroxide, prepared under NIOSH Contract No. 210-76-0123, available as Order No. PB 275-838 from NTIS, Springfield, VA 22161 (July 8, 1977).
- [2] NIOSH Manual of Analytical Methods, 2nd. ed., V. 4, S381, U.S. Department of Health, Education, and Welfare, Publ. (NIOSH) 78-175 (1978).
- [3] Ibid, V. 1, P&CAM 241, U.S. Department of Health, Education, and Welfare, Publ. (NIOSH) 77-157-A (1977).

- [4] Criteria for a Recommended Standard...Occupational Exposure to Sodium Hydroxide, U.S. Department of Health, Education, and Welfare, Publ. (NIOSH) 76-105 (1975).
- [5] Standard Methods for the Examination of Water and Wastewater, 15th Edition, American Public Health Association, American Water Works Association and Water Pollution Control Federation (1981).
- [6] NIOSH Research Report-Development and Validation of Methods for Sampling and Analysis of Workplace Toxic Substances, U.S. Department of Health and Human Services, Publ. (NIOSH) 80-133 (1980).

METHOD REVISED BY:

Mary Ellen Cassinelli, NIOSH/DPSE; S381 originally validated under NIOSH Contract No. 210-76-0123.