

Pall Bacterial Detection System (eBDS)  
-7 days storage (release use)

Pall Pre Storage Pooling Systems  
-5 and 7 Days of Storage

Stein Holme, Ph.D.  
Pall Life Sciences

## eBDS and 7 day storage (release use)

**Random Donor PCs (RD PCs) with the Pall Leukotrap RC-PL system have been approved (quality) for 7 days storage**

**eBDS approved Feb 2004 for QC use**

**510(k) for release use of eBDS:**

**Target date for submission February, 2005.**

- 1) Field Data with testing conducted under actual use conditions**
- 2) Post-marketing protocol**

**(Approval will allow for 7 day storage of single products already approved in terms of quality)**

# eBDS and 7 day storage (release use)

## Results of 118067 tests performed at 23 blood centers in the US - March to Nov 2004

118 (0.1%) “failed” results where:

23 confirmed true positives ( presence of bacteria in both the eBDS pouch and the platelet mother bag by culture)

76 false positives ( no presence of bacteria in neither the eBDS pouch nor the platelet mother bag by culture)

1 false negative (sepsis) with confirmed presence of bacteria (*S.epidermidis*) in the mother bag by culture)

18 were not confirmed false positives (no bacteria in the mother bag by culture, eBDS pouch not tested)

# Field Data with eBDS:

## 118067 tests performed at 23 blood centers

	<b>Tested Samples that "Failed"</b>	<b>Confirmed True Positives</b>	<b>Confirmed False Negatives</b>	<b>Confirmed False Positives*</b>	<b>Not Confirmed False Positives**</b>
<b>Number</b>	118	23	1	76	18
<b>Percentage of total tested</b>	0.100% (1/1001)	0.019% (1/5133)	0.001% (1/118056)	0.064% (1/1554)	0.015% (1/6559)


\* Pouches confirmed to be negative by culture

\*\* Pouches not tested by culture



# System for Pre Storage Leukoreduced Pooled Platelet Products

**SYSTEM (#1) FOR 5 DAYS STORAGE**  
**1.5L CLX®-HP bag, with eBDS in line.**





# System for Leukoreduced Pre Storage Pooled Platelet Products

## SYSTEM (#1) FOR 5 DAYS STORAGE

**1.5L CLX®-HP bag, with eBDS in line.**

**Suitable for pooled leukoreduced PCs from whole blood collected with the Pall Leukotrap(R) RC-PLsystem using CP2D as anticoagulant.**

**Store pools of 4-6 units (40-65 mL) of leukoreduced PCs in plasma for up to 5 days with total yields of 2.2 – 5.8 x10<sup>11</sup> plts.**

**To used with an approved bacterial detection system (eBDS)**



# Concerns/challenges with pre storage pooling of random donor PC

- Risk of elevated bacteria levels after storage
  - Sensitive bacterial detection system
- Lymphocyte activation, generation of harmful levels of cytokines, complement and clotting factors
  - Pre storage leukoreduction
- Platelet Storage Quality
  - Bag able to handle 4-6 RD PCs with a large variability in yields and volume


# Bacterial Issues with Pre storage Pooling

- Bacteria growth and final CFU levels may be higher on a pre storage pooled product as compared to a post storage pooled product at the time of transfusion
- At the time of pooling a potentially contaminated individual PC will be **diluted** in the pooled product – thus resulting in a lower CFU/mL level. This may thus challenge the sensitivity of detection with immediate sampling.



# Bacterial Growth

Generation (doubling) time of *S.epidermidis*  
in PC during storage (24-48 hrs)



# Effect of dilution with pooling

Probability of obtaining no organism in a 3 ml (eBDS) sample

<b>CFU/mL</b>	<b>1 PC</b>	<b>Pool 5 PC*</b>
5.000	0.000	0.049
4.000	0.000	0.089
3.000	0.000	0.164
2.000	0.002	0.300
1.000	0.049	0.550
0.800	0.089	0.620
0.600	0.164	0.700
0.400	0.300	0.790
0.200	0.550	0.890
0.100	0.740	0.940

\* 1 contaminated and 4 non contaminated PC




# Satisfactory In vitro and in vivo Quality and lymphocyte activation levels at 5 Days of Storage with pooled leukoreduced RDP stored in CLX- HP

## STUDY SITES:

- 1) **Dr. Joe Sweeney - Miriam Hospital, Providence –**  
*Prestorage pooled whole-blood-derived leukoreduced platelets stored for seven days, preserve acceptable quality and do not show evidence of a mixed lymphocyte reaction.*  
*Sweeney JD, Kouttab NM, Holme S et al.*  
*Transfusion. 2004 Aug;44(8):1212-9.*
- 2) **Dr. Scott Murphy - ARC Penn Jersey, Philadelphia**
- 3) **Nancy Heddle - McMaster University, Hamilton -**  
*Whole blood derived platelets stored as a pool: a randomized block noninferiority trial. Heddle NM et al. Vox sang 2004, 87(suppl.3).p 6*


# Bag Capacity

## Yield Distribution of individual processed PC (Pall RC-PL)




# Bag Capacity

## Expected Platelet Yield Distribution with a Pool of 4 PC




# Bag Capacity

## Expected Platelet Yield Distribution with a Pool of 5 PC



# PLATELET CAPACITY STUDIES

## CLX-HP – 5 DAY STORAGE



# Pall System for Leukoreduced Pre Storage Pooled Platelet Products

**SYSTEM (#1) FOR 5 DAYS STORAGE  
1.5L CLX®-HP bag, with eBDS in line.**


<i>Milestone</i>	<i>Status</i>
<b>Lymphocyte/Plasma Activation Studies</b>	<b>Completed</b>
<b>In vitro and In vivo Storage Quality</b>	<b>Completed</b>
<b>CLX –HP storage bag capacity studies (low/high yield, and low/ high volume)</b>	<b>Completed</b>
<b>Studies on eBDS bacteria testing in pooled PCs</b>	<b>February '05</b>
<b>Submit to FDA</b>	<b>March '05</b>





# Pooling Set # 2 with Filter for non leukoreduced PC - 7 Days of Storage

Better Blood, Better Care<sup>SM</sup>





# Pre storage pooling and with an inline Sensitive Bacteria Detection System

**What are the advantages in terms  
*Blood Safety and Availability?***

**The Simplicity and Affordability with  
1 bacterial test for 4-6 RD PC:**

**Enables continuous use of random  
donor platelets as an important  
source of platelets**

**Enables improved bacterial detection  
compared to current practice**