



## 6.6 ALKALINITY AND ACID NEUTRALIZING CAPACITY

**By D.B. Radtke, F.D. Wilde, J.V. Davis, and  
T.J. Popowski**

	Page
<b>Alkalinity and acid neutralizing capacity .....</b>	<b>ALK-3</b>
<b>6.6.1 Equipment and supplies .....</b>	<b>5</b>
<b>6.6.1.A Maintenance, cleaning, and storage .....</b>	<b>8</b>
<b>6.6.2 Calibration and standardization .....</b>	<b>9</b>
<b>6.6.3 Collection and processing .....</b>	<b>11</b>
<b>6.6.4 Measurement .....</b>	<b>13</b>
<b>6.6.4.A Titration procedures .....</b>	<b>14</b>
<b>Buret titrator .....</b>	<b>19</b>
<b>Digital titrator .....</b>	<b>22</b>
<b>6.6.4.B Inflection point titration method .....</b>	<b>25</b>
<b>6.6.4.C Gran function plot method .....</b>	<b>28</b>
<b>6.6.5 Troubleshooting .....</b>	<b>32</b>
<b>6.6.6 Reporting .....</b>	<b>33</b>
 <b>Illustrations</b>	
<b>6.6-1. Flow chart showing a summary of             alkalinity or ANC titration procedures .....</b>	<b>16</b>
<b>6.6-2. Graph showing example of an inflection             point titration using a buret .....</b>	<b>26</b>
<b>6.6-3. Graph showing example of inflection             point titration using digital titrator .....</b>	<b>27</b>
<b>6.6-4. Graph showing example of Gran function plot             of a bicarbonate titration using a digital titrator .....</b>	<b>30</b>

**Tables**

<b>6.6-1. Equipment and supplies used for alkalinity or ANC titrations .....</b>	<b>6</b>
<b>6.6-2. Suggested sample volume and titrant normality (<i>N</i>) for use with digital titrator at various ranges of alkalinity or ANC .....</b>	<b>13</b>
<b>6.6-3. Digital titration factors .....</b>	<b>24</b>
<b>6.6-4. Results of typical inflection point titration using a buret .....</b>	<b>26</b>
<b>6.6-5. Results of typical inflection point titration using a digital titrator .....</b>	<b>27</b>
<b>6.6-6. Example of information needed for a Gran function plot .....</b>	<b>31</b>

## ALKALINITY AND ACID NEUTRALIZING CAPACITY 6.6

Alkalinity and the acid neutralizing capacity (ANC) are determined using identical electrometric procedures involving the acidimetric titration of a sample; the only difference is that the alkalinity sample is filtered, but the ANC sample remains unfiltered. The terms alkalinity, ANC, and carbonate alkalinity are used in this manual as follows:

- ▶ **Alkalinity** applies to the acid neutralizing capacity of solutes in a water sample, reported in equivalents per liter (or milliequivalents or microequivalents). Alkalinity, thus, consists of the sum of titratable carbonate and noncarbonate chemical species in a **filtered water sample** (filter membrane of 0.45- $\mu\text{m}$  or less pore size).
- ▶ **ANC** applies to the acid-neutralizing capacity of solutes plus particulates in an **unfiltered water sample**, reported in equivalents per liter (or milliequivalents or microequivalents). (The ANC value would be equivalent to alkalinity for samples without titratable particulate matter.)
- ▶ **Carbonate alkalinity** represents the acid-neutralizing capacity of carbonate solutes ( $\text{HCO}_3^- + 2m\text{CO}_3^{2-}$ , where  $m$  is moles), reported either in equivalents per liter (or milliequivalents or microequivalents) or in milligrams per liter as a carbonate species, and titrated on a **filtered water sample**. In many aqueous systems, alkalinity is controlled by carbonate chemistry and most commonly is attributable to bicarbonate ( $\text{HCO}_3^-$ ) and less frequently to carbonate ( $\text{CO}_3^{2-}$ ).

Alkalinity: the capacity of solutes in an aqueous system to neutralize acid.

Acid Neutralizing Capacity (ANC): the equivalent sum of all bases or base-producing materials in an aqueous system that can be titrated with a strong acid to an equivalence point.

#### 4 — ALK

Alkalinity is used routinely in checking the charge balance of a solution and to gain insights on the evolution of aqueous systems. Alkalinity and ANC provide information on the suitability of water for uses such as irrigation, determining the efficiency of wastewater processes, determining the presence of contamination by anthropogenic wastes, and maintaining ecosystem health.

Alkalinity is determined on a filtered sample.

Any substance in the water sample that reacts with the strong titrant acid can contribute to the water's acid neutralizing capacity.

- ▶ Important noncarbonate contributors include organic ligands (especially acetate and propionate) and ions of hydroxide, phosphate, ammonium, silicate, sulfide, borate, and arsenate (Hem, 1985).
- ▶ Except for unusual natural waters and waters significantly affected by human activity, noncarbonate ionized contributors generally are not present in large enough quantities to affect alkalinity or ANC determinations.
- ▶ Particulate matter can be an important contributor, and must be filtered from samples to be titrated for the alkalinity determination.

TECHNICAL NOTE: Alkalinity and ANC are independent of exchange with carbon dioxide and other atmospheric gases. However, atmospheric gas exchange can alter concentrations of individual species, such as bicarbonate. Also, aeration of a sample during filtration can cause mineral precipitation on the filter—this alters alkalinity, especially in water systems closed to the atmosphere under ambient conditions.

---

## EQUIPMENT AND SUPPLIES 6.6.1

Equipment and supplies for the electrometric method to determine alkalinity and ANC are listed in table 6.6-1. The equipment must be tested before each field trip and cleaned immediately after each use.

### ***Buret, micrometer buret, digital titrator***

The buret provides good accuracy and precision when used by a trained operator.

- ▶ Micrometer burets provide better accuracy and precision than burets—they can deliver acid increments to 0.0001 mL and are available commercially (for example, Gilmont™ micrometer burets).
- ▶ The digital titrator is not as accurate as a buret because it was designed for endpoint titrations. The digital titrator is popular because it is more convenient and less fragile than a buret and keeps the acid in a virtually closed system. (The Hach brand is used as an example in this document.) **Delivery tubes of clear (instead of red) plastic are recommended. Empty titrant cartridges are available.**

### ***Volumetric pipet, graduated cylinder, digital balance***

The volumetric pipet is used for dispensing the correct volume of filtered sample only.

- ▶ Use only class A “TD” pipets. TD is a calibration designation meaning “To Deliver” (TD is distinguished from “TC” or “To Contain” pipets).
- ▶ Class A pipets should not be used to aspirate or dispense solutions containing suspended particles—the small bore of this pipet tends to reject particles during aspiration and retain them during delivery (C.J. Patton, U.S. Geological Survey, written commun., 1995).

The graduated cylinder and digital balance are used for dispensing the correct volume of unfiltered (ANC) sample only. The digital balance yields higher precision than the graduated cylinder.

Table 6.6–1. Equipment and supplies used for alkalinity or ANC titrations<sup>1</sup>  
 [mL, milliliters; ANC, acid neutralizing capacity; g, gram;  $\mu\text{S}/\text{cm}$ , microsiemens per centimeter at 25°Celsius; N, normal]

**Equipment and supplies common for using either a digital titrator or a buret**

- ✓ pH meter, preferably with automatic temperature compensator (see NFM 6.4 for selection and associated supplies)
- ✓ pH electrode, calibrated, combination or equivalent, and appropriate filling solution, if required
- ✓ Thermometer, calibrated (see NFM 6.1 for selection and calibration criteria)
- ✓ Stirrer, magnetic (battery operated) or glass stir rods
- ✓ Stirring bars, Teflon™ coated, smallest size (always carry spare bars)
- ✓ Volumetric pipets, class A—25 mL, 50 mL, and 100 mL (for alkalinity)
- ✓ Graduated cylinder (for ANC). For higher precision, use a digital balance, 0.1-g accuracy, 200-g capacity, pocket-sized (available from Acculab Company)
- ✓ Pipet squeeze bulb or pipet pump
- ✓ Sample bottle, 500 mL, acid rinsed or deionized-water rinsed
- ✓ Beakers, glass—50 mL, 100 mL, and 150 mL
- ✓ Beaker, Berzelius, 300 mL, tall form, spoutless, with two- or three-hole stopper
- ✓ Deionized water (DIW) (maximum conductivity of 1  $\mu\text{S}/\text{cm}$ )
- ✓ Dispenser bottle, squeeze, for deionized water
- ✓ Filtration unit, in-line capsule or pressure unit with inert gas (for alkalinity)
- ✓ Sodium carbonate standard solution
- ✓ Safety gloves, glasses, acid spill kit, and apron
- ✓ Titrant solution, sulfuric acid solution, 0.1600*N*, 1.600*N*, and 0.01639*N* (pre-filled cartridges for use with the Hach system are available for 0.1600*N* and 1.600*N* solutions)
- ✓ Paper tissues, disposable, soft and lint free

**Equipment and supplies for using a digital titrator**

- ✓ Digital titrator and mounting assembly
- ✓ Delivery tubes, bent-stem, colorless, transparent
- ✓ Titrant cartridges, empty or prefilled according to study needs (used with Hach system)

**Equipment and supplies for using a buret**

- ✓ Calibrated buret, 25-mL capacity with 0.05-mL graduations and Teflon™ stopcock
- ✓ Calibrated buret, 10-mL capacity with 0.02-mL graduations and Teflon™ stopcock
- ✓ Micrometer buret (alternative to standard burets, for greater accuracy)
- ✓ Titrant solution, sulfuric acid solution, 0.01639*N*
- ✓ Buret stand and clamp
- ✓ Wire pen cleaner (for cleaning buret tip)
- ✓ Buret cap
- ✓ Buret meniscus reader
- ✓ Acid bottle, pump (for filling buret)

<sup>1</sup>Modify this list to meet the specific needs of the field effort.

### ***Sulfuric acid titrant***

Sulfuric acid is the titrant used routinely by USGS. The normality of sulfuric acid titrant is subject to change with time.

- ▶ 0.1600*N* or 1.600*N* solutions in prefilled cartridges for the Hach digital titrator are available from the QWSU (USGS, Ocala, Fla.). Normality of sulfuric acid in prefilled Hach titrant cartridges is monitored by the QWSU for 1 year after date of purchase; the QWSU notifies USGS personnel if a change in normality occurs. The lot number indicates the date of preparation (QWSU will translate lot numbers upon request). **QWSU recommends that Hach titrant cartridges be used within 3 months of purchase.**
- ▶ 0.01639*N* solution for the buret is available from the QWSU. Normality of 0.01639*N* titrant is not monitored by the QWSU after preparation and standardization. **Check the normality of this titrant each month.**
- ▶ Acid solutions of other normalities may be needed, depending on the sample chemistry or the ionic strength. Prepare under a fume hood. Check the normality monthly.

TECHNICAL NOTE: For samples vulnerable to precipitation reactions, a Berzelius beaker can help minimize gas exchange. Select a size of Berzelius beaker that fits the sample volume and associated titrating equipment and yet minimizes headspace above the sample. Fit the Berzelius beaker with a two- or three-hole stopper to accommodate the electrode(s), the thermometer, and the digital or buret titrator. Another option is to work in a glove box filled with an inert gas atmosphere. Oceanographers use a closed cell with an expanding plunger to avoid gas exchange (Almgren and others, 1977).

**CAUTION:** Use the safety precautions outlined on the Material Safety Data Sheets (MSDS) when handling chemicals—wear safety glasses, gloves, and protective clothing.

### 6.6.1.A MAINTENANCE, CLEANING, AND STORAGE

Proper maintenance, cleaning, and storage of the pH instrument is critical for ensuring the accuracy of alkalinity or ANC determinations, and guidance is provided in NFM 6.4.

Clean the volumetric pipets, beakers, bottles, burets, and stirring bars with hot water and nonphosphate detergent; rinse them copiously with tap water followed by deionized water. If oily or chemical residues are difficult to remove, soak the glassware and nonmetal equipment in a mild (5 percent by volume) hydrochloric acid or nitric acid solution (see Horowitz and others, 1994), and repeat the detergent wash. Store cleaned equipment wrapped or bagged in plastic until ready for use.

Reagents must not exceed their shelf life. Store reagents, as appropriate, in a dust-free cabinet, desiccator, or refrigerator. When chemicals to be used for preparation of reagents are received, mark the dates of receipt and expiration on container. When a reagent is prepared, label the container with the contents, date of preparation, expiration date, and preparer's initials. Store the 0.01639*N* standard sulfuric acid solution and filled Hach titrant cartridges in a cool, dark place (a storage cabinet or frost-free refrigerator). Seal the filled cartridges in plastic bags to avoid moisture loss or gain.

Some of the procedures recommended herein for equipment operation may be out of date if the equipment being used is different from that described or incorporates more recent technological advances—follow the manufacturer's instructions.



## CALIBRATION AND STANDARDIZATION 6.6.2

Calibration is required for the pH instrument—follow instructions in NFM 6.4.

Alkalinity and ANC determinations require standardization if the 0.01639*N* sulfuric acid titrant (buret method) is used. In the procedure that follows (Fishman and Friedman, 1989), the reagent concentrations were selected so that the titrant volume would correspond exactly to the bicarbonate equivalent in milligrams per liter of sample.

***Prepare fresh carbon dioxide-free deionized water (DIW):***

1. Boil the DIW in a beaker for 15 minutes.
2. Cool the DIW to room temperature. Cover the beaker to minimize atmospheric contamination while the DIW cools.

An alternative to the method in steps 1 and 2 is to infuse the DIW with compressed helium or nitrogen gas for about 10 minutes.

***Prepare a fresh standard solution of sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>)(1.00 mL = 1.00 mg HCO<sub>3</sub>):***

1. Dry 1.0 g primary standard Na<sub>2</sub>CO<sub>3</sub> at 150 to 160°C for 2 hours.
2. Cool in a desiccator; weigh out 0.8686 g Na<sub>2</sub>CO<sub>3</sub>; add to a 1-L volumetric flask.
3. Dilute with carbon dioxide-free DIW to the 1-L mark.

***Prepare sulfuric acid titration solution:***

1. Add 0.5 mL concentrated H<sub>2</sub>SO<sub>4</sub> (specific gravity 1.84 g/mL) to 950 mL DIW.
2. Mix thoroughly.

10 — ALK

To standardize the sulfuric acid, follow steps 1 through 9 below. To check the normality of a sulfuric acid titration solution, follow steps 3 through 6 below.

1. Calibrate the pH system (follow the instructions in NFM 6.4).
2. Prepare a sodium carbonate standard solution as outlined above.
3. Pipet 25 mL of sodium carbonate standard solution into a 100-mL beaker.
4. Titrate with the sulfuric acid solution to pH 4.5.
5. Record the volume of titrant used.
6. Determine the normality of acid by use of the following equation:

$$N = (25) (0.01639) / mL_a$$

where

$N$  = normality, and

$mL_a$  = volume of sulfuric acid added in milliliters to reach pH 4.5.

7. Adjust the concentration of sulfuric acid to exactly 0.01639*N* by dilution with deionized water or by the addition of dilute acid, as indicated by the first titration.
8. Confirm the exact normality by restandardization.
9. Keep the solution in a tightly sealed 1-L glass bottle until used.
10. If the sulfuric acid titrant solution is not exactly 0.01639*N*, divide the actual normality of the H<sub>2</sub>SO<sub>4</sub> titrant by 0.01639*N* to obtain the correction factor.

**CAUTION:** Wear appropriate safety gloves, glasses, and apron when working with corrosive and oxidizing solutions.

---

## COLLECTION AND PROCESSING 6.6.3

Collect and process the samples in a manner that ensures that they represent environmental concentrations at the time of collection. Minimize the effects of wind, rain, cold, dust, and sun on the samples. Collect and process the samples in a chamber to protect them from airborne particulates.

Before collecting or processing the sample, clean the samplers, compositing and splitting devices, sample bottles, measurement vessels, and other equipment that contacts the sample (for detailed procedures refer to Horowitz and others, 1994).

- ▶ Prerinse the sample bottles with deionized water and store in sealable plastic bags until ready for field sampling (acid-rinsed sample bottles are recommended, especially for samples with low alkalinity or ANC). Field rinse the bottle(s) three times with sample (sample filtrate for alkalinity).
- ▶ Do not field rinse the measurement vessels. Volumetric pipets and graduated cylinders should be clean and dry before use.

### ***To collect and process the sample:***

1. Filter the samples along with other anion samples, if making alkalinity determination. The 0.45- $\mu$ m flowthrough disposable filter capsule is the standard unit used by USGS. Record on field forms if a different unit or membrane is used, as this can affect the determination.
2. Fill and securely cap two 250-mL sample bottles with the sample (filtrate for alkalinity) to ensure there is enough sample to repeat the titration, to preserve the integrity of the second aliquot after the first has been opened, and to avoid losing the volume of sample needed to spillage.

3. Prevent agitation of the sample or prolonged exposure to air in order to avoid oxidation of hydrogen sulfide, ferrous iron, manganese, and prevent precipitation of mineral phases.
  - Loss of carbon dioxide (CO<sub>2</sub>) from the sample will not change the alkalinity or ANC determination, but chemical or physical reactions can cause concentrations of bicarbonate and carbonate to change within minutes.
  - Begin the titration as soon as possible, as there is less chance of chemical precipitation once acidification begins. If the titration is delayed, maintain the samples at the temperature of their ambient environment.
  - If there is a tendency for mineral precipitation, collect and process the sample in an inert gas atmosphere.

### ***Surface water***

Collect and process a representative sample according to USGS-approved methods (see Horowitz and others, 1994).

- ▶ NASQAN, BENCHMARK, and NAWQA programs require filtration of alkalinity samples through a 0.45- $\mu$ m membrane.
- ▶ To collect and process samples from anoxic lake or reservoir depth intervals, adapt procedures described for ground water.

### ***Ground water***

Collect the sample as close to the source as possible; minimize aeration of the sample; take the precautions described in step 3 above.

- ▶ Purge the well (NFM 6.0) and connect the filter unit in-line with the pump.
- ▶ Flush and fill the lines and filter unit with sample water so as to exclude air.

## MEASUREMENT 6.6.4

Alkalinity, ANC, and concentrations of bicarbonate, carbonate, and hydroxide species are determined using either the inflection point titration (IPT) method or the Gran function plot (Gran) method to analyze the titration data. Be familiar with the information in 6.6.4.B (IPT method) and 6.6.4.C (Gran method) before selecting the method to be used and before starting the sample titration.

- ▶ **The inflection point titration (IPT) method** is adequate for most waters and study needs. Difficulty in identifying the inflection points using an IPT method increases as the ratio of organic acids to carbonate species increases.
- ▶ **The Gran method** is recommended for water in which the alkalinity or ANC is expected to be less than about 0.4 meq/L (20 mg/L as CaCO<sub>3</sub>), or in which conductivity is less than 100 μS/cm, or if there are appreciable noncarbonate contributors or measurable concentrations of organic acids.

The IPT and Gran methods require electrometric titration of a sample with incremental additions of H<sub>2</sub>SO<sub>4</sub> of specified normality. Suggested combinations of titrant normality and sample volume for various ranges of alkalinity or ANC values are given in table 6.6–2. Generally, 1.600*N* acid is too strong for most samples and is used at alkalinity or ANC greater than 4.0 meq/L (200 mg/L as CaCO<sub>3</sub>).

Table 6.6–2. Suggested sample volume and titrant normality (*N*) for use with digital titrator or at various ranges of alkalinity or ANC [meq/L, milliequivalents per liter; mg/L, milligrams per liter; CaCO<sub>3</sub>, calcium carbonate; mL, milliliter; >, greater than]

ANC or alkalinity (meq/L)	Alkalinity (mg/L as CaCO <sub>3</sub> )	Sample volume (mL)	Titration normality ( <i>N</i> )
0–1.0	0–50	100	0.1600 (or lower)
1.0–4.0	50–200	50	.1600
4.0–20	200–1,000	100	1.600
>20	>1,000	50	1.600

Titrate a filtered sample if you will check the charge balance or will report the results as carbonate alkalinity.

### 6.6.4.A TITRATION PROCEDURES

Titration procedures are identical for surface-water and ground-water determinations on filtered or unfiltered aliquots of fresh to saline water samples. Become familiar with the information and detailed instructions for the buret and digital titration systems and the IPT and Gran methods before proceeding with the titration.

TECHNICAL NOTE: Fixed endpoint titration to pH 4.5 is no longer used by the USGS for reported alkalinity values because it can be less accurate than the IPT and Gran methods, particularly at small concentrations of total carbonate species and in water with significant organic and other noncarbonate contributors to alkalinity or ANC.

Before beginning titration, select the titration system to be used.

- ▶ The digital system is convenient but tends to be less precise or accurate than the buret system because of mechanical inadequacies. Good technique is requisite to produce acceptable results.
- ▶ The buret system can be cumbersome and fragile in the field, and requires experience to execute with precision and accuracy.
- ▶ A micrometer buret can achieve accuracy to one-tenth of a mg/L (routine determinations are reported to whole numbers).

Select the size of the volumetric pipet (for alkalinity sample) or the graduated cylinder or digital balance (for unfiltered ANC sample) and the measurement vessel, according to the volume of sample needed.


- ▶ 50 mL of a sample in a 100-mL beaker is selected for most routine work.
- ▶ Use 100 mL of a sample in a 150-mL beaker for low concentration of alkalinity or ANC. If the total sample volume is small, you may need to use a sample volume of 25 mL or less, although concentration is low.
- ▶ Use 25 mL or less of a sample and a 50-mL beaker for high concentration of alkalinity or ANC.

To achieve greater accuracy, use lower normality titrant, decrease the volume of acid increments, and increase the number of increments. Figure 6.6-1 provides a general summary of the electrometric titration procedures for alkalinity or ANC.

When pipetting the sample, a small amount of contents remains in the tip of class A “TD” volumetric pipets—do not blow it out, but use the following procedure:

1. Suspend the pipet tip vertically in a beaker, touching neither the walls nor the contents of the receiving vessel.
2. Allow the sample to drain freely until the liquid it contains reaches the bottom of the bulb.
3. Touch the pipet tip to the beaker wall until the flow from the pipet stops—leave the tip in contact with the beaker wall for an additional 10 seconds after the flow stops.

The buret system yields better accuracy than the digital titrator.





When titrating, stirring helps to establish a uniform mixture of sample and titrant and an equilibrium between sensors and sample.

- ▶ **If using a magnetic stirrer, stir the sample slowly and continuously**, using the smallest stir bar; avoid creating a vortex and large streaming potentials. If using a digital titrator, keep the delivery tube immersed throughout the procedure but keep the aperture of the tube away from the stir bar to avoid bleeding acid from the tube to the sample between titrant additions.
- ▶ If swirling the sample to mix, make the pH measurement as the sample becomes quiescent, after each addition of titrant.
- ▶ Avoid splashing the sample out of the beaker or onto the beaker walls. Droplets on the beaker walls can be rinsed down with deionized water. If you splash the sample out of the beaker you must start over.

**If concentrations of contributing carbonate species will be determined**—Titrate to pH of about 8.1 carefully, in small increments. If concentrations of contributing carbonate species **will not** be determined—titrate rapidly at first, adding relatively large acid increments to bring pH to about 5.5; then titrate slowly in small increments.

- ▶ **IPT method.** Titrate to at least pH 4.0 (to pH 3.5 if the alkalinity or ANC range is unknown for the water sampled or if the sample contains high concentrations of noncarbonate contributors, such as organic acids).
  - Titrate cautiously at and beyond the expected equivalence points.
  - Let the pH value stabilize before the next addition of titrant.
- ▶ **Gran method.** Titrate to pH 3.5, or to pH 3.0 or less if the alkalinity or ANC range is unknown for the waters sampled.
  - It is not necessary to develop incremental points above about pH 5.5 for a Gran determination of the bicarbonate equivalence point.
  - A sufficient number of titration points beyond the equivalence are needed to ensure the accuracy of the calculation.

To obtain the most accurate results for carbonate species, titrate at the field site without delay.

**Quality control (QC)**—Verify your ability to reproduce the alkalinity or ANC determinations by repeating the titration periodically on duplicate or triplicate samples. The frequency and distribution of QC determinations are established by study requirements.

**Rule of thumb**—QC should be no less than every tenth sample. Determination on a filtered sample should be reproducible within  $\pm 5$  percent when titrating a duplicate aliquot from the same batch of sample filtrate.

- ▶ For filtered samples with less than 0.4 meq/L (20 mg/L as  $\text{CaCO}_3$ ), reproducibility should be between 5 and 10 percent.
- ▶ If the alkalinity is about 0.02 meq/L or less, differences between duplicate samples are likely to exceed 10 percent using the standard titration methods because of rounding error alone.

When interferences are absent, titration on an unfiltered sample often results in a determination identical to or within 5 percent of the filtered sample and can be used as the QC check. If filtered and unfiltered values fail the  $\pm 5$ -percent criterion, repeat the titration on a replicate aliquot of filtered sample.

Reproducibility of the ANC determination to within 5 percent on duplicate unfiltered samples can be problematical when the sample has large amounts of particulate matter—extend the quality-assurance criterion to  $\pm 10$  percent.

### Calculation

Use the following equation to calculate the alkalinity or ANC in milliequivalents per liter:

$$(\text{meq/L}) = \frac{mL_a \times N (\text{meq/mL}_{\text{acid}}) \times 1,000 (\text{mL/L})}{mL_{\text{sample}}}$$

where

$mL_a$  = total volume of the standard acid solution used to reach the equivalence point, in milliliters

$N$  = normality of the standard acid solution used, in milliequivalents per milliliter acid and

$mL_{\text{sample}}$  = volume of the sample titrated, in milliliters.

---

## Buret titrator

When using a buret, exercise caution to ensure that the acid does not evaporate or become contaminated by extrinsic matter or moisture. The titrant temperature should be equilibrated to the sample temperature before use. Always empty the buret after each use. Never reuse the titrant solution; dispose of the solution properly.

1. Fill a clean, dry buret with 0.01639*N* sulfuric acid titrant.
  - Use a 10-mL semimicroburet with 0.05-mL graduations and a Teflon™ stopcock for samples with alkalinity or ANC less than 4 meq/L (200 mg/L as CaCO<sub>3</sub>).
  - Use a 25-mL buret with 0.1-mL graduations and a Teflon™ stopcock for samples with alkalinity or ANC of 4 meq/L (200 mg/L as CaCO<sub>3</sub>) or greater and when the sample pH exceeds 8.1.
  - If greater accuracy is desired, use a Gilmont™-type micrometer buret.
2. Make sure that no air bubbles are trapped in buret or buret stopcock. Record the sulfuric acid normality and initial buret reading on field forms.
3. Pipet the selected volume of sample to a clean beaker. Do not pipet by mouth.
  - If a magnetic stirrer is used, place a clean, dry, small stir bar into the beaker before pipetting the sample to the beaker. Do not use a magnetic stirrer if sample conductivity is <100 μS/cm. Place beaker on stirrer.
4. Rinse the pH and temperature sensors with DIW. Blot water droplets adhering to the sensors with lint-free paper (residual DIW will not affect the determination).
5. Insert the sensors in the beaker.
  - Do not let the sensors touch the beaker bottom or wall.
  - Sample depth in the beaker must be sufficient to cover the junction of the reference electrode, the electrode bulb, and the temperature sensor.
6. Measure the initial pH and temperature while gently stirring or after gently swirling the sample.
  - Do not splash the sample onto the beaker wall or out of the beaker.

- Minimize the vortex caused by magnetic stirring.
  - Record on the field form the pH and temperature values, the counter reading (it should read “0000”), the titrant normality, the time, and the sample size.
7. Begin titration. If using a magnetic stirrer, stir the sample slowly and continuously. Measure pH after each addition of titrant, and after the acid and sample are mixed homogeneously. If a magnetic stirrer is not used, swirl to mix the sample and acid after each addition of titrant. Allow 15 to 20 seconds after each addition for equilibration, then record pH.
- **pH  $\geq$  8.1**—To determine the bicarbonate inflection point, add the titrant drop by drop in 0.01-mL increments with a 25-mL buret until pH is less than 8.0.
  - **pH  $<$  8.1**—Titrate rapidly to pH of about 5.0, using several large acid increments (to pH of 5.5 for sample alkalinity or ANC of 0.4 meq/L ( $<20$  mg/L as  $\text{CaCO}_3$ ) or conductivity of sample  $<100$   $\mu\text{S/cm}$ ).
  - **pH  $<$  5.0**—Cautiously add the titrant drop by drop in 0.01-mL increments from pH 5.0 to 4.0 or less (the most sensitive part of the titration curve is usually between pH 4.8 and 4.3). Extend titration to pH 3.5 or less if using the Gran method or for a sample high in organic acid concentration.

TECHNICAL NOTE: 0.01 mL of a standard 0.05-mL drop of titrant tends to remain on the buret tip. To dispense this 0.01-mL titrant drop, quickly rotate the stopcock through 180 degrees (one-half turn) and then rinse the titrant from the buret tip into the filtration beaker with a small quantity of DIW (C.J. Patton, U.S. Geological Survey, written commun., 1995).

### Calculation

Use the following equations to calculate alkalinity or ANC and carbonate species from inflection points with 0.01639N sulfuric acid:

$$\text{Alkalinity or ANC (meq/L)} = \frac{mL_a \times N \text{ (meq/mL}_{\text{acid}}) \times 1000 \text{ (mL/L)}}{mL_{\text{sample}}}$$

$$\text{Alkalinity (mg/L as CaCO}_3) = 1000/mL_s \times (0.8202 \times mL_a) \times CF$$

$$\text{Bicarbonate (mg/L as HCO}_3^-) = 1000/mL_s \times [(mL_a \text{ from ip near pH 8.3 to ip near pH 4.5}) - (mL_a \text{ to ip near pH 8.3})] \times CF$$

$$\text{Carbonate (mg/L as CO}_3^{2-}) = 1000/mL_s \times [(mL_a \text{ at ip near pH 8.3}) \times 0.9835] - (\text{mg/L OH}^- \times 3.527) \times CF$$

$$\text{Hydroxide (mg/L as OH}^-) = 1000/mL_s \times [(mL_a \text{ to ip near pH 8.3}) - (mL_a \text{ from ip near pH 8.3 to ip near pH 4.5})] \times 0.2788 \times CF$$

where:

$mL_{\text{sample}}$  or  $mL_s$  = volume of sample, in milliliters

$mL_a$  = volume of H<sub>2</sub>SO<sub>4</sub> added, in milliliters

*ip* = inflection point

*CF* = correction factor for H<sub>2</sub>SO<sub>4</sub> normality (to be used if normality of H<sub>2</sub>SO<sub>4</sub> is not exactly 0.01639N; *CF* is determined by dividing the actual normality of the H<sub>2</sub>SO<sub>4</sub> by 0.01639)

*N* = normality of acid solution used.

The presence of hydroxide is indicated when the carbonate titrant volume exceeds the bicarbonate titrant volume.

## Digital titrator

It is necessary to be thoroughly familiar with the operation of the digital titrator before field use. A plunger in the digital titrator forces acid in the titrant cartridge into the delivery tube. The plunger is controlled by a main-drive screw, which in turn is controlled by rotation of the delivery knob. The delivery knob controls the volume of titrant delivered through the delivery tube, as indicated by a digital counter.

1. Record the sample volume to be titrated and the titrant normality. Equilibrate titrant temperature to sample temperature.
2. Assemble the digital titrator.
  - a. Depress the plunger-release button and retract the plunger.
  - b. Insert the titrant cartridge into the titrator and twist the cartridge one-quarter turn to lock it into position.
  - c. Carefully depress the plunger-release button and push the plunger forward until it makes contact with the Teflon™ seal inside the cartridge.
3. Remove the vinyl cap from the cartridge (save the cap) and insert the straight end of the delivery tube into the cartridge.
  - Do not push the delivery tube beyond the cartridge tip.
  - Do not alter the delivery tube.
  - Tag the delivery tube to avoid cross contamination. Do not interchange delivery tubes between cartridges with different titrant normality.
4. Purge the titrant through the delivery tube to ensure that no air bubbles or water are in the tube by holding the titrator with the cartridge tip up and turning the delivery knob to force a few drops of titrant through the end of the delivery tube. Rinse tube exterior with DIW and blot off acid or water droplets before inserting it into the sample.
5. Set the digital counter to zero using the counter-reset knob, taking care not to turn the delivery knob.
6. Pipet the selected volume of the sample to a clean beaker. If a magnetic stirrer is used, place a clean, dry, small stir bar into the beaker before pipetting the sample to the beaker. Do not use a magnetic stirrer for sample conductivity <100  $\mu\text{S}/\text{cm}$ .

7. Insert sensors into the beaker.
  - Rinse the pH and temperature sensors with DIW. Blot off water droplets adhering to the sensors.
  - Do not let sensors touch the bottom or wall of the beaker.
  - The amount of sample in the beaker must be sufficient to cover the junction of the reference electrode, the electrode bulb, and the temperature sensor.
8. Measure the initial pH and temperature while gently stirring or after gently swirling the sample.
  - Do not splash sample onto beaker wall or out of the beaker.
  - Minimize the vortex caused by magnetic stirring.
  - Record the pH and temperature values, the counter reading (it should read "0000"), the titrant normality, the time, and the sample size on field forms.
9. Immerse the end of the titrant delivery tube in the sample. To prevent bleeding of the titrant from the delivery tube, keep the aperture of the delivery tube away from the stir bar.
10. **Begin titration.** If using a magnetic stirrer, stir the sample slowly and continuously. Measure pH after each addition of titrant, and after the acid and sample are mixed homogeneously. If a magnetic stirrer is not used, swirl to mix the sample and acid after each addition of the titrant. Allow 15 to 20 seconds after each addition for equilibration, then record pH.
  - **pH  $\geq$  8.1**—Slowly add the titrant in replicate increments no greater than two to three digital counts until pH of the sample is about 8.0, to determine the carbonate inflection point. Record the pH and digital counter reading after each addition of the titrant. Larger increments can be used for samples containing high carbonate concentrations.
  - **pH  $<$  8.1**—Titrate rapidly with several large acid increments to pH of about 5.0 (to pH 5.5 for sample alkalinity at 0.4 meq/L ( $<$ 20 mg/L as  $\text{CaCO}_3$ ) or sample conductivity  $<$ 100  $\mu\text{S}/\text{cm}$ ).
  - **pH  $<$  5.0**—Add titrant cautiously, in increments of one to three digital counts, from pH 5.0 to 4.0 or less. (The most sensitive part of the titration curve is between pH 4.8 and 4.3 for many natural waters.) If using the Gran method, extend the titration to pH 3.5 or less and to pH 3.0 or less for samples high in organic acids and other noncarbonate contributors, or when the alkalinity or ANC range is unknown.

11. After completing the titration, depress the plunger release, retract the plunger, and remove the titrant cartridge. Immediately replace the vinyl cap on the cartridge tip. Discard the delivery tube after each use.
12. Calculate alkalinity/ANC in the field—

$$\text{Alkalinity or ANC (meq/L)} = \frac{mL_a \times N(\text{meq/mL}_{\text{acid}}) \times 1000 \text{ (mL/L)}}{mL_{\text{sample}}}$$

- To determine carbonate alkalinity concentrations and concentrations of contributing species, plot change in pH divided by change in digital counts against digital counts of the titrant or tabulate change in pH divided by change in digital counts.
- Refer to table 6.6–3 for digital titration factors.
- The factors and equations used for the 0.1600*N* or 1.600*N* titrant cartridges are as follows (these calculations are based on the same equations described under the buret titration procedures, except that milliliters of acid used is shown as digital counts for the Hach™ titrator; 800 counts = 1 mL):

$$\text{Alkalinity or ANC (meq/L)} = B(D3)(C_a)/mL_s$$

$$\text{Carbonate (mg/L as CO}_3^{2-}) = A(D1)/mL_s$$

$$\text{Bicarbonate (mg/L as HCO}_3^-) = (B-2A)(D2)/mL_s$$

$$\text{Hydroxide (mg/L as OH}^-) = (A-C)(D4)/mL_s$$

$$\text{Alkalinity (mg/L as CaCO}_3) = B(D3)/mL_s$$

where:

$mL_s$  = volume of the sample, in milliliters

$A$  = digital count from the initial pH to the inflection point near 8.3

$B$  = digital count from the initial pH to the inflection point near 4.5

$C$  = digital count from the inflection point near 8.3 to the inflection point near 4.5

$D$  = digital titration factor

$C_a$  = concentration of acid.

Factors	0.01600 <i>N</i> titrant	0.1600 <i>N</i> titrant	1.600 <i>N</i> titrant
D1	1.20	12.0	120
D2	1.22	12.2	122
D3	1.00	10.0	100
D4	.34	3.4	34



### INFLECTION POINT TITRATION METHOD 6.6.4.B

The IPT method uses the inflection points of plotted values to select equivalence points instead of assuming equivalence points to be at pH 8.3 and 4.5 (as in the fixed endpoint method). Inflection points are points of maximum rate of change in pH per volume of titrant added. **Near equivalence points, rapid pH changes occur with small additions of titrant. For this reason, titration as you near and pass the expected equivalence points must be slow and cautious, using small incremental additions of titrant.** Relative error of the determinations can be within  $\pm 4$  percent if the equivalence point is recognizable within  $\pm 0.3$  pH unit of the true equivalence point.

**Use either the buret or digital titrator.** To determine the inflection point, you can either construct a titration curve by plotting the change in pH divided by the change in titrant volume against the incremental volumes of the titrant added to the sample, or tabulate a titration. Figures 6.6-2 and 6.6-3 represent the titration of a sample that has both a carbonate and a bicarbonate inflection point, whereas tables 6.6-4 and 6.6-5 represent a titration that resulted only in a bicarbonate inflection point.

- ▶ More than one inflection point in close proximity indicates that the true inflection point has been missed. If this occurs, titrate a duplicate sample using smaller acid increments near the inflection point or use a Gran plot.
- ▶ If no clear point or points can be determined easily, interferences from weak organic acids are likely—use the Gran method.

#### EXAMPLES:

**IPT method using the buret system.** Referring to table 6.6-4, pH 4.51 at a titrant volume of 8.95 mL is the point of maximum rate of change of pH per volume titrant. The actual inflection point is before the titrant volume corresponding to the maximum change in pH per unit volume of acid added; therefore, the correct value lies between 8.95 mL and the previous value, 8.90 mL. In this example, the calculated titrant volume 8.93 mL would be the correct inflection point.

**IPT method using the digital titrator.** Referring to table 6.6-5, pH 4.51 is the point of maximum rate of change of pH per volume of titrant. The actual inflection point, however, is between the digital-counter value (454) representing the maximum change in pH per unit volume of acid added and the previous digital-counter value (452). The correct (calculated) digital-counter value for the inflection point would be 453. The error in computing concentration from the digital-counter value 454 instead of the digital-counter value 453 is considered insignificant. Note, however, that the larger the increments used, the greater the significance of the error. Calculation of the correct inflection point is recommended.

Table 6.6-4. Results of typical inflection point titration using a buret  
[mL, milliliter]

pH	Change in pH	Titrant volume (mL)	Change in volume	Change in pH per change in volume
7.28	—	0.00	—	—
7.00	0.28	2.00	2.00	0.140
6.70	.30	4.00	2.00	.150
6.33	.37	6.00	2.00	.185
5.73	.60	8.00	2.00	.300
5.24	.49	8.50	.50	.980
4.93	.31	8.70	.20	1.55
4.81	.12	8.80	.10	1.20
4.72	.09	8.85	.05	1.80
4.62	.10	8.90	.05	2.00
<b>4.51<sup>a</sup></b>	<b>.11</b>	<b>8.95</b>	<b>.05</b>	<b>2.20</b>
4.42	.09	9.00	.05	1.80
4.34	.08	9.05	.05	1.60
4.22	.12	9.12	.07	1.71
3.92	.30	9.22	.20	1.50
3.62	.30	9.32	.20	1.50

<sup>a</sup>pH 4.51 is the point at which the maximum rate of change of pH per volume of titrant occurs.




Figure 6.6-2. Example of an inflection point titration using a buret.

Table 6.6–5. Results of typical inflection point titration using a digital titrator

pH	Change in pH	Counter reading (mL)	Change in counter	Change in pH per change in counter
7.28	—	0	—	—
7.00	0.28	100	100	0.0028
6.70	.30	200	100	.0030
6.33	.37	300	100	.0037
5.73	.60	400	100	.0060
5.24	.49	435	35	.0140
4.93	.31	445	10	.0310
4.81	.12	448	3	.0400
4.72	.09	450	2	.0450
4.62	.10	452	2	.0500
<b>4.51<sup>a</sup></b>	<b>.11</b>	<b>454</b>	<b>2</b>	<b>.0550</b>
4.42	.09	456	2	.0450
4.34	.08	458	2	.0400
4.22	.12	461	3	.0400
3.92	.30	471	10	.0300
3.62	.30	481	10	.0300

<sup>a</sup>pH 4.51 is the point of maximum rate of change of pH per digital count.




Figure 6.6–3. Example of inflection point titration using a digital titrator.

#### 6.6.4.C GRAN FUNCTION PLOT METHOD

Gran function plots commonly are used to determine alkalinity and ANC in sea water, low ionic-strength water, water with low carbonate concentrations, and water with measurable concentrations of organic compounds. The Gran function plot method also is used for calculations of the base neutralizing capacity in waters of low ionic strength such as atmospheric deposition.

The Gran function plot method uses functions that linearize titration curves, making it possible to determine alkalinity or ANC with a few points rather than relying on the inflection point from an entire titration curve (Baedecker and Cozzarelli, 1992). Four Gran functions can be calculated over the entire titration curve. The  $F_1$  function, described below, is the most commonly applicable Gran function.

- ▶ The  $F_1$  function requires titration data beyond the equivalence point, usually to pH between 3.0 and 3.5. For systems with measurable concentrations of organic acids, titrate to pH 2.5 (Baedecker and Cozzarelli, 1992).
- ▶ The buret titration system is recommended to achieve the accuracy desired when using the Gran method: substitute the equivalence points determined by the Gran function plot into the equations given in the discussion of the buret system. The digital titration system also is used under some circumstances (illustrated on table 6.6–6 and fig. 6.6–4): refer to table 6.6–3 to calculate concentrations.

**To construct a Gran function plot:**

1. Titrate incrementally to about pH 3.5. It is not necessary to develop incremental points above about pH 5.5 for a Gran determination of the bicarbonate equivalence point.
2. Plot  $(V_o + V_t) 10^{-\text{pH}}$  against the titrant volume (fig. 6.6–4 and table 6.6–6) using the  $F_1$  function for the bicarbonate equivalence point,

where:

$V_o$  = volume of the sample

$V_t$  = volume of the titrant added

$V_s$  = volume of the titrant needed to reach the bicarbonate equivalence point

$V_w$  = volume of the titrant needed to reach the carbonate equivalence point.

3. Extrapolate the straight line in the region beyond the equivalence point to  $(V_o + V_t) 10^{-\text{pH}} = 0$  or  $V_t = V_s$ .
- ▶ **Carbonate equivalence point.** Plot  $(V_t - V_s) 10^{-\text{pH}}$  against the titrant volume (in mL or digital counts). To locate the equivalence point, extrapolate the straight line in the region beyond the equivalence point to  $(V_t - V_s) 10^{-\text{pH}} = 0$  or  $V_t = V_w$ .
  - ▶ **Hydroxide equivalence point.** Plot  $(V_o - 2V_w + V_t) 10^{-\text{pH}}$  against the titrant volume (in mL or digital counts). To locate the equivalence point, extrapolate the straight line in the region beyond the equivalence point to  $(V_o + V_t) 10^{-\text{pH}} = 0$  or  $V_t = V_w$ .

TECHNICAL NOTE: During alkalinity titration (carbonate system), the hydrogen ions added convert carbonate to bicarbonate and then bicarbonate to carbonic acid. The titration continues until no more species are reacting. When this process is complete, additional hydrogen ions will be in excess in the solution. The Gran function plot identifies the point at which all alkalinity has been titrated and hydrogen ions begin to be in excess. Beyond the equivalence point, the shape of the curve is determined by hydrogen ions in excess of all hydronium ion acceptors in the sample. The curvature results from the logarithmic relation between pH and hydrogen-ion activity.

**EXAMPLE:****Gran function plot method using the digital titrator.**

Refer to figure 6.6-4 and table 6.6-6. In the region beyond the equivalence point in figure 6.6-4, a straight line results. Extrapolation of this straight line to  $(V_o=V_i)10^{-\text{pH}}=0$  or  $V_i=V_s$  locates the equivalence point. The extrapolated straight line intercept at  $(V_o+V_i)10^{-\text{pH}}=0$  on figure 6.6-4 is 0.566 mL of titrant added and corresponds to an equivalence point at approximately pH of 4.58.




Figure 6.6-4. Example of Gran function plot of a bicarbonate titration using a digital titrator.

Table 6.6–6. Example of information needed for a Gran function plot

[Data shown in columns 1, 2, 6, and 7 are results of an incremental titration using a digital titrator—remaining columns show calculations used for the Gran function plot (fig. 6.6–4); V, volume of sample; v, volume of titrant; mL, milliliters; cts, digital counts; e-, base 10 logarithm to the negative exponent]

pH	Change in pH	$10^{-\text{pH}}$	(V+v) in mL	$\frac{(V+v)}{10^{-\text{pH}}}$ in mL	Counter reading cts/mL	Change in counter	Change in pH/change in counter
7.28	0.00	5.248e-8	50.000	2.624e-6	0/0.00	0	0
7.00	.28	1.000e-7	50.125	5.012e-6	100/.125	100	.0028
6.70	.30	1.995e-7	50.250	1.003e-5	200/.250	100	.0030
6.33	.37	4.677e-7	50.375	2.356e-5	300/.375	100	.0037
5.73	.60	1.862e-6	50.500	9.404e-5	400/.500	100	.0060
5.24	.49	5.754e-6	50.544	2.908e-4	435/.544	35	.0140
4.93	.31	1.175e-5	50.556	5.940e-4	445/.556	10	.0310
4.81	.12	1.549e-5	50.560	7.831e-4	448/.560	3	.0400
4.72	.09	1.906e-5	50.563	9.635e-4	450/.563	2	.0450
4.62	.10	2.399e-5	50.565	1.213e-3	452/.565	2	.0500
4.51	.11	3.090e-5	50.568	1.563e-3	454/.568	2	.0550
4.42	.09	3.802e-5	50.570	1.923e-3	456/.570	2	.0450
4.34	.08	4.571e-5	50.573	2.312e-3	458/.573	2	.0400
4.22	.12	6.026e-5	50.576	3.048e-3	461/.575	3	.0400
3.92	.30	1.202e-4	50.589	6.082e-3	471/.586	10	.0300
3.75	.17	2.399e-4	50.601	8.998e-2	477/.596	6	.0280
3.59	.16	2.399e-4	50.610	1.301e-2	487/.610	10	.0160

---

## 6.6.5 TROUBLESHOOTING

The major difficulties with equipment for alkalinity or ANC are the same as for pH measurement—refer to NFM 6.4. Particulate materials, including algae or other biota, can interfere with the stability and reproducibility of pH readings. Such difficulties normally are eliminated by filtering the sample.

When the sample has low ionic strength, or when dissolved organic compounds or noncarbonate inorganic species are present that can interfere with the titration (note color, odor, or previous chemical analysis), follow the procedures for a Gran function determination.



## REPORTING 6.6.6

Report alkalinity and ANC to three significant figures. Only the value from filtered samples is to be published as alkalinity. Titration values from unfiltered samples are to be entered and published as ANC (the NWIS parameter code dictionary uses the term “alkalinity, unfiltered” instead of ANC).

Alkalinity and ANC should be reported in milliequivalents per liter (or microequivalents per liter), if possible. If this option is not available in the data base, calculate alkalinity and ANC in milligrams per liter, assigning all neutralizing capacity to the carbonate system.

If calculating ANC, alkalinity, bicarbonate, or carbonate in milligram units, then report: less than 1,000 mg/L, to whole numbers; 1,000 mg/L and above, to three significant figures. Carbonate alkalinity usually is reported in the data base in milligrams per liter as calcium carbonate.

Conversion factors listed below are taken from Hem (1985).

<b>Multiply</b>	<b>By</b>	<b>To obtain</b>
Alkalinity (mg/L as $\text{HCO}_3^-$ )	0.8202	Alkalinity (mg/L as $\text{CaCO}_3$ )
Alkalinity (mg/L as $\text{HCO}_3^-$ )	0.4917	Alkalinity (mg/L as $\text{CO}_3^{2-}$ )
Alkalinity (mg/L as $\text{CaCO}_3$ )	0.08332	Alkalinity (meq/L as $\text{CaCO}_3$ )
Bicarbonate (mg/L as $\text{HCO}_3^-$ )	0.01639	Bicarbonate (meq/L as $\text{HCO}_3^-$ )
Bicarbonate (mg/L as $\text{HCO}_3^-$ )	16.388	Bicarbonate ( $\mu\text{eq/L}$ as $\text{HCO}_3^-$ )
Carbonate (mg/L as $\text{CO}_3^{2-}$ )	0.03333	Carbonate (meq/L as $\text{CO}_3^{2-}$ )
Hydroxide (mg/L as $\text{OH}^-$ )	0.05880	Hydroxide (meq/L as $\text{OH}^-$ )

Report the average value for duplicate samples or the median when more than two replicate samples are used for quality control and the value falls within the appropriate quality-assurance criterion for variability.

Use the correct parameter code to indicate (1) the method of titration or calculation and (2) a filtered or unfiltered sample.