

## **Sodium and Potassium Analysis by Flow Cytometry**

### **Preparation of the dyes:**

PBFI-AM (Molecular Probes) 50 ug is dissolved in 8.5 ul of DMSO and 8.5 ul of Pluronic F-127. Working concentration 2.5 mM

SBFI-AM (Molecular Probes) 50 ug is dissolved in 8.9 ul of DMSO and 8.9 ul of Pluronic F-127. Working concentration 2.5 mM

Corona Green (Molecular Probes) 50 ug is dissolved in 15.25 ul of DMSO and 15.25 ul of Pluronic F-127. Working concentration of 2.5 mM

### **Loading of the dyes:**

1. Add 2 ul of PBFI-AM or SBFI AM or Corona Green per 1 ml of cells (approximately  $5 \times 10^5$  cells per ml) for a final concentration of 5  $\mu$ M and vortex gently.
2. Incubate cells are at 37°C for 60 minutes to load the dye.
3. Immediately prior of flow cytometric examination, add 1 ul of propidium iodide (PI; 10 mg/ml stock) for a final concentration of 10  $\mu$ g/ml.

### **Flow Cytometric Analysis:**

Cells are examined by exciting the PBFI-AM or SBFI-AM dyes with a UV laser. PI and Corona Green are excited using a 488 nm laser.

#### **A standard optical set-up is used on FACSVantage SE:**

PBFI-AM and SBFI-AM are detected using a 424/44 nm filter in front of the UV FL-4 PMT.

PI is detected using a 630/22 nm filter in front of the FL-3 PMT.

Corona Green is detected using a 530/30 in front of the FL-1 PMT.

#### **A standard optical set-up used on the LSR II:**

PBFI-AM and SBFI-AM are detected using a 440/40 nm filter in front of the UV-B PMT.

PI is detected using a 695/40 nm filter in front of the 488-B PMT.

Corona Green is detected using a 530/30 in front of the 488-D PMT.

Cells that are PI positive indicating a loss of membrane integrity are eliminated from further ion analysis.