

DEPARTMENT OF HEALTH AND HUMAN SERVICES PUBLIC HEALTH SERVICE FOOD AND DRUG ADMINISTRATION <b>MILK LABORATORY EVALUATION FORM</b>	LABORATORY	
	LOCATION	LAB #
	DATE	X = DEVIATION      U = UNDETERMINED O = NOT USED      NA = NOT APPLICABLE

**ELECTRONIC SOMATIC CELL COUNT**  
**Foss 250/300/360/400**  
[Unless otherwise stated all tolerances are ±5%]

- 1. Laboratory Requirements (see CP, item 33 & 34)** .....
- a. Unpreserved samples may be run up to 72 hours after initial collection .....
  - b. Samples may be run up to 7 days after initial collection if preserved with 0.02% 2-bromo-2-nitropropane-1, 3-diol (Bronopol™) or 0.05% potassium dichromate (K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub>) .....
  - c. Comparative test with DMSCC .....
  - 1. Performed by each analyst performing ESCC test .....
  - 2. Test 4 samples (100K - 200K; 300K - 500K; 600K - 800K and 900K - 1.2M) in triplicate for both DMSCC (three separate smears each) and ESCC (three separate sub-samples each, do not read same sample three times) .....
  - 3. Results must be shown to be acceptable prior to official testing by analyst performing comparison, i.e. analyst is not certified until found acceptable (**co-requisite for certification**) .....
  - 4. Copy of comparison and results in QC record (or easily accessible file in laboratory) .....
  - d. Analysts certified for DMSCC .....

**APPARATUS**

- 2. See Cultural Procedures, items 1 - 5** .....
- 3. Automated Electronic Somatic Cell Counters** .....
- a. **Fossomatic 250** .....
  - b. **Fossomatic 300** .....
  - c. **Fossomatic 360** .....
  - d. **Fossomatic 400** .....
- 4. Water Bath** .....
- a. Circulating and thermostatically controlled to 37 - 42C .....

**REAGENTS**

- 5. Stock Dye/Buffer Solution** .....
- a. Dissolve 2.5g (or number of tablets specified by manufacturer) ethidium bromide (C<sub>21</sub>H<sub>20</sub>BrN<sub>3</sub>) in 1 liter MS water (caution TOXIC, use gloves when handling and do not breath dust), heat to 40 - 60C and mix to dissolve .....
  - b. Add 400g tripotassium citrate monohydrate (C<sub>6</sub>H<sub>5</sub>O<sub>7</sub>K<sub>3</sub>·H<sub>2</sub>O), 14.5g citric acid monohydrate (C<sub>6</sub>H<sub>8</sub>O<sub>7</sub>·H<sub>2</sub>O), and 4 liters MS water, heat to 40 - 60C and mix to dissolve .....
  - c. Add dye and buffer solutions together and mix .....
  - d. Add 50 mL neutral detergent, Triton X-100 to mixture and stir until dissolved .....
  - e. Store refrigerated (0 - 4.4C) in airtight, light-proof container for no longer than 90 days .....
  - f. Date prep. \_\_\_\_\_ Exp. Date \_\_\_\_\_
- 6. Stock Detergent Solution** .....
- a. Dissolve 10 mL neutral detergent, Triton X-100 in 1 liter of MS water and heat 40 - 60C to complete solution .....

- b. Store refrigerated (0 - 4.4C) in airtight, container for no longer than 30 days .....
- c. Date prep. \_\_\_\_\_ Exp. Date \_\_\_\_\_

- 7. Ammonium Hydroxide (NH<sub>4</sub>OH) Solution, reagent Grade, 25%** ....
- 8. All stock dye/buffer and detergent solutions labeled with date prepared and expiration date** .....

**WORKING SOLUTIONS**

- 9. Dye/Buffer Solution** .....
- a. Dilute 1 L dye/buffer stock solution (item 5) with 9 L MS water .....
  - b. Protect from light and use within 21 days .....
  - c. Date prep. \_\_\_\_\_ Exp. Date \_\_\_\_\_
- 10. Rinsing Solution (use within 7 days)** .....
- a. Add 10 mL of stock neutral detergent stock solution (item 6) and 25 mL of ammonium hydroxide solution (item 7) and suspend to 10 L with MS water .....
  - b. Date prep. \_\_\_\_\_ Exp. Date \_\_\_\_\_
- 11. Optionally, use manufacturer's reagent kits and instructions specific for each instrument** .....
- 12. All working dye/buffer and rinsing solutions labeled with date prepared and expiration date** .....

**START UP**

- 13. Cell Counter** .....
- a. Check that the amount of dye/buffer solution (item 9) and rinsing (cleaning) solution (item 10) in instrument supply containers is of sufficient volume for the number of samples to be run .....
  - b. Solutions not used beyond expiration date(s) .....
  - c. Turn on power and cycle at least six times .....
  - d. Perform a zero check before starting any measurements, within acceptable limits, single counts up to 5 and mean up to 3 .....
  - e. **IF ANY ABOVE PARAMETERS ARE WRONG, CORRECT BEFORE PROCEEDING** .....
  - f. Records maintained on all parameters each time instrument is used .....
- 14. Milk Standards** .....
- a. Commercially prepared: \_\_\_\_\_  
Lot# \_\_\_\_\_ Date Rcd. \_\_\_\_\_
  - 1. Four standards in ranges 100K - 200K, 300K - 500K, 600K - 800K and 900K - 1.2M .....
  - 2. Do DMSCC in triplicate on each standard in set and average counts, records maintained .....
  - 3. DMSCC check performed in rotation by all certified analysts .....
  - 4. Standards used within one week .....

LABORATORY	LAB #	LOCATION	DATE
------------	-------	----------	------

**ELECTRONIC SOMATIC CELL COUNT  
Foss 250/300/360/400**

**[Unless otherwise stated all tolerances are ±5%]**

- b. Certified provider: \_\_\_\_\_  
 Lot# \_\_\_\_\_ Exp. Date \_\_\_\_\_  
 Date Rcd. \_\_\_\_\_
- 1. Four standards in ranges 100K - 200K, 300K - 500K, 600K - 800K and 900K - 1.2M \_\_\_\_\_
- 2. Maintain copies of all provided DMSCC values \_\_\_\_\_
- 3. Measure and maintain records of temperature (0 - 7.2C) of standards as received \_\_\_\_\_
- 4. Maintain copies of all correspondence regarding problems \_\_\_\_\_
- c. Laboratory prepared (weekly) \_\_\_\_\_
- 1. Prepare from raw milk > 18 hours old preserved with 0.05% potassium dichromate (K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub>) \_\_\_\_\_
- 2. Or, preserved with 0.02% 2-bromo-2-nitropropane-1, 3-diol (Bronopol™) \_\_\_\_\_
- 3. Standards cannot be preserved with formalin \_\_\_\_\_
- 4. Prepare 4 standards in ranges 100K - 200K, 300K - 500K, 600K - 800K and 900K - 1.2M, used within one week \_\_\_\_\_  
 Date prep. \_\_\_\_\_ Exp. Date \_\_\_\_\_
- 5. Do DMSCC in triplicate on each standard and average counts, records maintained \_\_\_\_\_
- 6. DMSCC check performed in rotation by all certified analysts \_\_\_\_\_
- d. Hourly Control Sample (instrument drift check) \_\_\_\_\_
- 1. Use one of the standards (items 14a or b) in the 500 - 800k range, run in triplicate and determine average \_\_\_\_\_
- 2. Optionally, prepare sufficient control/sample 500 - 800k range, run in triplicate and determine average \_\_\_\_\_

**PROCEDURE**

- 15. Testing Standards (each time instrument used)** \_\_\_\_\_
- a. Heat standards to 37 - 42C (using a temperature control) and read within 30 minutes of reaching temperature, used once and then discarded, i.e. do not re-use \_\_\_\_\_
- b. Mix by inverting at least 2x, place in rack and put onto automatic track, run within 10 minutes \_\_\_\_\_
- c. Run the standards six times and average the counts for each level, records maintained \_\_\_\_\_
- d. Each standard's average must be within 10% of the DMSCC (item 14) for that level, except within 15% for 100K - 200K standard, records maintained \_\_\_\_\_

- e. Repeatability — a standard in the 300K to 800K range must have a coefficient of variation (C<sub>v</sub>) of 5% or less on 10 replicates (**Refer to Operating Manual**), records maintained \_\_\_\_\_
- f. **THESE PARAMETERS MUST BE ACHIEVED BEFORE PROCEEDING** \_\_\_\_\_
- 16. Testing Samples** \_\_\_\_\_
- a. Heat samples to 37 - 42C (using a temperature control) and read within 30 minutes of reaching temperature; samples must not be reused and must be discarded after use \_\_\_\_\_
- b. Mix by inverting at least 2x, place in rack and put onto automatic track, run within 10 minutes of reaching the testing temperature \_\_\_\_\_
- 17. With continuous operation:** \_\_\_\_\_
- a. Run a standard or optionally a control/sample (item 14d) in the 500K to 800K range hourly, must be within 5% of the original established instrument average value (optionally, within 10% of original DMSCC average) \_\_\_\_\_
- b. Run control 6x \_\_\_\_\_
- c. Run zero control (as in item 13d) \_\_\_\_\_
- d. Maintain records \_\_\_\_\_
- 18. Routine maintenance** \_\_\_\_\_
- a. Perform as described in operating manual \_\_\_\_\_
- b. Maintain records \_\_\_\_\_

**REPORTS**

- 19. Computing and Reporting Counts** \_\_\_\_\_
- a. Count obtained x 1000 is the cell count/mL milk \_\_\_\_\_
- b. In reporting electronic somatic cell counts (ESCC/mL), record only first two left hand digits, raising second digit to next higher number when third digit is 6 or more \_\_\_\_\_
- c. Report the two left hand digits (rounded) \_\_\_\_\_
- 1. If the third digit is 5 the second digit is rounded by the following rule \_\_\_\_\_
- a. When the second digit is odd round up, raise the second digit by 1 (odd up, 235 = 240) \_\_\_\_\_
- b. When the second digit is even round down, delete the and report the second digit as is (even down, 225 to 220) \_\_\_\_\_
- d. If count on instrument is < 100 report as < 100,000 ESCC/mL \_\_\_\_\_