

UNITED STATES OF AMERICA
DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION

+ + +

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
MEDICAL DEVICES ADVISORY COMMITTEE

+ + +

GASTROENTEROLOGY AND UROLOGY DEVICES PANEL

+ + +

June 25, 2008
8:00 a.m.

Hilton Washington DC North
620 Perry Parkway
Gaithersburg, MD 20877

PANEL MEMBERS:

MARK TALAMINI, M.D.	Chair, Voting Member
MANOOP BHUTANI, M.D.	Voting Member
JASON CONNOR, Ph.D.	Voting Member
PHILLIP DAHM, M.D.	Voting Member
SUSAN J. KALOTA, M.D.	Voting Member
TERRY N. LAYTON, Ph.D.	Industry Representative
FRANCINE STOKES, J.D.	Consumer Representative
CRAIG DONATUCCI, M.D.	Temporary Voting Member
MARGUERITE LIPPERT, M.D.	Temporary Voting Member
ROBERT MARCOVICH, M.D.	Temporary Voting Member
BRUCE REDMAN, D.O.	Temporary Voting Member

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

FDA REPRESENTATIVES:

NANCY BROGDON, M.D.
Director, Division of Reproductive,
Abdominal and Radiological Devices

JEFFREY COOPER, D.V.M.
Veterinary Medical Officer,
Gastroenterology Renal Devices,
Executive Secretary to Gastroenterology
and Urology Devices Panel

MEGAN MICKAL
Biomedical Engineer, Gastroenterology and Renal
Devices, Executive Secretary in training to
Gastroenterology and Urology Devices Panel

DANICA MARINAC-DABIC, M.D., Ph.D.
Chief of Epidemiology,
CDRH, Office of Surveillance and Biometrics

NILSA LOYO-BERRIOS, Ph.D.
Epidemiologist,
Office of Surveillance and Biometrics

DONNA-BEA TILLMAN, Ph.D.
Director, Office of Device Evaluation

FDA PRESENTERS:

JOHN BAXLEY, M.S.
HECTOR HERRERA, M.D., MPH
ROBERT KANE, M.D., FACP
XUEFENG LI, Ph.D.
SHAOKUI WEI, M.D., MPH

SPONSOR REPRESENTATIVES:

YAGEL KOREN, M.D.
MICHAEL O'DONNELL, M.D.
AHAVA STEIN
FRED WITJES, M.D.
H. BARTON GROSSMAN, M.D.
LISA DEUTSCH, Ph.D.

OTHER PARTICIPANTS:

ROBERT LIPMAN

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

INDEX

	PAGE
CALL TO ORDER - Mark Talamini, M.D.	5
FDA Conflict of Interest Disclosure Statement - Jeffrey Cooper, D.V.M.	6
Appointment to Temporary Voting Status - Jeffrey Cooper, D.V.M.	9
General Announcements - Jeffrey Cooper, D.V.M.	11
Introductions	12
Office of Surveillance and Biometrics Post-Market Studies Update - Danica Marinac-Dabic, M.D., Ph.D.	15
OPEN PUBLIC HEARING	28
OPEN COMMITTEE DISCUSSION Sponsor Presentation	
Introduction - Yagel Koren, M.D.	34
Treatment Options - Michael O'Donnell, M.D.	37
Device Description and Preclinical Studies - Ahava Stein	51
Overview of Clinical Studies - Fred Witjes, M.D.	55
Overall Summary - H. Barton Grossman, M.D.	78
Post-Approval Study - Michael O'Donnell, M.D.	84
Questions and Answer Session	89
FDA Presentation	
Regulatory History/Preclinical Studies - John Baxley, M.S.	100

INDEX

	PAGE
FDA Presentation (Continued)	
Clinical Overview of the Disease - Hector Herrera, M.D., MPH	109
Clinical Considerations - Robert Kane, M.D.	113
Statistical Considerations - Xuefeng Li, Ph.D.	130
Post Market Review - Shaokui Wei, M.D., MPH	136
Questions on the Presentations and Discussion	144
Panel Discussion and FDA Questions	
General Discussion	163
Reading of Questions for the Panel and Discussion	199
Question 1	199
Question 2	205
Question 3	209
Question 4	241
Open Public Comments (None)	261
Final Comments	
FDA Comments (None)	262
Sponsor Comments	262
Panel Recommendation Options for Pre-market Approval Applications - Jeffrey Cooper, D.V.M.	263
Panel Deliberations and Vote	265
Adjournment	300

M E E T I N G

(8:00 a.m.)

1
2
3 DR. TALAMINI: Good morning, everybody. I
4 would like to call this meeting of the
5 Gastroenterology and Urology Devices Panel to order.
6 I'm Dr. Mark Talamini, Chairperson of this Panel.
7 I'm a gastrointestinal surgeon, the Chairman of the
8 Department of Surgery at University of California,
9 San Diego.

10 If you haven't already done so, please sign
11 the attendance sheets that are on the tables by the
12 doors. If you wish to address this Panel during one
13 of the open sessions, please provide your name to
14 Ms. AnnMarie Williams at the registration table.

15 If you are presenting in any of the open
16 public sessions today and have not previously
17 provided an electronic copy of your presentation to
18 FDA, please arrange to do so with Ms. Tobey Lowe.
19 Tobey, if you could stand. There's Tobey. Thank
20 you.

21 I note for the record that the voting
22 members present constitute a quorum as required by 21
23 C.F.R. Part 14. I would also like to add that the
24 Panel participating in the meeting today has received
25 training in FDA device law and regulations.

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 No one from the public or the press is
2 allowed into the Panel area up here at anytime during
3 the break or during the conduct of this meeting, and
4 I would remind everybody to please silence your cell
5 phones for the smooth conduct of the meeting.

6 Dr. Cooper, the Executive Secretary for the
7 Gastroenterology and Urology Devices Panel, will make
8 some introductory remarks. Dr. Cooper.

9 DR. COOPER: Thank you. I'm now going to
10 read the Conflict of Interest Statement, the FDA
11 Conflict of Interest Disclosure Statement, with a
12 particular matter involving specific parties.

13 The date of the meeting is June 25, 2008.

14 The Food and Drug Administration is
15 convening today's meeting of the Gastroenterology and
16 Urology Devices Panel of the Medical Devices Advisory
17 Committee under the authority of the Federal Advisory
18 Committee Act, FACA, of 1972. With the exception of
19 the industry representative, all members and
20 consultants of the Panel are special government
21 employees or regular federal employees from other
22 agencies and are subject to federal conflict of
23 interest laws and regulations.

24 The following information on the status of
25 this Panel's compliance with the federal ethics and

1 conflict of interest law is covered by, but not
2 limited to, those found at 18 U.S.C. 208 and 712 of
3 the federal Food, Drug and Cosmetic Act, the FD&C
4 Act, are being provided to participants in today's
5 meeting and to the public. FDA has determined that
6 members and consultants of this Panel are in
7 compliance with federal ethics and conflict of
8 interest laws.

9 Under 18 U.S.C. 208, Congress has
10 authorized FDA to grant waivers to special government
11 employees who have potential financial conflicts when
12 it is determined that the Agency's need for a
13 particular individual's services outweighs his or her
14 potential financial conflict of interest.

15 Under 712 of the FD&C Act, Congress has
16 authorized FDA to grant waivers to special government
17 employees and regular government employees with
18 potential financial conflicts when necessary to
19 afford the Committee essential expertise.

20 Related to the discussions of today's
21 meetings, members and consultants of this Panel or
22 special government employees have been screened for
23 potential financial conflicts of interest of their
24 own as well as those imputed to them, including those
25 of their spouses or minor children and, for purposes

1 of 18 U.S.C. 208, their employers. These interests
2 may include investments, consulting, expert witness
3 testimony, contracts, grants, CRADAs, teaching,
4 speaking, writing, patents, royalties and primary
5 employment.

6 Today's agenda involves the discussion of a
7 pre-market approval application, a PMA, for the
8 Synergo SB-TS 101.1 Device and Mitomycin C, sponsored
9 by Medical Enterprises, Ltd. This drug/device
10 combination product is designed to prevent the
11 recurrence of bladder cancer. Synergo SB-TS 101.1
12 Device with Mitomycin C is indicated for use for
13 prophylactic treatment of recurrence in patients
14 following endoscopic removal of Ta to T1 and G1 to 3
15 superficial transitional cell carcinoma of the
16 bladder, STCCB. Synergo and Mitomycin C treatment is
17 clinically indicated for STCCB patients of
18 intermediate and high risk.

19 This is a particulars matter meeting which
20 specific matters related to the PMA will be
21 discussed.

22 Based on the agenda for today's meeting and
23 all financial interests reported by the Panel members
24 and consultants, no conflict of interest waivers have
25 been issued in accordance with 18 U.S.C. 208 and 712

1 of the FD&C Act.

2 A copy of this statement will be available
3 for review at the registration table during this
4 meeting and will be included as part of the official
5 transcript.

6 Terry Layton, Ph.D., is serving as the
7 Industry Representative acting on behalf of all
8 related industry and is self-employed by Laytech,
9 Incorporated.

10 We would like to remind members and
11 consultants that if the discussions involve any other
12 products or firms not already on the agenda for which
13 a FDA participant has the personal or imputed
14 financial interest, their participants need to
15 exclude themselves from such involvement and their
16 exclusion will be noted for the record.

17 FDA encourages all other participants to
18 advise the Panel of any financial relationships that
19 they may have with any firms at issue. Thank you.

20 Now, I will read the Appointment to
21 Temporary Voting Status.

22 Pursuant to the authority granted under the
23 Medical Devices Advisory Committee Charter, dated
24 October 27, 1990, and amended August 18, 2006, I
25 appoint the following as voting members to the

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 Gastroenterology and Urology Devices Panel for the
2 duration of this meeting on June 25, 2008: Craig
3 Donatucci, M.D., Marguerite Lippert, M.D. and Robert
4 Marcovich, M.D.

5 For the record, these people are special
6 government employees and are consultants to this
7 Panel or another panel under the Medical Devices
8 Advisory Committee. They have undergone the
9 customary conflict of interest review and have
10 reviewed the material to be considered at this
11 meeting.

12 This was signed by Daniel G. Schultz, M.D.,
13 Director, Center for Devices and Radiological Health
14 and dated May 28, 2008.

15 I'll also read a second appointment to
16 temporary voting status, and that is pursuant to the
17 authority granted under the Medical Devices Advisory
18 Committee Charter of the Center for Devices and
19 Radiological Health, dated October 27, 1990, and
20 amended August 18, 2006, I appoint Bruce G. Redman,
21 M.D., as a temporary voting member of the
22 Gastroenterology and Urology Devices Panel for the
23 duration of the meeting on June 25, 2008.

24 For the record, Dr. Redman serves as a
25 consultant to the Oncologic Drugs Advisory Committee

1 of the Center for Drug Evaluation and Research. He
2 is a special government employee who has undergone
3 the customary conflict of interest review and has
4 reviewed the material to be considered at this
5 meeting.

6 This was signed by Randall W. Lutter,
7 Ph.D., Deputy Commissioner for Policy, on May 13,
8 2008.

9 I would also like to note the absence of
10 our Patient Representative, Col. James D. Schultz.
11 He was enthusiastic about attending the meeting to
12 offer his viewpoints as a patient. Unfortunately,
13 his family informed us that he passed away on May
14 22nd. His son, Jim Schultz, Jr., told us that, "I
15 know that he appreciated the opportunity to return to
16 the medical community some measure of his
17 appreciation for the wonderful care he received over
18 the years."

19 Before I turn the meeting back over to
20 Dr. Talamini, there are a few general announcements.

21 Transcripts of today's meeting will be
22 available from the Free State Court Reporting. Their
23 contact information is available as a handout at the
24 registration table outside.

25 Information on purchasing videos of today's

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 meeting can be found on the table outside the meeting
2 room also.

3 Presenters to the Panel who have not
4 already done so should provide FDA with a hard copy
5 of their remarks including any overheads.

6 I'd like to remind everyone that members of
7 the public and press are not permitted around the
8 Panel area, beyond the speaker's podium.

9 The press contact for today's meeting is
10 Peper Long.

11 And I request that the reporters wait to
12 speak to FDA officials until after the Panel meeting.

13 Thank you. Now, Dr. Talamini.

14 DR. TALAMINI: Again, good morning,
15 everyone. At this meeting, the Panel will be making
16 a recommendation to the Food and Drug Administration,
17 the FDA, on the pre-market approval application, PMA,
18 P010045 for the Synergo SB-TS 101.1 Device and
19 Mitomycin C for Medical Enterprises, Ltd.

20 Before we begin, I would like to ask our
21 Panel members and FDA staff seated at this table to
22 introduce themselves. Please state your name, your
23 area of expertise, your position and your
24 affiliation, and I would remind, a technical detail
25 for those of us at the table, when you push your

1 button, we can only have four buttons pushed at a
2 time. So when you're done speaking, please push it
3 and turn it off again, and if we could begin over to
4 the right with Ms. Brogdon.

5 MS. BROGDON: Good morning. I'm Nancy
6 Brogdon. I'm not a member of the Panel. I'm the
7 Director of FDA's Division of Reproductive, Abdominal
8 and Radiological Devices.

9 DR. MARCOVICH: Good morning. I'm Robert
10 Marcovich. I'm a urologist at the University of
11 Texas Health Science Center in San Antonio.

12 DR. DONATUCCI: Good morning. Craig
13 Donatucci. I'm a urologist at Duke University in
14 North Carolina.

15 DR. LIPPERT: Good morning. I'm Marguerite
16 Lippert. I'm a urologist at the University of
17 Virginia in Charlottesville, Virginia.

18 DR. BHUTANI: Good morning. I'm Manoop
19 Bhutani. I'm a gastroenterologist at MD Anderson
20 Cancer Center in Houston.

21 DR. CONNOR: Good morning. I'm Jason
22 Connor. I'm a biostatistician. Basically I design
23 clinical trials typically in the regulatory
24 environment, and I work for Berry Consulting in
25 Orlando, Florida.

1 MS. MICKAL: Good morning. I'm Megan
2 Mickal. I'm a biomedical engineer in the
3 Gastroenterology and Renal Devices Branch, and I'm
4 the Executive Secretary in training.

5 DR. COOPER: Good morning. I'm Jeff
6 Cooper, veterinary medical officer in the Gastro
7 Renal Devices Branch of the FDA and also the
8 Executive Secretary for the Gastroenterology and
9 Urology Devices Panel.

10 DR. TALAMINI: Again, my name is Mark
11 Talamini, Panel Chair, gastrointestinal surgeon
12 employed at University of California, San Diego, as
13 the Chair of the Department of Surgery.

14 DR. DAHM: Good morning. My name is
15 Phillip Dahm. I'm a urologist at the University of
16 Florida in Gainesville.

17 DR. KALOTA: Good morning. I'm Susan
18 Kalota, private practice, urology, in Arizona.

19 DR. REDMAN: Good morning. Bruce Redman,
20 Medical Oncologist, University of Michigan,
21 Comprehensive Cancer Center.

22 MS. STOKES: Good morning. I'm Francine
23 Stokes, Esquire. I'm an assistant to the President
24 of Morgan State University for Government Relations.

25 DR. LAYTON: Good morning. I'm Terry

1 Layton, a biomedical engineer. I'm the Industry
2 Panel member here today and also from Laytech located
3 in Chicagoland area.

4 DR. TALAMINI: Thank you, everyone. Next
5 Danica Marinac-Dabic, I hope I didn't hurt that too
6 badly, from the Office of Surveillance and
7 Biometrics, would like to provide the Panel with a
8 Post-Market Studies Update. Dr. Marinac-Dabic.

9 DR. MARINAC-DABIC: Thank you. Good
10 morning, ladies and gentlemen, Dr. Talamini,
11 Dr. Brogdon, distinguished members of the Panel.

12 My name is Danica Marinac-Dabic. I am the
13 Chief of the Epidemiology Branch at the CDRH's Office
14 of Surveillance and Biometrics, and the Epidemiology
15 Unit is in charge of the review, monitoring and
16 tracking of post-approval studies also known as
17 conditional approvals that is. We're also in charge
18 of post-market surveillance studies also known as
19 Section 522 studies, another way how the FDA can ask
20 for additional post-market data, and also our unit is
21 in charge of FDA funded epidemiologic research
22 studies, which are the studies that are independently
23 conducted by the FDA to obtain additional post-market
24 data on the approved products.

25 As members of the CDRH Expert Advisory

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 Panel, you play a crucial role in our decision making
2 in terms of approval of medical devices and also very
3 often you make recommendations for so-called
4 condition of approval studies or post-approval
5 studies that are imposed by the PMA Order, and I know
6 that this afternoon, you're going to be also engaging
7 in some of the discussions about the post-approval
8 study for the product that is under discussion today.

9 So having that in mind, I would like to
10 give you today an update on significant changes that
11 had occurred in the CDRH Post-Approval Studies
12 Program during the last two to three years and also
13 specifically to give you a brief snapshot of what
14 urology devices post-approval studies we currently
15 have in place and what is their status. I'm not
16 going to go into a lot of details but just to give
17 you an idea of the information that is publicly
18 available on the status of our studies.

19 As you know, FDA can impose post-approval
20 requirements at the time of the approval of the PMA,
21 and this slide just lists that authority, and just to
22 make sure that we all understand that these studies
23 are done for continuing evaluation and reporting of
24 the safety, effectiveness and reliability of the
25 device for its intended use.

1 And this is certainly very important
2 statement because we would like to know up front that
3 post-approval studies should not be used to address
4 any pre-market questions. Anything that is essential
5 for the establishment of the reasonable assurance of
6 safety and effectiveness of the product has to be
7 demonstrated by the pre-market data.

8 This is just a summary of the need for
9 post-approval studies and, you know, during your
10 review of the submitted package, you know that, you
11 know, the sponsors usually gather a lot of pre-market
12 information that help you, the Panel, decide what
13 kind of recommendations you're going to make.

14 However, there still could be some
15 unanswered post-market questions for which post-
16 approval study route may be suitable. And these are
17 some of the reasons why we need the post-approval
18 studies. For example, we need to learn longer term
19 performance including the facts of re-treatments and
20 product changes. This is something for which post-
21 approval study can be asked for.

22 Now, as the devices are moving from the
23 highly trained and best clinical centers as is
24 usually chosen by the sponsors during the pre-market
25 trial, we would like to know how these devices are

1 performing in the community type of practices. And
2 sometimes the effectiveness of training programs can
3 be assessed by the post-approval studies, or if the
4 pre-market data did not have sufficient information
5 of subgroup performance, we can ask for the post-
6 approval study to address that question.

7 And in our intent to reduce the burden on
8 the sponsors and provide the least burdensome
9 mechanism for them to demonstrate reasonable safety
10 and effectiveness, we also consider post-approval
11 studies to obtain the longer term data in the broader
12 population for the intended use.

13 And very often Panel members bring their
14 thoughts, based on their experience, to us that we
15 didn't think of as the peer review team, and we would
16 like to incorporate those in the post-approval
17 studies.

18 It's very important again to state that the
19 objective is to evaluate the device performance and
20 potential device-related problem in a broader
21 population over an extended period of time after the
22 pre-market establishment of the reasonable assurance
23 of the device safety and effectiveness, and post-
24 approval studies should not be used to evaluate
25 unresolved issues from the pre-market data.

1 And again, it's always the balance between
2 less burdensome evidence to support pre-market
3 approval and assurance of continuous product safety
4 and effectiveness.

5 Well, a couple of years ago, we looked into
6 our post-approval studies across the Center, and
7 after this evaluation, we, the Center, had decided to
8 start major post-market transformation in this area
9 with the following major goals. To enhance
10 scientific rigor of post-approval studies and also to
11 establish and maintain the accountability for the
12 post-approval study commitments. We also wanted to
13 build the post-approval studies information
14 management system and to build a bridge between the
15 pre-market information and the post-market
16 information and certainly to feedback what we learn
17 in the post-market arena to our pre-market colleagues
18 that are reviewing the new PMA submissions.

19 And, finally, we wanted to increase our
20 transparency with the public, to make sure that all
21 interest of stakeholders have timely access to the
22 publicly available information on the status of the
23 post-market studies.

24 So these are the major areas that we tackle
25 in this post-market transformation effort. We have

1 changed the oversight. We developed a tracking
2 system. We made major changes in the review of post-
3 approval studies. We issued a guidance document. We
4 developed a web posting page that lists all the post-
5 approval studies, and we started this Post-Advisory
6 Panel Updates, again with an intent to give the Panel
7 the most recent information about these changes. And
8 we undertook a comprehensive strategit to build
9 public health partnerships with the clinical
10 communities, with manufacturers, with CROs, academia,
11 to help us better understand the design and the
12 conduct of post-approval studies.

13 So historically, post-approval studies were
14 housed in Office of Device Evaluation in the Review
15 Divisions, and in 2005, the initial transfer had
16 occurred and completed in 2007, where the post-
17 approval studies review was integrated into one post-
18 approval studies program housed in the Office of
19 Surveillance and Biometrics. And again, the
20 epidemiology group was the one who was given the
21 charge for this review and monitoring.

22 In addition to the oversight changes, as I
23 said, we have developed and instituted automated
24 tracking system which is based on the timelines that
25 are agreed upon between the sponsors and the FDA at

1 the time of the approval. And, based on that
2 timeline, our tracking system is designed to make
3 sure that we keep track of all the submissions, when
4 they come, and make sure that we remind sponsors that
5 if they do not comply with our commitments,
6 everything is posted on our webpage, and I'm going to
7 talk about that in a little bit.

8 So these are the major changes that have
9 occurred in the pre-market review process.

10 Again, in the past, the epidemiology group
11 was not part of the review process, and we would have
12 received information about the approved product and
13 then designed the post-approval study, helped to work
14 with sponsors on the design of post-approval studies.

15 The major change had occurred when an
16 epidemiologist was added to each PMA review team. So
17 now we are part of the pre-market review process, and
18 our role is to identify those reasons for the post-
19 approval studies, the rationale for the post-approval
20 studies and then work interactively with the sponsor
21 to help them design good study that is based on good
22 solid hypothesis, that has good, clearly stated
23 objectives and certainly our goal is to increase the
24 scientific rigor of those studies. And you will hear
25 also, if the device goes to the Panel, the

1 epidemiologists will be part of the Panel
2 presentation as well.

3 The goal is certainly to finalize the
4 protocol before the product is approved. So once the
5 approval order had been issued, the sponsor is ready
6 to go and apply for the IRB approvals and can start
7 the study.

8 What happens when the product is approved?
9 Again, the epidemiology group has the lead on all
10 post-approval study reports and supplements that
11 involve changes to the protocol, but we make sure
12 that our pre-market colleagues stay informed and
13 engaged in this review process, and this strategy is
14 designed to actually couple the epidemiologist's
15 expertise and observational study design and the
16 technical expertise that resides in our pre-market
17 office.

18 We issued a guidance document in 2006 and
19 with a slight revision in August of 2007, and in that
20 guidance document, we clearly stated what are the
21 reporting status definitions. As you can see, they
22 can range from on time, overdue or report already
23 received, and also these are our study status
24 definitions, and all of these things are also
25 available in our webpage, with a goal to increase the

1 transparency with the public.

2 This is the webpage that went live in on
3 April 6 of 2007, and all the reporting schedule
4 status and post-approval study status is on that
5 webpage.

6 We currently have over 120 studies, but
7 only those that initiated post-2005 are on the
8 webpage, and this is how the webpage looks like. You
9 can see there is an application number and the name
10 of the applicant, device name, also the medical
11 specialty. We also have study commitments listed
12 there, and when the protocol was approved and what is
13 the study population under the study and what is the
14 current status.

15 This database, this is linked also the PMA
16 database. So you can also search for more
17 information out there. And this is constantly
18 updated based on the feedback from our stakeholders.

19 Now, as I said, it's very important to
20 close the role certainly with Panel members as they
21 are part of our decision making process, and we
22 instituted these two initiatives, to prove this
23 general post-approval study updates at the beginning
24 of every Panel meeting. So first update was
25 presented in November 2007, and since then, we have

1 these updates at every Panel meeting. Those are so-
2 called general updates.

3 If there is a specific issue that we would
4 like the Panel to discuss, then we have a different
5 strategy. We have so-called specific post-approval
6 study updates, and they're a little bit more formal.
7 We invite the sponsor to give their presentation, and
8 we are also giving FDA perspective on this, and then
9 we come up with questions that we would like to get
10 Panel's input, but certainly that's a more joint
11 effort between the sponsors and the FDA.

12 And as I said, we started huge effort and
13 devoted significant amount of resources to build
14 public health partnership. This is just an example
15 of first in series of conferences that we cosponsored
16 with the Food and Drug Law Institute last year, and
17 there are two that are being planned for this year
18 and next year, where we invited the prominent members
19 of a clinical community also, you know, contract
20 research organizations, certainly lawyers,
21 manufacturers, Panel members, to talk about post-
22 approval studies and give us their feedback.

23 I would like also to tell you that we are
24 looking into innovative approaches, how we can design
25 those studies better and how we can use existing

1 databases to satisfy some of the post-market study
2 commitments.

3 Now, very briefly, I'm not going to go into
4 a lot of details, this is, you know, how many PMAs,
5 original PMAs and Panel track supplements were
6 approved in the period from 2005 to 2008, and how
7 many of those, as you can see, there were 5 approved
8 original PMA supplements with the post-approval
9 study. This is for all gastroenterology and urology
10 devices, and this is how the picture looks like for
11 the urology devices only. We have four ongoing post-
12 approval studies issued, at the time of the three PMA
13 approvals, and this is the list of the urology
14 devices for which the PMA approval asked for the
15 post-approval study. I'm sure that many of you had
16 participated in the discussion and the
17 recommendations. So I'm not going to go into a lot
18 of details.

19 This is just a brief summary of what type
20 of study design we are using. Certainly, there are
21 various design strategies that we can use in the
22 post-approval study. Sometimes we use the registry
23 as a framework to -- the study, and this is one
24 example, and you can see here what the objectives of
25 this particular registry is and what type of end

1 points and what type of duration we're talking about
2 here, also to say that these studies currently listed
3 on time on our webpage.

4 For Macroplastique, we have two different
5 post-approval study commitments. One is the real
6 time observation of safety and effectiveness
7 registry, again with 275 subjects, to look at the
8 durability of treatment effect and impact of re-
9 treatment, 5 years again follow-up, and this study is
10 on time as well.

11 This is the second piece for
12 Macroplastique, which is enhanced surveillance
13 system. Again, we try to compliment the post-
14 approval studies with the enhanced surveillance
15 system to try to get gather as much as possible of
16 the post-market data.

17 And finally, the -- this study was listed
18 overdue until recently, and the protocol was just
19 approved last week. Again, this was an example when
20 there were some challenges on implementing and
21 starting and drawing patients into the study, and we
22 had revisited our original approval order and worked
23 with the sponsor to design the study that will
24 address your post-market questions but would be less
25 burdensome to conduct.

1 And again this is just a summary of what
2 type of study designs we use. As you can see, we
3 have registries or perspective one-arm studies or
4 enhanced surveillance system in this case, and as far
5 as how diligent the sponsors are in terms of sending
6 their reports to the FDA, there were two reports that
7 were overdue, but they were received and they're
8 marked as such on our website, and two other reports
9 are on time.

10 And studies are on time and, you know, for
11 three of those studies and the protocol recently
12 approved which means the study will start very soon.

13 And this is again just a recap and this is
14 my last slide. This is our vision. We would like
15 for the Panel to know that post-approval studies, our
16 vision is to have studies that answer only important
17 post-market questions, not just the questions that
18 the FDA staff may be curious about but really the
19 ones that are important public health questions
20 because we understand the burden and the cost and the
21 effort that the sponsor have to put into making these
22 studies success. We would like those studies to be
23 realistic and founded on good science, and timely,
24 accurate and provide use for results that we can then
25 incorporate into labeling changes.

1 We also would like reports to be clearly
2 identified, effectively tracked and we certainly are
3 committed to keep our stakeholders apprised. I would
4 like again to say that nothing can be accomplished in
5 terms of this vision if we do not continue
6 cooperating with our pre-market colleagues. That's a
7 key for our success, and if we continue with this
8 effort, we believe that the enforcement actions will
9 be rare. We ask for those when it is necessary but
10 by proactively addressing these issues, I think we
11 are going to have not that frequent cases when we
12 need to do the reinforcement. Thank you very much.

13 DR. TALAMINI: Thank you, Dr. Marinac-
14 Dabic. That was extremely helpful.

15 We'll now proceed with the open public
16 hearing portion of the meeting.

17 Both the Food and Drug Administration and
18 the public believe in a transparent process for
19 information gathering and decision making. To insure
20 such transparency at the open public hearing session
21 of the Advisory Committee meeting, the FDA believe
22 that it is important to understand the context of any
23 individual's presentation. For this reason, FDA
24 encourages you, the open public hearing or industry
25 speaker, at the beginning of your written or oral

1 statement, to advise the Committee of any financial
2 relationship that you may have with the sponsor, its
3 product, and if know, its direct competitors.

4 For example, this financial information may
5 include the sponsor's payment of your travel, lodging
6 or other expenses in connection with your attendance
7 at the meeting. Likewise, FDA encourages you at the
8 beginning of your statement to advise the Committee
9 if you do not have any such financial relationships.
10 If you choose not to address this issue of financial
11 relationships at the beginning of your statement, it
12 will not preclude you from speaking.

13 Prior to the meeting, we have received one
14 formal request to speak during today's open public
15 hearing sessions. Our speaker is Mr. Bob Lipman.
16 Please come forward to the microphone. We ask that
17 you speak clearly into the microphone to allow the
18 transcriptionist to provide an accurate record of
19 this meeting, and we have about five minutes to stay
20 on track. Mr. Lipman.

21 MR. LIPMAN: Good morning. My name is
22 Robert Lipman. I am representing the Bladder Cancer
23 Advocacy Network, and I have no financial
24 relationship with the sponsor.

25 Thank you for the opportunity to speak here

1 today and to share my experience as a bladder cancer
2 survivor and patient. Like many who are ultimately
3 diagnosed with this disease, I was initially
4 misdiagnosed in 2003 by my internist before my
5 urologist was able to confirm with a cystoscopy and
6 bladder biopsy that I indeed did have bladder cancer.

7 After having the cancer removed and some
8 complications with bleeding, I received BCG treatment
9 once a week for six weeks. I did not experience any
10 side effects from the BCG treatment although it is
11 never pleasant to get treatment through a catheter.

12 A subsequent bladder biopsy showed that the
13 cancer had returned. Since BCG was not effective in
14 preventing the cancer from returning, an alternative
15 treatment choice at the time was BCG+Interferon. A
16 study had recently been published indicating that BCG
17 and Interferon could be effective in patients who
18 failed with BCG alone. BCG+Interferon is not an FDA
19 approved treatment.

20 I had much more severe side effects with
21 BCG+Interferon including intense irritation of my
22 bladder and extreme exhaustion. After getting
23 treated with BCG+Interferon once a week for six
24 weeks, another bladder biopsy showed that the cancer
25 had returned and again it was removed.

1 After failing with BCG+Interferon, what
2 were my treatment choices? I received a second
3 opinion at Johns Hopkins, and we made the unusual
4 decision to repeat the BCG+Interferon even though it
5 hadn't been effective.

6 Fortunately, the next biopsy in October
7 2005 was clear, and I've been cancer free since then.
8 However, bladder cancer requires frequent monitoring
9 and treatment. Since then I have undergone many more
10 bladder biopsies and almost 30 more BCG+Interferon
11 treatments.

12 Another side effect of BCG with Interferon
13 was dealing with my health insurance company.
14 Interferon requires special approval by the insurance
15 company, and approval is good for one year. After
16 having been approved for several years, I was denied
17 approval because the insurance company said that it
18 was not a FDA approved treatment. As a federal
19 employee, I was able to file an appeal with the
20 Office of Personnel Management who overturned the
21 insurance company's decision.

22 What happens if my bladder cancer returns?
23 What are my treatment choices? Of course, I want to
24 avoid at all cost having my bladder removed. I and
25 many other bladder cancer survivors and patients need

1 more safe and effective treatment choices. I urge
2 the FDA to encourage the development of new
3 treatments for this disease and to work quickly to
4 approve those treatments that are proven to be
5 effective.

6 While every cancer patient's story is
7 unique, there are several common issues that many of
8 us in the bladder cancer community share. The
9 disease is quite prevalent, currently the fifth most
10 commonly diagnosed cancer in the United States,
11 fourth among men. When I was diagnosed, I had never
12 even heard of it even though it is so common.
13 Despite its prevalence, many of us go undiagnosed or
14 misdiagnosed before finding out that we have bladder
15 cancer. A late diagnosis increases the chance that
16 the cancer will have invaded the bladder muscle wall,
17 and unfortunately medical treatments for muscle
18 invasive disease are limited and for most patients,
19 the standard of care is removal of the bladder as
20 well as the prostate in men, a life altering surgery.

21 Bladder cancer has a very high recurrence
22 rate, and bladder cancer patients must have ongoing
23 rigorous checkups and testing with respect to follow-
24 up care and, if necessary, treatment for recurrence.
25 Patients and their families bear the psychological

1 burden of knowing that the bladder cancer often comes
2 back and that there are limited treatments for the
3 disease, so that if the treatment does not work,
4 major life altering surgery may be the only option.

5 I began volunteering for the Bladder Cancer
6 Advocacy Network, BCAN, in 2005, to help raise the
7 awareness about bladder cancer and to be part of the
8 push for more research so that early detection
9 becomes a reality and that new treatment options are
10 available to patients who have been diagnosed.

11 BCAN is the first and only national
12 organization dedicated to raising awareness,
13 educating patients and clinicians and advocating for
14 more research into treatments for this disease.
15 BCAN's Scientific Advisory Board has more than 35
16 bladder cancer specialists, all urologists,
17 oncologists, radiologists or pathologists,
18 representing major cancer centers in the United
19 States and Canada.

20 On behalf of BCAN, and for all who are
21 currently living with bladder cancer, and for those
22 who have yet to be diagnosed, again I urge the FDA to
23 encourage the development of new treatments for the
24 disease and to carefully evaluate this application to
25 determine whether to recommend to the FDA that

1 Synergo should be approved as another treatment for
2 bladder cancer. Thank you.

3 DR. TALAMINI: Thank you very much,
4 Mr. Lipman. Is there anyone else in the audience who
5 would like to address the Panel at this time? Please
6 raise your hand and come forward to the microphone.
7 I'm not seeing any hands.

8 We will proceed with the sponsor
9 presentation for the Synergo SB-TS 101.1 Device and
10 Mitomycin C.

11 I would like to remind the public observers
12 at this meeting that while this meeting is open for
13 public observation, public attendees may not
14 participate except at the specific request of the
15 Panel.

16 We will begin with the sponsor
17 presentation. The first presenter is Dr. Yagel
18 Koren, Medical Director for Medical Enterprises, Ltd.
19 Dr. Koren.

20 DR. KOREN: Thank you very much, and good
21 morning. I would like to use this opportunity to
22 thank the FDA for giving us this opportunity to show
23 here today and the distinguished members of the Panel
24 for giving the time and effort to be here with us
25 today.

1 I would present the introduction for this
2 PMA approval. Our presenters today would be
3 Professor Michael O'Donnell who is the Director of
4 Urologic Oncology at the University of Iowa Hospitals
5 and Clinics, Professor Fred Witjes from the
6 Department of Urology in Radboud University, Nijmegen
7 Medical Center, in The Netherlands, Professor Barton
8 Grossman who is from the Department of Urology at the
9 University of Texas, MD Anderson Cancer Center,
10 Ms. Ahava Stein who is a regulatory consultant for
11 Medical Enterprises, and myself.

12 The order of our presentation today will be
13 first the introduction which will be given as I said
14 by myself. Later on, a brief overview of the disease
15 and current treatment options will be given by
16 Professor Michael O'Donnell. Device description and
17 overview of the clinical studies will be presented by
18 Ms. Ahava Stein, and later on, the overview of the
19 clinical studies will be presented by Professor
20 Witjes. Finally, Professor Barton Grossman will give
21 an overall summary of our application, and just to
22 add a note, we have been asked by the FDA to prepare
23 a plan for possible post-approval study should the
24 Panel recommend such one, and this plan will be
25 presented again by Professor Michael O'Donnell.

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 In my introduction, I will give you a few
2 short notes on the device itself, the history of the
3 company and the history our previous study, Study
4 101.1.

5 About the device. As we know, intravesical
6 chemotherapy has been widely used for a decades, both
7 in the U.S. and outside the United States for the
8 treatment of non-muscle invasive bladder cancer. But
9 its limited frequency gave reasons to find methods to
10 improve its efficacy. This thing in turn led to the
11 development of the Synergo hyperthermia device in San
12 Raffaele Hospital, which is located in Milan, in
13 Italy, back in the 1990s.

14 Our hyperthermia device is designed to heat
15 the bladder walls with this to increase the effect of
16 mitomycin for the treatment of bladder cancer.

17 Our history, in 1994, our pivotal study has
18 begun. It begun as the collaboration of three
19 investigators as an investigator-initiated study in
20 three academic centers in Italy and Israel. Only
21 three years later, Medical Enterprises was formed and
22 acquired the technology and the sponsorship of this
23 academic, scientific research. CRFs were then formed
24 and all the previous data was transcribed into these
25 newly formed CRFs. From 1997 until the end of the

1 study in 2001, all the data was already prospectively
2 registered on the CRFs.

3 In 2000, about eight years ago, our Synergo
4 device received the CE mark in Europe and the
5 approval of the Israeli Ministry of Health and ever
6 since is routinely used to treat patients with
7 bladder cancer in Europe and in Israel.

8 What about our company? Medical
9 Enterprises is a small company, and it has Synergo as
10 its only product. Should this PMA be approved, we
11 plan on cooperating with a local American company to
12 introduce the device into the United States.

13 Thank you for your attention, and now for a
14 few words on the disease and current treatment
15 options, I would like to call on Professor Michael
16 O'Donnell. Thank you.

17 DR. O'DONNELL: Thank you, Yagel.
18 Mr. Chairman, members and guests, it's my pleasure to
19 talk to you about the disease of bladder cancer.
20 Just to tell you about myself, I am a urologist. I
21 specialize in urologic cancer, particularly in
22 bladder cancer. I've served on committees with
23 national trials groups such as CALGB, to advise them
24 on trials for bladder cancer, and I've also reviewed
25 the AUA Guidelines Panel as a peer reviewer. I have

1 been a consultant for Medical Enterprises for about
2 three years, but I have no financial stake in the
3 company per se.

4 Let me tell you a little bit about the
5 disease. The bladder is obviously a cavity, and it's
6 lined by a surface epithelium, and it's this
7 epithelium that becomes malignant. And bladder
8 cancer occurs in three basic forms, a surface
9 spreading disease known as carcinoma in situ. This
10 is the minority, about 5 to 10 percent, and this is a
11 disease that is not being handled by the Synergo
12 device.

13 The vast majority of cases are this
14 papillary exophytic growth in the bladder that
15 projects into the lumen of the bladder. This is at
16 least two-thirds to three-quarters of the cases of
17 bladder cancer.

18 And the third type is this nodular form
19 which is a more ominous form. It can be in the
20 mucosa, the submucosa and eventually invade into the
21 muscle. This is the type that essentially kills
22 patients.

23 There are two major distinguishing factors
24 when we talk about bladder cancer. One is the stage,
25 and that is the depth of invasion of the disease, and

1 second is the grade which is the degree of
2 aggressiveness of the cancer. The stage is given in
3 a TMN, tumor, node, metastasis staging system as a
4 group here that includes CIS or Tis, Ta disease
5 limited to the mucosa, T1 disease that goes into the
6 submucosa. And we distinguish this group as being
7 the superficial group mostly because we can resect
8 this from an endoscopic approach. Anything beyond
9 that, including invasion into the muscle or through
10 the muscle into the fat or into other adjacent
11 organs, is the muscle invasive category that occupies
12 about 25 percent of the presentations of bladder
13 cancer. The Synergo is limited to the superficial
14 group.

15 We also distinguish the grade as I
16 mentioned. I'm just showing you here a picture of
17 low grade, these well-formed cells, papillae versus
18 this more anaplastic example of a high grade cancer.
19 Grade does correlate greatly with prognosis as well
20 and with progression.

21 Now, bladder cancer is a disease that is
22 not limited certainly to the United States but rather
23 has a certain distribution related to environmental
24 and occupational exposures. And so it tends to be
25 the same in the westernized countries. This is

1 demonstrating the male instance of bladder cancer in
2 North America. As you can see, a similar color code
3 for Europe, for Israel and for Russia.

4 Now, the disease is also very similar
5 between the two continents, between the U.S. or North
6 America and Europe with regard to the disease
7 characteristics, that is the percentage of
8 superficial versus invasive, the stage distribution
9 and the grade distribution. Likewise, the treatments
10 are very similar, such that in terms of diagnosis,
11 it's universally diagnosed by a cystoscopy, often
12 office based, and it's universally treated by a
13 procedure of going into the bladder and resecting or
14 removing the tumor as the first initial treatment.

15 And finally, and probably more importantly,
16 is that the recognition by two of the major governing
17 bodies, first the United States, the American
18 Urological Association, and the corresponding group
19 in Europe, The European Association of Urology, have
20 drafted guidelines that have been recently published
21 in November 2007, for the AUA and in March of 2008,
22 updated for the EAU, that represent consensus of
23 management that are remarkably similar, both in
24 regard how to initially diagnose, biopsy, resect and
25 then apply adjuvant treatment such as intravesical

1 therapy.

2 This is the process by which we diagnose
3 bladder cancer. It involves obviously placing a
4 resectoscope through the urethra in this case,
5 through the prostate of a male, and imaging the
6 internal surface of the bladder. This is the Gold
7 Standard and is felt to accurately detect bladder
8 cancer in over 90 percent of cases. We use the same
9 technology then by applying different instruments to
10 actually physically remove the tumors and this is
11 illustrated here in this cartoon but here's a real-
12 time image of a resection occurring through a hot
13 loop, through a papillary projecting bladder cancer
14 here.

15 The problem with doing surgery though is
16 that most patients recur. In fact, this represents
17 the aggregate results of that analysis from the
18 European community of about 1200 patients with a
19 reasonably good follow-up of about five years,
20 representing that there's an exponential decay curve,
21 which is similar what you see for many cancer
22 treatments, or cancer recurrences, and it levels off
23 somewhere around five years.

24 If you look at the event rate, the number,
25 the percentage of those cases that occur, most of

1 them occur early. In fact, over 50 percent within
2 the first year and close to 90 percent in the 2
3 years, and this is the reason why a 2 year point is a
4 reasonable point to evaluate the efficacy of
5 treatments on this disease.

6 Now, the high recurrence rate remains a
7 very significant problems, and it's for this reason
8 that adjuvant therapy or additional treatment,
9 intravesical treatment, placing medication in the
10 bladder, have become advocated by both national
11 governing bodies, the American Urological Association
12 and the European Association of Urology, and they've
13 done so according to basically a risk adapted policy,
14 putting patients into categories such that one can
15 apply more universal treatment recommendations so to
16 speak. The AUA does this by giving index patients
17 and the European Association of Urology has
18 formulated actually these three categories here which
19 I put for convenience.

20 The low risk group, which is by the way not
21 the group that is being treated by Synergo and
22 represents about half of the patients with bladder
23 cancer, are the single, solitary papillary low grade
24 small tumors, the best actors, and they have a
25 recurrence rate of about 40 percent within 2 years,

1 as opposed to the high-risk group, which has a very
2 high recurrence rate, over 70 percent within 2 years.
3 These include any high grade disease, any grade 3 Ta
4 or T1, any CIS, and a subset of grade 2 T1 disease
5 that a multifocal. And then the intermediate risks
6 which represents about a third of all the patients,
7 and these patients are essentially those that are
8 between the extremes. Any patient that has had a
9 recurrent cancer, that's papillary, non-high grade,
10 any patient with multifocal disease falls into this
11 group or any very large tumors. These are the two
12 groups that are being targeted for the Synergo
13 device.

14 Now, the guidelines then reflect this risk
15 adapted policy, advocated both in very similar
16 formats by the American Urological Association and
17 the European Association of Urology. So for the low
18 risk group, those small papillary tumors, it's felt
19 that one single treatment of medication, and the most
20 common one in the United States is mitomycin, is
21 given immediately after, within six hours of the
22 urethral resection, transurethral resection, and
23 that's felt to be recommended and sufficient. It's
24 listed as an option and a recommendation respectively
25 for the AUA.

1 The intermediate group which represents
2 about a third of the patients, includes that first
3 single dose of chemotherapy, more advocated by the
4 EAU than AUA, but relatively recognizes as being an
5 important first step, followed by additional
6 treatment. And here, both the guidelines recommend
7 either mitomycin C or equivalent chemotherapy or BCG,
8 often with a maintenance program, about 1 year with
9 BCG and 6 to 12 months for chemotherapy.

10 And finally we come to the high-risk group
11 which represents a subset, 15 percent, nonetheless
12 the most important ones because these are the ones
13 that go to progression. And here, the options really
14 are only BCG with maintenance, or viable options, or
15 cystectomy, removing of the bladder. Now, it used to
16 be in the original guidelines in 1999 and 2002, that
17 mitomycin C was listed here. That was until we had
18 sufficient information to understand that the results
19 are significantly inferior to BCG unless only BCG is
20 the non-surgical alternative for these high-risk
21 patients.

22 Unfortunately, even the best treatment that
23 we have, have serious limitations especially in terms
24 of efficacy. This graph represents a large southwest
25 oncology group study that was performed in the early

1 nineties, and what it shows here is a randomized
2 trial of BCG versus mitomycin C. It includes a small
3 percentage of patients with CIS but nonetheless these
4 results are representative of what happens also in
5 the U.S., that is patients relapse early and then
6 there is a plateau somewhere around two years or so.

7 And in this case, although there is a
8 strong difference between BCG and mitomycin, in fact,
9 for this reason, the study it shows statistical
10 significance. It was actually even stopped early as
11 part of an interim analysis, and I bring this up
12 because this study, in fact, recapitulates some of
13 the issues and the points that you'll see in the
14 Synergo program as well.

15 And these points include, number one, the
16 study was not a FDA design study. It wasn't one
17 brought to the FDA for approval, BCG, though the data
18 was used retrospectively to submit an approval and
19 for which BCG gained an approval for use against
20 papillary cancer.

21 Number two, it was stopped prematurely
22 because the safety analysis indicated that there was
23 such a significant difference here that it was
24 unethical to continue it any further.

25 Number three, is mitomycin C, although it's

1 not a FDA approved, recognized drug for this
2 indication, was actually used as part of the
3 indication to FDA and accepted as such.

4 And four, these studies were not blinded,
5 in fact. Studies of this nature, intravesical
6 therapy, to my knowledge, have never really been
7 blinded, and I know that will become an issue that is
8 raised later on.

9 The final point is that about half of these
10 patients do represent the Ta grade 1 tumors, and if
11 you take those out, both curves drove by about 10
12 percent. So, if anything, these curves are a little
13 bit optimistic in terms of what you can expect from
14 this form of therapy.

15 The second major problems with the forms of
16 therapy that we have right now relates to toxicity.
17 You don't get something for nothing essentially.
18 I've made two major comparisons here between
19 mitomycin C and BCG because these are the most
20 relevant issues that we talk about today. As you can
21 see, the toxicity can be related to local toxicity
22 which is mostly frequency, dysuria, irritative
23 symptoms, so-called urinary tract symptoms including
24 hematuria, incontinence and so forth. Systemic
25 symptoms, fever, flu-like symptoms or arthralgias,

1 skin rash, infectious related issues including
2 bacterial cystitis or UTI, and the continuation rate
3 or treatment discontinuation rate, incomplete or
4 interruption.

5 First, notice a couple of things. Wide
6 ranges from multiple studies, this was a table done
7 as part of the comprehensive analysis by the 1990 AUA
8 Guidelines Panel, represented the most complete
9 analysis of the literature at the time, wide ranges
10 throughout the studies. In general, a high amount of
11 lower urinary tract symptoms for both mitomycin C and
12 for BCG but generally higher for BCG by about 50
13 percent, a significant rate of hematuria as well, 16
14 to 29 percent. Some of the ranges here are up to a
15 third of the patients. Bladder contracture which is
16 a loss of functional capacity of the bladder, 5
17 percent in the mitomycin C group alone, a small
18 amount in BCG. Generally more systemic fever-like
19 symptoms with BCG which we expect, you know, from
20 therapy but nonetheless some in the mitomycin group.
21 Skin rash more prevalent in the mitomycin group.

22 Here I've highlighted the incidence of what
23 I would call more serious side effects, and the
24 problem with BCG is that it's a live, though
25 attenuated microbacteria. We give 10 to the 8th

1 organisms in the bladder once a week for 6 weeks.
2 Some of it can get into the blood stream and cause
3 serious infection, and that serious infection rate is
4 approximately 5 percent and includes some systemic
5 organ manifestations as well as more specific ones
6 like rumentitis.

7 And finally as you can see with both
8 therapies, there is a real discontinuation or
9 interruption rate that hovers around 10 percent.
10 There are very rare cases, by the way, of lethal
11 consequences of BCG, an overwhelming TB-like
12 infection. Fortunately, they're relatively rare but
13 they remain in the back of everyone's mind that gives
14 BCG.

15 Now, bladder cancer, you've heard a little
16 bit from our public speaker, is a very significant
17 disease, and just to put this in perspective, there
18 are about 69,000 new cases expected this year in the
19 United States. That's just the tip of the iceberg
20 because for every patient with a new diagnosis, there
21 are 9 to 10 more patients living with the diagnosis.
22 So the estimate is over 500,000, half a million, in
23 the United States with this disease. About 14,000
24 cancer deaths.

25 Bladder cancer is also not appreciated that

1 it is extremely costly. Patients who live a long
2 time have multiple treatments as you've heard about
3 and the average cost is estimated to be somewhere
4 between 100 and \$200,000 per patient from the time of
5 diagnosis to death from whatever cause.

6 This results in or is a reflection of this,
7 many procedures in the U.S., 300,000 over, for just
8 the surgical procedure per year of scraping out the
9 tumors in the bladder and an estimated over 2.5
10 million instillations of various drugs, mitomycin,
11 BCG, Interferon, et cetera, for this disease.

12 So why do we need new treatments? Well,
13 for the intermediate risk patients, it's certainly
14 the case that neither BCG nor mitomycin C provide
15 reliable long-term disease relapse rates.

16 The same problem is that for high risk, but
17 in the high risk, we have even less options. All we
18 have is BCG or we have the removal of the bladder
19 which is not a very acceptable option. What about
20 the patients that can't get BCG because they have
21 immune related issues or they develop intolerance or
22 significant side effects? We really don't have
23 anything to offer this.

24 And just as a personal note, I've been
25 dealing with patients with bladder cancer for over 15

1 years. I have, you know, talked with them. I've
2 given lectures all over the country, talked with the
3 patients and the physicians. There's a real need
4 there for some alternatives for some advancement.
5 There have been no significant improvement in bladder
6 cancer treatment in the United States since the
7 approval of BCG in the early 1990s. And we've got
8 this disease, we've got a high prevalence. We have
9 patients still suffering from this.

10 I also can tell you that I've been looking
11 at the Synergo program for over 10 years, when I was
12 originally asked to do so as part of a due diligence
13 from a different company that was investigating
14 whether they should invest in this group. I've
15 stayed in contact. I have reviewed the papers for
16 the journals. I've seen their presentations at the
17 AUA. What really impresses me is that you get
18 consistency and the efficacy of the safety for this
19 device. I talked with the -- investigators. I've
20 been to Milan. I've seen the patients being treated.
21 I've talked with the patients, and it's the real
22 thing. And I hope you will stay open minded, just to
23 hear the results that we show you today, to show you
24 that we really need to have a new treatment for
25 bladder cancer in the United States. Thank you.

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 MS. STEIN: My name is Ahava Stein, and I'm
2 a regulatory consultant to Medical Enterprises. I do
3 not have a financial interest in the company, and I
4 appreciate the opportunity to address the Panel.

5 I'll begin with the indications for use for
6 the Synergo Device. The Synergo hyperthermia device,
7 in conjunction with mitomycin C is intended for use
8 for prophylactic treatment of patients with STCCB,
9 following endoscopic removal of Ta-T1 and G1-G3
10 tumors. The Synergo treatment is clinically

11 indicated for intermediate and high-risk patients

12 The Synergo treatment delivers heat to the
13 urinary bladder wall using RF energy and the Synergo
14 hyperthermia is delivered concomitantly with cooled
15 mitomycin C drug solution.

16 In this picture here, you see the Synergo
17 catheter system inserted into the bladder. You see
18 the catheter balloon fixating the catheter within the
19 bladder. The thermocouples here, when deployed,
20 extend into the bladder wall to monitor the
21 temperatures at the bladder wall. There are an
22 additional two thermocouples below the balloon which
23 monitor the temperatures within the urethra. The
24 antenna, along the catheter, emits RF energy to heat
25 the bladder wall.

1 The Synergo catheter performs basically
2 three functions. One, as I mentioned, uniform
3 heating of the bladder wall. The thermocouples
4 monitor the temperature of the bladder wall and
5 circulation of the mitomycin C drug solution into and
6 out of the bladder.

7 I'm going to discuss some of the
8 preclinical testing that was conducted on the Synergo
9 device, the catheter and the mitomycin drug. The
10 testing with the mitomycin C drug was submitted as
11 part of the NDAs, submitted by the drug
12 manufacturers, Bedford Laboratories and Bristol Myers
13 Squibb. Medical Enterprises has been granted letters
14 of authorization from the drug manufacturers to
15 reference the data contained in the mitomycin C NDAs.

16 Additionally, the company conducted their
17 own pharmacokinetic study, that is the Paroni Study,
18 which assessed the effective local hyperthermia on
19 the systemic absorption of mitomycin C during a
20 Synergo treatment. This study showed that
21 hyperthermia caused an increase in mitomycin C
22 penetration into the bladder wall and showed that
23 there were higher plasma concentrations of MMC, in
24 the group administered hyperthermia in conjunction
25 with MMC, compared to the group receiving MMC alone.

1 Despite the increase in the penetration of
2 the mitomycin C into the bladder wall and the higher
3 concentrations of MMC, the plasma levels of MMC at
4 twice the indicated dose, the highest MMC plasma
5 concentration was still 6 times less than the
6 critical toxic systemic dose that causes bone marrow
7 suppression.

8 So in summary, this test demonstrates that
9 the hyperthermia treatment enhances the mitomycin C
10 uptake into the bladder wall, while maintaining MMC
11 plasma concentrations that are still well below the
12 toxic levels.

13 Another study that the company conducted
14 evaluated the degradation of mitomycin C when heated.
15 Mitomycin C was dissolved in different IV fluids
16 including dextrose and NACL, when heated to 50
17 degrees, temperatures that are higher than
18 administered during a Synergo treatment. The results
19 demonstrated that the MMC degradents were below the
20 Gensia Sicor specification limits for these
21 impurities.

22 Testing of the Synergo device included
23 mechanical and electrical safety testing, and
24 electromechanical compatibility testing as well as
25 software validation according to international

1 standards. All of the tests passed and the results
2 of the tests met the requirements of these
3 international standards.

4 The catheter component of the device was
5 tested according to the ASTM standard for Foley
6 Catheters. The materials of the catheter were tested
7 for biocompatibility according to the ISO 10993
8 standard. Bench testing of the catheter included
9 mapping of the electromagnetic field in a liquid
10 phantom to show that there was minimal absorption of
11 the RF energy by the liquid.

12 And a second bench test in a simulated
13 tissue model demonstrated that the RF energy
14 generated by the antenna homogeneously heated the
15 bladder wall and then the temperature rapidly
16 decreased over four to six millimeters across the
17 bladder wall. This was done in a bench test and then
18 further validated in the animal study.

19 The animal study was conducted to
20 demonstrate that during normal treatment conditions,
21 there are no risks to the bladder tissue or to the
22 adjacent organs. The sheep model was chosen as the
23 sheep bladder is similar to the human bladder.
24 Temperature mapping was conducted of the internal
25 bladder wall as well as the external bladder wall and

1 the adjacent organs during a Synergo treatment.

2 At the end of the study, the animals were
3 sacrificed and pathological evaluation of the bladder
4 organs, the bladder tissue and the adjacent organs
5 were compared to control animals.

6 The results of the animal study
7 demonstrated that the Synergo thermocouples
8 temperature measurements were accurate when compared
9 to independent temperature monitoring system. The
10 results of the study also demonstrated that the
11 temperature was homogeneous over the bladder wall and
12 decreased by a magnitude of three to five degrees
13 over the bladder wall, internal to the external
14 bladder wall, and a temperature drop of five to seven
15 degrees was measure at the adjacent organs.

16 And finally, the animal study demonstrated
17 that there were not risks of irreversible damage to
18 the bladder or to the adjacent organs at temperatures
19 that were higher than administered during a normal
20 Synergo treatment.

21 I will turn this over now to Dr. Witjes for
22 a summary of the clinical studies.

23 DR. WITJES: Okay. Ahava, thank you very
24 much. Dr. Talamini, Panel members and guests, thank
25 you very much for the opportunity to speak here about

1 this technology. My name is Fred Witjes. I'm a
2 oncological urologist working in The Netherlands, and
3 I'm involved in the treatment of bladder cancer
4 patients, guidelines and things like that, as you can
5 see. As such, I am already treating patients with
6 this machine since 2001. So I have around seven
7 years of clinical experience with this device, and I
8 am a principal investigator of one of the studies
9 that I will talk about, namely Study 102.1, but that
10 will come later.

11 I'm very thrilled that we finally are able
12 to present this data for the American people since
13 I've been treating already for seven years, people in
14 Europe with the device, and I'm very impressed by the
15 results we've had so far and hope we can achieve
16 similar results in American patients.

17 These will be the studies that I will be
18 shortly addressing. The first 101 and 102 are
19 randomized controlled trials which will be used for
20 safety and efficacy. Then we have the European
21 prophylactic patients group which is an uncontrolled
22 commercial use dataset which we also will be using
23 for safety and efficacy. And then you can see that
24 there are three smaller trials, that are listed
25 below, uncontrolled and one controlled study which we

1 will only use for safety data.

2 I'll first start, of course, with pivotal
3 study, the 101.1. The objectives of the study were
4 to compare safety and efficacy of the Synergo
5 hyperthermia versus mitomycin C alone for
6 prophylactic treatment of superficial bladder cancer.

7 The primary endpoint was the comparison of
8 the recurrence rate at two years. You've heard from
9 Professor O'Donnell that this is a very common and
10 also very logic endpoint for superficial bladder
11 cancer trials.

12 Secondary endpoints obviously were
13 comparison of progression, stage and grade,
14 comparison of the occurrence of CIS, comparison of
15 the occurrence of upper urinary tract tumors or
16 tumors in the prostatic urethra, and finally
17 comparison of the occurrence of distant metastasis.

18 The sample size calculation was initially
19 based on the primary endpoint of the two-year
20 recurrence rate. The initial assumptions for
21 calculation were a 2 year recurrence rate in the
22 mitomycin C only group, the control group of 40
23 percent based on the scientific literature of the
24 '90s. The study was designed to detect a 50 percent
25 reduction in this recurrence rate in the Synergo

1 group versus the mitomycin C control group with a
2 power of 80 percent and a 5 percent level of
3 significance. So the initial sample size calculation
4 case up with 158 patients.

5 The protocol did call for interim analysis
6 when 80 patients completed the one-year follow-up.
7 It was done a little bit earlier than planned due to
8 ethical reasons, and you can see why, because the
9 interim analysis clearly showed a major difference
10 between the recurrence rate in the Synergo which was
11 11 percent versus the recurrence rate in the
12 mitomycin C which was 62 percent. That resulted in a
13 recalculated sample size which was now 84.

14 One of the things we have to address is
15 there have been some randomization errors in this
16 101.1 study. Five pairs of administrative
17 randomization errors were done at the central
18 randomization office. The clinical sites nor the
19 sponsor were aware of these errors until years after
20 the trial had closed. The total number of patients
21 in each group obviously remains unchanged because
22 they were paired randomization errors, and realizing
23 that there have been some randomization errors, we
24 reanalyzed the study. We have not only looked at how
25 they were treated but also we looked at the results

1 as they were randomized and then still you see that
2 there was a significant advantage for the Synergo
3 treatment in this case, smaller than 0.01.

4 There were some protocol deviations and
5 withdrawals from the study. We had five Synergo
6 patients withdrawn from the study. Three withdrew
7 consent prior to receiving any treatment. For
8 example, one patient from Sicily who lived too far
9 from the hospital. One physician withdrew the
10 patient from treatment due to deteriorating health
11 before starting treatment, and we have one patient in
12 the Synergo group with skin allergic to mitomycin C.
13 That is one of the things that does happen, and you
14 can also see that in the mitomycin C group, there's
15 one patient with skin allergy. Two additional
16 Synergo patients were not included in the Per
17 Protocol cohort due to major protocol deviations.

18 This results in the following table. You
19 see that the top line is all study patients
20 randomized as treated. It's 42 in Synergo arm and 41
21 in the mitomycin C arm.

22 The second line shows you the randomization
23 as intended. So that means that some of the patients
24 who did receive Synergo treatment were put in the
25 mitomycin C arm and the other way around. That

1 results in 36 patients in the Synergo arm and 41 in
2 the mitomycin C arm.

3 And the third line shows you randomized as
4 treated, and again you have 37 in the Synergo arm and
5 40 in the mitomycin C arm.

6 In the second and third column you see that
7 you missed sort of one patient, but one patient who
8 had a randomization error also had an allergy to
9 mitomycin C. So that's why you skipped in this small
10 table and that small table.

11 After two protocol deviations, you are left
12 with the protocol analysis which has 35 patients in
13 the Synergo arm and 40 patients in the mitomycin C
14 arm.

15 What did we do after randomization?
16 Patients were randomized for Synergo or mitomycin C
17 therapy which is 2 times 20 milligrams. We had eight
18 weekly sessions and after that, four monthly
19 treatment sessions. Follow-up was as stated, two
20 years, and we did that every three months in the
21 first two years. And the endpoint assessment was a
22 histologically proven, biopsy proven tumor
23 recurrence. So not only by visual cystoscopy but
24 biopsy proven.

25 Here you can see the three centers which

1 have been doing this trial. They are, of course, --
2 centers, one from Milan, one from Palermo and one
3 from Israel.

4 The clinical data of the study, we have
5 been monitoring them and 100 percent of those CRFs
6 according to the GCP requirements that was performed.
7 There was a FDA audit in 2005 of all the sites that
8 were in this study and they confirmed that the CRFs
9 were an adequate reflection of source documentation.
10 So we think that safety and efficacy data are
11 adequately captured on CRFs. Since we did
12 retrospective until '97 and prospective after '97, we
13 compared those two datasets, and we see that there is
14 a consistent reporting of the adverse events
15 throughout the study. So even before and after 1997.

16 Here you see the baseline characteristics
17 of our patients. You can see that they are well
18 balanced between the two Synergo groups, and despite
19 of the facts that we have for sample patients with --
20 result, for example patients with the -- and for
21 example also patients with -- lung tumors, according
22 to the EAU risk criteria, none of those patients was
23 in the low risk category, half of them approximately
24 in the intermediate risk category and half of them
25 approximately in the high-risk category.

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 Blinding has been an issue. Well, as
2 Professor O'Donnell already told you, investigator
3 blinding is not typically performed or actually not
4 performed at all in intravesical therapy trials
5 published in scientific literature. As Chairman of
6 the Non-Muscle Invasive Bladder Cancer Group, we've
7 done many trials in the -- to see. We've never done
8 on with blinding. The pivotal study submitted to the
9 FDA, for FDA approval also, were not blinded, for
10 example, the BCG trial or the Valrubicin trial.

11 Moreover, it's very difficult to blind for
12 Synergo treatment because patients are obviously
13 aware of the treatment they get. They feel the heat
14 during the treatment and also the urologist who
15 checks the patient is obviously aware of the
16 treatment that the patient had because of the thermal
17 effects that you see during cystoscopy.

18 Moreover, the long-term results of this
19 trial, the Synergo trial, actually confirm that the
20 study results that we have, for example, after two
21 years were not biased, but they are absolutely
22 consistent.

23 If you look at the efficacy results, you
24 see here three scenarios. The first line is
25 evaluated as treated. You see that the scenario

1 group has a two-year recurrence rate of approximately
2 19 percent versus 62 in the mitomycin C arm,
3 obviously statistically very significant.

4 Then we also looked, of course, at
5 randomized as intent so as they should have been
6 treated. That means that some of the patients who
7 had Synergo therapy had a low recurrence rate versus
8 skipped through the mitomycin C arm and the other way
9 around. That obviously results in a little bit
10 higher recurrence rate in the Synergo arm, 25 percent
11 and a little bit lower recurrence in the mitomycin C
12 arm, 55 percent, but still this is significant and
13 the protocol analysis also again shows you a very
14 significant advantage for the Synergo arm, around 17
15 percent versus 62 percent.

16 So Synergo treatment was consistently
17 significantly better than the mitomycin C in all
18 these patient analysis or patient populations.

19 Here you see the Kaplan-Meier curve for the
20 treatment as randomization is treated. You see that
21 again after approximately two years there is a sort
22 of leveling out of the results in the Synergo group
23 and that there's a constant drop in the patients who
24 do not have recurrence in the mitomycin C group, and
25 you see the final difference. This is up until two

1 and a half years is very clear and statistically
2 significant.

3 Here you have a similar curve for
4 randomized testing as intended. You see that the
5 difference is a little bit less, but still there is a
6 very large difference between the two treatment
7 groups.

8 And finally you have here the evaluation
9 per protocol again showing you a highly statistically
10 significant and clinically relevant advantage for the
11 Synergo treatment.

12 We also tried to come up with the worst-
13 case scenario. For both of these worst-case
14 scenarios we assumed that the one patient that
15 dropped out of the mitomycin C group would have been
16 recurrence free at two years which, of course, might
17 be possible. And we also assumed that the five
18 patients that had Synergo treatment would have had
19 disease recurrence at first follow-up. Of course,
20 it's not very realistic but it was the worst-case
21 scenario we could think of and even then you see that
22 in the scenario group, the recurrence rate is around
23 31 percent and the mitomycin C group is around 60
24 percent and that remains only clinically relevant but
25 also statistically significant.

1 If you then also like the FDA has done
2 reversed the treatments, so again put the mitomycin C
3 treated patients in the Synergo group and the other
4 way around, due to the randomization errors, you
5 still see that there is an advantage with the Synergo
6 treatment, but then the statistical significance is
7 lost. However, I think this is not a very realistic
8 scenario.

9 We had some secondary endpoint analysis,
10 and they didn't really reveal any difference. There
11 was no progression in tumor stage or grade. There
12 was no occurrence of CIS in the Synergo group. No
13 patients had carcinoma in the upper urinary tract or
14 the urethra in the Synergo group, and no patients had
15 occurrence of distant metastasis in the Synergo group
16 nor in the mitomycin C group.

17 However, if we look at longer-term follow-
18 up because this predominantly looks at the first two
19 years, in the long-term follow-up, we had three
20 patients who had distant metastasis in the control
21 group.

22 Here you see the long-term efficacy
23 analysis. You see that it goes beyond 12 years and
24 actually just like at the beginning of the curve, you
25 see that there is a consistent difference between the

1 two treatment arms in favor of the Synergo treatment.
2 You'll see here a 10 year recurrence rate of 48
3 percent in the Synergo group and 10 year recurrence
4 rate of 85 percent in the mitomycin C group. Again
5 illustrating by the way, that it is a very nasty
6 disease in patients and tends to come back very
7 often.

8 Long-term follow-up also shows some results
9 with regard to overall mortality or radical
10 cystectomy, two in the Synergo group and five in the
11 mitomycin C group. Overall mortality five in the
12 Synergo group, nine in the mitomycin C group. None
13 were treatment related. None were disease related.

14 We did a subgroup analysis with regard to
15 some of potential significant predictive factors.
16 You can see the list there. There was no significant
17 effect of any of those listed there like H gender,
18 number of previous occurrences, et cetera, et cetera,
19 of the efficacy analysis. However, we did find
20 significant effect of the history of recurrence. So
21 the first episode, recurrent, first is recurrent to a
22 high recurrence and also of the EAU risk category,
23 which is not surprising because in the EAU risk
24 category, one of the very important things is this
25 history of recurrence, and you'll see that even if

1 you adjust for these prognostic factors, that the
2 Synergo treatment remains better than the mitomycin C
3 group, giving the impression that the results are a
4 little bit more outspoken if you use that in higher
5 risk patients.

6 The next thing after talking about the
7 results is talking about the side effects. You see
8 here that many of the side effects, there was no
9 statistically significant difference. Like for the
10 dysuria, hematuria, tissue reaction, urethral
11 stenosis and skin allergy, urethra urinary tract
12 infection and bladder wall necrosis, two were in
13 favor of mitomycin C only which was pain and
14 posterior wall tissue reaction. The ones highlighted
15 in yellow, I will discuss in one of the next slides.

16 Other adverse events are very seldom.
17 You'll see, for example, anxiety, amnesia, and then
18 hypertonic bladder which has been noted in one
19 patient each in the Synergo group and fever and
20 urgency, general weakness and the false passage which
21 is only again noted in one patient, so that it is not
22 significantly different, and reduced bladder capacity
23 was also not very oftenly seen, but I will also
24 address it in one of the next slides.

25 With regard to pain, we included all forms

1 of pain being bladder spasms, intolerance to the
2 treatment, pain in general and urethral pain. This
3 only happened actually in a small number of patients
4 resulting in shortening of treatment, 10 out of 425,
5 or skipping the treatment for one week, delay of one
6 week of treatment due to pain which was 7 out of 425.
7 All these patients, all these pain problems, and I
8 recognize that from my own patient population, is
9 actually mild. It's easily manageable with
10 medication, and it is transient.

11 Posterior wall tissue reaction is typically
12 something that is caused by the hyperthermia. It is
13 asymptomatic, and you detect it only at the follow-up
14 cystoscopy. We made a visual scoring system which
15 was mild, moderate or severe, mild meaning some
16 redness, moderate meaning some mucosal damage and
17 severe you also see some necrosis. The severity was
18 not related to symptomatic yes or no. We had, for
19 example, 10 percent of severe reactions but they were
20 still asymptomatic.

21 As I told you, these were resolved with out
22 medical intervention. It is, of course, due to the
23 radio frequency antenna which is in the bladder. It
24 is superficial. You don't see any involvement of the
25 muscle, and it actually gives minor or no residual

1 effect. Sometimes you see some scars like you see
2 after TUR or you'll see some residual hyperemia.

3 Some other adverse events which were not
4 significantly different between the two groups, but
5 might be important to address is reduced bladder
6 capacity. We had two patients in the Synergo group
7 who had reduction up until 250 and 330 MLH. So it's
8 not a shrunken bladder but there is some reduction in
9 bladder capacity which is actually, of course, known
10 to be a problem after any kind of intravesical
11 therapy.

12 We had some patients with urethral stenosis
13 and stricture. Again also our urologists recognized
14 that. That is something that you do see in patients
15 which are treated for superficial bladder cancer
16 because of their multiple catheterizations, multiple
17 TURs and, of course, multiple cystoscopy procedures.
18 The fact that the catheter we use for Synergo is a
19 little bit larger, 20 French, than you normally use,
20 14 French, of course, will identify less significant
21 stenosis in the urethra earlier. And finally dysuria
22 which we found in some of the patients. The majority
23 of those patients did not require any treatment and
24 none of the patients had shortening or delay of one
25 of the treatments due to dysuria.

1 We have some serious adverse events in
2 Synergo group of bronchial bleeding, suspected MI and
3 nephrolithiasis. None of those were related to the
4 treatment, and also in the mitomycin C group, we had
5 one patient with hydronephrosis, one cerebrovascular
6 accident and one patient with leukemia, again not to
7 be considered to be treatment related.

8 So to conclude the 101.1 study, I think
9 with regard to the efficacy, I hope I have shown to
10 you that there is a highly significant reduction in
11 the two-year recurrence rate in Synergo group, which
12 to my opinion also is clinically very relevant for
13 these patients. There are compelling study results
14 even with relatively small sample size, and that the
15 results over time have shown to be durable.

16 Safety, Synergo is absolutely well
17 tolerated and the toxicity is comparable to the
18 literature with the wide ranges that Dr. O'Donnell
19 has shown to you with the wide ranges for
20 intravesical therapy.

21 The first supportive trial is the 102.1
22 study which is Synergo versus BCG. This is a
23 randomized controlled trial, comparing Synergo
24 treatment to BCG immunotherapy, again for the
25 prophylactic treatment of patients with intermediate

1 or high-risk papillary superficial bladder carcinoma.
2 It is anticipated to end around 2013, and we will use
3 it as supportive data. So not to statistically
4 compare the study endpoints but we will try to
5 demonstrate the consistency of the results for the
6 Synergo treatment in this other randomized controlled
7 clinical trial.

8 Primary endpoint was also a comparison of
9 the two-year recurrence rate between the two groups,
10 same as in the 101. Secondary endpoints obviously
11 also again were comparison of the progression rate to
12 state higher than T1 or comparison of metastatic
13 disease. And then an additional endpoint obviously
14 also was the local and systemic adverse events.

15 The treatment was a little bit different
16 from the treatment in the 101. The reason for that
17 was that we use in Europe the -- schedule for BCG
18 meaning 6 initial weekly instillations followed by 3
19 instillations at months 3, 6 and 12, and we came up
20 with the Synergo treatment which was more or less the
21 same, being 6 initial weekly instillations and then 6
22 monthly instillations. As you might remember in 101,
23 we had eight weekly instillations and four monthly
24 instillations, small difference. Follow-up again was
25 every three months for two years, and then endpoint

1 assessment again was biopsy or histologically proven
2 tumor recurrence.

3 The population that we are talking about is
4 104 patients for toxicity and evaluable for efficacy,
5 it's 92 patients, and you can see that there are well
6 balanced between the two treatment groups.

7 The baseline characteristics just like in
8 the 101 trial are actually nicely comparable between
9 the two study arms. You can see here all the
10 potential baseline characteristics that you can think
11 of. They are obviously the same in all kinds of
12 these trials. Here you see the Kaplan-Meier curves
13 for the time to recurrence, and you'll see again, if
14 you look at the approximately two-year endpoint, that
15 it is around 17 percent for Synergo which is very,
16 very close to the results in the 101, and it is 32
17 percent in the BCG arm.

18 Secondary endpoints, we didn't have any
19 progression in tumor stage or grade in Synergo
20 treated patients. Of course, this is an interim
21 analysis and based on a limited follow-up.

22 Here you see the side effects that did
23 significantly differ between the two study arms. In
24 the upper part of the graph, you'll see in red the
25 results or the adverse events that were more common,

1 statistically more common in the BCG arm. As you
2 might expect, things like fever, arthralgia, fatigue
3 and in the lower part of the graph you'll see those
4 results, those side effects that were statistically
5 significantly more common in the hyperthermia group,
6 such as obviously pain during treatment and the
7 posterior wall tissue reaction.

8 And then there were some other side effects
9 that did not differ between the two groups. So they
10 were similar in the BCG and the Synergo groups, and I
11 will not burden you with all these. You can read
12 what the results are.

13 There were some serious adverse events. We
14 had the urethral stricture in the Synergo group. One
15 patient had a contracted bladder in that Synergo
16 group. Professor O'Donnell showed you that that is
17 on average around 5 percent with mitomycin C
18 treatment. We had one in this study, although he was
19 recurrent free. So his cancer obviously was cured.
20 We had one patient with urethral bleeding with
21 withdraw consent, and one patient with dysuria,
22 urinary urgency and fever which was by the way
23 transient and he went on with his therapy.

24 We had two serious adverse events in the
25 BCG treatment arm. One patient with macrohematuria

1 and the clot retention which was treated with a
2 catheter and antibiotics and resolved. One patient
3 with fever, conjunctivitis and urinary tract
4 infection, again resolved after antibiotic treatment.

5 So the safety conclusions for 102 are again
6 that the expected adverse events are similar as in
7 the study of 101.1, like dysuria, hematuria, tissue
8 reaction, urinary tract infection, pain, posterior
9 wall tissue reaction and bladder wall necrosis and
10 the other adverse events are also similar in nature
11 to the Study 101.1. And in all, I think the Synergo
12 treatment currently is very well tolerated.

13 We combined these two studies with regard
14 to some of the results we have. We are talking then
15 about 93 Synergo treated patients at 12 unique sites,
16 3 sites in Study 101, 10 sites in Study 102. So one
17 site is entering patients in both studies. And we
18 see that there was a consistent two-year recurrence
19 rate, consistent results across sites and consistent
20 safety profile.

21 Here you see a summary of the results. You
22 see the two-year recurrence rate in the Synergo group
23 from the 101 study which is 19 percent. It is 17
24 percent in the 102 study. So with a similar patient
25 profile, it is very close to each other. We see a 63

1 percent recurrence rate in the mitomycin C group.
2 That is a little bit higher than you see from the
3 meta-analysis from literature. That's 42 percent but
4 obviously in the meta-analysis also, low risk
5 patients have been included and in the Study 101, we
6 only included intermediate to high-risk patients.

7 If you look at the BCG results, they're
8 very nicely comparable to the literature, 32 percent
9 recurrence rate in our study as compared to 35 in the
10 literature and obviously in the literature, you only
11 treat intermediate and high-risk patients with BCG.
12 So that's why you don't find the difference that we
13 find in the mitomycin C treatment.

14 Then shortly, some of the supportive
15 trials. The European Prophylactic Patient Trial,
16 it's a single-arm trial, it's uncontrolled, and it's
17 for commercial use of the Synergo device. The
18 patient selection, treatment sessions and follow-up
19 examinations are actually similar to the two studies
20 I already discussed. You'll see that it is, for
21 example, similar amount of high-risk patients, that
22 is 58 percent in the European Prophylactic Patient
23 Study, 55 in the 101 and also 55 in the 102, and we
24 are now talking about 186, it says here 68, but I'm
25 sorry for the mistake. It's 186 patients with close

1 to 1600 Synergo treatments.

2 We have a little bit higher recurrence rate
3 in this EPP results. The reason for that is that
4 although we have a similar amount of high-risk
5 patients, we have more patients with highly recurring
6 tumors in these EPP dataset. It is 60 percent in the
7 EPP with highly recurrent tumors as opposed to 36
8 percent in the 101 and 22 in the 102. And highly
9 recurrent is defined as at least three recurrences in
10 the last two years. And you know from one of my
11 previous slides that that was one of the most
12 important prognostic factors. So we come up with a
13 32.2 percent estimated recurrence rate at two years,
14 which is still far better than mitomycin C and at
15 least comparable to BCG treatment.

16 We have a group of patients who have had
17 bladder salvage treatment, as those are patients
18 again extremely high-risk patients, also at least 3
19 recurrences within the last 24 months, but also
20 patients who failed prior BCG treatments, and
21 actually they were candidates for cystectomy. It's a
22 group of 82 patients with over 800 Synergo
23 treatments, and we will present this only for safety.

24 The next supportive trial is the 101.4.
25 This is an older study from the nineties, controlled,

1 monitored, one arm clinical study. It's ablative
2 indication for the use of patients with transitional
3 cell carcinoma of the bladder. So patients for whom
4 TUR was not possible or was not recommended, for
5 example, very old patients, this is a dataset of 42
6 patients with close to 400 Synergo treatments, and we
7 will use these data for safety.

8 And a similar group of patients is the
9 European Ablation Patients -- which is after this
10 101.4 study. We went on with treating patients with
11 ablative intent. So again ablative indication for
12 the use of the Synergo machine, patients in whom the
13 TUR was not possible or not recommended. That's a
14 relatively large group of patients, 104, with close
15 to 800 Synergo treatments and again these patients
16 are presented for safety.

17 And if you look at these supportive
18 studies, the EPP, the Bladder Salvage group, the
19 101.4 and the European Ablation group of patients, we
20 see that again the expected adverse events are quite
21 similar to the ones reported in the 101.1 and 102.1
22 and also adverse events were similar to these two
23 studies, and that there were no serious adverse
24 events related to the Synergo device.

25 Having said this, I hope I've convinced you

1 or I've shown you that the Synergo treatment is
2 absolutely clinically very effective, not only
3 significantly better but also clinically very
4 relevant, and that the side effects are very
5 manageable, limited and usually self-limiting, and at
6 least in the same range as we see for other
7 intravesical treatments, and I hope that American
8 patients in the future will have some advantage of
9 this treatment. Thank you very much.

10 DR. GROSSMAN: I'm Barton Grossman, and I
11 am a consultant to Medical Enterprises. I have no
12 equity interest in this company.

13 The treatment of intermediate and high-risk
14 bladder cancer continues to be a significant problem
15 in the United States, and personally from my
16 practice, which is about 90 percent related to
17 bladder cancer. Both mitomycin C and BCG are
18 recommended by the American Urologic Association and
19 are commonly used for the treatment of these
20 diseases.

21 BCG is characterized by high initial
22 efficacy, but there is a significant deterioration in
23 the proportion of patients that remain disease free
24 over time. This continues to be a serious problem.

25 Furthermore, BCG has significant local and

1 systemic toxicity, and this is an ongoing problem,
2 and I continue to see patients referred to me with
3 serious systemic toxicity from BCG.

4 There is a need for more effective and less
5 toxic treatment particularly for intermediate and
6 high-risk bladder cancer.

7 There are safety data on over 4,500 Synergo
8 treatment sessions in 506 patients. Similar
9 toxicities have been reported in the pivotal study
10 and across all five supportive clinical studies.

11 The most common toxicity is that of
12 posterior wall tissue reaction and pain due to the
13 hypothermia. Importantly, the posterior wall tissue
14 reaction was found only at surveillance cystoscopy.
15 These patients were completely asymptomatic. If a
16 cystoscopy was not performed, patients wouldn't even
17 know they had these lesions in their bladder. They
18 resolve spontaneously. These lesions have also been
19 seen, non-healing ulcers have also been seen with
20 intravesical mitomycin C without hyperthermia, and I
21 regularly see patients like this in my own practice
22 who receive intravesical mitomycin C, and again
23 they're asymptomatic, self-limited and usually
24 resolve over time without specific therapy.

25 Pain was seen in these patients. The

1 proportion of patients who actually had pain is much
2 greater than the number of sessions which involved
3 pain. Only 4 percent of the Synergo treatments were
4 shortened or skipped due to transient pain during the
5 session, and again it's important to realize that
6 when sessions were stopped, that did not necessarily
7 prohibit future successful treatments with Synergo.

8 The adverse events that have been observed
9 in the Synergo studies commonly occur with other
10 forms of intravesical chemotherapy and intravesical
11 immunotherapy. There very few serious adverse events
12 that were treatment related, and overall the Synergo
13 therapy was well tolerated.

14 In the pivotal trial, Study 101.1,
15 considering the evaluable patients, randomized as
16 treated, there was an 80 percent reduction in the
17 rate of recurrence with Synergo compared to controls
18 with a hazard ratio of 0.23.

19 Again, importantly Synergo treatment was
20 consistently better than mitomycin in all patient
21 analyses.

22 This is the data for the 101.1 and 102.1
23 Synergo treatment arms, and you can see that the two
24 arms is the pivotal trial. This is the 102.1 trial.
25 The two arms did quite well and the curves are, in

1 fact, very, very similar.

2 If we compare that with the BCG, it
3 demonstrates the Synergo appears to be at least as
4 good as BCG and much better than the mitomycin C in
5 the control arm.

6 How does this compare with the overall
7 experience of BCG and mitomycin C because there is
8 considerable literature about that. If you look at
9 BCG in the literature, the results with 102.1 are, in
10 fact, very similar to what you could expect from
11 previously reported studies.

12 Mitomycin overall, as reported, appears
13 somewhat better than the 101.1 arm but that includes
14 low risk patients as well as intermediate and high-
15 risk patients. If you censor the low risk patients
16 and you just look at the literature for the
17 intermediate and high-risk patients, you find out
18 again the results that were attained in the pivotal
19 trial were very similar to what you could expect from
20 previously reported studies.

21 Long-term data has demonstrated the
22 durability of these responses, and you can see both
23 at 5 and 10 years, Synergo is much better than
24 mitomycin. This data is also important because it
25 demonstrates that there was no bias in the initial

1 early reports of efficacy.

2 The regulatory standard for valid
3 scientific evidence includes adequate and controlled
4 investigations, partially controlled studies, studies
5 without matched controls and well documented case
6 histories. Furthermore, significant human experience
7 with a market device is considered valid scientific
8 evidence.

9 We do have valid scientific evidence. The
10 pivotal trial, 101.1, Synergo versus BCG, 102.1 and
11 the Ablation Study, 101.4, were adequate and
12 controlled investigations. The European Prophylactic
13 Patients, the Bladder Salvage Patients and the
14 European Ablation Patients involved significant human
15 experience with a marketed device.

16 The three trials for efficacy involved
17 high-risk patients. More than 52 percent of these
18 patients enrolled in the Synergo arm were high-risk,
19 a very at risk population, and the total number of
20 these studies, these are patients that only received
21 Synergo was 201 patients, a significant group. The
22 safety, listed here, involves a total of 506
23 patients.

24 Synergo treatment was demonstrated in the
25 pivotal trial to be much better than mitomycin C, for

1 the prophylactic treatment of intermediate and high-
2 risk non-muscle invasive bladder cancer. The data
3 suggested that Synergo may be comparable if not
4 better than BCG. Synergo has low, acceptable and
5 predictable toxicity, without the potential life
6 threatening adverse events that have ultimately been
7 reported with BCG but still occur on a regular basis
8 throughout the United States.

9 Patients treated with Synergo are virtually
10 identical to the intermediate and high-risk patients
11 in the United States, and I must say in my patient
12 population.

13 The pivotal trial results are compelling
14 and furthermore, they're consistent across studies.
15 The long-term results show that there was no
16 assessment bias. Synergo therapy fills an important
17 need for treatment of intermediate and high-risk
18 patients in the United States. The pivotal study,
19 101.1 and the supportive data, provide reasonable
20 assurance of safety and effectiveness based on valid
21 scientific evidence.

22 I've had the opportunity of seeing Synergo
23 therapy. It's amazingly easy to give. It is very
24 well tolerated. I've reviewed the data and am very
25 impressed with the results. I hope you also agree

1 that this is an important new therapy which is safe,
2 effective and needed for our patients in the United
3 States.

4 DR. O'DONNELL: Thank you. I'll conclude
5 with a discussion about the post-approval study, and
6 I'd first like to begin by emphasizing that the
7 sponsor is committed to documenting the consistency
8 of the Synergo device in the treatment of the U.S.
9 patient population and update the labeling to
10 accurately reflect its performance and precautions
11 for its use in the U.S.

12 And while we strongly believe that the data
13 that we've presented so far, showing the efficacy and
14 safety, is compelling, we realize that it's important
15 that for the U.S. population to see the device in use
16 and to have data from our own patients. In fact, it
17 would be a poor marketing strategy to attempt even to
18 bring a new therapy into the U.S. without any U.S.
19 experience and try to sell it to the U.S. population.
20 So it's understood that even if the panel were to say
21 we don't need a post-approval study, we would
22 strongly do this as a company. I say we. I would
23 support, the sponsor would support this as a company
24 for bringing this experience into the U.S. to provide
25 exactly the type of information that we heard about

1 in the previous lecture about how it performs in the
2 subgroup of populations, that demonstrate that the
3 toxicity profile is consistent in our group of
4 patient as well.

5 So with that in mind, I'm going to give you
6 just some brief ideas of where we put the post-
7 approval study design with the caveat that the
8 company remains flexible and open to further input
9 from the Panel, from the FDA, to make this the best
10 study to provide the benefit for the patients and for
11 the physicians in the U.S.

12 So we would begin with the major objective
13 being to demonstrate and to substantiate the safety
14 of this Synergo system in the U.S., and for this, we
15 feel that a single-arm study would be the appropriate
16 study group to use. We would use the treatment
17 regime as demonstrated in the pivotal trial, namely
18 eight weekly sessions with four monthly maintenance
19 sessions. With a follow-up program that conforms to
20 the standard of care in the United States,
21 essentially every three months cystoscopy for the
22 first year, to obtain certainly all the safety data
23 but also to record and to provide the data for the
24 recurrence rate including the results from
25 cystoscopy, cytology and biopsies as appropriate.

1 The key eligibility criteria would be those
2 that have been used already in the studies that have
3 been mentioned, namely restricting this to stage Ta
4 and T1, grade 1 through grade 3 bladder cancer,
5 superficial bladder cancers that conform to the
6 intermediate and high-risk categories in the EAU
7 definition. And in all cases, complete eradication
8 of the tumors attempted ahead of time with the
9 transurethral resection.

10 The key exclusion criteria will remain not
11 to treat patients in the low risk group with a
12 single, low grade papillary tumors, not to treat
13 patients with muscle invasive disease, anything above
14 stage T1, and not to treat patients with carcinoma in
15 situ, CIS or Tis.

16 The appropriate endpoints certainly for the
17 safety would include those that have come out through
18 the previous safety studies, to indicate these are
19 the events that we would expect to see with the
20 Synergo treatment including, of course, the posterior
21 wall tissue reaction, pain, dysuria including other
22 urinary tract symptoms as frequency and urgency, the
23 incidence of stenosis and stricture, hematuria, false
24 passage, hypotonic bladders, reduced bladder
25 capacity, bladder contracture, urinary tract

1 infection and bladder wall necrosis, but certainly
2 also to include any other adverse events that would
3 emerge during this kind of post-approval analysis and
4 data collection.

5 I wish to apologize in advance to the FDA
6 that we originally put together a post-approval
7 study. It wasn't clear to us all the details of what
8 a post-approval study would be. We have evolving
9 thoughts on the subject and came to realize that the
10 original idea of putting together a non-inferiority
11 trial based on the set points for the eight different
12 adverse events occurring in the 101 and 102 study was
13 really unrealistic, unworkable, and didn't conform to
14 the spirit of what a post-approval study really is
15 meant to do. So that is not a clinical meaningful
16 post-approval study, and we really don't feel that
17 this is the appropriate kind of study that we should
18 do, and so we withdraw that formal study concept from
19 the field.

20 What we now feel is more appropriate is a
21 representative group from the U.S. populations, that
22 include at least about 120 subjects that would
23 represent a similar amount that you saw on the
24 combined 101 and 102 studies. It would represent
25 about a quarter more patients, about 20, 25 percent

1 more patients from the cumulative experience for the
2 safety and about another 50 percent increase in the
3 number that we already have for efficacy, and would
4 involve about 5 to 10 U.S. sites.

5 The type of analysis would be more
6 descriptive certainly of the adverse events, the
7 adverse event rate per treatment session and per
8 patient, and would provide point estimates at 95
9 percent confidence intervals to be reported to be
10 used to update the labeling so that we have a label
11 that actually reflects and confirms what we've
12 already seen in the European studies.

13 This would include as well a training
14 program. In fact, as you can imagine, we've seen a
15 lot of new technology come through in urology. We
16 began really first with extra corporeal shock wave
17 lithotripsy. It's evolved into microwave
18 hyperthermia for BPH, green light lasers,
19 cryosurgery. You know, urologists are a group that
20 tends to embrace new technology but with that comes
21 the incumbent need to have a training program to make
22 sure that physicians and their staff are properly
23 educated in the use. And so the company feels that
24 this is an important part of this process as well
25 which would include a training program with didactic

1 elements, a written set of format for teaching and
2 videos, a mentorship by physicians that are
3 experienced with the technique, to come and train
4 physicians and technical staff including on-site
5 training and then, of course, an assessment of
6 proficiency.

7 So I want to thank you, the entire Panel
8 and the public and FDA, for being open minded to
9 listen to this, about Synergo, which we feel is
10 really a novel, advanced and very efficient and new
11 breakthrough for bladder cancer. I hope you'll
12 conclude positively with us. Thank you.

13 DR. TALAMINI: I want to thank the sponsor
14 for their presentation and their punctuality. We're
15 right on time.

16 It's now time for the panel to ask
17 questions. For the Panel, please remember that you
18 may also ask the sponsor questions during the Panel
19 deliberations later on today. So if anyone on the
20 Panel has an extensive question for the sponsor, it
21 would be good to ask that now so that the sponsor
22 would have time to prepare an answer for later today.

23 In addition, it would be important to ask
24 clarifying questions at this time regarding the
25 sponsor's presentation.

1 So with that, I'll ask the Panel if there
2 are questions, and please indicate your desire to ask
3 a question by raising your hand. Yeah, Dr. Connor.

4 DR. CONNOR: My, I think my first question
5 would ask if you could describe a little bit more. I
6 understand how we design trials and we put a great
7 deal of thought into designing trials, but then when
8 trials are implemented, the people at the sites don't
9 do exactly what we hope they'd do. So I'd like to
10 understand more on your first trial about how the
11 randomizations just didn't work out the way you
12 thought. I think there were at least five patients
13 who randomizations got switched on and can someone
14 speak to how exactly that happened. I'd just like to
15 understand that better, why it occurred so often.

16 DR. TALAMINI: And I would ask the sponsor
17 in your response, if it appears that there needs to
18 be an extended response, let us know so we might do
19 that this afternoon. We only have 15 minutes for
20 this question and answer session.

21 MS. STEIN: Okay. The randomization errors
22 occurred at the site due to administrative/clerical
23 errors when the envelopes were pulled. The
24 randomization forms came in and the envelopes were
25 pulled. These were only discovered years later when

1 the FDA had asked us to send in the randomization
2 scheme. We pulled the randomization scheme from the
3 files and were reviewing it before we gave it to the
4 FDA, that's when these randomization errors were
5 discovered. So they were really discovered long
6 after the study was completed. The clinical sites
7 were unaware of it. The sponsor was unaware of these
8 randomization errors until much later. We had
9 mentioned that the numbers were, they were switched
10 in pairs. So both study groups ended up with the
11 same number of patients at the end and I think most
12 importantly is that the statistical analysis that was
13 at the end, where we took that into consideration,
14 and we did a worst-case scenario analysis, we still
15 found that there was a statistical significant
16 between the groups.

17 DR. TALAMINI: Dr. Connor, further
18 questions.

19 DR. CONNOR: So I think I'm still not clear
20 how this mixing up, especially since it sounded like
21 it was one to one, where a patient was randomized to
22 Synergo versus the other treatment, were individually
23 switched. It wasn't someone was randomized to
24 something and got the other treatment. It was that
25 this mixing or, you know, mismatching occurred which

1 I don't understand.

2 DR. O'DONNELL: A piece of paper went to
3 the wrong -- got switched between two groups of
4 patients, between the pairs. One thing that Ahava
5 didn't mention though that I think is important is
6 that of the characteristics of the patients that were
7 switched were looked to see if, well, did that result
8 in, you know, shifting of risk groups or shifting of
9 the higher stage tumors or something or one or the
10 other, and they were really indistinguishable.

11 DR. CONNOR: And in particular, I was just
12 wondering about age, since it looked like the age in
13 the controlled population in that trial was higher.
14 I wondered if there was any systemic issue of higher
15 age patients being involved in that mismatching.

16 DR. O'DONNELL: That I don't know.

17 DR. TALAMINI: If you're not certain, we
18 can certainly look into it.

19 DR. O'DONNELL: I think that whatever -- I
20 mean obviously it was small numbers, five patients
21 each. There was nothing that was obvious. It should
22 be noted that 65 is close to the median point where
23 patients present with bladder cancer. So it's
24 relatively an arbitrary thing. We probably picked it
25 because it's kind of a convenient Medicare related,

1 you know, endpoint, and it has no real clinical
2 significant. Mitomycin C doesn't appear to have a
3 different activity level based on age related
4 differences.

5 DR. CONNOR: Okay. Thank you.

6 DR. TALAMINI: Other topics? Dr. Redman.

7 DR. REDMAN: Yeah, just on slide 37 and
8 also 39. Was the 80 percent follow-up at one year a
9 predetermined, a priority, that that would be the
10 interim analysis? Usually accept events occurring
11 for a priority analysis. That's my one question on
12 that. And the other, why 80 -- and if that was, why
13 80 patients at 1 year follow-up instead of a 2 year
14 follow-up which was a primary endpoint which I
15 understand was not a priority.

16 DR. TALAMINI: And again, if the sponsor
17 needs time to process that question, we can, you
18 know, hold it for later this afternoon. Please just
19 indicate if that would be favorable.

20 MS. DEUTSCH: I haven't been formally
21 introduced. I'm Lisa Deutsch, a biostatistician. I
22 took over the analysis of this project after the
23 monitoring started, after 1997. The interim analysis
24 was conducted by the original statistician that
25 designed the study. It was planned that 80 patients,

1 that after one year follow-up but for ethical
2 reasons, the study was stopped earlier, and there was
3 a data monitoring committee that had decided to stop
4 the study and provide an interim analysis and see if
5 there was -- because they actually saw, because there
6 was no blinding, that the Synergo patients had a very
7 much better safety profile than, efficacy profile,
8 survival profile than the other patients, than the
9 mitomycin patients, and I assume that they had
10 decided to provide the interim analysis at that time
11 based on that information.

12 DR. TALAMINI: Dr. Redman.

13 DR. REDMAN: That's the first I've heard
14 that there was a data monitoring committee. So was
15 there an independent data safety monitoring committee
16 set up at the time this trial was done or that made
17 that decision?

18 MS. DEUTSCH: Yes, there was.

19 DR. TALAMINI: Other questions? Dr. Dahm.

20 DR. DAHM: I have a question with regards
21 to the case report forms. From my reading, it's my
22 understanding that there were in the initial phase no
23 case report forms, and that those were instituted
24 secondarily I believe in 1997. The trial was started in
25 1994, and the case reports were instituted in 1997,

1 and then retrospectively completed. I just wanted to
2 understand that a little better especially in the
3 context that I think most of your events occurred in
4 the first two years. So if the trial accrued over
5 three years and then you went back and did the case
6 report forms retrospectively, if that is correct, if
7 my understanding is correct, what the potential
8 impact of that may have been.

9 MS. STEIN: In 1997, that is correct. The
10 CRFs were prepared in 1997, and the information from
11 the patient's hospital records were transcribed from
12 the hospital records onto the CRF. From that point
13 on, all the patient data was prospectively completed
14 on the CRFs as well as on the hospital source
15 documentation. We had mentioned that there was 100
16 percent monitoring according to GCP requirements of
17 all the CRFs versus hospital source documentation
18 including those that were transcribed before 1997 and
19 throughout the study until the end of the two-year
20 follow-up in 2001. The FDA conducted a BIMO audit.
21 That's a bioresearch monitoring audit by a FDA
22 inspector in 2005 where he inspected also the CRFs
23 versus the source documentation. As we presented in
24 this slide, we mentioned that in his report he had
25 written that his finding was that the data on the

1 CRFs were an adequate reflection of the hospital
2 source documentation.

3 DR. WITJES: A small additional remark,
4 what comes out of the results of those CRFs before
5 '97, after '97 and in the supportive trials is
6 similar. So we didn't find any change in reporting
7 of, for example, adverse event because that might be
8 one of the things that you're afraid of.

9 DR. TALAMINI: So let me just ask the Panel
10 if there are other major issues that might require
11 the sponsor to come back to us this afternoon after
12 further analysis. If we could make sure that we get
13 those now. Are there any Panel members that would
14 have such a topic or issue? Dr. Redman.

15 DR. REDMAN: Again, just -- this is
16 regarding supportive data. On the 5 and 10 year
17 follow-up, do you have the number of patients that
18 were followed for that period of time? In other
19 words, I'm sure there are dropouts. I'm sure it
20 wasn't 100 percent. That's just because of the
21 supporting data that you're presenting, and also on
22 slide 71, just for clarification, you're claiming 90
23 patients. I think it was 48 and 42 for an arm in
24 support of that. Were those 90 patients followed up
25 beyond the two years or at two years of follow-up

1 because I didn't count the hash marks, but it looks
2 like a lot of patients never made it -- hadn't made
3 it yet to that two years. So is that 90 patients
4 that were followed at two years or beyond?

5 DR. TALAMINI: And again, if the sponsors
6 require more time to answer, that's fine.

7 MS. DEUTSCH: I assume slide 71 refers to
8 the 102.1 study?

9 DR. REDMAN: Yes.

10 MS. DEUTSCH: So in that case, well, that
11 interim analysis was provided and the data lock was
12 in 2007 and we have more follow-up to date but we're
13 not going to look at that right now. And then the
14 101 study, I can tell you that in the long-term
15 follow-up, in the mitomycin group, after 2 years,
16 there were 11 patients that were still alive without
17 recurrence after 23 of the patients had already had
18 recurrence prior to the 2 years endpoint. So -- and
19 in the Synergo, there were 28 patients that were
20 followed up after 2 years, between the 2 and 10 year
21 and the long-term analysis out of which 6 had already
22 recurred prior to the 2 year endpoint if that answers
23 your question.

24 DR. TALAMINI: Other --

25 DR. O'DONNELL: Yeah, just to clarify a

1 point that one reason the numbers are small in the
2 102 group in the long-term follow-up is that most of
3 the patients of the mitomycin C arm recurred. So
4 there weren't many patients at risk left to follow
5 them for a long period of time, but the numbers were
6 greater in the Synergo group.

7 DR. TALAMINI: Thank you. Other questions
8 from the Panel? Dr. Connor.

9 DR. CONNOR: I think this is a brief
10 clarifying question. It's my slide 39, but in your
11 101.1 study, you do this interim analysis when there
12 were 80 patients and you said it was on the slide
13 that you recalculated the sample size to be 84, and I
14 wanted to clarify, I assume that there wasn't the
15 sample size recalculation but rather by the time, you
16 know, if someone took a dataset with 80 patients,
17 looked at it and said, oh, there's a difference here,
18 let's stop this trial, but by the time you stopped
19 the trial, 4 more patients were enrolled. Is that
20 true versus a conceivable size recalculation?

21 MS. DEUTSCH: I'll answer that. The
22 interim analysis was called for when 80 patients had
23 completed 1 year follow-up but the data safety
24 monitoring committee alerted the company, the
25 sponsor, to provide interim data and the study

1 statistician at the time provided the analysis on 39
2 patients who had completed a 2 year follow-up.

3 DR. CONNOR: And how many were in the trial
4 at that point?

5 MS. DEUTSCH: There were 64 patients --

6 DR. CONNOR: Okay.

7 MS. DEUTSCH: -- in the trial at that
8 point, and he recalculated the sample size and
9 extended the study so that they would complete 84
10 patients with 2 year follow-up.

11 DR. CONNOR: Okay.

12 DR. TALAMINI: Dr. Donatucci.

13 DR. DONATUCCI: Yes. I would just like to
14 understand a little bit better the -- just looking at
15 the numbers, three centers from '94 to '99, if I
16 understand it, and we have a breakdown in the number
17 of patients per center, but what I don't see and
18 don't understand is how many patients were pre-1997
19 when the case report forms were generated and how
20 many were post-1997?

21 MS. STEIN: Approximately two-thirds of the
22 patients were already involved in the study by 1997.
23 That means that they were enrolled but they did not
24 necessarily have two-year follow-up by that time.

25 DR. TALAMINI: Okay. I think seeing no

1 further questions, our schedule now provides for a
2 15-minute break.

3 I remind Panel members that there should be
4 no discussion of the PMA during the break, amongst
5 themselves, with the sponsor or with the public.

6 It is now 10:15. So we will resume
7 promptly at 10:30. Thanks.

8 (Off the record.)

9 (On the record.)

10 DR. TALAMINI: It is now 10:30 by my watch.
11 I'd like to call the meeting back to order.

12 We will now hear FDA's presentation. The
13 first FDA presenter is Mr. John Baxley, the review
14 team leader for this PMA. Mr. Baxley.

15 MR. BAXLEY: Good morning. I'd like to
16 thank the Panel for your time and effort in reviewing
17 this PMA. My name is John Baxley, and I would like
18 to present the FDA review of the Medical Enterprises
19 Synergo SB-TS 101.1 Device and mitomycin C.

20 The Synergo SB-TS 101.1 Device and
21 mitomycin C, collectively referred to as the Synergo
22 system, is a device/drug combination product. The
23 lead review was conducted by the Center for Devices
24 and Radiological Health.

25 As presented earlier by the sponsor, the