

DEPARTMENT OF HEALTH AND HUMAN SERVICES PUBLIC HEALTH SERVICE FOOD AND DRUG ADMINISTRATION MILK LABORATORY EVALUATION FORM	LABORATORY	
	LOCATION	LAB #
	DATE	X = DEVIATION U = UNDETERMINED O = NOT USED NA = NOT APPLICABLE

PASTEURIZED MILK CONTAINERS
 [Unless otherwise stated all tolerances are ±5%]

1. Laboratory Requirements

- a. Record time and date when samples received _____
- b. Record time and date when samples examined _____

RINSE METHOD

APPARATUS

- 2. See Cultural Procedures, items 1-23** _____
- 3. Adding Rinse Solution to Containers** _____
 - a. Sterile hypodermic syringes (capacity 20 or 100 mL) and needles _____
 - b. Or, sterile pipets _____
 - c. Or, sterile automatic syringe _____
 - d. Or, sterile graduate cylinder _____
 - e. Or, pre-dispensed dilution bottles or tubes with rinse solution, volumes checked as in CP item 28i _____

MATERIALS

- 4. See Cultural Procedures, items 24-32** _____
- 5. Rinse Solutions, see CP items 27 i-k** _____
- 6. Ethyl Alcohol, 70%** _____
- 7. Plastic Tape** _____

PROCEDURE

- 8. Identify Plates (see SPC, item 4)** _____
- 9. Rinse Solution Volumes for Collection of Surface Rinse Samples** _____
 - a. 20 mL for ½ pints (236 mL), pints (473 mL), and quarts (946 mL) _____
 - b. 50 mL for ½ gallons (1892 mL) _____
 - c. 100 mL for gallons (3784 mL) or larger _____
 - d. Containers less than 100 mL use swab method, see items 19-34 _____
- 10. Collection of Surface Rinse Samples** _____
 - a. Firm walled paper containers, sealed on line _____
 - 1. Swab top of containers with 70% alcohol _____
 - 2. Add required amount of rinse solution to each container by injection and seal puncture with plastic tape _____
 - 3. Vigorously shake container length-wise 10 times, holding container horizontally _____
 - 4. Each shake a complete back and forth movement of about 8 in _____
 - 5. Turn container 90° and repeat horizontal shaking treatment _____
 - 6. Turn container 90° twice more and repeat horizontal shaking _____
 - 7. Grasp container and swirl 20 times in a small circle while upright (top up) _____
 - 8. Invert (top down) and repeat swirling of container 20 times _____
 - 9. Stand upright and allow to drain for 1 min _____

- b. Plastic capped containers (submitted with caps) _____
 - 1. Swab top of container with 70% alcohol when appropriate _____
 - 2. Add required amount of rinse solution by aseptically removing cap, pouring in solution without touching the top and replace cap _____
 - 3. Complete rinse procedure as described in 10a 3/9 above _____
- c. Flexible-walled containers/bags _____
 - 1. Add 100 mL aseptically by swabbing an area of tube adjacent to liner with 70% alcohol; introduce rinse by syringe and seal puncture with plastic tape _____
 - 2. Place container/bag on smooth, clean, firm horizontal surface as flat as its construction permits _____
 - 3. With hands or roller, move rinse solution back and forth 10 times, contacting all surfaces completely _____
 - 4. Lift liner and hang with "fill tube" down to permit rinse solution to collect for one minute _____
 - 5. Transfer rinse solution to sterile container by cutting "fill tube" with sterile scissors _____

11. Sample Measurements _____

- a. As described in SPC, item 6 and 7 except _____
- b. For Residual Bacterial Count (RBC), plate the following amounts per container _____
 - 1. ½ pint: 5 mL in two plates _____
 - 2. Pint and larger: 2 mL in a single plate _____
- c. For Residual Coliform Count (RCC), distribute 10 mL of remaining rinse solution among three plates _____

12. Plating (See SPC, item 13) _____

13. Controls (See SPC, item 14) _____

14. Incubation (See SPC, item 15) _____

15. Counting Colonies (See SPC, items 16-18) _____

16. Confirmation Test (See SPC, item 18e) _____

17. Recording Counts _____

- a. Count obtained from RBC plate recorded as colonies counted _____
- b. If no colonies on RBC plate, record as 0 _____
- c. Count obtained from RCC plate recorded as colonies counted _____
- d. If no colonies on RCC plate, record as 0 _____
- e. Values are recorded as number of colonies per container _____

REPORTS

18. Reporting Counts _____

- a. Report computed count as Residual Bacterial Count (RBC)/ specified container capacity in mL _____
 - 1. Containers rinsed with 20 mL _____
 - a. 5 mL plated for RBC, multiply colonies by 4 _____
 - b. 2 mL plated for RBC, multiply colonies by 10 _____
 - 2. Containers rinsed with 50 mL _____
 - a. 2 mL plated for RBC, multiply colonies by 25 _____

PASTEURIZED MILK CONTAINERS
[Unless otherwise stated all tolerances are ±5%]

- 3. Containers rinsed with 100 mL _____
 - a. 2 mL plated for RBC, multiply colonies by 50 _____
- b. Report computed coliform count as Residual Coliform Count/container _____
 - 1. Containers rinsed with 20 mL _____
 - a. 10 mL plated for RCC, multiply colonies by 2 _____
 - 2. Containers rinsed with 50 mL _____
 - a. 10 mL plated for RCC, multiply colonies by 5 _____
 - 3. Containers rinsed with 100 mL _____
 - a. 10 mL plated for RCC, multiply colonies by 10 _____
- c. If no colonies appear on plates, report as less than n/size of specified container in mL, substituting for n the number that would be reported if 1 colony had been counted from the volume of rinse solution plated and multiply by appropriate factor _____

SWAB METHOD
APPARATUS

- 19. See Cultural Procedures, items 1-23 _____
- 20. Screw-Capped Containers _____
 - a. 7 to 10 cm long to contain 5 mL solution for non-soluble swabs _____
 - b. Containing 4.5 mL solution for alginate swabs _____
 - c. Sterile _____
- 21. Swabs _____
 - a. Non-toxic _____
 - 1. Test each lot by swirling several swabs in 5 mL of sterile dilution buffer, test using *Bacillus stearothermophilus* type assay _____
 - 2. Maintain records _____
 - b. Packaged in convenient protective containers and sterilized _____
 - c. Cotton, non-absorbent (firmly twisted to about 5 mm diameter by 2 cm long over one end of wood applicator stick 12-15 cm long) _____
 - d. Or, calcium alginate fibers _____
 - e. Or, dacron or rayon fibers _____

MATERIALS

- 22. See items 4 and 5 _____
- 23. Sodium Hexa-metaphosphate Solution, 10% (if alginate swabs used), sterile _____
- 24. Shaking Machine, Optional _____

PROCEDURE

- 25. Identify Plates (See SPC, item 4) _____
- 26. Collection of Swab Samples _____

- a. Aseptically remove sterile swab from container _____
- b. Open vial of solution, wet swab and press out excess solution _____
- c. Holding swab at 30° angle to surface, rub over 50 sq. cm area (calculate or use sterile template) three times, reversing direction between successive strokes _____
- d. Rinse swab in solution and press out excess _____
- e. Swab four additional 50 sq. cm areas _____
- f. After fifth area has been swabbed, position swab head in vial and break stick, leaving swab head in vial _____

- 27. Sample Measurement _____
 - a. As described in SPC item 6 and 7 except _____
 - 1. If calcium alginate wool is used, add 0.5 mL of sterile Na Hexa-metaphosphate solution to 4.5 mL rinse solution in vial and shake until dissolved _____
 - 2. For all other fibers _____
 - a. Shake swab container 50 times _____
 - b. Each shake a complete back and forth movement of 15 cm _____
 - c. Strike palm of hand at end of each cycle _____
 - d. Complete shaking in about 10 sec _____
 - e. Or, use shaking machine to give equivalent mixing/disintegration, approximately 30 sec _____
 - b. Transfer 1 mL for Residual Bacterial Count (RBC) _____
 - c. Transfer 3 mL for Residual Coliform Count (RCC) _____

- 28. Plating (See SPC, item 13) _____
- 29. Controls (See SPC, item 14) _____
- 30. Incubation (See SPC, item 15) _____
- 31. Counting Colonies (See SPC, items 16-18) _____
- 32. Confirmation Test (See SPC, item 18e) _____
- 33. Recording Counts _____
 - a. Count obtained from RBC plate recorded as colonies counted _____
 - b. If no colonies on RBC plate, record as 0 _____
 - c. Count obtained from RCC plate recorded as colonies counted _____
 - d. If no colonies on RCC plate, record as 0 _____

REPORTS

- 34. Reporting Counts _____
 - a. Report RBC count as n/50 sq cm by dividing count in 33a ... _____
 - b. If no colonies on RBC plate, report as < 1/50 sq cm _____
 - c. Report RCC count as recorded in 33c _____
 - d. If no colonies on RCC plate, report as < 1/50 sq cm _____