

APPENDIX H
OHIO WATER MICROBIOLOGY LABORATORY
MEI AGAR PREPARATION

BASAL MEDIUM (mE AGAR)

*Don't need to mix—use manufactured mE agar

| | |
|-----------------|----------|
| Peptone | 10.0 g/L |
| Sodium chloride | 15.0 g/L |
| Yeast extract | 30.0 g/L |
| Esculin | 1.0 g/L |
| Actidione | 0.05 g/L |
| Sodium azide | 0.15 g/L |
| Agar | 15.0 g/L |

- Add 71.2 g of pre-mixed dehydrated basal medium plus 0.75 g of indoxyl- β -D glucoside to 1 L of reagent water in a Wheaton bottle.
- Heat to dissolve (to boiling).
- Autoclave the Wheaton bottle and 10 empty dilution bottles with caps for 15 minutes.
- Mix 0.24 g nalidixic acid in 5 mL sterile reagent water (add a few drops of 0.1N NaOH to dissolve if needed). Add this solution to the tempered medium.
- Add 0.02 g triphenyl tetrazolium chloride to the medium and mix.
- Pour medium into the sterilized 100-mL dilution bottles.
- Store at 4°C for up to 6 months.

PREPARATION OF AGAR PLATES

- Melt the basal medium using a beaker with water on a hot plate. **DO NOT AUTOCLAVE!**
- Pour the plates after the agar is tempered (50-60°C).
- Store the plates at 4°C for up to 2 weeks in a tightly sealed container.