F



Laboratory Component: Vitamin A, E and Carotenoids

Survey Years: 2003 to 2004

SAS Export File: L45VIT_C.XPT



First Published: October 2007 Last Revised: N/A

NHANES 2003–2004 Data Documentation

Laboratory Assessment: Lab 45 -	- Vitamin A (retinol, retinyl palmitate, retinyl
	stearate), Vitamin E (α -tocopherol, δ -tocopherol,
	γ-tocopherol), Carotenoids (lutein, zeaxanthin, α-
	cryptoxanthin, β -cryptoxanthin, trans-lycopene,
	cis-lycopene, α -carotene, trans- β -carotene, cis- β -carotene; phytoene and phytofluene)
First Published: October 2007	Last Revised: N/A

This file is updated to add vitamins A, E, and Carotenoids data and regression analysis comparing data from 2001-2002 to 2003-2004.

Component The objectives of this component are: 1) to provide data for monitoring Description secular trends in measures of nutritional status in the U.S. population; to evaluate the effect of people's habits and behaviors such as physical activity and the use of alcohol, tobacco, and dietary supplements on people's nutritional status; and 3) to evaluate the effect of changes in nutrition and public health policies including welfare reform legislation, food fortification policy, and child nutrition programs on the nutritional status of the U.S. population. These data will be used to estimate deficiencies and toxicities of specific nutrients in the population and subgroups, to provide population reference data, and to estimate the contribution of diet, supplements, and other factors to serum levels of nutrients. Data will be used for research to further define nutrient requirements as well as optimal levels for disease prevention and health promotion. Eligible Participants aged 6 years and older who do not meet any of the exclusion Sample criteria are eligible. **Description of** Vitamin A, Vitamin E, and Carotenoids

Laboratory Methodology Serum concentrations of vitamin A (retinol, retinyl palmitate, retinyl stearate), vitamin E (α -, δ -, and γ -tocopherols), and twelve carotenoids (lutein, zeaxanthin, α -cryptoxanthin, β -cryptoxanthin, *trans*-lycopene, cis-lycopene, α -carotene, trans- β -carotene, *cis*- β -carotene, *cis*-combined lutein/zeaxanthin,phytoene and phytofluene) are measured using high performance liquid chromatography with multiwavelength photodiode-array absorbance detection. A small volume (150 µL) of serum/plasma is mixed with an equal volume of buffer, and then mixed with 2 volumes of ethanol containing the internal standard (tocol). The analytes are extracted from the aqueous phase into hexane. The combined hexane extracts are dried under vacuum. The extract is redissolved in ethyl acetate and diluted in mobile phase. An aliquot is injected onto a C18 reversed phase column and eluted isocratically. The analytes all possess absorbance and/or fluorescence which are proportional to their concentration in solution; therefore these properties are used for quantitative analysis. The mode of detection is chosen to provide the highest sensitivity and selectivity. Carotenoids are measured by absorbance at 450 nm. Retinol, retinyl esters, phytoene and phytofluene are measured by UV absorbance near their absorption maxima of 325 nm, 280 nm and 340 nm. Tocopherols have absorption maxima between 292 and 300 nm. Chromatograms are recorded using a computer data system. Analytes are quantified by external standard quantitation using standards to calculate response factors based on the peak area of the analyte. The quantities of analytes are corrected for recovery post-run based upon tocol as an internal standard.

Laboratory testing changed from 2001-2002 and 2003-2004. The National Center for Environmental Health, Centers for Disease & Prevention performed the testing in 2001-2002. Craft Technologies performed the testing in 2003-2004. Craft Technologies measured more analytes in 2003-2004 than the National Center for Environmental Health, Centers for Disease & Prevention measured in 2001-2002. See the codebooks on the NHANES website for the different tests between the two laboratories.

The age group in which these measurements were made changed from 2001-2002 to 2003-2004. Serum concentrations of vitamin A, vitamin E and carotenoids were tested on participants 3 years and older in 2001-2002 and on participants 6 years and older in 2003-2004.

The variable names and/or label descriptors changed from 2001-2002 to 2003-2004 for vitamin E and lutein/zeaxanthin. Also, additional analytes were tested in 2003-2004. The following pairs of test names or codes are used synonymously; vitamin E and α -tocopherol; and vitamin A and retinol. The variable names for vitamin E/ α -tocopherol changed from LBXVIE/LBDVIESI in 2001-2002 to LBXATC/LBDATCSI in 2003-2004.

The variable name and label descriptor for lutein (LBXLUT/LBXLUTSI) in the 2002 and earlier data is equivalent to the variable name "combined lutein/zeaxanthin" (LBXLUZ/LBXLUZSI) in the 2003-2004 data. This is because the lutein peak in the CDC method included both lutein and zeaxanthin.

	Serum concentrations of α -cryptoxanthin (LBXACY/LBDACYSI), δ - tocopherol (LBXDTC/LBDDTCSI), cis-lycopene (LBXCLC/LBDCLCSI), total lycopene (LBXLCC/LBDLCCSI), lutein (LBXLUT/LBDLUTSI), phytofluene (LBXPHF/LBDPHFSI), phytoene (LBXPHE/LBDPHESI), and zeaxanthin (LBXZEA/LBDZEASI) were measured in 2003-2004, but not in 2001-2002.
	A detailed description of the laboratory method used can be found on the NHANES website.
Laboratory Quality Control and Monitoring	The NHANES quality assurance and quality control (QA/QC) protocols meet the 1988 Clinical Laboratory Improvement Act mandates. Detailed QA/QC instructions are discussed in the NHANES Laboratory/Medical Technologists Procedures Manual (LPM). Read the LABDOC file for detailed QA/QC protocols.
	A detailed description of the quality assurance and quality control procedures can be found on the NHANES website.
Data Processing	Data Processing and Editing Serum specimens were processed, stored, and shipped to CTI in Wilson, NC for analysis.
and Editing	Detailed specimen collection and processing instructions are discussed in the NHANES LPM. Vials were stored under appropriate frozen (–20°C) conditions until they were shipped to CTI.
	Twenty-one derived variables were created in this data file. The formula for its derivatization is as follows:
	The α -tocopherol (Vitamin E) results in ug/dL were converted into umol/L by multiplying by 0.02322.
	The α -carotene results in ug/dL were converted into umol/L by multiplying by 0.01863 .
	The α -cryptoxanthin results in ug/dL were converted into umol/L by multiplying by 0.01810.
	The trans-β-carotene results in ug/dL were converted into umol/L by multiplying by 0.01863.
	The total β -carotene results in ug/dL were converted into umol/L by multiplying by 0.01863 .
	The cis- β carotene results in ug/dL were converted into umol/L by

multiplying by 0.01863.

The cis-lycopene results in ug/dL were converted into umol/L by multiplying by 0.01863.

The cis-lutein/zeaxanthin results in ug/dL were converted into umol/L by multiplying by 0.01758.

The β -cryptoxanthin results in ug/dL were converted into umol/L by multiplying by 0.01810.

The $\delta\text{-tocopherol}$ results in ug/dL were converted into umol/L by multiplying by 0.02485

The γ -tocopherol results in ug/dL were converted into umol/L by multiplying by 0.02402

The total lycopene results in ug/dL were converted into umol/L by multiplying by 0.01863

The lutein results in ug/dL were converted into umol/L by multiplying by 0.01758.

The combined lutein/zeaxanthin results in ug/dL were converted into umol/L by multiplying by 0.01758.

The trans-lycopene results in ug/dL were converted into umol/L by multiplying by 0.01863.

The phytofluene results in ug/dL were converted into umol/L by multiplying by 0.01843.

The phytoene results in ug/dL were converted into umol/L by multiplying by 0.01837.

The retinyl palmitate results in ug/dL were converted into umol/L by multiplying by 0.03491.

The retinyl stearate results in ug/dL were converted into umol/L by multiplying by 0.03491.

The vitamin A (retinol) results in ug/dL were converted into umol/L by multiplying by 0.03491.

The zeaxanthin results in ug/dL were converted into umol/L by multiplying by 0.01758.

Detailed instructions on specimen collection and processing can be found on the NHANES website. Analytic Notes The analysis of NHANES 2003–2004 laboratory data must be conducted with the key survey design and basic demographic variables. The NHANES 2003–2004 Household Questionnaire Data Files contain demographic data, health indicators, and other related information collected during household interviews. The Household Questionnaire Data Files also contain all survey design variables and sample weights required to analyze these data. The Phlebotomy Examination file includes auxiliary information on duration of fasting, the time of day of the venipuncture, and the conditions precluding venipuncture. The Household Questionnaire and Phlebotomy Exam files may be linked to the laboratory data file using the unique survey participant identifier SEQN.

Regression Analysis for 2001-2002 and 2003-2004 for Vitamins A,

E, and Carotenoids:

In 2001-2002, Vitamins A, E and Carotenoids were analyzed using a HPLC method performed at CDC/NCEH [See Lab 6 in NHANES 2001-2002 for method details]. In 2003-2004, Vitamins A, E and Carotenoids were analyzed using a comparable HPLC method at Craft Technologies, Inc. (CTI). Crossover studies between CDC/NCEH and CTI were done in early 2003 and late 2004 and differences existed between the methods for some analytes. The following Deming regression analyses were performed to compare data from 2001-2002 to 2003-2004:

Early 2003 Deming Regression Analysis for A/E/Carotenoids Y(CDC)=X(CTI) in ug/dL :					
<u>Test</u>	<u>n</u>	<u>Y(CDC) = X (CTI)</u>	<u>r</u> ²		
Alpha-Carotene	100	y = 0.9734x + 0.1847	0.9944		
Trans-Beta-Carotene	100	y = 1.0764x – 0.0248	0.9971		
Cis-Beta-Carotene	100	y = 0.8017x - 0.0372	0.9685		
Beta-Cryptoxanthin	100	y = 1.1359x – 0.2116	0.9945		
Gamma-Tocopherol	100	y = 1.0346x – 8.5749	0.9818		
Lutein/Zeaxanthin	99	y = 0.9321x + 0.1231	0.9455		
Trans-Lycopene	100	y = 1.0719x – 0.1817	0.9947		
Retinyl Palmitate	100	y = 1.0245x + 0.5039	0.1988		
Retinyl Stearate	97	y = 0.803x - 0.0048	0.9807		
Vitamin A	100	y = 0.9406x – 0.1199	0.9887		
Vitamin E	100	y = 0.9541x – 22.231	0.9941		

	n Analy	sis for A/E/Carotenoids Y(CDC)	=X(CTI) in
ug/dL:			
Test	n	Y(CDC) = X (CTI)	r ²
Alpha-Carotene	99	y = 0.9989x + 0.2204	0.9768
Trans-Beta-Carotene	98	y = 1.0275x + 1.2976	0.9960
Cis-Beta-Carotene	99	y = 0.7154x + 0.1344	0.9826
Beta-Cryptoxanthin	98	y = 1.0824x - 0.056	0.9696
Gamma-Tocopherol	98	y = 1.0306x - 5.523	0.9922
Lutein/Zeaxanthin	99	y = 1.1082x - 1.602	0.9141
Trans-Lycopene	98	y = 1.0971x + 0.0645	0.9737
Retinyl Palmitate	98	y = 0.7989x + 1.2106	0.8845
Retinyl Stearate	98	y = 0.7681x + 0.1483	0.9621
Vitamin A	99	y = 0.9806x + 0.2901	0.9934
Vitamin E	98	y = 0.93x + 22.614	0.9905

Combined Deming Regression Analysis for A/E/Carotenoids Y(CDC)=X(CTI) in ug/dL					
Test	<u>n</u>	$\underline{Y(CDC)} = X(CTI)$	<u>r</u> ²		
Alpha-Carotene	199	y = 0.9837x + 0.2137	0.9877		
Trans-Beta-Carotene	198	y = 1.0339x + 0.8812	0.9960		
Cis-Beta-Carotene	199	y = 0.7232x + 0.0879	0.9815		
Beta-Cryptoxanthin	198	y = 1.0958x - 0.0542	0.9793		
Gamma-Tocopherol	198	y = 1.029x - 6.287	0.9893		
Lutein/Zeaxanthin	198	y = 1.0641x - 1.2402	0.9225		
Trans-Lycopene	198	y = 1.0852x - 0.0681	0.9849		
Retinyl Palmitate	198	y = 0.8227x + 0.99	0.7603		
Retinyl Stearate	195	y = 0.7881x + 0.0711	0.9675		
Vitamin A	199	y = 0.9805x - 1.0069	0.9871		
Vitamin E	198	y = 0.9397x + 1.649	0.9914		

The combined regression was performed using crossover data from early 2003 and late 2004. The data user may wish to apply regression formulas to compare results from 2001-2002 and 2003-2004.

References N/A

Locator Fields

Title: Serum concentrations of Vitamin A, E and Carotenoids

Contact Number: 1-866-441-NCHS

Years of Content: 2003–2004

First Published: October 2007

Revised: N/A

Access Constraints: None

Use Constraints: None

Geographic Coverage: National

Subject: fat-soluble micronutrients

Record Source: NHANES 2003–2004

Survey Methodology: NHANES 2003–2004 is a stratified multistage probability sample of the civilian non-institutionalized population of the U.S.

Medium: NHANES Web site; SAS transport files

National Health and Nutrition Examination Survey Codebook for Data Production (2003-2004)

Vitamin A, E and Carotenoids (L45VIT_C) Person Level Data

First Published: October 2007 Last Revised: N/A



SEQN	Target			
	B(6 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Respondent sequence number			
English Text: Respondent sequence number.				
English Instructions:				

LBXATC		Ta	rget		
		B(6 Yrs. to	o 150 Yrs.)		
Hard Edit	s	SAS Label			
		a-Tocophe	erol(ug/dL)		
English Text: a-Toco	opherol(ug/dL)	/dL)			
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
39 to 7939	Range of Values	7266	7266		
7	Fill Value of Limit of Detection	of Limit of Detection 2 7268			
•	Missing	714	7982		

LBDATCSI		Target			
	•		B(6 Yrs. to	o 150 Yrs.)	
Hard Edits	5	SAS Label			
			a-Tocophe	rol(umol/L)	
English Text: a-Toco	pherol(umo	nol/L)			
English Instructions	:				
Code or Value	I	Description Count Cumulative Skip to Item			
0.9056 to 184.3436	Ra	nge of Values	7266	7266	
0.1625	Fill Value	of Limit of Detection	2	7268	
•		Missing	714	7982	

LBXALC		Ta	arget		
			B(6 Yrs.	to 150 Yrs.)	
Hard Edit	S		SAS	Label	
			a-Carote	ene(ug/dL)	
English Text: a-Caro	otene(ug/dL)				
English Instructions	:				
Code or Value	D	Description Count Cumulative Skip to Item			Skip to Item
0.31 to 123.35	Rar	nge of Values	7198	7198	
0.21	Fill Value	of Limit of Detection	70	7268	
		Missing	714	7982	

LBDALCSI		Ta	rget			
		B(6 Yrs. to	o 150 Yrs.)			
Hard Edit	s	SAS Label				
		a-Caroten	e(umol/L)			
English Text: a-Caro	tene(umol/L)	/L)				
English Instructions	:					
Code or Value	Description	Description Count Cumulative Skip to Item				
0.0058 to 2.298	Range of Values	7198	7198			
0.0039	Fill Value of Limit of Detection	of Limit of Detection 70 7268				
•	Missing	714	7982			

LBXACY	,	Target			
		B(6 Yrs.	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		a-Cryptoxa	anthin(ug/dL)		
English Text: a-Cryp	otoxanthin(ug/dL)	(ug/dL)			
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.22 to 14.21	Range of Values	7250	7250		
0.14	Fill Value of Limit of Detection	of Limit of Detection 18 7268			
•	Missing	714	7982		

LBDACYSI		Ta	rget		
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	s	SAS Label			
		a-Cryptoxan	thin(umol/L)		
English Text: a-Cryp	toxanthin(umol/L)	umol/L)			
English Instructions	:				
Code or Value	Description	Description Count Cumulative Skip to Item			
0.004 to 0.2572	Range of Values	7250	7250		
0.0025	Fill Value of Limit of Detection	of Limit of Detection 18 7268			
•	Missing	714	7982		

LBXBEC		Target				
		B(6 Yrs. to 150 Yrs.)				
Hard Edits		SAS Label				
		trans-b-carotene(ug/dL)				
English Text: trans-b-ca	arotene(ug/dL)					
English Instructions:						
Code or Value	Description	Count	Cumulative	Skip to Item		
0.43 to 388.19	Range of Values	7268	7268			
•	Missing	714	7982			

LBDBECSI		Target			
			B(6 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label			
			trans-b-caro	tene(umol/L)	
English Text: trans-b-c	carotene(ui	nol/L)			
English Instructions:					
Code or Value	I	Description	Count	Cumulative	Skip to Item
0.008 to 7.232	Ra	nge of Values	7268	7268	
		Missing	714	7982	

LBXBCC		Target			
		B(6 Yrs. to 150 Yrs.)			
Hard Edits			SAS	Label	
		total b-Carotene(ug/dL)			
English Text: total b-Ca	English Text: total b-Carotene(ug/dL)				
English Instructions:					
Code or Value	D	escription	Count	Cumulative	Skip to Item
0.43 to 411.48	Ran	nge of Values 7268 7268			
		Missing	714	7982	

LBDBCCSI		Target				
		B(6 Yrs.	to 150 Yrs.)			
Hard Edits SAS Label						
		total b-Caro	otene(umol/L)			
English Text: total b-Ca	rotene(umol/L)					
English Instructions:						
Code or Value	Description	Count	Cumulative	Skip to Item		
0.008 to 7.6659	Range of Values	nge of Values 7268 7268				
· .	Missing	714	7982			

LBXCBC	۱	Target				
LDACD		B(6 Yrs.	to 150 Yrs.)			
Hard Edits SAS Label						
		cis-b-car	otene(ug/dL)			
English Text: cis-b-c	English Text: cis-b-carotene(ug/dL)					
English Instructions	5:					
Code or Value	Description	Count	Cumulative	Skip to Item		
0.3 to 23.28	Range of Values	6870	6870			
0.21	Fill Value of Limit of Detection	of Limit of Detection 398 7268				
	Missing	714	7982			

LBDCBCS	T	Target				
		B(6 Yrs. to 150 Yrs.)				
Hard Edits SAS Label						
		cis-b-carotene(umol/L)				
English Text: cis-b-c	ext: cis-b-carotene(umol/L)					
English Instructions	:					
Code or Value	Description	Count	Cumulative	Skip to Item		
0.0056 to 0.4337	Range of Values	6870	6870			
0.0039	Fill Value of Limit of Detection	of Limit of Detection 398 7268				
•	Missing	714	7982			

LBXCLC	1	Target			
		B(6 Yrs.	to 150 Yrs.)		
Hard Edits SAS Label			Label		
		cis-Lycop	pene(ug/dL)		
English Text: cis-Ly	copene(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.4 to 93.63	Range of Values	7265	7265		
0.21	Fill Value of Limit of Detection	of Limit of Detection 3 7268			
•	Missing	714	7982		

LBDCLCS	T	Target			
		B(6 Yrs. to	o 150 Yrs.)		
Hard Edits SAS Label					
		cis-Lycopene(umol/L)			
English Text: cis-Ly	ext: cis-Lycopene(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0075 to 1.7443	Range of Values	7265	7265		
0.0039	Fill Value of Limit of Detection	of Limit of Detection 3 7268			
•	Missing	714	7982		

LBXCLZ		Target				
			B(6 Yrs.	to 150 Yrs.)		
Hard Edit	s		SAS	5 Label		
			cis- Lutein/Ze	eaxanthin(ug/dL)		
English Text: cis- Lu	cis-Lutein/Zeaxanthin(ug/dL)					
English Instructions	:					
Code or Value	I	Description	Count	Cumulative	Skip to Item	
0.23 to 16.28	Ra	nge of Values	7262	7262		
0.14	Fill Value	of Limit of Detection	6	7268		
		Missing	714	7982		

LBDCLZSI		Target			
			B(6 Yrs. to	o 150 Yrs.)	
Hard Edits SAS Label					
		cis- Lutein/Zeaxanthin(umol/L)			
English Text: cis- Lu	English Text: cis- Lutein/Zeaxanthin(umol/L)				
English Instructions	:				
Code or Value	D	escription	Count	Cumulative	Skip to Item
0.004 to 0.2862	Ran	ge of Values	7262	7262	
0.0025	Fill Value o	f Limit of Detection	6	7268	
•		Missing	714	7982	

LBXCRY		TargetB(6 Yrs. to 150 Yrs.)				
Hard Edits		SAS Label				
		b-crypt	oxanthin(ug/dL)			
English Text: b-cryp	English Text: b-cryptoxanthin(ug/dL)					
English Instructions	:					
Code or Value	Description	Count	Cumulative	Skip to Item		
0.49 to 97.1	Range of Values	7252	7252			
0.14	Fill Value of Limit of Detec	ction 16	7268			
	Missing	714	7982			

LBDCRYS	T	Target			
LDDCKI		B(6 Yrs. t	o 150 Yrs.)		
Hard Edits SAS Label					
		b-cryptoxanthin(umol/L)			
English Text: b-cryp	xt: b-cryptoxanthin(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0089 to 1.7575	Range of Values	7252	7252		
0.0025	Fill Value of Limit of Detection	of Limit of Detection 16 7268			
•	Missing	714	7982		

LBXDTC	1	Target			
		B(6 Yrs.	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		d-Tocoph	erol(ug/dL)		
English Text: d-Toco	opherol(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
3 to 126	Range of Values	5146	5146		
2	Fill Value of Limit of Detection	of Limit of Detection 2122 7268			
•	Missing	714	7982		

LBDDTCS	T.	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	S	SAS	Label		
		d-Tocophe	rol(umol/L)		
English Text: d-Toco	opherol(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0746 to 3.1311	Range of Values	5146	5146		
0.0497	Fill Value of Limit of Detection	2122	7268		
•	Missing	714	7982		

LBXGTC	<u>, </u>	Target			
		B(6 Yrs.	to 150 Yrs.)		
Hard Edit	s	SAS	5 Label		
		g-tocoph	erol(ug/dL)		
English Text: g-toco	pherol(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
15 to 3122	Range of Values	7266	7266		
7	Fill Value of Limit of Detection	of Limit of Detection 2 7268			
•	Missing	714	7982		

LBDGTCS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	S	SAS	Label		
		g-tocophe	rol(umol/L)		
English Text: g-tocop	English Text: g-tocopherol(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.3603 to 74.9904	Range of Values	7266	7266		
0.1681	Fill Value of Limit of Detection	n 2	7268		
•	Missing	714	7982		

LBXLCC		Target			
		B(6 Yrs.	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		total Lyco	pene(ug/dL)		
English Text: total L	ycopene(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.76 to 175.62	Range of Values	7266	7266		
0.21	Fill Value of Limit of Detection	2	7268		
•	Missing	714	7982		

LBDLCCS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	S	SAS	Label		
		total Lycop	ene(umol/L)		
English Text: total L	ycopene(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0142 to 3.2718	Range of Values	7266	7266		
0.0039	Fill Value of Limit of Detectio	of Limit of Detection 2 7268			
•	Missing	714	7982		

LBXLUT		Target			
			B(6 Yrs	. to 150 Yrs.)	
Hard Edit	S		SA	S Label	
			Lute	in(ug/dL)	
English Text: Lutein	(ug/dL)				
English Instructions	:				
Code or Value	I	Description	Count	Cumulative	Skip to Item
0.38 to 88.18	Ra	nge of Values	7264	7264	
0.14	Fill Value	of Limit of Detection	4	7268	
•		Missing	714	7982	

LBDLUTSI		Target			
			B(6 Yrs.	to 150 Yrs.)	
Hard Edits			SAS	S Label	
			Lutei	n(umol/L)	
English Text: Lutein	(umol/L)				
English Instructions	:				
Code or Value	D	Description	Count	Cumulative	Skip to Item
0.0067 to 1.5502	Rar	nge of Values	7264	7264	
0.0025	Fill Value	of Limit of Detection	4	7268	
•		Missing	714	7982	

LBXLUZ		Target			
			B(6 Yrs.	to 150 Yrs.)	
Hard Edits SAS Label					
		С	ombined Lutein	/zeaxanthin (ug/dI	.)
English Text: Combi	ned Lutein/z	n/zeaxanthin (ug/dL)			
English Instructions	:				
Code or Value	D	Description	Count	Cumulative	Skip to Item
0.57 to 113.05	Range of Values		7267	7267	
0.14	Fill Value of Limit of Detection		1	7268	
		Missing	714	7982	

LBDLUZS	T	Target				
		B(6 Yrs. to 150 Yrs.)				
Hard Edits SAS Label						
	C	ombined Lutein/	zeaxanthin (umol/	L)		
English Text: Combi	ned Lutein/zeaxanthin (umol/L)	n/zeaxanthin (umol/L)				
English Instructions	:					
Code or Value	Description	Count	Cumulative	Skip to Item		
0.01 to 1.9874	Range of Values	7267	7267			
0.0025	Fill Value of Limit of Detection	1	7268			
	Missing	714	7982			

LBXLYC		Target			
			B(6 Yrs. 1	to 150 Yrs.)	
Hard Edit	S		SAS	Label	
			trans-lyco	pene(ug/dL)	
English Text: trans-l	ycopene(ug/a	dL)			
English Instructions	:				
Code or Value	E	Description	Count	Cumulative	Skip to Item
0.34 to 81.99	Rar	nge of Values	7265	7265	
0.21	Fill Value	of Limit of Detection	3	7268	
		Missing	714	7982	

LBDLYCS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	5	SAS	Label		
		trans-lycop	ene(umol/L)		
English Text: trans-ly	ycopene(umol/L)	umol/L)			
English Instructions					
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0063 to 1.5275	Range of Values	7265	7265		
0.0039	Fill Value of Limit of Detection	3	7268		
•	Missing	714	7982		

LBXPHF			Target B(6 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label			
	ເວ 				
			Phytoflue	ene(ug/dL)	
English Text: Phytof	English Text: Phytofluene(ug/dL)				
English Instructions	:				
Code or Value	D	Description	Count	Cumulative	Skip to Item
0.35 to 42.08	Rar	nge of Values	7209	7209	
0.21	0.21 Fill Value of Limit of Detection		59	7268	
		Missing	714	7982	

LBDPHFS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	s	SAS	Label		
		Phytofluene(umol/L)			
English Text: Phytof	luene(umol/L)	nol/L)			
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0065 to 0.7755	Range of Values	7209	7209		
0.0039	Fill Value of Limit of Detection	59	7268		
•	Missing	714	7982		

LBXPHE			TargetB(6 Yrs. to 150 Yrs.)		
Hard Edits			SAS Label		
		Phytoene(ug/dL)			
English Text: Phytoe	English Text: Phytoene(ug/dL)				
English Instructions	:				
Code or Value	D	escription	Count	Cumulative	Skip to Item
0.31 to 168.92	Rar	nge of Values	7174	7174	
0.21	Fill Value of Limit of Detection		94	7268	
		Missing	714	7982	

LBDPHES	T	Target			
		B(6 Yrs. to 150 Yrs.)			
Hard Edit	S	SAS	Label		
		Phytoene(umol/L)			
English Text: Phytoe	: Phytoene(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0057 to 3.1031	Range of Values	7174	7174		
0.0039	Fill Value of Limit of Detection	94	7268		
•	Missing	714	7982		

LBXRPL	,	Target			
		B(6 Yrs. 1	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		Retinyl palmitate(ug/dL)			
English Text: Retiny	l palmitate(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.13 to 58.52	Range of Values	7127	7127		
0.09	Fill Value of Limit of Detection	141	7268		
•	Missing	714	7982		

LBDRPLS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	S	SAS	Label		
		Retinyl palmitate(umol/L)			
English Text: Retiny	English Text: Retinyl palmitate(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0045 to 2.0429	Range of Values	7127	7127		
0.0031	Fill Value of Limit of Detection	141	7268		
•	Missing	714	7982		

LBXRST		Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	s	SAS	Label		
		Retinyl stearate(ug/dL)			
English Text: Retiny	l stearate(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.11 to 13.01	Range of Values	4761	4761		
0.08	Fill Value of Limit of Detection	2507	7268		
	Missing	714	7982		

LBDRSTS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	s	SAS Label			
		Retinyl stearate(umol/L)			
English Text: Retiny	inyl stearate(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.0038 to 0.4542	Range of Values	4761	4761		
0.0028	Fill Value of Limit of Detection	2507	7268		
•	Missing	714	7982		

LBXVIA		Target			
		B(6 Yrs. 1	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		Retinol(ug/dL)			
English Text: Retinc	ol(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
6.31 to 164.6	Range of Values	7267	7267		
0.21	Fill Value of Limit of Detection	1	7268		
•	Missing	714	7982		

LBDVIAS	T	Target			
		B(6 Yrs. t	o 150 Yrs.)		
Hard Edit	S	SAS	Label		
		Retinol(umol/L)			
English Text: Retino	n Text: Retinol(umol/L)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.2203 to 5.7462	Range of Values	7267	7267		
0.0073	Fill Value of Limit of Detectio	n 1	7268		
•	Missing	714	7982		

LBXZEA		Target			
		B(6 Yrs. t	to 150 Yrs.)		
Hard Edit	s	SAS	Label		
		Zeaxanthin(ug/dL)			
English Text: Zeaxa	nthin(ug/dL)				
English Instructions	:				
Code or Value	Description	Count	Cumulative	Skip to Item	
0.29 to 25.05	Range of Values	7266	7266		
0.14	0.14 Fill Value of Limit of Detection		7268		
•	Missing	714	7982		

LBDZEASI		Target			
			B(6 Yrs. to 150 Yrs.)		
Hard Edits			SAS Label		
		Zeaxanthin(umol/L)			
English Text: Zeaxa	nthin(umol/L	.)			
English Instructions	:				
Code or Value	Ι	Description	Count	Cumulative	Skip to Item
0.0051 to 0.4404	Rai	nge of Values	7266	7266	
0.0025	Fill Value	of Limit of Detection	2	7268	
		Missing	714	7982	