

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21

**LUMI-CELL[®] ER ASSAY
AGONIST PROTOCOL**

**National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
Toxicological Methods (NICEATM)**

**Developed by:
Xenobiotic Detection Systems, Inc.
1601 E. Geer St., Suite S
Durham, NC 27704
*17 October 2008***

22 **TABLE OF CONTENTS**

23 **LIST OF ACRONYMS AND ABBREVIATIONS vi**

24 **LIST OF FIGURES vii**

25 **LIST OF TABLES viii**

26 **1.0 Purpose 1**

27 **2.0 Sponsor 1**

28 **2.1 Substance Inventory and Distribution Management..... 3**

29 **3.0 Definitions..... 3**

30 **4.0 Testing Facility and Key Personnel 4**

31 4.1 Testing Facility 4

32 4.2 Key Personnel..... 4

33 **5.0 Identification of Test and Control Substances 4**

34 5.1 Test Substances..... 4

35 5.2 Controls..... 4

36 **6.0 Overview of General Procedures For Agonist Testing 5**

37 6.1 Range Finder Testing..... 6

38 6.2 Comprehensive Testing..... 6

39 **7.0 Materials for LUMI-CELL® ER Agonist Testing 7**

40 7.1 BG1Luc4E2 Cells..... 7

41 7.2 Technical Equipment 7

42 7.3 Reference Standard, Controls, and Tissue Culture Supplies..... 9

43 **8.0 Preparation of Tissue Culture Media and Solutions..... 10**

44 8.1 RPMI 1640 Growth Medium (RPMI)..... 10

45 8.2 Estrogen-free DMEM Medium 11

46	8.3	1X Trypsin Solution.....	11
47	8.4	1X Lysis Solution	12
48	8.5	Reconstituted Luciferase Reagent	12
49	9.0	Overview of Propagation and Plating of BG1Luc4E2 Cells	13
50	9.1	Procedures for Thawing Cells and Establishing Tissue Cultures.....	13
51	9.1.1	Thawing Cells.....	13
52	9.1.2	Establishing Tissue Cultures.....	14
53	9.2	Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free	
54		Medium, and Plating Cells for Experimentation.....	16
55	9.2.1	Ongoing Tissue Culture Maintenance.....	17
56	9.2.2	Conditioning in Estrogen-free Medium.....	17
57	9.2.3	Plating Cells Grown in Estrogen-free DMEM for	
58		Experimentation	18
59	10.0	Preparation of Test Substances.....	21
60	10.1	Determination of Test Substance Solubility.....	21
61	10.2	Preparation of Reference Standards, Control and Test Substances.....	22
62	10.2.1	Preparation of Reference Standard and Positive Control Stock	
63		Solutions	22
64	10.2.1.1	E2 Stock Solution	22
65	10.2.1.2	Methoxychlor Stock Solution.....	23
66	10.2.2	Preparation of Reference Standard and DMSO Control for	
67		Range Finder Testing	23
68	10.2.2.1	Preparation of E2 Reference Standard for	
69		Range Finder Testing.....	23
70	10.2.2.2	Preparation of DMSO Control for Range	
71		Finder Testing.....	24

72	10.2.3 Preparation of Test Substance Dosing Solutions for Range	
73	Finder Testing	24
74	10.2.4 Preparation of Reference Standard and Positive Control	
75	Dosing Solutions for Comprehensive Testing	25
76	10.2.4.1 Preparation of E2 Reference Standard for	
77	Comprehensive Testing	25
78	10.2.4.2 Preparation of Methoxychlor Control Dosing	
79	Solution for Comprehensive Testing	26
80	10.2.4.3 Preparation of DMSO Control Dosing	
81	Solution for Comprehensive Testing	26
82	10.2.5 Preparation of Test Substance Dosing Solutions for	
83	Comprehensive Testing	26
84	10.2.5.1 Preparation of Test Substance 1:2 Serial	
85	Dilutions for Comprehensive Testing	27
86	10.2.5.2 Preparation of Test Substance 1:5 Serial	
87	Dilutions for Comprehensive Testing	28
88	11.0 General Procedures for the Testing of Coded Substances	28
89	11.1 Application of Reference Standard, Controls, and Test Substances	29
90	11.1.1 Preparation of Excel® Data Analysis Template for Range Finder	
91	Testing	29
92	11.1.2 Preparation of Excel® Data Analysis Template for Comprehensive	
93	Testing	30
94	11.2 Visual Evaluation of Cell Viability	31
95	11.3 Lysis of Cells for LUMI-CELL® ER	32
96	11.4 Measurement of Luminescence	32
97	11.5 Data Analysis	32

98	11.5.1	Collection and Adjustment of Luminometer Data for	
99		Range Finder Testing	33
100	11.5.2	Collection and Adjustment of Luminometer Data for	
101		Comprehensive Testing	35
102	11.5.3	Determination of Outliers	37
103	11.5.4	Acceptance Criteria	38
104	11.5.4.1	Range Finder Testing	38
105	11.5.4.2	Comprehensive Testing.....	38
106	12.0	Range Finder Testing	39
107	13.0	Comprehensive Testing	45
108	14.0	Compilation of the Historical Quality Control Database	47
109	14.1	Compilation of Historical Quality Control Database for	
110		Phase IIa Testing.....	47
111	14.1.1	Methoxychlor Control	48
112	14.1.2	EC ₅₀	48
113	14.1.3	DMSO Control.....	48
114	14.2	Compilation of Historical Quality Control Database for	
115		Phase IIb Testing	49
116	14.3	Compilation of Historical Quality Control Database for Phase	
117		III and IV Testing	49
118	15.0	Quality Testing of Materials	49
119	15.1	Tissue Culture Media	49
120	15.2	G418.....	50
121	15.3	DMSO	51
122	15.4	Plastic Tissue Culture Materials	51

123	16.0	References.....	52
124			

LIST OF ACRONYMS AND ABBREVIATIONS		
124		
125	13 mm test tube	13 x 100 mm glass test tubes
126	DMEM	Dulbecco's Modification of Eagle's Medium
127	DMSO	Dimethyl Sulfoxide
128	DMSO control	1% v/v dilution of DMSO in tissue culture media used as a
129		vehicle control
130	E2	17 β -estradiol
131	E2 reference standard	11 Point Serial Dilution of 17 β -estradiol reference standard
132		for the LUMI-CELL® ER agonist assay
133	EC ₅₀ value	Concentration that produces a half-maximal response as
134		calculated using the four parameter Hill function.
135	ER	Estrogen Receptor
136	Estrogen-free DMEM	DMEM (phenol red free) supplemented with 1%
137		Penicillin/Streptomycin, 2% L-Glutamine, and 5%
138		Charcoal-dextran treated FBS
139	FBS	Fetal Bovine Serum
140	G418	Gentamycin
141	Methoxychlor	<i>p,p'</i> -Methoxychlor
142	Methoxychlor control	3.13 μ g/mL Methoxychlor Weak Positive Control for the
143		LUMI-CELL® ER Agonist Assay
144	RPMI	RPMI 1640 growth medium
145	TA	Transcriptional Activation
146	T25	25 cm ² tissue culture flask
147	T75	75 cm ² tissue culture flask
148	T150	150 cm ² tissue culture flask
149		

149

LIST OF FIGURES

150 Figure 7-1 pGudLuc7.ERE Plasmid 7

151 Figure 9-1 Hemocytometer Counting Grid 19

152 Figure 12-1 Agonist Range Finder Test Plate Layout 40

153 Figure 12-2 Agonist Range Finder (example 1)..... 42

154 Figure 12-3 Agonist Range Finder (example 2)..... 43

155 Figure 12-4 Agonist Range Finder (example 3)..... 43

156 Figure 12-5 Agonist Range Finder (example 4)..... 44

157 Figure 12-6 Agonist Range Finder (example 5)..... 44

158 Figure 12-7 Agonist Range Finder (example 6)..... 45

159 Figure 13-1 Agonist Comprehensive Test Plate Layout..... 46

160

161

161 **LIST OF TABLES**

162 Table 6-1 Concentrations of E2 Reference Standard Used in
 163 Comprehensive Testing..... 5

164 Table 10-1 Preparation of E2 Stock Solution 22

165 Table 10-2 Preparation of E2 Reference Standard Dosing Solution for Range
 166 Finder Testing..... 24

167 Table 10-3 Preparation of Test Substance Dosing Solutions for Range
 168 FinderTesting..... 25

169 Table 10-4 Preparation of E2 Reference Standard Dosing Solution for
 170 Comprehensive Testing..... 26

171 Table 10-5 Preparation of Test Substance 1:2 Serial Dilutions for
 172 Comprehensive Testing..... 27

173 Table 10-6 Preparation of Test Substance 1:5 Serial Dilutions for
 174 Comprehensive Testing..... 28

175 Table 11-1 Visual Observation Scoring..... 32

176 Table 11-2 Q Test Values 37

177

178

179 **1.0 PURPOSE**

180 This protocol is designed to evaluate coded test substances for potential estrogen receptor (ER)
181 agonist activity using the LUMI-CELL® ER assay.

182 **2.0 SPONSOR**

183 The National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
184 Toxicological Methods (NICEATM), P.O. Box 12233 Research Triangle Park, NC 27709

185

186 William S. Stokes, DVM, DACLAM
187 Rear Admiral, U.S. Public Health Service
188 Chief Veterinary Officer, USPHS
189 Director, NICEATM
190 National Institute of Environmental Health Sciences, NIH, DHHS
191 Bldg. 4401, Room 3129, MD EC-14
192 79 T.W. Alexander Drive
193 Research Triangle Park, NC 27709
194 Phone: 919-541-7997
195 Fax: 919-541-0947
196 Email: stokes@niehs.nih.gov

197

198 Raymond Tice, Ph.D.
199 Deputy Director, NICEATM
200 National Institute of Environmental Health Sciences
201 MD EC-17, P.O. Box 12233
202 Research Triangle Park, NC 27709
203 Phone: 919-541-4482
204 FAX: 919-541-0947
205 Email: tice@niehs.nih.gov

206

207

207 David Allen, Ph.D.
208 Principal Investigator
209 ILS, Inc./Contractor supporting NICEATM
210 National Institute of Environmental Health Sciences
211 MD EC-17, P.O. Box 12233
212 Research Triangle Park, NC 27709
213 Phone: 919-316-4664
214 FAX: 919-541-0947
215 Email: allen7@niehs.nih.gov

216

217 Frank Deal, M.S.
218 Staff Toxicologist
219 ILS, Inc./Contractor supporting NICEATM
220 National Institute of Environmental Health Sciences
221 MD EC-17, P.O. Box 12233
222 Research Triangle Park, NC 27709
223 Phone: 919-316-4587
224 FAX: 919-541-0947
225 Email: dealf@niehs.nih.gov

226

227 Patricia Ceger, M.S.
228 Project Coordinator/Technical Writer
229 ILS, Inc./Contractor supporting NICEATM
230 National Institute of Environmental Health Sciences
231 MD EC-17, P.O. Box 12233
232 Research Triangle Park, NC 27709
233 Phone: 919-316-4556
234 Fax: 919-541-0947
235 E-Mail: cegerp@niehs.nih.gov

236

237

237 **2.1 Substance Inventory and Distribution Management**

238 Cynthia Smith, Ph.D.

239 Chemistry Resources Group Leader

240 National Institute of Environmental Health Sciences

241 MD EC-06, P.O. Box 12233

242 Research Triangle Park, NC 27709

243 Phone: 919-541-3473

244

245 **3.0 DEFINITIONS**

246 • **Dosing Solution:** The test substance, control substance, or reference standard
247 solution, which is to be placed into the tissue culture wells for experimentation.

248 • **Raw Data:** Raw data includes information that has been collected but not
249 formatted or analyzed, and consists of the following:

250 ○ Data recorded in the Study Notebook

251 ○ Computer printout of initial luminometer data

252 ○ Other data collected as part of GLP compliance, e.g.:

253 ■ Equipment logs and calibration records

254 ■ Test substance and tissue culture media preparation logs

255 ■ Cryogenic freezer inventory logs

256 • **Soluble:** Test substance exists in a clear solution without visible cloudiness or
257 precipitate.

258 • **Study Notebook:** The study notebook contains recordings of all activities related
259 to the conduct of the LUMI-CELL® ER agonist assay.

260 • **Test Substances:** Substances supplied to the testing laboratories that are coded
261 and distributed such that only the Project Officer, Study Management Team
262 (SMT), and the Substance Inventory and Distribution Management have
263 knowledge of their true identity. The test substances will be purchased, aliquoted,

264 coded, and distributed by the Supplier under the guidance of the NIEHS/NTP
265 Project Officer and the SMT.

266 **4.0 TESTING FACILITY AND KEY PERSONNEL**

267 **4.1 Testing Facility**

268 Xenobiotic Detection Systems, Inc. (XDS), 1601 E. Geer St., Suite S, Durham, NC 27704

269 **4.2 Key Personnel**

- 270 • Study Director: John Gordon, Ph.D.
- 271 • Quality Assurance Director: Mr. Andrew

272 **5.0 IDENTIFICATION OF TEST AND CONTROL SUBSTANCES**

273 **5.1 Test Substances**

274 Test substances are coded and will be provided to participating laboratories by the Substance
275 Inventory and Distribution Management team.

276 **5.2 Controls**

277 Controls for the ER agonist protocol are as follows:

278 *Vehicle control (dimethyl sulfoxide [DMSO]):* 1% (v/v) DMSO (CASRN 67-68-5) diluted in
279 tissue culture media.

280 *Reference standard (17 β -estradiol [E2]):* Three concentrations of E2 (CASRN 50-28-2) in
281 duplicate for range finder testing and a serial dilution consisting of 11 concentrations of E2 in
282 duplicate for comprehensive testing

283 *Positive control (p,p'-Methoxychlor [methoxychlor]):* Methoxychlor (CASRN 72-43-5), 3.13
284 μ g/mL in tissue culture media, used as a weak positive control.

285 **6.0 OVERVIEW OF GENERAL PROCEDURES FOR AGONIST TESTING**

286 All experimental procedures are to be carried out under aseptic conditions and all solutions,
 287 glassware, plastic ware, pipettes, etc., shall be sterile. All methods and procedures shall be
 288 documented in the study notebook.

289 Agonist range finder testing is conducted on 96-well plates using four concentrations of E2
 290 (5.00×10^{-5} , 1.25×10^{-5} , 3.13×10^{-6} and 7.83×10^{-7} $\mu\text{g/mL}$) in duplicate as the reference standard
 291 and four replicate wells for the DMSO control. Range finder testing uses all wells of the 96-well
 292 plate to test six substances as seven point 1:10 serial dilutions in duplicate.

293 Comprehensive testing is conducted on 96-well plates using 11 concentrations of E2 in duplicate
 294 as the reference standard (**Table 6-1**). Four replicate wells for the DMSO control and four
 295 replicate wells for the methoxychlor control are included on each plate. Comprehensive testing
 296 uses all wells of the 96-well plate to test 2 substances as 11 point serial dilutions in triplicate.

297 **Table 6-1 Concentrations of E2 Reference Standard Used in Comprehensive Testing**

E2 Concentrations ¹		
1.00×10^{-4}	6.25×10^{-6}	3.92×10^{-7}
5.00×10^{-5}	3.13×10^{-6}	1.95×10^{-7}
2.50×10^{-5}	1.56×10^{-6}	9.78×10^{-8}
1.25×10^{-5}	7.83×10^{-7}	

298 ¹Concentrations are presented in $\mu\text{g/mL}$.

299 Visual observations for cell viability are conducted for all experimental plates just prior to
 300 luminescence measurements, as outlined in **Section 11.2**.

301 Luminescence data, measured in relative light units (RLUs), is corrected for background
 302 luminescence by subtracting the mean RLU value of the vehicle control (DMSO) wells from the
 303 RLU measurements for each of the other wells of the 96-well plate. Data is then transferred into
 304 Excel® data management spreadsheets and GraphPad PRISM® 4.0 statistical software, graphed,
 305 and evaluated as follows:

- 306 • A response is considered positive for agonist activity when the average adjusted
 307 RLU for a given concentration is greater than the mean RLU value plus three
 308 times the standard deviation for the vehicle control.

309 • Any response below this threshold is considered negative for agonist activity.

310 For substances that are positive at one or more concentrations, the concentration that causes a
311 half-maximal response (EC₅₀) is calculated using a Hill function analysis. The Hill function is a
312 four-parameter logistic mathematical model relating the substance concentration to the response
313 (typically following a sigmoidal curve) using the equation below:

$$314 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log \text{EC}_{50} - X) \text{HillSlope}}}$$

315 where Y = response (i.e., relative light units); X = the logarithm of concentration; Bottom = the
316 minimum response; Top = the maximum response; log EC₅₀ = the logarithm of X as the response
317 midway between Top and Bottom; and HillSlope describes the steepness of the curve. The model
318 calculates the best fit for the Top, Bottom, HillSlope, and EC₅₀ parameters. See **Section 11.6.5**
319 for more details.

320 Acceptance or rejection of a test is based on evaluation of reference standard and control results
321 from each experiment conducted on a 96-well plate. Results for these controls are compared to
322 historical results compiled in the historical database, as seen in **Section 14.0**.

323 **6.1 Range Finder Testing**

324 Agonist range finding for coded substances consists of a seven point, 1:10 serial dilution using
325 duplicate wells per concentration. Concentrations for comprehensive testing are selected based
326 on the response observed in range finder testing. If necessary, a second range finder test can be
327 conducted to clarify the optimal concentration range to test (see **Section 12.0**).

328 **6.2 Comprehensive Testing**

329 Comprehensive agonist testing for coded substances consists of 11 point, serial dilutions, with
330 each concentration tested in triplicate wells of the 96-well plate. Three separate experiments are
331 conducted for comprehensive testing on three separate days, except during Phases III and IV of
332 the validation effort, in which comprehensive testing experiments are conducted once (see
333 **Section 13.0**).


334 7.0 MATERIALS FOR LUMI-CELL® ER AGONIST TESTING

335 This section provides the materials needed to conduct LUMI-CELL® ER testing, with associated
336 brand names/vendors¹ in brackets.

337 7.1 BG1Luc4E2 Cells:

338 Human ovarian cancer cell line stably transfected with a plasmid containing an estrogen response
339 element pGudLuc7.0 (Figure 7-1) [XDS].

340 **Figure 7-1 pGudLuc7.ERE Plasmid.**



341

342 7.2 Technical Equipment:

343 All technical equipment may be obtained from Fisher Scientific International, Inc. (Liberty Lane
344 Hampton, NH, USA 03842). Equivalent technical equipment from another commercial source
345 can be used.

- 346
- Analytical balance (Cat. No. 01-910-320)
 - 347
 - Berthold Orion 1 Microplate Luminometer [Berthold CatNo.: Orion 1 MPL3] or
348 equivalent and dedicated computer
 - 349
 - Biological safety hood, class II, and stand (Cat. No. 16-108-99)

¹Brand names and vendors should not be considered an endorsement by the U.S. Government or any member of the U.S. Government; such information is provided as examples.

- 350 • Centrifuge (low speed, tabletop with swinging bucket rotor) (Cat. No. 04-978-50
- 351 centrifuge, and 05-103B rotor)
- 352 • Combustion test kit (CO₂ monitoring) (Cat. No. 10-884-1)
- 353 • Drummond diaphragm pipetter (Cat. No. 13-681-15)
- 354 • Freezers, -20°C (Cat. No. 13-986-150), and -70°C (Cat. No. 13-990-86)
- 355 • Hand tally counter (Cat. No. 07905-6)
- 356 • Hemocytometer, cell counter (Cat. No. 02-671-5)
- 357 • Light microscope, inverted (Cat. No. 12-561-INV)
- 358 • Light microscope, upright (Cat. No. 12-561-3M)
- 359 • Liquid nitrogen flask (Cat. No. 11-675-92)
- 360 • Micropipetter, repeating (Cat. No. 21-380-9)
- 361 • Pipettors, air displacement, single channel (0.5 –10µl (Cat. No. 21-377-191), 2 –
- 362 20 µl (Cat. No. 21-377-287), 20 – 200 µl (Cat. No. 21-377-298), 200 - 1000 µl
- 363 (Cat. No. 21-377-195))
- 364 • Refrigerator/freezer (Cat. No. 13-986-106A)
- 365 • Shaker for 96-well plates (Cat. No. 14-271-9)
- 366 • Sodium hydroxide (Cat. No. 5318-500)
- 367 • Sonicating water bath (Cat. No. 15-335-30)
- 368 • Tissue culture incubator with CO₂ and temperature control (Cat. No. 11-689-4)
- 369 • Vacuum pump with liquid trap (side arm Erlenmeyer) (Cat. No. 01-092-29)
- 370 • Vortex mixer (Cat. No. 12-814)

371 Equipment should be maintained and calibrated as per GLP guidelines and individual laboratory
372 SOPs.

373

373 7.3 Reference Standard, Controls, and Tissue Culture Supplies

374 All tissue culture reagents must be labeled to indicate source, identity, storage conditions and
375 expiration dates. Tissue culture solutions must be labeled to indicate concentration, stability
376 (where known), and preparation and expiration dates.

377 Equivalent tissue culture media and sera from another commercial source can be used, but must
378 first be tested as described in **Section 15.0** to determine suitability for use in this test method.

379 The following are the necessary tissue culture reagents and possible commercial sources (in
380 brackets) based on their use in the pre-validation studies:

- 381 • BackSeal-96/384, white adhesive bottom seal for 96-well and 384-well microplate
382 [Perkin-Elmer, Cat. No. 6005199]
- 383 • 17 β -estradiol (CAS RN: 50-28-2) [Sigma-Aldrich, Cat. No. E8875]
- 384 • Cryovial, 2 mL (Corning Costar) [Fisher Scientific Cat. No. 03-374-21]
- 385 • Culture tube 13 x 100mm (case) [Thomas Scientific Cat. No.: 10009186R38]²
- 386 • Culture tube, 50 mL conical (Corning Costar) [Fisher Scientific Cat. No. 05-
387 526C]
- 388 • DMSO, U.S.P. analytical grade. [Sigma-Aldrich, Cat. No. 34869-100ML]
- 389 • Dulbecco's Modification of Eagle's Medium (DMEM), containing 4.5 g/L
390 glucose, with sodium pyruvate, without phenol red or L-glutamine
391 [Mediatech/Cellgro, Cat. No. 17-205-CV]
- 392 • Fetal Bovine Serum [Mediatech/Cellgro Cat. No. MT 35-010-CV]
- 393 • Fetal Bovine Serum, charcoal/dextran treated, triple 0.1 μ m sterile filtered
394 [Hyclone, Cat. No. SH30068.03]
- 395 • Gentamycin Sulfate (G418), 50 mg/mL [Mediatech/Cellgro Cat. No. 30-234-CR]
- 396 • L-glutamine, 29.2 mg/mL [Cellgro, Cat. No. 25005-CI]

² If glass tubes can not be obtained from Thomas Scientific, the preference is for flint glass, then lime glass, then borosilicate glass.

- 397 • Luciferase Assay System (10-Pack) [Promega Cat. No. E1501]
- 398 • Lysis Solution 5X [Promega, Cat. No. E1531]
- 399 • Methoxychlor (CAS RN: 72-43-5) [Sigma-Aldrich, Cat. No. 49054]
- 400 • Penicillin/streptomycin solution, 5000 I.U. penicillin, 5000 µg/mL streptomycin
- 401 [Cellgro, Cat. No. 30-001-CI].
- 402 • Phosphate buffered saline (PBS, 1X) without calcium and magnesium [Cellgro,
- 403 Cat. No. 21-040-CV]
- 404 • Pipettes, serological: 2.0 mL [Sigma-Aldrich, Cat. No. P1736], 5.0 mL [Sigma-
- 405 Aldrich, Cat. No. P1986], 25 mL [Sigma-Aldrich, Cat. No. P2486]
- 406 • RPMI 1640 medium, containing L-glutamine [Mediatech, Cat. No. 10-040-CV]
- 407 • Tissue culture flasks (Corning-Costar): 25 cm² (T25) [Fisher Cat. No. 10-126-28];
- 408 75 cm² (T75) [Fisher Cat. No. 10-126-37]; and 150 cm² (T150) [Fisher Cat. No.
- 409 10-126-34]
- 410 • Tissue culture plates (Corning-Costar): 96-well [Thomas Scientific Cat. No.
- 411 6916A05]
- 412 • Trypsin (10X), 2.5% in Hank's balanced salt solution (HBSS), without calcium
- 413 and magnesium, without phenol red [Cellgro, Cat. No. 25-054-CI].

414 All reagent lot numbers and expiration dates must be recorded in the study notebook.

415 **8.0 PREPARATION OF TISSUE CULTURE MEDIA AND SOLUTIONS**

416 All tissue culture media and media supplements must be quality tested before use in experiment
417 (see Section 15.0).

418 **8.1 RPMI 1640 Growth Medium (RPMI)**

419 RPMI 1640 is supplemented with 0.9% Pen-Strep and 8.0% FBS to make RPMI growth medium
420 (RPMI).

421 Procedure for one 549 mL bottle:

- 422 1. Remove FBS from -70°C freezer, and Pen-Strep from -20°C freezer and allow to
423 equilibrate to room temperature.
- 424 2. Add 44 mL of FBS and 5 mL Pen-Strep to the bottle of RPMI 1640.
- 425 3. Label RPMI bottle as indicated in **Section 7.3**

426 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
427 *component.*

428 **8.2 Estrogen-Free DMEM Medium**

429 DMEM is supplemented to contain 4.5% charcoal/dextran treated FBS, 1.9% L-glutamine, 0.9%
430 Pen-Strep.

431 Procedure for one 539 mL bottle:

- 432 1. Remove charcoal/dextran treated FBS from -70°C freezer, and L-glutamine and
433 Pen-Strep from -20°C freezer and allow to equilibrate to room temperature.
- 434 2. Add 24 mL of charcoal/dextran treated FBS, 10 mL L-glutamine, and 5 mL Pen-
435 Strep to one 500 mL bottle of DMEM.
- 436 3. Label estrogen-free DMEM bottle as indicated in **Section 7.3**

437 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
438 *component..*

439 **8.3 1X Trypsin Solution**

440 1X Trypsin solution is prepared by dilution from a 10X premixed stock solution. The 10X stock
441 solution should be stored in 10 mL aliquots in a -20°C freezer.

442 Procedure for making 100 mL of 1X trypsin:

- 443 1. Remove a 10 mL aliquot of 10X trypsin from -20°C freezer and allow to
444 equilibrate to room temperature.
- 445 2. Aliquot 1 mL Trypsin (10X) along with 9 mL of 1X PBS into ten 15 mL sterile
446 centrifuge tubes.
- 447 3. Label 1X trypsin aliquots as indicated in **Section 7.3**

448 *1X Trypsin should be stored at -20°C.*

449 **8.4 1X Lysis Solution**

450 Lysis solution is prepared by dilution from a 5X premixed stock solution. Both the 5X and 1X
451 solutions can be repeatedly freeze-thawed.

452 The procedure for making 10 mL of 1X lysis solution:

- 453 1. Thaw the 5X Promega Lysis solution and allow it to reach room temperature.
- 454 2. Remove 2 mL of 5X solution and place it in a 15 mL conical centrifuge tube.
- 455 3. Add 8 mL of distilled, de-ionized water to the conical tube.
- 456 4. Cap and shake gently until solutions are mixed.

457 *Store at -20°C for no longer than 1 year from receipt.*

458 **8.5 Reconstituted Luciferase Reagent**

459 Luciferase reagent consists of two components, luciferase buffer and lyophilized luciferase
460 substrate.

461 For long term storage, unopened containers of the luciferase buffer and lyophilized luciferase
462 substrate can be stored at -70°C for up to one year.

463 To reconstitute luciferase reagent:

- 464 1. Remove luciferase buffer and luciferase substrate from -70°C freezer and allow
465 them to equilibrate to room temperature.
- 466 2. Add 10 mL of luciferase buffer solution to luciferase substrate container and swirl
467 or vortex gently to mix; the Luciferase substrate should readily go into solution.
- 468 3. After solutions are mixed, aliquot to a 15mL centrifuge tube.
- 469 4. Store complete solution at -20°C.

470 Reconstituted luciferase reagent is stable for up to 1 month at -20°C.

471 **9.0 OVERVIEW OF PROPAGATION AND EXPERIMENTAL PLATING OF**
472 **BG1Luc4E2 CELLS**

473 The BG1Luc4E2 (BG-1) cells are stored in liquid nitrogen in 2 mL cryovials. BG-1 cells are
474 grown as a monolayer in tissue culture flasks in a dedicated tissue culture incubator at 37°C ±
475 1°C, 90% ± 5% humidity, and 5.0% ± 1% CO₂/air. The cells should be examined, on a daily
476 basis during working days, under an inverted phase contrast microscope and any changes in
477 morphology and/or adhesive properties must be noted in the study notebook.

478 Two T150 flasks containing cells at 80 to 90% confluence will usually yield a sufficient number
479 of cells to fill three 96-well plates for use in experiments.

480 **9.1 Procedures for Thawing Cells and Establishing Tissue Cultures**

481 Warm all of the tissue culture media and solutions to room temperature by placing them under
482 the tissue culture hood several hours before use.

483 All tissue culture media, media supplements, and tissue culture plasticware must be quality
484 tested before use in experiments (**Section 15.0**).

485 9.1.1 Thawing Cells

- 486 1. Remove a cryovial of frozen BG-1 cells from the liquid nitrogen flask.
- 487 2. Facilitate rapid thawing by loosening the top slightly (do not remove top) to
488 release trapped gasses and retightening it. Roll vial between palms.
- 489 3. Use a micropipette to transfer cells to a 50 mL conical centrifuge tube.
- 490 4. Rinse cryovial twice with 1X PBS and add PBS rinse material to the conical tube.
- 491 5. Add 20 mL of RPMI to the conical tube.
- 492 6. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
493 for an additional 5 minutes.
- 494 7. Aspirate media from pellet and re-suspend it in 5 mL RPMI, drawing the pellet
495 repeatedly through a 1.0 mL serological pipette to break up any clumps of cells.
- 496 8. Transfer cells to a T25 flask, place them in an incubator (see conditions in
497 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

498 9.1.2 Establishing Tissue Cultures

499 Once cells have reached 80 to 90% confluence, transfer the cells to a T75 flask by performing,
500 for example, the following steps:

- 501 1. Remove the T25 flask from the incubator.
- 502 2. Aspirate the RPMI, then add 5 mL 1X PBS, making sure that the cells are coated
503 with PBS.
- 504 3. Aspirate 1X PBS, then add 1 to 2 mL 1X trypsin to the T25 flask, gently swirling
505 the flask to coat all cells with the trypsin.
- 506 4. Place the flask in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 507 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
508 hand.
- 509 6. Confirm cell detachment by examination under an inverted microscope. If cells
510 have not detached, return the flask to the incubator for an additional 2 minutes,
511 then hit the flask again.
- 512 7. After cells have detached, add 5 mL PBS, and transfer the suspended cells to a 50
513 mL centrifuge tube. Wash the flask one additional time with 5 mL PBS.
- 514 8. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
515 digestion by residual trypsin.
- 516 9. Pellet the cells by centrifugation, as described in **Section 9.1.1**, and re-suspend the
517 cells in 10 mL RPMI medium.
- 518 10. Draw the pellet repeatedly through a 25 mL serological pipette to break up
519 clumps of cells
- 520 11. Transfer cells to a T75 flask, then place the flask in an incubator (see conditions
521 in **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

522 When cells have reached 80% to 90% confluency, transfer them into a T150 flask by performing,
523 for example, the following steps:

- 524 12. Remove the T75 flask from the incubator, aspirate the old media and add 5 mL
525 1X PBS.
- 526 13. Aspirate 1X PBS, add 2 mL of 1X trypsin to the flask, and place it in an incubator
527 (see conditions in **Section 9.0**) for 5 to 10 min.
- 528 14. Repeat steps **5** through **11** in **Section 9.1.2**, re-suspending the pellet in 20 mL of
529 RPMI.
- 530 15. Transfer cells to a T150 flask and place it in the incubator (see conditions in
531 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 532 16. Remove the T150 flask from the incubator.
- 533 17. Aspirate the RPMI and add 5 mL 1X PBS.
- 534 18. Aspirate 1X PBS and add 3 mL 1X trypsin to the T150 flask, making sure that the
535 cells are coated with the trypsin.
- 536 19. Incubate cells in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 537 20. Detach cells by hitting the side of the flask sharply against the palm or heel of the
538 hand.
- 539 21. Confirm cell detachment by examination under an inverted microscope. If cells
540 have not detached, return the flask to the incubator for an additional 2 minutes,
541 then hit the flask again.
- 542 22. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
543 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
544 flask, swirl around the flask, and then transfer the PBS to the 50 mL conical tube.
- 545 23. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
546 digestion by residual trypsin.
- 547 24. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
548 for an additional 5 minutes.

549 25. Aspirate the media from the pellet and re-suspend it in 40 mL RPMI, drawing the
550 pellet repeatedly through a 25 mL serological pipette to break up any clumps of
551 cells.

552 26. Transfer 20 mL of cell suspension to each of two T150 flasks, place them in an
553 incubator (see conditions in **Section 9.0**) and grow to 80% to 90% confluence
554 (approximately 48 to 72 hrs).

555 **9.2 Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free Medium,** 556 **and Plating Cells for Experimentation**

557 The following procedure is used to condition the BG1Luc4E2 cells to an estrogen-free
558 environment prior to plating the cells in 96-well plates for analysis of estrogen dependent
559 induction of luciferase activity.

560 To start the tissue culture maintenance and estrogen-free conditioning, split the two T150 culture
561 flasks into four T150 flasks. Two of these flasks will be used for continuing tissue culture and
562 will use the RPMI media mentioned above. The other two flasks will be cultured in estrogen-free
563 DMEM for experimental use. Extra care must be taken to avoid contaminating the estrogen-free
564 cells with RPMI.

- 565 1. Remove both T150 flasks from the incubator.
- 566 2. Aspirate the medium and rinse the cells with 5 mL 1X PBS.
- 567 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
568 to coat all cells with the trypsin.
- 569 4. Incubate cells in the incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 570 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
571 hand.
- 572 6. Confirm cell detachment by examination under an inverted microscope. If cells
573 have not detached, return the flask to the incubator for an additional 2 minutes,
574 then hit the flask again.
- 575 7. After cells have detached, add 5 mL 1X PBS to the first T150 flask and transfer
576 the suspended cells to the second T150 flask.

- 577 8. Transfer the contents of both flasks to a 50 mL conical tube. Repeat step 7 with an
578 additional 5 mL 1X PBS and transfer to the 50 mL conical tube.
- 579 9. Immediately add 20 mL estrogen-free DMEM to the 50 mL conical tube to inhibit
580 further cellular digestion by residual trypsin.
- 581 10. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
582 for an additional 5 minutes.
- 583 11. Aspirate media from pellet and re-suspend it in 4 mL estrogen-free DMEM,
584 drawing the pellet repeatedly through a 1 mL serological pipette to break up
585 clumps of cells.

586 At this point, cells are ready to be divided into the ongoing tissue culture and estrogen-free
587 conditioning groups.

588 9.2.1 Ongoing Tissue Culture Maintenance

- 589 1. Add 20 mL RPMI to two T150 flasks.
- 590 2. Add 220 µl G418 to the RPMI in the T150 flasks
- 591 3. Add 1 mL of cell suspension from **9.2 step 11** to each flask.
- 592 4. Place T150 flasks in tissue culture incubator (see conditions in **Section 9.0**) and
593 grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 594 5. Tissue culture medium may need to be changed 24 hours after addition of G418 to
595 remove cells that have died because they do not express reporter plasmid.
- 596 6. G418 does not need to be added to the flasks a second time.
- 597 7. Repeat **Section 9.2 steps 1-11** for ongoing tissue culture maintenance.

598 9.2.2 Conditioning in Estrogen-free Medium

- 599 1. Add 20 mL estrogen-free DMEM to two T150 flasks.
- 600 2. Add 150 µL G418 to the estrogen-free DMEM in the T150 flasks.
- 601 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.


- 602 4. Tissue culture medium may need to be changed 24 hours after addition of G418 to
603 remove cells that have died because they do not express reporter plasmid.
- 604 5. G418 does not need to be added to the flasks a second time.
- 605 6. Place the T150 flasks in the incubator (see conditions in **Section 9.0**) and grow to
606 80% to 90% confluence (approximately 48 to 72 hrs).

607 9.2.3 Plating Cells Grown in Estrogen-free DMEM for Experimentation

- 608 1. Remove the T150 flasks that have been conditioned in estrogen-free DMEM for
609 48 to 72 hours from the incubator.
- 610 2. Aspirate the medium, then rinse the cells with 5 mL 1X PBS.
- 611 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
612 to coat all cells with the trypsin.
- 613 4. Place the flasks in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 614 5. Detach cells by hitting the side of the flask sharply against the palm or the heel of
615 the hand.
- 616 6. Confirm cell detachment by examination under an inverted microscope. If cells
617 have not detached, return the flask to the incubator for 2 additional minutes, then
618 hit the flask again.
- 619 7. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
620 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
621 flask, gently swirl around the flask, and then transfer to the 50 mL conical tube.
- 622 8. Immediately add 20 mL estrogen-free DMEM to each conical tube to inhibit
623 further cellular digestion by residual trypsin.
- 624 9. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
625 for an additional 5 minutes.
- 626 10. Aspirate the media from the pellet and re-suspend it in 20 mL DMEM, drawing
627 the pellet repeatedly through a 25 mL serological pipette to break up any clumps
628 of cells.

- 629 11. Pipette 15 μ L of the cell suspension into the “v” shaped slot on the
 630 hemocytometer. Ensure that the solution covers the entire surface area of the
 631 hemocytometer grid, and allow cells to settle before counting.
- 632 12. Using 100x magnification, view the counting grid.
- 633 13. The counting grid on the hemocytometer consists of nine sections, four of which
 634 are counted (upper left, upper right, lower left, and lower right, see **Figure 9-1**).
 635 Each section counted consists of four by four grids. Starting at the top left and
 636 moving clockwise, count all cells in each of the four by four grids. Some cells
 637 will be touching the outside borders of the square, but only count those that touch
 638 the top and right borders of the square. This value is then used in the calculation
 639 below to get to the desired concentration of 200,000 cells/mL.
- 640

Figure 9-1 Hemocytometer Counting Grid.



641 The volume of each square is 10^{-4} mL, therefore:

642
$$\text{Cells/mL} = (\text{average number per grid}) \times 10^{-4} \text{ mL} \times 1/(\text{starting dilution}).$$

643 Starting dilution: 20 mL (for T150 flasks)

644

645 Harvested cells for a T150 flask are suspended in 20 mL of estrogen-free DMEM and sampled
 646 for determination of concentration of cells/mL.

647

648 Example Calculation:

649 • Grids 1, 2, 3, and 4 are counted and provide the following data:

650 ○ 50, 51, 49, and 50: average number of cells per grid is equal to 50.

651 Cells/mL = 50 cells per grid ÷ 10^{-4} mL volume of grid = 50×10^4 cells/mL (or 500,000
652 cells/mL)

653 Total # of Cells Harvested = 500,000 cells/mL x 20 mL

654 Desired Concentration (or Concentration_{Final}) = 200,000 cells/mL

655 Formula: (Concentration_{Final} x Volume_{Final} = Concentration_{Initial} x Volume_{Initial})

656 Concentration_{Final} = 200,000 cells/mL

657 Concentration_{Initial} = 500,000 cells/mL

658 Volume_{Initial} = 20 mL

659 Volume_{Final} – to be solved for.

660 Therefore: 200,000 cells/mL x Volume_{Final} = 500,000 cells/mL x 20 mL

661 Solving for Volume_{Final} we find = 50 mL

662 Therefore, add 30 mL of estrogen-free DMEM to the cell suspension for a total volume of 50
663 mL, which will yield the desired concentration of 200,000 cells/mL for plating.

664 14. This dilution scheme will give a concentration of 200,000 cells/mL. 200 µL of
665 this cell suspension is used for each well of a 96-well plate (i.e., 40,000 cells per
666 well).

667 15. Remove a 96-well plate from its sterile packaging. Use a repeater pipetter to
668 pipette 200 µL of cell suspension into each well for to be used for the testing of
669 coded substances, reference standard and controls (**note:** add 200 µL of estrogen-
670 free DMEM only to any wells not being used for testing).

671 16. Incubate plate(s) in an incubator (see conditions in **Section 9.0**) for a minimum of
672 24 hours, but no longer than 48 hours before dosing.

673 Two T150 flasks containing cells at 80% to 90% confluence will typically yield sufficient cells
674 to fill four 96-well plates.

675 **10.0 PREPARATION OF TEST SUBSTANCES**

676 The solvent used for dissolution of test substances is 100% DMSO. All test substances should be
677 allowed to equilibrate to room temperature before being dissolved and diluted. Test substance
678 solutions (except for reference standards and controls) should not be prepared in bulk for use in
679 subsequent tests. Test substances are to be used within 24 hours of preparation. Solutions should
680 not have noticeable precipitate or cloudiness.

681 All information on weighing, solubility testing, and calculation of final concentrations for test
682 substances, reference standards and controls is to be recorded in the study notebook.

683 **10.1 Determination of Test Substance Solubility**

- 684 1. Prepare a 100 mg/mL solution of the test substance in 100% DMSO in a 4 mL
685 conical tube.
- 686 2. Vortex to mix.
- 687 3. If the test substance does not dissolve at 100 mg/mL, prepare a 10 mg/mL
688 solution and vortex as above.
- 689 4. If the test substance does not dissolve at 10 mg/mL solution, prepare a 1 mg/mL
690 solution in a 4 mL conical tube and vortex as above.
- 691 5. If the test substance does not dissolve at 1 mg/mL, prepare a 0.1 mg/mL solution
692 in a 4 mL conical tube and vortex as above.
- 693 6. Continue testing, using 1/10 less substance in each subsequent attempt until test
694 substance is solubilized in DMSO.

695 Once the test substance has fully dissolved in 100% DMSO, the solubility of the test substance
696 must be determined in the 1% DMSO/99% estrogen-free DMEM mixture used for LUMI-
697 CELL® ER testing.

- 698 7. Add 4 µL of the highest concentration of the test substance/DMSO solution to a
699 13 mm test tube.

- 700 8. Add 400 µL estrogen-free DMEM to the test tube and vortex gently,
 701 9. If cloudiness or precipitate develop, vortex for up to 10 minutes.
 702 10. If test substance has visible precipitate or is cloudy return to **10.1 step 7** to try the
 703 next lower concentration for the test substance.

704 The Testing Facility shall forward the results from the solubility tests assay to the SMT through
 705 the designated contacts in electronic format and hard copy upon completion of testing.

706 **10.2 Preparation of Reference Standards, Control and Test Substances**

707 All “dosing solutions” of test substance concentrations are to be expressed as µg/mL in the study
 708 notebook and in all laboratory reports.

709 All information on preparation of test substances, reference standards and controls is to be
 710 recorded in the study notebook.

711 10.2.1 Preparation of Reference Standard and Positive Control Stock Solutions

712 Stock solutions of E2 and methoxychlor are prepared in 100% DMSO and stored at room
 713 temperature for up to three years or until the expiration date listed in the certificate of analysis
 714 for that substance.

715 10.2.1.1 *E2 Stock Solution*

716 The final concentration of the E2 stock solution is 1.0×10^{-2} µg/mL. Prepare the E2 stock as
 717 shown in **Table 10-1**.

718 **Table 10-1 Preparation of E2 Stock Solution**

Step #	Action	DMSO	E2 Concentration
1	Make a 10 mg/mL stock solution in 100% DMSO in a 4mL vial.	-	10 mg/mL
2	Transfer 10 µL E2 solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 10 µL E2 solution from Step #2 to a new 4mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	1 µg/mL
4	Transfer 10 µL E2 solution from Step #3 to a 13 mm test tube to create the working solution.	Add 990 µL of 100% DMSO. Vortex to mix.	1.0×10^{-2} µg/mL

719

720 10.2.1.2 *Methoxychlor Stock Solution*

721 The final concentration of the methoxychlor stock solution is 313 µg/mL.

722 To prepare the methoxychlor stock solution, proceed as follows:

- 723 1. Make a 10 mg/mL stock solution of Methoxychlor in 100% DMSO in a 4 mL
- 724 vial.
- 725 2. Remove 94 µL of the methoxychlor solution and place it in a new 4 mL vial.
- 726 3. Add 2.906 mL of 100% DMSO to the 4mL vial and gently vortex to mix.

727 10.2.2 Preparation of Reference Standard and DMSO Control for Range Finder Testing

728 Range finder testing is conducted on 96-well plates using four concentrations of E2 in duplicate
729 as the reference standard. Four replicate wells are used for the DMSO control. All wells on the
730 96 well plate are used during range finder testing.

731 Store dosing solutions at room temperature. Use within 24 hours of preparation.

732 10.2.2.1 *Preparation of E2 Reference Standard for Range Finder Testing*

733 To make E2 dosing solutions:

- 734 1. label four 4 mL conical tubes with numbers 1 through 4 and place them in a tube
- 735 rack
- 736 2. label four 13 mm glass test tubes with numbers 1 through 4, place them in a tube
- 737 rack and add 600 µL of estrogen-free DMEM to each tube

738 Prepare dilutions to give final concentrations of the E2 as shown in **Table 10-2**.

739

739 **Table 10-2 Preparation of E2 Reference Standard Dosing Solution for Range Finder**
 740 **Testing**

Tube Number	100% DMSO	E2 ¹	Estrogen-free DMEM ²	Final Volume	E2 Concentration
1	6 µL	6 µl of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	5.00 x 10 ⁻⁵ µL
2	18 µL	6 µL of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	1.25 x 10 ⁻⁵ µL
3	18 µL	6 µL from conical tube #2	600 µL	606 µL	3.13 x 10 ⁻⁶ µL
4	18 µL	6 µL from conical tube #3	600 µL	606 µL	7.83 x 10 ⁻⁷ µL

741 ¹Add specified volume of 100% DMSO and 6 µl of the specified E2 solution to labeled 4 mL conical
 742 tubes, and vortex.

743 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
 744 DMEM and vortex.

745 10.2.2.2 *Preparation of DMSO Control for Range Finder Testing*

746 1. Add 10 µL of 100% DMSO to a 13 mm glass tube.

747 2. Add 1000 µL of estrogen-free DMEM to tube and vortex.

748 10.2.3 Preparation of Test Substance Dosing Solutions for Range Finder Testing

749 Range finder experiments are used to determine the concentrations of test substance to be used
 750 during comprehensive testing. Agonist range finding for coded substances consists of seven
 751 point, 1:10 serial dilutions run in duplicate.

752 To make dosing solutions for coded substances:

753 1. label seven 4 mL conical tubes with numbers 1 through 7 and place them in a tube
 754 rack

755 2. label seven 13 mm glass test tubes with numbers 1 through 7, place them in a tube
 756 rack and add 600 µL of estrogen-free DMEM to each tube

757 Prepare dilutions as shown in **Table 10-3**.

758 **Table 10-3 Preparation of Test Substance Dosing Solutions for Range Finder Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	6 µL of test substance solution from Section 10.1 step 10	6 µL	600 µL	606 µL

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
2	90 µL	10 µL of test substance solution from Section 10.1 step 10	6 µL	600 µL	606 µL
3	90 µL	10 µL from conical tube #2	6 µL	600 µL	606 µL
4	90 µL	10 µL from conical tube #3	6 µL	600 µL	606 µL
5	90 µL	10 µL from conical tube #4	6 µL	600 µL	606 µL
6	90 µL	10 µL from conical tube #5	6 µL	600 µL	606 µL
7	90 µL	10 µL from conical tube #6	6 µL	600 µL	606 µL

759 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
760 vortex.

761 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
762 DMEM and vortex.

763
764 Determination of whether a substance is positive in range finder testing and selection of starting
765 concentrations for comprehensive testing will be discussed in **Section 12.0**.

766 10.2.4 Preparation of Reference Standard and Positive Control Dosing Solutions for
767 Comprehensive Testing

768 Comprehensive testing is conducted on 96-well plates using 11 concentrations of E2 in duplicate
769 as the reference standard. Four replicate wells for the DMSO control and three replicate wells for
770 the methoxychlor control are included on each plate.

771 Store dosing solutions at room temperature. Use within 24 hours of preparation.

772 10.2.4.1 *Preparation of E2 Reference Standard for Comprehensive Testing*

773 To make E2 dosing solutions:

- 774 1. label 11 4 mL conical tubes with numbers 1 through 11 and place them in a tube
775 rack
- 776 2. label 11 13 mm glass test tubes with numbers 1 through 11, place them in a tube
777 rack and add 600 µL of DMEM to each tube

778 Prepare dilutions to give final concentrations of E2 as shown in **Table 10-4**.

779

779 **Table 10-4 Preparation of E2 Reference Standard Dosing Solution for**
 780 **Comprehensive Testing**

Tube Number	100% DMSO	E2 ¹	Estrogen-free DMEM ²	Final Volume	E2 Concentration
1	-	6 µl of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	1.00 x 10 ⁻⁴ µL
2	6 µL	6 µL of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	5.00 x 10 ⁻⁵ µL
3	6 µL	6 µL from conical tube #2	600 µL	606 µL	2.50 x 10 ⁻⁵ µL
4	6 µL	6 µL from conical tube #3	600 µL	606 µL	1.25 x 10 ⁻⁵ µL
5	6 µL	6 µL from conical tube #4	600 µL	606 µL	6.25 x 10 ⁻⁶ µL
6	6 µL	6 µL from conical tube #5	600 µL	606 µL	3.13 x 10 ⁻⁶ µL
7	6 µL	6 µL from conical tube #6	600 µL	606 µL	1.56 x 10 ⁻⁶ µL
8	6 µL	6 µL from conical tube #7	600 µL	606 µL	7.83 x 10 ⁻⁷ µL
9	6 µL	6 µL from conical tube #8	600 µL	606 µL	3.92 x 10 ⁻⁷ µL
10	6 µL	6 µL from conical tube #9	600 µL	606 µL	1.95 x 10 ⁻⁷ µL
11	6 µL	6 µL from conical tube #10	600 µL	606 µL	9.78 x 10 ⁻⁸ µL

781 ¹Add specified volume of 100% DMSO and 6 µl of the specified E2 solution to labeled 4 mL conical
 782 tubes, and vortex.

783 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
 784 DMEM and vortex.
 785

786 10.2.4.2 *Preparation of Methoxychlor Control Dosing Solution for Comprehensive Testing*

787 1. Add 10 µL of the 313 µg/mL methoxychlor to a 13 mm glass tube.

788 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex.

789 10.2.4.3 *Preparation of DMSO Control Dosing Solution for Comprehensive Testing*

790 1. Add 10 µL of 100% DMSO to four 13 mm tubes (solvent/negative controls).

791 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex.

792 10.2.5 Preparation of Test Substance Dosing Solutions for Comprehensive Testing

793 Comprehensive testing experiments are used to determine whether a substance possesses ER
 794 agonist activity in the LUMI-CELL® ER test method. Agonist comprehensive testing for coded
 795 substances consists of either an 11 point 1:2 serial dilution or an 11 point 1:5 serial dilution,

796 depending on the results from range finder testing (see **Section 12.0**) with each concentration
 797 tested in triplicate wells of the 96-well plate.

798 10.2.5.1 *Preparation of Test Substance 1:2 Serial Dilutions for*
 799 *Comprehensive Testing*

800 Start the 11-point serial dilution according to criteria in **Section 12.0**.

801 To make test substance 1:2 serial dilutions for comprehensive testing:

- 802 1. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
 803 tube rack
- 804 2. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
 805 tube rack and add 800 µL of estrogen-free DMEM to each tube

806 Prepare dilution of test substance as shown in **Table 10-6**.

807 **Table 10-5 Preparation of Test Substance 1:2 Serial Dilutions for Comprehensive**
 808 **Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
2	8 µL	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
3	8 µL	8 µL from conical tube #2	8 µL	800 µL	808 µL
4	8 µL	8 µL from conical tube #3	8 µL	800 µL	808 µL
5	8 µL	8 µL from conical tube #4	8 µL	800 µL	808 µL
6	8 µL	8 µL from conical tube #5	8 µL	800 µL	808 µL
7	8 µL	8 µL from conical tube #6	8 µL	800 µL	808 µL
8	8 µL	8 µL from conical tube #7	8 µL	800 µL	808 µL
9	8 µL	8 µL from conical tube #8	8 µL	800 µL	808 µL
10	8 µL	8 µL from conical tube #9	8 µL	800 µL	808 µL
11	8 µL	8 µL from conical tube #10	8 µL	800 µL	808 µL

809 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
 810 vortex.

811

812 10.2.5.2 *Preparation of Test Substance 1:5 Serial Dilutions for Comprehensive*
 813 *Testing*

814 Start the 11-point serial dilution according to criteria in **Section 12.0**.

815 To make test substance 1:5 serial dilutions for comprehensive testing:

816 3. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
 817 tube rack

818 4. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
 819 tube rack and add 800 µL of estrogen-free DMEM to each tube

820 Prepare dilution of test substance as shown in **Table 10-6**.

821 **Table 10-6 Preparation of Test Substance 1:5 Serial Dilutions for Comprehensive**
 822 **Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
2	16 µL	4 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
3	16 µL	4 µL from conical tube #2	8 µL	800 µL	808 µL
4	16 µL	4 µL from conical tube #3	8 µL	800 µL	808 µL
5	16 µL	4 µL from conical tube #4	8 µL	800 µL	808 µL
6	16 µL	4 µL from conical tube #5	8 µL	800 µL	808 µL
7	16 µL	4 µL from conical tube #6	8 µL	800 µL	808 µL
8	16 µL	4 µL from conical tube #7	8 µL	800 µL	808 µL
9	16 µL	4 µL from conical tube #8	8 µL	800 µL	808 µL
10	16 µL	4 µL from conical tube #9	8 µL	800 µL	808 µL
11	16 µL	4 µL from conical tube #10	8 µL	800 µL	808 µL

823 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
 824 vortex.
 825

826 **11.0 GENERAL PROCEDURES FOR THE TESTING OF CODED SUBSTANCES**

827 Range finder experiments are used to determine the concentrations of test substance to be used
 828 during comprehensive testing. Comprehensive testing experiments are used to determine whether
 829 a substance possesses ER agonist activity in the LUMI-CELL® ER assay.

830 General procedures for range finder and comprehensive are similar. For specific details (such as
831 plate layout) of range finder testing see **Section 12.0**. For specific details of comprehensive
832 testing, see **Section 13.0**.

833 **11.1 Application of Reference Standard, Controls, and Test Substances**

- 834 1. Remove seeded 96-well plates from the incubator, inspect them using an inverted
835 microscope. Only use plates in which the cells in all wells giving a score of 1
836 according to **Table 11-1**.
- 837 2. Remove medium by inverting the plate onto blotter paper. Gently tap plate against
838 the bench surface to remove residual liquid trapped in the wells.
- 839 3. Add 200 µL of reference standard, control, or test substance to each well (see
840 **Sections 12.0** and **13.0** for specific plate layouts).
- 841 4. Return plates to incubator and incubate (see **Section 9.0** for details) for 19 to 24
842 hours to allow maximal induction of luciferase activity in the cells.

843 **11.1.1 Preparation of Excel® Data Analysis Template For Range Finder Testing**

- 844 1. In Excel®, open a new “AgRFTemplate” and save it with the appropriate project
845 name as indicated in the NICEATM Style Guide.
- 846 2. Fill out the table at the top of the “Raw Data” worksheet with information
847 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
848 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
849 default template “AgRFTemplate” on a laboratory specific basis).
- 850 3. Add the following information regarding the assay to the “Compound Tracking”
851 worksheet.
 - 852 ▪ Plate # - Enter the experiment ID or plate number into cell E1
 - 853 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
854 experiment into cell B5
 - 855 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
856 Media in cells B6 and B7

- 857 ▪ Test Substance Code – Enter the test substance codes into cells C13 to
858 C18
- 859 ▪ Name: Enter the experimenter name into cell G6
- 860 ▪ Date: Enter the experiment date in the format day\month\year into cell
861 G10
- 862 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
863 plate contaminated)
- 864 4. Enter the following substance testing information to the “List” page:
- 865 ▪ Concentration – Type in the test substance concentration in µg/ml in
866 descending order.
- 867 ▪ Also add any replicate-specific comments on this page (e.g, spilled tube,
868 etc.), in the comments section
- 869 ▪ All of the remaining cells on the List tab should populate automatically.
- 870 ▪ The “Template”, “Compound Mixing” and “Visual Inspection” tabs
871 should automatically populate with the information entered into the
872 Compound Tracking and List tabs.
- 873 5. Save the newly named project file.
- 874 6. Print out either the “List” or “Template” page for help with dosing the 96-well
875 plate. Sign and date the print out and store in study notebook.

876 11.1.2 Preparation of Excel® Data Analysis Template for Comprehensive Testing

- 877 1. In Excel®, open a new “AgCTTemplate” and save it with the appropriate project
878 name as indicated in the NICEATM Style Guide.
- 879 2. Fill out the table at the top of the “Raw Data” worksheet with information
880 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
881 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
882 default template “AgCTTemplate” on a laboratory specific basis).

883

- 884 3. On the “Compound Tracking” tab, enter the following information:
- 885 ▪ Plate # - Enter the experiment ID or plate number into cell E1
- 886 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
- 887 experiment into cell C5
- 888 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
- 889 Media in cells C6 and C7
- 890 ▪ Test Substance Code – Enter the test substance codes into cells C15 and
- 891 C16. Enter the test substance dilution into cells E25 and E26.
- 892 ▪ Name: Enter the experimenter name into cell G6
- 893 ▪ Date: Enter the experiment date in the format day\month\year into cell
- 894 G10
- 895 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
- 896 plate contaminated)
- 897 4. Enter substance testing concentrations to the “List” page. Also add any replicate-
- 898 specific comments on this page (e.g, spilled tube, etc.).
- 899 5. Save the newly named project file.
- 900 6. Print out either the “List” or “Template” page for help with dosing the 96-well
- 901 plate. Sign and date the print out and store in study notebook.

902 11.2 **Visual Evaluation of Cell Viability**

- 903 1. 19 to 24 hours after dosing the plate, remove the plate from the incubator and
- 904 remove the media from the wells by inverting the plate onto blotter paper. Gently
- 905 tap plate against the bench surface to remove residual liquid trapped in the wells.
- 906 2. Use a repeat pipetter to add 50 µL 1X PBS to all wells. Immediately remove PBS
- 907 by inversion.
- 908 3. Using an inverted microscope, inspect all of the wells used in the 96-well plate
- 909 and record the visual observations using the scores in **Table 11-1**.

910 **Table 11-1 Visual Observation Scoring**

Viability Score	Brief Description ¹
1	Normal Cell Morphology and Cell Density
2	Altered Cell Morphology and/or Small Gaps between Cells
3	Altered Cell Morphology and/or Large Gaps between Cells
4	Few (or no) Visible Cells
P	Wells containing precipitation are to be noted with “P”

911 ¹Reference photomicrographs are provided in the LUMI-CELL® ER Validation Study “Visual Observation Cell
 912 Viability Manual.”
 913

914 **11.3 Lysis of Cells for LUMI-CELL® ER**

- 915 1. Apply the reflective white backing tape to the bottom of the 96-well plate (this
 916 will increase the effectiveness of the luminometer).
- 917 2. Add 30 µL 1X lysis reagent to the assay wells and place the 96-well plate on an
 918 orbital shaker for one minute.
- 919 3. Remove plate from shaker and measure luminescence (as described in **Section**
 920 **11.4**).

921 **11.4 Measurement of Luminescence**

922 Luminescence is measured in the range of 300 to 650 nm, using an injecting luminometer and
 923 with software that controls the injection volume and measurement interval. Light emission from
 924 each well is expressed as RLU per well. The luminometer output is saved as raw data in an
 925 Excel® spreadsheet. A hard copy of the luminometer raw data should be signed, dated and stored
 926 in the study notebook.

927 **11.5 Data Analysis**

928 LUMI-CELL® ER uses an Excel® spreadsheet to collect and adjust the RLU values obtained
 929 from the luminometer and a GraphPad Prism® template to analyze and graph data. The Excel®
 930 spreadsheet subtracts background luminescence (average DMSO solvent control RLU value)
 931 from test substance, reference standard and control RLU values. Plate induction is calculated
 932 using these corrected RLU values. Test substance, reference standard, and control RLU values
 933 are then adjusted relative to the highest E2 reference standard RLU value, which is set to 10,000.
 934 After adjustment, values are transferred to GraphPad Prism® for data analysis and graphing.

935 11.5.1 Collection and Adjustment of Luminometer Data for Range Finder Testing

936 The following steps describe the procedures required to populate the Excel® spreadsheet that has
937 been configured to collect and adjust the RLU values obtained from the luminometer.

- 938 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
939 from **Section 11.1.1**.
- 940 2. Copy the raw data using the Excel® copy function, then paste the copied data into
941 cell B19 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
942 the **Paste Special – Values** command. This position corresponds to position A1 in
943 the table labeled Table 1 in this tab.
- 944 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
945 whether there are any potential outliers. See **Section 11.6.2** for further explanation
946 of outlier determinations.
- 947 4. If an outlier is identified, perform the following steps to remove the outlier from
948 calculations:
- 949 ▪ correct the equation used to calculate DMSO background in Table 1 [e.g.,
950 if outlier is located in cell F26, adjust the calculation in cell H40 to read
951 =AVERAGE(G26:I26)]
 - 952 ▪ then correct the equation used to calculate the average DMSO value in
953 Table 2 [e.g., following the above example, adjust cell M42 to read
954 =AVERAGE(G26:I26)]
 - 955 ▪ then correct the equation used to calculate the standard deviation of the
956 DMSO value in Table 2 [e.g., following the above example, adjust cell
957 M43 to read =STDEV(G36:I36)]
- 958 5. Excel® will automatically subtract the background (the average DMSO control
959 value) from all of the RLU values in Table 1 and populate Table 2 with these
960 adjusted values.
- 961 6. To calculate plate induction, identify the cell containing the E2a replicate in Table
962 1, plate row H that has the highest RLU value (i.e., cell B26, C26, D26, or E26).

- 963 7. Click into cell D14 and enter the cell number from the previous step into the
964 numerator.
- 965 8. Identify the cell containing the E2b replicate in Table 1, plate row H that has the
966 highest RLU value (i.e., cell J26, K26, L26, or M26).
- 967 9. Click into cell E14 and enter the cell number from the previous step into the
968 numerator.
- 969 10. Click on the “ER Agonist Report” worksheet.
- 970 11. The data for the E2 reference standard, methoxychlor, and DMSO replicates
971 populate the left portion (columns A – F) of the spreadsheet. The data is
972 automatically placed in an Excel® graph.
- 973 12. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
974 D2 of “ER Agonist Report” tab and check the formula contained within that cell.
975 The divisor should be the cell number of the cell containing the highest Mean E2
976 RLU value ((i.e., cell A16, A17, A18, or A19).
- 977 13. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
978 scores for each well on the 96-well plate. This data will be linked to the “ER
979 Agonist Report” worksheet.
- 980 14. After the testing results have been evaluated and reviewed for quality control,
981 enter the following information into the Compound Tracking worksheet:
- 982 ▪ Enter pass/fail results for plate reference standard and control parameters
983 into the Plate Pass/Fail Table
 - 984 ▪ Enter information from the testing of coded substances into the Testing
985 Results Table
 - 986 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
987 data into cell A34
 - 988 ▪ Date – Enter the date on which the data was reviewed into cell D34
989

989 11.5.2 Collection and Adjustment of Luminometer Data for Comprehensive Testing

990 The following steps describe the procedures required to populate the Excel® spreadsheet that has
991 been configured to collect and adjust the RLU values obtained from the luminometer.

- 992 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
993 from **Section 11.1.2**.
- 994 2. Copy the raw data using the Excel® copy function, then paste the copied data into
995 cell B16 of the “RAW DATA” worksheet in the experimental Excel® spreadsheet
996 using the **Paste Special – Values** command. This position corresponds to position
997 A1 in the table labeled Table 1 in this worksheet.
- 998 3. Fill out the table at the top of the “Raw Data” worksheet with information
999 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
1000 Meas. Time/Well (s), etc. If desired, this information can be added to the
1001 Laboratory Template File.
- 1002 4. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
1003 whether there are any potential outliers. See **Section 11.6.2** for further explanation
1004 of outlier determinations.
- 1005 5. If an outlier is identified, perform the following steps to remove the outlier from
1006 calculations:
- 1007 ▪ correct the equation used to calculate DMSO background in Table 1[e.g.,
1008 if outlier is located in cell M17, adjust the calculation in cell H37 to read
1009 =AVERAGE(M16,M18:M19)]
 - 1010 ▪ then correct the equation used to calculate the DMSO mean and SD
1011 values [e.g., following the above example, adjust cell M39 to read
1012 =AVERAGE(M28,M30:M31), and adjust cell M40 to read
1013 =STDEV(M28,M30:M31)]
- 1014 6. Excel® will automatically subtract the background (the average DMSO control
1015 value) from all of the RLU values in Table 1 and populate Table 2 with these
1016 adjusted values.

- 1017 7. To calculate plate induction, identify the cell in containing the E2 replicate in
1018 Table 1, plate row G that has the highest RLU value.
- 1019 8. Click into cell D11 and enter the cell number from the previous step into the
1020 numerator.
- 1021 9. Identify the cell containing the E2 replicate in plate row H that has the highest
1022 RLU value.
- 1023 10. Click into cell E11 and enter the cell number from the previous step into the
1024 numerator.
- 1025 11. Open the “ER Agonist Report” worksheet.
- 1026 12. The data for the E2 reference standard, methoxychlor, and DMSO replicates
1027 populate the left portion (columns A – E) of the spreadsheet. The data is
1028 automatically placed in an Excel® graph.
- 1029 13. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
1030 E2 of “ER Agonist Report” tab and check the formula contained within that cell.
1031 The divisor should be the cell number of the cell containing the highest Avg E2
1032 RLU value (cells A16 through A26).
- 1033 14. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
1034 scores for each well on the 96-well plate. This data will be linked to the “ER
1035 Agonist Report” worksheet.
- 1036 15. Copy the data from the “ER Agonist Report” worksheet into GraphPad Prism® for
1037 the calculation of EC₅₀ values and to graph experimental results as indicated in the
1038 NICEATM Prism® Users Guide.
- 1039 16. After the testing results have been evaluated and reviewed for quality control,
1040 enter the following information into the Compound Tracking worksheet:
- 1041 ▪ Enter pass/fail results for plate reference standard and control parameters
1042 into the Plate Pass/Fail Table
 - 1043 ▪ Enter information from the testing of coded substances into the Testing
1044 Results Table

- 1045 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
- 1046 data into cell A32
- 1047 ▪ Date – Enter the date on which the data was reviewed into cell D32

1048 11.5.3 Determination of Outliers

1049 The Study Director will use good statistical judgment for determining “unusable” wells that will
 1050 be excluded from the data analysis and will provide an explanation in the study notebook for any
 1051 excluded data. This judgment for data acceptance will include Q-test analysis.

1052 The formula for the Q test is:

1053
$$\frac{\text{Outlier} - \text{Nearest Neighbor}}{\text{Range (Highest} - \text{Lowest)}}$$

1054 where the outlier is the value proposed for exclusion, the nearest neighbor is the value closest to
 1055 the outlier, and the range is the range of the three values (Q values for samples sizes from 3 to 10
 1056 are provided in **Table 11-2**). For example, if the value of this ratio is greater than 0.94 (the Q
 1057 value for the 90% confidence interval for a sample size of three) or 0.76 (the Q value for the 90%
 1058 confidence interval for a sample size of four), the outlier may be excluded from data analysis.

1059 **Table 11-2 Q Test Values**

Number Of Observations	Q Value
2	-
3	0.94
4	0.76
5	0.64
6	0.56
7	0.51
8	0.47
9	0.44
10	0.41

1060

1061 For E2 reference standard replicates (sample size of two), any adjusted RLU value for a replicate
 1062 at a given concentration of E2 is considered and outlier if its value is more than 20% above or
 1063 below the adjusted RLU value for that concentration in the historical database.

1064 11.5.4 Acceptance Criteria1065 11.5.4.1 *Range Finder Testing*

1066 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1067 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1068 (QC) for these parameters derived from the historical database, which are summarized below.

1069 • Induction: Plate induction, as measured by dividing the averaged highest E2
1070 reference standard RLU value by the averaged DMSO control RLU value, must
1071 be greater than three-fold.

1072 • DMSO control results: Solvent control RLU values must be within 2.5 times the
1073 standard deviation of the historical solvent control mean RLU value.

1074 An experiment that fails either acceptance criterion will be discarded and repeated.

1075 11.5.4.2 *Comprehensive testing*

1076 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1077 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1078 (QC) for these parameters derived from the historical database, which are summarized below.

1079 • Induction: Plate induction, as measured by dividing the averaged highest E2
1080 reference standard RLU value by the averaged DMSO control RLU value, must
1081 be greater than three-fold.

1082 • Reference standard results: The E2 reference standard concentration-response
1083 curve should be sigmoidal in shape and have at least three values within the linear
1084 portion of the concentration-response curve.

1085 • DMSO control results: DMSO control RLU values must be within 2.5 times the
1086 standard deviation of the historical solvent control mean RLU value.

1087 • Positive control results: Methoxychlor control RLU values must be above the line
1088 representing the DMSO mean plus three times the standard deviation from the
1089 DMSO mean.


1090 An experiment that fails any single acceptance criterion will be discarded and repeated.

1091

1092 **12.0 RANGE FINDER TESTING**

1093 Agonist range finding for coded substances consists of seven point, 1:10 serial dilutions, with
 1094 each concentration tested in duplicate wells of the 96-well plate. **Figure 12-1** contains a template
 1095 for the plate layout to be used in agonist range finder testing.

1096 **Figure 12-1 Agonist Range Finder Test Plate Layout**



- **Four Point E2 Reference Standard**
- **DMSO (Solvent Control)**
- **Range Finder for Sample #1**
- **Range Finder for Sample #2**
- **Range Finder for Sample #3**
- **Range Finder for Sample #4**
- **Range Finder for Sample #5**
- **Range Finder for Sample #6**

1097

1098 Evaluate whether range finder experiments have met the acceptance criteria
 1099 (see **Section 11.5.4.1**).


1100 To determine starting concentrations for comprehensive testing use the following criteria:

- 1101
- 1102
- 1103
- 1104
- 1105
- 1106
- 1107
- 1108
- 1109
- 1110
- 1111
- 1112
- 1113
- 1114
- 1115
- 1116
- 1117
- 1118
- 1119
- 1120
- 1121
- 1122
- 1123
- 1124
- 1125
- 1126
- 1127
- 1128
- If results in the range finder test suggest that the test substance is negative for agonist activity (i.e., if there are no points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control, see **Figure 12-2**), comprehensive testing will be conducted using an 11 point 1:2 serial dilution with the limit concentration as the starting concentration (i.e., the maximum soluble concentration in the range finder).
 - If results in the range finder test suggest that the test substance is negative for agonist activity (i.e., if there are no points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control), and the higher concentrations in the range finder are cytotoxic, comprehensive testing will be conducted using an 11 point 1:2 serial dilution with the lowest cytotoxic concentration as the starting concentration (see **Figure 12-3**).
 - If results in the range finder test suggest that the test substance is positive for agonist activity (i.e., if there are points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control), the starting concentration to be used for the 11-point dilution scheme in comprehensive testing should be one log higher than the concentration giving the highest adjusted RLU value in the range finder. The 11-point dilution scheme will be based on either 1:2 or 1:5 dilutions according to the following criteria:
 - An 11-point 1:2 serial dilution should be used if the resulting concentration range (note: an 11-point 1:2 serial dilution will cover a range of concentrations over approximately three orders of magnitude [three logs]) will encompass the full range of responses based on the concentration response curve generated in the range finder test (see **Figure 12-4**).
 - If the concentration range that would be generated with the 1:2 serial dilution will not encompass the full range of responses based on the

1129 concentration response curve in the range finder test (see **Figure 12-5** and
 1130 **12-6**), an 11-point 1:5 serial dilution should be used instead.

- 1131 • If a substance exhibits a biphasic concentration response curve in the range finder
 1132 test, both phases should also be resolved in comprehensive testing. In order to
 1133 resolve both curves, the starting concentration should be based on the peak
 1134 associated with the higher concentration and should be one log higher than the
 1135 concentration giving the highest adjusted RLU value in the range finder. As an
 1136 example, an 11-point 1:5 serial dilution should be used based on the range finder
 1137 results presented in **Figure12-7**.

1138 **Figure 12-2 Agonist Range Finder (example 1)**




1139 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.
 1140

1141

1142


1142 **Figure 12-3 Agonist Range Finder (example 2)**



1143

1144 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.


1145 **Figure 12-4 Agonist Range Finder (example 3)**



1146

1147 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.


1148 **Figure 12-5 Agonist Range Finder (example 4)**



1149

1150 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.

1151 **Figure 12-6 Agonist Range Finder (example 5)**




1152

1153 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.

1154

1154 **Figure 12-7 Agonist Range Finder (example 6)**



1155
 1156 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.
 1157

1158 **13.0 COMPREHENSIVE TESTING**

1159 Agonist comprehensive testing for coded substances consists of 11 point serial dilutions (either
 1160 1:2 or 1:5 serial dilutions based on the starting concentration for comprehensive testing criteria in
 1161 **Section 12.0**) with each concentration tested in triplicate wells of the 96-well plate. **Figure 13-1**
 1162 contains a template for the plate layout to be used in agonist comprehensive testing.
 1163

1163 **Figure 13-1 Agonist Comprehensive Test Plate Layout**

1164

1165

1166

1167

1168

1169

1170


1171

1172

1173

1174

1175



1176

 **11 Point Duplicate E2 Reference Standard**

1177

 **DMSO (Solvent Control)**

1178

 **Test Substance #1**

1179

1180

 **Test Substance #2**

1181

 **Methoxychlor Control**

1182

1183 Evaluate whether comprehensive experiments have met acceptance criteria (see **Section 11.6.4**)

1184 and graph the data as described in the NICEATM Prism® users guide.

1185 • If the substance has been tested up to the limit dose or the maximum soluble dose,

1186 without causing a significant decrease in cell viability, and there are no points on

1187 the concentration curve that are greater than the mean plus three times the

1188 standard deviation of the DMSO control, the substance is considered negative for

1189 agonism

1190 • If the substance has a positive response (See **Section 6.0**) at any concentration,

1191 the substance is considered positive for agonism.

1192 14.0 COMPILATION OF THE HISTORICAL QUALITY CONTROL DATABASE

1193 Historical databases are maintained in order to ensure that the assay is functioning properly.

1194 Historical databases are compiled using Excel® spreadsheets and are separate from the
1195 spreadsheets used to collect the data for individual test plates. Reference standard and control
1196 data are used to develop and maintain the historical database and are used as quality controls to
1197 determine acceptance of individual test plates.

1198 The sources of the data needed to compile the historical database for the DMSO control are the
1199 experiment specific Excel® data collection and analysis spreadsheets used for LUMI-CELL® ER
1200 agonist and antagonist testing (see **Section 11.5.2** of the LUMI-CELL® ER agonist protocol and
1201 **Section 13.5.2** in the LUMI-CELL® ER antagonist protocol).

1202 14.1 DMSO Control

1203 Open the combined agonist and antagonist LUMI-CELL® ER historical database Excel®
1204 spreadsheet (LUMI_AgandAntQC.xls) and save under a new name using the Excel® “Save As”
1205 function, adding the laboratory designator to the file name (e.g., for Laboratory H, the new name
1206 would be HLUMI_AgandAntQC.xls). Enter the date and experiment name into worksheet
1207 columns A and B respectively. Enter the experimental mean DMSO control value (from cell H37
1208 in the RAW DATA worksheet of the agonist and antagonist Excel® data collection and analysis
1209 spreadsheet) into worksheet column C. Acceptance or rejection of the plate DMSO control data
1210 for range finding and comprehensive testing is based on whether the mean plate DMSO RLU
1211 value falls within 2.5 times the standard deviation of the DMSO value in the historical database
1212 (columns G and H in the DMSO worksheet).

1213

1214 15.0 QUALITY TESTING OF MATERIALS

1215 All information pertaining to the preparation and testing of media, media supplements, and other
1216 materials should be recorded in the Study Notebook.

1217

1217 **15.1 Tissue Culture Media**

1218 Each lot of tissue culture medium must be tested in a single growth flask of cells before use in
1219 ongoing tissue culture or experimentation (**note:** each bottle within a given lot of
1220 Charcoal/Dextran treated FBS must be tested separately).

- 1221 1. Every new lot of media (RPMI and DMEM) and media components (FBS,
1222 Charcoal/Dextran treated FBS, and L-glutamine) must first be tested on the
1223 LUMI-CELL® ER assay prior to being used in any GLP acceptable assays.
- 1224 2. Add 4 µL of DMSO (previously tested) into four separate 13 mm tubes.
- 1225 3. Add 400 µL media (to be tested) to the same tubes.
- 1226 4. Dose an experimental plate as in **Section 12.0**, treating the media being tested as a
1227 test substance.
- 1228 5. Analyze 96-well plate as described in **Section 12.0**, comparing the data from the
1229 DMSO controls made using previously tested tissue culture media to the new
1230 media being tested.
- 1231 6. Use the agonist historical database to determine if the new media with DMSO lies
1232 within 2.5 standard deviations of the mean for the media. If the RLU values for
1233 the new media with DMSO lie within 2.5 standard deviation of the mean for the
1234 historical data on DMSO, the new lot of media is acceptable. If the RLU values
1235 for the new media with DMSO do not lie within 2.5 standard deviations of the
1236 DMSO mean from historical database, the new lot may not be used in the assay.
- 1237 7. Note date and lot number in study notebook.
- 1238 8. If the new bottle passes quality testing as described in **Section 15.1 step 6**, apply
1239 the media to a single flask of cells and observe cell growth and morphology over
1240 the following 2 – 3 days. If there is no change in growth or morphology, the new
1241 media is acceptable for use.

1242

1242 **15.2 G418:**

- 1243 1. New lots of G418 must first be tested on the LUMI-CELL® ER assay prior to
1244 being used in any GLP acceptable assays.
- 1245 2. Add 220 µL of G418 (previously tested) to a single flask containing cells growing
1246 in RPMI.
- 1247 3. Add 220 µL of G418 (to be tested) to a different flask containing cells growing in
1248 RPMI.
- 1249 4. Observe cellular growth and morphology in both tissue culture flasks over a 48 to
1250 72 hour period. If there are no differences in observed growth rate and
1251 morphology between the two flasks, the new G418 lot is acceptable.
- 1252 5. If cellular growth is decreased, or the cells exhibit abnormal morphology, the new
1253 lot of G418 is not acceptable.
- 1254 6. Note date and lot number in study book.

1255 **15.3 DMSO**

- 1256 1. Every new bottle of DMSO must be tested on the LUMI-CELL® ER assay prior
1257 to use in any GLP acceptable assays.
- 1258 2. Add 4 µL of DMSO (to be tested) into four separate 13 mm tubes.
- 1259 3. Add 400 µL media (previously tested) to the same tubes.
- 1260 4. Dose an experimental plate as in **Section 12.0**, treating the DMSO containing
1261 media being tested as a test substance.
- 1262 5. Analyze 96-well plate as described in **Section 12.0**, comparing the data from the
1263 DMSO controls made using previously tested tissue culture media.
- 1264 6. Use the agonist historical database to determine if media with new DMSO lies
1265 within 2.5 standard deviations of the DMSO mean from historical database. If the
1266 RLU values for the media with new DMSO lie within 2.5 standard deviations of
1267 the DMSO mean from the historical database, the new lot of DMSO is acceptable.
1268 If the RLU values for media with new DMSO do not lie within 2.5 standard

1269 deviations of the DMSO mean from historical database, the new lot may not be
1270 used in the assay.

1271 7. Note the date, lot number, and bottle number in study book.

1272 8. If no DMSO has been previously tested, test several bottles as described in
1273 **Section 15.3**, and determine whether any of the bottles of DMSO have a lower
1274 average RLU than the other bottle(s) tested. Use the DMSO with the lowest
1275 average RLU for official experiments.

1276 **15.4 Plastic Tissue Culture Materials**

1277 1. Grow one set of cells, plate them for experiments on plastic ware from the new lot
1278 and one set of cells in the plastic ware from a previous lot, and dose them with E2
1279 reference standard and controls.

1280 2. Perform the LUMI-CELL® ER experiment with both sets of cells.

1281 3. If all of the analysis falls within acceptable QC criteria, then the new
1282 manufacturer's products may be used.

1283

1283 **16.0 REFERENCES**

- 1284 Eli Lilly and Company and National Institutes of Health Chemical Genomics Center. 2005.
1285 Assay Guidance Manual Version 4.1. Bethesda, MD: National Institutes of Health. Available:
1286 http://www.ncgc.nih.gov/guidance/manual_toc.html [accessed 05 September 2006]
- 1287 ICCVAM. 2001. Guidance Document on Using *In Vitro* Data to Estimate *In Vivo* Starting Doses
1288 for Acute Toxicity. NIH Pub. No. 01-4500. Research Triangle Park, NC: National Institute of
1289 Environmental Health Sciences. Available: [http://iccvam.niehs.nih.gov/methods/invidocs/
1290 guidance/iv_guide.pdf](http://iccvam.niehs.nih.gov/methods/invidocs/guidance/iv_guide.pdf) [accessed 31 August 2006]