

# Screening group in Protein Chemistry Laboratory/NCI-Frederick

Robert J. Fisher, Head PCL

## Criteria for NCp7 Small Molecule Screen

- Novel high affinity, zinc finger dependent, nucleic acid binding activity of NCp7 provides an assay for a new class of inhibitors
- The objective of this work was to find reversible inhibitors of the high affinity interaction between NCp7 and nucleic acid


P7 screening assay  
(Andy Stephen 2-01)

Bind p7 to plate (costar high-bind polystyrene Cat# 3631 or 3601)  
Dilute p7 (200-1000nM) in 1xPBS/5mM b-mercaptoethanol  
Add 100ml to each well, allow to sit overnight at 4<sup>0</sup>C.  
Add 200ml of 2% BSA in 1xPBS/5mM b-mercaptoethanol, incubate at room temperature for 1-2hours.  
Aspirate off and wash with 2 x 200ml 1xPBS/5mM b-mercaptoethanol.  
The plates can be stored at -20<sup>0</sup>C, but more reproducible results are obtained if the plates are made fresh

Bind oligo to p7  
Dilute oligo (concentration depends on amount of p7 immobilized, approx 1mM) in 1xPBS/100mM b-mercaptoethanol/0.05% Tween/5-10% DMSO.  
Bind at room temperature for 1 hour, wash with 3 x 200ml 1xPBS/0.05% Tween

Detection.  
Add 100ml 1:20000 dilution of nutraavidin-HRP (stock 0.8mg/ml, from Pierce)  
Incubate at room temperature for 1 hour, wash with 3 x 200ml 1xPBS/0.05% Tween  
Add 100ml "supersignal" from Pierce and count in Wallac plate reader.

# The NCp7 primary assay



- p7NC immobilized to plate
- Block with BSA
- Wash, store at -20°C until use
- Add oligo (biotin-labeled)  $\pm$  10uM cpd
- Incubate for one hour
- Wash
- Incubate with streptavidin HRP
- Wash
- Develop
- read on Wallac Victor
- “hit threshold”=100% inhibition (yes/no)
- 8 positive, 8 negative control/plate
- background ~1% false positives