

**Non-radioactive Murine Local Lymph Node Assay: BrdU-ELISA Test Method Protocol
(LLNA: BrdU-ELISA)**

Draft Background Review Document

January 2008

[This Page Intentionally Left Blank]

Table of Contents

	Page Number
1	
2	
3	List of Tables iii
4	List of Abbreviations and Acronyms iv
5	Interagency Coordinating Committee on the Validation Of Alternative Methods
6	(ICCVAM) Designated Agency Representatives vi
7	Acknowledgements vii
8	Preface x
9	Executive Summary xii
10	1.0 Introduction 1
11	1.1 Historical Background 1
12	1.2 The LLNA: BrdU-ELISA 2
13	2.0 LLNA: BrdU-ELISA Test Method Protocol 2
14	3.0 LLNA: BrdU-ELISA Validation Database 3
15	4.0 Reference Data 6
16	5.0 Test Method Data and Results 7
17	6.0 Test Method Accuracy 7
18	6.1 Total LLNA: BrdU-ELISA Database Analysis Using $SI \geq 3$ Decision
19	Criteria 7
20	6.2 Accuracy Analysis ($SI \geq 3.0$) Based on Revised Draft ICCVAM
21	Performance Standards Reference Substances 10
22	6.3 Discordant Results for Accuracy Analysis of $SI \geq 3.0$ Decision Criterion..... 14
23	6.4 LLNA: BrdU-ELISA Accuracy Analysis Using Alternative Decision
24	Criteria 17

25	6.5	Discordant Results for Accuracy Analysis of Alternative Decision	
26		Criterion	18
27	7.0	Test Method Reliability	21
28	7.1	Intralaboratory Reproducibility – Qualitative Results	21
29	7.2	Intralaboratory Reproducibility – SI	22
30	7.3	Intralaboratory Reproducibility – EC3	23
31	8.0	Data Quality	25
32	9.0	Other Scientific Reports and Reviews	26
33	10.0	Animal Welfare Considerations	26
34	10.1	Rationale for the Need to Use Animals	27
35	10.2	Basis for Determining the Number of Animals Used	27
36	10.3	Reduction Considerations	27
37	11.0	Practical Considerations	27
38	11.1	Transferability of the LLNA: BrdU-ELISA	28
39	11.2	Facilities and Major Fixed Equipment Required to Conduct the	
40		LLNA: BrdU-ELISA	28
41	11.3	LLNA: BrdU-ELISA Training Considerations	28
42	12.0	References	29
43	Appendix A	LLNA: BrdU-ELISA Protocol	A-1
44	Appendix B	Analysis of Alternative Decision Criteria for Determination of	
45		Sensitizers	B-1
46	Appendix C	Physico-Chemical Properties of Substances Tested Using the	
47		LLNA: BrdU-ELISA	C-1
48	Appendix D	Comparative LLNA: BrdU-ELISA, Traditional LLNA,	
49		Guinea Pig Skin Sensitization, and Human Data	D-1

50	List of Tables		Page Number
51			
52	Table 3-1	Traditional LLNA EC3 Values and Chemical Classification of Substances	
53		Tested in the LLNA: BrdU-ELISA.....	5
54	Table 6-1	Evaluation of the Performance of the LLNA: BrdU-ELISA In Predicting	
55		Skin Sensitizing Potential Using Decision Criteria of $SI \geq 3.0$ to Identify	
56		Sensitizers.....	9
57	Table 6-2	Evaluation of the Performance of the LLNA: BrdU-ELISA ($SI \geq 3.0$)	
58		Using the Revised Draft ICCVAM Performance Standards Reference	
59		Substances	11
60	Table 6-3	Characteristics of the Substances Tested in the LLNA: BrdU-ELISA vs.	
61		the Revised Draft ICCVAM Performance Standards Substances List	14
62	Table 6-4	Discordant Results for LLNA: BrdU-ELISA (Using $SI \geq 3.0$ for Sensitizers)	
63		Compared to Traditional LLNA and Guinea Pig Reference Data.....	15
64	Table 6-5	Discordant Results for LLNA: BrdU-ELISA ($SI \geq 3.0$) when Compared	
65		to Traditional LLNA and Human Outcome Data	16
66	Table 6-6	Evaluation of the Performance of the LLNA: BrdU-ELISA In Predicting	
67		Skin Sensitizing Potential Using Alternative Decision Criteria to Identify	
68		Sensitizers.....	19
69	Table 6-7	Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision	
70		Criteria Compared to the Traditional LLNA.....	20
71	Table 7-1	Intralaboratory Reproducibility for the SI of Tested Substances in the	
72		LLNA: BrdU-ELISA - Coefficient of Variation	23
73	Table 7-2	Intralaboratory Reproducibility for the EC3 of Tested Substances in the	
74		LLNA: BrdU-ELISA - Coefficient of Variation	24
75	Table 7-3	Intralaboratory Reproducibility for the EC3 of Tested Substances in the	
76		Traditional LLNA.....	25
77			

78

List of Abbreviations and Acronyms

79	ACD	Allergic contact dermatitis
80	AOO	Acetone: olive oil
81	BRD	Background Review Document
82	BrdU	Bromodeoxyuridine
83	BT	Buehler Test
84	CASRN	Chemical Abstracts Service Registry Number
85	Conc.	Concentration tested
86	CPSC	U.S. Consumer Product Safety Commission
87	DMF	Dimethylformamide
88	DMSO	Dimethyl sulfoxide
89	EC3	Estimated concentration needed to produce a stimulation index
90		of three
91	ECVAM	European Centre for the Validation of Alternative Methods
92	ELISA	Enzyme-linked immunosorbent assay
93	EPA	U.S. Environmental Protection Agency
94	ESAC	ECVAM Scientific Advisory Committee
95	FDA	U.S. Food and Drug Administration
96	FR	<i>Federal Register</i>
97	GHS	United Nations Globally Harmonized System for the
98		Classification and Labelling of Chemicals
99	GLP	Good Laboratory Practice
100	GPMT	Guinea Pig Maximization Test
101	HCA	Hexyl cinnamic aldehyde
102	HMT	Human Maximization Test
103	HPTA	Human Patch Test Allergen
104	ICCVAM	Interagency Coordinating Committee on the Validation of
105		Alternative Methods
106	IR	Information requested
107	ISO	International Standards Organization
108	IWG	Immunotoxicity Working Group
109	Java	Japanese Center for the Validation of Alternative Methods
110	K _{ow}	Octanol-water partition coefficient
111	LNC	Lymph node cells
112	LLNA	Local Lymph Node Assay
113	LLNA: BrdU-ELISA	LLNA with enzyme-linked immunosorbent assay detection of
114		bromodeoxyuridine
115	MEK	Methyl ethyl ketone
116	MeSH	Medical Subject Headings
117	Min	Minimal
118	Mod	Moderate
119	NA	Not available
120	NC	Not calculated
121	NICEATM	National Toxicology Program Interagency Center for the
122		Evaluation of Alternative Toxicological Methods
123	NIEHS	National Institute of Environmental Health Sciences

124	NT	Not tested
125	NTP	National Toxicology Program
126	OECD	Organisation for Economic Co-operation and Development
127	OPPTS	Office of Prevention, Pesticides and Toxic Substances
128	Res	Result
129	SACATM	Scientific Advisory Committee on Alternative Toxicological
130		Methods
131	S.D.	Standard Deviation
132	SI	Stimulation Index
133	SLS	Sodium lauryl sulfate
134	TG	Test Guideline
135	U.S.	United States
136	Unk	Unknown
137	Veh.	Vehicle
138	vs.	Versus
139	w/v	Weight to volume ratio

140 **Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM)**
 141 **Designated Agency Representatives¹**

142	<p>Agency for Toxic Substances and Disease Registry</p> <ul style="list-style-type: none"> • Moiz Mumtaz, Ph.D. <p>Consumer Product Safety Commission</p> <ul style="list-style-type: none"> • Marilyn L. Wind, Ph.D. (Chair) ◊ Kristina Hatlelid, Ph.D. * Joanna Matheson, Ph.D. <p>Department of Agriculture</p> <ul style="list-style-type: none"> • Jodie Kulpa-Eddy, D.V.M. ◊ Elizabeth Goldentyer, D.V.M. <p>Department of Defense</p> <ul style="list-style-type: none"> • Robert E. Foster, Ph.D. ◊ Patty Decot * Peter J. Schultheiss, D.V.M., D.A.C.L.A.M. * Harry Salem, Ph.D. <p>Department of Energy</p> <ul style="list-style-type: none"> • Michael Kuperberg, Ph.D. ◊ Marvin Stodolsky, Ph.D. <p>Department of the Interior</p> <ul style="list-style-type: none"> • Barnett A. Rattner, Ph.D. ◊ Sarah Gerould, Ph.D. <p>Department of Transportation</p> <ul style="list-style-type: none"> • George Cushmac, Ph.D. ◊ Steve Hwang, Ph.D. <p>Environmental Protection Agency</p> <p>Office of Science Coordination and Policy</p> <ul style="list-style-type: none"> • Karen Hamernik, Ph.D. <p>Office of Research and Development</p> <ul style="list-style-type: none"> ◊ Julian Preston, Ph.D. * Suzanne McMaster, Ph.D. <p>OECD Test Guidelines Program</p> <ul style="list-style-type: none"> * Jerry Smrcek, Ph.D. <p>Office of Pesticides Programs</p> <ul style="list-style-type: none"> * Amy Rispin, Ph.D. * Deborah McCall 	<p>Food and Drug Administration</p> <p>Office of Science</p> <ul style="list-style-type: none"> • Suzanne Fitzpatrick, Ph.D., D.A.B.T. <p>Center for Drug Evaluation and Research</p> <ul style="list-style-type: none"> ◊ Abigail C. Jacobs, Ph.D. <p>Center for Devices and Radiological Health</p> <ul style="list-style-type: none"> * Melvin E. Stratmeyer, Ph.D. <p>Center for Biologics Evaluation and Research</p> <ul style="list-style-type: none"> * Richard McFarland, Ph.D., M.D. * Ying Huang, Ph.D. <p>Center for Food Safety and Nutrition</p> <ul style="list-style-type: none"> * David G. Hattan, Ph.D. * Robert L. Bronaugh, Ph.D. <p>Center for Veterinary Medicine</p> <ul style="list-style-type: none"> * Devaraya Jagannath, Ph.D. * M. Cecilia Aguila, D.V.M. <p>National Center for Toxicological Research</p> <ul style="list-style-type: none"> * William T. Allaben, Ph.D. * Paul Howard Ph.D. <p>Office of Regulatory Affairs</p> <ul style="list-style-type: none"> * Lawrence A. D’Hoostelaere, Ph.D. <p>National Cancer Institute</p> <ul style="list-style-type: none"> • Alan Poland, M.D. ◊ T. Kevin Howcroft, Ph.D. <p>National Institute of Environmental Health Sciences</p> <ul style="list-style-type: none"> • William S. Stokes, D.V.M., D.A.C.L.A.M. ◊ Raymond R. Tice, Ph.D. * Rajendra S. Chhabra, Ph.D., D.A.B.T * Jerrold J. Heindel, Ph.D. <p>National Institute for Occupational Safety and Health</p> <ul style="list-style-type: none"> • Paul Nicolaysen, V.M.D. ◊ K. Murali Rao, M.D., Ph.D. <p>National Institutes of Health</p> <ul style="list-style-type: none"> • Margaret D. Snyder, Ph.D. <p>National Library of Medicine</p> <ul style="list-style-type: none"> ◊ Jeanne Goshorn, M.S. <p>Occupational Safety and Health Administration</p> <ul style="list-style-type: none"> • Surender Ahir, Ph.D.
-----	---	--

-
- Principal Agency Representative
 - ◊ Alternate Principal Agency Representative
 - * Other Designated Agency Representative

¹ Roster as of January 2008.

142
143
144
145
146
147
148

Acknowledgements

The following individuals are acknowledged for their contributions to the LLNA test method review process

Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM) Immunotoxicity Working Group (IWG)

U.S. Consumer Product Safety Commission

Marilyn Wind, Ph.D.
Joanna Matheson, Ph.D. (IWG Co-Chair)

U.S. Environmental Protection Agency

Karen Hamernik, Ph.D.
Masih Hashim, Ph.D.
Marianne Lewis
Debbie McCall
Timothy McMahon, Ph.D.
Amy Rispin, Ph.D.
MaryJane Selgrade, Ph.D.
Marsha Ward, Ph.D.
Ronald E. Ward, Ph.D.

U.S. Food and Drug Administration

Ruth Barrett, Ph.D., D.V.M.
Paul C. Brown, Ph.D.
Abigail Jacobs, Ph.D. (IWG Co-Chair)
Dan Lyle, Ph.D.
Jiaqin Yao, Ph.D.

National Institute of Environmental Health Sciences

Dori Germolec, Ph.D.
William S. Stokes, D.V.M., D.A.C.L.A.M.
Raymond R. Tice, Ph.D.

National Institute for Occupational Health & Safety

Jean Meade, D.V.M., Ph.D.

ECVAM Liaison

Silvia Casati, Ph.D.

JaCVAM Liaison

Hajime Kojima, Ph.D.

149
150

150 **National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative**
151 **Toxicological Methods (NICEATM)**

152 **National Institute of Environmental Health Sciences**

William Stokes, D.V.M., D.A.C.L.A.M.	Director; Project Officer
Raymond Tice, Ph.D.	Deputy Director
Deborah McCarley	Special Assistant; Asst. Project Officer

153

154 **NICEATM Support Contract Staff (Integrated Laboratory Systems [ILS], Inc.)**

David Allen, Ph.D.	Principal Investigator
Douglas Winters, M.S.	Project Manager
Neepa Choksi, Ph.D.	Senior Staff Toxicologist
Judy Strickland, Ph.D., D.A.B.T.	Senior Staff Toxicologist
Frank Deal, M.S.	Staff Toxicologist
Elizabeth Lipscomb, Ph.D.	Staff Toxicologist
Eleni Salicru, Ph.D.	Staff Toxicologist
Thomas Burns, M.S.	Senior Project Coordinator/Technical Writer
Michael Paris	Senior Project Coordinator/Technical Writer
Patricia Ceger, M.S.	Project Coordinator/Technical Writer
James Truax, M.A.	Project Coordinator/Technical Writer
Catherine Sprankle	Communications Specialist/Web Developer
Linda Litchfield	Meeting Planner and Coordinator

155

Other Acknowledgements

156 *ICCVAM and NICEATM gratefully acknowledge Masahiro Takeyohsi, Ph.D., of the Chemicals*
157 *Evaluation and Research Institute (CERI-Japan), Saitama, Japan, for submitting data to*
158 *NICEATM used for the evaluation of the LLNA: BrdU-ELISA.*
159

160 Preface

161 In 1998, the Interagency Coordinating Committee on the Validation of Alternative Methods
162 (ICCVAM) in conjunction with the National Toxicology Program (NTP) Interagency Center for
163 the Evaluation of Alternative Toxicological Methods (NICEATM) evaluated the validation status
164 of the murine local lymph node assay (LLNA) as an alternative to guinea pig test methods for
165 assessing the allergic contact dermatitis (ACD) potential of substances. As described in the 1999
166 ICCVAM evaluation report², ICCVAM recommended that the LLNA could be used as a valid
167 substitute for the accepted guinea pig test methods, in most ACD testing situations.

168 Based on the ICCVAM recommendations, the ICCVAM member agencies that require the
169 regulatory submission of ACD data accepted the LLNA, with identified limitations, as an
170 alternative to guinea pig tests for assessing ACD. In 2002, the LLNA was adopted as Test
171 Guideline 429 by the 30-member countries of the Organisation for Economic Co-operation and
172 Development (OECD)³.

173 On January 10, 2007, the U.S. Consumer Product Safety Commission (CPSC) formally
174 nominated several activities related to the LLNA for evaluation by ICCVAM and NICEATM⁴.
175 One of the nominated activities was an assessment of the validation status of non-radioactive
176 alternatives to the current version of the LLNA, which uses radioactivity. After considering
177 comments from the public and the Scientific Advisory Committee on Alternative Toxicological
178 Methods (SACATM) on this nomination, ICCVAM assigned it a high priority, and directed
179 NICEATM and the ICCVAM Immunotoxicity Working Group (IWG) to conduct a review of the
180 current literature and an evaluation of the available data. The information described in this
181 background review document (BRD) was compiled by ICCVAM in response to this nomination.
182 ICCVAM and its IWG developed draft test method recommendations based on this evaluation.
183 An independent peer review panel (Panel) is being convened to peer review the BRD and to
184 evaluate the extent to which the information contained in the BRD support the draft

² ICCVAM 1999. The murine local lymph node assay: A test method for assessing the allergic contact dermatitis potential of chemical/compounds. NIH Publication No. 99-4494. Research Triangle Park, NC: National Toxicology Program (available at

http://iccvam.niehs.nih.gov/docs/immunotox_docs/llna/llnarep.pdf

³ OECD. 2002. Test guideline 429. Skin Sensitisation: Local Lymph Node Assay, adopted April 24, 2002. In: OECD Guidelines for Testing of Chemicals. Paris:OECD (available at http://www.oecd.org/document/55/0,2340,en_2649_34377_2349687_1_1_1_1,00.html)

⁴ Available at http://iccvam.niehs.nih.gov/methods/immunotox/llnadocs/CPSC_LLNA_nom.pdf

185 recommendations. ICCVAM will consider the conclusions and recommendations of the Panel,
186 along with comments received from the public and SACATM, when developing a final BRD and
187 final recommendations on the usefulness and limitations of each non-radioactive alternative
188 LLNA test method that is being considered.

189 We gratefully acknowledge the organizations and scientists who provided data and information
190 for this document. We would also like to recognize the efforts of the individuals who contributed
191 to the preparation of this BRD. These include David Allen, Ph.D., Thomas Burns, M.S., Neepa
192 Choksi, Ph.D., Michael Paris, Eleni Salicru, Ph.D., Catherine Sprankle, Judy Strickland, Ph.D.,
193 and Doug Winters, M.S., of Integrated Laboratory Systems, Inc., the NICEATM Support
194 Contractor, as well as the members of the ICCVAM IWG and the ICCVAM representatives who
195 subsequently reviewed and provided comments throughout the process leading to this final draft
196 version. We also want to thank Raymond Tice, Ph.D., Deputy Director of NICEATM, for his
197 contributions to this project. Finally, we want to recognize the excellent leadership of the IWG
198 Co-chairs, Abigail Jacobs, Ph.D. (FDA) and Joanna Matheson, Ph.D. (CPSC).

199 Marilyn Wind, Ph.D.
200 Deputy Associate Executive Director
201 Directorate for Health Sciences
202 U.S. Consumer Product Safety Commission
203 Chair, ICCVAM

204
205 William S. Stokes, D.V.M., D.A.C.L.A.M.
206 Rear Admiral, U.S. Public Health Service
207 Director, NICEATM
208 Executive Director, ICCVAM
209 January 7, 2008

210

211

212 **Executive Summary**

213 In 1999, the Interagency Coordinating Committee on the Validation of Alternative Methods
214 (ICCVAM) recommended to U.S. Federal agencies that the murine local lymph node assay
215 (LLNA) is a valid substitute for currently accepted guinea pig test methods to assess the allergic
216 contact dermatitis (ACD) potential of many, but not all, types of substances. The
217 recommendation was based on a comprehensive evaluation that included an independent
218 scientific peer review panel (Panel) assessment of the validation status of the LLNA. The Panel
219 report and the ICCVAM recommendations (ICCVAM 1999) are available at the National
220 Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative Toxicological
221 Methods (NICEATM)/ICCVAM website
222 (http://iccvam.niehs.nih.gov/docs/immunotox_docs/llna/llnarep.pdf). The LLNA was
223 subsequently incorporated into national and international test guidelines for the assessment of
224 skin sensitization (Organisation for Economic Co-operation and Development [OECD] Test
225 Guideline 429 [OECD 2002]; International Standards Organization [ISO] 10993-10: Tests for
226 Irritation and Sensitization [ISO 2002]; U.S. Environmental Protection Agency [EPA] Health
227 Effect Testing Guidelines on Skin Sensitization [EPA 2003]).

228 On January 10, 2007, the U.S. Consumer Product Safety Commission (CPSC) formally
229 nominated several activities related to the LLNA for evaluation by ICCVAM and NICEATM
230 (Available at http://iccvam.niehs.nih.gov/methods/immunotox/llnadocs/CPSC_LLNA_nom.pdf).
231 One of the nominated activities was an assessment of the validation status of non-radioactive
232 alternatives to the current version of the LLNA ([ICCVAM 1999; Dean et al. 2001] referred to
233 hereafter as the “traditional LLNA”), which uses radioactivity to detect sensitizers. The
234 information described in this background review document (BRD) was compiled by ICCVAM
235 and NICEATM in response to this nomination. The BRD provides a comprehensive review of
236 available data and information regarding the usefulness and limitations of one of these methods,
237 the LLNA with detection of bromodeoxyuridine (BrdU) incorporation by an enzyme-linked
238 immunosorbent assay (ELISA) (referred to hereafter as the “LLNA: BrdU-ELISA”).

239 The LLNA: BrdU-ELISA was developed by Takeyoshi et al. (2001). While the traditional
240 LLNA assesses cell proliferation by measuring the incorporation of radioactivity into the
241 deoxyribonucleic acid (DNA) of dividing lymph node cells, the LLNA: BrdU-ELISA assesses

242 cell proliferation by measuring the incorporation of the thymidine analog BrdU into the DNA of
243 dividing lymphocytes using ELISA. A Stimulation Index (SI), the ratio of the mean BrdU
244 incorporation into the lymph nodes of mice in the test substance group to the mean BrdU
245 incorporation into the lymph nodes of mice in the vehicle group, greater than three identifies a
246 substance as a sensitizer. Other than the procedure for measuring lymph node cell proliferation,
247 the protocol for the LLNA: BrdU-ELISA is similar to that of the traditional LLNA (Dean et al.
248 2001; ICCVAM 1999).

249 The accuracy and reliability of the LLNA: BrdU-ELISA were assessed using data for 24
250 substances generated from six published studies (Takeyoshi et al. 2003; 2004a; 2004b; 2005;
251 2006; 2007a) and one platform presentation (Takeyoshi 2007b). The reference test data for these
252 substances were obtained from the traditional LLNA, guinea pig (GP) skin sensitization tests,
253 and/or human skin sensitization tests. Of the 24 substances with traditional LLNA data, 16 were
254 classified by the traditional LLNA as skin sensitizers and eight were classified as non-sensitizers.

255 When the LLNA: BrdU-ELISA was compared to the traditional LLNA, accuracy was 75%
256 (18/24), sensitivity was 71% (12/17), specificity was 86% (6/7), the false positive rate was 14%
257 (1/7), and false negative rate was 29% (5/17). Using the traditional LLNA as the reference
258 classification, five non-sensitizers and one sensitizer were not classified correctly. No
259 commonalities were identified among these substances.

260 The LLNA: BrdU-ELISA results included eight of the 18 minimum substances proposed in the
261 *Revised Draft ICCVAM Murine Local Lymph Node Assay Performance Standards* (ICCVAM
262 2007); there were seven sensitizers and one non-sensitizer. The sensitizer/non-sensitizer outcome
263 of the LLNA: BrdU-ELISA was consistent with the outcome of the traditional LLNA with the
264 exception of one sensitizer.

265 When the decision criteria were altered to include SI values below three to identify a positive
266 response, improved performance was achieved. Best overall performance was achieved using an
267 $SI \geq 1.3$ with an accuracy of 96% (22/23), sensitivity of 100% (17/17), specificity of 83% (5/6),
268 a false positive rate of 17% (1/6), and false negative rate of 0% (0/17). Using an $SI \geq 1.3$ also
269 correctly classified all of the ICCVAM performance standards reference substances.

270 Intralaboratory reproducibility was assessed using a concordance analysis of sensitizer/non-
271 sensitizer results, and a coefficient of variation (CV) analysis of SI values and EC3 values

272 (calculated concentration corresponding to SI = 3). Four of six substances yielded 100%
273 concordance for sensitizer/non-sensitizer outcomes. Discordant test results were noted for two
274 substances: a sensitizer (isoeugenol) and a non-sensitizer (propylene glycol). Isoeugenol was
275 correctly identified as a sensitizer in 75% (3) of the four tests. Propylene glycol was correctly
276 identified as a non-sensitizer in 50% (1) of the two tests. The CVs for the SI values of five
277 substance/concentration combinations that were tested two times each ranged from 0.6% to
278 51.3%. The CVs for the EC3 values of four substances that were tested two to three times each
279 ranged from 10.1% to 47.1%.

280 Interlaboratory reproducibility could not be assessed because all LLNA: BrdU-ELISA results
281 were produced in one laboratory.

282 The LLNA: BrdU-ELISA will use the same number of animals when compared to the traditional
283 LLNA. However, since use of the traditional LLNA is restricted in some institutions because it
284 involves radioactivity, availability and use of the non-radioactive LLNA: BrdU-ELISA and the
285 test methods may lead to further reduction in use of the GP tests, which would provide for
286 reduced animal use and increased refinement due to the avoidance of pain and distress in the
287 LLNA procedure.

288 The transferability of the LLNA: BrdU-ELISA is expected to be similar to the traditional LLNA.
289 Compared to the traditional LLNA, the LLNA: BrdU-ELISA will not require facilities,
290 equipment, and licensing permits for handling radioactive materials. The level of training and
291 expertise needed to conduct the LLNA: BrdU-ELISA should be similar to the traditional LLNA
292 except that the understanding and use of ELISA is required.

293 ICCVAM has developed draft recommendations for the LLNA: BrdU-ELISA with regard to its
294 usefulness and limitations, test method protocol, and future studies to further characterize its
295 usefulness and limitations. These are provided in a separate document, *Draft ICCVAM Test*
296 *Method Recommendations, Non-radioactive Murine Local Lymph Node Assay: BrdU-ELISA Test*
297 *Method Protocol (LLNA: BrdU-ELISA)*.

298

299 1.0 Introduction

300 1.1 Historical Background

301 In 1999, the Interagency Coordinating Committee on the Validation of Alternative Methods
302 (ICCVAM) recommended that the murine local lymph node assay (LLNA) is a valid substitute
303 for currently accepted guinea pig (GP) test methods to assess the allergic contact dermatitis
304 (ACD) potential of many, but not all, types of substances. The recommendation was based on a
305 comprehensive evaluation that included an independent scientific peer review panel (Panel)
306 assessment of the validation status of the LLNA. The Panel report and the ICCVAM
307 recommendations (ICCVAM 1999) are available at the National Toxicology Program (NTP)
308 Interagency Center for the Evaluation of Alternative Toxicological Methods
309 (NICEATM)/ICCVAM website
310 (http://iccvam.niehs.nih.gov/docs/immunotox_docs/llna/llnarep.pdf).

311 ICCVAM forwarded recommendations to U.S. Federal agencies that the LLNA should be
312 considered for regulatory acceptance or other non-regulatory applications for assessing the ACD
313 potential of substances, while recognizing that some testing situations would still require the use
314 of traditional GP test methods (ICCVAM 1999, Sailstad et al. 2001). The LLNA was
315 subsequently incorporated into national and international test guidelines for the assessment of
316 skin sensitization (Organisation for Economic Co-operation and Development [OECD] Test
317 Guideline 429 [OECD 2002]; International Standards Organization [ISO] 10993-10: Tests for
318 Irritation and Sensitization [ISO 2002]; U.S. Environmental Protection Agency [EPA] Health
319 Effect Testing Guidelines on Skin Sensitization [EPA 2003]).

320 On January 10, 2007, the U.S. Consumer Product Safety Commission (CPSC) formally
321 nominated several activities related to the LLNA for evaluation by ICCVAM and NICEATM
322 (Available at http://iccvam.niehs.nih.gov/methods/immunotox/llndocs/CPSC_LLNA_nom.pdf).

323 One of the nominated activities was an assessment of the validation status of non-radioactive
324 alternatives to the current version of the LLNA ([ICCVAM 1999, Dean et al. 2001] referred to
325 hereafter as the “traditional LLNA”), which uses radioactivity to detect sensitizers. The
326 information described in this background review document (BRD) was compiled by ICCVAM
327 and NICEATM in response to this nomination. The BRD provides a comprehensive review of
328 available data and information regarding the usefulness and limitations of one of these methods,

329 the LLNA with detection of bromodeoxyuridine (BrdU) incorporation by enzyme-linked
330 immunosorbent assay (ELISA) (referred to hereafter as the “LLNA: BrdU-ELISA”).

331 **1.2 The LLNA: BrdU-ELISA**

332 The LLNA: BrdU-ELISA was developed by Takeyoshi et al. (2001) as a non-radioactive
333 alternative to the traditional LLNA. While the traditional LLNA assesses cellular proliferation by
334 measuring the incorporation of radioactivity into the deoxyribonucleic acid (DNA) of dividing
335 lymph node cells, the LLNA: BrdU-ELISA assesses the same endpoint by measuring the
336 incorporation of the thymidine analog BrdU, which is detected and quantified with an ELISA,
337 which is available as a kit commercially from several sources.

338 This document provides:

- 339 • A comprehensive summary of the LLNA: BrdU-ELISA test method protocol
- 340 • The substances used in the validation of the test method and the test results
- 341 • The performance characteristics (accuracy and reliability) of the test method
- 342 • Animal welfare considerations
- 343 • Other considerations relevant to the usefulness and limitations of this test method
344 (e.g., transferability, cost of the test method).

345 **2.0 LLNA: BrdU-ELISA Test Method Protocol**

346 The LLNA: BrdU-ELISA protocol (see **Appendix A**) is similar to the ICCVAM-recommended
347 protocol for the traditional LLNA (see
348 <http://iccvam.niehs.nih.gov/methods/immunotox/llnadocs/LLNAProt.pdf>), except for the method
349 used to assess lymphocyte proliferation. In both the LLNA: BrdU-ELISA and the traditional
350 LLNA, the test substance is administered on three consecutive days. In the traditional LLNA,
351 ³H- thymidine or ¹²⁵I-iododeoxyuridine (in phosphate buffered saline; 250 µL/mouse) is
352 administered via the tail vein two days after the final application of the test substance. In the
353 LLNA: BrdU-ELISA, 5 mg BrdU in a volume of 0.5 mL physiological saline (concentration of
354 10 mg/mL) is administered via intraperitoneal injection two days after the final application of the
355 test substance. Takeyoshi et al. (2001) reported that one injection of 5 mg BrdU was selected
356 over two injections to minimize the incorporation of BrdU in the control group. Injection of
357 BrdU two days after topical treatment with test substance yielded efficient incorporation of BrdU

358 in comparison to injection one day or three days after topical treatment with a test substance
359 (Takeyoshi et al. 2001). On the day following BrdU injection, lymph nodes are excised and a
360 single cell suspension is prepared from the lymph nodes of each animal. A standard aliquot of
361 the cell suspension is added in triplicate to the wells of a flat-bottom 96-well microplate and
362 centrifuged. Supernatants are then removed. Fix-Denat solution, which fixes the cells and
363 denatures the DNA, is added to each well, and the plate is incubated at room temperature. The
364 Fix-Denat solution is removed and the diluted anti-BrdU antibody solution is added to each well.
365 After each well is washed with phosphate-buffered saline, an aliquot of substrate solution
366 containing tetramethylbenzidine is added. After incubation at room temperature, the absorbance
367 is measured using a microplate reader.

368 2.1. Decision Criteria

369 Like the traditional LLNA, a stimulation index (SI) is used in the LLNA: BrdU-ELISA to
370 distinguish skin sensitizers from non-sensitizers. The SI is the ratio of the mean absorbance of
371 the incorporated BrdU in a lymph node suspension from individual mice in the test substance
372 group to the mean absorbance of the incorporated BrdU in a lymph node suspension from
373 individual mice in the vehicle control group as indicated by the formula below:

$$374 \quad SI = \frac{\text{Mean absorbance of the treatment group lymph nodes}}{\text{Mean absorbance of the vehicle control group lymph nodes}}$$

375 Consistent with the traditional LLNA, an $SI \geq 3$ was initially used as the threshold for labeling a
376 substance as a sensitizer. However, Takeyoshi et al. (2007b) also evaluated the use of
377 statistically significant differences in BrdU incorporation between treated and control groups,
378 and/or other SI values to distinguish sensitizers from non-sensitizers, and found that improved
379 accuracy resulted from using lower cutoff values for the SI as the decision criteria for whether a
380 substance was a sensitizer or non-sensitizer (see **Appendix B**).

381 **3.0 LLNA: BrdU-ELISA Validation Database**

382 To evaluate the validity of the LLNA: BrdU-ELISA, data were available for a total of 29
383 substances that had been tested in one laboratory (**Table 3-1**). Most of these substances
384 (24/29) had been previously tested in the traditional LLNA. No traditional LLNA data were
385 available for the remaining five substances, which are trans-cinnamaldehyde, two dimers of
386 eugenol (dihydroxyl-3,3'-dimethoxy-5,5'-diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-2,3'-

387 dimethoxyphenyl ether) and two dimers of isoeugenol (4-[1-Hydroxy-2-(2-methoxy-4-propenyl-
388 phenoxy)-propyl]-2-methoxy-phenol and 2-methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-
389 dihydro-benzofuran-2yl)-phenol).

390 Twenty-one of the 24 substances previously tested in the traditional LLNA were considered in
391 the original evaluation of the LLNA by ICCVAM (ICCVAM 1999). The traditional LLNA data
392 for the three remaining substances, cyclamen aldehyde, glutaraldehyde, and isopropyl myristate,
393 and for aniline were obtained from Basketter et al. (2005), Hilton et al. (1998), Ryan et al.
394 (2000), and Gerberick et al. (2005), respectively.

395 Of the 24 substances with traditional LLNA data, 16 were classified by the traditional LLNA as
396 skin sensitizers and eight were classified as non-sensitizers. As shown in **Table 3-1**, the
397 traditional LLNA EC3 values (i.e., estimated concentration needed to produce an SI =3) for the
398 16 sensitizers ranged from 0.005% to 44%.

399 **Appendix C** provides information on the physicochemical properties (e.g., physical form tested),
400 Chemical Abstracts Service Registry Number (CASRN), and chemical class for each substance
401 tested. When available, chemical classes for each substance were retrieved from the National
402 Library of Medicine's ChemID Plus database. If chemical classes were unavailable, they were
403 assigned to each test substance using a standard classification scheme based on the National
404 Library of Medicine Medical Subject Headings (MeSH) classification system (available at
405 <http://www.nlm.nih.gov/mesh/meshhome.html>). A substance could be assigned to more than one
406 chemical class; however, no substance was assigned to more than three classes. Chemical class
407 information is presented only to provide an indication of the variety of structural elements that
408 are present in the structures that were evaluated in this analysis. Classification of substances into
409 chemical classes is not intended to indicate the impact of structure on biological activity with
410 respect to sensitization potential. **Table 3-1** shows that 13 chemical classes are represented by
411 the substances tested in the LLNA: BrdU-ELISA. Three substances are classified in more than
412 one chemical class. The classes with the highest number of substances are carboxylic acids (9
413 substances) and aldehydes (5 substances).

414 **Table 3-1 Traditional LLNA EC3 Values and Chemical Classification of Substances**
 415 **Tested in the LLNA: BrdU-ELISA**

Substance Name	Chemical Class ¹	Traditional LLNA EC3 (%) ²
Diphenylcyclopropenone	Hydrocarbons, Cyclic	0.005
p-Benzoquinone	Quinones	0.01
2,4-Dinitrochlorobenzene*	Hydrocarbon, Halogenated; Nitro Compounds; Hydrocarbons, Cyclic	0.049
4-Phenylenediamine*	Amines	0.11
Glutaraldehyde	Aldehydes	0.14
Isoeugenol*	Carboxylic Acids	1.5
Cinnamic aldehyde	Aldehydes	2.4
3-Aminophenol	Amines; Phenols	3.2
4-Chloroaniline	Amines	6.5
2-Mercaptobenzothiazole*	Heterocyclic Compounds	9.8
Citral*	Hydrocarbons, Other	9.8
Hexyl cinnamic aldehyde*	Aldehydes	9.9
Eugenol*	Carboxylic Acids	10
Cyclamen aldehyde	Aldehydes	22.3
Hydroxycitronellal	Hydrocarbons, Other	23.8
Isopropyl myristate	Lipids	44
Aniline	Amines	63
Diethylphthalate	Carboxylic Acids	NA
Dimethyl isophthalate	Carboxylic Acids	NA
Glycerol	Alcohols; Carbohydrates	NA
Hexane	Hydrocarbons, Acyclic	NA
2-Hydroxypropyl methacrylate	Carboxylic Acids	NA
Isopropanol*	Alcohols	NA
Propylene glycol	Alcohols	NA
trans-Cinnamaldehyde	Aldehydes	NK
2,2'-Dihydroxyl-3,3'-dimethoxy-5,5'-diallyl-biphenyl	Carboxylic Acids	NK
2-Methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-dihydro-benzofuran-2yl)-phenol	Carboxylic Acids	NK
4,5'-Diallyl-2'-hydroxy-2,3'-dimethoxyphenyl ether	Carboxylic Acids	NK
4-[1-Hydroxy-2-(2-methoxy-4-propenyl-phenoxy)-propyl]-2-methoxy-phenol (Synonym: □-O-4-Dilignol)	Carboxylic Acids	NK

416 Abbreviations: LLNA: BrdU-ELISA= Local lymph node assay with enzyme-linked immunosorbent assay detection
 417 of bromodeoxyuridine; EC3 = Estimated concentration needed to produce a stimulation index (SI) = 3; NA = Not
 418 applicable since maximum SI < 3.0; NK = Not known (information requested but not yet obtained).

419 *Reference substance from ICCVAM (2007).

420 ¹Chemical classifications based on the Medical Subject Headings classification for chemicals and drugs, developed
421 by the National Library of Medicine (<http://www.nlm.nih.gov/mesh/meshhome.html>).

422 ²Mean EC3 values from the NICEATM database of traditional LLNA studies.

423

424 **4.0 Reference Data**

425 The traditional LLNA reference data used for the accuracy evaluation described in **Section 6.0**
426 were obtained from ICCVAM (1999), Basketter et al. (2005), Hilton et al. (1998), Ryan et al.
427 (2000), or Gerberick et al. (2005) (**Appendix D**). An independent quality assurance contractor
428 for the NTP audited the traditional LLNA data provided in ICCVAM (1999). Audit procedures
429 and findings are presented in the quality assurance report on file at the National Institute of
430 Environmental Health Sciences (NIEHS). The audit supports the conclusion that the transcribed
431 test data in the submission were accurate, consistent, and complete as compared to the original
432 study records. No traditional LLNA data could be located for trans-cinnamaldehyde.

433 The reference data for the GP tests (Guinea Pig Maximization Test [GPMT] or Buehler Test
434 [BT]) and human tests (Human Maximization Test [HMT], Human Patch Test Allergen [HPTA],
435 or other human data) were obtained from ICCVAM (1999), Basketter et al. (2000), Bjorkner
436 (1984), Hilton et al. (1998), Marzulli et al. (1974), Opdyke (1976), Takeyoshi et al. (2004a), or
437 Takeyoshi et al. (2007a). Although there were no traditional LLNA data available for the
438 eugenol dimers (dihydroxyl-3,3'-dimethoxy-5,5'-diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-
439 2,3'-dimethoxyphenyl ether) or the isoeugenol dimers (4-[1-Hydroxy-2-(2-methoxy-4-propenyl-
440 phenoxy)-propyl]-2-methoxy-phenol and 2-Methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-
441 dihydro-benzofuran-2yl)-phenol), Takeyoshi et al. (2004a and 2007a, respectively) provided
442 results from the GPMT for these compounds.

443 **5.0 Test Method Data and Results**

444 The LLNA: BrdU-ELISA data evaluated in this technical summary were obtained from six
445 published studies (Takeyoshi et al. 2003; 2004a; 2004b; 2005; 2006; 2007a) and one platform
446 presentation (Takeyoshi 2007b). All test results were obtained using the protocol in **Appendix**
447 **A. Appendix D** contains the LLNA: BrdU-ELISA data for the 29 substances tested in these
448 studies. The test substances were not coded to prevent the possibility of bias.

449

449 6.0 Test Method Accuracy

450 A critical component of a formal evaluation of the validation status of a test method is an
451 assessment of the accuracy of the proposed tested method when compared to the current
452 reference test method (ICCVAM 2003). Additional comparisons should also be made against
453 available human data, including experience from testing or accidental exposures. This aspect of
454 assay performance is typically evaluated by calculating:

- 455 • Accuracy (concordance): the proportion of correct outcomes (positive and
456 negative) of a test method
- 457 • Sensitivity: the proportion of all positive substances that are classified as positive
- 458 • Specificity: the proportion of all negative substances that are classified as
459 negative
- 460 • False positive rate: the proportion of all negative substances that are incorrectly
461 identified as positive
- 462 • False negative rate: the proportion of all positive substances that are incorrectly
463 identified as negative.

464 6.1 Total LLNA: BrdU-ELISA Database Analysis Using $SI \geq 3$ Decision Criteria

465 The performance characteristics of the LLNA: BrdU-ELISA were first evaluated using the
466 criterion of $SI \geq 3$ to identify sensitizers. Twenty-four of the 29 substances listed in **Table 3-1**
467 had sufficient LLNA: BrdU-ELISA and traditional LLNA data to conduct an accuracy analysis.
468 Of the remaining substances tested with the LLNA: BrdU-ELISA, 17 had LLNA: BrdU-ELISA,
469 traditional LLNA, and GP data; and 21 had LLNA: BrdU-ELISA, traditional LLNA, and human
470 data. 3-Aminophenol was excluded from the accuracy analyses for the dataset with LLNA:
471 BrdU-ELISA, traditional LLNA, and GP data since the available GP data were generated with a
472 nonstandard GPMT protocol. The nonstandard protocol did not include the 48-hour topical patch
473 induction that should follow induction by intradermal injections and it replaced the 24-hour skin
474 patch challenge (usually two weeks after topical induction) with a 6-hour skin patch challenge
475 (Basketter D, personal communication). Trans-cinnamaldehyde, the eugenol dimers (dihydroxyl-
476 3,3'-dimethoxy-5,5'-diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-2,3'-dimethoxyphenyl ether),
477 and the isoeugenol dimers (4-[1-Hydroxy-2-(2-methoxy-4-propenyl-phenyloxy)-propyl]-2-

478 methoxy-phenol and 2-methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-dihydro-benzofuran-
479 2yl)-phenol) were excluded from the accuracy analyses because traditional LLNA data for these
480 substances were not identified. The complete set of data for each substance is located in
481 **Appendix D.**

482 Discordant test results were noted for two substances. Four LLNA: BrdU-ELISA test results
483 were available for isoeugenol; three tests were positive for skin sensitization potential
484 (Takeyoshi et al. 2005; 2007a) and one was negative (Takeyoshi et al. 2006). Based on a weight
485 of evidence, a positive result was used for analysis of performance characteristics. Two
486 discordant test results were noted for propylene glycol. The most conservative result with respect
487 to prediction of adverse health effects, positive for skin sensitization, was used for the accuracy
488 analyses. The test result in Takeyoshi et al. (2005) produced a positive result as indicated by
489 individual animal data submitted by Dr. Takeyoshi to support the graphical data shown in the
490 paper. The test result in Takeyoshi et al. (2006) was negative.

491 6.1.1 *Accuracy vs. the Traditional LLNA*

492 When compared to the traditional LLNA and using a decision criteria of $SI \geq 3.0$ to identify
493 sensitizers, the LLNA: BrdU-ELISA had an accuracy of 75% (18/24), a sensitivity of 71%
494 (12/17), a specificity of 86% (6/7), a false positive rate of 14% (1/7), and a false negative rate of
495 29% (5/17) (**Table 6-1**).

496

497 **Table 6-1 Evaluation of the Performance of the LLNA: BrdU-ELISA In Predicting Skin Sensitizing Potential Using**
 498 **Decision Criteria of SI \geq 3.0 to Identify Sensitizers**

Comparison	n ¹	Accuracy		Sensitivity		Specificity		Positive Predictivity		Negative Predictivity		False Positive Rate		False Negative Rate	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
BrdU-ELISA vs. Traditional LLNA	24	75	18/24	71	12/17	86	6/7	92	12/13	55	6/11	14	1/7	29	5/17
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and GP Data</i>															
BrdU-ELISA vs. Traditional LLNA	17	71	12/17	67	8/12	80	4/5	89	8/9	50	4/8	20	1/5	33	4/12
LLNA: BrdU-ELISA vs. GP³	17	71	12/17	67	8/12	80	4/5	89	8/9	50	4/8	20	1/5	33	4/12
Traditional LLNA vs. GP³	17	100	17/17	100	12/12	100	5/5	100	12/12	100	5/5	0	0/5	0	0/12
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and Human Data</i>															
BrdU-ELISA vs. Traditional LLNA	21	76	16/21	73	11/15	83	5/6	92	11/12	56	5/9	17	1/6	27	4/15
LLNA: BrdU-ELISA vs. Human⁴	21	62	13/21	61	11/18	67	2/3	92	11/12	22	2/9	33	1/3	39	7/18
Traditional LLNA vs. Human⁴	21	76	16/21	78	14/18	67	2/3	93	13/14	33	2/6	33	1/3	22	4/18

499 Abbreviations: LLNA: BrdU-ELISA= Murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA) detection of
 500 bromodeoxyuridine (BrdU); GP = Guinea pig skin sensitization outcomes; ICCVAM = Interagency Coordinating Committee on the Validation of
 501 Alternative Methods; LLNA = Murine local lymph node assay; No. = Number.

502 ¹n = Number of substances included in this analysis.

503 ²The data on which the percentage calculation is based.

504 ³GP refers to outcomes obtained by studies conducted using either the Guinea Pig Maximization Test or the Buehler Test.

505 ⁴Human refers to outcomes obtained by studies conducting using the Human Maximization Test, inclusion of the test substance in a Human Patch Test
 506 Allergen Kit, and/or published clinical case studies/reports.

507 6.1.2 *Accuracy vs. Guinea Pig Data*

508 When the accuracy of the LLNA: BrdU-ELISA ($SI \geq 3.0$) and the traditional LLNA were
509 compared based on their performance relative to the GP test, the LLNA: BrdU-ELISA had a
510 lower accuracy rate (71% [12/17] vs. 100% [17/17]) and sensitivity (67% [8/12] vs. 100%
511 [12/12]), and higher false negative rate (33% [4/12] vs. 0% [0/12]) (**Table 6-1**). The specificity
512 was lower (80% [4/5] vs. 100% [5/5]) and the false positive rate was higher (20%, [1/5] vs. 0%
513 [0/5]) for the LLNA: BrdU-ELISA than that for the traditional LLNA.

514 6.1.3 *Accuracy vs. Human Data*

515 When the accuracy of the LLNA: BrdU-ELISA ($SI \geq 3.0$) and the traditional LLNA were
516 compared based on their performance relative to the available human data, the LLNA: BrdU-
517 ELISA had a lower accuracy (62% [13/21] vs. 76% [16/21]) and sensitivity (61% [11/18] vs.
518 78% [14/18]) and a higher false negative rate (39% [7/18] vs. 22% [4/18]) than the traditional
519 LLNA (**Table 6-1**). The specificity (67% [2/3]) and the false positive rate (33%, [1/3]) for the
520 LLNA: BrdU-ELISA was not different from that of the traditional LLNA.

521 6.2 Accuracy Analysis ($SI \geq 3.0$) Based on Revised Draft ICCVAM Performance
522 Standards Reference Substances

523 As shown in **Table 6-2**, eight of the 18 minimum reference substances included in the *Revised*
524 *Draft ICCVAM Performance Standards for the Local Lymph Node Assay* (ICCVAM 2007) have
525 been tested in the LLNA: BrdU-ELISA. Seven of the eight substances yielded the same
526 sensitizer/non-sensitizer outcome in the LLNA: BrdU-ELISA as in the traditional LLNA. 2-
527 Mercaptobenzothiazole, a sensitizer in the LLNA, tested negative in the LLNA: BrdU-ELISA.
528 One could suspect that testing in different vehicles produced the difference between the
529 ICCVAM (2007) EC3 values and those for the LLNA: BrdU-ELISA. The test results used by
530 ICCVAM (2007) used acetone:olive oil (4:1; AOO) for the vehicle while those for the LLNA:
531 BrdU-ELISA used dimethylformamide (DMF). However, the NICEATM database of traditional
532 LLNA studies indicates that 2-mercaptobenzothiazole has a higher EC3 when tested in AOO
533 (mean EC3 =9.8%) compared with DMF (mean EC3 =2.5%). Thus, the use of DMF as the
534 vehicle should have made 2-mercaptobenzothiazole more likely to be positive in the LLNA:
535 BrdU-ELISA.

536 **Table 6-2 Evaluation of the Performance of the LLNA: BrdU-ELISA (SI ≥ 3.0) Using the Revised Draft ICCVAM**
 537 **Performance Standards Reference Substances¹**

Substance	ICCVAM Draft LLNA Performance Standards ¹					LLNA: BrdU-ELISA ²		
	Vehicle	Result	EC3 (%) ¹	N ²	0.5x – 2.0x EC3 (%)	Vehicle	Result	EC3 (%)
5-Chloro-2-methyl-4-isothiazolin-3-one	DMF	+	0.009	1	0.0045 – 0.018	NA	NT	NT
2, 4-Dinitrochlorobenzene	AOO	+	0.049	15	0.025 – 0.099	AOO	+	0.15 (n=2)³
4-Phenylenediamine	AOO	+	0.11	10	0.055 – 0.22	AOO	+	NR (+)
4-Methylaminophenol sulfate	DMF	+	0.8	1	0.4 – 0.12	NA	NT	NT
Isoeugenol	AOO	+	1.5	49	0.77 – 3.1	AOO	+	9.3 (n=3)³
2-Mercaptobenzothiazole	AOO	+	2.5	2	1.25 – 5.0	DMF	-	NA (-)
Cobalt chloride	DMSO	+	4.8	1	2.4 – 9.6	NA	NT	NT
Citral	AOO	+	9.8	2	4.9 – 19.6	AOO	+	NR (+)
Hexyl cinnamic aldehyde	AOO	+	9.9	22	5.0 – 19.9	AOO	+	41.2 (n=3)³
Eugenol	AOO	+	10.1	11	5.05 – 20.2	AOO	+	29.5 (n=3)³
Phenyl benzoate	AOO	+	13.6	3	6.8 – 27.2	NA	NT	NT
Cinnamic alcohol	AOO	+	21	1	10.5 - 42	NA	NT	NT
Imidazolidinyl urea	DMF	+	24	1	12 - 36	NA	NT	NT
Chlorobenzene	AOO	-	NA	1	NA	NA	NT	NT
Isopropanol	AOO	-	NA	1	NA	AOO	-	NA (-)
Lactic acid	DMSO	-	NA	2	NA	NA	NT	NT
Methyl salicylate	AOO	-	NA	10	NA	NA	NT	NT
Salicylic acid	AOO	-	NA	1	NA	NA	NT	NT
Ethylene glycol dimethacrylate	MEK	False +	28 (FP)	1	14-56	NA	NT	NT
Sodium lauryl sulfate	DMF	False +	8.1 (FP)	5	4.05 – 16.2	NA	NT	NT
Nickel sulfate	DMF	False -	NA (FN)	2	NA	NA	NT	NT
Sulfanilamide	DMF	False -	NA (FN)	1	NA	NA	NT	NT

538 Bolded italics text highlights discordant LLNA:BrdU-ELISA vs. traditional LLNA test results.

539 Abbreviations: AOO = acetone:olive oil (4:1); LLNA: BrdU-ELISA= Murine local lymph node assay with enzyme-linked immunosorbent assay detection of
 540 bromodeoxyuridine; DMF = *N,N*-dimethylformamide; DMSO = Dimethyl sulfoxide; EC3 = Calculated concentration that corresponds to SI=3; FN = False
 541 negative in traditional LLNA when compared to guinea pig and/or human results; FP = False positive in traditional LLNA when compared to guinea pig and/or
 542 human results; LLNA = Murine local lymph node assay; MEK = Methyl ethyl ketone; NA = Not applicable (Stimulation Index < 3); NR = Not reported
 543 (information requested by NICEATM); NT = Not tested.

544 + = Sensitizer.

545 - = Non-sensitizer.

546 ¹Mean EC3 when more than one value was available. From Revised Draft ICCVAM Performance Standards for the LLNA (ICCVAM 2007; available:
547 http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm.

548 ²From Takeyoshi et al. (2003, 2004b, 2005, 2006, 2007a, 2007b). Substances for which EC3 values were not available include in parentheses the outcome of the
549 BrdU-ELISA test (+ = sensitizer; - = nonsensitizer>

550 ³Number of values used to derive the mean EC3.

551

552

553 The *Revised Draft ICCVAM Performance Standards for the Local Lymph Node Assay* (ICCVAM
554 2007) recommend a range of 0.5X to 2.0X the historical mean EC3 traditional LLNA as the
555 criteria for accuracy for nontraditional LLNA methods. EC3 values from the LLNA: BrdU-
556 ELISA were reported for four of the seven ICCVAM reference sensitizers tested. The EC3
557 values for all four substances were outside of the proposed acceptability range.

558 Selected characteristics of the substances tested using the LLNA: BrdU-ELISA were compared
559 with the characteristics of the 18 minimum reference substances included in the *Revised Draft*
560 *ICCVAM Performance Standards for the Local Lymph Node Assay* (ICCVAM 2007). **Table 6-3**
561 provides the range of substances tested in the LLNA: BrdU-ELISA based on the overall database
562 of 24 substances in comparison to the range of substances included on the revised draft
563 ICCVAM LLNA performance standards substances list. The table indicates that although not all
564 of the draft ICCVAM performance standards reference substances have been tested, the range of
565 the substances tested in the LLNA: BrdU-ELISA based on traditional LLNA EC3 values,
566 physical form of the test substance, and availability of reference data is similar to that included in
567 the draft performance standards list. In general, there are a proportionally increased number of
568 substances tested in the LLNA: BrdU-ELISA in each of the categories included in the table.

569

569 **Table 6-3 Characteristics of the Substances Tested in the LLNA: BrdU-ELISA vs. the**
 570 **Revised Draft ICCVAM Performance Standards Substances List¹**

EC3 range (%)	No. Chems	Solid/Liquid	Actual EC3 Range (%) ²	Human Data	Peptide Reactivity (High/Mod/Min/Unk) ³
<0.1	3	3/0	0.005 - 0.049	2	3/0/0/0
	2	1/1	0.009-0.05	2	0/1/0/1
≥ 0.1 to <1	2	1/1	0.11 - 0.14	2	1/0/0/1
	2	2/0	0.11-0.8	2	1/0/0/1
≥ 1 to <10	7	2/5	1.5 - 9.9	6	1/0/1/5
	5	2/3	1.6-9.9	5	1/0/1/3
≥ 10 to <100	5	0/5	10 - 63	4	0/0/3/2
	4	3/1	10.1-24	4	0/1/0/3
Negative	7	2/5	NC	4	0/0/6/1
	5	2/3	NC	3	0/0/2/3
Overall	24	8/16	0.005 - 63	18	5/0/10/9
	18	10/8	0.009-24	16	2/2/3/11

571 Bolded text represents characteristics of the LLNA: BrdU-ELISA database.

572 Abbreviations: Chems = Chemicals; EC3 = Estimated concentration needed to produce a stimulation index of three;
 573 NC = Not calculated because maximum SI <3.0; No. = Number; Min = Minimal; Mod = Moderate; SI = Stimulation
 574 Index; Unk = Unknown.

575 ¹From Revised Draft ICCVAM Performance Standards for the LLNA (available:

576 http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm). Includes the 18 "required" substances for
 577 testing.

578 ²Based on traditional LLNA studies for substances in the LLNA: BrdU-ELISA database (bold values) and the draft
 579 ICCVAM LLNA performance standards substances.

580 ³Data obtained from: Gerberick et al. (2007)

581 6.3 Discordant Results for Accuracy Analysis of SI ≥ 3.0 Decision Criterion

582 When the outcomes for the 24 substances tested in the LLNA: BrdU-ELISA (using SI ≥ 3.0) and
 583 the traditional LLNA were compared, the classifications for six substances were different. The
 584 LLNA: BrdU-ELISA classified aniline, 4-chloroaniline, cyclamen aldehyde, hydroxycitronellal,
 585 and 2-mercaptobenzothiazole as non-sensitizers while the traditional LLNA classified them as
 586 sensitizers (**Table 6-4**). No commonalities in chemical class, physical form, or EC3 value (based
 587 on the traditional LLNA) were noted among these substances. The LLNA: BrdU-ELISA
 588 classified propylene glycol as a sensitizer while the traditional LLNA classified it as a non-
 589 sensitizer. Two discordant LLNA: BrdU-ELISA results for propylene glycol were obtained. The
 590 most conservative result was used for the accuracy analyses to be conservative with respect to

590 prediction of adverse health effects. Additionally, the substances were tested in the same vehicle
 591 in both the LLNA: BrdU-ELISA and the traditional LLNA tests. Aniline, 4-chloroaniline,
 592 cyclamen aldehyde, propylene glycol, and hydroxycitronellal were tested in AOO, while 2-
 593 mercaptobenzothiazole was tested in DMF.

594 **Table 6-4 Discordant Results for LLNA: BrdU-ELISA (Using SI \geq 3.0 for Sensitizers)**
 595 **Compared to Traditional LLNA and Guinea Pig Reference Data**

Substance Name	LLNA: BrdU-ELISA ¹	Traditional LLNA ²	Guinea Pig Studies ²
Aniline	-	+ ³	+
4-Chloroaniline	-	+	+
Cyclamen aldehyde	-	+ ⁵	NA
Hydroxycitronellal	-	+	+
2-Mercaptobenzothiazole	-	+	+
Propylene glycol	+ ⁴	-	-

596 Abbreviations: LLNA: BrdU-ELISA= Murine local lymph node assay with enzyme-linked
 597 immunosorbent assay detection of bromodeoxyuridine; GP = Outcomes of guinea pig skin sensitization
 598 tests; LLNA = Murine local lymph node assay; NA = Not available.

599 + = Sensitizer.

600 - = Non-sensitizer.

601 ¹From Takeyoshi et al. (2005, 2006, 2007b).

602 ²From ICCVAM (1999) unless otherwise noted.

603 ³From Gerberick et al. (2005)

604 ⁴The test result in Takeyoshi et al. (2005) produced a positive result as indicated by individual animal data
 605 submitted by Dr. Takeyoshi to support the graphical data shown in the paper. The test result in Takeyoshi
 606 et al. (2006) produced a negative result. Both tests used a maximum concentration of 50%. The overall
 607 result was deemed to be positive (i.e., a conservative approach was used).

608 ⁵From Basketter et al. (2005).

609

610 For the 17 substances with LLNA: BrdU-ELISA, traditional LLNA, and GP test results, The
 611 results for aniline, 4-chloroaniline, hydroxycitronellal, 2-mercaptobenzothiazole, and propylene
 612 glycol were discordant with the GP test results. The LLNA: BrdU-ELISA results for aniline, 4-
 613 chloroaniline, hydroxycitronellal, and 2-mercaptobenzothiazole were negative, while the
 614 traditional LLNA and GP results were positive. The LLNA: BrdU-ELISA result for propylene
 615 glycol was positive, while the traditional LLNA and GP results were negative.

616 When analyses were restricted to the 21 substances with LLNA: BrdU-ELISA, traditional
 617 LLNA, and human outcomes, both LLNA methods misclassified three sensitizers (2-

618 hydroxypropylmethacrylate, isopropanol, and diethyl phthalate) as non-sensitizers (Table 6-5).
 619 The LLNA: BrdU-ELISA also misclassified four more sensitizers as non-sensitizers that were
 620 correctly classified by the traditional LLNA: aniline, hydroxycitronellal, cyclamen aldehyde, and
 621 2-mercaptobenzothiazole. No commonalities in chemical class, physical form, or EC3 range
 622 (based on the traditional LLNA) were noted among these substances. Both the LLNA: BrdU-
 623 ELISA and the traditional LLNA misclassified isopropyl myristate as a sensitizer.

624 **Table 6-5 Discordant Results for LLNA: BrdU-ELISA (SI \geq 3.0) When Compared to**
 625 **Traditional LLNA and Human Outcome Data**

Substance Name	LLNA: BrdU-ELISA ¹	Traditional LLNA ²	Human Outcome ²
Aniline	-	+ ⁷	+
Cyclamen aldehyde	-	+ ⁸	+ ⁸
Hydroxycitronellal	-	+	+
2-Hydroxypropylmethacrylate	-	-	+ ³
Isopropanol	-	-	+ ⁴
Isopropyl myristate	+	+ ⁵	- ⁶
2-Mercaptobenzothiazole	-	+	+
Diethyl phthalate	-	-	+

626 Abbreviations: LLNA: BrdU-ELISA= Murine local lymph node assay with enzyme-linked
 627 immunosorbent assay detection of bromodeoxyuridine; GP = outcomes of guinea pig skin sensitization
 628 tests; LLNA = Murine local lymph node assay.

629 + = Sensitizer.

630 - = Nonsensitizer.

631 ¹From Takeyoshi et al. (2005, 2006, 2007b).

632 ²From ICCVAM (1999) unless otherwise noted.

633 ³From Bjorkner (1984).

634 ⁴From Kwon et al. (2003).

635 ⁵From Ryan et al. (2000).

636 ⁶From Opdyke (1976).

637 ⁷From Gerberick et al. (2005).

638 ⁸From Basketter et al. (2005).

639

640 The accuracy analyses for the eight reference substances from the Revised Draft ICCVAM
 641 Performance Standards (ICCVAM 2007) tested in LLNA: BrdU-ELISA yielded one discordant
 642 substance, 2-mercaptobenzothiazole. The LLNA: BrdU-ELISA classified this substance as a
 643 non-sensitizer, while the traditional LLNA, GP, and human tests classified it as a sensitizer.
 644 While the vehicles for the historical results reported in ICCVAM (2007) were AOO and that

645 used by Takeyoshi et al. (2007b) was DMF, the different vehicles were not responsible for the
646 discordant results. Other reports of traditional LLNA tests using DMF have also classified 2-
647 mercaptobenzothiazole as a sensitizer (e.g., Ashby et al. 1995; Gerberick et al. 2005).

648 6.4 LLNA: BrdU-ELISA Accuracy Analysis Using Alternative Decision Criteria

649 Takeyoshi et al (2007b) evaluated the effect of decision criteria other than $SI \geq 3$ to determine
650 skin sensitization potential on test performance characteristics with the traditional LLNA serving
651 as the reference test. The performance characteristics for nine alternate decision criteria for
652 determining whether the skin sensitization potential for the substances were positive or negative
653 are reported in this section. **Appendix B** also reports results for intermediate SI cutoff values.
654 The substances evaluated were the same as those evaluated in **Sections 6.1** through **6.3** except
655 that hexane was not included. The decision criteria included:

- 656 1. SI values ≥ 1.3 , ≥ 1.5 , ≥ 2 , ≥ 2.5 , or ≥ 3
- 657 2. Mean BrdU labeling index is statistically different from control group
- 658 3. Mean BrdU labeling index $\geq 95\%$ confidence interval of the control group
- 659 4. Mean BrdU labeling index is ≥ 2 standard deviations (SD) or ≥ 3 SD from the
660 control group mean

661 Using a decision criteria of $SI \geq 3.0$ to identify sensitizers for these 23 substances, the LLNA:
662 BrdU-ELISA had an accuracy of 74% (17/23), a sensitivity of 71% (12/17), a specificity of 83%
663 (5/6), a false positive rate of 17% (1/6), and a false negative rate of 29% (5/17) (**Table 6-6**).

664 However, when the decision criteria are altered to include lower SI values, improved
665 performance was achieved. When the mean labeling index for the treatment group was outside
666 the 95% confidence interval of the control group or ≥ 2 SD from the index for the vehicle control
667 group, or when $SI \geq 1.5$, the LLNA: BrdU-ELISA accuracy improved to 91% (21/23), with a
668 sensitivity of 94% (16/17), and a false negative rate of 6% (1/17). The specificity (83% [5/6])
669 and false positive rate (17% [1/6]) were the same as that for $SI \geq 3$.

670 The best overall performance was achieved using an $SI \geq 1.3$ with an accuracy of 96% (22/23),
671 sensitivity of 100% (17/17), specificity of 83% (5/6), a false positive rate of 17% (1/6), and false

672 negative rate of 0% (0/17). Using an $SI \geq 1.3$ also correctly classified all of the ICCVAM
673 performance standards reference substances.

674 6.5 Discordant Results for Accuracy Analysis of Alternative Decision Criterion

675 Using the decision criteria of $SI \geq 3.0$ to identify sensitizers for the 23 substances used in the
676 analysis of alternative decision criteria, the six discordant substances (when compared to the
677 traditional LLNA) were propylene glycol, 4-chloroaniline, hydroxycitronellal, aniline, cyclamen
678 aldehyde, and 2-mercaptobenzothiazole (**Table 6-7**). As indicated in **Section 6.3**, 4-
679 chloroaniline, aniline, hydroxycitronellal, cyclamen aldehyde, and 2-mercaptobenzothiazole
680 were misclassified as non-sensitizers, and propylene glycol was misclassified as a sensitizer
681 when compared to the traditional LLNA.

682 **Table 6-7** shows how the number and identity of discordant substances changes with the
683 alternate decision criteria. Using $SI \geq 2.0$ or $SI \geq 2.5$ yielded the same five discordant substances:
684 propylene glycol, hydroxycitronellal, aniline, cyclamen aldehyde, and 2-mercaptobenzothiazole.
685 Three discordant substances, propylene glycol, aniline, and hydroxycitronellal, were noted when
686 a statistical test was used to determine a difference between the treatment and vehicle control
687 group means and when the treatment group mean >3 SD from the control group mean. Using SI
688 ≥ 1.5 , $\geq 95\%$ CI, or ≥ 2 SD yielded two discordant substances, propylene glycol and
689 hydroxycitronellal. Using $SI \geq 1.3$ to classify substances as sensitizers yielded only one
690 discordant substance (propylene glycol). As noted in **Section 6.4**, using the decision criterion of
691 $SI \geq 1.3$ would correctly classify the ICCVAM reference substance 2-mercaptobenzothiazole as
692 a sensitizer, which was incorrectly classified as a non-sensitizer using $SI \geq 3.0$ as the criterion
693 (see **Section 6.2**).

694

695 **Table 6-6 Evaluation of the Performance of the LLNA: BrdU-ELISA In Predicting Skin Sensitizing Potential Using**
 696 **Alternative Decision Criteria to Identify Sensitizers**

Alternate Criterion	N ¹	Accuracy		Sensitivity		Specificity		Positive Predictivity		Negative Predictivity		False Positive Rate		False Negative Rate	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
Statistics ³	23	87	20/23	88	15/17	83	5/6	94	15/16	71	5/7	17	1/6	12	2/17
≥ 95% CI	23	91	21/23	94	16/17	83	5/6	94	16/17	83	5/6	17	1/6	6	1/17
≥ 2 SD	23	91	21/23	94	16/17	83	5/6	94	16/17	83	5/6	17	1/6	6	1/17
≥ 3 SD	23	87	20/23	88	15/17	83	5/6	94	15/16	71	5/7	17	1/6	12	2/17
SI ≥ 3.0	23	74	17/23	71	12/17	83	5/6	92	12/13	50	5/10	17	1/6	29	5/17
SI ≥ 2.5	23	78	18/23	77	13/17	83	5/6	93	13/14	56	5/9	17	1/6	24	4/17
SI ≥ 2.0	23	78	18/23	77	13/17	83	5/6	93	13/14	56	5/9	17	1/6	24	4/17
SI ≥ 1.5	23	91	21/23	94	16/17	83	5/6	94	16/17	83	5/6	17	1/6	6	1/17
SI ≥ 1.3	23	96	22/23	100	17/17	83	5/6	94	17/18	100	5/5	17	1/6	0	0/17

697 Abbreviations: LLNA: BrdU-ELISA = Murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA) detection of bromodeoxyuridine (BrdU); CI =
 698 Confidence interval; No. = Number; SD = Standard deviation; SI = Stimulation index

699 ¹ n = Number of substances included in this analysis.

700 ² The data on which the percentage calculation is based.

701 ³ Statistical test for difference of group means.

702

703 **Table 6-7 Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision Criteria Compared to the Traditional**
 704 **LLNA**

Discordant Substance	Alternate Decision Criterion								
	Statistics ¹	≥ 95% CI	≥ 2 SD	≥ 3 SD	SI ≥ 3.0	SI ≥ 2.5	SI ≥ 2.0	SI ≥ 1.5	SI ≥ 1.3
Propylene glycol	X	X	X	X	X	X	X	X	X
4-Chloroaniline					X				
Hydroxycitronellal	X	X	X	X	X	X	X	X	
Aniline	X			X	X	X	X		
Cyclamen aldehyde					X	X	X		
2-Mercaptobenzothiazole					X	X	X		

705 Abbreviations: LLNA: BrdU-ELISA = Murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA) detection of
 706 bromodeoxyuridine (BrdU); CI = Confidence interval; SD = Standard deviation; SI = Stimulation index; X = Discordant result obtained in the
 707 LLNA: BrdU-ELISA when compared to the traditional LLNA

708 ¹ Statistical test for difference of group means.

709
 710

711 7.0 Test Method Reliability

712 An assessment of test method reliability (intra- and inter-laboratory reproducibility) is an
713 essential element of any evaluation of the performance of an alternative test method (ICCVAM
714 2003). Intralaboratory reproducibility refers to the extent to which qualified personnel within the
715 same laboratory can replicate results using a specific test protocol at different times.
716 Interlaboratory reproducibility refers to the extent to which different laboratories can replicate
717 results using the same protocol and test substances, and indicates the extent to which a test
718 method can be transferred successfully among laboratories.

719 Since several substances were tested in the LLNA: BrdU-ELISA multiple times, data were
720 available for an evaluation of intralaboratory reproducibility. However, interlaboratory
721 reproducibility could not be assessed because the test results were generated in one laboratory.
722 The test results for the LLNA: BrdU-ELISA are amenable to intralaboratory reproducibility
723 analyses for three endpoints: sensitizer (positive) or non-sensitizer (negative) classification, SI
724 values, and EC3 values. Analyses of intralaboratory reproducibility were performed using a
725 concordance analysis for the qualitative results (sensitizer vs. non-sensitizer) (**Section 7.1**) and a
726 coefficient of variation (CV) analysis for the quantitative results (SI values and EC3 values)
727 (**Sections 7.2** and **7.3**).

728 7.1 Intralaboratory Reproducibility – Qualitative Results

729 The dataset available for an intralaboratory concordance analysis of the qualitative test results for
730 the LLNA: BrdU-ELISA included six substances that were tested multiple times and classified
731 as sensitizers or non-sensitizers. Eugenol and isoeugenol were each tested four times, hexyl
732 cinnamic aldehyde and isoeugenol were each tested three times, and 2,4-dinitrochlorobenzene
733 and propylene glycol were each tested twice (Takeyoshi et al. 2003, 2004a, 2005, 2006, 2007a).
734 All substances were sensitizers in the traditional LLNA except for propylene glycol. The
735 multiple test results for eugenol, hexyl cinnamic aldehyde, and 2,4-dinitrochlorobenzene were
736 100% concordant.

737 Discordant test results were noted for isoeugenol and propylene glycol. Three of the four LLNA:
738 BrdU-ELISA results for isoeugenol were positive for skin sensitization potential (Takeyoshi et
739 al. 2005; 2007a) and one was negative (Takeyoshi et al. 2006), which yields a 75%
740 intralaboratory concordance. Two of the positive results (Takeyoshi et al. 2006) and one negative

741 result (Takeyoshi et al. 2006) were obtained at a maximum concentration of 10% isoeugenol.
742 The remaining positive result was obtained using a maximum concentration of 30% isoeugenol
743 (Takeyoshi et al. 2007a). Two discordant test results were noted for propylene glycol. A positive
744 result was indicated by individual animal data submitted by Dr. Takeyoshi to support the
745 graphical data shown in Takeyoshi et al. (2005) (although the graphical display indicated a
746 negative result). The test result from Takeyoshi et al. (2006) produced a negative result. Both
747 tests used a maximum concentration of 50%.

748 The qualitative intralaboratory concordance analysis for the traditional LLNA (ICCVAM 1999)
749 was based on a dataset of six substances that included six results each for benzocaine and hexyl
750 cinnamic aldehyde, five results for eugenol, four results each for isoeugenol and methyl
751 salicylate, and three results for 2,4-dinitrochlorobenzene. All intralaboratory results for each
752 substance were 100% concordant with the exception of one of the six benzocaine (5/6 or 83%
753 concordance) results that was reported as +/- (i.e., equivocal). An equivocal result was described
754 as one in which SI increases with dose, but does not reach the criterion of three for classification
755 as a sensitizer.

756 Thus, the intralaboratory concordance of qualitative results for the LLNA: BrdU-ELISA was
757 lower than that of the traditional LLNA, but it was based on a smaller dataset.

758 7.2 Intralaboratory Reproducibility – SI

759 Three of the Takeyoshi et al. (2003, 2005, 2007a) studies reported numerical SI values (i.e.,
760 values were reported in the text or tables rather than plotted on graphs) that allowed for an
761 assessment of intralaboratory reproducibility. The SI values reported for five
762 substance/concentration combinations that were tested twice (in separate experiments) were used
763 to calculate a CV for the assessment of intralaboratory variability. Hexyl cinnamic aldehyde was
764 tested twice at three different concentrations. Eugenol and isoeugenol were each tested twice at
765 one concentration. As shown by **Table 7-1**, the CVs ranged from 0.6% (hexyl cinnamic
766 aldehyde) to 51.3% (isoeugenol). The intralaboratory reproducibility of the traditional LLNA
767 was not assessed by CV analysis of SI values (ICCVAM 1999).

768

769

769 **Table 7-1 Intralaboratory Reproducibility for the SI of Tested Substances in LLNA:**
 770 **BrdU-ELISA - Coefficient of Variation**

Substance	Concentration Tested (%)	SI	Mean	SD	CV (%)	Takeyoshi et al. Reference
Eugenol	30	3.30	3.57	0.37	10.5%	2004a
Eugenol	30	3.83				2007
Hexyl cinnamic aldehyde	12.5	1.87	1.73	0.21	11.9%	2003
Hexyl cinnamic aldehyde	12.5	1.58				2003
Hexyl cinnamic aldehyde	25	2.42	2.41	0.01	0.6%	2003
Hexyl cinnamic aldehyde	25	2.40				2003
Hexyl cinnamic aldehyde	50	3.63	3.62	0.02	0.6%	2003
Hexyl cinnamic aldehyde	50	3.60				2005
Isoeugenol	10	5.20	3.82	1.96	51.3%	2005
Isoeugenol	10	2.43				2007

771 Abbreviations: CV = Coefficient of variation; SD = Standard deviation, SI = Stimulation index
 772

773 7.3 Intralaboratory Reproducibility – EC3

774 CV values were also calculated for the EC3 values of substances that were tested multiple times.
 775 Five Takeyoshi et al. (2003, 2004a, 2005, 2006, 2007a) studies reported multiple EC3 values, or
 776 SI values that could be used to interpolate EC3 values, for multiple tests of the same substances.
 777 Multiple EC3 values were available for four substances. Two EC3 values were reported for 2-
 778 dinitrochlorobenzene and three EC3 values each were reported for isoeugenol, hexyl cinnamic
 779 aldehyde, and isoeugenol. As shown by **Table 7-2**, the CVs ranged from 10.1% (hexyl cinnamic
 780 aldehyde) to 47.1% (2, 4-dinitrochlorobenzene).

781

781 **Table 7-2 Intralaboratory Reproducibility for the EC3 of Tested Substances in LLNA:**
 782 **BrdU-ELISA - Coefficient of Variation**

Substance	EC3	Mean	SD	CV (%)	Takeyoshi et al. Reference
2, 4-Dinitrochlorobenzene	0.2	0.15	0.07	47.1%	2005
2, 4-Dinitrochlorobenzene	0.1				2006
Isoeugenol	5.6	9.3	3.6	38.3%	2005
Isoeugenol	9.6				2006
Isoeugenol	12.7				2007b
Hexyl cinnamic aldehyde	40.8	41.2	4.2	10.1%	2005
Hexyl cinnamic aldehyde	45.5				2006
Hexyl cinnamic aldehyde	37.2				2003
Eugenol	25.1	29.5	9.7	32.8%	2004a
Eugenol	40.6				2006
Eugenol	22.8				2007b

783 Abbreviations: CV = Coefficient of variation; EC3 = Estimated concentration needed to produce a
 784 stimulation index of three; SD = Standard deviation.
 785

786 The intralaboratory reproducibility of the traditional LLNA was assessed by CV analysis of EC3
 787 values using a larger dataset (ICCVAM 1999) than that available for the LLNA: BrdU-ELISA
 788 analysis. Two EC3 values were reported by each of five laboratories for
 789 2, 4-dinitrochlorobenzene, five EC3 values were reported by one laboratory for isoeugenol, six
 790 EC3 values were reported for hexyl cinnamic aldehyde by two laboratories, and five EC3 values
 791 were reported for eugenol by one laboratory (**Table 7-3**).

792 Most intralaboratory CV values for the EC3 values from LLNA: BrdU-ELISA tests were higher
 793 than those reported in ICCVAM (1999) for the traditional LLNA. At 47.1%, the intralaboratory
 794 EC3 CV values from the LLNA: BrdU-ELISA tests of 2, 4-dinitrochlorobenzene (**Table 7-2**)

795 were at the top of the range cited in ICCVAM (1999) (**Table 7-3**). The intralaboratory EC3 CV
 796 from the LLNA: BrdU-ELISA tests of isoeugenol was greater than that from ICCVAM (1999)
 797 (38.3% vs. 26.1%). The intralaboratory EC3 CV from the LLNA: BrdU-ELISA tests of hexyl
 798 cinnamic aldehyde was lower than that reported by ICCVAM (1999) (10.1% vs. 18.7 to 26.7%).
 799 The intralaboratory EC3 CV from the LLNA: BrdU-ELISA tests of eugenol were higher than that
 800 reported by ICCVAM (1999) (32.8% vs. 18.4%).

801 **Table 7-3 Intralaboratory Reproducibility for the EC3 of Tested Substances in the**
 802 **Traditional LLNA¹**

Substance	Number of Laboratories	Number of Tests per Laboratory	CV (%)
2, 4-Dinitrochlorobenzene	5	2	12.9 – 47.1
Isoeugenol	1	5	26.1
Hexyl cinnamic aldehyde	2	6	18.7-26.7
Eugenol	1	5	18.4

803 Abbreviations: CV = Coefficient of variation; EC3 = Estimated concentration needed to produce a
 804 stimulation index of three.

805 ¹From ICCVAM (1999).

806 **8.0 Data Quality**

807 All of data were generated at the Hita Laboratory of the Chemicals Evaluation and Research
 808 Institute, Japan. Although the Hita Laboratory is a GLP-conforming facility, the studies on the
 809 LLNA: BrdU-ELISA did not conform fully with GLP guidelines since they were not intended
 810 for regulatory purposes. However, all systems employed for these studies (i.e., test facilities,
 811 study staff, reagents, and the other study elements) were reportedly the same as those employed
 812 in the fully GLP-compliant studies conducted in the laboratory. Although multiple staff
 813 members checked the reported data for consistency with the raw data, no audit report is available
 814 (Takeyoshi M, personal communication). The raw data are also not available for audit.

815

815 **9.0. Other Scientific Reports and Reviews**

816 A multi-laboratory validation study of the LLNA: BrdU-ELISA is underway in Japan (Kojima
817 H, personal communication). Seven laboratories are testing 10 substances (different from those
818 evaluated by Dr. Takeyoshi) using a revised version of Dr. Takeyoshi's protocol. The final tests
819 were scheduled for completion by the end of December 2007. The validation study management
820 team is scheduled to meet on February 15, 2008, to discuss the results. More information about
821 the validation study, including the protocol, will be added as it is received. NICEATM has
822 requested the identities of the substances tested, the number of laboratories participating, and the
823 number of times each substance was tested in each laboratory.

824 A set of studies were conducted by Yamano et al. using a similar LLNA: BrdU-ELISA based
825 method (Yamano et al. 2003, 2004, 2005, 2006, 2007). The test method protocol (e.g.,
826 application of test substance to ear of mouse) was similar to what was described in the Takeyoshi
827 et al. studies discussed above. Compared to the method Takeyoshi et al., which administered 5
828 mg BrdU/mouse, the concentration of BrdU administered (via intraperitoneal injection) was 150
829 mg/kg/15 mL saline, which would be approximately 3 mg BrdU/mouse (based on a 20 g mouse).
830 The studies discussed the use of a BrdU-ELISA based method to assess the skin sensitization
831 potential of a variety of substances including metal salts of naphthenic acid, methylated phenols,
832 industrial biocides, and preservatives.

833 The outcomes of these studies were not included in this evaluation since comparative traditional
834 LLNA data were not available for the substances evaluated. Therefore, a comparison of the
835 accuracy of the LLNA: BrdU-ELISA versus the traditional LLNA, when outcomes were
836 compared to guinea pig or human results, could not be conducted.

837 **10.0 Animal Welfare Considerations**

838 The LLNA: BrdU-ELISA will require the use of the same number of animals when compared to
839 the traditional LLNA. However, since the traditional LLNA uses radioactivity and as such its use
840 is restricted in some institutions, broader use of the non-radioactive LLNA: BrdU-ELISA
841 protocol in place of the GP test could further reduce the number of guinea pigs that are still being
842 used to assess skin sensitization.

843

843 10.1 Rationale for the Need to Use Animals

844 The rationale for the use of animals in the LLNA: BrdU-ELISA is the same as the rationale for
845 the traditional LLNA; there are no valid and accepted non-animal ways to determine the ACD
846 potential of substances and products, except for situations where human studies could be
847 conducted ethically and where such studies would meet regulatory safety assessment
848 requirements. The most detailed information about the induction and regulation of
849 immunological responses are available for mice (ICCVAM 1999).

850 10.2 Basis for Determining the Number of Animals Used

851 The number of animals used for the experimental, vehicle, and positive control groups is based
852 on the number of animals specified in the ICCVAM recommended traditional LLNA protocol
853 (ICCVAM 2001).

854 10.3 Reduction considerations

855 A further reduction of 40% (15 vs. 25) could be achieved by using a limit dose version of the
856 LLNA: BrdU-ELISA in cases where dose response information is not needed for hazard
857 identification purposes. In such an approach, only the highest soluble dose of the test article that
858 does not produce skin irritation or systemic toxicity would be administered, and the two lower
859 dose groups would not be used. Additional reductions could be achieved by testing more
860 substances concurrently, so that the same vehicle and positive control group could be used for
861 multiple substances, thus further reducing the number of animals for each additional substance
862 by 10 animals, or 40% (15 vs. 25).

863 **11.0 Practical Considerations**

864 Several issues are taken into account when assessing the practicality of using an alternative to an
865 existing test method. In addition to performance evaluations, assessments of the laboratory
866 equipment and supplies needed to conduct the alternative test method, level of personnel
867 training, labor costs, and the time required to complete the test method relative to the existing
868 test method are necessary. The time, personnel cost, and effort required to conduct the proposed
869 test method(s) must be considered to be reasonable when compared to the existing test method it
870 is intended to replace.

871

871 11.1 Transferability of the LLNA: BrdU-ELISA

872 Test method transferability addresses the ability of a method to be accurately and reliably
873 performed by multiple laboratories (ICCVAM 2003), including those experienced in the
874 particular type of procedure as well as laboratories with less or no experience in the particular
875 procedure. It would be expected that the transferability of the LLNA: BrdU-ELISA would
876 similar to the traditional LLNA, since the protocols of the two methods (except for the detection
877 of lymphocyte proliferation) are similar.

878 11.2 Facilities and Major Fixed Equipment Required to Conduct the LLNA: BrdU-ELISA

879 Compared to the traditional LLNA, the LLNA: BrdU-ELISA will not require facilities,
880 equipment, and licensing permits for handling radioactive materials. The remaining facilities
881 (e.g., animal care facilities) are the same between the two methods.

882 11.3 LLNA: BrdU-ELISA Training Considerations

883 The level of training and expertise needed to conduct the LLNA: BrdU-ELISA should be similar
884 to the traditional LLNA. Additionally, individuals will need to understand and perform ELISAs.

885 **12.0 References**

- 886 Ashby J, Basketter DA, Patona D, Kimber I. 1995. Structure activity relationships in skin
887 sensitization using the murine local lymph node assay. *Toxicology* 103:177-194.
- 888 Basketter DA, Blaikie L, Dearman RJ, Kimber I, Ryan CA, Gerberick GF. 2000. Use of the local
889 lymph node assay for the estimation of relative contact allergenic potency. *Contact Dermatitis*
890 42:344-348.
- 891 Basketter DA, Clapp C, Jefferies D, Safford B, Ryan CA, Gerberick F, Dearman RJ, Kimber I.
892 2005. Predictive identification of human skin sensitization thresholds. *Contact Dermatitis*
893 53:260-267.
- 894 Basketter DA, Kimber I. 2006. Predictive tests for irritants and allergens and their use in
895 quantitative risk assessment. In “Contact Dermatitis”, 4th Edition. Eds, Frosch P J, Menné T,
896 Lepoittevin J-P, Springer Verlag, Heidelberg, pp 179 – 188.
- 897 Bjorkner B. 1984. Contact allergy to 2-hydroxypropyl methacrylate (2-HPMA) in an ultraviolet
898 curable ink. *Acta Derm Venereol* 64: 264-267.
- 899 Dearman DJ, Kimber I. 2001. Development of the Local Lymph Node Assay for risk assessment
900 of chemical and formulations. HSE contract research report 399. Syngenta Central Toxicology
901 Laboratory, UK.
- 902 EPA. 2003. Health Effects Test Guideline, OPPTS 870.2600. Skin Sensitization EPA 712–C–
903 03–197. Washington, DC: U.S. Environmental Protection Agency.
- 904 Gerberick GF, Ryan CA, Kern PS, Schlatter H, Dearman RJ, Kimber I, Patlewicz GY, Basketter
905 DA. 2005. Compilation of historical local lymph node data for evaluation of skin sensitization
906 alternative methods. *Dermatitis* 16:157-202.
- 907 Gerberick GF, Vassallo JD, Foertsch LM, Price BB, Chaney JG, Lepoittevin JP. 2007.
908 Quantification of chemical peptide reactivity for screening contact allergens: a classification tree
909 model approach. *Toxicol Sci* 97:417-27.
- 910 Hilton J, Dearman RJ, Harvey P, Evans P, Basketter DA, Kimber I. 1998. Estimation of relative
911 skin sensitizing potency using the local lymph node assay: A comparison of formaldehyde with
912 glutaraldehyde. *Am J Contact Dermat* 9:29-33.

- 913 ICCVAM. 1997. Validation and Regulatory Acceptance of Toxicological Test Methods: A
914 Report of the ad hoc Interagency Coordinating Committee on the Validation of Alternative
915 Methods. NIH Publication No.: 97-3981. Research Triangle Park:National Toxicology Program.
- 916 ICCVAM 1999. The murine local lymph node assay: A test method for assessing the allergic
917 contact dermatitis potential of chemical/compounds. NIH Publication No. 99-4494. Research
918 Triangle Park, NC: National Toxicology Program.
- 919 ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised, and
920 Alternative Test Methods. NIH Publication No: 03-4508. Research Triangle Park:National
921 Toxicology Program.
- 922 ICCVAM. 2007. Revised Draft ICCVAM Murine Local Lymph Node Assay Performance
923 Standards. Interagency Coordinating Committee on the Validation of Alternative Methods.
924 Research Triangle Park:National Toxicology Program. Available at
925 <http://iccvam.niehs.nih.gov/methods/immunotox/>.
- 926 ISO. 2002. Biological evaluation of medical devices -- Part 10: Tests for irritation and delayed-
927 type hypersensitivity. Available for purchase at: <http://www.iso.org/iso/home.htm>
- 928 Kimber I, Dearman RJ, Scholes EW, Basketter DA. 1994. The local lymph node assay:
929 developments and applications. *Toxicology* 93:13-31.
- 930 Kimber I, Basketter DA, Butler M, Gamer A, Garrigue J-L, Gerberick GF, Newsome C, Steiling
931 W, Vohr H-W. 2003. Classification of contact allergens according to potency: proposals. *Food*
932 *Chem Toxicol* 41:1799-1809.
- 933 Kwon JA, Lee MS, Kim MI, Park YM, Kim HO, Kim CW. 2003. Allergic contact dermatitis
934 from dodecyl diaminoethylglycine and isopropyl alcohol in a commercial disinfectant swab.
935 *Contact Dermatitis* 48:339-340.
- 936 Marzulli FN, Maibach HI. 1974. The use of graded concentrations in studying skin sensitizers:
937 experimental contact sensitization in man. *Food Cosmet Toxicol* 12:219-227.
- 938 OECD. 2002. Test guideline 429. Skin Sensitisation: Local Lymph Node Assay, adopted April
939 24, 2002. In: OECD Guidelines for Testing of Chemicals. Paris:OECD.

- 940 Opdyke DLJ. 1976. Fragrance raw materials monographs. Isopropyl myristate. Food Cosmet
941 Toxicol 14:323-325.
- 942 Ryan CA, Gerberick G F, Cruse LW, Basketter DA, Lea L, Blaikie L, Dearman RJ, Warbrick
943 EV, Kimber I. 2000. Activity of human contact allergens in the murine local lymph node assay
944 Contact Dermatitis 43:95-102.
- 945 Sailstad DM, Hattan D, Hill RN, Stokes WS. 2001. ICCVAM Evaluation of the Murine Local
946 Lymph Node Assay (LLNA) I: The ICCVAM Review Process. Regulatory Toxicology and
947 Pharmacology 34:249-257.
- 948 Takeyoshi M, Yamasaki K, Yakabe Y, Takatsuki M, Kimber I. 2001. Development of non-radio
949 isotopic endpoint of murine locallymph node assay based on 5-bromo-2%-deoxyuridine (BrdU)
950 Incorporation. Toxicol Lett 119:203–208.
- 951 Takeyoshi M, Sawaki M, Yamasaki K, Kimber I. 2003. Assessment of statistic analysis in non-
952 radioisotopic local lymph node assay (non-RI-LLNA) with alpha-hexylcinnamic aldehyde as an
953 example. Toxicology 191:259-63.
- 954 Takeyoshi M, Noda S, Yamazaki S, Kakishima H, Yamasaki K, Kimber I. 2004a. Assessment of
955 the skin sensitization potency of eugenol and its dimers using a non-radioisotopic modification
956 of the local lymph node assay. J Appl Toxicol 24:77-81.
- 957 Takeyoshi M, Noda S, Yamasaki K. 2004b. Differences in responsiveness of mouse strain
958 against p-benzoquinone as assessed by non-radioisotopic murine local lymph node assay. Exp
959 Anim 53:171-3.
- 960 Takeyoshi M, Iida K, Shiraishi K, Hoshuyama S. 2005. Novel approach for classifying
961 chemicals according to skin sensitizing potency by non-radioisotopic modification of the local
962 lymph node assay. J Appl Toxicol 25:129-34.
- 963 Takeyoshi M, Noda S, Yamasaki K, Kimber I. 2006. Advantage of using CBA/N strain mice in a
964 non-radioisotopic modification of the local lymph node assay. J Appl Toxicol 26:5-9.
- 965 Takeyoshi M, Iida K, Suzuki K, Yamazaki S. 2007a. Skin sensitization potency of isoeugenol
966 and its dimers evaluated by a non-radioisotopic modification of the local lymph node assay and
967 guinea pig maximization test. J Appl Toxicol DOI: 10.1002/jat.1305 [Online 10 October 2007].

- 968 Takeyoshi M, Minobe Y, Takatsuki M. 2007b. Promising endpoint and assay performance of
969 non-radioisotopic Local Lymph Node Assay based on BrdU incorporation. Presentation at the 6th
970 World Congress on Alternatives and Animal Use in the Life Sciences. Tokyo, Japan.
- 971 Yamano T, Shimizu M, Noda T. 2003. Allergenicity evaluation of p-chloro-m-cresol and p-
972 chloro-m-xyleneol by non-radioactive murine local lymph-node assay and multiple-dose guinea
973 pig maximization test. *Toxicology* 190:259-66.
- 974 Yamano T, Shimizu M, Noda T. 2004. Allergenicity evaluation of Bioban CS-1135 in
975 experimental animals. *Contact Dermatitis* 50:339-43.
- 976 Yamano T, Shimizu M, Noda T. 2005. Quantitative comparison of the results obtained by the
977 multiple-dose guinea pig maximization test and the non-radioactive murine local lymph-node
978 assay for various biocides. *Toxicology* 211:165-75.
- 979 Yamano T, Shimizu M, Noda T. 2006. Allergenicity and cross-reactivity of naphthenic acid and
980 its metallic salts in experimental animals. *Contact Dermatitis* 54:25-8.
- 981 Yamano T, Ichihara M, Shimizu M, Noda T, Tsujimoto Y. 2007. Immunomodulatory effects of
982 mono-, di-, and trimethylphenols in mice. *Toxicology* 232:132-7.