

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25

LUMI-CELL[®] ER ASSAY

ANTAGONIST PROTOCOL

**National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
Toxicological Methods (NICEATM)**

Developed by:

Xenobiotic Detection Systems, Inc.

1601 E. Geer St., Suite S

Durham, NC 2770413

17 October 2008

26	TABLE OF CONTENTS	
27	LIST OF ACRONYMS AND ABBREVIATIONS	vi
28	LIST OF FIGURES	viii
29	LIST OF TABLES	ix
30	1.0 Purpose	1
31	2.0 Sponsor	1
32	2.1 Substance Inventory and Distribution Management.....	3
33	3.0 Definitions.....	3
34	4.0 Testing Facility and Key Personnel	4
35	4.1 Testing Facility	4
36	4.2 Key Personnel.....	4
37	5.0 Identification of Test and Control Substances	4
38	5.1 Test Substances.....	4
39	5.2 Controls.....	4
40	6.0 Overview of General Procedures For Antagonist Testing.....	5
41	6.1 Range Finder Testing.....	6
42	6.2 Comprehensive Testing.....	7
43	7.0 Materials for LUMI-CELL® ER Antagonist Testing.....	7
44	7.1 BG1Luc4E2 Cells.....	7
45	7.2 Technical Equipment	8
46	7.3 Reference Standard, Controls, and Tissue Culture Supplies.....	9
47	8.0 Preparation of Tissue Culture Media and Solutions.....	11
48	8.1 RPMI 1640 Growth Medium (RPMI).....	11
49	8.2 Estrogen-Free DMEM.....	11

50	8.3	1X Trypsin Solution.....	12
51	8.4	1X Lysis Solution	12
52	8.5	Reconstituted Luciferase Reagent	12
53	9.0	Overview of Propagation and Experimental Plating of	
54		BG1Luc4E2 Cells	13
55	9.1	Procedures for Thawing Cells and Establishing Tissue Cultures.....	13
56	9.1.1	Thawing Cells	13
57	9.1.2	Establishing Tissue Cultures.....	14
58	9.2	Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free	
59		Medium, and Plating Cells for Experimentation.....	16
60	9.2.1	Ongoing Tissue Culture Maintenance	17
61	9.2.2	Conditioning in Estrogen-free Medium.....	18
62	9.2.3	Plating Cells Grown in Estrogen-free DMEM for	
63		Experimentation	18
64	10.0	Preparation of Test Substances.....	21
65	10.1	Determination of Test Substance Solubility.....	21
66	11.0	Preparation of Reference Standard, Control, and Test Substance	
67		Stock Solutions for Range Finder and Comprehensive Testing	22
68	11.1	Preparation of Ral/E2 Stock Solutions	22
69	11.1.1	E2 Stock Solution.....	22
70	11.1.2	Raloxifene Stock Solution	23
71	11.2	Ral/E2 Range Finder Testing Stock.....	23
72	11.2.1	Raloxifene Dilutions.....	23
73	11.2.2	Preparation of Ral/E2 Range Finder Working Stocks.....	24
74	11.3	Ral/E2 Comprehensive Testing Stock	24

75	11.3.1 Raloxifene Dilutions.....	24
76	11.3.2 Preparation of Ral/E2 Comprehensive Testing Working	
77	Stocks	25
78	11.4 Flavone/E2 Stock Solution.....	26
79	12.0 Preparation of Reference Standard, Control, and Test Substance	
80	Dosing Solutions for Range Finder and Comprehensive Testing	26
81	12.1 Preparation of Reference Standard and Control Dosing Solutions	
82	for Range Finder Testing	26
83	12.1.1 Preparation of Ral/E2 Reference Standard Range Finder	
84	Dosing Solutions	26
85	12.1.2 Preparation of DMSO Control Range Finder Dosing	
86	Solution.....	27
87	12.1.3 Preparation of E2 Control Range Finder Dosing	
88	Solution.....	27
89	12.2 Preparation of Test Substance Dosing Solutions for Range	
90	Finder Testing.....	27
91	12.3 Preparation of Reference Standard and Control Dosing	
92	Solutions for Comprehensive Testing.....	29
93	12.3.1 Preparation of Ral/E2 Reference Standard Dosing	
94	Solutions for Comprehensive Testing	29
95	12.3.2 Preparation of DMSO Control Comprehensive	
96	Testing Dosing Solution	30
97	12.3.3 Preparation of E2 Control Comprehensive Testing	
98	Dosing Solution.....	30
99	12.3.4 Preparation of Flavone/E2 Control Comprehensive	
100	Testing Dosing Solution	30
101		

101	12.4	Preparation of Test Substance Dosing Solutions for	
102		Comprehensive Testing.....	30
103	12.4.1	Preparation of Test Substance 1:2 Dilution for	
104		Comprehensive Testing.....	30
105	12.4.2	Preparation of Test Substance 1:5 Dilution for	
106		Comprehensive Testing.....	31
107	13.0	General Procedures for the Testing of Coded Substances	32
108	13.1	Application of Reference Standard, Control, and Test Substances	32
109	13.1.1	Preparation of Excel® Data Analysis Template for	
110		Range Finder Testing	33
111	13.1.2	Preparation of Excel® Data Analysis Template for	
112		Comprehensive Testing	34
113	13.2	Visual Evaluation of Cell Viability	35
114	13.3	Lysis of Cells for LUMI-CELL® ER.....	36
115	13.4	Measurement of Luminescence.....	36
116	13.5	Data Analysis.....	36
117	13.5.1	Collection and Adjustment of Luminometer Data for Range	
118		Finder Testing.....	37
119	13.5.2	Collection and Adjustment of Luminometer Data for	
120		Comprehensive Testing	39
121	13.5.3	Determination of Outliers.....	41
122	13.5.4	Acceptance Criteria.....	42
123	13.5.4.1	Range Finder Testing.....	42
124	13.5.4.2	Comprehensive Testing.....	42
125	14.0	Range Finder Testing	43
126	15.0	Comprehensive Testing.....	48

127	16.0	Compilation of the Historical Quality Control Database	50
128	16.1	E2 Control	51
129	16.2	DMSO Control	51
130	17.0	Quality Testing of Materials	51
131	17.1	Tissue Culture Media.....	52
132	17.2	G418.....	53
133	17.3	DMSO	53
134	17.4	Plastic Tissue Culture Materials.....	54
135	18.0	References.....	55
136			
137			

137	LIST OF ACRONYMS AND ABBREVIATIONS	
138	13 mm test tube	13 x 100 mm glass test tubes
139	DMEM	Dulbecco's Modification of Eagle's Medium
140	DMSO	Dimethyl Sulfoxide
141	DMSO Control	1% v/v dilution of DMSO in tissue culture media
142		used as a vehicle control
143	E2	17β-estradiol
144	E2 Control	2.5 x 10 ⁻⁵ μg/mL E2 used as a control.
145	IC ₅₀ Value	Concentration that produces a half-maximal response as
146		calculated using the four parameter Hill function.
147	ER	Estrogen Receptor
148	Estrogen-free DMEM	DMEM (phenol red free), supplemented with 1 %
149		Penicillin/Streptomycin, 2 % L-Glutamine, and 5%
150		Charcoal-dextran treated FBS
151	FBS	Fetal Bovine Serum
152	Flavone/E2 Control	25 μg/mL flavone + 2.5 x 10 ⁻⁵ μg/mL E2,
153		used as a weak positive control.
154	G418	Gentamycin
155	Ral/E2 Reference Standard	Nine point dilution of raloxifene HCl + 2.5 x 10 ⁻⁵ 17β-
156		estradiol reference standard for the LUMI-CELL® ER
157		antagonist assay
158	RPMI	RPMI 1640 growth medium
159	TA	Transcriptional Activation

160	T25	25 cm ² tissue culture flask
161	T75	75 cm ² tissue culture flask
162	T150	150 cm ² tissue culture flask
163		

163

LIST OF FIGURES

164 Figure 7-1 pGudLuc7.ERE Plasmid 7

165 Figure 9-1 Hemocytometer Counting Grid 19

166 Figure 14-1 Antagonist Range Finder Test Plate Layout 44

167 Figure 14-2 Antagonist Range Finder (example 1) 46

168 Figure 14-3 Antagonist Range Finder (example 2) 47

169 Figure 14-4 Antagonist Range Finder (example 3) 47

170 Figure 14-5 Antagonist Range Finder (example 4) 48

171 Figure 15-1 Antagonist Comprehensive Test Plate Layout 49

172

172	LIST OF TABLES	
173	Table 6-1	Concentration of Ral/E2 Reference Standard Used for
174		Comprehensive Testing..... 5
175	Table 11-1	Preparation of E2 Stock Solution 23
176	Table 11-2	Preparation of Raloxifene Stock Solution..... 23
177	Table 11-3	Preparation of Raloxifene Dilutions for Range Finder Testing..... 24
178	Table 11-4	Concentrations of Raloxifene and E2 in the Ral/E2 Range
179		Finder Working Stocks 24
180	Table 11-5	Preparation of Raloxifene Dilutions for Comprehensive Testing 25
181	Table 11-6	Concentrations of Raloxifene and E2 in the Ral/E2 Working
182		Stock 26
183	Table 12-1	Preparation of Test Substance Serial Dilution for Range Finder
184		Testing..... 27
185	Table 12-2	Addition of E2 to Test Substance Serial Dilution for Range
186		Testing..... 28
187	Table 12-3	Preparation of Ral/E2 Reference Standard Dosing Solution for
188		Comprehensive Testing..... 29
189	Table 12-4	Preparation of Test Substance 1:2 Dilutions for Comprehensive
190		Testing..... 31
191	Table 12-5	Preparation of Test Substance 1:5 Dilutions for Comprehensive
192		Testing..... 31
193	Table 13-1	Visual Observation Scoring..... 36
194	Table 13-2	Q Test Values 42
195		

196 **1.0 PURPOSE**

197 This protocol is designed to evaluate coded test substances for potential estrogen receptor (ER)
198 antagonist activity using the LUMI-CELL® ER assay.

199 **2.0 SPONSOR**

200 The National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
201 Toxicological Methods (NICEATM), P.O. Box 12233 Research Triangle Park, NC 27709

202 William S. Stokes, DVM, DACLAM

203 Rear Admiral, U.S. Public Health Service

204 Chief Veterinary Officer, USPHS

205 Director, NICEATM

206 National Institute of Environmental Health Sciences, NIH, DHHS

207 Bldg. 4401, Room 3129, MD EC-14

208 79 T.W. Alexander Drive

209 Research Triangle Park, NC 27709

210 Phone: 919-541-7997

211 Fax: 919-541-0947

212 Email: stokes@niehs.nih.gov

213

214 Raymond Tice, Ph.D.

215 Deputy Director, NICEATM

216 National Institute of Environmental Health Sciences

217 MD EC-17, P.O. Box 12233

218 Research Triangle Park, NC 27709

219 Phone: 919-541-4482

220 FAX: 919-541-0947

221 Email: tice@niehs.nih.gov

222

223

223 David Allen, Ph.D.
224 Principal Investigator
225 ILS, Inc./Contractor supporting NICEATM
226 National Institute of Environmental Health Sciences
227 MD EC-17, P.O. Box 12233
228 Research Triangle Park, NC 27709
229 Phone: 919-316-4587
230 FAX: 919-541-0947
231 Email: allen7@niehs.nih.gov
232
233 Frank Deal, M.S.
234 Staff Toxicologist
235 ILS, Inc./Contractor supporting NICEATM
236 National Institute of Environmental Health Sciences
237 MD EC-17, P.O. Box 12233
238 Research Triangle Park, NC 27709
239 Phone: 919-316-4587
240 FAX: 919-541-0947
241 Email: dealf@niehs.nih.gov
242
243 Patricia Ceger, M.S.
244 Project Coordinator/Technical Writer
245 ILS, Inc./Contractor supporting NICEATM
246 National Institute of Environmental Health Sciences
247 MD EC-17, P.O. Box 12233
248 Research Triangle Park, NC 27709
249 Phone: 919-316-4556
250 Fax: 919-541-0947
251 E-Mail: cegerp@niehs.nih.gov
252
253

253 **2.1 Substance Inventory and Distribution Management**

254 Cynthia Smith, Ph.D.

255 Chemistry Resources Group Leader

256 National Institute of Environmental Health Sciences

257 MD EC-06, P.O. Box 12233

258 Research Triangle Park, NC 27709

259 Phone: 919-541-3473

260

261 **3.0 DEFINITIONS**

262 • **Dosing Solution:** The test substance, control substance, or reference standard
263 solution which is to be placed into the tissue culture wells for experimentation.

264 • **Raw Data:** Raw data includes information that has been collected but not
265 formatted or analyzed, and consists of the following:

266 ○ Data recorded in the Study Notebook

267 ○ Computer printout of initial luminometer data

268 ○ Other data collected as part of GLP compliance, e.g.:

269 ■ Equipment logs and calibration records

270 ■ Test substance and tissue culture media preparation logs

271 ■ Cryogenic freezer inventory logs

272 • **Soluble:** Test substance exists in a clear solution without visible cloudiness or
273 precipitate.

274 • **Study Notebook:** The study notebook contains recordings of all activities related
275 to the conduct of the LUMI-CELL[®] ER TA antagonist assay.

276 • **Test Substances:** Substances supplied to the testing laboratories that are coded
277 and distributed such that only the Project Officer, Study Management Team
278 (SMT), and the Substance Inventory and Distribution Management have
279 knowledge of their true identity. The test substances will be purchased, aliquoted,

280 coded, and distributed by the Supplier under the guidance of the NIEHS/NTP
281 Project Officer and the SMT.

282 **4.0 TESTING FACILITY AND KEY PERSONNEL**

283 **4.1 Testing Facility**

284 Xenobiotic Detection Systems, Inc. (XDS), 1601 E. Geer St., Durham, NC 27704

285 **4.2 Key Personnel**

- 286 • Study Director: John Gordon, Ph.D.
- 287 • Quality Assurance Director: Mr. Carlos Daniel

288 **5.0 IDENTIFICATION OF TEST AND CONTROL SUBSTANCES**

289 **5.1 Test Substances**

290 Test substances are coded and will be provided to participating laboratories by the Substance
291 Inventory and Distribution Management team.

292 **5.2 Controls**

293 Controls for the ER antagonist protocol are as follows:

294 *Vehicle control (dimethyl sulfoxide [DMSO]):* 1% v/v dilution of DMSO (CASRN 67-68-5)
295 diluted in tissue culture media.

296 *Ral/E2 reference standard for range finder testing:* Three concentrations (1.56×10^{-3} ,
297 3.91×10^{-4} , and 9.77×10^{-5} $\mu\text{g/mL}$) of raloxifene HCl (Ral), CASRN 84449-90-1, plus a fixed
298 concentration (2.5×10^{-5} $\mu\text{g/mL}$) of 17β -estradiol (E2), CASRN: 50-28-2, in duplicate wells.

299 *Ral/E2 reference standard for comprehensive testing:* A serial dilution of Ral plus a fixed
300 concentration (2.5×10^{-5} $\mu\text{g/mL}$) of E2 consisting of nine concentrations of Ral/E2 in duplicate
301 wells.

302 *E2 control:* 17β -estradiol, 2.5×10^{-5} $\mu\text{g/mL}$ E2 in tissue culture media used as a base line
303 negative control.

304 *Flavone/E2 Control*: Flavone, CASRN 525-82-6, 25 µg/mL, with 2.5×10^{-5} µg/mL E2 in tissue
 305 culture media used as a weak positive control.

306 **6.0 OVERVIEW OF GENERAL PROCEDURES FOR ANTAGONIST TESTING**

307 All experimental procedures are to be carried out under aseptic conditions and all solutions,
 308 glassware, plastic ware, pipettes, etc., shall be sterile. All methods and procedures shall be
 309 documented in the study notebook.

310 Antagonist range finder testing is conducted on 96-well plates using three concentrations of
 311 Ral/E2 (1.56×10^{-3} , 3.91×10^{-4} , and 9.77×10^{-5} µg/mL Ral) with 2.50×10^{-5} µg/mL E2) in
 312 duplicate as the reference standard, with three replicate wells for the E2 and DMSO controls.

313 Comprehensive testing is conducted on 96-well plates using nine concentrations of Ral/E2 in
 314 duplicate as the reference standard (**Table 6-1**). Four replicate wells for the DMSO control,
 315 Flavone/E2 and E2 controls are included on each plate.

316 **Table 6-1 Concentrations of Ral/E2 Reference Standard**
 317 **Used for Comprehensive Testing**

Raloxifene Concentrations ¹	E2 Concentrations
1.25×10^{-2}	2.5×10^{-5}
6.25×10^{-3}	2.5×10^{-5}
3.13×10^{-3}	2.5×10^{-5}
1.56×10^{-3}	2.5×10^{-5}
7.81×10^{-4}	2.5×10^{-5}
3.91×10^{-4}	2.5×10^{-5}
1.95×10^{-4}	2.5×10^{-5}
9.77×10^{-5}	2.5×10^{-5}
4.88×10^{-5}	2.5×10^{-5}

318 ¹Concentrations are presented in µg/mL.

319
 320 Visual observations for cell viability are conducted for all experimental plates just prior to
 321 LUMI-CELL® ER evaluation, as outlined in **Section 11.4**.

322 Luminescence data, measured in relative light units (RLUs), is corrected for background
 323 luminescence by subtracting the mean RLU value of the vehicle control (DMSO) wells from the

324 RLU measurements for each of the other wells of the 96-well plate. Data is then transferred into
325 Excel® data management spreadsheets and GraphPad PRISM® 4.0 statistical software, graphed,
326 and evaluated for a positive or negative response as follows:

- 327 • A response is considered positive for antagonist activity when the average
328 adjusted RLU for a given concentration is less than the mean RLU value minus
329 three times the standard deviation for the E2 control.
- 330 • Any luminescence at or above this threshold is considered a negative response.

331 For substances that are positive at one or more concentrations, the concentration of test substance
332 that causes a half-maximal response (the relative IC₅₀) is calculated using a Hill function
333 analysis. The Hill function is a four-parameter logistic mathematical model relating the
334 substance concentration to the response (typically following a sigmoidal curve) using the
335 equation below

$$336 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log \text{IC}_{50} - X) \text{HillSlope}}}$$

337 where Y = response (i.e., relative light units); X = the logarithm of concentration; Bottom = the
338 minimum response; Top = the maximum response; log IC₅₀ = the logarithm of X as the response
339 midway between Top and Bottom; and HillSlope describes the steepness of the curve. The model
340 calculates the best fit for the Top, Bottom, HillSlope, and IC₅₀ parameters. See **Section 13.6.5** for
341 more details.

342 Acceptance or rejection of a test is based on evaluation of reference standard and control results
343 from each experiment conducted on a 96-well plate. Results for these controls are compared to
344 historical results compiled in the historical database, as seen in **Section 16.0**.

345 **6.1 Range Finder Testing**

346 Antagonist range finding for coded substances consists of a seven-point 1:10 serial dilution using
347 duplicate wells per concentration. Concentrations for comprehensive testing are selected based
348 on the response observed in range finder testing. If necessary, a second range finder test can be
349 conducted to clarify the optimal concentration range to test (see **Section 14.0**).

350

350 6.2 Comprehensive Testing

351 Comprehensive antagonist testing for coded substances consists of 11-point serial dilutions, with
 352 each concentration tested in triplicate wells of the 96-well plate. Three separate experiments are
 353 conducted for comprehensive testing on three separate days, except during Phases III and IV of
 354 the validation effort, in which comprehensive testing experiments are conducted once (see
 355 **Section 15.0**).


356 7.0 MATERIALS FOR LUMI-CELL® ER ANTAGONIST TESTING

357 This section provides the materials needed to conduct LUMI-CELL® ER testing, with associated
 358 brand names/vendors¹ in brackets.

359 7.1 BG1Luc4E2 Cells:

360 Human ovarian cancer cell line stably transfected with a plasmid containing an estrogen response
 361 element (**Figure 7-1**) [XDS].

362 **Figure 7-1 pGudLuc7.ERE Plasmid.**



363

364

¹Brand names and vendors should not be considered an endorsement by the U.S. Government or any member of the U.S. Government; such information is provided as examples.

364 7.2 Technical Equipment:

365 All technical equipment may be obtained from Fisher Scientific International, Inc. (Liberty Lane
366 Hampton, NH, USA 03842). Equivalent technical equipment from another commercial source
367 can be used.

- 368 • Analytical balance (Cat. No. 01-910-320)
- 369 • Berthold Orion 1 Microplate Luminometer [Berthold CatNo.: Orion 1 MPL3] or
370 equivalent and dedicated computer
- 371 • Biological safety hood, class II, and stand (Cat. No. 16-108-99)
- 372 • Centrifuge (low speed, tabletop with swinging bucket rotor) (Cat. No. 04-978-50
373 centrifuge, and 05-103B rotor)
- 374 • Combustion test kit (CO₂ monitoring) (Cat. No. 10-884-1)
- 375 • Drummond diaphragm pipetter (Cat. No. 13-681-15)
- 376 • Freezers, -20°C (Cat. No. 13-986-150), and -70°C (Cat. No. 13-990-86)
- 377 • Hand tally counter (Cat. No. 07905-6)
- 378 • Hemocytometer, cell counter (Cat. No. 02-671-5)
- 379 • Light microscope, inverted (Cat. No. 12-561-INV)
- 380 • Light microscope, upright (Cat. No. 12-561-3M)
- 381 • Liquid nitrogen flask (Cat. No. 11-675-92)
- 382 • Micropipetter, repeating (Cat. No. 21-380-9)
- 383 • Pipetters, air displacement, single channel (0.5 –10µl (Cat. No. 21-377-191), 2 –
384 20 µl (Cat. No. 21-377-287), 20 – 200 µl (Cat. No. 21-377-298), 200 - 1000 µl
385 (Cat. No. 21-377-195))
- 386 • Refrigerator/freezer (Cat. No. 13-986-106A)
- 387 • Shaker for 96-well plates (Cat. No. 14-271-9)
- 388 • Sodium hydroxide (Cat. No. 5318-500)

- 389 • Sonicating water bath (Cat. No. 15-335-30)
 - 390 • Tissue culture incubator with CO₂ and temperature control (Cat. No. 11-689-4)
 - 391 • Vacuum pump with liquid trap (side arm Erlenmeyer) (Cat. No. 01-092-29)
 - 392 • Vortex mixer (Cat. No. 12-814)
- 393 Equipment should be maintained and calibrated as per GLP guidelines and individual laboratory
394 SOPs.

395 **7.3 Reference Standard, Controls, and Tissue Culture Supplies**

396 All tissue culture reagents must be labeled to indicate source, identity, storage conditions and
397 expiration dates. Tissue culture solutions must be labeled to indicate concentration, stability
398 (where known), and preparation and expiration dates.

399 Equivalent tissue culture media and sera from another commercial source can be used, but must
400 first be tested as described in **Section 17.0** to determine suitability for use in this test method.

401 The following are the necessary tissue culture reagents and possible sources based on their use in
402 the pre-validation studies:

- 403 • BackSeal-96/384, white adhesive bottom seal for 96-well and 384-well microplate
404 [Perkin-Elmer, Cat. No. 6005199]
- 405 • 17 β -estradiol (CAS RN: 50-28-2) [Sigma-Aldrich, Cat. No. E8875]
- 406 • CellTiter-Glo® Luminescent Cell Viability Assay [Promega Cat. No. G7572]
- 407 • Cryovial, 2 mL (Corning Costar) [Fisher Scientific Cat. No. 03-374-21]
- 408 • Culture tube 13 x 100mm (case) [Thomas Scientific Cat. No.: 10009186R38]²
- 409 • Culture tube, 50 mL conical (Corning Costar) [Fisher Scientific Cat. No. 05-
410 526C]
- 411 • DMSO, U.S.P. analytical grade. [Sigma-Aldrich, Cat. No. 34869-100ML]

²If glass tubes can not be obtained from Thomas Scientific, the preference is for flint glass, then lime glass, then borosilicate glass.

- 412 • Dulbecco's Modification of Eagle's Medium (DMEM), containing 4.5 g/L
413 glucose, with sodium pyruvate, without phenol red or L-glutamine
414 [Mediatech/Cellgro, Cat. No. 17-205-CV]
- 415 • Fetal Bovine Serum [Mediatech/Cellgro Cat. No. MT 35-010-CV]
- 416 • Fetal Bovine Serum, charcoal/dextran treated, triple 0.1 µm sterile filtered
417 [Hyclone, Cat. No. SH30068.03]
- 418 • Flavone (CASRN: 525-82-6) [Sigma-Aldrich, Cat. No. F2003]
- 419 • Gentamycin Sulfate (G418), 50 mg/mL [Mediatech/Cellgro Cat. No. 30-234-CR]
- 420 • L-glutamine, 29.2 mg/mL [Cellgro, Cat. No. 25005-CI]
- 421 • Luciferase Assay System (10-Pack) [Promega Cat. No. E1501]
- 422 • Lysis Solution 5X [Promega, Cat. No. E1531]
- 423 • Penicillin/streptomycin solution, 5000 I.U. penicillin, 5000 µg/mL streptomycin
424 [Cellgro, Cat. No. 30-001-CI].
- 425 • Phosphate buffered saline (PBS, 1X) without calcium and magnesium [Cellgro,
426 Cat. No. 21-040-CV]
- 427 • Pipettes, serological: 2.0 mL [Sigma-Aldrich, Cat. No. P1736], 5.0 mL [Sigma-
428 Aldrich, Cat. No. P1986], 25 mL [Sigma-Aldrich, Cat. No. P2486]
- 429 • Raloxifene (CASRN 84449-90-1) [Sigma-Aldrich Cat. No. R1402]
- 430 • RPMI 1640 medium, containing L-glutamine [Mediatech, Cat. No. 10-040-CV]
- 431 • Tissue culture flasks (Corning-Costar): 25 cm² (T25) [Fisher Cat. No. 10-126-28];
432 75 cm² (T75) [Fisher Cat. No. 10-126-37]; and 150 cm² (T150) [Fisher Cat. No.
433 10-126-34]
- 434 • Tissue culture plates (Corning-Costar): 96-well [Thomas Scientific Cat. No.
435 6916A05]
- 436 • Trypsin (10X), 2.5% in Hank's balanced salt solution (HBSS), without calcium
437 and magnesium, without phenol red [Cellgro, Cat. No. 25-054-CI].

438 All reagent lot numbers and expiration dates must be recorded in the study notebook.

439 **8.0 PREPARATION OF TISSUE CULTURE MEDIA AND SOLUTIONS**

440 All tissue culture media and media supplements must be quality tested before use in experiments
441 (see **Section 15.0**).

442 **8.1 RPMI 1640 Growth Medium (RPMI)**

443 RPMI 1640 is supplemented with 0.9% Pen-Strep and 8.0% FBS to make RPMI growth medium
444 (RPMI).

445 Procedure for one 549 mL bottle:

- 446 1. Remove FBS from -70°C freezer, and Pen-Strep from -20°C freezer and allow to
447 equilibrate to room temperature.
- 448 2. Add 44 mL of FBS and 5 mL Pen-Strep to the bottle of RPMI 1640.
- 449 3. Label RPMI bottle as indicated in **Section 7.3**

450 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
451 *component.*

452 **8.2 Estrogen-Free DMEM Medium**

453 DMEM is supplemented to contain 4.5% charcoal/dextran treated FBS, 1.9% L-glutamine, 0.9%
454 Pen-Strep.

455 Procedure for one 539 mL bottle:

- 456 1. Remove charcoal/dextran treated FBS from -70°C freezer, and L-glutamine and
457 Pen-Strep from -20°C freezer and allow to equilibrate to room temperature.
- 458 2. Add 24 mL of charcoal/dextran treated FBS, 10 mL L-glutamine, and 5 mL Pen-
459 Strep to one 500 mL bottle of DMEM.
- 460 3. Label estrogen-free DMEM bottle as indicated in **Section 7.3**

461 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
462 *component.*

463 **8.3 1X Trypsin Solution**

464 1X Trypsin solution is prepared by dilution from a 10X premixed stock solution. The 10X stock
465 solution should be stored in 10 mL aliquots in a -20°C freezer.

466 Procedure for making 100 mL of 1X trypsin:

- 467 1. Remove a 10mL aliquot of 10X trypsin from -20°C freezer and allow to
468 equilibrate to room temperature.
- 469 2. Aliquot 1 mL Trypsin (10X) along with 9 mL of 1X PBS into ten 15 mL
470 centrifuge tubes.
- 471 3. Label 1X trypsin aliquots as indicated in **Section 7.3**

472 *1X Trypsin should be stored at -20°C.*

473 **8.4 1X Lysis Solution**

474 Lysis solution is prepared by dilution from a 5X premixed stock solution. Both the 5X and 1X
475 solutions can be repeatedly freeze-thawed.

476 The procedure for making 10 mL of 1X lysis solution:

- 477 1. Thaw the 5X Promega Lysis solution and allow it to reach room temperature.
- 478 2. Remove 2 mL of 5X solution and place it in a 15 mL conical centrifuge tube.
- 479 3. Add 8 mL of distilled, de-ionized water to the conical tube.
- 480 4. Cap and shake gently until solutions are mixed.

481 *Store at -20°C for no longer than 1 year from receipt.*

482 **8.5 Reconstituted Luciferase Reagent**

483 Luciferase reagent consists of two components, luciferase buffer and lyophilized luciferase
484 substrate.

485 For long-term storage, unopened containers of the luciferase buffer and lyophilized luciferase
486 substrate can be stored at -70°C for up to six months.

487 To reconstitute luciferase reagent:

- 488 1. Remove luciferase buffer and luciferase substrate from -70°C freezer and allow
489 them to equilibrate to room temperature.
- 490 2. Add 10 mL of luciferase buffer solution to luciferase substrate container and swirl
491 or vortex to mix, the Luciferase substrate should readily go into solution.
- 492 3. Luciferase substrate should readily go into solution.
- 493 4. After solutions are mixed aliquot to a 15mL centrifuge tube.
- 494 5. Store complete solution at -20°C.

495 Reconstituted luciferase reagent is stable for 1 month at -20°C.

496 **9.0 OVERVIEW OF PROPAGATION AND EXPERIMENTAL PLATING OF** 497 **BG1Luc4E2 CELLS**

498 The BG1Luc4E2 (BG-1) cells are stored in liquid nitrogen in 2 mL cryovials. BG-1 cells are
499 grown as a monolayer in tissue culture flasks in a dedicated tissue culture incubator at 37°C ±
500 1°C, 90% ± 5% humidity, and 5.0% ± 1% CO₂/air. The cells should be examined on a daily basis
501 during working days under an inverted phase contrast microscope, and any changes in
502 morphology and adhesive properties must be noted in the study notebook.

503 Two T150 flasks containing cells at 80% to 90% confluence will usually yield a sufficient
504 number of cells to fill three 96-well plates for use in experiments.

505 **9.1 Procedures for Thawing Cells and Establishing Tissue Cultures**

506 Warm all tissue culture media and solutions to room temperature by placing them under the
507 tissue culture hood several hours before use.

508 All tissue culture media, media supplements, and tissue culture plasticware must be quality
509 tested before use in experiments (**Section 17.0**).

510 9.1.1 Thawing Cells

- 511 1. Remove a cryovial of frozen BG-1 cells from the liquid nitrogen flask.
- 512 2. Facilitate rapid thawing by loosening the top slightly (do not remove top) to
513 release trapped gasses and retightening it. Roll vial between palms.

- 514 3. Use a micropipette to transfer cells to a 50 mL conical centrifuge tube.
- 515 4. Rinse cryovial twice with 1X PBS and add PBS rinse material to the conical tube.
- 516 5. Add 20 mL of RPMI to the conical tube.
- 517 6. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
518 for an additional 5 minutes.
- 519 7. Aspirate media from pellet and re-suspend it in 5 mL RPMI, drawing the pellet
520 repeatedly through a 1.0 mL serological pipette to break up any clumps of cells.
- 521 8. Transfer cells to a T25 flask, place them in an incubator (see conditions in
522 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

523 9.1.2 Establishing Tissue Cultures

524 Once cells have reached 80% to 90% confluence, transfer the cells to a T75 flask by performing,
525 for example, the following steps:

- 526 1. Remove the T25 flask from the incubator.
- 527 2. Aspirate the RPMI, then add 5 mL 1X PBS, making sure that the cells are coated
528 with PBS.
- 529 3. Aspirate 1X PBS, then add 1 to 2 mL 1X trypsin to the T25 flask, gently swirling
530 the flask to coat all cells with the trypsin.
- 531 4. Place the flask in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 532 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
533 hand.
- 534 6. Confirm cell detachment by examination under an inverted microscope. If cells
535 have not detached, return the flask to the incubator for an additional 2 minutes,
536 then hit the flask again.
- 537 7. After cells have detached, add 5 mL PBS, and transfer the suspended cells to a 50
538 mL centrifuge tube. Wash the flask one additional time with 5 mL PBS.
- 539 8. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
540 digestion by residual trypsin.

- 541 9. Pellet the cells by centrifugation, as described in **Section 9.1.1**, and re-suspend the
542 cells in 10 mL RPMI medium.
- 543 10. Draw the pellet repeatedly through a 25 mL serological pipette to break up
544 clumps of cells
- 545 11. Transfer cells to a T75 flask, then place the flask in an incubator (see conditions
546 in **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 547 When cells have reached 80% to 90% confluency, transfer them into a T150 flask by performing,
548 for example, the following steps:
- 549 12. Remove the T75 flask from the incubator, aspirate the old media and add 5 mL
550 1X PBS.
- 551 13. Aspirate 1X PBS, add 2 mL of 1X trypsin to the flask, and place it in an incubator
552 (see conditions in **Section 9.0**) for 5 to 10 min.
- 553 14. Repeat **steps 5** through **11** in **Section 9.1.2**, re-suspending the pellet in 20 mL of
554 RPMI.
- 555 15. Transfer cells to a T150 flask and place it in the incubator (see conditions in
556 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 557 16. Remove the T150 flask from the incubator.
- 558 17. Aspirate the RPMI and add 5 mL 1X PBS.
- 559 18. Aspirate 1X PBS and add 3 mL 1X trypsin to the T150 flask, making sure that the
560 cells are coated with the trypsin.
- 561 19. Incubate cells in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 562 20. Detach cells by hitting the side of the flask sharply against the palm or heel of the
563 hand.
- 564 21. Confirm cell detachment by examination under an inverted microscope. If cells
565 have not detached, return the flask to the incubator for an additional 2 minutes,
566 then hit the flask again.

- 567 22. After cells have detached, add 5mL 1X PBS and transfer the suspended cells from
568 the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the flask,
569 then transfer to the 50 mL conical tube.
- 570 23. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
571 digestion by residual trypsin.
- 572 24. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
573 centrifuge for an additional 5 minutes.
- 574 25. Aspirate the media from the pellet and re-suspend it in 40 mL RPMI, drawing the
575 pellet repeatedly through a 25 mL serological pipette to break up any clumps of
576 cells.
- 577 26. Transfer 20 mL of cell suspension to each of two T150 flasks, place them in an
578 incubator (see conditions in **Section 9.0**) and grow to 80% to 90% confluence
579 (approximately 48 to 72 hrs).

580 **9.2 Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free Medium,**
581 **and Plating Cells for Experimentation**

582 The following procedure is used to condition the BG1Luc4E2 cells to an estrogen-free
583 environment prior to plating the cells in 96-well plates for analysis of estrogen dependent
584 induction of luciferase activity.

585
586 To start the tissue culture maintenance and estrogen-free conditioning, split the two T150 culture
587 flasks into four T150 flasks. Two of these flasks will be used for continuing tissue culture and
588 will use the RPMI media mentioned above. The other two flasks will be cultured in estrogen-free
589 DMEM for experimental use. Extra care must be taken to avoid contaminating the estrogen-free
590 cells with RPMI.

- 591 1. Remove both T150 flasks from the incubator.
- 592 2. Aspirate the medium and rinse the cells with 5 mL 1X PBS.
- 593 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
594 to coat all cells with the trypsin.

- 595 4. Incubate cells in the incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 596 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
- 597 hand.
- 598 6. Confirm cell detachment by examination under an inverted microscope. If cells
- 599 have not detached, return the flask to the incubator for an additional 2 minutes,
- 600 then hit the flask again.
- 601 7. After cells have detached, add 5 mL 1X PBS to the first T150 flask and transfer
- 602 the suspended cells to the second T150 flask.
- 603 8. Transfer the contents of both flasks to a 50 mL conical tube. Repeat step 7 with an
- 604 additional 5 mL 1X PBS and transfer to the 50 mL conical tube.
- 605 9. Immediately add 20 mL estrogen-free DMEM to the 50 mL conical tube to inhibit
- 606 further cellular digestion by residual trypsin.
- 607 10. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
- 608 centrifuge for an additional 5 minutes.
- 609 11. Aspirate media from pellet and re-suspend it in 4 mL estrogen-free DMEM,
- 610 drawing the pellet repeatedly through a 1 mL serological pipette to break up
- 611 clumps of cells.

612 At this point, cells are ready to be divided into the ongoing tissue culture and estrogen-free

613 conditioning groups.


614 9.2.1 Ongoing Tissue Culture Maintenance

- 615 1. Add 20 mL RPMI to two T150 flasks.
- 616 2. Add 220 µL G418 to the RPMI in the T150 flasks
- 617 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.
- 618 4. Place T150 flasks in tissue culture incubator (see conditions in **Section 9.0**) and
- 619 grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 620 5. Tissue culture medium may need to be changed 24 hours after addition of G418 to
- 621 remove cells that have died because they do not express reporter plasmid.

- 622 6. G418 does not need to be added to the flasks a second time.
- 623 7. Repeat **Section 9.2 steps 1-11** for ongoing tissue culture maintenance.
- 624 9.2.2 Conditioning in Estrogen-free Medium
- 625 1. Add 20 mL estrogen-free DMEM to two T150 flasks.
- 626 2. Add 150 µL G418 to the estrogen-free DMEM in the T150 flasks.
- 627 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.
- 628 4. Tissue culture medium may need to be changed 24 hours after addition of G418 to
- 629 remove cells that have died because they do not express reporter plasmid.
- 630 5. G418 does not need to be added to the flasks a second time.
- 631 6. Place the T150 flasks in the incubator (see conditions in **Section 9.0**) and grow to
- 632 80% to 90% confluence (approximately 48 to 72 hrs).
- 633 9.2.3 Plating Cells Grown in Estrogen-free DMEM for Experimentation
- 634 1. Remove the T150 flasks that have been conditioned in estrogen-free DMEM for
- 635 48 to 72 hours from the incubator.
- 636 2. Aspirate the medium, then rinse the cells with 5 mL 1X PBS.
- 637 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
- 638 to coat all cells with the trypsin.
- 639 4. Place the flasks in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 640 5. Detach cells by hitting the side of the flask sharply against the palm or the heel of
- 641 the hand.
- 642 6. Confirm cell detachment by examination under an inverted microscope. If cells
- 643 have not detached, return the flask to the incubator for 2 additional minutes, then
- 644 hit the flask again.
- 645 7. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
- 646 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
- 647 flask, then transfer to the 50 mL conical tube.

- 648 8. Immediately add 20 mL estrogen-free DMEM to each conical tube to inhibit
649 further cellular digestion by residual trypsin.
- 650 9. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
651 centrifuge for an additional 5 minutes.
- 652 10. Aspirate off the media from the pellet and re-suspend it in 20 mL DMEM,
653 drawing the pellet repeatedly through a 25 mL serological pipette to break up any
654 clumps of cells.
- 655 11. Pipette 15 μ L of the cell suspension into the “v” shaped slot on the
656 hemocytometer. Ensure that the solution covers the entire surface area of the
657 hemocytometer grid, and allow cells to settle before counting.
- 658 12. Using 100x magnification, view the counting grid.
- 659 13. The counting grid on the hemocytometer consists of nine sections, four of which
660 are counted (upper left, upper right, lower left, and lower right, see **Figure 9-1**).
661 Each section counted consists of four by four grids. Starting at the top left and
662 moving clockwise, count all cells in each of the four by four grids. Some cells
663 will be touching the outside borders of the square, but only count those that touch
664 the top and right borders of the square. This value is then used in the calculation
665 below to get to the desired concentration of 200,000 cells/mL.

Figure 9-1 Hemocytometer Counting Grid.



666

667 **The volume of each square is 10^{-4} mL, therefore:**
668 **Cells/mL = (average number per grid) x 10^{-4} mL. x 1/(starting dilution).**
669 **Starting dilution: 20mL (for T150 flasks)**
670
671 Harvested cells for a T150 flask are suspended in 20 mL of estrogen-free DMEM and sampled
672 for determination of concentration of cells/mL.

673

674 Example Calculation:

- 675 • Grids 1, 2, 3, and 4 are counted and provide the following data:
- 676 ○ 50, 51, 49, and 50: average number of cells per grid is equal to 50.

677 Cells/mL = 50 cells per grid ÷ 10^{-4} mL volume of grid = 50×10^4 cells/mL (or 500,000
678 cells/mL)

679 Total # of Cells Harvested = 500,000 cells/mL x 20 mL

680 Desired Concentration (or Concentration_{Final}) = 200,000 cells/mL

681 Formula: (Concentration_{Final} x Volume_{Final} = Concentration_{Initial} x Volume_{Initial})

682 Concentration_{Final} = 200,000 cells/mL

683 Concentration_{Initial} = 500,000 cells/mL

684 Volume_{Initial} = 20 mL

685 Volume_{Final} – to be solved for.

686 Therefore: 200,000 cells/mL x Volume_{Final} = 500,000 cells/mL x 20 mL

687 Solving for Volume_{Final} we find = 50 mL

688 Therefore, add 30 mL of estrogen-free DMEM to the cell suspension for a total volume of 50
689 mL, which will yield the desired concentration of 200,000 cells/mL for plating.

690 14. This dilution scheme will give a concentration of 200,000 cells/mL. 200 µL of
691 this cell suspension is used for each well of a 96-well plate (i.e., 40,000 cells per
692 well).

693 15. Remove a 96-well plate from its sterile packaging. Use a repeater pipetter to
694 pipette 200 µL of cell suspension into each well to be used for the testing of
695 coded substances, reference standard and controls (**note**: add 200 µL of estrogen-
696 free DMEM only to any wells not being used for testing).

697 16. Incubate plate(s) in an incubator (see conditions in **Section 9.0**) for a minimum of
698 24 hours, but no longer than 48 hours before dosing.

699 Two T150 flasks containing cells at 80% to 90% confluence will typically yield sufficient cells
700 to fill four 96-well plates (not including the perimeter wells).

701 **10.0 PREPARATION OF TEST SUBSTANCES**

702 The solvent used for dissolution of test substances is 100% DMSO. All test substances should be
703 allowed to equilibrate to room temperature before being dissolved and diluted. Test substance
704 solutions (except for reference standards and controls) should not be prepared in bulk for use in
705 subsequent tests. Test substances are to be used within 24 hours of preparation. Solutions should
706 not have noticeable precipitate or cloudiness.

707 All information on weighing, solubility testing, and calculation of final concentrations for test
708 substances, reference standards and controls is to be recorded in the study notebook.

709 **10.1 Determination of Test Substance Solubility**

- 710 1. Prepare a 200 mg/mL solution of the test substance in 100% DMSO in a 4 mL
711 conical tube.
- 712 2. Vortex to mix.
- 713 3. If the test substance does not dissolve at 200 mg/mL, prepare a 20 mg/mL
714 solution and vortex as above.
- 715 4. If the test substance does not dissolve at 20 mg/mL solution, prepare a 2 mg/mL
716 solution in a 4 mL conical tube and vortex as above.
- 717 5. If the test substance does not dissolve at 2 mg/mL, prepare a 0.2 mg/mL solution
718 in a 4 mL conical tube and vortex as above.

719 6. Continue testing, using 1/10 less substance in each subsequent attempt until test
720 substance is solubilized in DMSO.

721 Once a solution of test substance has been obtained that does not have any visible precipitate or
722 cloudiness in 100% DMSO, the solubility of the test substance must be determined in the 1%
723 DMSO/99% estrogen-free DMEM mixture used for LUMI-CELL® ER testing.

724 7. Add 2 µL of the highest concentration of the test substance/DMSO solution to a
725 13 mm test tube.

726 8. Add 400 µL estrogen-free DMEM to the test tube and vortex gently,

727 9. If cloudiness or precipitate develop, vortex for up to 10 minutes.

728 10. If vortexing does not dissolve test substance, sonicate test substance for up to 10
729 minutes.

730 11. If test substance has visible precipitate or is cloudy return to **Section 10.1 step 7**
731 to try the next lower concentration for the test substance.

732 The Testing Facility shall forward the results from the solubility tests assay to the SMT through
733 the designated contacts in electronic format and hard copy upon completion of testing.

734 **11.0 PREPARATION OF REFERENCE STANDARD, CONTROL AND TEST**
735 **SUBSTANCE STOCK SOLUTIONS FOR RANGE FINDER AND**
736 **COMPREHENSIVE TESTING**

737 All information on preparation of test substances, reference standards and controls is to be
738 recorded in the study notebook.

739 **11.1 Preparation of Ral/E2 Stock Solutions**

740 E2 and raloxifene stocks are prepared separately and then combined into Ral/E2 stocks, which
741 are then used to prepare dosing solutions in **Section 12**.

742 11.1.1 E2 Stock Solution

743 The final concentration of the E2 stock solution is 5.0×10^{-3} µg/mL. Prepare the E2 stock as
744 shown in **Table 11-1**.

745 **Table 11-1 Preparation of E2 Stock Solution**

Step #	Action	DMSO	E2 Concentration
1	Make a 10 mg/mL stock solution in 100% DMSO in a 4mL vial.	-	10 mg/mL
2	Transfer 10 µL E2 solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 10 µL E2 solution from Step #2 to a new 4mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	1 µg/mL
4	Transfer 100 µL E2 solution from Step #3 to a new glass container large enough to hold 15 mL.	Add 9.90 mL of 100% DMSO. Vortex to mix.	1.0×10^{-2} µg/mL
5	Transfer 5 mL E2 solution from Step #4 to a new glass container large enough to hold 15 mL	Add 5 mL of 100% DMSO. Vortex to mix.	5.0×10^{-3} µg/mL

746

747 11.1.2 Raloxifene Stock Solution748 Prepare a 2.5 µg/mL raloxifene working stock solution as shown in **Table 11-2**.749 **Table 11-2 Preparation of Raloxifene Stock Solution**

Step #	Action	DMSO	Raloxifene Concentration
1	Make a 10 mg/mL solution of raloxifene in a 4 mL glass vial.	-	1.0×10^4 µg/mL
2	Transfer 10 µL raloxifene solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 150 µL raloxifene solution from Step #2 to a new 4 mL vial.	Add 2.850 mL of 100% DMSO. Vortex to mix.	5 µg/mL
4	Transfer 1.5 mL raloxifene solution from Step #3 to a new 13 mm test tube.	Add 1.5 mL of 100% DMSO. Vortex to mix.	2.5 µg/mL

750

751 **11.2 Ral/E2 Range Finder Testing Stock**752 11.2.1 Raloxifene Dilutions

753 Number three 4 mL vials with the numbers 1 to 3 and use the raloxifene solution prepared in

754 **Section 11.1.2** to make raloxifene dilutions as shown **Table 11-3**.

755

755 **Table 11-3 Preparation of Raloxifene Dilutions for Range Finder Testing**

Step #	Action	DMSO	Raloxifene Concentration
1	Transfer 250 µL of the 2.5 µg/mL raloxifene working stock solution to a 4 mL tube	Add 750 µL of 100% DMSO and vortex	6.25×10^{-1} µg/mL
2	Transfer 500 µL of the 6.25×10^{-1} µg/mL raloxifene solution to a 4 mL tube	Add 500 µL of 100% DMSO and vortex	3.13×10^{-1} µg/mL
3	Transfer 250 µL of the 3.13×10^{-1} µg/mL raloxifene solution to a 4 mL tube	Add 750 µL of 100% DMSO and vortex	7.81×10^{-2} µg/mL
4	Transfer 125 µL of the 7.81×10^{-2} µg/mL raloxifene solution to a 4 mL tube	Add 375 µL of 100% DMSO and vortex	1.95×10^{-2} µg/mL

756

757 11.2.2 Preparation of Ral/E2 Range Finder Working Stocks:

758 Label three 4 mL conical tubes with numbers 1 through 3 and add 500 µL of the 5×10^{-3} µg/mL
 759 E2 solution prepared in **Section 11.1.1** to each tube. Add 500 µL of the 3.13×10^{-1} , 7.81×10^{-2} ,
 760 and 1.95×10^{-2} µg/mL raloxifene solutions prepared in **Section 11.2.1** to tubes 1, 2, and 3
 761 respectively. Vortex each tube to mix. The final concentrations for raloxifene and E2 are listed in
 762 **Table 11-4.**

763 **Table 11-4 Concentrations of Raloxifene and E2 in the**
764 **Ral/E2 Range Finder Working Stocks**

Tube #	Raloxifene (µg/ml)	E2 (µg/ml)
1	1.56×10^{-1}	2.5×10^{-3}
2	3.91×10^{-2}	2.5×10^{-3}
3	9.77×10^{-3}	2.5×10^{-3}

765

766 **11.3 Ral/E2 Comprehensive Testing Stock**767 11.3.1 Raloxifene Dilutions

768 Use the raloxifene solution prepared in **Section 11.1.2** to make a nine-point serial dilution of
 769 raloxifene as shown **Table 11-5.**

770

770 **Table 11-5 Preparation of Raloxifene Dilutions for Comprehensive Testing**

Step #	Action	DMSO	Discard	Raloxifene Concentration
1	Transfer 500 µL of the raloxifene working stock solution to a new 4 mL vial.	-	-	2.5 µg/mL
2	Transfer 500 µL of the raloxifene working stock solution to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.25 µg/mL
3	Transfer 500 µL raloxifene solution from Step #2 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	6.25×10^{-1} µg/mL
4	Transfer 500 µL raloxifene solution from Step #3 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	3.13×10^{-1} µg/mL
5	Transfer 500 µL raloxifene solution from Step #4 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.56×10^{-1} µg/mL
6	Transfer 500 µL raloxifene solution from Step #5 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	7.81×10^{-2} µg/mL
7	Transfer 500 µL raloxifene solution from Step #6 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	3.91×10^{-2} µg/mL
8	Transfer 500 µL raloxifene solution from Step #7 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.95×10^{-2} µg/mL
9	Transfer 500 µL raloxifene solution from Step #8 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	Discard 500 µL from Tube #9	9.77×10^{-3} µg/mL

771

772 11.3.2 Preparation of Ral/E2 Comprehensive Testing Working Stocks:

773 Add 500 µL of the 5×10^{-3} µg/mL E2 solution prepared in **Section 11.1.1** to each of the 9
774 raloxifene dilution vials (including the working stock solution in Tube #1). Vortex each tube to
775 mix. The final concentrations for raloxifene and E2 are listed in **Table 11-6**.

776

776 **Table 11-6 Concentrations of Raloxifene and E2 in the Ral/E2 Working Stocks**

Tube #	Raloxifene (µg/mL)	E2 (µg/mL)
1	1.25	2.5×10^{-3}
2	6.25×10^{-1}	2.5×10^{-3}
3	3.13×10^1	2.5×10^{-3}
4	1.56×10^{-1}	2.5×10^{-3}
5	7.81×10^2	2.5×10^{-3}
6	3.91×10^{-2}	2.5×10^{-3}
7	1.95×10^{-2}	2.5×10^{-3}
8	9.77×10^{-3}	2.5×10^{-3}
9	4.88×10^{-3}	2.5×10^{-3}

777

778 **11.4 Flavone/E2 Stock Solution**

779 To prepare the flavone/E2 stock solution, proceed as follows:

- 780 1. Prepare 1 mL of 5 mg/mL flavone
- 781 2. Add 1 mL of the 5×10^{-3} µg/mL E2 (prepared as in **Section 11.1.1**) to the 10
- 782 mg/mL flavone. This will make a working solution of 2.5 mg/mL flavone with
- 783 2.5×10^{-3} µg/mL E2.

784 **12.0 PREPARATION OF REFERENCE STANDARD, CONTROL AND TEST**

785 **SUBSTANCE DOSING SOLUTIONS FOR RANGE FINDER AND**

786 **COMPREHENSIVE TESTING**

787 **12.1 Preparation of Reference Standard and Control Dosing Solutions for Range**

788 **Finder Testing**

789 Range finder testing is conducted on 96-well plates using three concentrations of Ral/E2 in

790 duplicate as the reference standard. Three replicate wells for the DMSO, and E2 controls are

791 included on each plate.

792 All “dosing solutions” of test substance concentrations are to be expressed as µg/mL in the study

793 notebook and in all laboratory reports.

794 Dosing solutions are to be used within 24 hours of preparation.

795 **12.1.1 Preparation of Ral/E2 Reference Standard Range Finder Dosing Solutions**

- 796 1. Label three 13 mm glass tubes with the numbers 1 to 3.

- 797 2. Add 6 µL of Ral/E2 stock from tube #1 from **Section 11.2.2** to the 13 mm glass
798 test tube labeled #1.
- 799 3. Add 6 µL of Ral/E2 stock from tube #2 from **Section 11.2.2** to the 13 mm glass
800 test tube labeled #2. Repeat for tube #3.
- 801 4. Add 600 µL of estrogen-free DMEM to each tube and vortex.

802 12.1.2 Preparation of DMSO Control Range Finder Dosing Solution

- 803 1. Add 8 µL of 100% DMSO to a 13 mm glass test tube.
- 804 2. Add 800 µL of estrogen-free DMEM to each tube and vortex.

805 12.1.3 Preparation of E2 Control Range Finder Dosing Solution

- 806 1. Add 4 µL of the E2 stock from **Section 11.1.1** to a 13 mm glass test tube.
- 807 2. Add 4 µL of 100% DMSO to the tube.
- 808 3. Add 800 µL of estrogen-free DMEM to the tube and vortex to mix.

809 **12.2 Preparation of Test Substance Dosing Solutions for Range Finder Testing**

810 Range finder experiments are used to determine the concentrations of test substance to be used
811 during comprehensive testing. Antagonist range finding for coded substances consists of seven-
812 point 1:10 serial dilutions in duplicate.

813 To prepare test substance dosing solutions:

- 814 1. Label two sets of seven glass 13 mm test tubes with the numbers 1 through 7 and
815 place them in a test tube rack. Perform a serial dilution of test substance as shown
816 in **Table 12-1** using one set of tubes.

817 **Table 12-1 Preparation of Test Substance Serial Dilution for Range Finder Testing**

Tube #	100% DMSO	Test Substance ¹	Final Volume
1	-	100 µL of test substance solution from Section 10.1	100 µL
2	90 µL	10 µL of test substance solution from Section 10.1	100 µL
3	90 µL	10 µL from Tube #2	100 µL

Tube #	100% DMSO	Test Substance ¹	Final Volume
4	90 µL	10 µL from Tube #3	100 µL
5	90 µL	10 µL from Tube #4	100 µL
6	90 µL	10 µL from Tube #5	100 µL
7	90 µL	10 µL from Tube #6	100 µL

¹Vortex tubes #2 through 6 before removing test substance/DMSO solution to place in the next tube in the series.

818
819
820
821

- Transfer test substance/DMSO solutions to the second set of labeled tubes and add E2 as shown in **Table 12-2**.

822 **Table 12-2 Addition of E2 to Test Substance Serial Dilution for Range Finder Testing**

Tube Number	Test Substance	E2	Estrogen-free DMEM ³	Final Volume
1	Transfer 4 µL of test substance from Tube #1 in Section 12.2 step 1 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
2	Transfer 4 µL of test substance from Tube #2 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
3	Transfer 4 µL of test substance from Tube #3 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
4	Transfer 4 µL of test substance from Tube #4 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
5	Transfer 4 µL of test substance from Tube #5 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
6	Transfer 4 µL of test substance from Tube #6 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL
7	Transfer 4 µL of test substance from Tube #7 to a new tube	Add 4 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 . Vortex to mix.	800 µL	808 µL

823

824 Determination of whether a substance is positive in range finder testing and selection of starting
825 concentrations for comprehensive testing will be discussed in **Section 14.0**.

826 **12.3 Preparation of Reference Standard and Control Dosing Solutions for**
827 **Comprehensive Testing**

828 Comprehensive testing is conducted on 96-well plates using nine concentrations of Ral/E2 in
829 duplicate as the reference standard. Four replicate wells for the DMSO, E2 and flavone/E2
830 controls are included on each plate.

831 All “dosing solutions” of test substance concentrations are to be expressed as µg/mL in the study
832 notebook and in all laboratory reports.

833 Store dosing solutions at room temperature. Use within 24 hours of preparation.

834 12.3.1 Preparation of Ral/E2 Reference Standard Dosing Solutions for Comprehensive
835 Testing

836 In preparation for making Ral/E2 1:2 serial dilutions, label two sets of nine glass 13 mm test
837 tubes with the numbers 1 through 9 and place them in a test tube rack. Tube number 1 will
838 contain the highest concentration of raloxifene (**Table 12-3**).

839 **Table 12-3 Preparation of Ral/E2 Reference Standard Dosing Solution**
840 **for Comprehensive Testing**

Tube Number	Ral/E2 Stock	Estrogen-free DMEM	Final Volume
1	6 µL of Tube #1 from Section 11.3.2	600 µL	606 µL
2	6 µL of Tube #2 from Section 11.3.2	600 µL	606 µL
3	6 µL of Tube #3 from Section 11.3.2	600 µL	606 µL
4	6 µL of Tube #4 from Section 11.3.2	600 µL	606 µL
5	6 µL of Tube #5 from Section 11.3.2	600 µL	606 µL
6	6 µL of Tube #6 from Section 11.3.2	600 µL	606 µL
7	6 µL of Tube #7 from Section 11.3.2	600 µL	606 µL
8	6 µL of Tube #8 from Section 11.3.2	600 µL	606 µL
9	6 µL of Tube #9 from Section 11.3.2	600 µL	606 µL

841

842 12.3.2 Preparation of DMSO Control Comprehensive Testing Dosing Solution

- 843 1. Add 10 µL of 100% DMSO to a 13 mm glass test tube.
- 844 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex to mix.

845 12.3.3 Preparation of E2 Control Comprehensive Testing Dosing Solution

- 846 1. Add 5 µL of the E2 stock from **Section 11.1.1** to a 13 mm glass test tube.
- 847 2. Add 5 µL of 100% DMSO to the tube.
- 848 3. Add 1000 µL of estrogen-free DMEM to the tube and vortex to mix.

849 12.3.4 Preparation of Flavone/E2 Control Comprehensive Dosing Solution

- 850 1. Add 10 µL of flavone/E2 from **Section 11.4** to a 13 mm glass test tube.
- 851 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex to mix.

852 **12.4 Preparation of Test Substance Dosing Solutions for Comprehensive Testing**

853 Comprehensive testing experiments are used to determine whether a substance possesses ER
854 antagonist activity in the LUMI-CELL® ER test method. Antagonist comprehensive testing for
855 coded substances consists of either an 11-point 1:2 serial dilution, or an 11-point 1:5 serial
856 dilution with each concentration tested in triplicate wells of the 96-well plate.

857 12.4.1 *Preparation of Test Substance 1:2 Serial Dilutions for*
858 *Comprehensive Testing*

859 Start the 11-point serial dilution according to criteria in **Section 14.0**.

860 To make test substance 1:2 serial dilutions for comprehensive testing:

- 861 1. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
862 tube rack
- 863 2. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
864 tube rack and add 800 µL of estrogen-free DMEM to each tube

865 Prepare dilution of test substance as shown in **Table 12-4**.

866

866 **Table 12-4 Preparation of Test Substance 1:2 Serial Dilutions for Comprehensive Testing**

Tube Number	100% DMSO	Test Substance ¹	Discard	E2 Testing Stock	Estrogen-free DMEM ²	Final Volume
1	-	4 µL of test substance solution from Section 10.2.4 step 1	-	4 µL	800 µL	808 µL
2	4 µL	4 µL of test substance solution from Section 10.2.4 step 1	-	4 µL	800 µL	808 µL
3	4 µL	4 µL from Tube #2	-	4 µL	800 µL	808 µL
4	4 µL	4 µL from Tube #3	-	4 µL	800 µL	808 µL
5	4 µL	4 µL from Tube #4	-	4 µL	800 µL	808 µL
6	4 µL	4 µL from Tube #5	-	4 µL	800 µL	808 µL
7	4 µL	4 µL from Tube #6	-	4 µL	800 µL	808 µL
8	4 µL	4 µL from Tube #7	-	4 µL	800 µL	808 µL
9	4 µL	4 µL from Tube #8	-	4 µL	800 µL	808 µL
10	4 µL	4 µL from Tube #9	-	4 µL	800 µL	808 µL
11	4 µL	4 µL from Tube #10	4 µL	4 µL	800 µL	808 µL

867 ¹Vortex tubes #2 through 10 before removing test substance/DMSO solution to place in the next tube in the series.868 ²Vortex all tubes to mix media, test substance, and E2.

869

870 12.4.2 *Preparation of Test Substance 1:5 Serial Dilutions for*
871 *Comprehensive Testing*872 Start the 11-point serial dilution according to criteria in **Section 14.0**.

873 To make test substance 1:5 serial dilutions for comprehensive testing:

874 1. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
875 tube rack876 2. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
877 tube rack and add 800 µL of estrogen-free DMEM to each tube878 Prepare dilution of test substance as shown in **Table 12-5**.879 **Table 12-5 Preparation of Test Substance 1:5 Dilutions for Comprehensive Testing**

Tube Number	100% DMSO	Test Substance ¹	Discard	E2 Testing Stock	Estrogen-free DMEM ²	Final Volume
1	-	4 µL of test substance solution	-	4 µL	800 µL	808 µL

Tube Number	100% DMSO	Test Substance ¹	Discard	E2 Testing Stock	Estrogen-free DMEM ²	Final Volume
		from Section 10.2.4 step 1				
2	16 µL	4 µL of test substance solution from Section 10.2.4 step 1	-	4 µL	800 µL	808 µL
3	16 µL	4 µL from Tube #2	-	4 µL	800 µL	808 µL
4	16 µL	4 µL from Tube #3	-	4 µL	800 µL	808 µL
5	16 µL	4 µL from Tube #4	-	4 µL	800 µL	808 µL
6	16 µL	4 µL from Tube #5	-	4 µL	800 µL	808 µL
7	16 µL	4 µL from Tube #6	-	4 µL	800 µL	808 µL
8	16 µL	4 µL from Tube #7	-	4 µL	800 µL	808 µL
9	16 µL	4 µL from Tube #8	-	4 µL	800 µL	808 µL
10	16 µL	4 µL from Tube #9	-	4 µL	800 µL	808 µL
11	16 µL	4 µL from Tube #10	20 µL	4 µL	800 µL	808 µL

¹Vortex tubes #2 through 10 before removing test substance/DMSO solution to place in the next tube in the series.

²Vortex all tubes to mix media, test substance, and E2.

880
881
882

883 13.0 GENERAL PROCEDURES FOR THE TESTING OF CODED SUBSTANCES

884 Range finder experiments are used to determine the concentrations of test substance to be used
885 during comprehensive testing. Comprehensive testing experiments are used to determine whether
886 a substance possesses ER antagonist activity in the LUMI-CELL® ER test method.

887 General procedures for range finder and comprehensive testing are nearly identical. For specific
888 details (such as plate layout) of range finder testing see Section 14.0. For specific details of
889 comprehensive testing, see Section 15.0.

890 13.1 Application of Reference Standard, Control and Test Substances

- 891 1. Remove the 96-well plates (from Section 9.2.3 step 18) from the incubator;
892 inspect them using an inverted microscope. Only use plates in which the cells in
893 all wells receive a score of 1 according to Table 11-1.
- 894 2. Remove medium by inverting the plate onto blotter paper. Gently tap plate against
895 the bench surface to remove residual liquid trapped in the wells.

896 3. Add 200 µL of medium, reference standard, control or test substance to each well
897 (see **Sections 14.0** and **15.0** for specific plate layouts).

898 4. Return plates to incubator (see **Section 9.0** for details) for 19 to 24 hours to allow
899 maximal induction of luciferase activity in the cells.

900 13.1.1 Preparation of Excel® Data Analysis Template For Range Finder Testing

901 1. In Excel®, open a new “AntRFTemplate” and save it with the appropriate project
902 name as indicated in the NICEATM Style Guide.

903 2. Fill out the table at the top of the “Raw Data” worksheet with information
904 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
905 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
906 default template “AntRFTemplate” on a laboratory specific basis).

907 3. Add the following information regarding the assay to the “Compound Tracking”
908 worksheet.

909 ▪ Plate # - Enter the experiment ID or plate number into cell E1

910 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
911 experiment into cell B5

912 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
913 Media in cells B6 and B7

914 ▪ Test Substance Code – Enter the test substance codes into cells C14 to
915 C19

916 ▪ Name: Enter the experimenter name into cell G6

917 ▪ Date: Enter the experiment date in the format day\month\year into cell
918 G10

919 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
920 plate contaminated)

921 4. Enter the following substance testing information to the “List” worksheet:

- 922 ▪ Concentration – Type in the test substance concentration in µg/ml in
923 descending order.
- 924 ▪ Any specific comments about the test substance or condition of the wells
925 should be entered into this sheet, in the comments section
- 926 ▪ All of the remaining cells on the “List” worksheet should populate
927 automatically.
- 928 ▪ The “Template”, “Compound Mixing” and “Visual Inspection”
929 worksheet should automatically populate with the information entered
930 into the “Compound Tracking” and “List” worksheet.
- 931 5. Save the newly named project file.
- 932 6. Print out either the “List” or “Template” worksheet for help with dosing the 96-
933 well plate. Sign and date the print out and store in study notebook.
- 934 13.1.2 Preparation of Excel® Data Analysis Template for Comprehensive Testing
- 935 1. In Excel®, open a new “AntCTTemplate” and save it with the appropriate project
936 name as indicated in the NICEATM Style Guide.
- 937 2. Fill out the table at the top of the “Raw Data” worksheet with information
938 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
939 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
940 default template “AntCTTemplate” on a laboratory specific basis).
- 941 3. On the “Compound Tracking” worksheet, enter the following information:
- 942 ▪ Plate # - Enter the experiment ID or plate number into cell E1
- 943 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
944 experiment into cell C5
- 945 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
946 Media in cells C6 and C7
- 947 ▪ Test Substance Code – Enter the test substance codes into cells C15 and
948 C16. Enter the test substance dilution into cells D15 and D16.

- 949 ▪ Name: Enter the experimenter name into cell F6
- 950 ▪ Date: Enter the experiment date in the format day\month\year into cell
- 951 G10
- 952 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
- 953 plate contaminated)
- 954 4. Enter the following substance testing information to the “List” worksheet:
- 955 ▪ Concentration – Type in the test substance concentration in µg/ml in
- 956 descending order.
- 957 ▪ Any specific comments about the test substance or condition of the wells
- 958 should be entered into this sheet, in the comments section
- 959 ▪ All of the remaining cells on the “List” worksheet should populate
- 960 automatically.
- 961 ▪ The “Template”, “Compound Mixing” and “Visual Inspection”
- 962 worksheet should automatically populate with the information entered
- 963 into the “Compound Tracking” and “List” worksheet.
- 964 5. Save the newly named project file.
- 965 6. Print out either the “List” or “Template” worksheet for help with dosing the 96-
- 966 well plate. Sign and date the print out and store in study notebook.

967 **13.2 Visual Evaluation of Cell Viability**

- 968 1. 19 to 24 hours after dosing the plate, remove the plate from the incubator and
- 969 remove the media from the wells by inverting the plate onto blotter paper. Gently
- 970 tap plate against the bench surface to remove residual liquid trapped in the wells.
- 971 2. Use a repeat pipetter to add 50 µL 1X PBS to all wells. Immediately remove PBS
- 972 by inversion.
- 973 3. Using an inverted microscope, inspect all of the wells used in the 96-well plate
- 974 and record the visual observations using the scores in **Table 13-1**.
- 975

975 **Table 13-1 Visual Observation Scoring**

Viability Score	Brief Description ¹
1	Normal Cell Morphology and Cell Density
2	Altered Cell Morphology and/or Small Gaps between Cells
3	Altered Cell Morphology and/or Large Gaps between Cells
4	Few (or no) Visible Cells
P	Wells containing precipitation are to be noted with “P”

976 ¹Reference photomicrographs are provided in the LUMI-CELL® ER Validation Study “Visual Observation Cell
977 Viability Manual.”

978

979 **13.3 Lysis of Cells for LUMI-CELL® ER**

- 980 1. Apply the reflective white backing tape to the bottom of the 96-well plate (this
981 will increase the effectiveness of the luminometer).
- 982 2. Add 30µL 1X lysis reagent to the assay wells and place the 96-well plate on an
983 orbital shaker for one minute.
- 984 3. Remove plate from shaker and measure luminescence (as described in **Section**
985 **13.4**).

986 **13.4 Measurement of Luminescence**

987 Luminescence is measured in the range of 300 to 650 nm, using an injecting luminometer and
988 with software that controls the injection volume and measurement interval. Light emission from
989 each well is expressed as relative light units (RLU) per well. The luminometer output is saved as
990 raw data in an Excel® spreadsheet. A hard copy of the luminometer raw data should be signed,
991 dated and stored in the study notebook.

992 **13.5 Data Analysis**

993 LUMI-CELL® ER uses an Excel® spreadsheet to collect and adjust the RLU values obtained
994 from the luminometer and a GraphPad Prism® template to analyze and graph data. Plate
995 reduction is calculated using unadjusted RLU values.

996 The Excel® spreadsheet subtracts background luminescence (average DMSO solvent control
997 RLU value) from test substance, reference standard and control RLU values. Test substance,
998 reference standard, and control RLU values are then adjusted relative to the highest Ral/E2

999 reference standard RLU value, which is set to 10,000. After adjustment, values are transferred to
1000 GraphPad Prism® for data analysis and graphing.

1001 13.5.1 Collection and Adjustment of Luminometer Data for Range Finder Testing

1002 The following steps describe the procedures required to populate the Excel® spreadsheet that has
1003 been configured to collect and adjust the RLU values obtained from the luminometer.

- 1004 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
1005 from **Section 13.1.1**.
- 1006 2. Copy the raw data using the Excel® copy function, then paste the copied data into
1007 cell B19 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
1008 the **Paste Special – Values** command. This position corresponds to position A1 in
1009 the table labeled Table 1 in this tab.
- 1010 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
1011 whether there are any potential outliers. See **Section 13.5.3** for further explanation
1012 of outlier determinations.
- 1013 4. If an outlier is identified, perform the following steps to remove the outlier from
1014 calculations:
 - 1015 ▪ correct the equation used to calculate DMSO background in Table 1
1016 [e.g., if outlier is located in cell F26, adjust the calculation in cell H40 to
1017 read =AVERAGE(E26,G26)]
 - 1018 ▪ then correct the equation used to calculate the average DMSO value in
1019 Table 2 [e.g., following the above example, adjust cell M42 to read
1020 =AVERAGE(E38,G38)]
 - 1021 ▪ then correct the equation used to calculate the standard deviation of the
1022 DMSO value in Table 2 [e.g., following the above example, adjust cell
1023 M43 to read =STDEV(E38,G38)]
- 1024 5. Excel® will automatically subtract the background (the average DMSO control
1025 value) from all of the RLU values in Table 1 and populate Table 2 with these
1026 adjusted values.

- 1027 6. To calculate plate reduction, identify the cell containing the Ral/E2a replicate in
1028 Table 1, plate row H that has the lowest RLU value (i.e., cell B26, C26, or D26).
- 1029 7. Identify the cell containing the Ral/E2a replicate in Table 1, plate row H that has
1030 the highest RLU value (i.e., cell B26, C26, or D26).
- 1031 8. Click into cell D14 and enter the cell number from **Section 13.5.1 step 7** into the
1032 numerator and the cell number from **step 6** into the denominator.
- 1033 9. Identify the cell containing the Ral/E2b replicate in Table 1, plate row H that has
1034 the lowest RLU value (i.e., cell K26, L26, or M26).
- 1035 10. Identify the cell containing the Ral/E2b replicate in Table 1, plate row H that has
1036 the highest RLU value (i.e., cell K26, L26, or M26).
- 1037 11. Click into cell E14 and enter the cell number from **Section 13.5.1 step 10** into the
1038 numerator and the cell number from **step 9** into the denominator.
- 1039 12. Click on the “ER Antagonist Report” worksheet.
- 1040 13. The data for the Ral/E2 reference standard, DMSO, and E2, replicates populate
1041 the left portion (columns A-F) of the spreadsheet. The data is automatically
1042 placed into an Excel® graph.
- 1043 14. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
1044 C2 of “ER Antagonist Report” worksheet and check the formula contained within
1045 that cell. The divisor should be the cell number of the cell containing the highest
1046 averaged Ral/E2 RLU value (column A).
- 1047 15. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
1048 scores for each well on the 96-well plate. This data will be linked to the “ER
1049 Antagonist Report” worksheet.
- 1050 16. After the testing results have been evaluated and reviewed for quality control,
1051 enter the following information into the Compound Tracking worksheet:
- 1052 ▪ Enter pass/fail results for plate reference standard and control parameters
1053 into the Plate Pass/Fail Table

- 1054 ▪ Enter information from the testing of coded substances into the Testing
1055 Results Table
- 1056 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
1057 data into cell A34
- 1058 ▪ Date – Enter the date on which the data was reviewed into cell D34

1059

1060 13.5.2 Collection and Adjustment of Luminometer Data for Comprehensive Testing

1061 The following steps describe the procedures required to populate the Excel® spreadsheet that has
1062 been configured to collect and adjust the RLU values obtained from the luminometer.

- 1063 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
1064 from **Section 13.1.2**.
- 1065 2. Copy the raw data using the Excel® copy function, then paste the copied data into
1066 cell B14 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
1067 the **Paste Special – Values** command. This position corresponds to position A1 in
1068 the table labeled Table 1 in this tab.
- 1069 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
1070 whether there are any potential outliers. See **Section 13.5.3** for further explanation
1071 of outlier determinations.
- 1072 4. If an outlier is identified, perform the following steps to remove the outlier from
1073 calculations:
- 1074 ▪ correct the equation used to calculate DMSO background in Table 1 [e.g.,
1075 if outlier is located in cell M14, adjust the calculation in cell H40 to read
1076 =AVERAGE(M15:M17)]
- 1077 ▪ then correct the equation used to calculate the average DMSO value in
1078 Table 2 [e.g., following the above example, adjust cell M35 to read
1079 =AVERAGE(M25:M27)]

- 1080 ▪ then correct the equation used to calculate the standard deviation of the
1081 DMSO value in Table 2 [e.g., following the above example, adjust cell
1082 M36 to read =STDEV(M25:M27)]
- 1083 5. Excel® will automatically subtract the background (the average DMSO control
1084 value) from all of the RLU values in Table 1 and populate Table 2 with these
1085 adjusted values.
- 1086 6. To calculate plate reduction, identify the cell containing the Ral/E2 replicate in
1087 plate row G that has the lowest RLU value.
- 1088 7. Identify the cell containing the Ral/E2 replicate in plate row G that has the highest
1089 RLU value.
- 1090 8. Click into cell D14 and enter the cell number from **Section 13.5.2 step 7** into the
1091 numerator and the cell number from **step 6** into the denominator.
- 1092 9. Identify the cell containing the Ral/E2 replicate in plate row H that has the lowest
1093 RLU value.
- 1094 10. Identify the cell containing the Ral/E2 replicate in plate row H that has the highest
1095 RLU value.
- 1096 11. Click into cell E14 and enter the cell number from **Section 13.5.2 step 10** into the
1097 numerator and the cell number from **step 9** into the denominator.
- 1098 12. Click on the “ER Antagonist Report” worksheet.
- 1099 13. The data for the Ral/E2 reference standard, DMSO, E2, and Flavone/E2 replicates
1100 populate the left portion (columns A-E) of the spreadsheet. The data is
1101 automatically placed into an Excel® graph.
- 1102 14. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
1103 D2 of “ER Antagonist Report” worksheet and check the formula contained within
1104 that cell. The divisor should be the cell number of the cell containing the highest
1105 averaged Ral/E2 RLU value (column A).

- 1106 15. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
1107 scores for each well on the 96-well plate. This data will be linked to the “ER
1108 Antagonist Report” worksheet.
- 1109 16. Copy the data into GraphPad Prism® for the calculation of IC₅₀ values and to
1110 graph experimental results as indicated in the NICEATM Prism® Users Guide.
- 1111 17. After the testing results have been evaluated and reviewed for quality control,
1112 enter the following information into the Compound Tracking worksheet:
- 1113 ▪ Enter pass/fail results for plate reference standard and control parameters
1114 into the Plate Pass/Fail Table
 - 1115 ▪ Enter information from the testing of coded substances into the Testing
1116 Results Table
 - 1117 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
1118 data into cell A34
 - 1119 ▪ Date – Enter the date on which the data was reviewed into cell D32

1120 13.5.3 Determination of Outliers

1121 The Study Director will use good statistical judgment for determining “unusable” wells that will
1122 be excluded from the data analysis and will provide an explanation in the study notebook for any
1123 excluded data. This judgment for data acceptance will include Q-test analysis.

1124 The formula for the Q test is:

$$1125 \frac{\text{Outlier} - \text{Nearest Neighbor}}{\text{Range (Highest} - \text{Lowest)}}$$

1126 where the outlier is the value proposed for exclusion, the nearest neighbor is the value closest to
1127 the outlier, and the range is the range of the three values (Q values for samples sizes from 3 to 10
1128 are provided in Table 13-2). For example, if the value of this ratio is greater than 0.94 (the Q
1129 value for the 90% confidence interval for a sample size of three) or 0.76 (the Q value for the 90%
1130 confidence interval for a sample size of four), the outlier may be excluded from data analysis.

1131

1131 **Table 13-2 Q Test Values**

Number Of Observations	Q Value
2	-
3	0.94
4	0.76
5	0.64
6	0.56
7	0.51
8	0.47
9	0.44
10	0.41

1132

1133 For E2 reference standard replicates (sample size of two), any adjusted RLU value for a replicate
1134 at a given concentration of E2 is considered and outlier if its value is more than 20% above or
1135 below the adjusted RLU value for that concentration in the historical database.

1136 13.5.4 Acceptance Criteria1137 13.5.4.1 *Range Finder Testing*

1138 Acceptance or rejection of a range finder test is based on reference standard and solvent control
1139 results from each experiment conducted on a 96-well plate.

- 1140 • Reduction: Plate reduction, as measured by dividing the averaged highest Ra/E2
1141 reference standard RLU value by the averaged DMSO control RLU value, must
1142 be greater than three-fold.
- 1143 • DMSO control results: DMSO control RLU values must be within 2.5 times the
1144 standard deviation of the historical solvent control mean RLU value (see **Section**
1145 **16.5**).

1146 An experiment that fails either acceptance criterion will be discarded and repeated.

1147 13.5.4.2 *Comprehensive Testing*

1148 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1149 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1150 (QC) for these parameters derived from the historical database (see **Section 16.5**), which are
1151 summarized below.

- 1152 • Reduction: Plate reduction, as measured by dividing the averaged highest Ral/E2
1153 reference standard RLU value by the averaged lowest Ral/E2 control RLU value,
1154 must be greater than three-fold.
- 1155 • DMSO control results: DMSO control RLU values must be within 2.5 times the
1156 standard deviation of the historical solvent control mean RLU value (see **Section**
1157 **16.5**).
- 1158 • Reference standard results: The Ral/E2 reference standard concentration-response
1159 curve should be sigmoidal in shape and have at least three values within the linear
1160 portion of the concentration-response curve.
- 1161 • E2 control results: E2 control RLU values must be within 2.5 times the standard
1162 deviation of the historical E2 control mean RLU value.
- 1163 • Positive control results: Flavone/E2 control RLU values must be less than the E2
1164 control mean minus three times the standard deviation from the E2 control mean.


1165 An experiment that fails any single acceptance criterion will be discarded and repeated.










1166 **14.0 RANGE FINDER TESTING**

1167 Antagonist range finding for coded substances consists of seven point, 1:10 serial dilutions tested
1168 in duplicate wells of the 96-well plate. **Figure 14-1** contains a template for the plate layout used
1169 in antagonist range finder testing.

1170

1170 **Figure 14-1 Antagonist Range Finder Plate Layout**



-  **Three Point Ral/E2 Reference Standard**
-  **DMSO (Solvent Control)**
-  **Range Finder for Sample #1**
-  **Range Finder for Sample #2**
-  **Range Finder for Sample #3**
-  **Range Finder for Sample #4**
-  **Range Finder for Sample #5**
-  **Range Finder for Sample #6**
-  **E2 Control**

1171

1172


1172 Evaluate whether range finder experiments have met acceptance criteria (see **Section 13.6.3**).

1173 To determine starting concentrations for comprehensive testing use the following criteria:

- 1174
- 1175
- 1176
- 1177
- 1178
- 1179
- If results in the range finder test suggest that the test substance is negative for antagonist activity (i.e., if there are no points on the test substance concentration curve that are less than the mean minus three times the standard deviation of the E2 control, see **Figure 14-2**), comprehensive testing will be conducted using an 11-point 1:2 serial dilution with the limit dose as the starting concentration (i.e., the maximum soluble concentration in the range finder).
 - If results in the range finder test suggest that the test substance is positive for antagonist activity (i.e., if there are points on the test substance concentration curve that are less than the mean minus three times the standard deviation of the E2 control), the top concentration to be used for the 11-point dilution scheme in comprehensive testing should be one log concentration higher than the concentration giving the lowest adjusted RLU value in the range finder or the limit dose. The 11-point dilution scheme will be based on either a 1:2 or 1:5 serial or dilution according to the following criteria:
 - An 11-point 1:2 serial dilution should be used if the resulting concentration range (note: an 11-point 1:2 serial dilution will cover a range of concentrations over approximately three orders of magnitude [three logs]) will encompass the full range of responses based on the concentration response curve generated in the range finder test (see **Figure 14-3**).
 - If the concentration range that would be generated with the 1:2 serial dilution will not encompass the full range of responses based on the concentration response curve in the range finder test (see **Figure 14-4**), an 11-point 1:5 serial dilution should be used instead.
 - If a substance exhibits a biphasic concentration response curve in the range finder test (see **Figure 14-5**), both phases should also be resolved in comprehensive testing. In this case, two peaks could potentially be used to identify the top
- 1198
- 1199
- 1200


1201 concentration to be used for the 11-point dilution scheme in comprehensive
1202 testing. In order to resolve both curves, the top concentration should be based on
1203 the peak associated with the higher concentration and the top dose one log
1204 concentration higher than the concentration giving the lowest adjusted RLU value
1205 in the range finder. An 11-point 1:5 serial dilution should be used.

1206 **Figure 14-2 Antagonist Range Finder (example 1)**



1207
1208 The solid horizontal line represents the mean minus three times the standard deviation of the E2 control.
1209


1209 **Figure 14-3 Antagonist Range Finder (example 2)**



1210

1211 The solid horizontal line represents the mean minus three times the standard deviation of the E2 control.


1212 **Figure 14-4 Antagonist Range Finder (example 3)**



1213

1214 The solid horizontal line represents the mean minus three times the standard deviation of the E2 control.

1215 **Figure 14-5 Antagonist Range Finder (example 4)**



1216

1217 The solid horizontal line represents the mean minus three times the standard deviation of the E2 control.

1218

1219 **15.0 COMPREHENSIVE TESTING**


1220 Antagonist comprehensive testing for coded substances consists of 11 point, 1:2 serial dilutions,
 1221 with each concentration tested in triplicate wells of the 96-well plate. **Figure 15-1** contains a
 1222 template for the plate layout to be used in antagonist comprehensive testing.

1223

1224

1224 **Figure 15-1 Antagonist Comprehensive Test Plate Layout**

1225
1226
1227
1228
1229
1230
1231
1232
1233
1234
1235
1236
1237
1238
1239
1240
1241
1242
1243
1244
1245
1246
1247
1248
1249
1250
1251
1252
1253
1254



- 9 Point Duplicate Ral/E2 Reference Standard
- DMSO (Solvent Control)
- Test Substance #1
- Test Substance #2
- E2 Control
- Flavone Control

1254 Evaluate whether comprehensive experiments have met acceptance criteria (see **Section 13.6.3**)
1255 and graph the data as described in the NICEATM Prism® users guide.

1256 • If the substance has been tested up to the limit dose or the maximum soluble dose
1257 without causing a significant decrease in cell viability, and there are no points on
1258 the concentration curve that are less than the mean minus three times the standard
1259 deviation of the E2 control, the substance is considered negative for antagonism.

1260 • If the substance has been tested up to the limit dose and there are points on the
1261 concentration curve that are less than the mean minus three times the standard
1262 deviation of the E2 control, but cell viability has a visual inspection score of 2 or
1263 greater, at all points falling below the E2 line, the substance is considered
1264 negative for antagonism.

1265 • If there are points on the test substance concentration curve that are less than the
1266 mean minus three times the standard deviation of the E2 control that do not cause
1267 a visual inspection score of 2 or greater, the substance is positive for antagonism.

1268 – Points in the test substance concentration curve that cause a visual
1269 inspection score of 2 or greater, are not included in data analyses.

1270 **16.0 COMPILATION OF THE HISTORICAL QUALITY CONTROL DATABASE**

1271 Historical databases are maintained in order to ensure that the assay is functioning properly.

1272 Historical databases are compiled using Excel® spreadsheets and are separate from the
1273 spreadsheets used to collect the data for individual test plates. Reference standard and control
1274 data is used to develop and maintain the historical database and are used as quality controls to
1275 determine acceptance of individual test plates.

1276 The sources of data needed to compile the historical database for the E2 control and flavone/E2
1277 control values are the experiment specific Excel® data collection and analysis spreadsheets (see
1278 **Section 13.5.2**) used for LUMI-CELL® ER antagonist testing. The sources of the data needed to
1279 compile the historical database for the DMSO control are the experiment specific Excel® data
1280 collection and analysis spreadsheets used for LUMI-CELL® ER antagonist and agonist testing

1281 (see **Section 13.5.2** of the LUMI-CELL® ER antagonist protocol and **Section 11.5.2** in the
1282 LUMI-CELL® ER agonist protocol).

1283 **16.1 E2 Control**

1284 Open the LUMI-CELL® ER antagonist specific historical database Excel® spreadsheet
1285 (LUMI_AgandAntQC.xls) and save under a new name using the Excel® “Save As” function,
1286 adding the laboratory designator to the file name (e.g., for Laboratory H, the new name would be
1287 HLUMI_AgandAntQC.xls). Open the E2 Control worksheet and enter the date and experiment
1288 name into worksheet columns A and B respectively. Enter the experimental mean adjusted E2
1289 control value (from cell D37 in the ER Antagonist Report worksheet of the Excel® data
1290 collection and analysis spreadsheet) into the Antagonist E2 control worksheet, column C.
1291 Acceptance or rejection of plate E2 control data for comprehensive testing is based on whether
1292 the mean plate E2 RLU value falls within 2.5 times the standard deviation of the E2 value in the
1293 historical database (columns G and H in the E2 Control worksheet).

1294 **16.2 DMSO**

1295 Open the combined agonist and antagonist LUMI-CELL® ER historical database Excel®
1296 spreadsheet (LUMI_AgandAntQC.xls) and save under a new name using the Excel® “Save As”
1297 function, adding the laboratory designator to the file name (e.g., for Laboratory H, the new name
1298 would be HLUMI_AgandAntQC.xls). Enter the date and experiment name into worksheet
1299 columns A and B respectively. Enter the experimental mean DMSO control value (from cell H37
1300 in the RAW DATA worksheet of the agonist and antagonist Excel® data collection and analysis
1301 spreadsheet) into worksheet column C. Acceptance or rejection of the plate DMSO control data
1302 for range finding and comprehensive testing is based on whether the mean plate DMSO RLU
1303 value falls within 2.5 times the standard deviation of the DMSO value in the historical database
1304 (columns G and H in the DMSO worksheet).

1305

1306 **17.0 QUALITY TESTING OF MATERIALS**

1307 All information pertaining to the preparation and testing of media, media supplements, and other
1308 materials should be recorded in the Study Notebook.

1309 **17.1 Tissue Culture Media**

1310 Each lot of tissue culture medium must be tested in a single growth flask of cells before use in
1311 ongoing tissue culture or experimentation (**note:** each bottle within a given lot of
1312 Charcoal/Dextran treated FBS must be tested separately).

- 1313 1. Every new lot of media (RPMI and DMEM) and media components (FBS,
1314 Charcoal/Dextran treated FBS, and L-glutamine) must first be tested on the
1315 LUMI-CELL® ER assay prior to being used in any GLP acceptable assays.
- 1316 2. Add 4 µL of DMSO (previously tested) into four separate 13 mm tubes.
- 1317 3. Add 400 µL media (to be tested) to 13 mm tube.
- 1318 4. Dose an experimental plate as in **Section 12.0**, treating the media being tested as a
1319 test substance.
- 1320 5. Analyze 96-well plate as described in **Section 12.0**, comparing the data from the
1321 DMSO controls made using previously tested tissue culture media to the new
1322 media being tested.
- 1323 6. Use the agonist historical database to determine if the new media with DMSO lies
1324 within 2.5 standard deviations of the mean for the media. If the RLU values for
1325 the new media with DMSO lie within 2.5 standard deviations of the DMSO mean
1326 from the historical database, the new lot of media is acceptable. If the RLU values
1327 for the new media with DMSO do not lie within 2.5 standard deviations of the
1328 DMSO mean from the historical database, the new lot may not be used in the
1329 assay.
- 1330 7. Note date and lot number in study notebook.
- 1331 8. If the new bottle passes quality testing as described in **Section 15.1 step 6**, apply
1332 the media to a single flask cells and observe the cells growth and morphology
1333 over the following 2 to 3 days. If there is no change in growth or morphology, the
1334 new media is acceptable for use.

1335

1335 **17.2 G418**

- 1336 1. New lots of G418 must first be tested on the LUMI-CELL® ER assay prior to
1337 being used in any GLP acceptable assays.
- 1338 2. Add 220 µL of G418 (previously tested) to a single flask containing cells growing
1339 in RPMI.
- 1340 3. Add 220 µL of G418 (to be tested) to a different flask containing cells growing in
1341 RPMI.
- 1342 4. Observe cellular growth and morphology in both tissue culture flasks over a 48 to
1343 72 hour period. If there are no differences in observed growth rate and
1344 morphology between the two flasks, the new G418 lot is acceptable.
- 1345 5. If cellular growth is decreased, or the cells exhibit abnormal morphology, the new
1346 lot of G418 is not acceptable.
- 1347 6. Note date and lot number in study book.

1348 **17.3 DMSO**

- 1349 1. Every new bottle of DMSO must be tested on the LUMI-CELL® ER assay prior
1350 to use in any GLP acceptable assays.
- 1351 2. Add 4 µL of DMSO (to be tested) into four separate 13 mm tubes.
- 1352 3. Add 400 µL media (previously tested) the same tubes.
- 1353 4. Dose an experimental plate as in **Section 15.0**, treating the media being tested as a
1354 test substance.
- 1355 5. Analyze 96-well plate as described in **Section 15.0**, comparing the data from the
1356 DMSO controls made using previously tested tissue culture media to the new
1357 media being tested.
- 1358 6. Use the agonist historical database to determine if media with new DMSO lies
1359 within 2.5 standard deviations of the DMSO mean from historical database. If the
1360 RLU values for the media with new DMSO lie within 2.5 standard deviations of
1361 the DMSO mean from the historical database, the new lot of DMSO is acceptable.

1362 If the RLU values for media with new DMSO do not lie within 2.5 standard
1363 deviations of the DMSO mean from historical database, the new lot may not be
1364 used in the assay.

1365 7. Note the date, lot number, and bottle number in study book.

1366 8. If no DMSO has been previously tested, test several bottles as described in
1367 **Section 15.3**, and determine whether any of the bottles of DMSO have a higher
1368 average RLU than the other bottle(s) tested. Use the DMSO with the lowest
1369 average RLU for official experiments.

1370 **17.4 Plastic Tissue Culture Materials**

1371 1. Grow one set of cells, plate them for experiments on plastic ware from the new lot
1372 and one set of cells in the plastic ware from a previous lot, and dose them with E2
1373 reference standard and controls.

1374 2. Perform the LUMI-CELL® ER experiment with both sets of cells.

1375 3. If all of the analysis falls within acceptable QC criteria, then the new
1376 manufacturer's products may be used.

1377

1377 **18.0 REFERENCES**

- 1378 Eli Lilly and Company and National Institutes of Health Chemical Genomics Center. 2005.
1379 Assay Guidance Manual Version 4.1. Bethesda, MD: National Institutes of Health. Available:
1380 http://www.ncgc.nih.gov/guidance/manual_toc.html [accessed 05 September 2006]
- 1381 ICCVAM. 2001. Guidance Document on Using *In Vitro* Data to Estimate *In Vivo* Starting Doses
1382 for Acute Toxicity. NIH Pub. No. 01-4500. Research Triangle Park, NC: National Institute of
1383 Environmental Health Sciences. Available: [http://iccvam.niehs.nih.gov/methods/invidocs/
1384 guidance/iv_guide.pdf](http://iccvam.niehs.nih.gov/methods/invidocs/guidance/iv_guide.pdf) [accessed 31 August 2006]