




Paul Langan
 (505) 665-8125
 langan_paul@lanl.gov

Benno Schoenborn
 (505) 665-2033
 schoenborn@lanl.gov



Protein Crystallography Station (PCS)

The Bioscience Division has built a Protein Crystallography Station (PCS) at the Lujan Neutron Scattering Center for the international structural biology community to investigate the structure and dynamics of proteins, biological polymers, and membranes. Emerging efforts in structural genomics, to solve the three-dimensional structures of thousands of proteins sequenced in genome projects, largely use synchrotron x-rays. However, certain unique types of information relating to structurally and functionally important water molecules and hydrogen atoms can only be obtained using neutrons. Despite the importance of neutrons, the number of studies are limited by access to neutron sources and appropriate instrumentation. In particular, the neutron flux and detector technologies on presently available


instruments have led to prohibitively large sample size requirements. Large crystals of proteins are difficult to grow. The PCS, funded by the Department of Energy Office of Biological and Environmental Research, is the only resource of its kind in North America, and the first to be built at a spallation neutron source. A number of technological innovations have been incorporated, including a partially coupled water moderator and a large cylindrical detector that, along with the improved LANSCE neutron source, make data collection from smaller crystals feasible. These innovations will greatly broaden the application of neutrons to structural biology and make the PCS a benchmark for similar projects worldwide. Approximately 20 experiments are accommodated each year during an 8 month run cycle.




The PCS cylindrical neutron detector consists of an electrode structure contained in an aluminum pressure vessel. The pressure vessel is filled with a mixture of ^3He and propane. The ^3He has an extremely high cross-section for thermal neutron absorption. Neutrons diffracted by the sample are absorbed by the ^3He . This interaction results in the creation of a proton and triton. These primary ionization products drift toward the nearest electrode-anode wire where they multiply in the high electric field near the wire surface. The charges induced over different electrodes allow for the spatial detection of an event.



Neutron diffraction from fibers of biological polymers such as DNA and polysaccharides reveals details of hydrogen bonding patterns and hydration.



Details of the 2.2-Å neutron map clearly reveal D positions.



Details from the structure of myoglobin (a) Scattering density calculated from 1.8-Å x-ray data. Only C atoms are covered. (b) Negative (blue) and positive (red) scattering density calculated from 2-Å neutron data. H with its strong negative scattering

length is covered. (c) Positive scattering density calculated from predeuterated myoglobin. Deuterium (D) with its strong positive scattering density is covered. Neutron diffraction can be used to locate H or D with data extending to less than 2 Å.

Specifications	
Unit cell size	150 Å
Sample size	1 mm
d-spacing	1 Å–250 Å
Wavelength shaping	T ₀ and T ₁ choppers and tail cutter
Wavelength range	1 Å–5 Å
Sample-to-detector distance	700 mm
Beam size	5 mm on detector
Resolution	1.2 mm
Detector area	2930 cm ²
Counting rate	> 1,000,000 c/s