



Respirator Fitting Instructions

1. Hold the respirator in hand with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.
2. Press the respirator firmly against your face with the nosepiece on the bridge of your nose.
3. Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.
4. Using both hands, mold the metal nosepiece to the shape of your nose.
- 5.(a) To test fit; a) cup both hands over the respirator being careful not to disturb position, b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.
- 5.(b) To test fit; a) cup both hands over the respirator being careful not to disturb position, b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.
6. Remember, careful observance of these fitting instructions is respirator use.