

Detection of Myosin Heavy Chain in Formalin-Fixed, Paraffin-Embedded Rodent Tissue

Reagents:

[1X Automation Buffer](#)
[3% Hydrogen Peroxide](#)
[Antibody Diluent](#)
[Citrate Buffer](#)
[DAB Chromagen](#)
[Hematoxylin](#)

Antibody Information:

Kit: Vector M.O.M. Kit
Burlingame, CA 94010
www.vectorlabs.com
1-800-227-6666
Catalog #PK-2200

Avidin/Biotin Blocking Kit
Vector Laboratories, Inc.
Burlingame, CA 94010
www.vectorlabs.com
1-800-227-6666
Catalog #SP-2001

Primary antibody: Mouse anti-Cardiac Myosin Heavy Chain
Abcam Inc.
Cambridge, MA 02139
www.abcam.com
1-888-772-2226
Catalog# ab15

Negative control: Normal Mouse Serum
Jackson Immunoresearch Laboratories, Inc.
West Grove, PA 19390
www.jacksonimmuno.com
1-800-367-5296
Catalog #015-000-001

Label antibody: Vector Elite Label

Vector Laboratories, Inc.
Burlingame, CA 94010
www.vectorlabs.com
1-800-227-6666
Catalog #PK-6100

Staining Procedure

- Positive Control Tissue: heart
- Stain Localization: cytoplasm

Deparaffinize and hydrate slides through the following solutions.

| | | |
|----------------------|---------|-----------|
| Xylene | 2 times | 5 minutes |
| 100% EtOH | 2 times | 3 minutes |
| 95% EtOH | 2 times | 3 minutes |
| 1X Automation Buffer | 2 times | 3 minutes |

1. Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.

2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

3. Unmasking Techniques using the decloaker.

Add 500 ml D/W to the pan of the decloaker.

Place full rack of slides in 200 ml of 1X citrate buffer and place in the decloaker.

Decloak for 5 minutes. Pressure_____

Depressurize for 10 minutes.

Remove pan top and cool for 10 min. Temp_____

Rinse in D/W, 2x for 3 min each

4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

5. Incubate sections for 1 HOUR in M.O.M. specific IgG blocking reagent

Made via 2.5 mls 1x PBS plus 2 drops of Mouse IgG blocking reagent

Kit Lot#_____ Exp Date_____ New Kit: yes / no

6. Apply Avidin/Biotin block

Lot#_____ Exp Date_____ New Kit: yes / no

Apply avidin block - 15 min at RT.

Quick rinse in 1X AB.

Apply biotin block - 15 min at RT.

Wipe excess block

DO NOT RINSE SECTIONS WITH BUFFER.

Prepare Vector M.O.M .diluent: 600ul of protein concentrate stock in 7.5 mls of 1X PBS.

Make primary secondary, and label antibody dilution in Vector M.O.M. diluent.

7. Apply myosin primary antibody at 1:40 and incubate for one hour:

Lot# _____ Exp Date _____

For negative control slides, normalize the protein concentration of normal mouse serum to the protein concentration of the primary antibody and use this to make the 1:40 dilution. Apply to slides and Incubate for one hour.

Lot# _____ Reconstituted Date _____

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply M.O.M. biotinylated anti-mouse IgG and incubate for 10 minutes
Made via 10ul of antibody in 2.5mls of Vector M.O.M. diluent.

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply Vector Elite label for 5 minutes. (Prepare 30 minutes before use)
Made via 2 drops of Reagent A plus 2 drops of Reagent B in 2.5mls M.O.M. diluent
Exp Date _____ New Kit: yes / no

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.
(Add 1 drop of DAB per ml of substrate)

Lot # _____ Exp Date _____ New Kit: yes / no

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30 seconds.

16. Rinse in tap water until water is clear.

17. Rinse slides in 1x automation buffer for 1 min.

18. Dehydrate through the following solutions.

| | | |
|-----------|---------|-----------|
| 95% EtOH | 1 times | 3 minutes |
| 100% EtOH | 3 times | 3 minutes |
| Xylene | 2 times | 5 minutes |

19. Coverslip
updated 10/25/04