DENTAL RADIOGRAPHIC SYSTEMS

FIELD TEST RECORD EDIT CHECKS

(Test Procedure DRA - Form FDA 2785)

Verify that:

- 1. The certification status of each component is indicated in item 2.
- 2. For systems manufactured before May 1994, if a value of tube current is not entered in item 32, a value of mAs is entered at item 7 or the system is identified as fixed tube current in Remarks. For systems manufactured after May 1994, data item 32 is the mAs product.
- 3. The appropriate series of aluminum filters has been used in the Beam Quality determination:

Above 70kVp: 4.5 mm, 3.5 mm, 2.5 mm, 1.5 mm

50 - 70 kVp: 3.0 mm, 2.0 mm, 1.5 mm, 1.0 mm

- 4. If only four values are entered for reproducibility, no two exposures differ by more than ten percent of the highest value.
- 5. For controls manufactured after May 1994, item 32 is mAs product and not just mA.