

UNITED STATES OF AMERICA
 DEPARTMENT OF HEALTH AND HUMAN SERVICES
 FOOD AND DRUG ADMINISTRATION

+ + +

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
 MEDICAL DEVICES ADVISORY COMMITTEE

+ + +

CIRCULATORY SYSTEM DEVICES PANEL

+ + +

November 20, 2008
 8:00 a.m.

Holiday Inn
 Gaithersburg, Maryland

PANEL MEMBERS:

JEFFERY BORER, M.D.	Chairperson
VALLUVAN JEEVANANDAM, M.D.	Voting Member
DAVID NAFTEL, Ph.D.	Voting Member
JUDAH WEINBERGER, M.D.	Voting Member
SETH BILAZARIAN, M.D.	Consultant
PAMELA E. KARASIK, M.D.	Consultant
PATRICIA KELLEY, M.D.	Consultant
DAVID SLOTWINER, M.D.	Consultant
JOHN C. SOMBERG, M.D.	Consultant
MICHAEL HALPIN	Industry Representative
MIKE FLEMING, D.D.S., P.A.	Consumer Representative
JAMES P. SWINK	Executive Secretary

Free State Reporting, Inc.
 1378 Cape Saint Claire Road
 Annapolis, MD 21409
 (410) 974-0947

FDA REPRESENTATIVES:

BRIAN ZUCKERMAN, Director
Division of Cardiovascular Devices

SPONSOR PRESENTERS:

MARCIA S. YAROSS, Ph.D.
DONALD A. BERRY, Ph.D.
DAVID J. WILBER, M.D.
ALBERT L. WALDO, M.D.

FDA PRESENTERS:

BENJAMIN ELOFF, Ph.D.
LAURA THOMPSON, Ph.D.
RANDALL BROCKMAN, M.D.
ELLEN PINNOW, M.S.

REVIEWERS:

DAVID SLOTWINTER, M.D.
DAVID NAFTEL, Ph.D.

PRESENT ON BEHALF OF SPONSOR:

MATTHEW REYNOLDS, M.D.
HUGH CALKINS, M.D.

INDEX

	PAGE
CALL TO ORDER - Jeffery Borer, M.D.	5
CONFLICT OF INTEREST - James Swink	6
APPOINTMENT TO TEMPORARY VOTING STATUS - James Swink	9
PANEL INTRODUCTION	11
FIRST OPEN PUBLIC HEARING (No Speakers)	13
SPONSOR PRESENTATION	
Marcia S. Yaross, Ph.D.	15
Donald A. Berry, Ph.D.	24
David J. Wilber, M.D.	35
Albert L. Waldo, M.D.	55
Marcia S. Yaross, Ph.D.	62
SPONSOR Q&A	68
FDA PRESENTATION	
Benjamin Eloff, Ph.D.	97
Laura Thompson, Ph.D.	104
Randall Brockman, M.D.	129
Ellen Pinnow, M.S.	146
FDA Q&A	154
REVIEWER PRESENTATIONS	
David Slotwiner, M.D. - Panel Reviewer	170
David Naftel, Ph.D. - Statistical Reviewer	179
REVIEWER Q&A	188

INDEX

	PAGE
GENERAL PANEL DISCUSSION	222
PANEL DELIBERATIONS AND FDA QUESTIONS	
Question 1	241
Question 2	249
Question 3	261
Question 4	274
Question 5-A	285
Question 5-B	299
Question 5-C	300
Question 5-D	309
Question 5-E	313
Question 6-A	321
SECOND OPEN PUBLIC HEARING (No Speakers)	337
FDA AND SPONSOR SUMMATIONS	338
PANEL VOTE	340
PANEL RECOMMENDATION	359
CLOSING STATEMENTS	
Jeffrey Borer, M.D.	362
Brian Zuckerman, M.D.	362
ADJOURNMENT	

M E E T I N G

(8:02 a.m.)

DR. BORER: I would like to call the
Circulatory System Devices Panel to order.

I am Dr. Jeffery Borer, the Chairperson of
the Panel. I am the Professor-in-Chief of the
Division of Cardiovascular Medicine at State
University of New York, Downstate Medical Center, and
Director of the Cardiovascular Translational Research
Institute as well at that Institution.

If you haven't already done so, please sign
the attendance sheets that are on the tables next to
the doors. If you want to address this Panel during
one of the open sessions, please provide your name to
Ms. AnnMarie Williams at the registration table.

If you're presenting in any of the open
public sessions today and have not previously
provided an electronic copy of your presentation to
the FDA, please arrange to do so with Ms. Williams.

For the record, note that the voting
members present constitute a quorum as required by
C.F.R., Part 14. I'd also like to add that the Panel
participating in the meeting today has received
training in FDA device law and regulations.

No one from the public or press is allowed

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 into the Panel area at any time during the breaks or
2 during the conduct of this meeting.

3 In addition, please remember to put your
4 cell phones on vibrate or silent or something because
5 it would be good if they don't go off during the
6 meeting.

7 Mr. Swink, the Executive Secretary for the
8 Circulatory System Devices Panel will make some
9 introductory remarks. Mr. Swink.

10 MR. SWINK: I will now read the Conflict of
11 Interest Statement.

12 The Food and Drug Administration is
13 convening today's meeting of the Circulatory System
14 Devices Panel of the Medical Devices Advisory
15 Committee under the authority of the Federal Advisory
16 Committee Act of 1972. With the exception of the
17 industry representative, all members and consultants
18 of the Panel are special Government employees or
19 regular Federal employees from other agencies and are
20 subject to Federal conflict of interest laws and
21 regulations.

22 The following information on the status of
23 this Panel's compliance with Federal ethics and
24 conflict of interest laws covered by, but not limited
25 to, those found at 18 U.S.C. Section 208 and Section

1 712 of the Federal Food, Drug and Cosmetic Act are
2 being provided to participants in today's meeting and
3 to the public.

4 FDA has determined that members and
5 consultants of this Panel are in compliance with
6 Federal ethics and conflict of interest laws. Under
7 18 U.S.C. Section 208, Congress has authorized FDA to
8 grant waivers to special Government employees who
9 have potential financial conflicts when it is
10 determined that the Agency's need for that particular
11 individual's services outweighs his or her potential
12 financial conflict of interest. Under Section 712 of
13 the FD&C Act, Congress has authorized FDA to grant
14 waivers to special Government employees and regular
15 Government employees with potential financial
16 conflicts when necessary to afford the committee
17 essential expertise.

18 Related to the discussions of today's
19 meeting, members and consultants of this Panel who
20 are special Government employees have been screened
21 for potential financial conflicts of interest of
22 their own as well as those imputed to them, including
23 those of their spouses or minor children and, for
24 purpose of 18 U.S.C. Section 208, their employers.
25 These interests may include investments, consulting,

1 expert witness testimony, contracts, grants, CRADAs,
2 teaching, speaking, writing, patents and royalties,
3 and primary employment.

4 Today's agenda involves the discussion of a
5 premarket approval application sponsored by Biosense
6 Webster, a Johnson & Johnson Company, for the
7 NaviStar ThermoCool irrigated RF ablation catheter.
8 This device, an open lumen, irrigated tip, steerable
9 radiofrequency cardiac ablation catheter, is inserted
10 through the venous circulation to the heart, across
11 the intra-atrial septum to the left atrium to ablate
12 cardiac tissue for the purposes of creating lines of
13 block in the atria to eliminate conduction patterns
14 that theoretically generate or allow propagation of
15 electrical waves responsible for paroxysmal atrial
16 fibrillation. This is a particular matters meeting
17 during which specific matters related to this PMA
18 will be discussed.

19 Based on the agenda for today's meeting and
20 all financial interest reports by the Panel members
21 and consultants, no conflict of interest waivers have
22 been issued in accordance with 18 U.S.C. Section 208
23 and Section 712 of the FD&C Act. A copy of this
24 statement will be available for review at the
25 registration table during this meeting and will be

1 included as part of the official transcript.

2 Michael Halpin is serving as the industry
3 representative, acting on behalf of all related
4 industry, and is employed by Genzyme Corporation.

5 We would like to remind members and
6 consults that if the discussions involve any other
7 products or firms not already on the agenda for which
8 the FDA participant has a personal or imputed
9 financial interest, the participants need to exclude
10 themselves from such involvement and their exclusion
11 will be noted for the record.

12 FDA encourages all other participants to
13 advise the Panel of any financial relationships that
14 they may have with any firms at issue.

15 Thank you.

16 I will now read the appointment of
17 temporary voting status.

18 Pursuant to the authority granted under the
19 Medical Devices Advisory Committee Charter of the
20 Center for Devices and Radiological Health, dated
21 October 27, 1990, and as amended August 18, 2006, I
22 appoint the following individuals as voting members
23 of the Circulatory System Devices Panel for the
24 duration of this meeting on November 20, 2008.

25 Drs. Valluvan Jeevanandam, Patricia Kelley,

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 Seth Bilazarian, David Slotwiner, Pamela Karasik,
2 John Somberg.

3 In addition, I appoint Jeffery S. Borer,
4 M.D., to act as temporary Chair for the duration of
5 this meeting.

6 For the record, these individuals are
7 special Government employee who have undergone the
8 customary conflict of interest review and have
9 reviewed the material to be considered at this
10 meeting.

11 This was signed by Daniel G. Schultz, M.D.,
12 Director for the Center of Devices and Radiological
13 Health and dated November 14, 2008.

14 Before I turn the meeting back over to
15 Dr. Borer, I have a few general announcements.

16 Transcripts of today's meeting will be
17 available from Free State Reporting, Incorporated.

18 Information on purchasing videos of today's
19 meeting can be found on the table outside the meeting
20 room.

21 Presenters to the Panel who have already
22 done so should provide FDA with a hard copy of their
23 remarks, including overheads.

24 I would like to remind everyone that
25 members of the public and the press are not permitted

1 around the Panel area, beyond the speaker's podium.

2 The press contact for today's meeting is
3 Siobhan DeLancey and Scott McFarland.

4 I request that reporters wait to speak with
5 FDA officials until after the Panel meeting.

6 Thank you.

7 DR. BORER: Good morning, everyone. At
8 this meeting, the Panel will develop recommendations
9 to the FDA, to the Food and Drug Administration, on
10 the Premarket Approval Application, the PMA, P030031,
11 Biosense Webster NaviStar ThermoCool irrigated
12 radiofrequency, RF, ablation catheter.

13 Before we begin, however, I would like to
14 ask our Panel members who are generously giving their
15 time today, and other FDA staff seated at this table,
16 to introduce themselves. As you do so, please state
17 your name, your area of expertise, your position and
18 your affiliation. Why don't we start with
19 Mr. Halpin.

20 MR. HALPIN: Good morning. My name is
21 Michael Halpin. I am the Industry Representative for
22 today's Panel meeting. I'm the Vice President of
23 Regulatory Affairs at Genzyme Corporation. Thank
24 you.

25 DR. FLEMING: Good mooring. I'm Dr. Mike

1 Fleming. I'm the Consumer Representative on the
2 Panel. My area of expertise is dental materials
3 science.

4 DR. JEEVANANDAM: Good morning. My name is
5 Valluvan Jeevanandam. I'm the Chief of Cardiac and
6 Thoracic Surgery at the University of Chicago.

7 DR. KARASIK: I'm Pamela Karasik. I'm the
8 Assistant Chief of Cardiology at the VA here in
9 Washington, D.C., and I'm the Director of Clinical
10 EP.

11 DR. SOMBERG: Good morning. I'm John
12 Somberg. I'm a Professor of Medicine and
13 Pharmacology at Rush University in Chicago, Illinois.

14 DR. NAFTEL: Good morning. I'm David
15 Naftel. I'm a Professor of Surgery and Professor of
16 Biostatistics at the University of Alabama at
17 Birmingham, and my area of expertise is
18 biostatistics.

19 DR. BORER: I'm Jeff Borer as I mentioned
20 earlier.

21 DR. WEINBERGER: I'm Judah Weinberger. I'm
22 an Associate Professor of Medicine and Pharmacology
23 and interventional cardiologist at Columbia, New
24 York.

25 DR. SLOTWINER: I'm David Slotwiner. I'm a

1 cardiac electrophysiologist practicing at Long Island
2 Jewish Medical Center, Assistant Professor of
3 Medicine, Albert Einstein College of Medicine.

4 DR. KELLEY: Good morning. Patricia
5 Kelley. I'm a cardiac electrophysiologist at Montana
6 Heart Center in Missoula, Montana.

7 DR. BILAZARIAN: I'm Seth Bilazarian. I'm
8 a clinical and interventional cardiologist in private
9 practice in Haverhill, Massachusetts.

10 DR. ZUCKERMAN: Brian Zuckerman, Director,
11 FDA, Division of Cardiovascular Devices.

12 DR. BORER: Thank you very much.

13 We'll now proceed with the open public
14 hearing.

15 Both the Food and Drug Administration and
16 the public believe in a transparent process for
17 information gathering and for decision making. To
18 ensure this kind of transparency at the open public
19 hearing session of the Advisory Committee meeting,
20 the FDA believes that it is important to understand
21 the context of any individual's presentation. For
22 this reason, FDA encourages you, the open public
23 hearing or industry speaker, at the beginning of your
24 written or oral statement, to advise the Committee of
25 any financial relationship you may have with the

1 sponsor, its product, and if known, its direct
2 competitors.

3 For example, this financial information may
4 include the sponsor's payment of your travel,
5 lodging, or other expenses in connection with your
6 attendance at the meeting. Likewise, FDA encourages
7 you at the beginning of your statement to advise the
8 Committee if you do not have any such financial
9 relationships. If you choose not to address this
10 issue of financial relationships at the beginning of
11 your statement, it will not preclude you from
12 speaking.

13 Does anyone wish to address the Panel at
14 this time?

15 (No response.)

16 DR. BORER: There do not seem to be any.
17 Do we have anybody listed?

18 MR. SWINK: No.

19 DR. BORER: Okay. We will then proceed
20 with today's agenda. Please note there will be a
21 second open public session in the afternoon if there
22 are any comments to be made by anyone, any member of
23 the public.

24 We'll begin with the sponsor presentation
25 for the NaviStar ThermoCool Irrigated Radiofrequency

1 Ablation Catheter.

2 I'd like to remind public observers at this
3 meeting that while this meeting is open for public
4 observation, public attendees may not participate
5 except at the specific request of the Panel. That's
6 true for the Sponsor as well. Once you give your
7 presentation, any further comments really have to be
8 made in the context of a request from the Panel or
9 another specific opportunity in the program.

10 We'll then begin with the sponsor
11 presentation.

12 DR. YAROSS: Thank you, Dr. Borer. Good
13 morning. I'm Marcia Yaross, and I'm an employee of
14 Biosense Webster. I thank the Panel and the FDA for
15 the opportunity to present this morning on the
16 clinical evidence in support of the safety and
17 effectiveness of the ThermoCool catheter for
18 radiofrequency ablation of systematic paroxysmal
19 atrial fibrillation.

20 In this morning's presentation, I will
21 briefly introduce you to the clinical trial with an
22 overview of its history and design.

23 Dr. Donald Berry will then provide an
24 overview of the Bayesian statistics at the FDA's
25 request and then speak to the statistical rationale

1 for declaring early success following a planned
2 interim analysis.

3 Dr. David Wilber, the study's primary
4 investigator, will present the subject demographics
5 and effectiveness results for the trial.

6 Dr. Waldo, who chaired the study's clinical
7 events committee will present the safety data for
8 this trial. I'll then conclude our presentation.

9 Atrial fibrillation is an important public
10 health issue. AF represents the most prevalent
11 arrhythmia encountered in clinical practice. Recent
12 Miyasaka data estimate that between 2.3 and 5 million
13 U.S. adults have atrial fibrillation. AF is a
14 debilitating disease, particularly the symptomatic
15 paroxysmal form.

16 Patients with atrial fibrillation are at
17 increased risk of stroke, of heart failure, and other
18 significant comorbidities. AF has also been found to
19 be an independent marker of risk of death.

20 Patients with atrial fibrillation have a
21 significantly reduced quality of life.

22 While pharmacological therapy is available,
23 it has been proven to be ineffective in many AF
24 patients, estimated up to 50 percent in some series.

25 Surgical techniques such as the Cox-Maze

1 procedure can be effective, but they are highly
2 invasive with the associated morbidity and mortality
3 one would expect from an open chest cardiovascular
4 procedure.

5 In response to this need, radiofrequency
6 catheter ablation has become increasingly important
7 as a tool in the toolkit of electrophysiologists. RF
8 ablation is today a standard of care for simple
9 arrhythmias such as Wolff-Parkinson-White Syndrome,
10 atrial flutters, and AV node reentry and
11 tachycardias.

12 It's also increasingly being used for more
13 complex arrhythmias. This has been recognized by a
14 number of FDA approved indications for VT procedures.

15 Catheter ablation is also increasingly used
16 for atrial fibrillation, which is the subject of this
17 morning's discussion.

18 Catheter ablation was formally recognized
19 in 2006 as a second line therapy for atrial
20 fibrillation in the American College of Cardiology,
21 AHA, ESC practice guidelines, and we'll discuss this
22 some more this morning.

23 HRS and several European EP societies also
24 affirmed the importance of ablation in treating AF in
25 their expert consensus released in 2007. This

1 consensus statement identified catheter ablation as
2 indicated for symptomatic AF where it is refractory
3 or intolerant to at least one class I or class III
4 antiarrhythmic drug. The task force also recognized
5 in some rare circumstances ablation may be
6 appropriate as first line therapy. The document also
7 stated that it's appropriate for selected patients
8 with heart failure and/or reduced ejection fraction.

9 Despite this growth of support of the
10 recognition of the importance of AF, however, no
11 ablation catheter is today approved for treatment of
12 AF in the United States.

13 The ThermoCool catheters, which are the
14 subject of this morning's discussion, are well
15 established in the practice of electrophysiology.
16 They have been used in nearly 40 countries across the
17 world over the past decade. A quarter of a million
18 catheters have been distributed worldwide since their
19 initial introduction overseas, and about a quarter of
20 that experience has been in the United States since
21 they were first FDA approved for the treatment of
22 atrial flutter in 2004.

23 A ThermoCool catheter is a steerable,
24 multi-electrode, deflectable electrophysiology
25 catheter with saline irrigation provided through six

1 ports in the tip. These ports irrigate and cool the
2 catheter tip to maintain low temperature during RF
3 application. This thereby reduces the risk of char
4 or thrombus. There is a temperature sensor for
5 feedback to assure that the irrigation is cooling the
6 tip as intended. However, the catheters are intended
7 for use in power control and not temperature control
8 mode.

9 The NaviStar version used in this clinical
10 trial also contains a location sensor for
11 visualization with the CARTO electroanatomical
12 mapping system.

13 ThermoCool catheters are today approved by
14 the FDA for two indications: Type 1 atrial flutter
15 as well as post-MI ventricular tachycardia. Please
16 note that the atrial flutter indication applies to
17 both the location sensor enabled NaviStar models as
18 well as well as the Celsius catheters which lack such
19 sensors. The VT indication is specific to the
20 NaviStar Catheters.

21 We are now asking FDA to extend these
22 indications to include treatment of symptomatic AF as
23 detailed in your Panel packages.

24 The clinical study we're presenting to you
25 was a randomized clinical trial. Randomization was 2

1 to 1 with more subjects randomized to the ThermoCool
2 treatment arm. The population consisted of
3 symptomatic paroxysmal patients refractory to at
4 least one antiarrhythmic drug. The study was multi-
5 centered, and the final protocol included planned
6 interim analyses beginning with 150 subjects.

7 The primary study goals were to demonstrate
8 superior chronic effectiveness of the ThermoCool
9 catheter versus antiarrhythmic drug treatment in the
10 prevention of symptomatic AF recurrence. This was
11 measured during a nine-month effectiveness window in
12 each arm. The safety goal was to have an acceptable
13 safety profile versus a prespecified performance
14 goal.

15 Chronic effectiveness was evaluated during
16 comparable nine-month evaluation periods. The
17 ThermoCool group had a three-month blanking period
18 during which repeat ablation could be performed up to
19 80 days. Freedom from symptomatic AF recurrence was
20 then measured in accordance with protocol criteria
21 from days 91 to 361.

22 The control group had a two week dosing
23 window to titrate, the drug regimen to maximum
24 effectiveness. Chronic success was then evaluated
25 from days 15 to 285.

1 Freedom from symptomatic AF recurrence was
2 monitored throughout the effectiveness evaluation
3 period through regular transtelephonic ECG
4 monitoring.

5 The primary safety endpoint was comprised
6 of 18 serious adverse events occurring during the
7 first 7 days after an ablation procedure, and this
8 included new or prolonged hospitalization for any
9 reason. The primary safety performance goal was
10 prospectively established based on a literature
11 review. Incidence of pulmonary vein stenosis was
12 also deemed a primary safety measure, although no
13 quantitative hypothesis was established.

14 The study enrolled patients who had been
15 diagnosed with symptomatic AF and who had experienced
16 at least three episodes in the proceeding six months.
17 All had failed at least one antiarrhythmic drug.
18 Either rate or rhythm control agents were considered
19 AADs for study inclusion purposes.

20 Exclusion criteria were typical of
21 paroxysmal atrial fibrillation trials with exclusion
22 for amiodarone therapy in the six months prior to
23 enrollment in accordance with the FDA guidance
24 document.

25 I'd now like to take you through a number

1 of the key milestones in this study. Enrollment in
2 this trial was extremely challenging with nearly a
3 year elapsing between IDE approval and enrollment of
4 the first subject. Only 10 additional subjects were
5 enrolled in the next 12 months despite aggressive
6 advertising and subject recruitment efforts.

7 Over the next year, we managed to get up to
8 53 total subjects with a combination of sites both
9 inside and outside of the United States. We continue
10 to add high volume ablation sites, and as you can
11 see, the pace of enrollment increased substantially.
12 Nonetheless, we were still looking at a prolonged
13 time span to reach the original population of 230
14 subjects.

15 We therefore approached the FDA to discuss
16 the use of a Bayesian statistical analysis approach
17 to facilitate early completion of the trial should
18 the results so warrant.

19 During this period of time, the FDA also
20 convened a meeting of this Advisory Panel during
21 which the barriers to enrollment in AF IDE studies
22 were discussed in detail, and I note that a number of
23 you were here for that Panel meeting which was 14
24 months ago today.

25 FDA approved our protocol amendment for a

Free State Reporting, Inc.
1378 Cape Saint Claire Road
Annapolis, MD 21409
(410) 974-0947

1 Bayesian analysis plan in early September of 2007.
2 This allowed us to conduct our first interim analysis
3 in a prospective manner. This interim analysis was
4 done independently by our outside statisticians,
5 Berry and Associates. Based on this interim analysis
6 of a December 2007 dataset, we were able to declare
7 early success, and therefore subject enrollment was
8 terminated in early October of 2007.

9 Our subject database was locked in June
10 2008, after thorough site monitoring, and so it's the
11 June 2008 dataset that you will see for most of the
12 analyses presented today.

13 To further illustrate the challenges to
14 enrollment, about 5500 patients were screened in
15 order to enroll 167 subjects in the trial. While the
16 trial only enrolled 3 percent of total screened
17 candidates over 3 years, about 1/3 of the subjects
18 actually did meet inclusion criteria as shown in this
19 graph presented at last year's September Panel
20 meeting. You'll note that 62 percent were excluded
21 by protocol requirements. There were others that
22 were unable to return for personal reasons, refused
23 randomization, et cetera.

24 Final enrollment was therefore comprised of
25 167 subjects at 19 centers. While the largest two

1 centers were outside the United States, no center
2 exceeded 30 percent of the total population.

3 Poolability among the 19 centers is based
4 on a common protocol, identical data collection
5 instruments and rigorous monitoring proportional to
6 the number of subjects enrolled at each site.

7 With that, I'd like to introduce Dr. Donald
8 Berry. First, he will provide an overview of
9 Bayesian statistics at the request of the FDA. He
10 will then present the statistical analysis plan for
11 the study as well as the interim analysis that
12 provided the early stopping decision and declaration
13 of success.

14 DR. BERRY: Thank you, Dr. Yaross. My name
15 is Donald Berry. I'm a statistician at the M. D.
16 Anderson Cancer Center and a member of Berry
17 Consultants, a consultant to the company. Berry
18 Consultants designed the Bayesian aspect of the
19 ThermoCool trial and carried out the interim
20 analysis.

21 I want to give you a primer on Bayesian
22 methods. Dr. Laura Thompson of the FDA is going to
23 give you the FDA perspective. I think these two
24 presentations are supplementary and complementary.
25 There's a bit of redundancy, but as usual, if you

1 hear something twice, it means it's important.

2 In 1997, the FDA Modernization Act used the
3 term least burdensome. The Secretary shall consider
4 in consultation with the applicant the least
5 burdensome, appropriate means of evaluating device
6 effectiveness that would have a reasonable likelihood
7 of resulting in approval. The first Bayesian
8 approval at CDRH was in 1997. Draft Guidance was
9 presented in 2006.

10 This is the set of table of contents of the
11 Draft Guidance and the website where it's available,
12 and we chose the use of the least burdensome approach
13 and its relationship to Bayesian statistics.

14 Current use of Bayesian adaptive designs at
15 my home institution, M. D. Anderson Cancer Center,
16 we've designed and have run or are running over 300
17 trials since I arrived there nine years ago. Many
18 device companies are using the Bayesian approach.
19 There have been over 20 PMAs and many 510(k)'s. In
20 the last few years, virtually all of the top drug
21 companies have been taking the Bayesian adaptive
22 approach in at least some of their clinical trials
23 and many biotechs.

24 Some of the Bayesian device applications
25 that have been submitted and some approved by the FDA

1 are listed here.

2 Bayes' Rule is the basis of the Bayesian
3 approach and many have been exposed to the use of
4 Bayes' Rule. Many M.D.'s, for example, in the
5 context of diagnoses, diagnostic tests, many are
6 familiar with sensitivity of a diagnostic test.
7 Specificity, positive predictive value.

8 The bottom formula here shows the
9 relationship between sensitivity, the probability
10 that a test is positive if an individual has a
11 disease, and the positive predictive value. This is
12 what we really want, what is the probability of the
13 individual has the disease if he or she tests
14 positive and what we have is sensitivity and
15 specificity, not listed here.

16 This crazy looking means proportional.
17 There's a constant that's not included here.

18 The relationship between the inverse
19 relationship as given by Bayes' Rule, an aspect which
20 is crucial, is the so-called prevalence of the
21 disease, the proportion of individuals who have the
22 disease among those in the population that we are
23 testing.

24 Same relationship applies up here. This is
25 called the likelihood function, the probability of

1 the data that you've observed given a hypothesis or
2 the value of a parameter, an unknown parameter, and
3 this is the posterior probability, the probability of
4 the hypothesis given the data. The analog of the
5 prevalence is the probability, the prior probability
6 of the hypothesis, which is a fundamental part in the
7 use of the Bayesian approach.

8 The Bayesian approach is a formalism for
9 learning under uncertainty. Anything that's not
10 known has a probability distribution that includes
11 hypotheses and all hypotheses, parameter values, and
12 future data. Hypothesis test is the posterior
13 probability of no treatment effect. The analog of
14 the confidence interval is an interval that contains
15 the parameter with a particular posterior probability
16 such as 95 percent.

17 The Bayesian approach is inherently
18 synthetic. The probability of a hypothesis given
19 data means probability of a hypothesis given
20 everything that you know, which includes not only the
21 study at hand but other studies that are going on or
22 have been conducted.

23 Advantages of the Bayesian approach is
24 naturally adaptive, and I want to focus a bit more on
25 that. It leads to the ability to calculate

1 predictive probabilities. It uses early by-patient
2 information and, for example, in the ThermoCool
3 trial, we modeled the hazard as being potentially
4 different than the early time periods, and you see
5 that indeed in Dr. Wilber's presentation, that the
6 failure rate vary over the course of time and
7 remodeled that aspect uses historical data --
8 hierarchical modeling. That's an aspect that we, in
9 fact, did not use in ThermoCool trial.

10 A Bayesian update, I want to take you
11 through a very simple example. Suppose there are
12 paired observations within each pair. One of the
13 members gets a treatment and the other gets control.
14 The probability that the treatment wins the pair is
15 P_S here, and the null hypothesis is that that
16 probability is a half. If there's no difference
17 between the two, half of the pairs would be won by T
18 and half by C.

19 This is the first 17 observations in this
20 particular example. The first pair was won by the
21 treatment of success. The second pair won by the
22 treatment again of success, then a failure of the
23 first 10, you see there were 7 successes and 3
24 failures.

25 The Bayesian approach, as we said, starts

1 with a prior probability distribution, and this is a
2 so-called non-informative prior or plat prior. It
3 uses a minimal of historical information. It regards
4 each value of P , the probability of success, as being
5 as likely as every other value.

6 It also has the characteristic that it
7 leads to at least approximately the same type of
8 measures as the more standard, frequentist approach,
9 such as the confidence interval is actually a
10 Bayesian posterior probability interval.

11 After the first observation, which you
12 remember was a success for the treatment, this
13 distribution changes, and it changes by Bayes' Rule.
14 The technical aspect is that you multiply by the
15 likelihood the probability of what you actually
16 observe. You observe the success, and so the
17 probability changes by multiplying by the probability
18 of success, which is P itself. After the next
19 observation, there's another success, and you get P
20 squared as the new posterior distribution. After the
21 next, you multiply it by the probability of failure,
22 which was the third observation. So you get P
23 squared times $1 - P$.

24 The technical aspects are not too important
25 here. What I'm interested in demonstrating is that

1 you can do this. You make an observation and you
2 update what you know on the basis of that previous
3 observation, and you can do that daily or monthly or
4 periodically.

5 So just showing you after the next success,
6 after the next success, another failure, it shifts a
7 little bit to the left. After 10 observations, we're
8 at this point here. This is final after the 10.
9 It's the current distribution.

10 The other thing I want to show you is what
11 happens tomorrow. So the next pair of observations
12 will either give you a success in which case this
13 moves up a little bit or failure in which case this
14 moves down a little bit. The Bayesian approach is
15 unique in allowing for the calculation of the
16 probability of those two outcomes, and that is a
17 well-known Laplace's rule of succession. In this
18 instance, it says that the probability of its success
19 is $8/12$ or $2/3$. The probability of a failure is $1/3$.

20 Predicted probabilities are essential for
21 monitoring trials. It's a critical component of
22 experimental design and were a critical component in
23 the ThermoCool trial, and I give a quote here from
24 the famous Bayesian clinical trialist, we must ask,
25 we must ask where we are and whether we are tending.

1 This is the current distribution after 17
2 observations and, you know, the most likely values
3 are in the .7, .8 range. The Bayesian approach
4 allows for calculating the probability of
5 superiority, which means that the treatment success
6 is bigger than a half, and that's the area under this
7 posterior distribution which indicates again is .985.

8 Remember, we have 17 observations, 13
9 successes. If you say, well, how useful will it be
10 to increase the sample size to go to the next 17
11 observations, let's say, this is the predictive
12 probability distribution of the next 17. So ranging
13 from 0 to 17, and it incorporates two types of
14 uncertainty. One is the future, the variability, the
15 sampling variability, but also the variability in P
16 and the variability in P, what we know today is
17 critical to incorporate into this distribution. If
18 you didn't incorporate it, if you used what's called
19 a binomial distribution, a typical coin tossing, with
20 the most likely value of P, you'd get a distribution
21 which is much more concentrated. Notice the tail of
22 this distribution is greater than here and the same
23 on this side. This is the correct distribution.
24 This is artificially assuming that the variability is
25 less than it truly is in the future observations.

1 So an instance of that is supposedly to
2 calculate the probability of statistical
3 significance. After the next 17 observations, that
4 would require 10 of the 17 to be successes or 23 of
5 the 34 observations, the actual probability is 88
6 percent. If you assume that you knew the value of P,
7 which, of course, you don't, it's 96 percent. This
8 artificially inflates what we know and actually
9 explains some of the causes of failure of phase 3
10 trials in the drug world.

11 The Bayesian approach is rather more
12 conservative than this fixed P approach.

13 Important issues to discuss, the Bayesian
14 approach must have something fundamental in clinical
15 research as to have a prospective design and so even
16 though the Bayesian approach is adaptive, uses
17 information to possibly modify the trial design, that
18 has to be specified completely up front in order to
19 calculate operating characteristics like false
20 positive rate, controlling type on error, and in the
21 case at hand, require an agreement with the FDA as to
22 that this was an appropriate design before anything
23 happened.

24 Changing from frequentist to Bayes, this is
25 something that's addressed in the Bayesian guidance,

1 the FDA Draft Guidance, and it's discouraged. In the
2 instance of the ThermoCool trial, we discussed this
3 with the FDA very early on, and they agreed that
4 since most of the data was yet to accrue, they okayed
5 the switch to the Bayesian approach.

6 Why do people do Bayesian things? The hook
7 is smaller trials, and I say usually. Sometimes it's
8 actually bigger. Sometimes you get to the end of the
9 trial and you wished the trial were not over, and so
10 you can build that in prospectively to allow for
11 increasing the sample size. When I say you wish the
12 trial weren't over, I don't mean that you haven't
13 answered the question in the way that you want to. I
14 mean that you are on the cusp. You don't know
15 whether you have a device which is safe and
16 effective.

17 That leads to more accurate conclusions and
18 in many examples, and the thing that is so appealing
19 to physicians of patients in cancer is potentially
20 better treatment of patients in trials.

21 So now I'll turn to the ThermoCool AF
22 trial. The amended design of the Bayesian sample
23 size was prospective -- rate was controlled with
24 planned interim analyses, the technical aspect of
25 this, is that we simulated many, many trials under

1 many different scenarios, including the scenario
2 where the device has no benefit in comparison to the
3 control, AAD. The Bayesian adaptive sample size kept
4 the maximum, the previous maximum of 230 but allowed
5 for stopping on the basis of predictive probabilities
6 when we get to 150. And as Dr. Yaross indicated,
7 that occurred at the time that the FDA agreed to the
8 Bayesian approach in September 2007. We performed
9 the interim analysis, and it showed two things: one,
10 that the accrual could stop, and the other was that
11 the data were sufficiently compelling that, in fact,
12 we could claim early success and submit an
13 application.

14 So these were the results, and this is in
15 contrast to it, not contrast to but anticipates the
16 presentation of Dr. Wilber. These were the interim
17 results, and he will be presenting you the final
18 results or the current results, and his presentation
19 is the right one in the sense that the Bayesian
20 approach, even though it asks where we're going, when
21 you eventually get there, what matters is what you
22 have at the end of the day.

23 So there are 148 patients eligible for the
24 interim analysis. The predictive probability of
25 eventual success which was defined to be at the end

1 of the day, probability of superiority of at least
2 .98, the probability of that had to be at least 90
3 percent in order to stop accrual. So we could stop
4 accrual at the sample size, 150, if we had a
5 predictive probability of success of at least 29, and
6 if we had a predictive probability of success of at
7 least 29.9, then we could declare a success at that
8 time. We would, of course, continue to follow and
9 the company did that.

10 The result of the interim analysis was a
11 predictive probability of, in fact, greater than
12 .999, and therefore early success was declared.

13 I'll now introduce Dr. David Wilber who
14 will give you the present results of the trial.

15 DR. WILBER: Good morning. My name is Dave
16 Wilber. I'm Director of Cardiology at Loyola
17 University Medical Center. I'm the primary
18 investigator for the study, and I'm a consultant for
19 Biosense Webster. What I'll be presenting are the
20 final results of the study and demographics as a
21 follow-up in June of 2008.

22 Overall, there were 167 patients who met
23 the study criteria, were consented and randomized,
24 106 in the ThermoCool group and 61 in the
25 antiarrhythmic drug group.

1 There were seven patients that were
2 excluded from the study. Five of these were because
3 of withdrawal of consent after randomization. One
4 patient who was assigned to ablation, the insurance
5 company didn't approve the ablation procedure and
6 reimbursement for it, and so that patient was not in
7 the study. And then one patient was found
8 subsequently not to meet the enrollment criteria and
9 was excluded.

10 One additional subject in the
11 antiarrhythmic drug group actually underwent initial
12 treatment with the assigned drug and then decided to
13 withdraw consent for follow-up.

14 So the final cohort is 103 in the ablation
15 group, 56 in the drug group, for a total of 159
16 patients who comprised the effectiveness analysis
17 cohort.

18 These are the demographics. The two groups
19 were well matched in terms of gender and age.
20 Approximately a third of the patients enrolled were
21 female, and the mean age was approximately 56 years.

22 Patients in this trial were highly
23 symptomatic, and there was a mean of approximately 63
24 episodes of symptomatic afib that were reported in
25 the 6 months prior to randomization, and these were

1 equally distributed between the two groups.

2 This slide demonstrates the cardiac
3 comorbidities at baseline in the two treatment
4 groups, red being antiarrhythmic drugs, blue being
5 ThermoCool. As you can see, there were no
6 significant differences between these two groups.
7 Approximately half of the group had hypertension.
8 Approximately 15 percent structural heart disease.
9 These are very typical comorbidities in a population
10 of patients with paroxysmal atrial fibrillation.

11 This slide summarizes prior antiarrhythmic
12 drug experience before randomization. There was an
13 average of 2.2 antiarrhythmic drugs failed, and you
14 can see there's no difference between the ThermoCool
15 and antiarrhythmic drug group. This included a mean
16 of 1.5 class I or III antiarrhythmic drugs and a mean
17 of 1.3 class II or IV antiarrhythmic drugs.

18 Patients who were assigned to
19 antiarrhythmic therapy received a new, not previously
20 administered class I or III drug. The minimum
21 recommended dosing was based on the ACC, AHA, ESC
22 guidelines published in 2001. These were all drugs
23 that were approved at that time of study onset for
24 treatment of atrial fibrillation, and they included
25 sotalol, dofetilide, flecainide, propafenone and

1 quinidine.

2 The prescribed antiarrhythmic drug was
3 adjusted in dose for maximum efficacy during the 14-
4 day dosing period, and then the antiarrhythmic drug
5 and dose were fixed at day 15 and remained on the
6 same drug and dose for the chronic follow-up.

7 Amiodarone therapy was not an option by protocol
8 definition.

9 In the ablation group, it was required that
10 all patients have circumferential isolation of the
11 pulmonary veins and that there be electrophysiologic
12 confirmation of entrance block into the pulmonary
13 veins as the acute procedural endpoint.

14 CARTO electroanatomical mapping was used in
15 all patients.

16 At the discretion of the investigator, and
17 depending on the outcome of the procedure after
18 pulmonary vein isolation, patients could undergo
19 isolation of the superior vena cava if that was an
20 initiator of atrial fibrillation. They could undergo
21 ablation of non-PV -- that initiated atrial
22 fibrillation. They could undergo left atrial linear
23 lesions if atrial fibrillation could be induced after
24 pulmonary vein isolation. They could undergo a left
25 inferior pulmonary vein to mitral isthmus line if

1 left atrial flutter was induced. And finally, they
2 could undergo cavotricuspid isthmus ablation if
3 isthmus-dependent right atrial flutter was induced.

4 This slide addresses the time from
5 randomization to initial treatment in each group. In
6 the ablation group, a mean of 28 days and a mean of
7 43 days between randomization and treatment with
8 catheter ablation. This largely reflected the issues
9 with simply getting patients on the schedule and
10 having the procedure performed.

11 In the drug group, there was a mean of 10
12 days and a mean of 16 days from randomization to the
13 initiation of antiarrhythmic drug therapy. While
14 there was a slightly longer delay in initiating
15 treatment in the ablation group, given that atrial
16 fibrillation is in general a progressive disease, we
17 feel that it is unlikely that this biased the outcome
18 in favor of the ablation group.

19 This slide summarizes acute effectiveness
20 outcome for the protocol definition. 103 patients
21 underwent the ablation procedure, and entrance block
22 was confirmed in all 103 patients. However, 2
23 patients underwent a second ablation procedure
24 between day 80 and 90 of the blanking period. That
25 was outside the protocol-defined window. That

1 resulted in an acute effectiveness success of 101 of
2 103 patients or 98 percent.

3 Chronic success in the ThermoCool group was
4 defined by protocol as freedom from the following:
5 documented symptomatic afib recurrence during the
6 follow-up period; freedom from acute procedural
7 failure, irrespective of afib recurrence, and that
8 included failure to confirm entrance block into the
9 targeted pulmonary veins; EP afib ablation procedure
10 after 80 days; and freedom from protocol adjudicated
11 antiarrhythmic drug failure, again irrespective of
12 afib recurrence. And these drugs included class I
13 and III antiarrhythmic drugs but also beta blockers,
14 calcium channel blockers, digitalis, ARBs and ACE
15 inhibitors.

16 In the drug group, chronic success was
17 defined by protocol as freedom from the following:
18 documented symptomatic afib recurrence during the
19 efficacy evaluation period; protocol adjudicated
20 antiarrhythmic drug failure defined similarly to that
21 in the ablation group; and also safety failure which
22 required discontinuation of the assigned
23 antiarrhythmic drug during the efficacy evaluation.

24 A standardized transtelephonic monitoring
25 protocol was followed within both groups. Subjects

1 were instructed to transmit once a week for the
2 initial eight weeks and monthly for the remaining
3 seven months. They were also instructed to transmit
4 during any cardiac symptoms. These recordings were
5 initially reviewed by two laboratory technicians that
6 provided the initial interpretation and then finally
7 all recordings were reviewed by an independent
8 cardiologist who was blinded and adjudicated the
9 final outcome of the recording.

10 This slide summarizes the compliance to the
11 standard TTM protocol in both treatment groups. As
12 you can see, there was high compliance over time.
13 This was very similar between both treatment groups
14 and the average was 89 percent compliance with the
15 prespecified time periods for transmission of data.

16 Bayesian analysis was then used to examine
17 the significance of the treatment differences in the
18 two groups, and the critical results of this analysis
19 are the predicted probability of study success for
20 230 patients and the posterior probability of
21 superiority for the ThermoCool group. The posterior
22 probability that the ThermoCool group was superior to
23 the antiarrhythmic group is essentially 1. The
24 probability of success for a subject in the
25 ThermoCool group is 62.7 plus or minus 4.8 percent,

1 and the probability of success for a subject in
2 antiarrhythmic drug group is 17.2 plus or minus 4.9
3 percent.

4 This slide graphically illustrates that
5 outcome and shows the distribution of probability for
6 the antiarrhythmic drug and ThermoCool group, and as
7 you can see, there's overlap in the probability
8 distribution between those two groups.

9 This slide demonstrates the Kaplan-Meier
10 curve of time to first chronic failure, protocol
11 defined, by randomization group. At the end of
12 follow-up, 64 percent of ThermoCool ablation patients
13 and only 16 percent of antiarrhythmic drug patients
14 were free of any chronic failure as defined by the
15 protocol.

16 These circles in the graph represent the 14
17 censored ThermoCool subjects who at the time of the
18 June database had not yet completed follow-up.

19 This slide summarizes again the chronic
20 effectiveness failures in the ThermoCool group.
21 Twenty-four patients failed because of recurrent
22 symptomatic atrial fibrillation or 23 percent of the
23 ablation group. The remaining failures were not due
24 to symptomatic atrial fibrillation recurrence but
25 were due to protocol differences, the first one being

1 2 patients who had a re-ablation within 80 to 90
2 days, and 10 patients who had protocol adjudicated
3 antiarrhythmic failures, again not due to recurrence
4 but because of the initiation, in a small number of
5 patients, a class I or III antiarrhythmic drug, but
6 in the majority because of a new beta blocker,
7 calcium channel blocker, ACE inhibitor or ARB.

8 Similarly, in the chronic effectiveness
9 failures, there were 47 failures in the
10 antiarrhythmic drug group or 71 percent of the
11 control subjects; 40 of these were due to symptomatic
12 afib recurrence. In 7 patients, there was no
13 symptomatic recurrence, but all 7 failed due to
14 intolerance or serious adverse events related to the
15 prescribed antiarrhythmic drug.

16 Several additional considerations we should
17 review. In the ThermoCool group, 50 percent of the
18 patients underwent PV isolation alone. In the
19 remaining patients, other procedures more than one in
20 a single subject included cavotricuspid isthmus
21 ablation for flutter, 34 percent, superior vena cava
22 isolation in 16, other focal drivers in 17 percent, a
23 mitral isthmus line in 21 percent, and other left
24 atrial lines in 20 percent.

25 As allowed by the protocol, 13 subjects or

1 12 percent underwent a second procedure within the
2 first 80 days of the blanking period. The protocol
3 did allow for the use of previously failed
4 antiarrhythmic drugs during follow-up, but this was
5 limited to only 7 percent of subjects classified as
6 success during the 6 months of follow-up.

7 In the antiarrhythmic drug group, the vast
8 majority of patients were treated with either
9 flecainide or propafenone. Sotalol and dofetilide
10 accounted for the remaining patients in terms of
11 their assigned treatment drug.

12 As allowed by protocol, 64 percent of the
13 antiarrhythmic drug group ultimately had an ablation
14 procedure after symptomatic recurrence and
15 classification as a treatment failure.

16 To better characterize the effectiveness
17 results, we conducted several post-hoc Kaplan-Meier
18 analyses. These included freedom from symptomatic
19 atrial fibrillation recurrence and freedom from any
20 atrial fibrillation recurrence, either symptomatic or
21 asymptomatic.

22 This slide demonstrates a Kaplan-Meier
23 curve of time to first symptomatic afib recurrence.
24 As you can see, there's a dramatic difference between
25 the two groups. At the end of the follow-up period,

1 75 percent of the ThermoCool ablation patients, and
2 only 21 percent of the antiarrhythmic drug patients,
3 were free of any symptomatic atrial fibrillation
4 recurrence.

5 Similarly, for total AF recurrence,
6 including both symptomatic and asymptomatic events,
7 at the end of the treatment period or rather at the
8 end of the follow-up period, 72 percent of the
9 patients in the ThermoCool group and only 21 percent
10 of patients in the antiarrhythmic drug group were
11 free of any recurrence of atrial fibrillation.

12 Eleven of the antiarrhythmic drug subjects
13 were prescribed the same or a higher dose of a
14 previously failed antiarrhythmic drug. We did
15 perform sensitivity analysis removing these 11
16 subjects from the antiarrhythmic drug control group,
17 and the results were consistent with the primary
18 analysis to continue to show the superiority for the
19 ThermoCool group at the highly significant P value.

20 In addition, four antiarrhythmic drug
21 subjects received less than the protocol recommended
22 antiarrhythmic dosage. That included one subject
23 that was also in the above group. When Bayesian
24 multiple -- analysis was conducted for these 14
25 subjects receiving less than the protocol specified

1 minimum antiarrhythmic dosage, superiority was still
2 demonstrated.

3 We also analyzed effectiveness outcomes by
4 site and region, and this included OUS 1 versus
5 remaining sites and non-US versus US sites.

6 This slide is a Kaplan-Meier curve of time
7 to chronic failure per protocol at OUS site 1. There
8 were no protocol chronic failures and there were, and
9 this was the recurrence in the antiarrhythmic drug
10 group. This is the comparison of that one. This is
11 the comparison in excluding that site. Even in this
12 group of patients, there remained a significant
13 difference between the ThermoCool and antiarrhythmic
14 drug treated patients although the less, a somewhat
15 smaller magnitude.

16 But there's several potential reasons why
17 there were some differences between that site and the
18 remaining sites. First of all, this was one of the
19 highest volume AF ablation centers worldwide. It's
20 had access to the ThermoCool catheter since 1999,
21 whereas many of the, particularly in the United
22 States, the first experience with the catheter was at
23 the onset of the clinical trial since the catheter
24 has not been available commercially released until
25 two years ago.

1 There were minor differences in baseline
2 demographics in favor of the OUS site 1 in that there
3 were smaller atrial size, less hypertension, and
4 somewhat younger subjects.

5 There were also differences in procedures
6 between the OUS site 1 and the remaining sites.
7 Cavo-tricuspid isthmus ablation was performed more
8 often, 23 out of 31 compared to 13 out of 72. In
9 addition, at OUS site 1, left atrial linear lesions
10 were performed more frequently, 20 of 31 at OUS 1
11 versus 9 of 72 at the remaining sites.

12 In the patients assigned to ablation, there
13 were four early recurrences at OUS site 1, and all of
14 these patients went re-ablation within 80 days. In
15 contrast, there were only 9 re-ablations within the
16 protocol specified period in all of the remaining
17 sites. So the OUS 1 was very good at getting their
18 patients back and having them re-ablated within the
19 protocol if they had an early recurrence.

20 And finally, there were some differences in
21 medical management. One of them was that since it
22 was possible to continue a previously failed
23 antiarrhythmic drug post-ablation, this was done
24 somewhat more often at this site. However, there
25 were typically continued for only three to six

1 months, and as we had stated earlier, only a small
2 number of patients remained on these drugs at the end
3 of follow-up.

4 And probably and very important as we point
5 out in subsequent analyses, this site was very good
6 at strict protocol compliance so that non-compliance
7 with other drugs, not necessarily antiarrhythmic
8 drugs, but beta blockers, calcium channel blockers
9 and ARBs that were the source of protocol adjudicated
10 failures in U.S. and other non-OUS 1 sites, never
11 occurred at this particular site. So for the primary
12 effectiveness criterion, this made a substantial
13 impact on differences and outcome.

14 To examine this further, a Bayesian
15 analysis was conducted excluding subjects from OUS 1.
16 The resulting posterior mean probability of success
17 in all of the remaining sites, again excluding OUS 1,
18 was a 46 percent chronic success in the ThermoCool
19 group versus the 20 percent success in the
20 antiarrhythmic drug group, and the posterior
21 probability that the ThermoCool group is superior to
22 the antiarrhythmic drug group remains very high at
23 0.9975.

24 Finally, there were a variety of
25 sensitivity analyses conducted, varying strengths of

1 borrowing of OUS 1 and remaining OUS site data. Even
2 if one heavily discounts the OUS 1 and remaining OUS
3 sites, the result is still very compelling. For
4 example, if one borrows only 20 percent or discounts
5 by 80 percent, the OUS 1 data and remaining OUS
6 sites, the probability of superiority remains at
7 0.991.

8 This impact, I think, is made more clear by
9 additional effectiveness outcomes that excluded OUS
10 1. If we look at time to symptomatic afib recurrence
11 on the Kaplan-Meier curve, and time to any observed
12 AF recurrence in the bottom curve, again these are
13 from all of the remaining centers excluding OUS 1.
14 There was still a highly meaningful and strong
15 difference between the two groups. At the end of the
16 follow-up period, 64 percent of the ThermoCool group
17 and 26 percent of the antiarrhythmic drug group
18 remained free of symptomatic atrial fibrillation.
19 And for any recurrence of atrial fibrillation, 60
20 percent of the ThermoCool and 26 percent of the
21 antiarrhythmic drug patients remained free.

22 So again, this particularly suggests that
23 many of the differences in protocol defined success
24 versus the actual recurrence of atrial fibrillation
25 involved these differences in the use of non-

1 antiarrhythmic drugs such as beta blockers, calcium
2 channel blockers, and ARBs in sites other than OUS 1.

3 Finally we looked at the outcome of the
4 United States population alone. This is the Kaplan-
5 Meier curve from time of first chronic failure per
6 protocol at U.S. sites. There remains a difference
7 between the ThermoCool group and the antiarrhythmic
8 drug group of almost twofold, 44 percent of the
9 patients in the ThermoCool group, 18 percent in the
10 antiarrhythmic drug group, remained free of chronic
11 failure as defined by the protocol.

12 But if you look again at the curves of time
13 to symptomatic afib recurrence, and any observed afib
14 recurrence, there remain very clinically meaningful
15 differences between the two groups. At the end of
16 the follow-up period, 61 percent of the ThermoCool
17 patients and only 28 percent of the antiarrhythmic
18 drug treated patients were free of symptomatic atrial
19 fibrillation recurrence. Similar results for any
20 observed afib recurrence.

21 We also looked at quality of life, and
22 quality of life was significantly improved in the
23 ThermoCool ablation patients. Quality of life was
24 assessed at baseline and at three, six, and nine
25 months of the follow-up period. Both the SF-36 and

1 the atrial fibrillation symptom frequency and
2 severity checklist were used.

3 This slide demonstrates that SF-36 results.
4 This is the mental component on the top and the
5 physical component on the bottom. Baseline values
6 for both groups were similar and below the population
7 norm of 50, and this implies that, in fact, this was
8 a significantly symptomatic group at the onset. The
9 scores were similar for both the ThermoCool and
10 antiarrhythmic drug groups.

11 Those patients in the antiarrhythmic drug
12 group who subsequently underwent ablation were
13 excluded from the then subsequent quality of life
14 analyses. Overall for the SF-36, a three to five
15 unit change is considered clinically significant and
16 meaningful, and as you can see, in both of sets
17 scales, mental and physical, there was a five to
18 eight point difference improvement in quality of life
19 that was maintained over the trial for both scales
20 and again, very little improvement, really change in
21 the quality of life for patients in the
22 antiarrhythmic drug group.

23 Similar results were seen for the symptom
24 checklist outcomes, both in terms of symptom
25 frequency and symptom severity. In this case,

1 there's a decrease in score. That correlates with a
2 decrease in symptoms. There was a greater than 50
3 percent decrease in symptom frequency and severity
4 scores from baseline in the ablation group at all
5 timeframes.

6 However, there was very little change in
7 the antiarrhythmic drug group, the only exception
8 being toward the end of the trial, symptomatic
9 severity in the remaining patients that had not yet
10 undergone ablation was somewhat less, and they think
11 simply reflects the fact that most of the patients
12 with severe symptoms in the antiarrhythmic drug group
13 by that time had elected to have a catheter ablation
14 procedure performed.

15 In terms of the relevance of this trial to
16 heart failure, only New York Heart Association Class
17 1 and 2 subjects were eligible for study inclusion.
18 There were five subjects, three in the ThermoCool arm
19 and two in the antiarrhythmic drug arm, that were
20 enrolled with a history of heart failure at baseline.
21 There were no heart failure related primary adverse
22 events reported in any of the three ThermoCool
23 subjects. However, the safety and effectiveness
24 inference based on these small numbers I think is
25 extremely difficult to make.

1 There are, however, some data outside of
2 the trial, relating to both catheter ablation and the
3 use of the ThermoCool catheter in heart failure
4 patients. The safety of the ThermoCool catheter has
5 been adequately characterized in the VT population by
6 a previous PMA, and in that PMA more than 56 percent
7 of patients that were enrolled and treated with the
8 ablation catheter had heart failure.

9 Finally, there is data not relating
10 specifically to this catheter, that restoration of
11 sinus rhythm by ablation in subjects with heart
12 failure and afib significantly improves cardiac
13 function, symptoms, exercise capacity and quality of
14 life with a low complication rate.

15 We did look at the small subgroup of
16 patients who were enrolled based on the failure of
17 only a class II or IV antiarrhythmic drug. So only
18 16 percent of all enrolled subjects, 20 were in the
19 ThermoCool group, 7 in the antiarrhythmic drug group.
20 These were the protocol defined chronic
21 effectiveness. As you can see here, because these
22 numbers are very small, I think it's difficult to
23 make inferences about the impact of this particular
24 subset of patients on chronic effectiveness.

25 We feel that the results of this trial are

1 generalizable to the U.S. population. There were 15
2 U.S. sites that contributed to the study population.
3 The statistical results of the trial were insensitive
4 to the exclusion of OUS 1 and to discounting of all
5 OUS sites. Analysis of time to symptomatic afib
6 recurrence and time of any observed afib recurrence,
7 demonstrates substantial -- effects in the U.S.
8 population alone.

9 And finally while amiodarone in this study
10 was excluded by protocol, it is considered an
11 unacceptable option by many patients and
12 practitioners for paroxysmal atrial fibrillation due
13 to potential long-term side effects.

14 In this study, electroanatomical mapping
15 was incorporated as part of the ablation procedure.
16 However, alternative mapping guides for AF ablation
17 including fluoroscopy, intra-cardiac,
18 echocardiological, and circulatory mapping catheters,
19 and these have been documented in the literature. We
20 feel that this study does not address whether
21 electroanatomical mapping is superior to these
22 alternative approaches.

23 So we conclude that the superiority for
24 ThermoCool ablation versus antiarrhythmic drug
25 therapy is demonstrated in achieving the primary

1 effectiveness endpoint. This was a randomized
2 control trial. There was a conservative
3 effectiveness endpoint definition that included the
4 addition of other non-antiarrhythmic drugs for
5 example. There was excellent transtelephonic
6 monitoring and compliance and rigorous adjudication
7 of these outcomes. The statistical conclusions were
8 robust to some protocol deviations, and the
9 directionality of treatment was robust across many
10 subsets and many different analyses.

11 We also feel that there were clinical
12 meetings that treatment effects in favor of the
13 ThermoCool arm, in terms of other secondary endpoints
14 including freedom from symptomatic afib or any
15 observed afib recurrence and in quality of life.
16 Thank you.

17 I'd like to now introduce Dr. Al Waldo who
18 will present the results of the safety analysis.

19 DR. WALDO: I'm Al Waldo, a Professor of
20 Medicine at Case Western Reserve University and I'm a
21 consultant for Biosense Webster.

22 So we start now with the primarily safety
23 analysis and the primary safety endpoint for this
24 study was to find as the incidence of early onset,
25 that is within seven days of the ablation procedure,

1 primary adverse events included in this long list
2 which I think you will see in your handout.

3 This looks at the accountability for
4 primary safety analysis, that is all subjects
5 undergoing ablation. So in the ThermoCool group
6 there were 106 patients, 3 of whom were excluded.
7 That meant the overall safety cohort was 103. So the
8 primary safety cohort in this group was 103 patients.
9 In the antiarrhythmic drug group, that is the control
10 group, there were 61 patients, 4 of whom were
11 excluded, leaving 57, and the overall safety cohort.
12 One subject was then discounted, 20 subjects did not
13 undergo ablation. In other words, the remaining
14 patients crossed over to the ablation arm. So there
15 were 36 of them. So we get the final number of 139
16 patients.

17 Now, this looks at the primary safety
18 endpoint. That is the early onset, within 7 days,
19 and you can see that of the 139 patients, there were
20 15 subjects who had 16 events. So the percent was
21 10.8 percent with primary events which made the 95
22 percent upper confidence bound of 16.1 percent which
23 missed by .1 percent of the primary adverse event
24 performance goal of 16 percent.

25 But this breaks down what the nature of the

1 adverse events were. This list is actually a
2 decreasing instance of severity beginning with death,
3 atrioesophageal fistula, et cetera, and you'll notice
4 that there are no events until we get to one
5 pulmonary edema, one pericarditis, seven
6 hospitalizations, one pericardial effusion, and five
7 vascular access complications. So the great
8 majority, that is seven hospitalizations and five
9 vascular access complications made up the most of
10 what this was about. The vascular access
11 complications, of course, are something that can be
12 seen in any kind of catheterization.

13 And this further breaks down the primary
14 safety events for seven days. And of all of these
15 events, we either satisfactorily resolved or
16 improved. So that was pulmonary edema, pericarditis,
17 eight hospitalizations in seven patients, one
18 pericardial effusion and five vascular access
19 complications, and just to again give you some
20 perspective, for instance, the pericardial effusion
21 patient was asymptomatic. It was picked up on the
22 mandated echocardiogram later on, and that was
23 completely resolved.

24 This looks at the hospitalizations in
25 detail. In an extended study, one subject for

1 decrease in hemoglobin level, one subject for
2 hematuria related to traumatic Foley catheter
3 insertion, one subject for atrial flutter. Regarding
4 readmission during the first week, three subjects
5 developed an atrial fibrillation recurrence. One
6 subject for pneumonia and one subject for shortness
7 of breath. These two patients were the same.

8 So again I think you'll find that when you
9 look at these events, they were really a clinically
10 acceptable category.

11 So if we look at the primary adverse events
12 by causality. There were 15 subjects, 16 events,
13 only 1 was possibly device-related. This was a
14 patient who developed heart failure two days out
15 after going home for two days, and this was
16 successfully resolved. There were nine procedure-
17 related, and there were five that were unrelated to
18 device or procedure.

19 Now, very important is to look at what
20 happens in pulmonary vein stenosis. Pulmonary vein
21 stenosis was defined in this study as 70 percent or
22 greater reduction in the diameter of the pulmonary
23 vein compared to baseline. The study cohort included
24 all subjects undergoing an ablation procedure with
25 follow-up CT or MRA. So it included ThermoCool

1 patients and also the antiarrhythmic drug patients
2 who crossed over.

3 This looks at the instance of pulmonary
4 vein stenosis that is greater than 70 percent or
5 greater from the baseline in any targeted vein per
6 subject. You notice there were none. So pulmonary
7 vein stenosis at three months, none. In the
8 ThermoCool group, pulmonary vein stenosis at 12
9 months, there were none, and the same for those who
10 crossed over, none at 3 months and none at 12 months.
11 So of the 82 patients at 3 months and no pulmonary
12 vein stenosis clinically, and of the 29 patients at
13 12 month period with no pulmonary vein stenosis.

14 This slides look at the pulmonary vein
15 diameter changes at these 3 and 12 month by pulmonary
16 vein. So we're looking here at the total number of
17 pulmonary veins that were examined, and you can see
18 that some of them decreased and some of them
19 increased. The histograms in yellow are those that
20 were looked at 3 months into the bio CT or MRA, and
21 the green were those who were looked at in 12 months.
22 So you can see most of the patients were in the mild
23 category with a decrease, in fact, the vast majority
24 were either 20 percent to less. The remainder were
25 the 50 percent at 3 months. At 12 months, there was

1 one patient whose diameter had decreased by 51.9
2 percent. And interestingly enough, some of these
3 veins increased in size, again some 20 percent and
4 some a little more than 20 percent. And really what
5 I think this reflects is two things. Just look at
6 this. It looks like a normal bell-shaped curve
7 suggesting normal variability. The other thing is
8 that when you make atria remodel and you shorten,
9 which will decrease some of the diameter, and also
10 there's variability from time to time, the diameter
11 of these veins are physiologic dependent depending on
12 your HO volume and what stage during filling or
13 contraction, et cetera. So I think basically the
14 important thing is to look at this. There were no
15 patients who had over seventy percent, and not even
16 close.

17 Now, this looks at the secondary safety
18 analysis. So Biosense Webster developed a
19 hierarchical classification of adverse events based
20 on two categories of level of severity. Category one
21 was adverse events that resulted in permanent injury
22 or impairment, including death, cerebral vascular
23 accident, myocardial infarction, pulmonary vein
24 stenosis, diaphragmatic paralysis, or atriopharyngeal
25 fistula. Category two was a long list of temporary

1 or reversible causes.

2 So when we look at this thing, if we look
3 at the secondary safety endpoint, that is early onset
4 of serious adverse events within 90 days of initial
5 treatment, there were no category one events in
6 either the ThermoCool group or in the antiarrhythmic
7 drug group. There were 19 events of category two in
8 the ThermoCool group which was 18.4 percent of the
9 cohort, and 20 events or 35.1 percent of the cohort
10 in the antiarrhythmic drug group, so it was basically
11 twice as many in the antiarrhythmic drug group or if
12 you will, half as many in the ThermoCool group with a
13 P value of .0221.

14 Now, here we look at the secondary safety
15 analysis of the late onset of serious adverse events,
16 that is beyond 90 days of initial treatment, and here
17 in the ThermoCool group there was one death. This
18 death occurred 284 days after the ablation. It was a
19 gentleman who had a known history of coronary artery
20 disease. He had chest pain at night, went to bed
21 with the chest pain, and unhappily never woke up.

22 Then if we look at the category two, there
23 were 7.8 percent in the ThermoCool group with 14
24 percent in the antiarrhythmic group and the category
25 other that is these are some events that were not on

1 that list that I just showed you before. So we
2 looked at the total. The incidence was 10.7 percent
3 with the ThermoCool and 15.8 percent of the
4 antiarrhythmic drug group. Again, the ThermoCool
5 group did better considerably.

6 So we conclude then that there is an
7 excellent safety profile from the ThermoCool atrial
8 fibrillation catheter ablation patients. The primary
9 adverse event incidence, the performance goal was 16
10 percent. The observed missed it by .1 percent.
11 There was one possibly device-related event. There
12 were no deaths, myocardial infarctions, strokes,
13 cerebral vascular accidents, heart block, atrial
14 perforation or the like within seven days. There is
15 no clinically significant pulmonary vein stenosis.

16 Regarding the early onset serious adverse
17 events, there was a lower incidence in the ThermoCool
18 group of 18.4 percent compared to the antiarrhythmic
19 drug group of 35.1 percent. Late onset serious
20 adverse events, again a lower incidence in the
21 ThermoCool group of 10.7 percent compared to the
22 antiarrhythmic drug group of 15.8 percent.

23 So now I'm pleased to turn the podium back
24 to Dr. Yaross.

25 DR. YAROSS: Thank you, Dr. Waldo. To wrap

1 up, we have presented valid scientific evidence as
2 defined in FDA regulations for your consideration
3 today.

4 The NaviStar AF study had randomized
5 control design. The Bayesian analytical approach was
6 used to permit efficient study completion following
7 the challenges to enrollment that I described to you
8 earlier.

9 The ThermoCool AF study was rigorously
10 conducted. We had excellent TTM adjudication, and it
11 was a thoroughly vetted dataset.

12 The sponsor continuously monitored the
13 sites throughout the trial, and in addition, we had
14 underwent FDA fire research monitoring audits of the
15 two highest enrolling sites and at the sponsor.

16 There were no Form 483 inspectional
17 observations at those sites, and thus this
18 corroborates the validity of the dataset.

19 This therefore likely represents the most
20 rigorously and thoroughly vetted AF ablation dataset
21 to date.

22 The results therefore demonstrate the
23 safety and effectiveness of the ThermoCool catheter
24 for the treatment of AF. The primary trial
25 objectives I outlined earlier were met. We have

1 shown the superior chronic effectiveness of the
2 ThermoCool catheter ablation versus antiarrhythmic
3 drug treatment by a rigorous protocol definition.

4 While we acknowledge on the safety side
5 that we missed the other confidence limit at the
6 safety endpoint by 0.1 percent, based on the lack of
7 severe events, we also submit that this provides an
8 acceptable safety profile as Dr. Waldo has outlined.

9 We have also shown to you additional
10 important results for the ThermoCool ablation
11 subjects. They, versus the AAD control subjects,
12 were more likely to be free of any observed atrial
13 fibrillation occurrence. They had improved quality
14 of life and have fewer severe side effects. In fact,
15 the secondary, the prespecified analysis that we
16 provided, we had hypothesized non-inferiority for
17 that analysis and, in fact, the results show
18 superiority.

19 We'd also like to point out that both
20 genders were well represented in the ThermoCool AF
21 trial. Women represented one-third of the
22 population, and logistic analysis showed that gender
23 was not a predictor of product success outcome or
24 primary adverse events in this study. We therefore
25 conclude that the product is effective and safe in

1 both men and women.

2 Biosense Webster believes that the results
3 of this trial should be incorporated with the
4 ThermoCool catheter instructions for use. Public
5 interest is best served by rapidly communicating the
6 results of this trial so that information on AF
7 ablation risks and benefits are communicated in an AF
8 approved package insert with the device label.

9 Currently we, as a company, are unable to
10 train physicians specifically on atrial fibrillation
11 treatment. Given the growing performance of AF
12 ablation off-label in the community, it is in the
13 public interest for us to be able to conduct formal
14 training on safe and efficient use of this catheter
15 in the AF population.

16 Biosense Webster is committed to a formal
17 training program. We currently require clinical
18 training prior to first shipment of the catheter to a
19 hospital, and we will commit to training on atrial
20 fibrillation as well as part of this program post-
21 approval.

22 As I showed you earlier, the current
23 approved indications for use are Type 1 atrial
24 flutter and recurrent drug and device refractory VT
25 in the post-MI population.

1 Our proposed change to this indication
2 produced is simply to add drug refractory,
3 symptomatic paroxysmal atrial fibrillation. Draft
4 information for use is in your handout.

5 We have also agreed to conduct a post-
6 approval study to confirm that the results of the
7 study may be generalized in a postmarket setting. In
8 the study, we also propose to look at long-term
9 safety and effectiveness.

10 Our proposal currently under review by the
11 FDA includes 5 year follow-up and at least 50 percent
12 more investigational sites within the United States
13 to demonstrate that atrial fibrillation ablation with
14 the ThermoCool catheter is safe and effective in
15 other sites in addition to those in the clinical
16 trial.

17 We look forward to the Panel's input on the
18 design of this study later today.

19 In closing, we conclude that we have met
20 the statutory burden for premarket approval. We have
21 provided a clinical study that meets the highest tier
22 of the criteria for valid scientific evidence as
23 outlined in Federal Regulations, and this is a
24 randomized control trial.

25 The probable benefits from AF ablation with

1 the ThermoCool catheter have been proven. The
2 probable risks when used as directed are clinically
3 acceptable in symptomatic paroxysmal AF population.

4 Biosense Webster therefore respectfully
5 requests that the Panel recommend this application
6 for approval.

7 With that, I thank the Panel for its
8 attention and the FDA for a very -- review.

9 DR. BORER: Thank you very much, Dr. Yaross
10 and colleagues, for a very effective and efficient
11 presentation.

12 I'm going to ask the Panel if anyone on the
13 Panel has any questions, but I'd like to set some
14 ground rules if I may.

15 At this point, I'd like to limit the
16 questions to clarifications of the data. We're going
17 to have the opportunity to ask questions of the Panel
18 again later, and we're going to hear a FDA
19 presentation, and I think that some of the questions
20 that we might be asking going beyond clarification of
21 the data might best be handled when we discuss the
22 questions after you've heard everything.

23 In addition, we've heard a wonderful
24 presentation by Dr. Berry about the statistics, and
25 I'm sure that for us as cardiovascular investigators

1 and people focused on cardiovascular diseases in some
2 but not all cases sitting around the table, the
3 Bayesian approach is not one that we're very much
4 accustomed to in clinical trials, although it's
5 commonly used in oncologic trials. We're going to
6 hear from Laura Thompson of the FDA about this as
7 well, and I would suggest that we hold the questions
8 about statistics until after we've heard the FDA
9 presentation, and then Dr. Berry will be welcome to
10 come up and answer some of the questions as well,
11 even though it doesn't say that on the program.

12 And finally, please wait until the Chair
13 recognizes you before asking a question so that we
14 can maintain some order. I promise, we'll sit here
15 all night if we have to, to get all the questions
16 answered that anyone may have.

17 Having said all those things, I see
18 Dr. Somberg's hand up already. So, Dr. Somberg.

19 DR. SOMBERG: Well, number one, I promise
20 you we won't have to be here all night because I'm
21 not going to be here.

22 Number two, congratulations to the company.
23 It's a very nice study, and it certainly sets an
24 example for many other device companies.

25 I have a couple of very small questions

1 about the data, that you present so much of it,
2 therefore there's always some little tidbits here and
3 there that I don't understand, and maybe you've said
4 this.

5 But number one, the people who received the
6 ThermoCool versus the antiarrhythmic drugs, the
7 patients in the ThermoCool group were continued on
8 drugs that were previously considered ineffective.
9 Then it was said that at the end of the follow-up
10 period, most of these people were off of it. Can you
11 be a little bit more specific? How many people were
12 off of it? And what is the -- and what I'm trying to
13 drive at is, is this an ablation without drugs at the
14 end or is it an ablation that the clinician should be
15 encouraged to continue antiarrhythmic therapy? I
16 think that's an important implication. So if you
17 could be a little bit more specific, I'd appreciate
18 it.

19 DR. YAROSS: I'll be happy to -- actually,
20 I'll ask Dr. Wilber to respond to your question.

21 DR. WILBER: In general, when patients were
22 allowed to have previously ineffective antiarrhythmic
23 drugs, they were most often stopped at the end of the
24 blanking period. However, they were allowed to
25 continue by protocol through the chronic efficacy

1 period. However, this was largely at the first three
2 months of the chronic efficacy evaluation period, and
3 at six months and nine months of the chronic efficacy
4 period, only seven percent of those successes were on
5 antiarrhythmic drugs. So the protocol did not
6 specify when they were stopped, but by the last two
7 follow-ups at six and nine months, only seven
8 percent.

9 DR. SOMBERG: That's very helpful. Maybe
10 you should stay there for a minute because I think
11 you're going to be asked to stand up again unless you
12 want the exercise and all that.

13 There were patients who were initially or
14 discontinued or considered inappropriate because you
15 said beta blockers, calcium channel blocker, ACE
16 inhibitors were added to that. Is there any
17 information in the materials that was added because
18 of arrhythmia or instability, et cetera, or was that
19 for other things like blood pressure and what have
20 you?

21 And the other question that sort of ties in
22 with that is that there's been recently a lot of
23 reports with statins being useful. Is there any
24 imbalance in the statin use between the two
25 populations?

1 DR. WILBER: With respect to your first
2 question, investigators were asked as to the reason
3 why a patient was prescribed with a specific drug if
4 it was changed. There were times when a specific
5 classification wasn't made. So the assumption was
6 always by default that it was then used potentially
7 as a antiarrhythmic drug. So the protocol by default
8 took a very strict definition of changes in drug
9 therapy, and as we pointed out, this tended to be
10 more typical of U.S. sites than non-U.S. sites.

11 With regard to your second question about
12 the statins, I don't have that information at this
13 time.

14 DR. SOMBERG: And my last quick question, I
15 guess to you, Dr. Wilber, there were 14 subjects I
16 understand that were not, and this doesn't seem --
17 when the FDA's review, which we'll hear and then your
18 presentation, which didn't seem to fit together the
19 number of people who were excluded or lost to follow-
20 up for some reason or not in the group. Was there 14
21 of the people, I don't know the exact number offhand,
22 but of the people on the ThermoCool study, was there
23 that number lost to follow-up?

24 DR. WILBER: It's not lost to follow-up.
25 Let's bring up the Kaplan-Meier curve if we can for

1 the first chronic efficacy evaluation by protocol.

2 There were 14 patients that as of June had
3 not yet completed their follow-up. So they weren't
4 lost to follow-up. They were followed until the time
5 of -- that's when the dataset was sealed or the last
6 update. So we have no further data. They were
7 censored at that time, and you can actually see very
8 specifically when those times were, and that's what
9 these circles are. It represents those -- these are
10 the times the patients were followed up until.

11 So as with any other Kaplan-Meier curve,
12 they're censored when follow-up is done. So the
13 follow-up is not complete in 14 of those patients,
14 and this data represents what was available when we
15 submitted the data in June.

16 DR. SOMBERG: Can you comment on, you know,
17 updating us on that? I mean that is a number of
18 people.

19 DR. ZUCKERMAN: Dr. Somberg, perhaps as
20 Dr. Borer said, Dr. Thompson will get into the
21 predictive modeling of how you account for those 14
22 patients, and then the question and answer will be
23 richer.

24 DR. SOMBERG: Okay.

25 DR. ZUCKERMAN: So we'll cover that.

1 DR. SOMBERG: I'll ask the question again
2 later. That's what I've been told.

3 DR. BORER: Okay. Dr. Kelley.

4 DR. KELLEY: I have two questions for
5 Dr. Wilber. One's simple, and I think I know the
6 answer but just to clarify. On slide 61, you talked
7 about freedom from symptomatic AF recurrence, and
8 then any AF recurrence, and I'm assuming that's any
9 documented on the --

10 DR. WILBER: Any observed afib recurrence
11 using the scheduled transmissions during which the
12 patient was not necessarily, in fact, usually was not
13 symptomatic.

14 DR. KELLEY: Okay. And then the second
15 question is the rationale where to -- for the
16 inclusion in the study that patients had to have had
17 three episodes in six months, only one of which was
18 documented. The others it seems like less than 30
19 seconds of symptoms would be counted in the absence
20 of documentation. But afterwards, it had to be
21 documented. What was the rationale for that
22 difference?

23 DR. WILBER: I think that the protocol
24 which is there needed to be electrocardiographic
25 documentation of the episode in at least one of the

1 three episodes that they had, part of that is simply
2 the issues about enrollment. If you, and for folks
3 who see atrial fibrillation patients, it's often very
4 hard to get that documentation. This is the protocol
5 as it was adjudicated when it was worked out with the
6 FDA, that at least one of them had to be documented
7 by atrial fibrillation. Otherwise, the collection of
8 data would have been much more difficult in terms of
9 finding patients suitable for enrollment, and I think
10 that was one of our barriers.

11 DR. KELLEY: So post-ablation, suppose they
12 had the same symptoms 30 seconds of what was called
13 afib before but didn't have their monitor with them.
14 It was still counted as atrial fibrillation or no?

15 DR. WILBER: If they had symptoms without
16 an electrocardiographic definition, they would not --
17 there would be no way to count it as afib because we
18 wouldn't know that, although I think much of that
19 information would be picked up in the quality of life
20 if there were patients who had continued symptomatic
21 afib but didn't use their monitor during that
22 recording, but it wouldn't capture all.

23 DR. KELLEY: Okay. Thank you.

24 DR. BORER: Dr. Slotwiner.

25 DR. SLOTWINER: Okay. Thanks. I just

1 wanted to clarify the catheter used. The ThermoCool
2 technology from what I've seen is incorporated into
3 five different catheters that are manufactured,
4 unidirectional, bidirectional, with and without nav
5 and the remote magnetic steering. But all these
6 patients who underwent ablation had the
7 unidirectional navigation capable catheter used?

8 DR. YAROSS: That's correct. The clinical
9 trial was conducted using the NaviStar ThermoCool
10 catheter, the original, not the bidirectional model.

11 DR. SLOTWINER: But the request for
12 approval is for all five of the flavors, correct?

13 DR. YAROSS: That's correct.

14 DR. BORER: Dr. Jeevanandam.

15 DR. JEEVANANDAM: I have a simple question.
16 Your transtelephonic monitoring, was that via an
17 implanted device that was capturing everything and
18 then downloaded, or was this something that the
19 patient wore just when they were transmitting or when
20 they were symptomatic?

21 DR. YAROSS: It's an external device.
22 Dr. Wilber, do you want to comment any further? It
23 was an external device that they applied, you know,
24 that they used at the time of transmission.

25 DR. JEEVANANDAM: So if somebody had afib

1 that wasn't asymptomatic a week before they
2 transmitted, it would not be picked up.

3 DR. YAROSS: They transmitted weekly during
4 the first two months of the effectiveness evaluation
5 period, monthly thereafter, but they did scheduled
6 transmissions plus were instructed to transmit at any
7 time they had cardiac symptoms.

8 DR. JEEVANANDAM: So within that week, if
9 they didn't have symptoms but they were in afib, you
10 wouldn't pick it up. You'd only pick it up when they
11 actually transmitted?

12 DR. YAROSS: That's correct. The
13 prospective endpoint was symptomatic paroxysmal
14 recurrence.

15 DR. BORER: I have a series of questions
16 but will interrupt them as people come up with
17 others. First of all, I didn't understand from
18 protocol, the presentation, how anticoagulants were
19 managed. It appeared that people who got the
20 ThermoCool therapy were on anticoagulation generally
21 for a period of a couple of months afterwards. It
22 wasn't clear what the AAD, the people on AAD were
23 getting, and that might have some impact on the AEs
24 that were or weren't recorded. There may be some
25 difference of opinion about how anticoagulants should

1 be used, but I'd like to know how they were used.

2 DR. YAROSS: I'll ask Dr. Wilber to respond
3 to that question.

4 DR. WILBER: So all anticoagulation therapy
5 was recommended in the ablation group for a minimum
6 of the three months after the procedure. If patients
7 had continuing atrial fibrillation, it would be
8 continued thereafter. If they did not have atrial
9 fibrillation in follow-up, it was possible and
10 generally with a low CHADS score to allow
11 discontinuation of anticoagulation in a limited
12 number of patients. Immediately after the ablation,
13 patients were treated with anticoagulation obviously
14 until the NR reached two or three.

15 For the control group, it was recommended
16 that patients be on oral anticoagulation in general
17 throughout the course of the study.

18 DR. BORER: Okay. But some were and some
19 weren't, I assume. Was there any relation to outcome
20 of anticoagulation or it just didn't show anything?

21 DR. WILBER: We do not have any data to
22 suggest that any of the outcome were related to, and
23 since there were, in the sense of thromboembolic
24 events, it was a small study. There weren't
25 thromboembolic events --

1 DR. BORER: Or bleeding.

2 DR. WILBER: -- or significant bleeding.

3 DR. BORER: Judah.

4 DR. WEINBERGER: Yeah, I just want to ask
5 at this point, there seems to be this issue about the
6 intensity of AF monitoring. Since you're only
7 looking for symptomatic AF, and you're not
8 particularly interested in asymptomatic AF, I mean
9 that's not a piece of data that you capture, unless
10 it happens to be there during the transtelephonic
11 transmission. Nevertheless, you allow
12 discontinuation of anticoagulation if there's no
13 symptomatic AF. So were there any Holters that were
14 mandated during this time?

15 DR. WILBER: Just to correct one statement,
16 in general, in accordance with the current guidelines
17 we would not recommend discontinuation of
18 antiarrhythmic therapy unless patients have a low
19 CHADS score.

20 DR. WEINBERGER: I said anticoagulation.

21 DR. WILBER: Anticoagulation, I'm sorry.
22 So that presumably to some extent is independent of
23 the occurrence of atrial fibrillation. Your second
24 question.

25 DR. WEINBERGER: They both tied into

1 intensity of monitoring.

2 DR. WILBER: And I think the other question
3 is that it's always problematic to decide how much
4 monitoring you need to do to pick up asymptomatic
5 atrial fibrillation and there's probably -- and there
6 are some data to suggest that if you monitored them
7 forever during the entire trial, you would obviously
8 pick up more than taking 15 isolated times, and at
9 some point you have to make a decision about what's
10 practical. For this study, the decision was 15 time
11 points. I would agree with you that that certainly
12 doesn't capture all potentially asymptomatic
13 episodes.

14 DR. BORER: Dr. Kelley.

15 DR. KELLEY: Still for Dr. Wilber. This is
16 about the antiarrhythmic drug issue. Certainly the
17 difference between OUS 1 and the other sites is
18 pretty dramatic and one of the theories is that those
19 patients were on -- drugs. Also if you look at the
20 difference between the results of the patients who
21 had just failed class II and IV drugs versus the
22 other patients, those patients had more afib, and one
23 wonders if it's because by protocol they couldn't be
24 on another antiarrhythmic drug. So I wondered if
25 there were any analyses looking at results by

1 patients on antiarrhythmic drugs versus not.

2 DR. WILBER: Is the question that you're
3 asking about the chronic effectiveness evaluation.

4 DR. KELLEY: For chronic efficacy.

5 DR. WILBER: Yeah. For patients who were
6 on anti -- again, because there were so few patients
7 on antiarrhythmic drug therapy that were successes at
8 the end of the trial that, in fact, most patients who
9 were on any antiarrhythmic drugs that had not failed,
10 in other words, had not had a recurrence, were really
11 only on them for the first three or four months of
12 the chronic efficacy evaluation.

13 DR. KELLEY: But not at OUS 1.

14 DR. WILBER: Even at OUS 1, yes, that
15 includes OUS 1. So in general, U.S. sites tended to
16 have no antiarrhythmic therapy at all, even a failed
17 drug, whereas OUS 1 did use the protocol that allowed
18 one for the first three or four months of the chronic
19 efficacy period, and then tended to stop it. And
20 there are only seven percent of those for the six and
21 nine-month follow-up.

22 So in a preliminary analysis of just
23 looking at who was and who wasn't, there's no impact
24 on outcome. We did not sort of do some time
25 dependent modeling to account for the fact that the

1 drugs were withdrawn during that, and that is an
2 analysis that could be done.

3 DR. KELLEY: In other words, did we see the
4 failure when those drugs were withdrawn or --

5 DR. WILBER: No, no. In other words, they
6 didn't suddenly pop up a month after that drug was
7 stopped. So the ultimate impact of doing that I
8 think is still a question that's unknown, and is it
9 about remodeling? Is the remodeling more than three
10 months? Is the remodeling less than three months?
11 As you know, clinically there's still a debate about
12 that in terms of when does the benefit of healing
13 end, and so I'm not sure that you can say much of
14 anything except there's no obvious signal that the
15 continued use of a prior failure antiarrhythmic drug
16 had an impact on chronic success.

17 DR. KELLEY: Okay. Thank you.

18 DR. BORER: Let me follow-up on that with a
19 question about the amiodarone. I understand the
20 protocol, but you had one patient, I believe that it
21 was a patient in OUS 1, who had an initial ablation,
22 a recurrence, was given intravenous amiodarone,
23 converted, I think it was intravenous, converted, and
24 then had a second ablation about two months later
25 which was within the protocol and then was fine. It

1 sounds as if there may be some inconsistency in the
2 application of the amiodarone rule, or maybe I don't
3 understand the amiodarone rule. Was this patient
4 continued in this study after the amiodarone was
5 given? How was that handled? And were other
6 patients given amiodarone for recurrences?

7 DR. WILBER: That was the only patient that
8 received intravenous amiodarone. Because it was an
9 intravenous administration, was short-lived, and
10 there are quite a bit of data that much of the
11 efficacy of acute amiodarone intravenously may have
12 more to do with not as class III but its class I
13 effects, and it's unlikely that that single dose,
14 that it persisted two months later to have an impact
15 on the trial, so that that patient was continued in
16 the group and that it was not considered a failure
17 because this was all in the blanking period. That
18 patient did not receive chronic amiodarone therapy,
19 and there are no patients that for purposes of
20 chronic success had received amiodarone. There was
21 some amiodarone use after patients had recurrences
22 and had met the failure criteria.

23 DR. BORER: Yes. Dr. Fleming.

24 DR. FLEMING: I wanted to follow-up on
25 Dr. Somberg's question. What's not clear to me from

1 reading this is whether the established protocol
2 post-ablation is to put all your patients on
3 antiarrhythmic drugs. Some centers do, some don't.
4 So I think it's a legitimate question from a
5 consumer's point of view, and I would also tell the
6 Panel that I have PAF. So I'm very interested in not
7 taking those drugs. And as a matter of fact, I would
8 not take them if offered. So what I'm asking you is
9 what is the general sense of antiarrhythmic drugs. I
10 think that was part of what Dr. Somberg was asking
11 earlier, and it's a legitimate question.

12 DR. WILBER: It's absolutely a legitimate
13 question. I'm not sure that the data from this study
14 can address that in that they're not, similar to the
15 specific ablation set that was done, there wasn't a
16 mandated criteria that you either must do this or
17 must do that for antiarrhythmic drugs. So leeway
18 allowed the investigator, and as one might expect,
19 there were differences in practices across group.
20 I'm not sure that even in the larger world of afib
21 ablation there's currently agreement among all
22 investigators on how this is handled, and so I wish I
23 would provide that result. We can't do it from this
24 study, but one thing we did try to assure is that the
25 use of those drugs didn't somehow have an impact on

1 the endpoint, which for that I'm satisfied they
2 don't, but whether or not it's a benefit otherwise,
3 unfortunately and where there are solid data, people
4 feel very passionately on different sides of the
5 issue about whether they should or shouldn't have it
6 after three months. Probably the most common
7 practice is that the drugs are stopped within three
8 months, certainly within the U.S. That is not
9 necessarily a European practice.

10 DR. BORER: John.

11 DR. SOMBERG: Yeah, this is to Dr. Waldo.
12 I'm going to give you a break for a minute, but I
13 didn't want to leave Professor Waldo who has come
14 this far without a question. I think it was very
15 nice to break it up into category 1 and category 2
16 safety actions, and it's really testimony to the
17 device that it didn't have category 1 or the drugs
18 for that matter.

19 But, you know, then you went on to compare
20 category 2 effects between the drug and the
21 ThermoCool, and it favors the ThermoCool. But isn't
22 one of the criteria prolonged QT and the drugs
23 prolong the QT. So the type 3 drugs are acting as
24 part of the efficacy by -- in the QT, but that's
25 considered a safety problem. Did you -- I mean what

1 is prolonged QT? Did it have to be more than 5 -- 40
2 milliseconds or something? And if you leave out the
3 QT, was that the driver of the difference between
4 ThermoCool and the drug?

5 DR. WALDO: No, it was not. First, most of
6 the patients were on class Ic drugs. So the class
7 III drugs were really not very, very important in the
8 trial. None of the patients had a prolonged QT
9 interval. There's no -- or nothing even close. So
10 it really wasn't an issue.

11 DR. SOMBERG: Okay. And then the other
12 question I would have is do you have any comment on
13 why most, you know, it just struck me as interesting
14 that they would use mostly Ic's, like -- and do you
15 think that is routine practice not to use the sotalol
16 and dofetilide?

17 DR. WALDO: Well, dofetilide, as you know
18 does not have an indication for paroxysmal atrial
19 fibrillation, and as you also probably know, Pfizer
20 did six different trials trying to see if dofetilide
21 would be effective in paroxysmal atrial fib, but it
22 didn't show that. There's not a lot of data about
23 sotalol in paroxysmal atrial fibrillation, but the
24 class Ic's have that indication. Not only that,
25 remember the patient population, the mean age was

1 about 55, 56, and young patients and they didn't have
2 underlying structural heart disease for the most
3 part. So I think it was very appropriate, very
4 understandable as far as I could see.

5 DR. BORER: Dr. Karasik.

6 DR. KARASIK: This is question is probably
7 best for Dr. Wilber, but perhaps you can address it.
8 So we weren't allowed to use amiodarone. You weren't
9 allowed to use amiodarone in this particular
10 protocol, but I'm wondering what percentage of the
11 patients had actually been exposed to amiodarone
12 prior to the six-month exclusionary procedure? If
13 patients had failed that drug in the past, it would
14 suggest that it was unlikely that any other drug
15 would be effective. So was there potentially a bias
16 there?

17 DR. YAROSS: The protocol was only specific
18 to the six months prior to enrollment. I'll have our
19 team check and see if we have that information and
20 provide that to you a little bit later.

21 DR. BORER: Okay. A couple of more -- oh,
22 Dr. Weinberger.

23 DR. WEINBERGER: This is for Dr. Wilber.
24 The OUS site number 1 did a lot better than everybody
25 else as you pointed out, and you presented several

1 hypotheses as to why that could be, including
2 previous experience with a catheter or previous AF
3 ablation experience. And I was wondering whether or
4 not you tested those hypotheses against the rest of
5 the group. In other words, did people improve, did
6 outcomes improve as more experience was garnered with
7 the NaviStar catheter?

8 DR. WILBER: It's an excellent question,
9 and it was an analysis we hoped to do, and the reason
10 why we couldn't do it is that for most of those
11 centers, there weren't enough volume that we could
12 rely on their self-reported afib volume. The problem
13 is since it's not verifiable and you'd have no hard
14 data, I'm not sure of the value of the observation.
15 So the question you raise is an important one, but
16 within the context of the study, we couldn't answer
17 it.

18 DR. WEINBERGER: And then one thing I'll
19 ask you, Dr. Waldo listed as possible hypothesis for
20 why there was some pulmonary vein shrinkage as
21 improvement and shrinkage of the left atrial size,
22 and I wondered whether, although this wasn't
23 presented, whether the echocardiographic information
24 in fact showed return of left atrial function in
25 these patients. Did they all have LA function and

1 shrinkage of the LA size?

2 DR. WILBER: We have not yet finalized the
3 analysis of that data. So I don't have that
4 specifically at this time.

5 DR. YAROSS: I can go back and answer the
6 previous question on amiodarone use. In the
7 ThermoCool group, there were seven subjects who had
8 had prior exposure to amiodarone, and in the control
9 group there were six subjects. So 6.7 percent and
10 10.2 percent, and there was no significant difference
11 between those two values.

12 DR. KARASIK: Thank you.

13 DR. BORER: I'd like to ask something about
14 the protocol itself in the pre-randomization period
15 and actually after randomization. You allowed a 14-
16 day dosing interval in the AAD group. At the end of
17 the day, the results are what they are, and I have no
18 major concern about this, but is it reasonable to
19 allow such a relatively short period to add multiple
20 drugs together and titrate them, et cetera, et
21 cetera? I mean what was the basis for selecting 14
22 days?

23 DR. YAROSS: I'll have Dr. Wilber return to
24 the podium.

25 DR. WILBER: The protocol was designed as a

1 single drug. So a possible alternative might have
2 been to have had multiple drugs in which case a
3 longer window might have been appropriate, but the
4 design of the study was to compare ablation to a
5 single antiarrhythmic drug, and if that drug didn't
6 work out in the two week dosing period, then it was a
7 failure.

8 DR. BORER: Okay. As a follow-up to that,
9 everyone had to have three episodes of atrial
10 fibrillation prior to entry, and I understand that
11 you were not capturing data on these because you
12 couldn't. It would be retrospective. However, maybe
13 you have some information.

14 Do you have any idea how long these
15 episodes lasted and how they were aborted? This
16 would be anecdotal, but I'd be interested to know.

17 DR. WILBER: We did not collect specific
18 data on the duration of those episodes. So I could
19 only sort of extrapolate from our own experience in
20 the trial in those patients, but typically since the
21 protocol definition would allow up to 30 days, the
22 vast majority of patients had symptoms from hours to
23 a day or two in duration and very rarely beyond that.

24 DR. BORER: And how were these episodes
25 aborted generally?

1 DR. WILBER: Generally spontaneously. In
2 fact, it would be rare. There was an occasional
3 patient that we all see that after two hours of
4 atrial fibrillation comes to the emergency room and
5 insists on being cardioverted, but in general, most
6 of these patients, since they've been long-time
7 atrial fibrillation patients, know if they wait it
8 out, their episode will stop and typically don't come
9 to the emergency room.

10 DR. BORER: Okay. Also just for
11 clarification, Dr. Waldo, you referred to the
12 pulmonary edema, the one patient who had pulmonary
13 edema, as having had a device-related complication.

14 DR. WALDO: Possibly device-related.

15 DR. BORER: My recollection of that case is
16 that the three and a half liters of fluid were given
17 to the patient during the procedure. I would have
18 thought that one would consider that procedure-
19 related rather than device-related, but can you just,
20 how would it? Why was it classified as device-
21 related?

22 DR. WALDO: Well, it was possibly device-
23 related because it happened so late. The patient was
24 discharged, don't forget, and came back, but that's
25 the reason we called it possibly device-related

1 because the ThermoCool injects saline as part of its
2 mechanism, and that could be an issue.

3 DR. BORER: So the injection of saline is
4 automatic. It's not --

5 DR. WALDO: It's part and parcel of the
6 procedure. And David could talk about it better than
7 I could.

8 DR. BORER: But is the saline injected
9 automatically or is it a decision by the operator
10 to --

11 DR. WALDO: No, the amount of saline
12 injected can be controlled, and there are
13 recommendations for that in the protocol.

14 DR. BORER: Okay. It's not a big deal
15 obviously. I just wanted to clarify it for myself.

16 DR. WILBER: The flow rate of the saline
17 from the tip is dependent on the power, and there's
18 certain recommendations. If the power is 30 watts or
19 lower, it's 17 ccs per minute, and if it's greater
20 than 30 watts, it's generally 33 ccs per minute. So
21 that can be quite a bit of fluid over the course of a
22 long ablation.

23 MS. YAROSS: Dr. Borer, if I can, I'd also
24 like to just go back to a previous question. We were
25 asked about the duration of the AF episodes. There

1 was a baseline Holter, and in the ThermoCool group,
2 the mean AF duration in a 24-hour Holter was 8.3
3 hours, and this was versus the AAD group, which was
4 10.9. So we did have patients with significant
5 duration.

6 DR. BORER: Dr. Somberg.

7 DR. SOMBERG: The blanking period was 90
8 days for the ThermoCool, but you only have 14 days
9 for the drug. Was an analysis looked at if you -- if
10 you looked at 90 days for both was -- if the drug had
11 prolonged time for action and then took that period
12 and, you know, made 90 days for the drug comparable
13 to 90 days for that, would there be any change? You
14 didn't do that analysis. Okay.

15 DR. WILBER: That specific analysis was not
16 done, although because all of the drugs that we used
17 had a half-life of hours, as opposed to days, at
18 least from the standpoint of achieving adequate
19 efficacy, the 14-day window should have been adequate.
20 Obviously if amiodarone had been one of the drugs,
21 then the trial design might have needed to be
22 different than it was.

23 DR. BORER: Dr. Weinberger.

24 DR. WEINBERGER: I'm sorry to go back to
25 OUS 1. It's just such a striking outcome. There was

1 reference made by somebody that maybe the patients
2 were a little less -- had somewhat smaller atria,
3 maybe the patients were a little bit younger, maybe
4 there was a little bit less hypertension, and I'm
5 wondering whether or not in that particular group of
6 patients who are younger and healthier, AF is
7 frequently triggered by another supraventricular
8 tachycardia. And the ablation of a bypass tract
9 might be sufficient to inhibit subsequent events of
10 atrial fibrillation. So my question was whether
11 specifically at other sites, there was a search for
12 other initiators of atrial fibrillation at the time
13 of electrophysiological study?

14 DR. WILBER: All patients had atrial pacing
15 as part of their evaluation procedure, and so it
16 would be -- they have no data from any of the sites
17 that we did collect data on other arrhythmias that
18 the patients had and whether any other arrhythmias
19 were ablated. None of them had, other than flutter,
20 there weren't other arrhythmias documented other than
21 triggering afib is another possible one or atrial --

22 DR. BORER: I'd like to ask again only for
23 clarification because this is, you know, it's not
24 nitpicking. I just want to know the answer. You
25 used a SF-36 to assess quality of life. I'm a strong

1 proponent of assessing disease burden as QOL. So I
2 like that you did that.

3 However, has the SF-36 ever been formally
4 validated for arrhythmias?

5 DR. YAROSS: I'd like to ask Dr. Matthew
6 Reynolds to come and address this issue please.

7 DR. REYNOLDS: Thanks. I'm Dr. Matthew
8 Reynolds, an electrophysiologist from Boston, the
9 Director of the Economics and Quality of Life
10 Assessment group with the Harvard Clinical Research
11 Institute and a consultant to the company.

12 Thanks for your question. To my knowledge,
13 the SF-36 instrument, which, as you know, is a
14 generic quality of life tool, has not been
15 specifically validated in terms of a validation study
16 in AF population. It has, however, very commonly
17 been used in previous AF studies, including NIH
18 funded studies and other antiarrhythmic drug studies.

19 DR. BORER: Okay.

20 DR. REYNOLDS: So there's a good body of
21 experience with it in AF populations. In terms of a
22 specific validation, the answer is no.

23 DR. BORER: And that's probably good enough
24 but, you know, these days with several formal
25 statistical tests used to validate QOL instruments, I

1 would have liked to have known, but you've answered
2 the question. It really hasn't been.

3 Are there any other issues around the
4 table?

5 (No response.)

6 DR. BORER: I have a couple of more if
7 you'll just bear with me for a second.

8 On slide CP76, it says no HF-related
9 primarily AEs reported in any of the three ThermoCool
10 subjects. These were the patients who had heart
11 failure, rather mild, at inclusion. So that did not
12 include the one person who developed heart failure
13 because that person presumably didn't have any
14 symptoms before. Is that right?

15 DR. YAROSS: That's our understanding, yes.

16 DR. BORER: Okay. That's fine. That's
17 great. Let me just see. Finally, the difference
18 between OUS 1 and other centers, which is going to be
19 I think an important topic for future discussion, the
20 difference seemed to me to be in the effectiveness of
21 ThermoCool which may have been due to any number of
22 factors, all of which you mentioned. There was also
23 a difference in the effectiveness of drugs. As I
24 recall, for the OUS 1, at the end of the day, there
25 was 11 percent effectiveness, and for everybody else,

1 17 percent. It doesn't seem like a big difference to
2 me. But I wonder if you did some formal statistical
3 assessment to see if it was likely that the
4 difference was significant, you know, if there was an
5 interaction by sight on the AAD side?

6 DR. YAROSS: I'll ask Dr. Berry to speak to
7 the statistical analyses.

8 DR. BERRY: We can provide a formal answer
9 after lunch but --

10 DR. BORER: Okay.

11 DR. BERRY: -- but in my quick and dirty
12 assessment of those numbers, is it's not
13 statistically different.

14 DR. BORER: Okay. If you want to check and
15 tell us after lunch, that's fine, but that was my
16 impression.

17 Now, the issue of mapping, I infer from the
18 data that we were sent that every site mapped every
19 patient. Is that right or wrong.

20 DR. YAROSS: The trial did call for CARTO
21 mapping in all cases, yes.

22 DR. BORER: And it was done.

23 DR. YAROSS: Yes.

24 DR. BORER: Okay. So the OUS 1 success
25 could not possibly have been related to the fact that

1 they mapped and maybe not everybody else did. Okay.
2 That's good.

3 Are there any other non-statistical
4 questions from the Panel?

5 (No response.)

6 DR. BORER: If not, I want to thank you
7 very much. If we're allowed to, Mr. Swink, you'll
8 have to tell me, if we're allowed to, can we take the
9 break now instead of waiting until 10:30?

10 MR. SWINK: Absolutely.

11 DR. BORER: Okay. You know, at Government
12 meetings, you have to be sure you're allowed to do
13 things. Okay. We will take a 15-minute break and
14 then we'll go onto the FDA presentation.

15 (Off the record.)

16 (On the record.)

17 DR. BORER: We will now have the FDA
18 presentation. And according to my slide here, this
19 will begin with Dr. Benjamin Eloff of the Division of
20 Cardiovascular Devices.

21 DR. ELOFF: Good morning. My name is Ben
22 Eloff. On behalf of the FDA review team, I would
23 like to extend my thanks to the distinguished members
24 of this Panel for their participation today.

25 I'm serving as the leader of the review

1 team for the application being discussed today from
2 Biosense Webster, to indicate the NaviStar ThermoCool
3 catheter for radiofrequency ablation of symptomatic
4 paroxysmal atrial fibrillation. This application is
5 a Panel-Track supplement to the original PMA Number
6 P030031.

7 The FDA review team consists of Randall
8 Brockman, a clinical officer; Laura Thompson, who
9 reviewed the statistical aspects of this application;
10 Ellen Pinnow, a postmarket epidemiologist; Martin
11 Hamilton, from our Division of Bioresearch
12 Monitoring; and myself, a biomedical engineer and
13 lead reviewer.

14 This device has a regulatory history dating
15 back to December 2003 when FDA approved the original
16 IDE study to investigate the use of the NaviStar
17 ThermoCool catheter in the treatment of symptomatic
18 paroxysmal atrial fibrillation or AF.

19 Subsequently, in 2004, FDA approved the
20 NaviStar and Celsius ThermoCool catheters for
21 treatment of type I atrial flutter.

22 In August 2006, FDA approved another PMA
23 for the NaviStar variant of the ThermoCool family
24 only, for the treatment of recurrent drug/device
25 refractory sustained monomorphic ventricular

1 tachycardia due to prior myocardial infarction in
2 adults.

3 The sponsor notified FDA in October 2007,
4 that they had completed an interim analysis and met
5 criteria for closure of the atrial fibrillation study
6 to enrollment.

7 In August of this year, the sponsor
8 submitted the present Panel-Track PMA supplement to
9 add an indication for treatment of drug-refractory
10 symptomatic paroxysmal AF to the approved indications
11 for the ThermoCool family.

12 FDA granted this application expedited
13 status based on the novelty of the device to address
14 a serious unmet public health need.

15 The existing indications for the ThermoCool
16 family of devices are as follows: All variants are
17 approved for the treatment of type 1 atrial flutter.
18 The NaviStar variants are approved for the treatment
19 of recurrent drug device refractory sustained
20 monomorphic ventricular tachycardia due to prior MI
21 in adults.

22 For all variants of the device, the sponsor
23 has proposed the new indication of treatment of drug-
24 refractory symptomatic paroxysmal atrial
25 fibrillation, which is the primary point of

1 discussion today.

2 The ThermoCool family of devices have in
3 common that they are deflectable catheters capable of
4 performing electrophysiologic mapping and ablation
5 procedures. A maximum of 50 watts of radiofrequency
6 ablation energy is delivered through a 3.5 millimeter
7 tip electrode that is actively cooled by an open loop
8 irrigation system.

9 The ThermoCool family of devices consists
10 of five distinct variants. The NaviStar ThermoCool,
11 the EZ Steer ThermoCool Nav, NaviStar ThermoCool RMT,
12 Celsius ThermoCool and EZ Steer ThermoCool.

13 The NaviStar ThermoCool, the EZ Steer
14 ThermoCool Nav, and NaviStar ThermoCool RMT are
15 approved for both treatment of atrial flutter and
16 ventricular tachycardia, as opposed to the Celsius
17 and EZ Steer ThermoCool catheters which are approved
18 for atrial flutter only.

19 This distinction is due to the presence of
20 a location sensor in the navigational variants which
21 allows for generation of advanced electroanatomic
22 maps of the heart which aid the physician in
23 diagnosing arrhythmias and identifying specific sites
24 for placing ablation lesions. This technology is
25 independent of the deflection mechanism which include