

Preparation of o-nitrophenyl- β -D-galactopyranoside

after Seidman, M., and Link, K. P.

See Babers & Goebel, JBC 105:473, 1934.; Glaser & Wulwek, Bioch. Z. 145,514, '21
42 g. o-nitrophenol in sol. 16.8 g NaOH/420 ml H₂O. Add sol. 88 g acetobrom-
galactose in 620 ml acetone. Keep at r.t. 5h. Remove solvent under reduced
pressure-- filter off long needles, and concentrate further until crystals
no longer come out. Wash ppt. water, air dry, recr. 95% EtOH. Weight: 56g.
mp 172.5. (onpg-~~xxxx~~ AcO₄)

To deacetylate, ~~add~~ 1g. onpg-OAc₄ in 50 ml MeOH and 1 ml. .4N Ba(OMe)₂ added.
Refrigerate, shake periodically. After 4 h. , clear sol., then crystals as long
hairlike needles. After 24 h. conc. under red. pr. to 10 ml., and a quant.
yield of onpg is obtained. Recr. 2x EtOH. MP 196-7.

For Acetobrom galactose:

~~mix~~ Gal-OAc₅ (accdg. to USBS c440) , 50 g. with 50 ml HBr/HAcO (EK)
Stand 2h. Add 200 cc CHCl₃, pour mix into 1 l. icewater, separate, wash CHCl₃
layer with NaHCO₃, then water (cold). Dry over Na₂SO₄. Evap. Take up syrup
in Et₂O to crystallize. If necessary, add Petether to turb. to ~~initiate~~
initiate cryst.