


# ROTENONE

5007



MW: 394.43

CAS: 83-79-4

RTECS: DJ2800000

METHOD: 5007, Issue 2

EVALUATION: FULL

Issue 1: 15 February 1984

Issue 2: 15 August 1994

OSHA : 5 mg/m<sup>3</sup>  
 NIOSH: 5 mg/m<sup>3</sup>; Group II Pesticide  
 ACGIH: 5 mg/m<sup>3</sup>

PROPERTIES: solid, MP 163 °C or 181 °C;  
 BP 220 °C @ 0.5 mm Hg;  
 d ca. 1 g/cm<sup>3</sup>; VP not significant

SYNONYMS: tubatoxin; cube

SAMPLING		MEASUREMENT	
<b>SAMPLER:</b>	FILTER (1-mm PTFE membrane)	<b>TECHNIQUE:</b>	HPLC; UV DETECTION
<b>FLOW RATE:</b>	1 to 4 L/min	<b>ANALYTE:</b>	Rotenone
<b>VOL-MIN:</b>	8 L	<b>EXTRACTION:</b>	4 mL acetonitrile; 30 min
<b>-MAX:</b>	400 L	<b>INJECTION VOLUME:</b>	10 µL
<b>SHIPMENT:</b>	routine	<b>MOBILE PHASE:</b>	60% methanol/40% water, 2 mL/min
<b>SAMPLE STABILITY:</b>	at least 7 days @ 25 °C in dark	<b>DETECTOR:</b>	UV @ 290 nm; 0.1A full-scale; 1-cm cell
<b>BLANKS:</b>	2 to 10 field blanks per set	<b>COLUMN:</b>	µ-Bondapak C <sub>18</sub> (30 cm x 3.9-mm ID stainless steel); ambient temperature
<b>BULK SAMPLE:</b>	desirable; 1 g	<b>CALIBRATION:</b>	solutions of Rotenone in acetonitrile
<b>ACCURACY</b>		<b>RANGE:</b>	0.04 to 1 mg per sample
<b>RANGE STUDIED:</b>	1 to 11 mg/m <sup>3</sup> [1] (100-L sample)	<b>ESTIMATED LOD:</b>	4 µg per sample [1,2]
<b>BIAS:</b>	- 0.6%	<b>PRECISION (S<sub>r</sub>):</b>	0.024 [1]
<b>OVERALL PRECISION (S<sub>r,T</sub>):</b>	0.079		
<b>ACCURACY:</b>	± 13.5%		

**APPLICABILITY:** The working range is 0.4 to 10 mg/m<sup>3</sup> for a 100-L air sample and the method is applicable to commercial formulations.

**INTERFERENCES:** None known. Rotenone, a naturally occurring insecticide, is adequately separated by HPLC from other compounds (e.g., sumatrol, α-toxicarol, deguelin, elliptone, malaccol, and tephrosin [3]) present in commercial cube root extracts [4]. Rotenone is sensitive to photodecomposition.

**OTHER METHODS:** This is Method S300 [2] in a revised format.

**REAGENTS:**

1. Acetonitrile, HPLC grade.\*
2. Methanol, HPLC grade.
3. Rotenone, 97% purity.
4. Water, distilled, HPLC grade.
5. Calibration stock solution, 3 mg/mL. Dissolve 0.075 g Rotenone in 25 mL acetonitrile. Prepare fresh daily in duplicate.
6. Recovery stock solution, 50 mg/mL. Dissolve 0.500 g Rotenone in acetone. Dilute to 10 mL. Prepare fresh daily.

\* See SPECIAL PRECAUTIONS.

**EQUIPMENT:**

1. Sampler: 37-mm, two-piece cassette containing 1- $\mu$ m PTFE membrane filter with backup pad.  
NOTE: Use an opaque cassette or otherwise shield the filter from light to minimize photodecomposition of Rotenone during and after sampling.
2. Personal sampling pump, 1 to 4 L/min, with flexible connecting tubing.
3. HPLC, UV detector, integrator and column (page 5007-1).
4. Jars, ointment, 60-mL, with PTFE-lined caps.
5. Vials, 4-mL, with PTFE-lined caps.
6. Syringes, 5-mL.
7. Filtration device, Swinney 13-mm with 1- $\mu$ m PTFE filters, or PTFE syringe filters.
8. Volumetrics, 10- and 25-mL.
9. Syringes, microliter, for sample injection and standard preparation.
10. Pipet, 4-mL, with pipet bulb.

---

**SPECIAL PRECAUTIONS:** Avoid breathing acetonitrile vapors; may cause skin irritation.

---

**SAMPLING:**

1. Calibrate each personal sampling pump with a representative sampler in line.
2. Sample at an accurately known flow rate between 1 and 4 L/min for a total sample size of 8 to 400 L. Do not exceed 2 mg total dust loading on the filter.
3. Collect a bulk sample (1 g) in a glass vial with PTFE-lined cap; ship separately from filters.

**SAMPLE PREPARATION:**

4. Open filter cassette; transfer filter to ointment jar.
5. Add 4.0 mL acetonitrile; gently swirl for 30 min.
6. Filter each sample using a 5-mL syringe with PTFE syringe filter or Swinney filtration device. Deliver filtrate to a 4-mL vial.

**CALIBRATION AND QUALITY CONTROL:**

7. Prepare at least six working standards daily in the range 0.01 to 1 mg Rotenone per sample.
  - a. Add known amounts of calibration stock solution to acetonitrile in 10-mL volumetric flasks and dilute to the mark.
  - b. Analyze together with samples and blanks (steps 9 and 10).
  - c. Prepare calibration graph (peak area vs. mg Rotenone).
8. Check recovery (R) with at least three spiked media blanks per sample set in the calibration range (step 7).
  - a. Add aliquots of recovery stock solution to blank filters with a microliter syringe. Air dry.
  - b. Analyze together with working standards (steps 4 through 6, 9 and 10).
  - c. Calculate recovery [(mg recovered - mg blank)/mg added].
  - d. Prepare recovery graph (R vs. mg Rotenone).

**MEASUREMENT:**

9. Set HPLC system according to manufacturer's recommendations and to conditions given on page 5007-1. Inject 10- $\mu$ L sample.  
NOTE: If peak area is above linear range of calibration graph, dilute, reanalyze, and apply appropriate dilution factor in calculations.
10. Measure peak area.

**CALCULATIONS:**

11. Read the mass, mg (corrected for recovery) of Rotenone found on the filter (W) and average media blank (B) from the calibration graph.
12. Calculate the concentration, C (mg/m<sup>3</sup>), of Rotenone in the air volume sampled, V (L):

$$C = \frac{(W - B) \cdot 10^3}{V}, \text{ mg/m}^3.$$

**EVALUATION OF METHOD:**

Method S300 [2] was issued on May 11, 1979, and validated over the range 1.16 to 11.1 mg/m<sup>3</sup> at 25°C and 760 mm, using 100-L samples [1, 5]. Overall precision,  $\hat{S}_{rT}$ , was 0.079 with average recovery 100.4%, representing a non-significant bias. The concentration of Rotenone (generated by Wright dust feeder using Ortho Rotenone Dust [1%; Chevron Chemical Co.] enriched to 10% Rotenone with analytical grade Rotenone [Aldrich Chemical Co.]) was independently verified by collection in dioxane and HPLC analysis. Recovery was 0.98 in the range 250 to 1000  $\mu$ g Rotenone per sample. Collection efficiency of the PTFE filter was found to be greater than 99% and no detectable Rotenone (LOD = 4  $\mu$ g) was found on Chromosorb 102 tubes placed behind the PTFE filters at 11.8 mg/m<sup>3</sup>. No loss was seen from spiked filters stored in the dark at room temperature for seven days.

**REFERENCES:**

- [1] Backup Data Report prepared under NIOSH Contract 210-76-0123, available as Order No. PB 82-114729 from NTIS, Springfield, VA 22161.
- [2] NIOSH Manual of Analytical Methods, 2nd. ed., V. 5, S300, U.S. Department of Health, Education, and Welfare, Publ. (NIOSH) 79-141 (1979).
- [3] Gunther, F. A., and R. G. Blinn. Analysis of Insecticides and Acaricides, 419-420, Interscience, NY (1955).
- [4] Bushway, R. J., B. S. Engdahl, B. M. Colvin, and A. R. Hanks. J. Assoc. Official Anal. Chemists, 58, 965 (1975).
- [5] NIOSH Research Report-Development and Validation of Methods for Sampling and Analysis of Workplace Toxic Substances, U.S. Department of Health and Human Services, Publ. (NIOSH) 80-133 (1980).

**METHOD REVISED BY:**

Jerome Smith, Ph.D., NIOSH/DPSE; S300 originally validated under NIOSH Contract 210-76-0123.