

# TUNGSTEN (soluble and insoluble)

7074

W      MW: 183.85      CAS: 7440-33-7      RTECS: Y07175000

**METHOD:** 7074, Issue 2

**EVALUATION:** FULL

**Issue 1:** 15 February 1984  
**Issue 2:** 15 August 1994

**OSHA :** no PEL

**PROPERTIES:** brittle metal; MP 3410 °C; valence 2,3,4,5,6

**NIOSH:** 1 mg/m<sup>3</sup> solubles; STEL 3 mg/m<sup>3</sup>  
5 mg/m<sup>3</sup> insolubles; STEL 10 mg/m<sup>3</sup>

**ACGIH:** 1 mg/m<sup>3</sup> solubles, STEL 3 mg/m<sup>3</sup>  
5 mg/m<sup>3</sup> insolubles, STEL 10 mg/m<sup>3</sup>

**SYNONYMS:** wolfram.

SAMPLING		MEASUREMENT	
<b>SAMPLER:</b>	FILTER (0.8-µm cellulose ester membrane)	<b>TECHNIQUE:</b>	ATOMIC ABSORPTION, FLAME
<b>FLOW RATE:</b>	1 to 4 L/min	<b>ANALYTE:</b>	tungsten
<b>VOL-MIN:</b>	200 L @ 1 mg/m <sup>3</sup>	<b>EXTRACTION:</b>	H <sub>2</sub> O; 6 min, 25 °C (soluble)
<b>-MAX:</b>	1000 L	<b>ASHING:</b>	1:1 HF:HNO <sub>3</sub> , 10 mL; 6 h, 150 °C (insoluble)
<b>SHIPMENT:</b>	routine	<b>FINAL SOLUTION:</b>	0.05 M NaOH/2% Na <sub>2</sub> SO <sub>4</sub> ; 25 mL
<b>SAMPLE STABILITY:</b>	stable at least 2 weeks @ 25 °C	<b>FLAME:</b>	nitrous oxide-acetylene, reducing
<b>BLANKS:</b>	2 to 10 field blanks per set	<b>WAVELENGTH:</b>	255.1 nm
<b>ACCURACY</b>		<b>BACKGROUND CORRECTION:</b>	none used
<b>RANGE STUDIED:</b>	0.42 to 2.0 mg/m <sup>3</sup> (soluble) [1]; 0.84 to 19.7 mg/m <sup>3</sup> (insoluble) [1]	<b>CALIBRATION:</b>	W in 0.05 M NaOH and 2% Na <sub>2</sub> SO <sub>4</sub>
<b>BIAS:</b>	+ 0.015	<b>RANGE:</b>	0.1 to 0.5 mg per sample (soluble); 0.25 to 12 mg per sample (insoluble)
<b>OVERALL PRECISION (S<sub>r,T</sub>):</b>	0.055 (soluble) [1]; 0.056 (insoluble) [1]	<b>ESTIMATED LOD:</b>	50 µg soluble W per sample; 125 µg insoluble W per sample [2]
<b>ACCURACY:</b>	± 12.5%	<b>PRECISION (S<sub>r</sub>):</b>	0.029 [1]

**APPLICABILITY:** The working range of this method is 0.25 to 5 mg/m<sup>3</sup> (soluble) and 0.6 to 5 mg/m<sup>3</sup> (insoluble) for a 400-L air sample. This is an elemental analysis and not compound-specific. For cemented tungsten carbide samples, cobalt can be quantitatively determined from the HCl extract used to remove interferences.

**INTERFERENCES:** None known. Ni, Mo, V, Mn, Cr, Co and Fe do not interfere at concentrations up to 50 times the tungsten concentration [3].

**OTHER METHODS:** This revises P&CAM 271 [1,3]. Method 7300 (ICP-AES) is an alternate measurement method; however, the extraction and ashing procedures described in this method must be followed.

**REAGENTS:**

1. Nitric acid, conc.\*
2. Hydrofluoric acid, conc.\*
3. Hydrochloric acid, conc.\*
4. Hydrochloric acid, 1% (v/v). Add 26.5 mL conc. HCl to 500 mL water; dilute to 1 L.
5. Sodium hydroxide, 0.5 M. Dissolve 20 g NaOH in 500 mL water; dilute to 1 L.
6. Sodium sulfate, 20% (w/v). Dissolve 20 g Na<sub>2</sub>SO<sub>4</sub> in 80 mL water; dilute to 100 mL.
7. Calibration stock solution, 10 mg W/mL. Commercially available or dissolve 1.5985 g Na<sub>2</sub>WO<sub>4</sub> (dried at 125 °C) in water; add 10 mL 0.5 M NaOH and dilute to 100 mL.
8. Water, distilled, deionized.

\* See SPECIAL PRECAUTIONS.

**EQUIPMENT:**

1. Sampler: cellulose ester membrane filter, 0.8-µm pore size, 37-mm diameter; in cassette filter holder.
2. Personal sampling pump, 1 to 4 L/min, with flexible connecting tubing.
3. Atomic absorption spectrophotometer with a nitrous oxide-acetylene burner head.
4. Tungsten hollow cathode lamp.
5. Two-stage regulators for N<sub>2</sub>O and acetylene.
6. Beakers, PTFE, 100-mL (with covers).\*
7. Volumetric flasks, 10-, 25- and 100-mL.\*
8. Assorted volumetric pipets, as needed, with pipet bulb.\*
9. Hotplate for use at 150 °C.
10. Steam bath or 100 °C hotplate.
11. Filtering apparatus (funnel, clamp, frit, holder) with 47-mm, 0.45-µm pore size cellulose ester membrane filter and collection vessel.

\* Clean all labware with conc. nitric acid and rinse thoroughly with distilled water before use.

---

**SPECIAL PRECAUTIONS:** Wear protective equipment (gloves, labcoats, and safety glasses) when working with HF. Do not use HF with glass labware. Perform all acid digestions in a fume hood.

---

**SAMPLING:**

1. Calibrate each personal sampling pump with a representative sampler in line.
2. Sample at an accurately known flow rate between 1 and 4 L/min for a total sample size of 200 to 1000 L. Do not exceed 2 mg total dust loading on the filter.

**SAMPLE PREPARATION:**

3. Open the cassette and transfer the filter to the filtering apparatus on top of a 47-mm filter.
4. Add 3 mL deionized water to the sample filter, allow it to stand 3 min, apply vacuum to transfer the extract to the holding vessel, and repeat with an additional 3 mL water.
5. Transfer the extracts to a 10-mL volumetric flask, add 1.0 mL 20% Na<sub>2</sub>SO<sub>4</sub>, dilute to volume, and analyze for soluble tungsten (steps 17 through 19).
6. Digest the filters and residue in a covered PTFE beaker with 5.0 mL conc. HNO<sub>3</sub> and 5.0 mL conc. HF at 150 °C.
7. Remove the beaker cover and reduce the volume to 2 mL (150 °C). Take to dryness at 100 °C.
8. Remove from heat, cool, add 10 mL 1% HCl, and agitate manually for 5 min.
9. Filter through an additional 47-mm filter (the filtrate may be analyzed for cobalt).
10. Digest the filter and residue in a PTFE beaker with 5.0 mL conc. HNO<sub>3</sub> and 5.0 mL conc. HF at 150 °C.
11. Remove the beaker cover and reduce the volume to 1 mL (150 °C). If filter residue (dark, charred material) is visible, add additional HNO<sub>3</sub> and HF (2 mL each).
12. Take to near-dryness at 100 °C.
13. Dissolve residue in beaker with 2.5 mL 0.5 M NaOH and 2.5 mL 20% Na<sub>2</sub>SO<sub>4</sub> at 100 °C (15 min).
14. Transfer to a 25-mL volumetric flask and dilute to volume.

### **CALIBRATION AND QUALITY CONTROL:**

15. Calibrate with at least six working standards over the range 0.05 to 12 mg W per sample.
  - a. Add known amounts of calibration stock solution to 100-mL volumetric flasks containing 10 mL 0.5 M NaOH and 10 mL 20% Na<sub>2</sub>SO<sub>4</sub>; dilute to volume with distilled, deionized water.
  - b. Analyze together with samples and blanks (steps 17 through 19).
  - c. Prepare calibration graph (absorbance vs. solution concentration, µg/mL).
16. Analyze three quality control blind spikes and three analyst spikes to ensure that the recovery and calibration graph are in control (steps 6 through 14 and 17 through 19).

### **MEASUREMENT:**

17. Set spectrophotometer to conditions on page 7074-1.
18. Aspirate standards and samples. Record absorbance readings.
19. If the absorbance values for the samples are outside of the range of the standards, dilute the solutions with 0.05 M NaOH/2% Na<sub>2</sub>SO<sub>4</sub> solution, reanalyze, and use the appropriate dilution factor in calculations.

### **CALCULATIONS:**

20. Determine the solution concentrations of tungsten in the sample, C<sub>s</sub> (µg/mL), and average media blank, C<sub>b</sub> (µg/mL), from the calibration graph.
21. Using the solution volumes, mL of the sample (V<sub>s</sub>) and media blank (V<sub>b</sub>), calculate the concentration of tungsten, C (mg/m<sup>3</sup>), in the air volume sampled, V (L):

$$C = \frac{(C_s V_s - C_b V_b)}{V}, \text{ mg/m}^3.$$

### **EVALUATION OF METHOD:**

This method is based on P&CAM 271 [3,4] and was further evaluated and ruggedized under contract [1]. Laboratory testing with spiked samples over the ranges of 0.17 to 0.8 mg per sample (soluble) and 0.34 to 7.9 mg per sample (insoluble) gave precisions (  $\bar{S}_r$  ) of 0.079 (soluble) and 0.076 (insoluble) [1]. Both soluble and insoluble tungsten species collected on filters were stable for two weeks [1].

### **REFERENCES:**

- [1] Carlin, L. M., G. Colovos, D. Garland, M. Jamin, M. Klenck, T. Long, and C. Nelson. Analytical Methods Evaluation and Validation: Arsenic, Nickel, Tungsten, Vanadium, Talc and Wood Dust, NIOSH Contract No. 210-79-0060 (1981), available as Order No. PB 83-155325 from NTIS, Springfield, VA 22161.
- [2] User check, UBTL, Inc., NIOSH Sequence #4213-M (unpublished, May 23, 1984).
- [3] NIOSH Manual of Analytical Methods, 2nd ed., Vol. 4, P&CAM 271, U.S. Department of Health, Education, and Welfare, Publ. (NIOSH) 78-175 (1978).
- [4] Hull, R. D. Analysis of Tungsten by Atomic Absorption Spectrophotometry: A Feasibility Study, NIOSH, DPSE, MRB, IMDS Technical Report (unpublished, December, 1977).

### **METHOD WRITTEN BY:**

R. DeLon Hull, Ph.D., NIOSH/DBBS; additional data obtained under NIOSH Contract 210-79-0060.