Appendix B1

Protocol for HepG2 Cells + Receptor + Reporter and/or β -gal plasmids for Use in Steroid Hormone Receptor Assays

(Provided by Dr. Kevin Gaido, CIIT Centers for Health Research, Research Triangle Park, NC, USA) [This page intentionally left blank]

TransIT Transfection Method of HepG2 Cells for Use in Steroid Hormone Receptor Assays

1. MATERIALS AND SOURCES:

- a. TransIT-LT1 Transfection Reagent, supplier: Mirus Corporation, CAT. #: MIR 2300.
- b. 1X Phosphate Buffered Saline Solution.
- c. Plasmid DNA's of choice: i.e., receptor, reporter, and/or gal plasmids.
- d. Phenol red-free Minimum Essential Medium (MEM).
- e. Complete phenol red-free Minimum Essential Medium (MEM), with stripped (or charcoal/dextran treated) fetal bovine serum.
- f. 0.02% EDTA.
- g. Trypsin, 2.5%.
- h. Dimethyl sulfoxide.
- i. 1M Sodium pyruvate.
- j. L-glutamine (100X).

2. EQUIPMENT AND SUPPLIES:

- a. Incubator with 5% CO₂/air, 37°C
- b. Vortexer
- c. 10 µl, 100 µl, 200 µl, and 1000 µl Eppendorf pipettor or equivalent
- d. pipet tips
- e. 1, 2, 5, 10, 25, and 50 ml pipets
- f. 500 ml screw cap glass bottles, sterile
- g. 24 well tissue culture plates
- h. 15 and 50 ml centrifuge tubes, sterile, polypropylene
- i. 17x100, polypropylene snap-cap tubes, sterile, round bottom
- j. 1.5 ml siliconized polypropylene screw-cap vials

3. PREPARATION:

a. 0.12% Trypsin/0.02% EDTA.

In 500 ml sterile screw cap glass bottle, sterilely transfer 190 ml 0.02% EDTA. Add 10 ml of 2.5% trypsin. Store at 4°C.

b. Complete phenol red-free MEM.

To 500 ml of phenol red-free MEM, add 0.5 ml 1M sodium pyruvate solution, 10.0 ml glutamine, and 50 ml resin-stripped (or charcoal dextran treated) fetal bovine serum. Store 4°C.

c. Chemicals.

Dissolve chosen chemical to make a 0.1M stock solution using appropriate vehicle. Make serial dilutions in 1.5 ml polypropylene screw-cap vials to yield a standard curve of concentrations varying from 10^{-5} M to 10^{-11} M (may be changed as necessary).

4. **PROCEDURE:**

Plating Cells.

- a. Aspirate medium from 150 mm plate of 75-80% confluent HepG2 cells and rinse with 10 ml of 0.02% EDTA.
- b. Place 10 ml of 0.12% trypsin/0.02% EDTA on plate.
- c. Place in incubator until cells begin to detach (~5 min).
- d. After cells have detached, pipette vigorously to remove the cells and transfer to 50 ml polypropylene centrifuge tube containing complete phenol red-free MEM.
- e. Rinse plate with complete phenol red-free MEM and add to tube.
- f. Centrifuge at 1000 RPM for 5 min at 4°C.
- g. Carefully aspirate supernatant and resuspend the pellet in phenol red-free complete MEM.
- h. Take cell count. Plate cells in 24-well tissue culture dishes at 10⁵ cells/0.5 ml complete phenol red-free MEM. Swirl the plate gently to spread cells evenly in wells.
- i. Place cells in 37°C incubator with 5% CO₂/air for 18 hours.

<u>Transfecting Cells</u>.

In a 17x100 ml round bottom, polypropylene, snap cap tube, add the following reagents: (For transfection of a 24-well tissue culture plate)

- a. 0.65 ml of phenol red-free MEM without any additives.
- b. Appropriate amount of TransIT LT1 reagent. For every μg of DNA plasmid, add 2 μ1 of TransIT LT1 reagent. (11 μl of TransIT LT-1 reagent is needed for the suggested amounts of plasmid listed in 3. Below.) Mix <u>very gently</u> and let sit at RT for at least 5 min.
- Carefully add appropriate amounts of receptor, promoter, and gal plasmids.
 This may vary depending on the application. A suggestion for amounts is as follows:

Estrogen Assay Example

Receptor Plasmid: 7 ng/well pCMV Plasmid (-gal): 30 ng/well Promoter Plasmid: 200 ng/well

- a. Mix <u>very gently</u> and let sit at RT for at least 5 min.
- b. To each well of the 24 well plate containing HepG2 cells, carefully add 25 µl of the TransIT/DNA complex.
- c. Place plate in incubator and allow to incubate for 3 hr at 37°C.

Treating cells.

- Dilute chosen chemicals 1:1000 in complete phenol red-free MEM, to create final a. concentrations ranging from 10⁻⁵ to 10⁻¹¹ M (this may vary as necessary).

 After the 3 hr incubation, aspirate the media and add 0.5 ml/well of the chemical
- b. diluted in media.
- Return plate to incubator and incubate for 24 hr. Collect cell lysate for -gal and c. luciferase assays.

Lysis Procedure

1. MATERIALS:

- a. Phosphate Buffered Saline (1X PBS).
- b. Tris base.
- c. Trans-1, 2-diaminocyclohexane-N, N, N', N'-tetraacetic acid (CDTA).
- d. Glycerol.
- e. Phosphoric Acid.
- f. Triton X-100.
- g. 1M Dithiothreitol (DTT).
- h. Transfected cells plated in 24-well plate.

2. EQUIPMENT AND SUPPLIES:

- a. $1-200 \mu l$ Pipettor
- b. Multi-channel pipettor, $1-100 \mu l$
- c. $1-200 \mu l$ pipette tips
- d. Pipette aid
- e. Vacuum system with hazardous waste flask attached
- f. pH meter
- g. 5 3/4" Pasteur pipette
- h. 500 ml squeeze water bottle
- i. 96 well ELISA plate
- j. 96 well Plate, white
- k. 250 ml Glass beakers
- 1. 100 and 200 ml Graduated cylinders
- m. Stirrer and stir bars

3. PREPARATION:

- a. 5X Lysis Solution.
- a. Weigh out 3.03g Tris Base and 0.695g CDTA and place in 250 ml beaker.
- b. Dissolve completely in 60 ml of dH₂O.
- c. Measure 100 ml glycerol in 100 ml graduated cylinder, pour into fresh 250 ml beaker.
- d. Rinse 100 ml cylinder with Tris base/CDTA. Add to glycerol in 250 ml beaker. Mix well.
- e. pH to 7.8 with phosphoric acid (H_3PO_4) if necessary.
- f. Add dH₂O to 200 ml.
- g. Add 5 ml of 100% Triton X-100 (solution will look cloudy/milky). Store room temperature.

b. 1X Lysis Solution.

In 50 ml centrifuge tube, dilute 5X Lysis Solution to 1X by diluting 1 ml 5X lysis solution into 4 ml dH₂O. Add **30 \mul** 1M DTT per **10 ml** 1X lysis solution. Make fresh each time. Make up enough 1X lysis solution to dispense 65 μ l per well.

4. PROCEDURE:

- a. Aspirate media from wells and rinse with 0.5ml of PBS per well.
- b. Aspirate PBS from wells and with multi-channel pipettor; dispense 65 μl of 1X lysis solution per well.
- c. Let sit at room temperature for 20 min, rocking occasionally.
- d. Transfer 30 μl of cell lysate to 96 well ELISA plate. This will be used for the galactosidase assay.
- e. Transfer 20 µl of cell lysate to a 96 well white plate. This will be used for the luciferase assay.

ß-Galactosidase Assay Using Chlorophenol Red-β-D-galactopyranoside

1. MATERIALS AND SOURCES:

- a. Chlorophenol red-β-D-galactopyranoside (CPRG).
- b. Disodium phosphate (Na₂HPO₄•7H₂0).
- c. Monosodium phosphate (NaH₂PO₄•H₂0).
- d. Potassium chloride (KCl).
- e. Magnesium sulfate (MgSO₄•7H₂O).
- f. -Mercaptoethanol (2-ME).

2. EQUIPMENT AND SUPPLIES:

- a. Spectrophotmetric microplate reader, with a 575 nm filter and kinetics capability
- b. Multi-channel pipettor
- c. Graduated cylinder, 1000 ml
- d. Balance
- e. Stir plate
- f. Magnetic stir bar
- g. 1-100 μl pipettor
- h. 1-100 μl pipet tips
- i. Pipettor reservoirs
- j. 0.2 µ Filter unit
- k. 96 well ELISA plate
- 1. 1 L beaker
- m. 50 ml centrifuge tube, polypropylene, sterile

3. CPRG BUFFER PREPARATION:

- a. Weigh out in 1 L beaker:
- b. 16.1 g $Na_2HPO_4 \cdot 7H_2O$
- c. 5.5 g NaH₂PO₄•H₂0
- d. 0.75 g KCl
- e. 0.25 g MgSO₄•7H₂O
- f. Dissolve in 800 ml of distilled water with stirring.
- g. Adjust pH to 7.8.
- h. Transfer to 1000 ml graduated cylinder. Bring up to 1000 ml with distilled water.
- i. Filter sterilize. Store at room temperature.

4. ASSAY PROCEDURE:

a. Pipet 30 µl of cell lysate into a 96 well plate (usually done in triplicate).

- b. **PER WELL OF 96 WELL PLATE**, add 170 µl of CPRG reagent made up as follows: 80 µg CPRG dissolved in 20 µl distilled water, 150 µl of CPRG buffer, and 0.84 µl 2-ME (1/200 dilution).
- c. Using multi-channel pipettor, dispense 170 μ l of CPRG reagent into each well containing lysate. For plate blank, use 30 μ l of lysis solution and add 170 μ l of CPRG reagent.

Set spectrophotometer microplate reader to kinetic endpoint and read the plate at 575 nm at 1 min intervals for 30 min to obtain Vmax. Samples will change from yellow to dark red as reaction occurs.

[This page intentionally left blank]