

EMBRYO CRYOPRESERVATION PROTOCOL

Superovulate females with PMS (5 IU per female);
48 hours later, inject females with HCG (5 IU per female) & mate;
Embryo collection must be scheduled in 72 hours.

Embryo Collection and Cryopreservation:

- 1- Dissect the oviducts.
- 2- Flush the embryos using a dull needle by injecting Modified D-PBS (Specialty Media CAT# MR-006-C) through the infundibulum and Oviduct.
- 3- Collect the 8-cell embryos into cryovials containing 0.1 ml of M-DPBS, 25-30 embryos per vial.
- 4- Set on ice bath at 0°C.
- 5- Add 0.1 ml CPA (containing DMSO) and let the vials sit on ice for 30 min.
- 6- Place vials in -6°C ice bath for 2 min.
- 7- Seed vials (media in Pasteur pipettes for seeding must be -10°C).
- 8- Transfer vials to the controlled-rate freezer already set at -6°C & freeze to -80°C (approximate freezing time 2 ½ hours, -1°C every 2 min).
- 9- Transfer vials to LN2 (liquid phase) for long-term storage.

Thawing Procedure:

- 1- Thaw each vial at room temperature for 10 min.
- 2- Slowly add 0.8 ml of M-DPBS to the vial containing the embryos (this process must be done gently, drop by drop). Mix the volume in the cryovial by pipetting the sample.
- 3- Draw up total volume containing the embryos and place in a 60mm. organ culture dish.
- 4- Place the dish on ice and transfer the embryos immediately.
- 5- Embryos should be transferred into the oviduct of the pseudopregnant recipients. We recommend that 10-12 embryos be transferred per female recipient. Our laboratory uses B6D2F1 females for this purpose.