

Detection of MMP-9 in Formalin-Fixed Paraffin-Embedded Mouse Tissue

Reagents:

[1X Automation Buffer](#)
[3% Hydrogen Peroxide](#)
[Antibody Diluent](#)
[Citrate Buffer](#)
[DAB Chromagen](#)
[Hematoxylin](#)

Antibody Dilutions:

Blocking Serum: Normal Donkey Serum
Jackson Immunoresearch Laboratories, Inc.
West Grove, PA 19390
www.jacksonimmuno.com
1-800-367-5296
Catalog# 017-000-001

Avidin Biotin Blocking Kit
Vector Laboratories, Inc.
Burlingame, CA 94010
www.vectorlabs.com
1-800-227-6666
Catalog #SP-2001

Primary antibody: Goat anti-mouse MMP-9 Antibody
R&D Systems, Inc
Minneapolis, MN 55413
www.rndsystems.com
1-800-343-7475
Catalog # AF909

Negative Serum Control: Normal Goat Serum
Jackson Immunoresearch Laboratories, Inc.
West Grove, PA 19390
www.jacksonimmuno.com
1-800-367-5296
Catalog # 005-000-121

Secondary antibody: Biotin-SP-conjugated Donkey anti-goat IgG
Jackson Immunoresearch Laboratories, Inc.
West Grove, PA 19390
www.jacksonimmuno.com
1-800-367-5296
Catalog # 705-065-147

Label antibody: Super Sensitive Peroxidase-Conjugated Streptavidin Label
Biogenex Laboratories
San Ramon, CA 94583
www.biogenex.com
1-800-421-4149
Catalog# HK330-9K

Staining Procedure

- Positive Control Tissue: LPS-treated liver
- Stain localization: Cytoplasmic

Deparaffinize and hydrate slides through the following solutions.

Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Quench endogenous peroxidase by placing slides in 3% Hydrogen Peroxide for 15 minutes.
2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
3. Unmasking Techniques using the decloaker.
Add 500 ml D/W to the pan of the decloaker.
Place full rack of slides in 200 ml of 1X citrate buffer and place in the decloaker.
Decloak for 5 min. Pressure _____
Depressurize for 10 min.
Remove pan top and cool for 10 min. Temp _____
Rinse in D/W, 2x for 3 min each
4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes.
5. Apply 5% Normal Donkey Serum for 20 minutes at room temperature.
Lot# _____ Reconstituted Date _____

DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.

6. Apply Avidin/Biotin block

Lot# _____ Exp Date _____ New Kit yes / no

Apply avidin block - 15 min @ RT.

Quick rinse in 1X AB.

Apply biotin block - 15 min @ RT.

Wipe excess block.

DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.

7. Apply primary antibody (MMP-9) at a 1:1000 dilution and incubate for one hour at room temperature.

Lot# _____ Exp Date _____

For negative control slides, normalize the protein concentration of normal goat serum to the protein concentration of the primary antibody and use this to make a 1:1000 dilution. Apply to the slides and incubate for one hour at room temperature.

Lot# _____ Reconstituted Date _____

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply Donkey anti-goat secondary antibody at a 1:200 dilution and incubate for 30 minutes at room temperature.

Lot# _____ Reconstituted Date _____

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply Biogenex Streptavidin Label and incubate for 30 minutes at room temperature.

Lot # _____ Exp. Date _____

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.

(Add 1 drop of DAB per ml of substrate)

Lot# _____ Exp Date _____ new kit yes / no

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30 seconds.

16. Rinse in tap water until water is clear.

17. Place slides in 1X Automation buffer for 1 minute with gentle agitation to blue slides.

18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% Ethanol	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip

Updated 04/06/06