

1
2
3
4
5
6
7
8
9
10
11
12
13

Appendix A

LLNA: BrdU-ELISA Protocol

13
14
15
16
17
18
19
20
21
22
23
24
25
26
27

[This Page Intentionally Left Blank]

28 **1.0 Introduction/Principle**

29 This document describes the recommended standard operating procedure for the non-
30 radioisotopic modification of the LLNA, which is based on BrdU incorporation in place of
31 tritiated thymidine to measure lymph node cell proliferation. This document is based on the
32 protocol used in the JSAAE multi-laboratory validation study of the LLNA: BrdU-ELISA,
33 *Recommended Standard Operating Procedure for the Non-Radioisotopic Local Lymph Node*
34 *Assay using BrdU-ELISA (Non-RI LLNA), version 1.20, July 31, 2008*, by Masahiro
35 Takeyoshi, Ph.D., Chemicals Evaluation and Research Institute, Japan. This

36 **2.0 Description of the Method**

37 The method is practically identical to the standard LLNA methodology excluding the use of
38 BrdU and colorimetric detection. A single intraperitoneal injection (5 mg/mouse per
39 injection) of BrdU is made on day 4. This administration schedule was decided as the most
40 effective labeling protocol to yield maximum SI values based on preliminary study data with
41 several different protocols. Approximately 24 h after the BrdU injection, the auricular lymph
42 nodes are removed, weighed, and stored at -20°C until analysis using an enzyme-linked
43 immunosorbent assay to measure the level of BrdU incorporation (BrdU-ELISA).

44 The cell proliferation response is measured by a commercial BrdU detection kit (i.e., Roche
45 Diagnostics GmbH, Roche Applied Science, 68298 Mannheim, Germany; Cat. No. 11 647
46 229 001). To perform the BrdU-ELISA, the lymph nodes are crushed, passed through a #70
47 nylon mesh. The lymph node cells (LNC) from individual animals are suspended in 15 ml of
48 physiological saline. The cell suspension is added to the wells of a flatbottom microplate in
49 triplicate. After fixation and denaturation of the LNC, anti-BrdU antibody is added to each
50 well, and after rinsing, substrate solution containing tetramethylbenzidine (TMB) is added
51 and allowed to produce chromogen. Absorbance at 370 nm with a reference wavelength of
52 492 nm is defined as the BrdU labeling index.

53

54

54 **2.1 Animals**

55 *2.1.1 Animal source*

56 Young adult female mice (nulliparous and non-pregnant) of the CBA/JN or other
57 recommended mouse strains, such as CBA/Ca or CBA/J strain, should be used at age 8-12
58 weeks. All animals should be age matched (preferably within a one-week time frame).

59 *2.1.2 Quarantine and Acclimation*

60 Healthy animals in good general condition on arrival should be quarantined for more than
61 five days. During the quarantine and acclimation period, clinical signs, body weights and
62 excrement of the animals should be observed.

63 *2.1.3 Grouping*

64 Animals confirmed to be in good health with favorable body weight gains during the
65 quarantine and acclimation period should be allocated to groups by a stratified randomization
66 or other appropriate methods before the start of the study.

67 *2.1.4 Identification*

68 Animals should be identified by colored marks on the tails, ear tags, or other appropriate
69 methods.

70 *2.1.5 Animal Husbandry*

71 The animals should be housed in an animal room maintained at a temperature of $22\pm 3^{\circ}\text{C}$ and
72 a relative humidity of 30-70%. The rooms should be artificially lighted for 12 h daily, and the
73 animals should be given free access to conventional laboratory diet and drinking water.

74 **2.2 Chemicals and Vehicle**

75 *2.2.1 Vehicle*

76 The solvent/vehicle should be selected on the basis of maximizing the test concentrations
77 while producing a solution/suspension suitable for application of the test substance. In order
78 of preference, recommended solvents/vehicles are acetone/olive oil (4:1 v/v), DMF, MEK,
79 propylene glycol, and DMSO, but others may be used.

80 2.2.2 *Test Chemicals*

81 Solid test substances should be dissolved in appropriate solvents or vehicles and diluted, if
82 appropriate, prior to dosing of the animals. Liquid test substances may be dosed directly or
83 diluted prior to dosing. Fresh preparations of the test substance should be prepared daily
84 unless stability data demonstrate the acceptability of storage.

85 2.2.3 *Controls*

86 Concurrent negative (vehicle) and positive controls should be included in each test. Positive
87 control (50% HCA, CAS RN. 101-86-0) should be used to ensure the appropriate
88 performance of the assay. The positive control should produce a positive LLNA response at
89 an exposure level expected to give an increase in the stimulation index (SI) >2 over the
90 negative (vehicle) control group.

91 2.2.4 *Dose selection*

92 Doses are selected from the concentration series 100%, 50%, 25%, 10%, 5%, 2.5%, 1%,
93 0.5%, etc. The maximum concentration tested should be the highest achievable level while
94 avoiding overt systemic toxicity and excessive local irritation. All test solutions should be
95 prepared in a day of application unless the stability is confirmed in advance.

96 2.2.5 *Preparation of BrdU*

97 BrdU should be accurately weighed and dissolved in physiological saline for injection) to
98 make 10 mg/ml solution. The BrdU solution should be sterilized by a commercial filtration
99 system (i.e. MILLEX®-HV, MILLIPORE etc.). The BrdU solution can be prepared before
100 administration and stored in a freezer below - 20°C until use.

101 **2.3 Animal Experiment**

102 2.3.1 *Grouping*

103 A minimum of four successfully treated animals is used per dose group, with a minimum of
104 three consecutive concentrations of the test substance plus a negative (vehicle) control and a
105 positive control group.

106

107

107 **Table 1 Structure of LLNA: BrdU-ELISA Test Groups**

| Group | Number of Animals |
|----------------------------|-------------------|
| Negative (vehicle) control | 4 |
| Positive control (50% HCA) | 4 |
| Test substance-low dose | 4 |
| Test substance-middle dose | 4 |
| Test substance-high dose | 4 |

108

109 **2.3.2 Sensitization Procedure**

110 Apply 25µl of test solution to the dorsum of both ears of the mice using micro volume pipette
111 daily for three consecutive days.

112 **2.3.3 BrdU Administration**

113 A single intraperitoneal injection of 0.5 ml of BrdU solution (5 mg/mouse/injection) should
114 be given to the mice 48 hours (h) after the final sensitization.

115 **2.3.4 General Condition**

116 Clinical signs should be observed at least once a day.

117 **2.3.5 Body Weights**

118 Body weights should be measured on the day of the first test substance application and on the
119 lymph nodes are collected.

120 **2.3.6 Collection of Lymph Nodes And Measurement of Lymph Node Weight**

121 Approximately 24 h after BrdU injection, the auricular lymph nodes should be removed. The
122 lymph nodes should be carefully dissected and trimmed of fascia and fat, weighed, and stored
123 individually in a 1.5 ml centrifuge tube at -20°C until the ELISA is performed.

124 **2.4 BrdU-ELISA**

125 The incorporation of BrdU into lymph node cells should be determined using a commercial
126 cell proliferation assay kit (Roche Diagnostics GmbH, Roche Applied Science, 68298
127 Mannheim, Germany; Cat. No. 11 647 229 001) after they are crushed and suspended in

128 physiological saline. The absorbance is defined as the BrdU labeling index. Follow the
129 instructions in the assay kit.

130 **2.5 Preparation of Reagents in the BrdU-ELISA Kit**

131 The assay method should be according to the instruction manual in the assay kit excluding
132 preparation of the BrdU labeling solution.

133 *2.5.1 Peroxidase Conjugated Anti-BrdU Antibody (Anti-BrdU-POD) Stock Solution*

134 Dissolve Anti-BrdU-POD (bottle 3) in 1.1 ml double distilled water for 10 minutes and mix
135 thoroughly. This solution can be stored at 2-8°C for several months. For long-term storage it
136 is recommended to store the solution in aliquots at -15 to -25°C.

137 *2.5.2 Anti-BrdU-POD Working Solution*

138 Dilute Anti-BrdU-POD stock solution 1:100 with antibody dilution solution (bottle 4). For
139 one 96-well microtiter plate dilute 100 µl Anti-BrdU-POD stock solution in 10 ml antibody
140 dilution solution (bottle 4). Prepare shortly before use.

141 *2.5.3 Washing Solution*

142 Dilute washing buffer concentrate (bottle 5) 1:10 with double distilled water. For one 96-well
143 microtiter plate, dilute 10 ml washing buffer concentrate (bottle 5) with 90 ml double
144 distilled water. This solution can be stored at 2-8°C for several weeks.

145 **2.6 Preparation of Cell Suspension of Lymph Nodes**

146 The procedure for preparing the LNC suspension is a critical step of this assay; it is most
147 important to crush the lymph node and suspend the LNC completely. Every technician
148 should establish the skill in advance. The lymph nodes in negative control animals are very
149 small, so careful operation is required to avoid an artificial effect on SI values.

150 *2.6.1 Optimizing Assay Condition*

151 Mean absorbance of negative (vehicle) control group should be within 0.1-0.2. Because the
152 absorbance depends on the combination of assay apparatus and the target volume of the LNC
153 suspension, every laboratory should decide their own optimal target volume of LNC
154 suspension in advance so that the absorbance of the negative control is within 0.1-0.2. The
155 volume is expected to be approximately 15 µl. The volume of the LNC suspension for all

156 test animals should be adjusted to the optimized volume.

157 2.6.2 Preparation of LNC Suspension

158 A small amount (approximately 0.3 ml) of physiological saline should be added to the
159 centrifuge tube that contains the collected lymph node. The lymph node should be crushed
160 with a disposable plastic pestle to make the LNC suspension. The LNC suspension should be
161 passed through a #70 nylon mesh and adjusted to the optimal target volume in a 50 ml Falcon
162 tube.

163 [Note: Although a crushing apparatus other than a plastic pestle can be used to prepare the
164 LNC, the target volume of the LNC suspension should be adjusted to the optimized volume.]

165 2.7 Assay Flow (BrdU-ELISA)

- 166 1. The cell suspension (100 μ l) is added to the wells of a flat-bottom microplate
167 (three wells per sample) after mixing thoroughly with a Vortex.
168 Simultaneously, three blank wells should be prepared by adding 100 μ l of
169 physiological saline.
- 170 2. After filling all sample wells and blank wells, the plate should be centrifuged
171 at 300 x g for 10 minutes.
- 172 3. Remove 3/4 of the supernatant volume. Great care should be taken so that the
173 LNC are not aspirated.
- 174 4. The assay plate should be dried completely in a hot-air oven.
- 175 5. Add 200 μ l of Fix-Denat solution and allow plate to stand for 30 minutes at
176 room temperature.
- 177 6. Remove the Fix-Denat solution completely.
- 178 7. Add 100 μ l of anti-BrdU-POD antibody working solution and allow it to react
179 for 1 h.
- 180 8. Remove the anti-BrdU-POD antibody solution completely.
- 181 9. Add 200 μ l of wash solution into each well, and wash the well by pipetting

- 182 10 times. Discard the wash solution completely.
- 183 10. The wash step (Step 9) should be repeated twice (three times in total).
- 184 11. Add 100 μ l of TMB substrate solution and let it stand for 15 minutes at room
185 temperature in a dark place.
- 186 12. Measure an absorbance (ABS) at 370 nm with a reference wavelength of 492
187 nm. When using stop solution (1M sulfuric acid, 25 μ l/well), measure ABS at
188 450 nm with a reference wavelength of 690 nm.

189 **3.0 Calculation of Results**

190 BrdU labeling index and Stimulation Index (SI) are defined as follows:

191 **3.1 Without Stop Solution**

192 BrdU labeling index = $(ABS_{370} - ABS_{\text{blank}370}) - (ABS_{490} - ABS_{\text{blank}490})$

193 **3.2 With Stop Solution**

194 BrdU labeling index = $(ABS_{450} - ABS_{\text{blank}450}) - (ABS_{650} - ABS_{\text{blank}650})$

195 **3.3 Stimulation Index**

196 Stimulation Index (SI) =
$$\frac{\text{BrdU labeling index for each test animal}}{\text{Mean BrdU labeling index for concurrent vehicle control group}}$$

197

198 **4.0 Evaluation of Results**

199 **4.1 Success Criteria for Each Experiment**

200 Employing the optimized assay condition described previously, the mean SI for the positive
201 control group (50% HCA) should be equal to or greater than 2. If not, all data derived from
202 the experiment should not be used for evaluation.

203 **4.2 Evaluation of the Results**

204 The mean BrdU labeling index for each animal should be calculated based on the results of
205 BrdU ELISA. The SI for each animal should be calculated by dividing of the mean BrdU
206 labeling index for each treated animal by the mean BrdU labeling index of the concurrent
207 vehicle control group. A positive response is defined as mean SI of the test group ≥ 2 .

208 5.0 References

- 209 Kimber I, Dearman RJ, Scholes EW, Basketter DA. 1994. The local lymph node assay:
210 developments and applications. *Toxicology* 93:13-31.
- 211 OECD. 2002. Guideline For Testing of Chemicals - Test Guideline 429: Skin Sensitisation:
212 Local Lymph Node Assay. Paris:OECD. Available:
213 [http://puck.sourceoecd.org/vl=4134815/cl=32/nw=1/rpsv/cw/vhosts/oecdjournals/1607310x/
214 v1n4/contp1-1.htm](http://puck.sourceoecd.org/vl=4134815/cl=32/nw=1/rpsv/cw/vhosts/oecdjournals/1607310x/v1n4/contp1-1.htm)
- 215 EPA. 2003. Health Effects Test Guidelines: OPPTS 870.2600 - Skin Sensitization. EPA 712–
216 C–98–197. Washington, DC: U.S. Environmental Protection Agency. Available:
217 [http://www.epa.gov/opptsfrs/publications/OPPTS_Harmonized/870_Health_Effects_Test_G
218 uidelines/Series/870-2600.pdf](http://www.epa.gov/opptsfrs/publications/OPPTS_Harmonized/870_Health_Effects_Test_Guidelines/Series/870-2600.pdf)
- 219 Takeyoshi M, Yamasaki K, Yakabe Y, Takatsuki M, and Kimber I. 2001. Development of
220 non-radio isotopic endpoint of murine local lymph node assay based on 5-bromo-2'-
221 deoxyuridine (BrdU) incorporation. *Toxicology Letters* 119:203-208.