

Detection of Marck's Protein in Formalin-Fixed, Paraffin Embedded Rodent Tissue

Reagents:

[1X Automation Buffer](#)

[3% Hydrogen Peroxide](#)

[Antibody Diluent](#)

[DAB Chromagen](#)

Hematoxylin

Antibody Information:

Blocking Serum: Normal Horse Serum

Jackson Immunoresearch Laboratories, Inc.

West Grove, PA 19390

www.jacksonimmuno.com

1-800-367-5296

Catalog #008-000-001

Avidin Biotin Blocking Kit

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Catalog #SP-2001

Primary antibody : Marck's Protein

(2F12)

Provided by Dr. Perry Blackshear

NIEHS

Suggested dilution: 1:100

Negative control serum: Pre-immune Mouse Serum

(6F6)

Supplied by Dr. Perry Blackshear

National Institute of Environmental Health Sciences

Research Triangle Park, NC

Suggested dilution 1:100

Secondary antibody: Biotinylated anti-Mouse IgG (made in Horse)

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Suggested dilution: 1:800

Catalog #BA-2001

Label antibody: Vector EliteVectastain® ABC

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Catalog #PK-6101

Staining Procedure

-Positive Control Tissue: Bovine Retina

-Stain Localization: Membrane

Deparaffinize and hydrate slides through the following solutions.

Xylene	2 times	5 minutes
100% EtOH	2 times	3 minutes
95% EtOH	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.

2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

3. Perform Heat Induced Epitope Retrieval using Microwave Oven

Place a full rack of slides in distilled water and place in the microwave oven.

When using a tissue tek™ container, add 250 mls of distilled water.

Microwave slides at 80% power for 2 minutes

Rest for one minute.

Repeat microwave/rest cycle at 80% power 2 times

Following the third microwave cycle, transfer slides to fresh distilled water and allow slides to cool for 15 minutes. Temp. _____

4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

5. Block in 5% Normal Horse Serum for 20 minutes.

Lot# _____ Reconstituted Date _____

6. Apply Avidin/Biotin block

New kit yes / no Lot# _____ Exp Date _____

Apply avidin block - 15 min @ RT.

Quick rinse in 1X AB.

Apply biotin block - 15 min @ RT.

No wash, wipe excess block and apply primary antibody

7. Apply primary antibody (Marck's protein) at 1:100 dilution and incubate for one hour.

Lot# _____ Aliquoted yes / no Date Aliquoted _____

For the negative control slides, normalize the protein concentration of the Pre-immune mouse serum to the protein concentration of the primary antibody (Marck's protein) and use this to make the 1:100 dilution. Apply to the slides and incubate for one hour.

Lot# _____

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply secondary antibody (Biotinylated Horse anti-mouse) @ 1:800 dilution and incubate for 30 minutes.

Lot# _____ Reconstituted Date _____

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply Label antibody and incubate for 30 minutes. (Prepare 30 minutes prior to use)

2 drops Reagent A + 5 ml diluent -> Mix and then add 2 drops Reagent B

Exp. Date _____ New kit yes / no

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.

(Add 1 drop of DAB per ml of substrate) Exp. Date _____

New kit yes / no Lot# _____

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30-45 seconds.

16. Rinse in tap water until water is clear.

17. Place slides in 1X Automation Buffer for 1 minute with gentle agitation to blue slides.

18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% EtOH	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip

updated 1/22/2003