

# Detection of PCNA in Formalin-Fixed, Paraffin-Embedded Human Tissue

## Reagent and Antibody Information

[1X Wash Buffer](#)

[3% Hydrogen Peroxide](#)

[1% BSA Diluent](#)

Distilled Water

[DAB Chromagen](#)

[Hematoxylin](#)

[1% Dry Milk](#)

Blocking Solution: 1% non-fat dry milk is used to dilute the primary antibody. This serves as the blocking solution for this protocol.

Primary Antibody: Monoclonal Anti-Mouse PCNA (19A2) Antibody

Chemicon International, Inc.

Temecula, CA 92590

[www.chemicon.com](http://www.chemicon.com)

1-800-437-7500

Catalog # MAB4078

Secondary Antibody: Goat Anti-Mouse IgM mu (Chain-Specific Biotin Conjugated)

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Catalog # 115-065-020

Label Complex: Peroxidase-Conjugated Streptavidin SS Label

Biogenex Laboratories

San Ramon, CA 94583

[www.biogenex.com](http://www.biogenex.com)

1-800-421-4149

Catalog # HK330-9K

## Staining Procedure

Positive Control Tissue: Tissue with a high cell turnover rate (breast carcinoma)

Stain Localization: The localization of the stain is dependent upon the cell cycle stage. (Foley, J et al.)

G0 = no staining

G1 = nuclear staining, 1+ just above background

S = nuclear, intense, dark brown staining

G2 = nuclear and cytoplasmic, 2+ distinct brown staining

M = cytoplasmic, 2+ distinct granular brown staining

1. Deparaffinize and hydrate slides through the following solutions:

<b>Solution</b>	<b>Repetitions</b>	<b>Time</b>
Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Wash Buffer	2 times	5 minutes

2. Quench endogenous peroxidase by placing the slides in 3% hydrogen peroxide for 15 minutes.

3. Rinse slides in 2 changes of 1X Wash Buffer for 5 minutes each.

4. Heat-Induced Epitope Retrieval Using The Microwave

Place a full rack of slides into a Tissue Tek® container with 200 ml of distilled water

(Insert blank slides into any empty slots in the rack to ensure even heating of slides)

Microwave for 5 minutes at power level 5.

Cool for 1 minute. (Add more distilled, if necessary.)

Microwave again for 5 minutes at power level 5. *Temperature Before Cooling Slides* \_\_\_\_\_

Cool 20 minutes at room temperature.

Rinse the slides in 2 changes of distilled water for 3 minutes each time.

5. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

6. Prepare the primary antibody at a 1:1200 dilution, and incubate for 30 minutes at room temperature.

(The diluent should consist of equal amounts of 1% non-fat dry milk and 1% BSA diluent.)

Lot # \_\_\_\_\_ Date Aliquoted \_\_\_\_\_

7. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each time.

8. Apply the goat anti-mouse IgM secondary antibody at a 1:400 dilution, and incubate for 30 minutes at room temperature. (This reagent is made in BSA diluent only.)

Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

9. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

10. Apply the Streptavidin SS Label and incubate for 30 minutes at room temperature.

Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_

11. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

12. Apply the DAB chromagen, and incubate in the dark for 6 minutes at room temperature.  
(Add 1 drop of DAB per ml of substrate)

Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_ New Kit: yes / no

13. Rinse the slides in tap water 3 minutes.

14. Counterstain with Harris Hematoxylin for 20 seconds.

15. Rinse the slides in tap water until water is clear.

16. Gently agitate slides in 1X Wash Buffer until they turn blue.

17. Dehydrate through the following solutions:

<b>Solution</b>	<b>Repetitions</b>	<b>Time</b>
95% Ethanol	1 time	3 minutes
100% Ethanol	3 times	3 minutes
Xylene	2 times	5 minutes

18. Coverslip

*Updated 02/07/06*