NIOSH ALERT

Controlling Exposures to Nitrous Oxide During Anesthetic Administration

WARNING!

Workers exposed to nitrous oxide (N2O) may suffer harmful effects.

Take the following steps to protect yourself from waste N2O in the workplace:

- Be aware that N₂O may cause the following health effects:
 - Decreases in mental performance, audiovisual ability, and manual dexterity
 - · Adverse reproductive effects
- Make sure that the following monitoring procedures are performed when the anesthetic equipment is installed and every 3 months thereafter:
 - Leak testing of equipment
 - Monitoring of air in the worker's personal breathing zone
 - Environmental (room air) monitoring
- Prevent leakage from the anesthetic delivery system through proper maintenance and inspection of equipment. Eliminate or replace the following:
 - Loose-fitting connections
 - Loosely assembled or deformed slip joints and threaded connections
 - Defective or worn seals, gaskets, breathing bags, and hoses
- Control waste N₂O with a well-designed scavenging system that includes the following:

- Securely fitting masks
- Sufficient flow rates for the exhaust system
- Properly vented vacuum pumps
- 5. Make sure that the room ventilation is effectively removing waste N₂O. If concentrations of N₂O are above 25 ppm, take the following steps:
 - Increase the airflow into the room.
 - Use supplemental local ventilation to capture N₂O at the source.

For additional information, see NIOSH Alert: Request for Assistance in Controlling Exposures to Nitrous Oxide During Anesthetic Administration [DHHS (NIOSH) 94–100], or call 1–800–35–NIOSH. Single copies of the Alert are available free from the following:

Publications Dissemination, DSDTT National Institute for Occupational Safety and Health

4676 Columbia Parkway Cincinnati, OH 45226

Fax number: (513) 533-8573



U.S. Department of Health and Human Services
Public Health Service
Centers for Disease Control and Prevention
National Institute for Occupational Safety
and Health