SOP HW-34 Revision 1 August 2007

SOP NO. HW-34/Trace VOA USEPA Contract Laboratory Program Statement of Work for Organic Analysis of Trace Concentration of Volatile Organic Compounds SOM01.2 Data Validation




TABLE OF CONTENTS

INTRODUC	-
Sco	ope and Applicability
Sur	mmary
Dat	ta Qualifiers
Lal	o Qualifiers
Res	viewer Qualifications
PACKAGE	COMPLETENESS AND DELIVERABLES
1.	Chain of Custody and Sampling Trip Reports
2.	<pre>Data Completeness and Deliverables</pre>
3.	Cover Letter SDG Narrative
	<u>Data Validation Checklist</u>
PART A:	VOA ANALYSES
1.	Sample Conditions/Problems
	Holding Times
3.	Deuterated Monitoring Compound (DMC) Recovery (Form II)
	Matrix Spike/Matrix Spike Duplicate Recovery (Form III)1
	Method Blanks (Form IV)
	<u>Contamination</u>
7.	GC/MS Instrument Performance Check (Form V) 1
8.	Target Compound List (TCL) Analytes (Form I) 1
9.	Tentatively Identified Compounds (TIC)
10	. Compound Quantitation and Reported Detection Limits 2
11	. Standards Data (GC/MS)
	. Standards Data (GC/MS)
	. GC/MS Continuing Calibration Verification (CCV) (Form
	<u>VII)</u>
14	. Internal Standard (Form VIII)
	. <u>Field Duplicates</u>
	. <u></u>
De	finitions
Do	5

INTRODUCTION

Scope and Applicability

This SOP offers detailed guidance in evaluating laboratory data generated according to the method in the "USEPA Contract Laboratory Program Statement of Work for Organics Analysis Multi-Media, Multi-Concentration, SOM01.2, May 2005". The method is based on EPA Volatile Method 524.2. The validation procedures and actions discussed in this document are based on the requirements set forth in the "USEPA Contract Laboratory Program National Functional Guidelines for Superfund Organic Methods Data Review, January 2005". This document attempts to cover technical problems specific to trace concentration of volatile compounds. Situations may arise where data limitations must be assessed based on the reviewer's own professional judgement.

In addition to technical requirements, contractual requirements may also be covered in this document. While it is important that instances of contract non-compliance be addressed in the Data Assessment, the technical criteria are always used to qualify the analytical data.

Summary

To ensure a thorough evaluation of each result in a data case, the reviewer must complete the checklist within this SOP, answering specific questions while performing the prescribed "ACTIONS" in each section. Qualifiers (or flags) are applied to questionable or unusable results as instructed. The data qualifiers discussed in this document are as follows:

Data Qualifiers

- U The analyte was analyzed for, but was not detected above the reported sample quantitation limit.
- J The analyte was positively identified; the associated numerical value is the approximate concentration of the analyte in the sample.
- N The analysis indicates the presence of an analyte for which there is presumptive evidence to make a "tentative identification."
- JN The analysis indicates the presence of an analyte that has been "tentatively identified" and the associated numerical value represents its approximate concentration.

- UJ The analyte was not detected above the reported sample quantitation limit. However, the reported quantitation limit is approximate and may or may not represent the actual limit of quantitation necessary to accurately and precisely measure the analyte in the sample.
- R The sample results are rejected due to serious deficiencies in the ability to analyze the sample and meet quality control criteria. The presence or absence of the analyte cannot be verified.

Lab Qualifiers:

- D The positive value is the result of an analysis at a secondary dilution factor.
- B The analyte is present in the associated method blank as well as in the sample. This qualifier has a different meaning when validating inorganic data.
- E The concentration of this analyte exceeds the calibration range of the instrument.

The reviewer must prepare a detailed data assessment to be submitted along with the completed SOP checklist. The Data Assessment must list all data qualifications, reasons for qualifications, instances of missing data and contract noncompliance.

Reviewer Qualifications:

Data reviewers must possess a working knowledge of the USEPA Statement of Work SOM01.1 and National Functional Guidelines mentioned above.

USEPA Region Method: CLP/		Date: 2 SOP HW-34			
		7	YES	NO	N/A
	PACKAGE COMPLETENESS AND DELIVERABLE	ls			
CASE NUMBER	C: LAB:				
SITE NAME:	SDG No(s).:				
1.0 Chain of	Custody and Sampling Trip Reports				
	re the Traffic Reports/Chain-of-Custody Rec resent for all samples?	ords -	<u>[]</u>		
ACTION:	If no, contact RSCC, or the TOPO to obtaireplacement of missing or illegible copie from the lab.				
	s the Sampling Trip Report present for all amples?				
ACTION:	If no, contact either RSCC or ask the TOPO obtain the necessary information from the contractor.				
2.0 Data Com	pleteness and Deliverables				
	ave any missing deliverables been received nd added to the data package?			<u>[]</u>	
ACTION:	Contact the TOPO to obtain an explanation resubmittal of any missing deliverables from If lab cannot provide them, note the effect review of the data package in the Contract Problems/Non-compliance section of the Data Assessment.	com the I ct on the			
	as CLASS CCS checklist included with the ackage?				

USEPA Reg Method: C		Date: Augu P HW-34, 1		
		YES	NO	N/A
2.3	Are there any discrepancies between the Traff Reports/Chain-of-Custody Records, Sampling Tr Report and Sample Tags?		_ [_]	
ACTI	CON: If yes, contact the TOPO to obtain an explain resubmittal of any missing deliverables from laboratory.			
3.0 Cover	Letter SDG Narrative			
3.1	Is the SDG Narrative or Cover Letter Present?	, []		
3.2	Are case number, SDG number and contract number contained in the SDG Narrative or cover letter (see SOW, Exhibit B, section 2.5.1)? EPA sample numbers in the SDG, detailed documentation of any quality control, sample, shipment, and/or analytical problems encounted in processing the samples? Corrective action taken?	er		
3.3	Does the Narrative contain description of coland trap used(see SOM, page B-12, section 2.5			
3.4	Does the narrative, VOA section, contain a li of all TICs identified as alkanes and their estimated concentrations?	.st <u>[_]</u>		
3.5	Did the contractor record the temperature of cooler on the Form DC-1, Item 9 - Cooler Temperature, and in the SDG Narrative?	the [_]		
3.6	Does the narrative contain a list of the pH values determined for each water sample submifor volatiles analysis (SOW, page B-13, secti 2.5.1.2)?			
	Does the Case Narrative contain the "verbatim" statement (page B-12, section 2.5.1 of the SOM)	ý []		

USEPA Region Method: CLP		Date:	_		
			YES	NO	N/A
the unav	"No", to any question in this section, cont TOPO to obtain necessary resubmittals. If vailable, document under the Contract Probl n-Compliance section of the Data Assessment	ems/			
4.0 Data Val	lidation Checklist				
	Check the package for the following (see SO requirements, section 2.1, page B-10):	M report	ing		
ā	a. Is the package paginated in ascending or starting from the SDG narrative?	der			
k	o. Are all forms and copies legible?		[]		
C	c. Assembled in the order set forth in the	SOW?	[_]		
C	d. Trace Concentration Volatiles Data prese	nt?	[]		
	PART A: Trace VOA ANALYSES	_			
1.0 Sample C	Conditions/Problems				
<u> </u>	Oo the Traffic Reports/Chain-of-Custody Rec Sampling Trip Report or Lab Narrative indic any problems with sample receipt, condition samples, analytical problems or special circumstances affecting the quality of the	ate of			
ACTION:	If samples were not iced or the ice was marrival at the laboratory and the tempera cooler was > 10°C, then flag all positive with a "J" and all non-detects "UJ".	ture of	the		
ACTION:	If both VOA vials for a sample have air b VOA vial analyzed had air bubbles, flag a results "J" and all non-detects "R".			è	

2.0 Holding Times

				YES NO
2.1 Hav	e any VOA technio	cal holding ti	mes. determined	l
fro	m date of collect eeded?	_		
water sam	Holding Times: T	from sample co	llection provid	led that samp
_	preserved to pH 2 Review the SDG		_	
	and arrived at			_
_	ication in the SI	=	- -	
	hat there was a p	-	- '	-
=	an be assumed to	-	-	-
	$\alpha \kappa Lii$ $\alpha \alpha \alpha L \alpha \alpha L \alpha Lii +$			
	erly cooled, but um holding time :			
the maxim	um holding time :	is 7 days from	sample collect	ion.
the maxim	-	is 7 days from SR, analysis d	sample collect	ion.
the maxim ACTION: L	um holding time :	is 7 days from SR, analysis d	sample collect	ion.
the maxim ACTION: L	um holding time ist sampling, VTS or samples which elow.	is 7 days from SR, analysis d	sample collect ates and preser g time in the t	ion.
the maxim ACTION: L	um holding time ist sampling, VTS or samples which elow.	is 7 days from SR, analysis d missed holdin	sample collect ates and preser g time in the t e Violations	ion.
the maxim ACTION: L fo	um holding time ist sampling, VTS or samples which elow.	is 7 days from SR, analysis d missed holdin of Holding Tim	sample collect ates and preser g time in the t e Violations	ion.
the maxim ACTION: L fo	um holding time ist sampling, VTS or samples which elow. Table of See (is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo	sample collect ates and preser g time in the t e Violations dy Records)	vation able
the maxim ACTION: L for book Sample	um holding time ist sampling, VTS or samples which elow. Table of See of Was Sample	is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo Date	sample collect ates and preser g time in the t e Violations dy Records) Date Lab	vation. able Date
the maxim	um holding time ist sampling, VTS or samples which elow. Table of See of Was Sample	is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo Date	sample collect ates and preser g time in the t e Violations dy Records) Date Lab	vation. able Date
the maxim	um holding time ist sampling, VTS or samples which elow. Table of See of Was Sample	is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo Date	sample collect ates and preser g time in the t e Violations dy Records) Date Lab	vation. able Date
the maxim	um holding time ist sampling, VTS or samples which elow. Table of See of Was Sample	is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo Date	sample collect ates and preser g time in the t e Violations dy Records) Date Lab	vation. able Date
the maxim ACTION: L for book Sample	um holding time ist sampling, VTS or samples which elow. Table of See of Was Sample	is 7 days from SR, analysis d missed holdin of Holding Tim Chain-of-Custo Date	sample collect ates and preser g time in the t e Violations dy Records) Date Lab	vation. able Date

ACTION: Qualify sample results using preservation and technical holding time information as follows:

a. If there is no evidence that the samples were properly preserved (acid and ice), but were analyzed within the technical holding time (7 days from sample collection), no qualification of the data is required.

	egion II CLP/SOW,	SOM01.2/Trace Volatiles		ite: Augu HW-34, R		
				YES	NO	N/A
	pr ou co	there is no evidence that the same served (acid and ice), and the same side of the technical holding the collection), qualify detects for all all non-detects "R".	samples we ime (7 day	ere analy s from s	zed ample	
	the s	the samples were properly presentantly the tamples were analyzed within the table 4 days from sample collection), retails is required.	technical	holding	time	
	we fro	the samples were properly presentere analyzed outside of the technic om sample collection), qualify det 2".	ical holdi	ng time	(14 d	ays
3.0 <u>De</u> t	terated M	Monitoring Compound (DMC) Recovery	y (Form II	<u>:)</u>		
3	1 Are t	the Volatile DMC Recovery Summarient?	es (Form I	:I [_]		
A	fro una	ntact the TOPO to obtain an explan om the lab. If missing deliverable available, document the effect in sessment.	les are	submittal		
3	2 Were	outliers marked correctly with ar	n asterisk	:? []		
AG	TION: Cir	ccle all outliers in red.				
3	Deute	more than three of the fourteen erated Monitoring Compounds (DMC's reries outside their corresponding	3)			
	If ye	es, were samples re-analyzed?				
	Were	method blanks re-analyzed?		<u>[]</u>		
A	bel	any DMC is outside the required low), qualify their associated tar				

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

VOLATILE DMC AND THEIR ASSOCIATED TARGET COMPOUNDS

Chloroethane-d5	1,2-Dichloropropane-d6	1,2-Dichlorobenzene-d4
Dichlorodifluoromethane Chloromethane Bromomethane Chloroethane Carbon Disulfide	Cyclohexane Methylcyclohexane 1,2-Dichloropropane Bromodichloromethane	Chlorobenzene 1,3-Dichlorobenzene 1,4-Dichlorobenzene 1,2-Dichlorobenzene 1,2,4-Trichlorobenzene 1,2,3-Trichlorobenzene
	trans-1,3- Dichloropropene-d4 cis-1,3-Dichloropropene trans-1,3- Dichloropropene 1,1,2-Trichloroethane	Chloroform-d 1,1-Dichloroethane Bromochloromethane Chloroform Dibromochloromethane Bromoform
2-Butanone-d5 Acetone 2-butanone	1,1-dichloroethene-d2 1,1-dichloroethene trans-1,2- Dichloroethene cis-1,2-Dichloroethene	2-Hexanone-d5 4-Methyl-2-pentanone 2-Hexanone
Vinyl Chloride-d3 Vinyl Chloride	Benzene-d6 Benzene	1,1,2,2- Tetrachloroethane- d2 1,1,2,2- Tetrachloroethane 1,2-Dibromo-3- chloropropane

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

1,2-Dichloroethane-d4	Toluene-d8
Trichlorofluoromethane 1,1,2-Trichloro-1,2,2- trifluoroethane Methyl Acetate Methylene Chloride Methyl tert-Butyl Ether Carbon Tetrachloride 1,2-Dichloroethane 1,1,1-Trichloroethane 1,2-Dibromoethane	Trichloroethene Toluene Tetrachloroethene Ethylbenzene o-Xylenes m,p-Xylene Styrene Isopropylbenzene

VOLATILE DEUTERATED MONITORING COMPOUND RECOVERY LIMITS

DMC	%RECOVERY LIMITS	DMC	%RECOVERY LIMITS
Vinyl Chloride-d3	65-131	1,2- Dichloropropane- d6	79-124
Chloroethane-d5	71-131	Toluene-d8	77-121
DMC	%RECOVERY LIMITS	DMC	%RECOVERY LIMITS
1,1- Dichloroethene-d2	55-104	trans-1,3- Dichloropropane-d4	73-121
2-Butanone-d5	49-155	2-Hexanone-d5	28-135
Chloroform-d	78-121		
1,2- Dichloroethane-d4	78-129	1,1,2,2- Tetrachloroethane-d2	73-125
Benzene-d6	77-124	1,2- Dichlorobenzene-d4	80-131

USEPA Region II		Date: August 2007
Method: CLP/SOW,	SOM01.2/Trace Volatiles	SOP HW-34, Revision 1

YES NO N/A

- 1. For any recovery greater than the upper limit:
 - a. Qualify "J" all positive associated target compounds.
 - b. Do not qualify associated non-detects.
- 2. For any recovery greater than or equal to 20%, but less than the lower limit:
 - a. Qualify "J" all positive associated target compounds.
 - b. Qualify "UJ" associated non-detects.
- 3. For any recovery less than 20%:
 - a. Qualify "J" all positive associated target compounds.
 - b. Qualify "R" all associated non-detects.
- NOTE: Up to three (3) DMC's per sample, and SIM analysis may fail to meet the recovery limits. (SOM, sec. 11.4.4, pg. D-36/Trace VOA).

As per SOM, any sample which has more than 3 DMC's outside the limits, it must be reanalyzed (sec. 11.5.3 pg. D-37/Trace VOA).

ACTION: Note in the Data Assessment under Contract Problems/ Non-Compliance if the Lab did not perform reanalysis.

- 3.4 Are there any transcription/calculation errors between raw data and form II? ____ [] ____
- ACTION: If large errors exist, ask the TOPO to obtain an explanation/resubmittal from the lab, make any necessary corrections and note errors in the data assessment.
- Note: DMC recovery limits criteria and qualification apply to samples diluted 5X and less. For samples diluted greater than 5X, recovery criteria does not apply Because it is assumed DMC is diluted below the quantitation range.

4.0 Matrix Spike/Matrix Spike Duplicate Recovery (Form III)

USEPA Region II Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

Note: Data for MS/MSD will not be present unless requested.

- 4.1 Are the MS/MSD Recovery Forms (Form III

 Trace VOA) present? [] _____
- 4.2 Was the MS/MSD analyzed at the required frequency (once per SDG, or every 20 samples, whichever is more frequent)?

 [] ____

ACTION: If any MS/MSD data are missing, take action as specified in section 3.1 above.

ACTION: No action is taken on MS/MSD data <u>alone</u>. However, using professional judgement, the validator may use the MS and MSD results in conjunction with other QC criteria and determine the need for some qualification of the data. If any MS/MSD % recovery or RPD is out of specification, qualify data to include the consideration of the existence of interference in the raw data. Consideration include, but not limited to the following "Action":

	Action		
Criteria	Detected Spiked Compounds	Non-detected Spiked Compounds	
%R or RPD > Upper acceptance Limits	J	No qualification	
20% < %R < Lower Acceptance Limits	J	UJ	
%R < 20%	J	Use Professional Judgement	
Lower Acceptance Limit < %R; RPD < Upper Acceptance Limit	No quali	fication	

5.0 Method Blanks (Form IV)

5.1	Is the Volatile	Method Blank	Summary (Fo	rm IV	
	Trace VOA) pres	ent?		1	<u> </u>

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1 YES NO N/A5.2 Frequency of Analysis: For the analysis of Trace Concentration VOA TCL compounds, has a method blank been analyzed for each SDG or every 20 samples, whichever is more frequent? [] 5.3 Has a VOA method blank been analyzed after the calibration standards and once every 12 hours time period for each GC/MS instrument used? 5.4 Was a VOA instrument blank analyzed after each sample/dilution that contains a target compound exceeding the initial calibration range (see SOM, page D-39/Trace VOA, section 12.1.1.3)? ACTION: If any method/instrument blank data are missing, notify the TOPO to obtain resubmittals or an explanation from the lab. If method blank data are unavailable, the reviewer may use professional judgement, or substitute field blank or trip blank data for missing method blank data. If an instrument blank was not analyzed after a sample containing a target analyte exceeding the initial calibration standards, inspect the sample chromatogram acquired immediately after this sample for possible carryover. The system is considered uncontaminated if the target analyte is below CRQL. Use professional judgement to determine if carryover occurred and qualify analyte(s) accordingly. 5.5 Was a storage blank analyzed once per SDG after all the samples were analyzed? [] ACTION: If storage blank data is missing, contact the TOPO to obtain any missing deliverables from the laboratory. If unavailable, note in the Contract Problems/Non-Compliance section of the Data Assessment. 5.6 The validator should verify that the correct identification scheme for EPA blanks was used. (See SOM page B-39, section 3.3.7.3 for more information.)

d: CL			∍: Augus √-34, Re		
			YES	NO	N/A
	Was the correct identification scheme used all Trace VOA blanks?	for	1_1		
ACTIO	N: Contact the TOPO to obtain corrections from or make the necessary corrections. Docume "Contract Problems/Non-Compliance section Assessment all corrections made by the value."	ent ir of th	n the ne Data		
5.7	<pre>Chromatography: review the blank raw data - (RICs), quant. reports, data system printout</pre>		_		
	Also compare the storage blank raw data with blank. Determine if contamination in the stalso present in the method blank.			is	
	Is the chromatographic performance (baseline stability) for each instrument acceptable for Trace VOAs?				
ACTIO	N: Use professional judgement to determine the data.	he eff	Eect on		
5.8	Are all detected hits for target compounds method, and storage blanks less than the CRO		<u>[]</u>		
	<pre>Exception: Methylene Chloride, Acetone and 2 be less than 2X times their respective CRQLs</pre>		anone mu	ıst	
ACTIOI	N: If no, an explanation and laboratory's conactions must be addressed in the case narrothe narrative contains no explanation, the note in the Contract Problems/Non-Complian of the Data Assessment.	rative en mak	e. If ke a		

6.0 Contamination

NOTE: "Water blanks", "drill blanks", and "distilled water blanks" are validated like any other sample, and are not used to qualify data. Do not confuse them with the other QC blanks discussed below.

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1 YES NO N/A6.1 Does the storage blank contain positive results (TCL and/or TICs) for Trace Concentration VOAs? 6.2 Do any method/reagent/instrument blanks contain positive results (including TICs) for Trace Concentration VOAs? <u>[]</u> NOTE: Contaminated instrument blanks are unacceptable under this SOW (see page D-41/Trace VOA, section 12.1.6.3). ACTION: Document in the Data Assessment under Contract Problems/Non-Compliance if a contaminated instrument blank was submitted. ACTION: Sample analysis results after the high concentration sample must be evaluated for carryover. Sample must meet the maximum carryover criteria as listed in SOM sec. 11.4.8.1, p. D-37/VOA.("the sample must not contain a concentration above the CROL for the target compounds that exceeded the limit in the contaminated sample.") 6.3 Do any field/trip/rinse blanks have positive Trace Concentration VOA results (including TICs)? [] ACTION: Prepare a list of the samples associated with each of the contaminated blanks. (Attach a separate sheet.) NOTE: All field blank results associated with a particular group of samples (may exceed one per case) must be used to qualify data. Trip blanks are used to qualify only those samples with which they were shipped. Blanks may not be qualified because of contamination in another blank. Field blanks & trip blanks must be qualified for system monitoring compound, instrument performance criteria, spectral or calibration QC problems.

TCL results due to contamination. Use the largest

ACTION: Follow the directions in the table below to qualify

STANDARD OPERATING PROCEDURE

USEPA Region II Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

value from all the associated blanks. If any blanks are grossly contaminated, all associated sample data should be qualified unusable (R).

Blank Type	Blank Result	Sample Result	Action for Samples
	Detects	Not detected	No qualification required
		< CRQL*	Report CRQL value with a U
	< CRQL *	<pre> CRQL and <2x the CRQL **</pre>	Report concentration of sample with a U
		≥ 2X CRQL **	No qualification required
	= CRQL *	< CRQL*	Report CRQL value with a U
Method, Field,		≥ CRQL*	No qualification required
Trip, Storage,		< CRQL*	Report CRQL value with a U
Instrument ***	> CRQL *	<pre></pre>	Report for sample concentration with a U
		<pre></pre>	No qualification required
	Gross contamination	Detects	Qualify results as unusable R
	TIC > 2ug/L	Detects	See "Action" below

^{* 2}x the CRQL for methylene chloride, 2-butanone and acetone

NOTE: Analytes qualified "U" for blank contamination are treated as "hits" when qualifying for calibration criteria.

Note: When applied as described in the table above, the contaminant concentration in the blank are multiplied by the sample dilution factor.

^{** 4}x the CRQL for methylene chloride, 2-butanone and acetone

^{***} Qualifications based on instrument blank results affect only the sample analyzed immediately after the sample that has target compounds that exceed the calibration range or non-target compounds that exceed 100 ug/L.

USEPA Region II Method: CLP/SOW, SOM01.2/Trace Volatiles	Date: August 2007 SOP HW-34, Revision 1
	YES NO N/A
Note: Gross contamination: greater than 2x the C the CRQL for methylene chloride, 2-butanon	_
ACTION: For TIC compounds, if the concentration sample is less than five times the conc the most contaminated associated blank, analyte "R" (unusable).	entration in
6.4 Are there field/rinse/equipment blanks ass with every sample?	ociated <u>[]</u>
ACTION: Note in data assessment that there is no field/rinse/equipment blank.	associated
<pre>Exception: samples taken from a drinking not have associated field blanks. 7.0 GC/MS Instrument Performance Check (Form V)</pre>	water tap do
7.1 Are the GC/MS Instrument Performance Check (Form V) present for Bromofluorobenzene (B	
7.2 Are the enhanced bar graph spectrum and mass/charge (m/z) listing for the BFB prov for each twelve hour shift?	ided <u>[]</u>
7.3 Did the 12-hour clock begin with either th injection of BFB, or in cases where a clos continuing calibration (CCV) was used as a opening CCV?	ing
Listed below are some, but not necessarily all, analytical sequences incorporating the use of th Use these examples as a guide for possible analy can be expected.	e opening/closing CCV.
Conditions for When Acceptable Criteria No Example Sequence is Appropriate:	otes:

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

If time remains on the 12 hour clock after initial calibration sequence	 BFB tunes meet instrument performance criteria. The five initial calibration standards meet initial calibration criteria. CCV A meets both opening and closing CCV criteria CCV B meets closing CCV criteria. 	The requirement of starting the new 12-hr clock for Analytical Sequence 2 with a new BFB tune is waived if CCV A meets opening CCV criteria. If CCV B meets opening CCV criteria, a method blank and subsequent samples may be analyzed immediately after CCV B.
If time remains on the 12 hour clock after initial calibration sequence	 BFB tunes meet instrument performance criteria. The five initial calibration standards meet initial calibration criteria. CCV A meets closing CCV criteria (but does not meet opening CCV criteria). CCV B meets opening CCV criteria. CCV C meets closing CCV Criteria. 	CCV A does not meet opening criteria, therefore a new BFB tune must be performed, immediately followed by CCV B before a method blank and any samples may be analyzed. In this case, the new 12 hr clock and Analytical Sequence 2 begins with the injection of the new BFB tune.
If more than 12 hrs have elapsed since the most recent initial calibration or closing CCV. OR If the most recent closing CCV was not or could not be used as an opening CCV.	 BFB tunes meet instrument performance criteria. CCV A meets opening CCV criteria. CCV B meets both opening and closing CCV criteria. CCV C meets both opening and closing CCV criteria. 	The requirement of starting the new 12 hour clock for Analytical Sequence 2 with a new BFB tune is waived if CCV B meets opening CCV criteria. If CCV C meets opening CCV criteria, a method blank and subsequent samples may be analyzed immediately after CCV B.

USEPA Region II Method: CLP/SOW, SOM01.2/Trace Volatiles			Date: August 2007 SOP HW-34, Revision 1
			YES NO N/A
elapsed since the most recent initial calibration or closing CCV OR If the most recent closing CCV was not or could not be used as an		 BFB tunes meet instrument performance criteria. CCV A meets opening CCV criteria. CCV B meets closing CCV criteria (but does not meet opening CCV criteria). CCV C meets opening CCV Criteria. CCV D meets both opening and closing CCV criteria. 	CCV B does not meet opening CCV criteria, therefore a new BFB tune must be performed, immediately followed by CCV B before a method blank and any samples may be analyzed. In this case, the new 12 hr clock and Analytical Sequence 2 begins with the injection of the new BFB tune. The requirement of starting the new 12 hr clock for Analytical Sequence 3 with a new BFB tune is waived if CCV D meets opening CCV criteria. If CCV D meets opening criteria, a method blank and subsequent samples may be analyzed after CCV B.
	All ion abunda	bundances been normalized to m/z since ratios must be normalized to more though the ion abundance of m/z	m/z 95, the nominal
ACTION:	If mass ass: unusable (R	ignment is in error, qualify all a	ssociated data as
7.5	Have the ion a instrument use	bundance criteria been met for eadd?	ch <u>Ll</u>
ACTION:	List all dat separate she	ca which do not meet ion abundance	criteria (attach a
ACTION:	Judgement ma	dance criteria are not met, profes by be applied to determine to what y be utilized.	
7.6	mass lists and	transcription/calculation errors be Form Vs? (Check at least two value)	
7.7	relative abund	of significant figures for the repances consistent with the number once criteria column on Form V?	

USEPA Regi Method: CL	Date: August 2007 SOP HW-34, Revision 1			
	P/SOW, SOM01.2/Trace Volatiles	YES		N/A
ACTION	: If large errors exist, take action as specified i above.	in section 3.1		
7.8	Is the spectrum of the mass calibration compound acceptable?	<u>[]</u>		
ACTION	: Use professional judgement to determine whether a should be accepted, qualified, or rejected.	associated dat	a	
8.0 <u>Target C</u>	ompound List (TCL) Analytes (Form I)			
8.1	Are the Organic Analysis Data Sheets (Form I) prese header information on each page, for each of the fo		red	
	a. Samples and/or fractions as appropriate?	[_]		
	b. Regional Control/MS/MSD samples?	<u>[]</u>		
	c. Blanks (method, trip, etc)?	П		
8.2	Are the VOA Reconstructed Ion Chromatograms, the maidentified compounds, and the data system printouts included in the sample package for each of the following	(Quant Report		
	a. Samples and/or fractions as appropriate?	[]		
	b. Regional Control/MS/MSD samples?	[]		
	c. Blanks (method, trip, etc)?	[]		
ACTION	: If any data are missing, take action specified in	n 3.1 above.		
8. 3	Is chromatographic performance acceptable with resp	ect to:		
	Baseline stability?			
	Resolution?	[_]		
	Peak shape?	<u>[]</u>		
	Full-scale graph (attenuation)?	<u>[]</u>		
	Other:?	[]		

ACTION: Use professional judgement to determine the acceptability of the data.

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1 YES NO N/A8.4 Are lab-generated standard mass spectra of the identified VOA compounds present for each sample? [] If any mass spectra are missing, take action as specified in 3.1 ACTION: above. If lab does not generate their own standard spectra, make note under the "Contract Problems/Non-Compliance" section of the Data Assessment. If spectra are unavailable reject "R" the reported results. Is the RRT of each reported compound within \pm 0.06 RRT 8.5 units of the standard RRT in the continuing calibration? [] 8.6 Are all ions present in the standard mass spectrum at a relative intensity greater than 10% also present in the sample mass spectrum? 8.7 Do sample and standard relative ion intensities agree to within ± 20%? [] ACTION: Use professional judgement to determine acceptability of data. If it is determined that incorrect identifications were made, all such data should be rejected (R) or changed to non-detected (U) at the calculated detection limit. In order to be positively identified, the data must comply with the criteria listed in sections 8.4-8.7 above. ACTION: When sample carry-over is suspected, review section 6.2/Action #2 above before determining if instrument cross-contamination has affected positive compound identifications. 9.0 Tentatively Identified Compounds (TIC) 9.1 Are all Tentatively Identified Compound Forms (Form I VOA-TIC) present? Do listed TICs include scan number or retention time, as well as the estimated "J" and/or "JN" qualifier? [] 9.2 Are the mass spectra for the tentatively identified compounds and associated "best match" spectra included in the sample package for each of the following: a. Samples and/or fractions as appropriate? b. Blanks? b. Are Alkanes listed in/or part of the Case Narrative?

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1 YES NO N/AACTION: If any TIC data are missing, take action specified in 3.1 above. ACTION: Verify "JN" qualifier is present for all chemically named TICs having a percent match of greater than or equal 85%. labeled "unknown" are qualified with a "J" qualifier. 9.3 Are any target compounds (from any fraction) listed as TICs? (Example: 1,2-dimethylbenzene is xylene - a VOA target analyte - and should not be reported as a TIC.) [] _ Flag with "R" only target compound detected in another fraction. (except blank contamination) 9.4 Are all ions present in the reference mass spectrum with a relative intensity greater than 10% also present in the sample mass spectrum? 9.5 Do TICs and "best match" reference spectra relative ion intensities agree within \pm 20%? Use professional judgement to determine the acceptability of TIC identifications. If it is determined that an incorrect identification was made, change its identification to "unknown" or to some less specific identification (example: "C3 substituted benzene") as appropriate. Action: When a compound is not found in any blank, but is detected in a sample and is a suspected artifact of a common laboratory contaminant, solvent preservatives or Aldo condensation, the result should be qualified as unusable (R). (i.e., common lab contaminants such as $CO_2(m/e 44)$, Siloxanes (m/e 73), diethyl ether, hexane, certain freons. Aldol condensation products: 4-hydroxy-4-methyl-2-pentanone, 4-methyl-2-penten-2one and 5,5-dimethyl-2(H)-furanone. Solvent preservatives: cyclohexene, and related by-products: cyclohexanone, cyclohexanol, cyclohexenone, chlorocyclohexene, and chlorocyclohexanol.). 10.0 Compound Quantitation and Reported Detection Limits 10.1 Are there any transcription/calculation errors in Form I results? (Check at least two positive values. Verify that the correct internal standards, quantitation ions, and RRFs were used to calculate Form I results.) <u>[]</u> __ [] 10.2 Are the CRQLs adjusted to reflect sample dilutions? ACTION: If errors are large, take action as specified in section 3.1

above.

	PA Regionod: CLI	on II P/SOW, SOM01.2/Trace Volatiles	Date: A			
			Y	ES	NO	N/A
	ACTION:	When a sample is analyzed at more than one diluting CRQLs are used (unless a QC exceedance dictates thigher CRQLs data from the diluted sample). Replaconcentrations that exceed the calibration range analysis by crossing out the "E" and its correspontate original Form I and substituting the data from sample. Specify which Form I is to be used, then across the entire page of all Form I's not to be any in the data summary package.	he use of the ace in the original noting value methe diluth draw a red	he inal on ed "X"		
11.0	Standard	ds Data (GC/MS)				
	11.1	Are the reconstructed ion chromatograms, and data syprintouts (quant. reports) present for each initial continuing calibration?	l and	<u>l</u> _		
	ACTION:	If any calibration standard data are missing, tak specified in section 3.1 above.	e action			
12.0	GC/MS Ir	nitial Calibration (Form VI)				
	12.1	Are the Initial Calibration Forms (Form VI LCV) present and complete for the volatile fraction at concentrate of 0.5, 1, 5, 10, and 25 $\mu g/\ell$ for non-ketones, 5, 10 100, and 200 ug/L for ketones.	tions	<u>l</u> _		
	Note:	The initial calibration standards for by Selected technique are 0.05, 0.1, 0.5, 1.0, and 2.0 ug/L.	. Ion Monito	ring	(SIM)	
	ACTION:	If any Initial Calibration forms are missing, tak specified in section 3.1 above.	e action as			
	12.2	Are the relative standard deviation (RSD) stable for over the concentration range of the calibration (i.e $RSD \le 30\%$, $\le 40\%$ for poor performers (see table below)	e.,	<u>l</u> _		
	ACTION:	Circle all outliers in red.				
	NOTE:	The twenty two (22) poor performers compounds and as	ssociated DN	MCs ar	re	

Volatile Compounds Exhibiting Poor Response

listed below. The relative response factor (RRF) for these compounds must

be greater than or equal to 0.010. All DMC must meet RRF \geq 0.010.

USEPA Region II	Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles	SOP HW-34, Revision 1
	YES NO N/A

Volatile Compounds		
Acetone	1,2-Dibromo-3-chloropropane	
2-Butanone	Isopropylbenzene	
Carbon disulfide	Methyl acetate	
Chloroethane	Methylene chloride	
Chloromethane	Methylcyclohexane	
Cyclohexane	Methyl tert-butyl ether	
1,4-Dioxane	trans-1,2-Dichloroethene	
1,2-Dibromoethane	4-Methyl-2-pentanone	
Dichlorodifluoromethane	2-Hexanone	
cis-1,2-dichloroethene	Trichlorofluoromethane	
1,2-Dichloropropane	1,1,2-Trichloro-1,2,2-trifluoroethane	

ACTION: If %RSD > 30.0%, (> 40.0% for the poor performers, qualify associated positive results for that analyte "J" (estimated). If %RSD is > 90, flag all non-detects for that analyte "R" (unusable) and positive hits "J".
NOTE: Analytes previously qualified "U" for blank contamination are still treated as "hits" when qualifying for initial calibration criteria.
12.3 Are any RRFs < 0.050 (< 0.010 for poor performers)? [_]
ACTION: Circle all outliers in red.
ACTION: If any $\overline{\text{RRF}}$ values are < 0.05 or < 0.01 for poor performers, qualify associated non-detects unusable (R) and associated positive results estimated (J).
ACTION: Document in the Data Assessment under Contract Problems/Non-Compliance the analytes that fail %RSD and/or RRF criteria.
12.4 Are there any transcription/calculation errors in the reporting of RRFs, RRFs or %RSD values? (Check at least 2 values, but if errors are found, check more.)
ACTION: Circle errors in red.

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1 YES NO N/AACTION: If errors are large, contact the TOPO to obtain an explanation/resubmittal from the lab, document in the Data Assessment under Contract Problems/Non-Compliance. 13.0 GC/MS Continuing Calibration Verification (CCV)(Form VII) 13.1 Are the Continuing Calibration Forms (Form VII) present and complete for the volatile fraction? [] Did the 12 hour clock begin with either the injection of 13.2 BFB or in cases where a closing CCV can be used as an opening CCV for each instrument? [] If any forms are missing or no continuing calibration standard ACTION: has been analyzed within twelve hours of every sample analysis, ask the TOPO to obtain explanation/resubmittal from the laboratory. If continuing calibration data are unavailable, flag all associated sample data as unusable (R). 13.3 Do any volatile compounds have a % Difference (% D) between the initial RRF and CCV RRF exceeding ± 50% for 1,4-Dioxane, ± 40% for the poor performers or ± 30% for the remaining compounds? [] ACTION: Circle all outliers in red. Do any volatile compounds have a RRF < 0.05 or < 0.01 for the poor performers? ACTION: Circle all outliers in red. Verify that the CCV was run at the required frequency (an opening and closing CCV must be run within 12-hour period) and the CCV was compared to the correct initial calibration. If the mid-point standard from the initial calibration is used as an opening CCV, verify that the result (RRF) of the mid-point standard was compared to the average RRF from the correct initial calibration. The closing CCV used to bracket the end of a 12-hour analytical sequence may be used as the opening CCV for the new 12-hour analyical sequence, provided that all the technical acceptance criteria are met for an opening CCV (see table below). If the closing CCV does not meet the technical acceptance criteria for an opening CCV, then a BFB tune followed by an opening CCV is required and the next 12-hour time period begins with the BFB tune. Use the following table to qualify data based on the technical Action: acceptance criteria for the opening CCV and closing CCV.

STANDARD OPERATING PROCEDURE .

USEPA Region II Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

YES NO N/A

Continuing Calibration Verification (CCV) Actions for Trace Volatiles Analyses

Criteria for	Criteria for	Action	
Opening CCV	Closing CCV	Detected Associated Compounds	Non-Detected Associated Compounds
RRF < 0.010 (poor responders) RRF < 0.050 (all other volatile target compounds)	RRF < 0.010 (for all volatile target compounds)	J	R
RRF \geq 0.010 (poor responders) RRF \geq 0.050 (for all other compounds)	RRF > 0.010 (for all target volatile compounds)	No	Action
%D > 40.0 or < -40.0 (poor responders) %D > 30.0 or < -30.0 (all other volatile target compounds)	%D > 50.0 or < -50.0 (for all volatile target compounds)	J	ŪJ
%D \leq 40.0 or \geq -40.0 (poor responders) %D \leq 30.0 or \geq -30.0 (all other volatile target compounds)	%D ≤ 50.0 or \geq -50.0 (for all volatile target compounds)	No	Action
Opening CCV not performed at required frequency *	Closing CCV not performed at required frequency *	R	

•	See	section	⊥3.2	apove

ACTION: Document in the Data Assessment under Contract Problems/Non-Compliance if more than two of the required analytes failed the above acceptance criteria.

13.5 Are there any transcription/calculation errors for the reporting of RRFs, or %D between initial RRFs and continuing RRFs? (Check at least two values but if errors are found, check more.)

ACTION: Circle errors with red pencil.

ACTION: If errors are large, notify the TOPO to obtain explanation/resubmittals from the lab. Document errors in the Contract Problems/Non-Compliance section of the Data Assessment.

Note: All DMCs must meet RRF \geq 0.010. No qualification of the data is necessary on the DMC RRF and %RSD/% Diff data alone. However, use professional judgment to evaluate the DMC RRF and %RSD/% Diff data in conjunction with the DMC recoveries to determine the need for qualification of data.

USEPA Region II Da			ite: August 2007			
Method: CI	P/SOW, SOM01.2/Trace Volatiles	OP HW-	34, R	evisi	on 1	
			YES	NO	N/A	
14.0 <u>Interna</u>	al Standard (Form VIII)					
14.1	Were the internal standard area counts for every same and blank within the range of 60.0% and 140.0% of its response in the most recent opening CCV standard calibration?	•	<u>[]</u>			
	If no, were affected sample reanalyzed?		[]			
ACTION	1: 1. Circle all outliers with red pencil.					
14.2	Are the retention times of the internal standards in sample or blanks within ±20 seconds from the RT of the internal standard in the 12-hour associated calibratistandard (opening CCV or mid-point standard from initiality calibration)?	ne ion	[]			

INTERNAL STANDARDS ACTIONS FOR TRACE VOLATILES

Action: Use the following table to qualify the data

	ACTION			
Criteria	Detected Associated Compounds *	Non-detected Associated Compounds *		
Area counts > 140% of 12-hour standard (opening CCV or mid-point standard from initial calibration)	J	No Action		
Area counts < 60% of 12-hour standard (opening CCV or mid-point standard from initial calibration)	J	R		
Area counts ≥ 60% but ≤ 140% of 12-hour standard (Opening CCV or mid-point standard from initial calibration)	No Action			
RT difference > 20.0 seconds between samples and 12-hour standard (Opening CCV or mid-point standard from initial calibration)	R **			
RT difference < 20.0 seconds between samples and 12-hour standard (Opening CCV or mid-point standard from initial calibration)	No Action			

For volatile compounds associated to each internal standard, see Table 3 - Trace Volatile Target Compounds and Deuterated Monitoring Compounds with Corresponding Internal Standards for Quantitation in SOM01.1, Exhibit D, available at:

Http://www.epa.gov/superfund/programs/clp/soml.htm

USEPA Region II Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

.

** Examine the chromatographic profile for that sample to determine if any false positives or negatives exist. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for that sample fraction. Detects should not need to be qualified as unusable "R" if the mass spectral are met.

NOTE: <u>Contract Requirements</u>: The SOM (section 11.5.1 page D-37/Trace VOA) states that any sample which fails the acceptance criteria for IS response must be reanalyzed.

ACTION: Document in the Data Assessment under Contract Problems/Non-Compliance any sample(s) which failed the above IS acceptance criteria.

15.0 Field Duplicates

15.1 Were any field duplicates submitted for Trace

Concentration VOA analysis?

[] ____ ____

ACTION: Compare the reported results for field duplicates and calculate the relative percent difference.

ACTION: Any gross variation between duplicate results must be addressed in the reviewer narrative. If large differences exist, contact the TOPO to confirm identification of field duplicates with the sampler.

USEPA Region II Date: August 2007
Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

.

Definitions

BFB - bromofluorobenzene CCS - contract compliance screening CLASS - Contract Laboratory Analytical Services Support CLP - Contract Laboratory Program CRQL - Contract Required Quantitation Limit GC/MS - gas chromatography/mass spectroscopy kg - kilogram uq - microgram ℓ - liter mℓ - milliliter QC - quality control RAS - Routine Analytical Services RIC - reconstructed ion chromatogram RPD - relative percent difference RRF - relative response factor RRF - average relative response factor (from initial calibration) RRT - relative retention time RSD - relative standard deviation RT - retention time RSCC - Regional Sample Control Center SDG - sample delivery group SOP - standard operating procedure SOW - Statement of Work TCL - Target Compound List TCLP - Toxicity Characteristics Leachate Procedure TIC - tentatively identified compound TPO - technical project officer VOA - volatile organic acid VTSR - validated time of sample receipt

TOPO - Task Order Project Officer

USEPA Region II Date: August 2007 Method: CLP/SOW, SOM01.2/Trace Volatiles SOP HW-34, Revision 1

method: CLP/SOW, SOMUL.2/ITate volatiles SOP HW-34, Revision

.

References

- 1. USEPA Contract Laboratory Program of Work for Organic Analysis Multi-Media, Multi-Concentration, SOW/CLPSOM01.1, October 2004
- 2. National Functional Guidelines for Superfund Organic Methods Data Review January 2005