

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

I. General Information

Device Generic Name: *In vitro* nucleic acid amplification test for the quantitation of Human Immunodeficiency Virus Type 1 (HIV-1) RNA in human plasma

Device Tradename: NucliSens[®] HIV-1 QT

Applicant's Name and Address: bioMérieux, Inc.
100 Akzo Avenue
Durham, North Carolina 27712

Premarket Approval Application (PMA) Number: BP010001

II. Intended use

The NucliSens[®] HIV-1 QT is an *in vitro* nucleic acid amplification test for the quantitation of Human Immunodeficiency Virus Type 1 (HIV-1) RNA in human plasma. The test can quantitate HIV-1 RNA over the range of 176 to 3.47X10⁶ copies/mL.

The test is intended for use in conjunction with clinical presentation and other laboratory markers of disease progression for prognostic assessment of HIV-1 infected patients, and for monitoring the effects of anti-retroviral therapy by serial measurements of plasma HIV-1 RNA for pediatric and adult patients with baseline viral loads greater than 93,000 and 28,000 copies of HIV-1 viral RNA/mL respectively.

The NucliSens[®] HIV-1 QT assay is not intended to be used as a screening test for HIV-1 nor is it to be used as a diagnostic test to confirm the presence of HIV-1 infection.

III. Device Description

(a) Test Principles

The NucliSens[®] HIV-1 QT assay comprises five separate stages as described below.

Nucleic acid release

The sample is added to NucliSens[®] Lysis Buffer containing guanidine thiocyanate and Triton X-100, which causes the lysis of viral particles and cells in the sample and the inactivation of RNases and DNases. Nucleic acid is thus released.

Nucleic Acid Isolation


Three synthetic RNAs (Qa, Qb, Qc) of known high, medium and low concentration, respectively, are added to the Lysis Buffer containing the released nucleic acid. These RNAs serve as internal calibrators, each differing from the HIV-1 wild-type (WT) RNA by only a short length of sequence. Under high salt conditions, all nucleic acid in the buffer, including the calibrators, binds to silicon dioxide particles.⁽³⁾ These particles, acting as the solid phase, are washed several times. Finally, the nucleic acid is eluted from the solid phase.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Nucleic Acid Amplification (see Figure 1 for illustration)

WT HIV-1 RNA present in the eluted nucleic acid is co-amplified along with the three internal calibrators. As illustrated in figure 1, amplification is based on repeated cycles of transcription reactions. Multiple copies of each WT and calibrator RNA target sequence are synthesized by T7-RNA polymerase by means of an intermediate DNA molecule that contains a double-stranded T7-RNA polymerase promoter. Each transcribed RNA molecule enters a new amplification cycle. The DNA intermediate is generated through a process that involves the binding of a primer to the RNA template, the extension of primer by AMV-RT (Avian Myeloblastosis Virus Reverse Transcriptase) to form an RNA-DNA duplex, the degradation of the RNA strand of the duplex by RNase H, the binding of a second primer to the remaining DNA strand and, finally, the extension of the second primer to form the double-stranded T7-RNA polymerase promoter needed for transcription. Once transcription is initiated, the RNA transcripts which are 'negatives' of the original RNA present in the sample will be subject to the same process, only in this case extension is not restricted to the second primer, since the extension product of the first primer will also be extended. The primers (one of which contains the sequence of the T7-RNA polymerase promoter) are complementary to two different parts of the HIV-1 RNA. This pair of primers defines, and allows the amplification of, a sequence within the HIV-1 gag region. Since the Nucleic Acid Sequence-Based Amplification (NASBA) process requires no strand separation, amplification is isothermal and continuous.

NucliSens[®] HIV-1 QT
 Summary of Safety and Effectiveness



Legend







- | | |
|---|---|
| A ss RNA (sense) |  sense RNA |
| B Primer 1 |  Primer 1 |
| C AMV-RT |  Primer 2 |
| D RNase H |  antisense DNA |
| E Primer 2 |  sense DNA |
| F T7 RNA pol |  antisense RNA |
| G ss RNA (antisense) amplificate | |
| H Isothermal amplification | |

Figure 1: Nucleic acid amplification principle

NucliSens[®] HIV-1 QT Summary of Safety and Effectiveness

Nucleic Acid Detection

The quantitation of HIV-1 RNA in a sample is based on the measurement of electrochemiluminescence (ECL) with the NucliSens[®] Reader. To differentiate among the amplicons (WT, Qa, Qb and Qc), aliquots of the amplified sample are added to four hybridization solutions, each specific for one type of amplicon. Here, the respective amplicons are hybridized with a bead-oligo (i.e. an oligo bound to streptavidin coated paramagnetic beads acting as the solid phase) and a ruthenium-labeled probe. The paramagnetic beads carrying the hybridized amplicon/probe complex are captured on the surface of an electrode by means of a magnet. Voltage applied to this electrode triggers the (ECL) reaction. The light emitted by the hybridized ruthenium-labeled probes is proportional to the amount of amplicons, which in turn are proportional to the HIV-RNA in the input samples. Calculation based on the relative amounts of the four amplicons gives an estimation of the amounts of HIV-1 RNA in the sample.

Quantitation of HIV-1 RNA

The NucliSens[®] HIV-1 QT system uses three internal calibrators for quantitation. These internal calibrator RNAs are produced *in vitro* from plasmid constructs encoding the wild-type *gag* region sequence. These three calibrator RNAs are identical to the wild-type sequence except for a 20 base sequence recognized specifically by the detector probe. In the case of calibrator RNAs (designated Qa, Qb, and Qc), the detector probe site was mutagenized to generate three distinct sequences, but similar in A, C, G, and T content to the wild type sequence. Thus, the calibrator RNAs are amplified with the same kinetics as the wild-type target, each reaction product is captured onto the magnetic bead with the same efficiency, and each calibrator reaction product can be distinguished from the other and from the wild-type product by the detector probes with distinct sequences. Because the A, C, G, and T content in each detector site is maintained, the melting temperature of all four probes is similar.

The internal calibrators for quantitation are introduced into the specimen immediately after the initial lysis step. Importantly, the three calibrators are included at known copy numbers (Qa ~ 10^5 , Qb ~ $10^{4.3}$, Qc ~ $10^{3.6}$). These calibrators are co-extracted and co-amplified with the wild-type nucleic acid in a single tube. The reaction product is then divided into four independent detection assays, each with one of the calibrator detector probes or with the wild-type detector probe. At the completion of the hybridization reactions, the four products are loaded into the NucliSens[®] Reader. The ECL signal for each of the hybridization reactions is determined and the input copy number for wild-type HIV-1 RNA, relative to the input quantities of calibrator RNAs is determined by calculating the ratio of ECL signals for WT to Qa, Qb, and Qc.

The actual computation of the input quantity of wild-type RNA is achieved through the application of a specific curve fitting program.⁽¹²⁾ Importantly, the inclusion of the calibrators in the analysis of the independent samples permits individual determinations of the acceptance of each wild-type calculation through a validated algorithm designed to evaluate the relationship of the calibrator signals to each other. The three calibrator signals (Qa, Qb, Qc) are corrected for background noise and checked against a number of parameters: upper detection limit, fixed minimum values, variance, and correlation. If the Qc calibrator is discarded as a result of these checks, the remaining two calibrators (Qa and Qb) are used to calculate the wild-type RNA concentration.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

(b) HIV-1 Target Description

The region of HIV-1 RNA targeted for amplification by the NucliSens[®] HIV-1 QT is a highly conserved region of the *gag* gene. The oligonucleotide primers were designed to amplify a region spanning nucleotides 1358 – 1499 of HXB2, AC K03455. Primer 1 is comprised of a T7 promoter sequence and a sequence complementary to nucleotides 1471 - 1499 of HIV-1. Primer 2 is comprised of a sense sequence corresponding to nucleotides 1358 - 1386 of HIV-1.

(c) HIV-1 RNA Quantitation

The relationship between the calibrator signals and the calibration input amounts is used for the determination of the unknown amount of Wild Type (WT) RNA in a sample. Within the linear range of the assay, the ECL signal for each target or standard is proportional to the amount of HIV-1 RNA, Qa RNA, Qb RNA, or Qc RNA present in the sample tube and therefore, proportional to the total amount of individual RNA present in the original sample. This relationship is given by the following equation:

$$\text{ECL signal} = \text{constant} \times \text{RNA input.}$$

Since the three calibrators are added to the sample prior to isolation, amplification, and detection the proportionality between each RNA and its corresponding signal is the same. So for each sample there are four equations:

$$\text{ECL signal Qa} = \text{constant} \times \text{Qa RNA input} \quad (1)$$

$$\text{ECL signal Qb} = \text{constant} \times \text{Qb RNA input} \quad (2)$$

$$\text{ECL signal Qc} = \text{constant} \times \text{Qc RNA input} \quad (3)$$

$$\text{ECL signal WT} = \text{constant} \times \text{WT RNA input} \quad (4)$$

As the ECL signal and the Q input amount are known for each level of the calibrators, three estimates of the constant can be determined. That is:

$$\text{Constant Qa} = \text{ECL signal Qa/Qa RNA input}$$

$$\text{Constant Qb} = \text{ECL signal Qb/Qb RNA input}$$

$$\text{Constant Qc} = \text{ECL signal Qc/Qc RNA input}$$

The average constant then is:

$$\frac{\text{Constant Qa} + \text{Constant Qb} + \text{Constant Qc}}{3}$$

This average constant can then be used to estimate the WT RNA input amount by the formula:

$$\text{WT RNA input} = \text{ECL signal WT/average constant.}$$

As mentioned above, each of the calibrator signals are measured against a number of parameters. This is to ensure that a calibrator that gives an unreasonable estimate of the constant, which can skew the calculation of the initial WT RNA input amount, is not included in the calculations. If two calibrators fail to meet the parameters above, the test is considered invalid and marked as such.

(d) Primers and Probes

NASBA amplification uses two primers to amplify the viral target as well as the three internal calibrators. The P1 primer has a length of 54 nucleotides, 29 nucleotides on the 3'-end code are complementary to the HIV sequence while the other nucleotides facilitate transcription by T7 RNA polymerase. Essential is a 20 nucleotide consensus sequence (5'-TAATACGAC-

NucliSens[®] HIV-1 QT Summary of Safety and Effectiveness

TCACTATAGGG) for the T7 RNA polymerase promoter. In addition, an extra sequence of 5 nucleotides (5'-AATTC) flanks the consensus sequence. This extra sequence serves as a "buffer" for the core sequence that is required for T7 RNA polymerase activity. Because of this design choice, partial length oligos n-1 through n-5 are expected to show equivalent functionality compared to the full length oligo since neither (1) the core consensus sequence for the T7 RNA polymerase nor (2) the 3'-end with the HIV specific sequence where chain elongation takes place is affected.

The P2 primer has a length of 29 nucleotides, which is completely HIV specific. A typical minimal oligo length needed for sufficient hybridization efficiency and to ensure sequence specificity is 20 nucleotides. Partial length oligos lack nucleotides on the 5'-end since DNA synthesis takes place in the 3' to 5' direction. This means that hybridization characteristics on the 3'-end of the oligo are not affected which is important since chain elongation during the amplification reaction takes place in the 5' to 3' direction. Therefore, it is expected that partial length oligos missing up to approximately 5 nucleotides will show equivalent performance as the full-length oligo.

The generic capture probe is an oligonucleotide comprised of a sense sequence that corresponds to nucleotides 1388 – 1411 of the HIV -1 RNA *gag* sequence. Four separate detection probes are used to detect the amount of amplified WT RNA and each of the calibrator RNAs present in the hybridized mixture. The detection probe for the WT is an oligonucleotide comprised of a sense sequence that corresponds to nucleotides 1419 – 1442 of the HIV -1 sequence. In order to have a unique detection probe for each of the calibrators a unique replacement is made in this position of the sequence for each calibrator RNA. Each replacement consists of a randomized sequence of the same length as present in the WT RNA sequence (24 nucleotides). Each sequence contains one less 'G' and one more 'A' than the WT sequence. The number of 'T' and 'C' bases is the same in all detection probes. Therefore, the efficiency of amplification for each of the calibrators is the same as that of the WT target.

(e) Quality Control

The integrity of each individual result can be monitored by reference to the performance of the three internal calibrators. However, it is recommended that a High Positive Control, Low Positive Control, and a Negative Control be included with the first run of each kit lot to verify product performance. Furthermore, the inclusion of a Low Positive and a Negative Control with each subsequent run of the same kit lot is also recommended. Refer to the control manufacturer's package insert for instructions.

As with any new laboratory procedure, new operators should consider the use of additional controls to establish a high degree of confidence in performing the assay.

If using the NucliSens[®] HIV-1 RNA Controls, the expected range for the Positive Controls is stated in the package insert for the controls. The copy number per ml for each Positive Control should fall within the range indicated in the package insert. The Negative Control should give a less than lower detection limit result. If using commercially available controls or internally prepared controls, a typical acceptance range is approximately ± 3 standard deviations from the established copy level. If controls are not as expected, the run should be re-evaluated and a determination made as to whether the sample results are acceptable or should be repeated.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

(f) Kit Components

Nucleic Acid Release

Component	Size	Volume	Composition
Lysis Buffer	50 Tubes	9.0 mL	5 mol/L guanidine thiocyanate, Triton X-100, TRIS/HCl.

Nucleic Acid Isolation

Components in each Isolation Reagent Module

Component	Size	Volume	Composition
Wash Buffer	5 tubes	22 mL	5 mol/L guanidine thiocyanate, TRIS/HCl.
Silica	5 tubes	0.8 mL	Hydrochloric acid-activated silicon dioxide particles.
Elution Buffer	5 tubes	1.5 mL	TRIS / HCl. Color code: white capped tube

Nucleic Acid Amplification

Components in each Amplification Reagent Module

Component	Size	Volume	Composition
Calibrators	5 tubes	6 mg	Lyophilized synthetic RNA (Qa, Qb, Qc) sphere, each tube contained in a foil pack with desiccant. Color code: yellow capped tube
Enzymes	5 tubes	6.5 mg	Lyophilized sphere containing AMV-RT, RNase H, T7-RNA polymerase and BSA, each tube contained in a foil pack with desiccant. Color code: red capped tube
Enzyme Diluent	5 tubes	0.5 mL	TRIS/HCl. Color code: red capped tube
Primers	5 tubes	10 mg	Lyophilized sphere with synthetic primers, nucleotides, and stabilizers; each tube contained in a foil pack with silica gel desiccant. Color code: blue capped tube
Primer Diluent	5 vials	0.5 mL	TRIS/HCl, 30% DMSO. Color code: blue capped tube

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Nucleic Acid Detection

Components in each Detection Reagent Module

Component	Size	Volume	Composition
Bead-Oligo	2 tubes	1.68 mL	DNA oligo bound to streptavidin-coated paramagnetic beads with preservative. Color code: pink capped tube
WT probe	1 vial	0.84 mL	Ruthenium-labeled DNA oligo with preservative. Color code: white capped tube
Qa probe	1 vial	0.84 mL	Ruthenium-labeled DNA oligo with preservative. Color code: red capped tube
Qb probe	1 vial	0.84 mL	Ruthenium-labeled DNA oligo with preservative. Color code: yellow capped tube
Qc probe	1 vial	0.84 mL	Ruthenium-labeled DNA oligo with preservative. Color code: blue capped tube
Detection Diluent	2 vials	15 mL	TRIS/HCl with preservative.
Instrument Reference Solution	1 vial	1.7 mL	Streptavidin-coated paramagnetic beads.

IV. Warnings and Limitations

- Monitoring the effects of anti-retroviral therapy by serial measurement of plasma HIV-1 RNA has only been validated for patients with baseline viral loads = 93,000 copies/mL for pediatric patients and = 28,000 copies/mL for adult patients.
- The performance of the NucliSens[®] HIV-1 QT assay has only been adequately validated with HIV-1 subtype B specimens.
- NucliSens[®] HIV-1 QT assay is generally insensitive for detection of HIV-2 RNA. However, samples from individuals infected with HIV-2 may exhibit cross-reactivity in this assay.
- When testing specimens with viral load < 300 copies/mL, the user should consider the use of well-characterized reference materials titrated from 150 copies/mL to 300 copies/mL, inclusive.

The NucliSens[®] HIV-1 QT is an *in vitro* nucleic acid amplification test for the quantitation of Human Immunodeficiency Virus Type 1 (HIV-1) RNA in human plasma. The test can quantitate HIV-1 RNA over the range of 176 to 3.47X10⁶ copies/mL. The NucliSens HIV-1 QT assay can accurately detect a 0.3 log₁₀ (2-fold) or greater change in HIV-1 RNA for patients whose viral load is between 3500 copies/mL and 3,500,000 copies/mL. The NucliSens HIV-1 QT can accurately detect a 0.9 log₁₀ (9-fold) change in HIV-1 RNA for patients whose viral load is between 180 and 850 copies/ml.

Use of this product should be limited to personnel who have been trained in all aspects and techniques of the NucliSens[®] HIV-1 QT assay. Attention should be paid to adequate specimen collection, transporting, storage and processing procedures. Good laboratory techniques should be used in all procedures to ensure proper performance of this assay. Because of the high analytical sensitivity of this assay and the need to avoid all possible sources of contamination,

NucliSens® HIV-1 QT
Summary of Safety and Effectiveness

extreme care should be taken to maintain the purity of the reagents and mixtures. All reagents should be monitored for purity. Discard any reagents that may be suspect or have exceeded their recommended shelf life. As with any diagnostic test, results from the NucliSens® HIV-1 QT assay should be interpreted with consideration of all clinical and laboratory findings.

A negative test result does not exclude the possibility of exposure to, or infection with, HIV-1.

The assay must be performed in strict accordance with the instructions in this package insert (especially those regarding contamination risks) and the NucliSens® Reader Operator Manual to obtain accurate, reproducible results.

Plasma should be used for determination of HIV-1 RNA levels. EDTA, citrate, or heparin may be used as anticoagulants.

In addition to quantitative HIV-1 RNA load, other virological or immunopathological factors may contribute to variable rates of CD4⁺ T-cell counts and clinical outcome of the disease.⁽⁸⁾

V. Alternative Practices and Procedures

There are various types of direct and indirect methods for detection and quantitation of human immunodeficiency virus in clinical samples. These methods allow the measurement of disease progression as a result of HIV infection and monitoring of patients' response to anti-retroviral therapy. Examples of these methods are listed below.

- a. HIV-1 Antigen assays including p24 core Antigen and Immune Complex Dissociation p24 Antigen tests
- b. Nucleic acid probe technologies for direct detection and quantitation of circulating viral particles
- c. Surrogate markers including CD4 and CD8 cell surface receptors
- d. Quantitative Cell Culture.

VI. Potential Adverse Effects of the Device on Health

Incorrect test results are possible in all *in vitro* diagnostic products. An erroneous low test result may lead to inappropriate treatment, or instill a false sense of security in a patient, which could lead to deterioration of the patient's condition. The possibility of incorrect results can happen with assignable causes such as a technician's error in following the procedures in the package insert or a device malfunction. An erroneous high test result on the other hand may contribute to unnecessary treatment or create anxiety or trauma to the patient. However, if the appropriate direction is followed as stated in the package insert, the likelihood of erroneous results are minimal from the use of this device.

The performance of the product in the clinical studies indicates that the benefit to the patient far outweighs any potential risk of adverse effect to the patient as a result of its use.

VII. Marketing History

The first generation NASBA HIV-1 RNA QT assay for quantitative determination of HIV-1 RNA "viral load" was introduced at the Biotechnology Against AIDS meeting in Florence, Italy in April of 1994. In 1997, the second generation NucliSens® HIV-1 QT assay was introduced. Organon Teknika markets this assay in more than 30 countries throughout North America, South America, Europe, Africa, Asia and Australia.

Since these initial placements in 1994, the NASBA/NucliSens® has been adopted by many users throughout the world. The device has also been introduced to investigators studying a wide

NucliSens® HIV-1 QT Summary of Safety and Effectiveness

variety of parameters related to AIDS. HIV-1 RNA “viral load” is now widely accepted as an important surrogate marker in research on the pathogenesis of HIV infection and AIDS. Research studies utilizing the NASBA assay include antiviral effects on progression of disease, mother to infant transmission, early detection of transmission, relationship of plasma viral load to HIV-1 RNA quantitation in other body compartments (e.g., CSF, seminal plasma, tissues), natural history cohorts, and others.

The NucliSens® HIV-1 QT assay has proven to be a valuable tool for studying the AIDS epidemic worldwide. Since introduction of the first NASBA assay and the subsequent NucliSens® HIV-1 QT, the product has not been withdrawn from marketing in any country for reasons related to safety or effectiveness or for any other reasons.

VIII. Summary of Studies

Potentially Interfering Substances

Elevated levels of lipids, bilirubin, and hemoglobin in specimens do not interfere with the quantitation of HIV-1 RNA by the NucliSens® HIV-1 QT assay. The presence of antinuclear antibodies in specimens had no detectable detrimental effect on the performance of the assay. Specimens from individuals known to be positive for rheumatoid factor, multiparous women and pregnant women showed no detectable detrimental effect on the quantitation of HIV-1 RNA. In addition, the presence of platelets in plasma did not appear to interfere with the performance of the NucliSens® HIV-1 QT assay.⁽³⁰⁾

The following compounds were found not to interfere with the quantitation of viral load by this assay: AZT, ddI, d4T, ddC, 3TC, indinavir, ritonavir, saquinavir, ganciclovir, acyclovir, zithromax-azithromycin, biaxin-clarithromycin, clofazamine, ethionamide, pentamidine, bactrin-trimethoprim sulfamethoxazole, dapsone, and diflucan.

Three anticoagulants, (EDTA, citrate, and heparin) have been evaluated for use with the assay. None of them were found to exhibit any significant interference on the quantitation of HIV-1 RNA load by the NucliSens® HIV-1 QT assay.

Linearity and Limit of Detection/Quantitation


The linear range (linearity) and limits of detection/quantitation of NucliSens® HIV-1 QT was assessed with specimens derived from a well-characterized HIV-1 RNA stock.

Linearity

Testing of diluted specimens derived from the HIV-1 RNA stock indicated a direct proportional relationship between the dilutions tested and the number of HIV-1 RNA copies reported. The linear range of the assay was determined using the data combined from three kit lots and linear regression analysis. Based on this study, the quantitation of HIV-1 RNA using this assay was determined to be linear over a range of 51 to 5,390,000 HIV-1 RNA copies.

NucliSens® HIV-1 QT
Summary of Safety and Effectiveness

Figure 2: NucliSens® HIV-1 QT Linearity



Limits of Detection and Quantitation.

- The limit of detection (LOD) is defined as the lowest HIV-1 RNA input level where at least 95% of the tests produce a result indicative of reactivity of the input sample for HIV-1.
- The limit of quantitation (LOQ) is defined as the lowest HIV-1 RNA input level where at least 95% of the tests produce quantifiable (within linear range) results with reasonable accuracy and precision.

The limits of detection and quantitation for NucliSens® HIV-1 QT were determined by testing at three sites of specimens containing HIV-1 RNA concentrations ranging from 659 to 41 copies/ml. Five specimens representing a range of HIV-1 RNA concentrations were selected and tested 72 times each using three lots of reagents.

Logistic regression was used to determine the relationship between the proportion detected and the log nominal input. The limit of detection at a 95% rate (LOD) was calculated to be 176 HIV-1 RNA copies/ml. This HIV-1 RNA input level for the Limit of Detection was verified by the results of panel member testing (Table 1). Panel member 16 (nominal input = 164 HIV-1 RNA copies) gave a positivity rate of 95.8% while lower viral concentrations resulted in positivity rates of less than 95% (Table 1).

**NucliSens® HIV-1 QT
Summary of Safety and Effectiveness**

Table 1: Analysis of NucliSens® HIV-1 QT Detection Rates

Panel Member	HIV-1 RNA input	HIV-1 RNA log input	Total Tested	Number Tested Positive	Observed fraction positive	% CV
14	659	2.8189	72	72	1.0000	72
15	329	2.5172	72	72	1.0000	53
16	164	2.2148	72	69	0.9583	94
17	82	1.9138	72	51	0.7083	98
18	41	1.6128	72	38	0.5278	100

The 95% detection rate was determined by logistic regression analysis to be 176 copies/ml. The coefficients of variation (CV) for panel members 15 and 16 with input concentrations of 329 and 164 copies/ml, respectively, were 53% and 94% (Table 9-1). Based on these results the Limit of Quantitation was determined to be equal to the Limit of Detection, i.e. 176 copies/ml.

Terminal Dilution of Clinical Specimens

To evaluate the relative sensitivity of NucliSens® HIV-1 QT, dilutions of six clinical specimens from HIV-1 infected individuals were made with human plasma negative for HIV-1. In each case, the response in copy number was a function of the dilution factor, with a diminishing HIV-1 RNA copy number as dilution factor increased. The detection of HIV-1 RNA was similar for each of the three NucliSens® HIV-1 QT lots evaluated.

Comparison to FDA Approved Device (Roche AMPLICOR HIV-1 MONITOR ä)

Dilution panel members used for the linearity study, clinical specimens, and panels provided by CBER (FDA) and Boston Biomedica (BBI panel 709)) were tested with NucliSens® HIV-1 QT and the Roche AMPLICOR HIV-1 MONITOR™ Test, standard procedure. If available, samples with results reported as <400 copies/ml with the Standard MONITOR™ procedure were re-tested with the AMPLICOR UltraSensitive procedure. The results are shown in Table 2 below.

Table 2: Comparison of NucliSens® HIV-1 QT to Roche AMPLICOR HIV-1 MONITOR ä (Standard)

Group	Number	Number of Specimens Reported with Copy Number	
		NucliSens® HIV-1 QT	Roche AMPLICOR HIV-1 MONITOR ä (Standard)
VQA Dilutions	14	14	9
Clinical Specimens	76	76	63
BBI & CBER Panels	13	12	8
Totals (%)	103	102 (99%)	80 (78%)*

*Three of these not reported were above the Monitor™ upper limit of >750,000 copies/ml

Seventeen (17) specimens that were <400 in the Standard Roche AMPLICOR HIV-1 MONITOR™ Test were re-tested with the UltraSensitive Roche AMPLICOR HIV-1 MONITOR™ procedure. Of the 17 specimens, NucliSens® HIV-1 QT detected 16 (94%) and the UltraSensitive Roche AMPLICOR HIV-1 MONITOR™ procedure detected 15 (88%). One specimen was not detected by either test (CBER-B2). Results are presented in Table 3 below

**NucliSens® HIV-1 QT
Summary of Safety and Effectiveness**

Table 3: Comparison of NucliSens® HIV-1 QT to Roche AMPLICOR HIV-1 MONITOR Test (Standard and UltraSensitive)


Specimen No.	NucliSens® HIV-1 QT		Roche AMPLICOR HIV-1 MONITOR®		
	Cp/ml	Log10	Standard	UltraSensitive	
			Cp/ml	Cp/ml	Log10
1	34	1.5315	<400	195	2.2900
2	920	2.9638	<400	744	2.8716
3	27	1.4314	<400	108	2.0334
4	72	1.8573	<400	95	1.9777
5	110	2.0414	<400	119	2.0755
6	260	2.4150	<400	92	1.9638
7	610	2.7853	<400	718	2.8561
8	670	2.8261	<400	563	2.7505
9	110	2.0414	<400	153	2.1847
10	970	2.9868	<400	1478	3.1697
11	630	2.7993	<400	402	2.6042
12	550	2.7404	<400	423	2.6263
13	71	1.8513	<400	108	2.0334
14	<LDL		<400	<50	
15	52	1.7160	<400	69	1.8388
16	160	2.2041	<400	<50	
17	340	2.5315	<400	207	2.3160

When NucliSens® HIV-1 QT and Roche AMPLICOR HIV-1 MONITOR Standard results were combined, an estimate of the Pearson product moment correlation, on a log copy basis, was 0.96. The 95% confidence interval for the correlation is (0.945, 0.977). The results reported by both tests are highly correlated.

An X-Y plot was created to show the log base 10 HIV-1 RNA copy numbers for all data obtained with NucliSens® HIV-1 QT and Roche AMPLICOR HIV-1 MONITOR™ Test. The slope and interception for this plot were 0.982 and 0.058, respectively. This further indicates the agreement of results obtained with both tests.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 3: The Log Results for the NucliSens[®] HIV-1 QT and Roche AMPLICOR HIV-1 MONITOR[®] Test are Plotted Compared to the Line of Perfect Agreement



Reproducibility

Reproducibility was determined by testing a panel (5 specimens) with varying levels of HIV-1 RNA with three lots at three sites by multiple technicians over four days. The variance components associated with run and replicate within run were used to estimate inter-assay variability and intra-assay variability, respectively. Inter-assay variability includes variability due to site, lot, and day. Intra-assay variability is replicate testing variability. The following table shows total inter-assay and intra-assay precision by specimen on an HIV-1 RNA copy basis.

Table 4: Per Specimen Precision Estimates - Copy Units

SID	N	MEAN	TOTAL		INTER-ASSAY		INTRA-ASSAY	
			SD	CV	SD	CV	SD	CV
SIII-1	72	3473611	1018421	29.3	832075	24.0	595469	17.1
SIII-2	72	356806	73806	20.7	52050	14.6	52691	14.8
SIII-3	72	47181	12729	27	10801	22.9	6857	14.5
SIII-4	72	3529	1257	35.6	998	28.3	774	21.9
SIII-5	71	373	383	102.8	234	62.8	305	81.7

NucliSens® HIV-1 QT
Summary of Safety and Effectiveness

As shown in Table 4, the NucliSens® HIV-1 QT was found to have relatively low variability. However, higher variability was observed for specimens with lower HIV-1 RNA concentrations.

Analytical Specificity

Analytical specificity was evaluated by testing a panel of specimens from individuals exhibiting the following infections or diseases, but serologically negative for HIV-1 infection. The majority of these specimens were tested with or without spiked HIV-1 RNA to assess false positive or false negative reactivity. Of the specimens tested, there was no evidence of false reactivity as shown in Table 5. The specificity was 100% across all sixteen-specimen categories listed below with an overall 95% confidence interval of 97.9% - 100%. Overall recovery of spiked HIV-1 RNA ranged from 97% to 101%.

Table 5: Possible Interfering Specimens Tested

Specimen Category	Number of Specimens Tested	Unspiked Results ¹	Mean Spiked RNA Copies Recovered	Mean Recovery Relative to Expected ²
HIV-2 positive	10	Not detected	24,500	97%
HTLV antibody positive	10	Not detected	33,500	100%
Herpes Simplex antibody positive	10	Not detected	28,100	100%
Icteric	10	Not detected	25,600	98%
Hemolyzed	10	Not detected	28,200	98%
Antinuclear Antibody	10	Not detected	24,900	99%
Lupus	10	Not detected	26,100	99%
Rheumatoid factor Positive	10	Not detected	35,000	100%
Specimens from Multiparous Woman	10	Not detected	26,500	98%
Specimens from Pregnant Woman	10	Not detected	24,900	97%
Syphilis (RPR) antibody positive	10	Not detected	51,500	Not Determined
Mycobacterium tuberculosis	10	Not detected	30,600	101%
Bacterially Contaminated	10	Not detected	69,300	Not Determined
Antibiotic/ Anti-Viral therapy	10	Not detected	60,410	Not Determined
Platelet	10	Not detected	26,400	101%
Lipemic	10	Not detected	30,600	99%
Summary Means	160	Not detected	34,132	99%

1 All specimens with no spiked HIV-1 virus were tested nonreactive (less than lower limit of detection).

2 Mean percent recovery was calculated by dividing the test result in log copy number with expected log copy number of the spiked specimens.

To evaluate the performance of NucliSens® HIV-1 QT assay on non-clade B HIV-1 subtypes, a panel of specimens with non-clade B subtypes obtained from BBI was tested using the NucliSens® HIV-1 QT assay and Roche AMPLICOR HIV-1 MONITOR Test (Standard). The results are presented in Table 6. Except for clade C, the NucliSens® HIV-1 QT assay showed improved detection efficiency as compared to the Roche AMPLICOR HIV-1 MONITOR Test. However, the difference of detection efficiency may not be statistically significant since only one specimen was tested for each subtype.

**NucliSens® HIV-1 QT
Summary of Safety and Effectiveness**

Table 6: NucliSens® HIV-1 QT and Roche AMPLICOR HIV-1 MONITOR® Test Results with BBI Panel PRD201

	NucliSens® Results	Roche AMPLICOR HIV-1 MONITOR® Results
Specimen	Observed copy/ml (log)	Observed copy/ml (log)
PRD201-01 Clade A	2,400 (3.3802)	Not detected
PRD201-02 Clade B	140,000 (5.1461)	149,360 (5.1742)
PRD201-03 Clade C	9,400 (3.9731)	23,861 (4.3776)
PRD201-04 Clade D	210,000 (5.3222)	394,638 (5.5961)
PRD201-05 Clade E	55,000 (4.7403)	12,387 (4.0929)
PRD201-06 Clade F	130,000 (5.1139)	50,423 (4.7026)
PRD201-07 Clade G	110,000 (5.0413)	8,401 (3.9243)
PRD201-08 Clade H	77,000 (4.8864)	1,508 (3.1784)
PRD201-09 (Diluent)	<LDL	Not Detected

<LDL = Less than Detectable Limit.

Clinical specificity

The clinical specificity of NucliSens® HIV-1 QT was assessed by testing a total of 502 plasma specimens collected from a low risk whole blood donor population. These specimens were non-reactive for antibodies to HIV-1 and HIV-2 as determined with an FDA licensed assay. The specificity of NucliSens® HIV-1 QT in this population was 100% (95% confidence limits of 99.27% to 100%).

Clinical Performance – Patient Prognosis and Drug Monitoring

Study Descriptions

The use of the assay was evaluated for patient prognosis and drug monitoring in conjunction with two clinical studies: 1) ACTG (AIDS Clinical Trial Group) Study 0152⁽³¹⁾ and 2) the Delta Trial.⁽³²⁾ ACTG 0152 evaluated the efficacy, safety, and long-term tolerance of zidovudine (AZT, also abbreviated as ZDV) or didanosine (ddl) as monotherapy or AZT and ddl in combination in a pediatric population. Children who were three months to eighteen years of age and who had received no or less than six weeks of previous anti-retroviral therapy and who had laboratory evidence of HIV infection were eligible for enrollment. Multiple specimens from 295 patients were tested with the assay. Two hundred eighty-seven (287) patients had reportable subsequent measurements for the calculation of proportion of progression to endpoint (Table 7). However, all 295 patients had reportable baseline HIV-1 RNA measurements prior to randomization into intent-to-treat groups (Figures 4 and 5).

The Delta Trial was a multi-center study evaluating the efficacy and safety of combinations of the anti-retroviral drug AZT plus ddl or AZT plus zalcitabine (ddC) compared to AZT alone in HIV infected adults. The Delta Trial began in March, 1992 and continued through 1995. Participants in the Delta Trial included one hundred seventy-five (175) centers from Europe, Australia and New Zealand. All testing was performed at two sites on repository specimens from this study. In order to be eligible for the extended virology study, samples must have been stored frozen at – 70°C or colder and must have been taken from individuals with a baseline sample and at least one later sample. Multiple samples from 1280 HIV infected participants in the Delta Trial were evaluated. Twelve hundred fifty-six (1256) had reportable subsequent measurement for the calculation of proportion of progression to endpoint (Table 8), of which eleven hundred thirty-three

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

(1133) also had reportable baseline HIV-1 RNA measurement prior to randomization into intent-to-treat groups (Figures 6 and 7).

Study Objectives

In each virology sub-study tied to these two investigations, the clinical utility was evaluated by 1) examining HIV-1 RNA measurement as a prognostic indicator of disease progression and 2) assessing the association, if any, of decreased HIV-1 RNA level with a change in patient therapy.

Patient Prognosis

To assess the association of RNA as measured by the assay with HIV-1 disease progression, a logistic function was used to model the probability of progressing to endpoint as a function of logarithm baseline HIV-1 RNA level. This model fits well, as measured by the Pearson Chi-Square Goodness-of-Fit Test for both the ACTG 152 study ($p=0.4304$) and Delta study ($p=0.7863$). For both studies, an overall association across the entire range of HIV-1 RNA concentrations was demonstrated with the test of slope coefficient ($p<0.001$). The estimate of slope was 0.976 for the ACTG 152 study and 0.999 for the Delta study. The estimates of slope for both studies are very similar and close to one. Using the logistic model, the odds ratio was defined and used to calculate how much difference the assay can distinguish in terms of risk of clinical endpoints. For an individual with a viral load of 10,000 copies/ml, the risk of disease progression is 1.35 times as high as for an individual with a viral load of 5,000 copies/ml. These data suggest that a higher viral load is associated with an increased risk in both studies.

For each study, the distributions of actual baseline HIV-1 RNA levels for patients progressing and not progressing to endpoint were plotted. Figures 4 and 5 show the distribution of HIV-1 RNA levels for patients progressing and not progressing to endpoint, respectively, for the ACTG 152 study. Similarly, Figures 6 and 7 show the distribution of HIV-1 RNA levels for patients progressing and not progressing to endpoint, respectively, for the Delta Study.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 4: Distribution of Baseline HIV-1 RNA Levels for Patients (N=90) Progressing to Endpoint – ACTG 152 Study





Figure 5: Distribution of Baseline HIV-1 RNA Levels for Patients (N=205) Not Progressing to Endpoint – ACTG 152 Study



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 6: Distribution of Baseline HIV-1 RNA Levels for Patients (N=144) Progressing to Endpoint – Delta Study





Figure 7: Distribution of Baseline HIV-1 RNA Levels for Patients (N=989) Not Progressing to Endpoint – Delta Study



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

The relationship between baseline RNA levels grouped by decile and the percent of subjects in each group progressing to endpoint in each decile was evaluated. Analyses were performed to look at differences in the probability of disease progression by sequentially testing combinations of data from different groups of deciles. With this approach, the proportion of subjects responding in decile 1 was contrasted with the proportion responding in the group formed by combining deciles 2-10. Next, the group formed by combining deciles 1 and 2 was contrasted against combined deciles 3-10, and so on. Tables 7 and 8 show decile, HIV-1 RNA copy number representation, and the numbers of individuals from the ACTG 152 and Delta Studies, respectively, progressing to endpoint and total number in each decile across all three treatment groups. The tables also show the contrast statement results of the test for differences between the proportions in sequential groupings of deciles. P-values associated with the contrasts for each study are not adjusted for multiple testing. These results show that the biggest increase in significance between successive decile groupings occurs between deciles two and three in both studies. Thus, in the ACTG 152 study, this jump occurs around 93,000 copies/mL; in the Delta around 28,000 copies/mL.

Table 7: Proportion to Endpoint by RNA Decile – ACTG 152 Study

Decile	Midpoint	Range	Progressing	N	Proportion	Cumulative Proportion	Sequential Testing	p value
1	<37000	<37000	3	29	0.10	0.10	1 vs 2-10	0.0183
2	65500	37000-93000	1	29	0.03	0.07	1-2 vs 3-10	0.0001
3	151500	93000-210000	7	31	0.23	0.12	1-3 vs 4-10	0.0001
4	275000	210000-340000	5	26	0.19	0.14	1-4 vs 5-10	0.0001
5	475000	340000-610000	8	29	0.28	0.17	1-5 vs 6-10	0.0001
6	855000	610000-1100000	11	33	0.33	0.20	1-6 vs 7-10	0.0001
7	1300000	1100000-1500000	12	25	0.48	0.23	1-7 vs 8-10	0.0001
8	1900000	1500000-2300000	11	28	0.39	0.25	1-8 vs 9-10	0.0001
9	3250000	2300000-4200000	14	29	0.48	0.28	1-9 vs 10	0.0001
10	>4200000	>4200000	18	28	0.64	0.31		
	Total		90	287	0.31			

Table 8: Proportion to Endpoint by RNA Decile – Delta Study

Decile	Midpoint	Range	Progressing	N	Proportion	Cumulative Proportion	Sequential Testing	p value
1	<12600	<12600	6	126	0.05	0.05	1 vs 2-10	0.0103
2	20300	12600-28000	6	126	0.05	0.05	1-2 vs 3-10	0.0001
3	38000	28000-48000	10	128	0.08	0.06	1-3 vs 4-10	0.0001
4	64000	48000-80000	10	134	0.07	0.06	1-4 vs 5-10	0.0001
5	96000	80000-112000	7	122	0.06	0.06	1-5 vs 6-10	0.0001
6	132000	112000-152000	15	118	0.13	0.07	1-6 vs 7-10	0.0001
7	186000	152000-220000	20	126	0.16	0.08	1-7 vs 8-10	0.0001
8	290000	220000-360000	22	125	0.18	0.10	1-8 vs 9-10	0.0001
9	480000	360000-600000	30	127	0.24	0.11	1-9 vs 10	0.0001
10	>600000	>600000	30	124	0.24	0.12		
	Total		156	1256	0.12			

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Cox proportional hazards regression models, stratified by anti-retroviral treatment group, were used to estimate the capability of reductions in RNA to predict lowering the risk of disease progression at various reference weeks throughout the studies. Separate analyses were performed to include baseline HIV-1 RNA and later weeks' HIV-1 RNA and to include baseline CD4⁺ count and later weeks' CD4⁺ count as predictors. Only patients surviving without disease progression up to the reference week were included in each analysis. If no progression occurred within the time on study, the observation was treated as right censored. HIV-1 RNA copies and CD4⁺ counts were log transformed. The hazard ratio in each model was defined as a 10-fold decrease in HIV-1 RNA or a 10-fold increase in CD4⁺ count.

For each target week, a proportional hazards (PH) model was fit to the data. For the ACTG 152 study, a target week of 24 weeks was analyzed. Table 9 shows the results for the model including baseline HIV-1 RNA and change from baseline HIV-1 RNA to the target week HIV-1 RNA as predictors. A change in HIV-1 RNA from baseline to 24 weeks was significant (p=0.0042) in predicting a reduced risk of disease progression after the baseline HIV-1 RNA adjustment. For every 10-fold decrease in HIV-1 RNA, the hazard of disease progression decreases by 61%.

Table 9: ACTG 152 Study
Cox Proportional Hazards Model Results – Pediatric Population
Risk of Disease Progression vs log₁₀(baseline RNA) and change from baseline RNA to 24 Weeks
Response: Time from Week to Endpoint for Patients Surviving Without Progression Through
Given Week

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
24	-Log ₁₀ (baseline RNA)	-0.9463	0.1845	26.3121	<0.0001	0.388	0.270	0.557
24	Change from baseline RNA to 24 wks	-0.6059	0.2118	8.1875	0.0042	0.546	0.360	0.826

Table 10 shows the results for the model including baseline CD4⁺ count and change from baseline CD4⁺ count to the target week CD4⁺ count as predictors. A change in CD4⁺ count from baseline to 24 weeks was significant (p=0.0001) in predicting a reduced risk of disease progression after the baseline CD4⁺ count adjustment. For every 10-fold increase in CD4⁺ count, the hazard of disease progression decreases by 84%.

Table 10: ACTG 152 Study
Cox Proportional Hazards Model Results – Stratified by Treatment Therapy – Pediatric Population
Risk of Disease Progression vs log₁₀(baseline CD4⁺) and change from baseline CD4⁺ to 24 Weeks
Response: Time from Week to Endpoint for Patients Surviving Without Progression Through
Given Week

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
24	Log ₁₀ (baseline CD4)	-1.8300	0.2925	39.1409	<0.0001	0.160	0.090	0.285
24	Change from baseline CD4 to 24 wks	-1.5180	0.3930	14.9230	0.0001	0.219	0.101	0.473

For the Delta study, target weeks of 8 and 16 weeks were analyzed. Table 11 shows the results for the model including baseline HIV-1 RNA and change from baseline HIV-1 RNA to the target

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

week HIV-1 RNA as predictors. A change in HIV-1 RNA from baseline to 8 or 16 weeks was significant ($p < 0.0001$) in predicting a reduced risk of disease progression after the baseline HIV-1 RNA adjustment. For every 10-fold decrease in baseline HIV-1 RNA, the hazard of disease progression decreases by 75%.

Table 11: Delta Study
Cox Proportional Hazards Model Results – Stratified by Treatment Therapy – Adult Population
Risk of Disease Progression vs log₁₀ (baseline RNA) and change from baseline RNA to Given
Weeks
Response: Time from Week to Endpoint for Patients Surviving Without Progression Through
Given Week

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
8	-Log ₁₀ (baseline RNA)	-1.4004	0.2019	48.1301	<0.0001	0.247	0.166	0.366
8	Change from baseline RNA to 8 wks	-1.5727	0.1894	68.9804	<0.0001	0.207	0.143	0.301

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
16	-Log ₁₀ (baseline RNA)	-1.3996	0.2462	32.3143	<0.0001	0.247	0.152	0.400
16	Change from baseline RNA to 16 wks	-1.2657	0.2116	35.7955	<0.0001	0.282	0.186	0.427

Table 12 shows the results for the model including baseline CD4⁺ count and change from baseline CD4⁺ count to the target week CD4⁺ count as predictors. A change in CD4⁺ count from baseline to 8 or 16 weeks was significant ($p=0.0304$ and $p=0.0185$, respectively) in predicting a reduced risk of disease progression after the baseline CD4⁺ count adjustment. For every 10-fold increase in CD4⁺ count, the hazard of disease progression decreases by 96-97%.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Table 12: Delta Study

**Cox Proportional Hazards Model Results – Stratified by Treatment Therapy – Adult Population
Risk of Disease Progression vs log10 (baseline CD4⁺) and change from baseline CD4⁺ to Given
Weeks**

**Response: Time from Week to Endpoint for Patients Surviving Without Progression Through
Given Week**

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
8	Log10(baseline CD4)	-3.2414	0.3054	112.6384	<0.0001	0.039	0.021	0.071
8	Change from baseline CD4 to 8 wks	-1.0266	0.4741	4.6882	0.0304	0.358	0.141	0.907

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
16	Log10(baseline CD4)	-3.7019	0.3852	92.3740	<0.0001	0.025	0.012	0.052
16	Change from baseline CD4 to 16 wks	-1.1009	0.4673	5.5507	0.0185	0.333	0.133	0.831

For each target week, another proportional hazards (PH) model was fit to the data. This model included baseline HIV-1 RNA, change from baseline HIV-1 RNA to the target week HIV-1 RNA, baseline CD4⁺ count, and change from baseline CD4⁺ count to the target week CD4⁺ count as predictors. For the ACTG 152 study, a target week of 24 weeks was analyzed. Table 13 shows the results for the model including baseline HIV-1 RNA, change from baseline HIV-1 RNA to the target week HIV-1 RNA, baseline CD4⁺ count, and change from baseline CD4⁺ count to the target week CD4⁺ count as predictors. A change in HIV-1 RNA from baseline to 24 weeks was marginally significant (p=0.0703) in predicting a reduced risk of disease progression after the baseline HIV-1 RNA adjustment and adjustment for CD4⁺ count and 24 week change in CD4⁺ count. For every 10-fold decrease in baseline HIV-1 RNA, the hazard of disease progression decreases by 52%.

Table 13: ACTG 152 Study

**Cox Proportional Hazards Model Results – Stratified by Treatment Therapy – Pediatric Population
Risk of Disease Progression vs log10 baseline RNA), change from baseline RNA to 24 Weeks,
log10(baseline CD4⁺), and change from baseline CD4⁺ to 24 Weeks**

**Response: Time from Week to Endpoint for Patients Surviving Without Progression Through 24
Weeks**

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
24	-Log10(baseline RNA)	-0.7332	0.1934	14.3665	0.0002	0.480	0.329	0.702
24	Change from baseline RNA to 24 wks	-0.4092	0.2261	3.2760	0.0703	0.664	0.426	1.034
24	Log10(baseline CD4)	-1.4861	0.2991	24.6853	<0.0001	0.226	0.126	0.407
24	Change from baseline CD4 to 24 wks	-1.3443	0.3954	11.5596	0.0007	0.261	0.120	0.566

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

For the Delta study, target weeks of 8 and 16 weeks were analyzed. Tables 14 and 15 show the results for the model including baseline HIV-1 RNA, change from baseline HIV-1 RNA to the target week HIV-1 RNA, baseline CD4⁺ count, and change from baseline CD4⁺ count to the target week CD4⁺ count as predictors. A change in HIV-1 RNA from baseline to 8 or 16 weeks was significant (p=0.0001) in predicting a reduced risk of disease progression after the baseline HIV-1 RNA adjustment and adjustment for CD4⁺ count and 8 or 16 week change in CD4⁺ count. For every 10-fold decrease in baseline HIV-1 RNA, the hazard of disease progression decreases by 64-69%.

Table 14: Delta Study
Cox Proportional Hazards Model Results – Stratified by Treatment Therapy – Adult Population
Risk of Disease Progression vs log10 (baseline RNA), change from baseline RNA to 8 Weeks,
log10(baseline CD4⁺), and change from baseline CD4⁺ to 8 Weeks
Response: Time from Week to Endpoint for Patients Surviving Without Progression Through 8 Weeks

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
8	-Log10 (baseline RNA)	-1.1808	0.2134	30.6045	0.0001	0.307	0.202	0.467
8	Change from baseline RNA to 8 wks	-1.4035	0.2038	47.4392	0.0001	0.246	0.165	0.366
8	Log10 (baseline CD4)	-2.8880	0.3509	67.7321	0.0001	0.056	0.028	0.111
8	Change from baseline CD4 to 8 wks	-0.0940	0.5236	0.0322	0.8575	0.910	0.326	2.540

Table 15: Delta Study
Cox Proportional Hazards Model Results – Adult Population
Risk of Disease Progression vs log10 (baseline RNA), change from baseline RNA to 16 Weeks,
log10(baseline CD4⁺), and change from baseline CD4⁺ to 16 Weeks
Response: Time from Week to Endpoint for Patients Surviving Without Progression Through 16 Weeks

Week	Model Covariate	Parameter Estimate	Standard Error	Wald Chi-Square	P-Value	Risk Ratio	Lower 95% Conf. Lim	Upper 95% Conf. Lim
16	-Log10 (baseline RNA)	-1.0224	0.2589	15.5925	0.0001	0.360	0.217	0.598
16	Change from baseline RNA to 16 wks	-0.8819	0.2242	15.4682	0.0001	0.414	0.267	0.642
16	Log10 (baseline CD4)	-3.1842	0.4185	57.9052	0.0001	0.041	0.018	0.094
16	Change from baseline CD4 to 16 wks	-0.7114	0.5046	1.9872	0.1586	0.491	0.183	1.320

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness


Measuring Response to Anti-retroviral Therapy

The association between HIV-1 RNA levels and survival time during the study without disease progression was assessed using Kaplan-Meier product limit survival curves. For this analysis the population was nominally divided into lower, middle, and upper thirds based on RNA levels. In the pediatric population, percent survival without progression referenced from week 24 was plotted as function of time in study based on the most recent RNA level prior to week 24. Only subjects who had not experienced disease progression prior to week 24 were included in the analysis. A similar analysis was performed with the adult population with a reference point of 8 weeks. A separate set of curves was generated for each treatment group.

Figures 8, 9, and 10 show the Kaplan-Meier curves for survival without progression to endpoint for the 3 anti-retroviral treatments studied with the pediatric population. Figures 11, 12, and 13 show similar curves for the 3 anti-retroviral treatments studied with the adult population. These figures demonstrate that, in each treatment group, subjects in the lowest RNA third survived longest without progression while subjects in the highest third progressed to endpoint most rapidly.

In another analysis, percent median change from baseline RNA by drug treatment was plotted versus weeks on study. Figure 14 plots percent change for the pediatric population study; Figure 15 shows a similar plot for the adult population study. In each study, the greatest and most enduring changes are seen in the treatment arms reflecting combination therapies.

Figure 8: ACTG 152
Kaplan-Meier Estimated Survival by Third of RNA Level at Week 24



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 9: ACTG 152 Kaplan-Meier Estimated Survival by Third of RNA Level at Week 24

Patients Surviving without Disease Progression Through Week 24
Pediatric Population Study - Treatment Therapy: ddl





Figure 10: ACTG 152 Kaplan-Meier Estimated Survival by Third of RNA Level at Week 24

Patients Surviving without Disease Progression Through Week 24
Pediatric Population Study - Treatment Therapy: ZDV+ ddl



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 11: Delta Trial Kaplan-Meier Estimated Survival by Third of RNA Level at Week 8

Patients Surviving without Disease Progression Through Week 8
 Adult Population Study - Treatment Therapy: AZT





Figure 12: Delta Trial Kaplan-Meier Estimated Survival by Third of RNA Level at Week 8

Patients Surviving without Disease Progression Through Week 8
 Adult Population Study - Treatment Therapy: AZT+ddI



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 13: Delta Trial Kaplan-Meier Estimated Survival by Third of RNA Level at Week 8

Patients Surviving without Disease Progression Through Week 8
 Adult Population Study - Treatment Therapy: AZT+ddC






Figure 14: ACTG 152
Percent Change in Median RNA Level from Baseline RNA Level as A Function of Weeks on Study
by Treatment Therapy



NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Figure 15: Delta Trial
Percent Change in Median RNA Level from Baseline RNA Level as A Function of Weeks on Study
by Treatment Therapy



Validation of Blood Collection and Processing

- Plasma should be used as the sample input. The standard sample volume is 1.0 ml (in the 9.0 ml Lysis Buffer tube).
- Blood should be collected in sterile tubes using EDTA, citrate or heparin as an anticoagulant. Clinical results indicate that there is no significant difference in the quantitation of HIV-1 RNA with the NucliSens[®] HIV-1 QT system for these three anticoagulants.
- Whole blood specimens collected using EDTA as an anticoagulant can be stored at room temperature (15-30°C) for up to 24 hours before processing without detectable loss of HIV-1 RNA. EDTA plasma can be stored at 2-8° C for up to 14 days.
- Plasma collected using EDTA, citrate or heparin is stable for ≥ 1 year at -70°C and can be frozen and thawed up to three times with no significant loss of HIV-1 RNA reported by the NucliSens[®] HIV-1 RNA QT system. Specimens repeatedly frozen and thawed or those containing particulate matter may give erroneous results.
- In Lysis Buffer, EDTA plasma specimens can be stored:
 - up to one year at -70°C **or**
 - for a maximum of 14 days at 2-8°C **or**
 - for a maximum of 24 hours at room temperature (15-30°C).

Note: Do not store specimens in Lysis Buffer at -20°C .

- Purified RNA eluate (post Isolation) can be stored for 12 months at -20 or -70°C or for 14 days at 2-8°C.
- Amplified material may be stored for 12 months at -20°C .

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

Specimen transport

If the shipment can be accomplished within 24 hours of collection, the specimen may be shipped at room temperature (15-30°C). Otherwise, the specimen should be shipped on dry ice. Specimens may be transported to the laboratory via courier, air freight, or regular mail in accordance with applicable Federal, state and local regulations that apply to the transportation of diagnostic specimens.

Note: Handle tubes stored at -70°C (or on dry ice) with care to avoid breakage associated with low temperature storage.

Note: Specimens in Lysis Buffer should be shipped at a temperature at or below -30°C (i.e, on dry ice).

IX Conclusions Drawn from the Studies

The studies demonstrate that the NucliSens[®] HIV-1 QT assay is safe and effective for use in measuring HIV-1 RNA within a range between 176 to 3.47 X 10⁶ copies/mL.

The data from the clinical studies contained in this application support the prognostic value of HIV-1 RNA measurements and demonstrate its clinical utility as a tool for monitoring the effects of antiretroviral therapy. For example, for an individual with a viral load of 10,000 copies/ml, the risk of disease progression is 1.35 times as high as for an individual with a viral load of 5,000 copies/ml. These data suggest that a higher HIV-1 RNA viral load is associated with an increased risk of disease progression in both the ACTG 152 and Delta studies. In both the ACTG and Delta studies, the greatest and most enduring changes for better clinical outcome were seen in the treatment arms with patients under combination therapies, which suggest that anti-retroviral therapies may reduce the measurable HIV-1 RNA and the progression of disease.

X. Benefit Analysis

The NucliSens[®] HIV-1 QT is an *in vitro* nucleic acid amplification test for the quantitation of HIV-1 RNA in human plasma. The test is intended for use in conjunction with clinical presentation and other laboratory markers as an indicator of disease prognosis by measuring baseline HIV-1 RNA levels, and for monitoring the effects of anti-retroviral therapy by serial measurements of plasma HIV-1 RNA. The NucliSens[®] HIV-1 QT is not intended to be used as a screening test for the HIV-1 virus nor is it to be used as a test to confirm the presence of an HIV-1 infection.

Quantitative measurements of HIV-1 viremia in the peripheral blood showed that higher levels of viral load correlated with a higher risk of clinical progression toward AIDS or death. Therefore, determination of viral load is a valuable marker for the prediction of disease progression. In addition, monitoring the efficacy of anti-viral therapy through serial measurement of viral loads may improve the clinical outcome of the infected patients under therapy^(19,20,21). The performance of the product in the clinical studies indicates that, when used properly, this assay could be used for the measurement of baseline HIV-1 RNA levels for disease prognosis and for the monitoring of effectiveness of anti-viral drug therapy, both of which would significantly benefit the patient while inflicting little, if any, adverse effect to the patient as a result of its use.

XI. FDA Decision

The FDA concluded that the results of the clinical studies demonstrated that the NucliSens[®] HIV-1 QT is a safe and effective test for the intended use. An approval letter was thus issued to bioMérieux Inc. for the intended use of this device.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

XII. References

- 1 Kievits T, van Gemen B, Lens P, et al: NASBATM isothermal enzymatic in vitro nucleic acid amplification optimized for the diagnosis of HIV-1 infection, *J Vir Meth* 1991; 35:273-286.
- 2 van Gemen B, Kievits T, Lens P, et al: Qualitative and quantitative detection of HIV-1 RNA by nucleic acid sequence-based amplification, *AIDS* 1993, 7 (suppl 2):S107-S110.
- 3 Boom R, Sol CJA, van der Noordaa J, et al: Rapid and simple method for purification of nucleic acids, *J Clin Microbiol* 1990; 28:495-503.
- 4 Bruisten S, van Gemen B, Huisman H, et al: Detection of HIV-1 Distribution in Different Blood Fractions of HIV-1 Seropositive Persons by Two Nucleic Acid Amplification Assays, *AIDS Research and Human Retroviruses* 1993; (9):259-265.
- 5 Bagnarelli P, Menzo S, Clementi M, et al: Detection of Human Immunodeficiency Virus Type 1 Genomic RNA in Plasma Samples by Reverse-Transcription Polymerase Chain Reaction, *J Med Vir* 1991;34:89-95.
- 6 Ottmann M, Innocenti P, Seigneurin JM, et al: The polymerase chain reaction for the detection of HIV-1 genomic RNA in plasma from infected individuals, *J Vir Meth* 1991; 31:273-284.
- 7 van Gemen B, Kievits T, Lens P, et al: Quantification of HIV-1 RNA in plasma using NASBATM during HIV-1 primary infection, *J Vir Meth* 1993; 43:177-188.
- 8 Piatak M, Saag MS, Lifson JD, et al: High Levels of HIV-1 in Plasma During All Stages of Infection Determined by Competitive PCR, *Science* 1993; 259:1749-1754.
- 9 Kessler HA, Blaauw B, Spear J, et al: Diagnosis of human immunodeficiency virus infection in seronegative homosexuals presenting with an acute viral syndrome, *JAMA* 1987; 258:1196-1199.
- 10 Lu W, Andrieu JM: Early Identification of Human Immunodeficiency Virus-Infected Asymptomatic Subjects Susceptible to Zidovudine by Quantitative Viral Coculture and Reverse Transcription-Linked Polymerase Chain Reaction, *J Inf Dis* 1993; 167:1014-1020.
- 11 Semple M, Loveday C, Weller I, Tedder R: Direct Measurement of Viraemia in Patients Infected With HIV-1 and Its Relationship to Disease Progression and Zidovudine Therapy, *J Med Vir* 1991; 35:38-45.
- 12 van Gemen B, van Beuningen R, Kievits T, et al: A one-tube quantitative HIV-1 RNA NASBA nucleic acid amplification assay using electrochemiluminescent (ECL) labelled probes, *J Vir Meth*, 1994; 49:157-168.
- 13 Nedelman J, Heagerty P, Lawrence C: Quantitative PCR with internal controls, *CABIOS* 1992; 8:65-70.
- 14 Nedelman J, Heagerty P, Lawrence C: Quantitative PCR: procedures and precisions, *Bull Math Biol* 1992; 54:477-502.
- 15 Blackburn GF, Shah HP, Kenten JH, Leland J, Kamin RA, Link J, Peterman J, Powell MJ, Shah A, Talley DB et al: Electrochemiluminescence detection for development of immunoassays and DNA probe assays for clinical diagnostics, *Clin Chem* 1991; 37:1534-1539.
- 16 Kenten JH, Gudibande S, Link J, Willey JJ, Curfman B, Major EO, Massey RJ: Improved electrochemiluminescent label for DNA probe assays: rapid quantitative assays of HIV-1 polymerase chain reaction products, *Clin Chem* 1992; 38:873-879.
- 17 Ho DD, Neumann AU, Markowitz M, et al: Rapid turnover of plasma virions and CD4 lymphocytes in HIV-1 infection. *Nature* 1995; 373:123-126.
- 18 Perelson AS, Neumann AU, Ho DD, et al: HIV-1 dynamics in vivo: virion clearance rate, infected cell life-span and viral generation time. *Science* 1996; 271:1582-1586.
- 19 Mellors JW, Kingsley LA, Kokka RP et al: Quantitation of HIV-1 RNA in plasma predicts outcome after seroconversion. *Ann. Intern. Med.* 1995; 122:573-579.
- 20 Jurriaans S, van Gemen B, Goudsmit J, et al: The natural history of HIV-1 infection: viral load and virus phenotype independent determinants of clinical course? *Virology* 1994; 204:223-233.
- 21 Mellors JW, Rinaldo CR, Kingsley LA, et al: Prognosis in HIV-1 infection predicted by the quantity of virus in plasma. *Science* 1996; 272:1167-1170.
- 22 de Wolf F, Spijkerman I, Schellekens PT, et al: Reciprocal predictive value of the quantity of HIV-1 RNA versus the number and function of CD4⁺ T-cells in time since seroconversion, *Published yet?*
- 23 O'Brien WA, Hartigan PM, Martin D, Esinhart D; Changes in plasma HIV-1 RNA and CD4⁺ lymphocyte count relative to treatment and progression to AIDS. *N. Engl. J. Med.* 1996; 334:426-431.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

- 24 Carpenter CJ, Fischl MA, Volberding PA, et al: Antiretroviral therapy for HIV infection in 1996. JAMA 1996; 146-154.
- 25 Hammer SM for the ACTG 175 Virology Substudy Team. Virological markers and outcome in ACTG 175. In: Third conference on retroviruses and opportunistic infections, Washington DC. 1996, abstract S24.
- 26 Gazzard B for the International Coordinating Committee. Further results from European/Australian Delta Trial. In: Third conference on retroviruses and opportunistic infections, Washington DC, 1996, abstract LB5A.
- 27 Gulick R, Mellors J, Havlir D: Potent and sustained antiretroviral activity of indinavir in combination with zidovudine and lamivudine. In: Third conference on retroviruses and opportunistic infections. Washington DC, 1996, abstract LB7.
- 28 Bruisten SM, Oudshoorn P, Cuypers HTM, et al: Stability of HIV-1 RNA in whole blood, plasma and serum during specimen handling and storage prior to amplification by NASBA-QT. published yet?
- 29 Holodniy M, Mole L, Jackson JB, et al: Comparative stabilities of human immunodeficiency virus RNA in plasma from samples collected in VACUTAINER CPT, VACUTAINER PPT, and standard VACUTAINER tubes. J. Clin. microbiol. 1995;33:1562-1566.
- 30 Witt, DJ, Kemper, M: Techniques for the evaluation of nucleic acid amplification technology performance with specimens containing interfering substances: efficacy of Boom methodology for extraction of HIV-1 RNA, J Virological Methods 1999, 79:97-111.
- 31 Palumbo, PE, Raskino, C, Fiscus, S, Pawha, S, Fowler, MG, Spector, SA, Englund, JA, Baker, CJ: Predictive Value of Quantitative Plasma HIV RNA and CD4⁺ Lymphocyte Count in HIV-Infected Infants and Children, JAMA 1998, 279:756-761.
- 32 Delta Coordinating Committee and Delta Virology Committee: HIV-1 RNA response to antiretroviral treatment in 1280 participants in the Delta Trial: an extended virology study, AIDS 1999, 13:57-65.

NucliSens[®] HIV-1 QT
Summary of Safety and Effectiveness

PMA REVIEW COMMITTEE

Signature

Date

<First name> <Last Name>, Chairperson

<First name> <Last Name>, Member

<First name> <Last Name>, Member

<First name> <Last Name>, Member

<First name> <Last Name>, Member

<First name> <Last Name>, Member

<First name> <Last Name>, Member