

UNITED STATES OF AMERICA
 DEPARTMENT OF HEALTH AND HUMAN SERVICES
 FOOD AND DRUG ADMINISTRATION
 CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
 DENTAL PRODUCTS PANEL
 MEDICAL DEVICES ADVISORY COMMITTEE

+ + + + +

MEETING

+ + + + +

THURSDAY, NOVEMBER 9, 2006

The meeting came to order at 8:30 a.m. in the ballroom of the Gaithersburg Holiday Inn, Two Montgomery Village Avenue, Gaithersburg, Maryland, RICHARD G. BURTON, Chairman, presiding.

PRESENT:

RICHARD G. BURTON, D.D.S., Chairman

SALOMON AMAR, D.D.S., Ph.D., Voting Member

WILLIAM J. O'BRIEN, M.S., Ph.D., Voting
Member

YIMING LI, Ph.D., Non-Voting Member

MASON DIAMOND, D.D.S., Industry
Representative

KURT C. GUNTER, M.D., Industry Representative

MICHAEL FLEMING, D.D.S., Consumer
Representative

JANINE E. JANOSKY, Ph.D., Temporary Voting
Member

MARK R. PATTERS, D.D.S., Ph.D., Temporary
Voting Member

JOHN R. ZUNIGA, Ph.D., D.M.D., Temporary
Voting Member

MICHAEL J. RYAN, Executive Secretary

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

I-N-D-E-X

AGENDA ITEM	PAGE
OPEN SESSION -- Welcome and Introductory....	3
Remarks	
Dr. Richard G. Burton, Chairman	
Mr. Michael J. Ryan,	
Executive Secretary	
Open Public Hearing.....	12
 Presentation by the Sponsor -	
Fuse Bone Graft.....	17
(P050053)	
 Presentation by the FDA - InFuse Bone Graft	
(P050053)	124
Dr. Peter L. Hudson, Biologist, Division of	
General, Restorative, and Neurological	
Devices	
 Dr. Zhiwei Zhang, Statistician, Office....	137
Of Surveillance and Biometrics	
 Dr. Robert S. Betz, Periodontist,.....	155
Dental Devices Branch	
 Open Public Session.....	277
 Panel Deliberations.....	289
 Summation	
FDA.....	278
Sponsor.....	278
 Panel Recommendation and Vote.....	325

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22

P-R-O-C-E-E-D-I-N-G-S

(8:33 a.m.)

CALL TO ORDER

OPEN SESSION -- WELCOME AND

INTRODUCTORY REMARKS

CHAIRMAN BURTON: Good morning. I am Dr. Richard Burton from the University of Iowa. I would like to welcome all of you to this meeting of the Dental Products Panel and to the CDRH Medical Devices Advisory Committee. I am the Chairman of the Dental Products Panel at this time, and I would like to call this meeting to order.

We are gathered here today to discuss the premarket approval application for the InFuse bone graft sponsored by Medtronic Sofamor Danek. This device consists of recombinant bone morphogenic protein, rhBMP-2, combined with a bovine collagen sponge.

I would like to go around the table, starting over here on the left, and have each of the members introduce themselves.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 MEMBER GUNTER: My name is Kurt
2 Gunter. And I'm the non-voting industry
3 representative.

4 MEMBER FLEMING: My name is Mike
5 Fleming. I am the non-voting consumer
6 representative on the Dental Products Panel.

7 MEMBER DIAMOND: My name is Mason
8 Diamond. I am the industry representative to
9 the Dental Products Panel.

10 MEMBER AMAR: Good morning. My
11 name is Salomon Amar. I am professor of
12 periodontology at Boston University. I am a
13 voting member.

14 MEMBER O'BRIEN: Bill O'Brien,
15 professor of biologic and material sciences at
16 the University of Michigan School of
17 Dentistry. And I am a voting member of the
18 panel.

19 MEMBER LI: I am Yiming Li,
20 professor of restorative dentistry at Loma
21 Linda University. I also serve as the
22 Director for Center for Dental Research. I am

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 a member of the Dental Products Panel.

2 MEMBER ZUNIGA: Good morning. My
3 name is John Zuniga. I am a professor of oral
4 and maxillofacial surgery at the University of
5 Texas Southwestern Medical Center at Dallas.
6 And I am a voting member of the panel.

7 MEMBER JANOSKY: Janine Janosky,
8 an associate professor at the University of
9 Pittsburgh School of Medicine. And I am a
10 consultant.

11 MEMBER PATTERS: Mark Patters.
12 I'm the Associate Dean for Academic Affairs
13 and professor of periodontology at the
14 University of Tennessee.

15 DR. LIN: Good morning. My name
16 is Chu Lin. I am the Director of the Division
17 of Anesthesiology, General Hospital Infection
18 Control and Dental Devices in the Office of
19 Device Evaluation, CDRH, FDA.

20 CHAIRMAN BURTON: Thank you.

21 The Executive Secretary will make
22 some introductory remarks at this time.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 EXECUTIVE DIRECTOR RYAN: Thank
2 you, Chairman Burton.

3 My name is Michael Ryan. I am the
4 Executive Secretary of the panel. I'll now
5 read into the record the deputization of
6 temporary voting member statement and the
7 conflict of interest statement.

8 "Pursuant to the authority granted
9 under the Medical Devices Advisory Committee
10 charter dated October 27, 1990, as amended on
11 April 20, 1995, I appoint the following
12 consultants as voting members of the Dental
13 Products Panel for the joint meeting to be
14 held on November 9th: Janine E. Janosky, Mark
15 R. Patters, John R. Zuniga.

16 "For the record, these individuals
17 are special government employees and are
18 consultants to this panel under the Medical
19 Advisory Committee. They have undergone the
20 customary conflict of interest review. They
21 have reviewed the material to be considered
22 for the meeting." This memo was signed by

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Daniel G. Schultz, M.D., Director of Center
2 for Devices and Radiological Health, FDA.

3 The conflict of interest statement
4 is as follows, "Food and Drug Administration
5 is convening today's meeting of the Dental
6 Products Panel of the Medical Devices Advisory
7 Committee under the authority of the Federal
8 Advisory Committee Act of 1972.

9 "With the exception of the
10 industry representative, all members and
11 consultants of the panel are special
12 government employees or regular federal
13 employees from other agencies and are subject
14 to federal conflict of interest laws and
15 regulations.

16 "Following information on the
17 status of this panel's compliance with federal
18 ethics and conflict of interest laws covered
19 by but not limited to those found at 18 USC
20 section 208 are being provided to participants
21 in today's meeting and to the public.

22 "FDA has determined that members

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 and consultants of this panel are in
2 compliance with federal ethics and conflict of
3 interest laws. Under 18 USC section 208,
4 Congress has authorized FDA to grant waivers
5 to special government employees who have
6 financial conflicts when it is determined that
7 the agency's need for particular individual
8 services outweighs his or her potential
9 financial conflict of interest.

10 "Members and consultants of this
11 panel who are special government employees at
12 today's meeting have been screened for
13 potential financial conflicts of interest of
14 their own as well as those imputed to them,
15 including those of their employer, spouse, or
16 minor child related to the discussion of
17 today's meeting.

18 "These interests may include
19 investments, consulting, expert witness
20 testimony, contracts, grants, CRADAs,
21 teaching, speaking, writing, patents and
22 royalties, and primary employment.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 "Today's agenda involves the
2 review of a premarket approval application for
3 the InFuse bone graft. This device is a
4 combination product which features a collagen
5 sponge that incorporates a recombinant bone
6 morphogenetic protein.

7 "The device is indicated for the
8 following oral maxillofacial bone grafting
9 procedures as an alternative to autogenous
10 bone graft for oral maxillofacial bone
11 grafting procedures, sinus augmentation, and
12 ridge augmentation at extraction socket sites.

13 "Particular matters during the
14 meeting or specific matters related to PMA
15 will be discussed. Based on the agenda for
16 today's meeting and all financial interests
17 reported by the panel members and consultants,
18 no conflict of interest waivers have been
19 issued in connection with this meeting.

20 "A copy of the statement will be
21 available for review at the registration table
22 during this meeting and will be included as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 part of the official transcript.

2 "Dr. Mason Diamond is serving as
3 the dental device industry representative,
4 acting on behalf of all related industry, and
5 is employed by TyRx Pharma, Incorporated.

6 "Dr. Kurt Gunter is serving as the
7 biologics industry representative, acting on
8 behalf of all related industry, and is
9 employed by Hospira, Incorporated.

10 "We would like to remind members
11 and consultants that if the discussions
12 involve any other products or firms not
13 already on the agenda for which an FDA
14 participant has a personal or imputed
15 financial interest, participants need to
16 exclude themselves from such involvement. And
17 their exclusion will be noted for the record.

18 "FDA encourages all other
19 participants to advise the panel of any
20 financial relationships that they may have
21 with any firms at issue. Thank you."

22 If you have not done so already, I

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 would like to request that everyone in
2 attendance please take the opportunity to sign
3 the attendance sheet that's available at the
4 door. I would also like to request that
5 everyone turn off their cell phone ringers.

6 Transcripts of today's meeting
7 will be available from Neal Gross and Company,
8 Incorporated. Information on purchasing
9 videos of today's meeting can be found on the
10 table outside the meeting room.

11 Presenters to the panel who have
12 not already done so should provide FDA with a
13 hard copy of their remarks, including
14 overheads. Ms. Annemarie Williams will
15 collect these for me at the podium.

16 With that, I will turn the meeting
17 over to Chairman Burton. Chairman?

18 CHAIRMAN BURTON: Thank you.

19 I would like to note for the
20 record that the voting members present
21 constitute a quorum, as required by 21 CFR
22 part 14. We will now proceed with the agenda.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 OPEN PUBLIC HEARING

2 CHAIRMAN BURTON: This is the
3 first of two open public hearing sessions for
4 this meeting. The second open public hearing
5 session will follow the panel discussion this
6 afternoon.

7 At these times, public attendees
8 are given the opportunity to address the
9 panel, to present data or views relevant to
10 the panel's activities.

11 I would like to remind public
12 observers at this meeting that while this
13 portion of the meeting is open to the public
14 for observation, public attendees may not
15 participate except at the specific request of
16 the Chair. You will be given no more than ten
17 minutes for your presentation.

18 Both the Food and Drug
19 Administration and the public believe in a
20 transparent process for the
21 information-gathering and decision-making
22 process. To ensure such transparency at the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 open public session of the Advisory Committee
2 meeting, the FDA believes that it is important
3 to understand the context of each individual's
4 presentation.

5 For this reason, the FDA
6 encourages you, the open public hearing
7 speaker, at the beginning of your written or
8 oral statement to advise the Committee of any
9 financial relationship that you may have with
10 the sponsor; its product; and, if known, its
11 direct competitors. For example, this
12 financial information may include the
13 sponsor's payment of your travel, lodging, or
14 other expenses in connection with your
15 attendance at the meeting.

16 Likewise, FDA encourages you at
17 the beginning of our statement to advise the
18 Committee if you do not have any such
19 financial relationships. If you choose not to
20 address this issue of financial relationships
21 at the beginning of your statement, it will
22 not preclude you from speaking.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 I would ask at this time that
2 persons addressing the panel come forward to
3 the microphone and speak clearly as the
4 transcriptionist is dependent upon this as a
5 means of providing an accurate transcription
6 of the proceedings of the meeting. If you
7 have a hard copy of your presentation, please
8 provide it to the FDA staff for use by the
9 transcriptionist to help provide an accurate
10 record of the proceedings.

11 Okay. The first speaker is Vivian
12 Roblin.

13 MS. ROBLIN: My name is Vivian
14 Roblin, and I am speaking on behalf of
15 Medtronic. Ten years ago, at the age of 62,
16 I had no teeth, no upper teeth, and I have no
17 bone.

18 If I laughed, I sneezed, or I
19 coughed, the denture fell out. No amount of
20 sticky stuff would hold the denture in. I was
21 limited to soft food. And it was a very
22 depressing time.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 I was referred to Dr. Spagnoli in
2 October of '96. He put the InFuse in my
3 mouth. Six months later, they were able to
4 put in eight implants. I had had that much
5 bone growth.

6 I have minimal discomfort from the
7 surgery. Would I do it again? Yes, I would.

8 Fortunately, I don't have to. I hope this
9 product will be available worldwide for people
10 that have my problem because, really, with
11 people living longer, it gives you a quality
12 of life that I did not have ten years ago.

13 That's my story. If you have any
14 questions, I would be happy to answer them.

15 MEMBER O'BRIEN: Where was the
16 surgery performed: Dr. Spagnoli's office?

17 MS. ROBLIN: Yes.

18 MEMBER O'BRIEN: Where is he
19 located?

20 MS. ROBLIN: Charlotte, North
21 Carolina.

22 MEMBER O'BRIEN: Thank you.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 CHAIRMAN BURTON: Do you have
2 other than your overall results any comments
3 you would like to make about your clinical
4 course in terms of problems you had or didn't
5 have during the course of treatment?

6 MS. ROBLIN: Yes. I did not have
7 any problems. Everything went just as Dr.
8 Spagnoli thought it would. I never dreamed
9 that I would be able to eat anything I want,
10 but I can now. It's a fabulous product.

11 CHAIRMAN BURTON: Thank you very
12 much for your input.

13 MS. ROBLIN: Thank you.

14 CHAIRMAN BURTON: Thank you for
15 coming.

16 That was the only preregistered
17 speaker that we had at this time. Are there
18 any others who wish to speak during this time
19 frame?

20 (No verbal response.)

21 CHAIRMAN BURTON: Hearing none,
22 we'll move on to the presentation by the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 sponsor. Medtronic Sofamor Danek will now
2 give their presentations on this PMA. And we
3 have three listed speakers that are Dr. Chin,
4 Dr. Marx, and Dr. Cochran. I don't know if
5 you care to stay in that order. Is that
6 correct? Okay. Dr. Edward Chin?

7 PRESENTATION BY THE SPONSOR -

8 INFUSE BONE GRAFT (P050053)

9 DR. CHIN: Good morning, members
10 of the panel, the Dental Products Advisory
11 Panel. My name is Ed Chin. And I am the
12 Group Director of Regulatory Affairs of
13 Medtronic Spinal and Biologics in Memphis,
14 Tennessee.

15 We have the pleasure to present to
16 you the results of decades of research and
17 development of rhBMP-2 for use in oral and
18 maxillofacial procedures. The InFuse bone
19 graft product is the combination of work of
20 hundreds of scientists and clinicians who have
21 worked over the years. And I would like to
22 acknowledge their efforts to make this product

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 available to the surgeons and their patients.

2 Today you will hear from two
3 investigators who participated in IDE clinical
4 trials. Dr Robert Marx of the University of
5 Miami will present the clinical problem that
6 patients face and the clinical data of InFuse
7 in sinus augmentation surgery.

8 Dr. Marx will be followed by Dr.
9 David Cochran of the University of Texas
10 Health Sciences Center in San Antonio, Texas.

11 Dr. Cochran will present the clinical data of
12 InFuse in extraction socket augmentation
13 surgery as well as the overall safety data
14 developed in our clinical trials. I will then
15 return for closing remarks.

16 We have also assembled here today
17 many of the scientists who performed their
18 preclinical research and some of the
19 investigators who participated in the clinical
20 trials as well as scientists and experts,
21 members of the clinical and regulatory staff
22 of Medtronic, Wyeth, and Allquest, who are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 available to answer questions from the panel.

2 I would especially like to
3 acknowledge Dr. Philip Boyne from Loma Linda
4 University, who is one of the recognized
5 pioneers of rhBMP-2 research in oral and
6 maxillofacial surgery and who wrote the
7 seminal paper on sinus lift procedures.

8 Discovery of osteoinductivity of
9 BMP was first made by Dr. Marshall Urist in
10 1965. In his landmark research, Dr. Urist
11 found that certain proteins, which he later
12 termed "bone morphogenetic proteins," BMPs,
13 stimulated the formation of new bone when
14 placed into a non-bony site of a rat. Thus,
15 the term "osteoinductivity" was coined to
16 describe this phenomenon.

17 Only BMPs have been demonstrated
18 to be osteoinductive. In the 1980s,
19 researchers of Wyeth Bioforma developed a
20 method to synthesize the osteoinductive bone
21 morphogenetic proteins commonly referred to as
22 rhBMP-2 using recombinant methods. We are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 also fortunate to have here today scientists
2 who cloned the BMP-2 and performed this work,
3 Dr. John Wozney from Wyeth.

4 As shown on this slide, rhBMP-2
5 production cells are grown in a bioreactor
6 that contains a well-defined nutrient media
7 free of human or animal-derived components.
8 The protein of interest is separated from
9 process stream components by a streamered series
10 of three chromatography steps resulting in
11 rhBMP-2 of very high quality and purity.

12 For added assurance of viral
13 safety, each batch is processed through a
14 nanofilter. Throughout the production
15 process, quality control testing is performed
16 to assess the consistency of the sample, the
17 processing and safety, purity and activity of
18 the resulting rhBMP-2 protein.

19 rhBMP-2 that has met the
20 established quality criteria is sterile
21 filtered; freeze dried in vials; and then
22 further tested for consistency, safety, and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 activity, involves meeting specification of
2 subsequently assembled into InFuse kits.

3 InFuse is commercially available
4 in four kit configurations containing either
5 4.2 milligrams or 12 milligrams of rhBMP-2.
6 The vials contain a free stripe powder that
7 has been reconstituted at the time of surgery
8 with sterile water to a final concentration of
9 1.5 milligrams per ml. The solution is then
10 applied to a type I bovine absorbable collagen
11 sponge referred to as ACS in this
12 presentation.

13 The ACS localizes the activity of
14 rhBMP-2 and that provides the scaffolding for
15 bone formation. The absorbable collagen
16 sponge is a commercially available product
17 that is manufactured by Integra LifeSciences.

18 FDA approved this hemostatic sponge in a PMA
19 application in 1981. The product for which we
20 are seeking approval is the same product
21 currently on the market.

22 RhBMP-2 is a specific

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 concentration, combined with ACS, is the
2 commercial product called InFuse bone graft,
3 which will from this point forward be referred
4 to as InFuse, and this is the product we are
5 discussing today.

6 The safety and effectiveness of
7 InFuse has already been demonstrated in two
8 previous PMA approvals. The first PMA
9 approval was in 2002 for interior lumbar
10 spinal fusion. The second PMA approval was
11 granted in 2004 for open tibia fractures.

12 There were 437 patients who
13 received InFuse in IDE clinical trials for
14 these indications. In addition, over 1,200
15 patients received InFuse or rhBMP-2 on other
16 carriers in clinical trials that are in
17 various stages of completion. Thus, our
18 clinical experience under rigorously
19 controlled, FDA-approved clinical trials is
20 very extensive.

21 Over the years, research sought to
22 find a bone grafting agent that is truly

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 osteoinductive. That search ended with the
2 approval of InFuse, providing surgeons with
3 the long-sought osteoinductive product to help
4 their patients.

5 The therapeutic benefits have been
6 available for patients requiring spinal
7 fusions or tibia repair. Today we are here to
8 seek a third therapeutic benefit for our
9 patients who need oral and maxillofacial
10 treatments to replace teeth.

11 Similar to other PMAs, we are
12 seeking a third indication for InFuse bone
13 graft. This is an oral indication where
14 InFuse bone graft again induces bone formation
15 that leads to a patient therapeutic benefits,
16 in this case to replace teeth.

17 The models studied in this PMA are
18 sinus augmentation supported by three studies
19 and extraction socket with buccal wall defects
20 augmentation supported by two studies.

21 BMP is one of the most studied of
22 all bone-forming agents. This graph

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 illustrates the large body of knowledge that
2 exists for BMP. The red line graphs the
3 cumulative number of publications over the
4 last 30 years.

5 Over 5,000 articles have been
6 published, and research continues. The blue
7 line shows 31 regulated clinical studies of
8 rhBMP-2 products conducted over the last 13
9 years, in which over 1,700 patients have been
10 enrolled.

11 In the oral and maxillofacial
12 space, early preclinical safety studies
13 provided the foundation for rhBMP-2
14 development. Preclinical studies were
15 conducted in lower to higher animal species,
16 as shown here, enabling human clinical trials
17 to begin in 1994.

18 Human experiences from five
19 prospective clinical studies provide the
20 evidence to unquestionably support an approval
21 recommendation for InFuse in oral and
22 maxillofacial bone grafting procedures.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Medtronic is seeking approval for
2 the following indications. InFuse bone graft
3 is indicated as an alternative to autogenous
4 bone graft for sinus augmentations and
5 localized alveolar ridge augmentations for
6 defects associated with extraction sockets.

7 InFuse has already been proven to
8 be safe and effective for two orthopedic uses.

9 We will present evidence from multi-centered,
10 prospective controlled clinical trials that
11 provide valid scientific evidence to support
12 that InFuse is safe and effective to, one,
13 regenerate bone; two, that that bone supports
14 dental implant placement; and, three, that the
15 restoration is stable over time.

16 This research has been recognized
17 by professional societies as outstanding high
18 quality work. The American Academy of Oral
19 and Maxillofacial Surgeons Journal editorial
20 board awarded the 2005 Daniel M. Laskin Award
21 for the most outstanding article published in
22 the Journal of Oral and Maxillofacial

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Surgeries to Dr. Boyne and others for their
2 article, "De Novo Bone Induction by rhBMP-2 in
3 Maxillofacial Sinus for Augmentation."

4 The American Academy of
5 Periodontology Foundation bestowed the 2005
6 Tarrson research award in oral plastic
7 surgeries to Drs. Florellini and others for
8 their paper, "Randomized Study Evaluating
9 rhBMP-2 for Extraction Socket Augmentation."

10 Today we have the privilege of
11 having several of those authors present. Dr.
12 Robert Marx and Dr. David Cochran of those
13 award-winning research papers will present
14 these clinical results, which is the basis for
15 this PMA.

16 First, I would like to introduce
17 Dr. Robert Marx, who will present the clinical
18 need for InFuse in oral and maxillofacial
19 surgery and the first part of the clinical
20 data that is the foundation of this PMA.

21 Thank you.

22 DR. MARX: Good morning, panel

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 members. My name is Dr. Robert Marx. I am an
2 oral and maxillofacial surgeon and Chief of
3 the Department of Oral and Maxillofacial
4 Surgery at the University of Miami, Miller
5 School of Medicine.

6 I have no direct financial
7 interest in the product under review today. I
8 am a consultant for Medtronic, which is
9 covering my expenses for attending this
10 meeting. I participated in the IDE clinical
11 trials of this device as a clinical
12 investigator since its inception in 1994.

13 My colleague, Dr. David Cochran,
14 and I have been asked to present the data from
15 the clinical studies of InFuse as an
16 alternative to autogenous grafts for
17 maxillofacial conditions, specifically sinus
18 augmentation and localized alveolar ridge
19 augmentations for defects associated with
20 extraction sockets. Autogenous bone grafts
21 will be referred to in this presentation as
22 bone graft.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 I would first like to describe the
2 clinical needs that require a bone graft. The
3 essential need occurs when individuals lose
4 teeth and thereby lose bone. There is a need
5 to replace bone loss due to disease, such as
6 this, due to trauma, and due to congenital
7 absence of bone.

8 We will show you that the use of
9 InFuse will provide bone support to replace
10 missing teeth and in doing so restore
11 structure and function as well as the
12 appearance of the individual.

13 These photographs show an extreme
14 example of an individual who lost significant
15 amounts of bone, which has resulted now in the
16 loosening of her dentures. She is
17 representative of a totally dentureless
18 patient enrolled in the sinus augmentation
19 study.

20 She did not have enough bone to
21 comfortably wear dentures or to have dental
22 implants placed. As you can see, this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 correlates to a facial collapse that affects
2 her speech and her eating abilities. It also
3 affects her appearance.

4 To regenerate enough bone to
5 support facial contours and to allow either
6 the placement of dental implants or the
7 wearing of dentures, this patient would
8 require an extensive bone harvest and a
9 grafting procedure.

10 The current standard of care is
11 autogenous bone grafting. It has certain
12 advantages. It is the patient's own bone.
13 And, therefore, it does not have any risk
14 related to transmissible diseases. It has
15 proven effectiveness as well.

16 However, autogenous bone grafts
17 such as these have also significant
18 disadvantages, mainly donor site morbidity of
19 pain, blood loss, and permanent scars. It
20 also extends the surgical time and the
21 anesthesia time. And its availability in some
22 patients is very limited.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Because of the risks associated
2 with extending operating time, postoperative
3 pain, infection, and even sensory bone loss,
4 many clinicians would prefer not to harvest
5 bone from their patients. And patients prefer
6 not to undergo this additional painful
7 procedure.

8 There are a number of
9 disadvantages to the bone graft. Today we
10 will show you that InFuse overcomes these
11 disadvantages. There is certainly a clear
12 need to grow bone with a product that doesn't
13 have the risks and morbidity associated with
14 such a bone harvest.

15 Bone grafting is also not a benign
16 procedure. Patients often will continue to
17 suffer from pain or numbness or sensory nerve
18 loss at the donor site long after the oral
19 surgery has healed. Patients will often have
20 permanent numbness at this site. There is
21 also significant blond loss and a significant
22 risk of donor site infection.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Because of the risk associated
2 with extended operating time, postoperative
3 pain, infection, or sensory loss, many
4 clinicians would prefer not harvest bone from
5 their patients. And many patients avoid
6 needed procedures due to their fear of bone
7 graft harvest and the pain associated with it,
8 essentially denying themselves access to care.

9 Ullman, et al., reported the
10 complication rates associated with iliac crest
11 bone grafts in this Journal of Bone and Joint
12 Surgery publication. They found that patients
13 had an average of over 200 milliliters of
14 blood loss, a 3 percent instance of hematomas,
15 an 8 percent sensory loss, of which 5 percent
16 was a permanent sensory loss with associated
17 numbness. They also had two percent with
18 chronic pain.

19 The proposed indications for
20 InFuse as an alternative to the autogenous
21 bone graft, there are a number of
22 disadvantages to bone grafts, as you can see.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Today we will show you that InFuse overcomes
2 these disadvantages.

3 The oral and maxillofacial surgery
4 program objectives were to demonstrate
5 effectiveness of rhBMP-2/ACS in the following:
6 one, regenerate or grow normal physiologic
7 bone; two, to provide an adequate amount of
8 good quality bone to support dental implants
9 and dental restorations; three, to produce
10 bone that remains stable under long-term
11 functional loading, providing a durable
12 result; and, finally, four, to demonstrate a
13 safety profile in the maxillofacial
14 indications.

15 The clinical studies were designed
16 to collect evidence to prove these. Bone
17 density measurements, bone biopsy, histologic
18 studies were accomplished. And CT scans were
19 performed to accomplish all of these
20 objectives.

21 The evidence for this PMA is
22 derived from two clinical models. The first

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 that I will present is for the sinus
2 augmentation indication illustrated here,
3 where an opening is made into the lateral wall
4 of the maxillary sinus, the membrane is
5 elevated, and bone or InFuse is placed.

6 The second will be the extraction
7 socket defect augmentation, in which this
8 lateral or buccal wall is lost and represents
9 a true critical-sized defect.

10 The science augmentation studies
11 were prospective, controlled, clinical trials.

12 These data provide a high level of clinical
13 evidence. After completion of pre-clinical
14 studies, a pilot study was initially performed
15 to assess the feasibility of using rhBMP-2/ACS
16 in sinus augmentation procedures utilizing a
17 concentration of 0.43 milligrams per
18 milliliter. That's this one.

19 Although bone formed in this
20 study, it was not optimal for dental implant
21 placement. Therefore, a randomized dosing
22 study was then conducted which evaluated 0.75

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 milligrams and 1.5 milligrams per milliliter
2 concentrations, which were selected based on
3 data from the preclinical, pharmacokinetic,
4 local bioavailability, and pharmacologic
5 studies.

6 Data from the dosing study
7 demonstrated the ability of rhBMP-2/ACS to
8 successfully induce an adequate amount of bone
9 in this surgical procedure for dental implant
10 placement and found that 1.5 milligrams per
11 milliliter concentration to be the most
12 effective concentration used.

13 To confirm these observations, a
14 randomized pivotal study was conducted with
15 the 1.5 milligrams per milliliter
16 concentration. The data demonstrated that
17 InFuse could safely and effectively induce new
18 bone, which could receive dental implants that
19 could then be functionally loaded and maintain
20 their functional loading over a long period of
21 time.

22 The data to be presented today and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 which was submitted in the PMA is an analysis
2 of subjects receiving 1.5 milligrams per
3 milliliter rhBMP-2/ACS in the sinus dosing
4 study and the sinus pivotal study. These data
5 were pooled following statistical
6 justification that they were homogeneous
7 populations with respect to demographics,
8 baseline characteristics, and clinical
9 outcomes.

10 By agreement with the FDA and with
11 similar justification, the autogenous bone
12 graft subjects in the sinus dosing study and
13 the sinus pivotal study were also pooled. The
14 primary efficacy endpoint for these analyses
15 was the same as that approved for the pivotal
16 study, which was the rate of functional
17 loading of the implant-borne restoration at
18 six months.

19 The target success rate calculated
20 for this study was 73 percent. This
21 predetermined target was selected based upon
22 data from the dosing study and a review of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 literature at the time.

2 The investigators determined that
3 70 percent was the minimum clinically
4 acceptable and meaningful success rate for
5 rhBMP-2/ACS success. And this target value
6 was adopted in the study design. The target
7 was then set at 73 percent for this
8 statistical consideration and sample size
9 requirements. This success rate was submitted
10 in the IDE protocol and approved by the
11 reviewers within the FDA.

12 The primary objectives of the
13 sinus augmentation study were to evaluate the
14 effectiveness of InFuse to induce adequate
15 bone to successfully support implant-borne
16 restorations after six months of functional
17 loading and to evaluate the safety of
18 rhBMP-2/ACS compared to a bone graft.

19 The secondary objectives of the
20 sinus augmentation studies were to evaluate
21 the overall quality of the bone from CT scans,
22 histology, and functional loading results of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 InFuse compared to those of the bone graft.

2 This slide illustrates the
3 sequence of events in our sinus augmentation
4 studies. Testing performed in each study
5 period is shown across the bottom of this
6 line. Following the baseline period with the
7 initial CT scans and other studies, patients
8 underwent the sinus lift procedure in which
9 InFuse was implanted under the sinus membrane
10 to induce new bone.

11 At approximately four to six
12 months, each patient was evaluated by CT scans
13 to determine whether or not a dental implant
14 could be placed. After implant placement, a
15 periapical radiograph was also taken. Core
16 biopsies of the bone for also taken for
17 histologic assessment.

18 If the clinician determined that
19 there was adequate quality and quantity of
20 bone, the dental implant was placed.
21 Following osseointegration of the dental
22 implants, the implants were then uncovered. A

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 dental prosthesis was fabricated and then
2 placed. Then functional loading of the
3 prosthesis began.

4 A third CT scan was taken at
5 either six months post-dental implant
6 placement or at six months post-functional
7 loading depending upon the study. The
8 patients were then assessed at six-month
9 intervals through 24 months post-functional
10 loading, which was approximately 36 months
11 following the original placement of
12 rhBMP-2/ACS.

13 Now, the initial cohort of the
14 sinus dosing study patients were randomized to
15 0.75 milligrams per milliliter rhBMP-2/ACS or
16 the bone graft.

17 After acute safety was
18 established, patients were randomized between
19 1.5 milligrams per milliliter rhBMP-2/ACS or
20 the bone graft itself as the second cohort.
21 In the sinus pivotal study, patients were
22 randomized to receive either InFuse or bone

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 graft.

2 The clinical data were collected
3 in three separate prospective controlled IDE
4 clinical trials. Patients were treated with
5 various concentrations of rhBMP-2/ACS or a
6 bone graft. In these studies, bone graft is
7 defined as either autogenous bone alone or a
8 combination of autogenous bone and allogenic
9 bone, consistent with the current standard of
10 care.

11 Two hundred twenty patients were
12 enrolled in the sinus augmentation studies at
13 21 different study sites. The effectiveness
14 data subset consists of 82 patients from the
15 pivotal study and 17 patients from the dosing
16 study, for a total of 99 patients, all of whom
17 were randomized to be implanted with the 1.5
18 milligram per milliliter concentration of
19 rhBMP-2/ACS plus collagen response, which is
20 the product of InFuse.

21 I will now present the results.
22 This representative slide is a preoperative

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 panoramic CT scan view of the lower portion of
2 the sinus cavity. As you can see here, the
3 sinus is hyperpneumatized and there is
4 insufficient bone to place an implant.

5 The lower CT scan now, taken from
6 the patient at 16 weeks post-InFuse placement,
7 shows a large amount of de novo bone
8 formation, new bone formation formed beneath
9 the sinus membrane on each side of the
10 maxilla. Indeed, InFuse was shown to induce
11 new bone in this indication.

12 On this next slide, CT scans from
13 another sinus augmentation patient are shown.

14 The height of bone is measured from the level
15 of the alveolar crest to the floor of the
16 maxillary sinus, as illustrated by this thin
17 yellow line. This patient had only 3.9
18 millimeters of bone at baseline, certainly an
19 insufficient amount to place a dental implant.

20 On this CT scan, you can see first
21 the opening of the lateral sinus window used
22 to place the InFuse, which is right here. At

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 6 months now, 16 millimeters of bone was
2 present, a fourfold increase in the bone
3 height, which is suitable now for a dental
4 implant placement.

5 On this bar graph, the bone height
6 gained for 98 patients who received InFuse in
7 the sinus augmentation study is displayed.
8 Nearly all patients great significant amounts
9 of bone. And most patients grew a substantial
10 amount of bone.

11 The clinical trial data show that
12 InFuse induced a substantial amount of new
13 bone in sinus augmentation procedures as a
14 conclusion. InFuse averaged a gain of 8.2
15 millimeters of bone in the sinus, quite
16 comparable to the current standard of
17 treatment of a bone graft, which had a
18 slightly higher average in bone gain at 9.7
19 millimeters.

20 At six months following InFuse
21 implantation, patients were eligible to
22 receive dental implants. Per the protocol,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 only those patients who received a dental
2 implant without further augmentation were
3 considered successful and allowed to be
4 followed for functional loading. If they were
5 not, they were considered treatment failures.

6 82.8 percent of the InFuse
7 patients were considered successful in this
8 study. 79.8 percent of patients went on to
9 receive a prosthesis and were evaluated for
10 the primary objective.

11 The combined results from the
12 sinus augmentation studies exceeded the target
13 success rate of 73 percent, achieving 79.6
14 percent success at 6 months of functional
15 loading. This was also seen in the separate
16 by study analyses.

17 One of the trial's secondary
18 objectives was to compare functional loading
19 success over time between InFuse and the bone
20 graft. In this slide, the by patient
21 functional loading success rates of the bone
22 graft and InFuse groups are compared. Note

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 that both treatments are highly successful at
2 the patient level and the results are
3 maintained all the way out to 24 months.

4 In this slide, the by implant
5 functional loading success rates of the bone
6 graft and InFuse groups are compared as well.

7 Note that both treatments are also highly
8 successful at the implant level, with 87
9 percent and 86 percent of implant target sites
10 receiving implants and a prosthesis.

11 After 6 months of functional
12 loading, 81 percent of the InFuse target sites
13 remained functionally loaded compared to 84
14 percent in the bone graft. This difference is
15 not statistically significant. As shown here,
16 the results are maintained once again out to
17 24 months with no statistical difference.

18 In summary, once the dental
19 implant is placed, almost all patients
20 continued to have a successful prosthesis
21 placement and long-term functional loading in
22 both groups.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Bone density was measured on CT
2 scans in a subset of sinus augmentation
3 patients treated with bone graft and InFuse.
4 At four months, dense mature enough bone to
5 receive dental implants in both groups
6 developed.

7 The higher bone density in the
8 bone graft group is probably due to the
9 residual mineral density of the bone graft,
10 rather than new bone, which is the nature of
11 autogenous bone.

12 From the time of surgery to
13 six-month post-functional loading, which is 10
14 to 12 months from the first density
15 measurement, the bone induced by InFuse
16 becomes much more dense with loading of the
17 dental implant.

18 The bone densities now are
19 comparable in each group. These results
20 demonstrate that bone induced by InFuse
21 responds as normal physiologic bone and
22 increases in density when loaded as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 anticipated by physiologically normal bone.

2 The histologic assessments
3 demonstrated that there were no clinically
4 significant differences in histologic
5 parameters between InFuse and the bone
6 graft-induced bone. Patients in these studies
7 had a core biopsy taken, which is represented
8 by this at the dental implant placement site.

9 That's why it is cylindrical in shape.

10 These core biopsies were used to
11 make qualitative and quantitative histologic
12 assessments. This representative specimen, as
13 you see here, is taken from a patient who
14 received InFuse. It has been prepared with
15 the Goldner stain.

16 Native bone is seen at the base of
17 the longitudinal section of the core biopsy,
18 which is here. This is the native bone of the
19 maxilla that's not induced by either a bone
20 graft or InFuse. And new bone is seen above
21 this level.

22 The histologic assessment

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 demonstrates that there are no significant
2 differences in histologic parameters between
3 InFuse and the bone graft-induced bone. Both
4 resulted in significant formation of new
5 trabecular bone comparable in density and
6 structure to native bone. Sufficient bone was
7 generated for osseointegration with the
8 implant.

9 Essentially this Goldner stain
10 shows the green to be bone, the red here to be
11 marrow spaces. You see end osteo,
12 osteoblasts. You see thick trabecular
13 connectivity here, which is ideal for dental
14 implant placement.

15 By everything we measured,
16 radiographs, histology, bone density, implant
17 placements, we have demonstrated that InFuse
18 induces growth of normal physiologic bone.

19 On this next slide, we show
20 representative samples from a patient who
21 received a bone graft and another who received
22 InFuse. The bone graft is here. Again, these

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 are core biopsies. And, therefore, they're
2 cylindrical and the InFuse-induced bone on
3 this side.

4 Autogenous bone and InFuse grafted
5 sites resulted in significant formation of new
6 trabecular bone comparable in density and
7 structure to the host site. Both show similar
8 trabecular bone volume as the amount of bone
9 in this total volume space, accounting for
10 narrow spaces, and bone thickness. And both
11 had a 90 to 95 percent lamellar bone
12 architecture, indicative of mature bone. Only
13 a small amount of residual immature bone was
14 present in each group.

15 Although statistically somewhat
16 different, perhaps due to the residual and
17 lamellar bone fragments in bone grafts, the
18 difference did not affect clinical outcomes.

19 Patients were administered
20 subtherapeutic doses of tetracycline and
21 doxycycline in order to label the new bone
22 that was formed. The early osseoinductive

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 event of InFuse is easily demonstrated on the
2 far right in this pair of micrographs. This
3 is the InFuse. This is a bone graft at ten
4 days. The early yellow/green fluorochrome
5 label shows all the new bone at ten days
6 post-grafting.

7 In comparison, the fluorochrome
8 label on the pair on the left shows 80 percent
9 residual allograft and autograft fragments
10 incorporated into only about 20 percent new de
11 novo bone. Essentially the yellow
12 fluorochrome, as you see here, is indicative
13 of new bone.

14 The remaining darker bone
15 particles are nonviable, residual bone
16 particles from the graft itself; whereas, the
17 rhBMP-2/ACS-produced bone at ten days, shows a
18 remarkable amount of new de novo bone
19 formation, as illustrated by the fluorochrome
20 labeling.

21 Dr. Stephen Cook, a bioengineer
22 and professor of orthopedics at Tulane

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 University with over three decades of
2 experience in bone histology, independently
3 reviewed the data from this PMA. He
4 concluded, and I quote, "Autogenous bone and
5 rhBMP-2/ACS grafted sites resulted in
6 significant formation of new trabecular bone
7 comparable in density and structure to the
8 host site. The bone that formed was
9 biologically and structurally normal."

10 It is well-established that InFuse
11 leads to bone growth where it is surgically
12 implanted. Logically it should follow that
13 the treatment sites have sufficient bone with
14 InFuse will lead to bone growth; and, in turn,
15 will allow for dental implants and successful
16 functional loading after a prosthesis is
17 placed.

18 The data presented here were
19 derived from prospective, randomized clinical
20 trials to confirm this is the oral and
21 maxillofacial population with a protocol to
22 find success definition with an agreed-upon

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 predetermined success rate.

2 The data collected from the sinus
3 augmentation studies indicate that InFuse
4 induced new bone growth. The 79.6 percent
5 success rate in the InFuse group exceeded the
6 73 percent target success rate for
7 implantation and long-term functional loading.

8 The evidence we presented here supports the
9 efficacy of InFuse in the sinus augmentation
10 procedure.

11 This concludes the sinus
12 augmentation presentation. I thank you for
13 your attention. I now would like to introduce
14 my colleague, Dr. David Cochran, who will
15 present the extraction socket data and the
16 safety profile.

17 Thank you.

18 DR. COCHRAN: Good morning, panel
19 members. My name is Dr. David Cochran. And I
20 am a periodontist and the Chairman of the
21 Department of Periodontics at the University
22 of Texas Health Science Center in San Antonio.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 I have no direct financial
2 interest in the product under review. I am a
3 consultant for Medtronic, who is covering my
4 expenses for attending this meeting. I
5 participated in the IDE clinical trials of
6 this device as a clinical investigator.

7 It is my pleasure to present why
8 there is a clinical need for extraction defect
9 augmentation with InFuse bone graft, the
10 clinical data from this IDE clinical trial and
11 the safety data that supports the entire PMA.

12 When a patient becomes edentulous,
13 the alveolus or alveolar ridge can collapse or
14 remodel through the resorption of bone to such
15 an extent that the patient is unable to have
16 dental implants placed or receive other dental
17 restoration.

18 Dentists want to augment these
19 extraction defects with a bone graft that can
20 induce new bone formation and preserve or
21 restore the height and width of the extraction
22 socket or ridge.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 If you look in this picture on the
2 left over here, as Dr. Marx pointed out, this
3 is an extraction socket area, where half the
4 buccal plate is missing. So the defect
5 exists, particularly on the facial side here,
6 where this bone is gone.

7 On the right side here is the
8 collagen sponge, or ACS as we refer to it in
9 the trial here. And what that is, is, in this
10 case either the sponge alone or the sponge
11 with the BMP placed in that sponge. We were
12 blinded as investigators. And I'm not sure
13 which case this is.

14 When we place a dental implant, we
15 need to have a sufficient volume of bone so
16 that the bone can completely surround the
17 cylindrical implant. The implant is normally
18 3.5 millimeters in diameter. So we have to
19 have a sufficient volume in here.

20 And when you take these teeth out,
21 this bone tends to collapse here. And we
22 don't have enough width, particularly from the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 palatal to the buccal side, to place our
2 implants. So we normally need some sort of
3 grafting material to bulk out this area and
4 form new bone.

5 The clinical data was derived from
6 a prospective randomized controlled 80-patient
7 human clinical trial. This trial design is
8 considered a high level of clinical evidence.

9 Patients were eligible for
10 inclusion if they had a buccal wall defect at
11 least 50 percent of the extraction socket
12 depth for maxillary teeth from the bicuspids
13 forward.

14 Other criteria were similar to the
15 sinus program, including no active nicotine
16 use or disease or medications that affected
17 bone metabolism. The efficacy endpoint was
18 the formation of adequate alveolar bone
19 formation similar to the sinus augmentation
20 study. Safety was also evaluated similar to
21 the sinus augmentation program.

22 This is the sequence of events in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 our extraction defect augmentation study. You
2 will notice that it's very similar to the
3 sinus augmentation studies. Testing performed
4 in each study period is shown again across the
5 bottom of this line here.

6 Following the baseline period
7 extraction of the tooth and collection of the
8 initial CAT scan or CT, rhBMP-2/ACS was
9 implanted in the extraction defect to induce
10 bone.

11 Approximately four months after
12 surgery, each patient was evaluated by CT
13 scans to determine whether sufficient bone
14 formation had occurred.

15 The surgeons then reentered the
16 surgical site. We took the core biopsy here
17 for histological evaluation. And then we
18 placed our dental implants. There was then a
19 period of time to allow for osseointegration
20 of the implant.

21 After the dental implant was
22 integrated, an abutment and prosthesis was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 placed on the implant and loading of the
2 prosthesis began. The patients were assessed
3 at six-month intervals through 24 months,
4 which was approximately 36 months following
5 rhBMP-2/ACS placement.

6 Two clinical trials were conducted
7 under the extraction socket IDE. And first I
8 want to mention about the pilot study: at the
9 same time as the sinus pilot study, a
10 two-center pilot study was conducted to assess
11 the feasibility of rhBMP-2/ACS in horizontal
12 ridge augmentation in six patients and
13 extraction socket augmentation also in another
14 six patients utilizing this concentration of
15 0.43 milligrams per ml.

16 Filling of the extraction sockets
17 was seen in this study, but it was not optimal
18 for dental implant placement. So we
19 progressed, then, to a dosing study, which was
20 then conducted similar to the sinus dosing
21 study.

22 In the first cohort on the bottom

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 down here, patients were randomized to receive
2 implantation of 0.75 milligrams per ml
3 concentration of rhBMP-2, the ACS sponge alone
4 down here, or unfilled extraction defects.

5 In the second cohort, which is up
6 on the right here, patients were randomized
7 again to receive implantation of the 1.5
8 milligram per ml concentration of rhBMP-2, the
9 ACS sponge alone, or unfilled extraction
10 defects.

11 The clinical data were collected
12 in these two separate prospective clinical
13 trials. Ninety-two patients were enrolled,
14 including the randomized dosing study of
15 localized alveolar ridge augmentation with
16 buccal wall defects, referred to throughout
17 this PMA as extraction defect augmentation
18 studies. And this was conducted at eight
19 different clinical study sites.

20 The effectiveness data set
21 consists of 21 patients from one study who
22 were treated with InFuse bone graft. The

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 other patients in these studies received lower
2 concentrations of rhBMP-2/ACS; no treatment,
3 which is referred to as the unfill control; or
4 the ACS alone. Let's now review some
5 representative CAT scan data.

6 This image is a pre-implant CAT
7 scan showing an extraction socket after the
8 tooth has been removed. This is where the
9 tooth was located right here. Here you see
10 the palatal wall of the extraction socket
11 that's fairly prominent and a buccal wall
12 that's not very prominent.

13 In fact, the buccal wall is
14 missing in this area here. Normally the
15 extraction socket would be down here. So
16 you're missing the buccal wall and all of the
17 space where the tooth was removed.

18 You will see that this
19 radiolucency makes this a much more
20 challenging defect because we're missing this
21 buccal wall. So we need bone fill in a
22 vertical direction as well as in a horizontal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 direction here.

2 This next CT scan shows bone
3 growth after 16 weeks following InFuse
4 placement. De novo bone exists between the
5 missing buccal wall and the palatal wall. So
6 all of this is new bone formation here. So
7 you see the palatal wall and then the new
8 buccal wall here. Thus, horizontal and
9 vertical ridge augmentation has occurred.

10 Now, this is really an exciting
11 picture when you think about it because we
12 have really never seen pictures like this
13 before. And we don't see a lot of residual
14 graft particles or anything like that, that we
15 have had to use in the past. This is all
16 induced bone by an osseointuctive protein.

17 And when this amount of new bone
18 is formed, particularly at the coronal aspect
19 of this, the surgeon can be very confident
20 that a dental implant can be placed in this
21 patient.

22 This slide is another

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 representative CAT scan for an InFuse patient,
2 both at baseline and 16 weeks post-placement.

3 Note again the nearly absent buccal wall.
4 You just see a wisp of this bone. This bone
5 is always very thin in this area, but we're
6 completely missing the remainder of the buccal
7 wall, both in a vertical and horizontal
8 direction. And then after InFuse placement,
9 we see bone has grown and completely filled
10 this area.

11 In contrast, on this next slide,
12 we show a set of CT scans from a patient with
13 the critical size defect that was left
14 unfilled. So the tooth was in this area here.

15 This is just a radiographic marker where the
16 implant -- we would like to place it in this
17 site here. There is just a wisp of buccal
18 plate here. Here's palatal wall.

19 Once again, the buccal wall is
20 barely present. And with no treatment,
21 continued resorption takes place over the 16
22 weeks such that a dental implant cannot be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 placed. There's just no fill in this defect,
2 just a little residual amount in this area,
3 certainly not enough bone to place an implant.

4 So in this slide, we compare the
5 unfilled and the InFuse-treated CT scans. The
6 InFuse treatment on the right provided
7 clinically significant results, allowing for
8 dental implant placement, which is shown on
9 this next slide for this same patient. So the
10 patient therapeutic benefit is demonstrated
11 here by showing this dental implant that's now
12 placed in this vertically enhanced bone growth
13 with the InFuse product.

14 Now, one of the challenges we
15 faced designing this trial was that no one had
16 evaluated extraction socket defects in the
17 healing over time from the radiographic point
18 of view.

19 So we collaborated with leading
20 radiologists and determined before the trial
21 began how to evaluate the fill of these
22 defects. We took serial section CT scans and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 established where the baseline was and where
2 the most coronal extension of the bone
3 occurred.

4 So here you can see in this
5 schematic the tooth was in this area. The
6 base of the extraction socket is down here.
7 And the most coronal aspect is right here.

8 Following the bone augmentation
9 surgery, we repeated this process. And then
10 we evaluated the change in vertical height,
11 which is represented by this green line, which
12 is D-1 plus D-2. So it's this vertical line
13 going through here.

14 And we also evaluated the width of
15 the defect at the one-quarter, one-half, and
16 three-quarter positions to allow us to look at
17 the width of the bone that's grown. It's
18 better to look at it over here. So there's
19 missing bone here we could evaluate what the
20 dimension of the new bone growth was.

21 On the next three slides, I will
22 show you the results of our extraction defect

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 augmentation study for the unfilled, the ACS
2 only, and the InFuse bone graft treatment.

3 In this slide, we are looking at
4 change in bone height. And we demonstrated
5 that in an unfilled defect, which is this
6 green bar here, the alveolar ridge height in
7 the extraction defect was lost. And with the
8 ACS only, the height of the extraction socket
9 was comparable to the unfilled defect.

10 What you are looking at here is
11 change in bone height. So here you see a
12 1.17-millimeter change in bone height, which
13 means that you started here, but then you went
14 down about 1.17 millimeters. When you had the
15 collagen only or the ACS, we still lost about
16 a millimeter of bone.

17 So in both these situations, we
18 resulted in loss of bone, which is important
19 for the model system, demonstrating that these
20 are indeed critical size defects that do not
21 heal when left alone.

22 With InFuse bone graft, the height

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 of the extraction socket was preserved. In
2 other words, there's no change, really, in
3 height of the extraction socket, which means
4 that essentially it filled in all the way and
5 preserved the height of that ridge. These
6 differences confirm that InFuse bone graft
7 leads to highly significant improvements in
8 bone height.

9 Now, on this next slide, we are
10 switching from bone height to the bone width.

11 The change in the width at the one-quarter
12 position is shown for each of the treatments
13 in this slide.

14 Here one sees another significant
15 gain in bone growth with InFuse bone graft
16 versus the unfilled and the collagen only-
17 treated patients. Bone growth at the socket
18 crest is significantly greater for InFuse than
19 with either the unfilled or the ACS only
20 treatment. So you see here a very significant
21 difference, which at the one-quarter point is
22 the most important point for us clinically to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 place an implant.

2 If we don't have the width up here
3 at the coronal area, then we can't place the
4 implant because the implants are generally
5 about 3.5 millimeters in width. So we're
6 looking at this area right here.

7 Patients with InFuse bone graft
8 experience an average of 2.7 millimeters
9 additional width gain where it counts the most
10 compared to patients with unfilled defects and
11 an average of 2.45 millimeters additional
12 width gain compared to patients treated with
13 ACS only.

14 This is not only significant, but
15 it is clinically relevant because bone is
16 needed here to place the implants and/or
17 support aesthetic restorations for our
18 patients; similarly, if we see another
19 significant gain in bone width at the one-half
20 position of the extraction defect, again with
21 the InFuse bone graft compared to the unfilled
22 or the ACS only.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Patients with InFuse experience an
2 average of 2.35 millimeters additional width
3 gain compared to patients with unfilled and an
4 average of 2.18 millimeters additional width
5 gain compared to patients treated with the
6 sponge only. As mentioned previously, this
7 bone is required to clinically support the
8 dental implant restorations.

9 In summary, InFuse is
10 significantly more effective than both
11 unfilled and ACS only in terms of the change
12 in alveolar ridge height and in width at the
13 one-half and one-quarter measurement
14 positions. This is clinically relevant for
15 our patients because the more bone, the better
16 the chance of implant placement and/or
17 prosthesis success in long-term function of
18 the dental restoration.

19 So the take-home message is that
20 over time InFuse induced bone growth in
21 significant and sufficient quantities to place
22 dental implants in this study.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Now, the previous slide showed the
2 effectiveness of InFuse for promoting bone
3 growth. This slide shows the relationship
4 between bone growth and implant success. So
5 we're looking at the relationship between that
6 bone growth that occurred and the success of
7 implant placement.

8 The data demonstrate that the
9 amount of bone growth is strongly associated
10 with successful dental implant placement.
11 Bone width gained at the one-half and
12 one-quarter positions and the increases in
13 bone height are all significantly associated
14 with implant success. Only at the
15 three-quarter position is there no association
16 in bone width in either treatment with implant
17 success.

18 In summary, greater bone growth is
19 associated with greater implant placement, a
20 logical conclusion, but this was proven in
21 this controlled clinical trial.

22 This slide shows the comparison of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 dental implant, prosthesis placement, and
2 functional loading by patient for the
3 unfilled, ACS, and InFuse patients. The
4 important time points are the dental implant
5 placement along this line here and the
6 six-month evaluation point along that line.
7 On the next two slides, we will look more
8 closely at these two time points.

9 When we compare dental implant
10 placement, prosthesis placement, and
11 functional loading by patient for unfilled,
12 ACS, and InFuse patients, we find significant
13 differences between the groups. Of the
14 patients that reached this phase of the study
15 without the need for further augmentation, 86
16 percent of the InFuse grafted patients
17 received a dental implant versus only 59
18 percent in the ACS only patients and 47
19 percent of the unfilled patients.

20 The greater bone growth achieved
21 with the use of InFuse bone graft led to a
22 significantly greater number of patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 receiving dental implants without an
2 additional augmentation procedure.

3 Focusing on the functional loading
4 at six months, functional loading at six
5 months, we find that a significantly greater
6 percentage of patients in the InFuse group, 74
7 percent, remained functionally loaded than
8 those in the unfilled group at 38 percent.

9 Similarly, a higher percentage of
10 the InFuse patients remained functionally
11 loaded at six months compared to those in the
12 ACS only group at 50 percent, although this
13 difference is not statistically significant.

14 The conclusion from these
15 effectiveness data is that InFuse bone graft
16 successfully outperformed both ACS only and
17 unfilled groups.

18 In this next slide are
19 representative core biopsies taken at the time
20 of dental implant placement from both the
21 extraction socket augmentation study, which is
22 on the left here; and the sinus augmentation

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 studies on the right.

2 This shows the new bone induced by
3 InFuse is nearly identical in every measurable
4 parameter to the native bone in terms of
5 trabecular volume, thickness, and number. And
6 this is the same, whether it's in the sinus
7 augmentation or in the extraction defect
8 augmentation studies.

9 This higher magnification section
10 shows the normal mix of lamellar and immature
11 bone pattern, which is indicative of maturing
12 bone produced by InFuse bone graft.

13 On this slide, the density of the
14 induced bone is compared in the extraction,
15 augmentation, and sinus augmentation studies,
16 the sinus on the right, extraction on the
17 left.

18 For the infused patients, the
19 purple bars, which we're looking at here, the
20 purple bars, on the right side of this slide,
21 in the sinus study, the density was 137
22 milligrams per cc at 4 months post-grafting

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 and 508 milligrams per cc 6 months
2 post-functional loading.

3 On the left, this graft shows a
4 mean density of 343 milligrams per cc at 4
5 months post-grafting in the extraction socket
6 augmentation study. The density achieved in
7 the extraction socket study was well on its
8 way to being comparable to the sinus study
9 results.

10 In summary, the clinical data from
11 the extraction defect augmentation studies
12 demonstrate that InFuse bone graft induces new
13 bone growth that leads to successful dental
14 implantation and long-term functional loading.

15 InFuse was most clinically effective
16 following tooth extraction for augmentation of
17 the alveolar ridge and dental restoration.

18 InFuse provides a new treatment
19 modality and a treatment alternative since
20 bone growth is stimulated by this
21 osteoinductive protein and really gives us
22 another option for our patients.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Now let's turn our attention to
2 the overall safety of InFuse bone graft. One
3 of the things that is prominent when you look
4 at the data related to the InFuse product is
5 that there is an extensive safety profile.
6 There are already two approved PMAs. And over
7 300,000 InFuse bone graft kits have been
8 distributed. Four hundred thirty-seven
9 patients support InFuse safety and
10 effectiveness in the two PMAs.

11 InFuse has more level I clinical
12 evidence than any other bone grafting agent.
13 There are 1,070 patients enrolled in
14 rigorously controlled FDA clinical trials.
15 InFuse has an established safety profile.

16 The oral and maxillofacial safety
17 data set consist of patients who were
18 implanted with any concentration of rhBMP-2
19 plus ACS sponge and sinus augmentation or
20 extraction defect studies.

21 The population includes 129
22 patients from the sinus augmentation studies

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 and 55 patients from the extraction defect
2 studies, for a total of 184 patients at all
3 concentrations of rhBMP-2/ACS. Of these, 120
4 represent the InFuse patients.

5 The table on the next two slides
6 shows the adverse events which occurred in
7 more than ten percent of the patients who
8 received InFuse or bone graft in the three
9 studies.

10 The majority of the events were
11 expected in oral surgical patients or patients
12 who undergo bone harvest procedures. They
13 included oral, facial, and general edema,
14 infection, mouth pain, arthralgia, and
15 abnormal gait. These events resolved in short
16 order.

17 When all of these categories of
18 adverse events in these two groups were
19 compared, only the AEs presented on this slide
20 were significantly different in the two
21 patient populations. These included
22 arthralgia, abnormal gait, hypoesthesia,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 erythema, general edema, infection, pain, and
2 hyperglycemia. All were significantly less
3 frequent in the patients treated with InFuse
4 bone graft.

5 This table shows the number of
6 adverse events reported in the IDE studies,
7 the percentages of patients who experienced at
8 least one adverse event, the relatedness of
9 the events to InFuse bone graft, plus the
10 number of grade 3 and grade 4 events.

11 The collection of adverse events
12 was very conservative and documented every
13 conceivable AE that patients experienced. As
14 you can see, virtually every patient reported
15 at least one adverse event. And there was no
16 difference between the InFuse and bone graft
17 groups.

18 The incidence of adverse events
19 related to InFuse was 17 percent and 24
20 percent in the 2 indications. These were
21 predominantly facial edema, oral edema, mouth
22 pain, and oral erythema, which rapidly

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 resolved.

2 The adverse events related to
3 autogenous bone graft were not recorded as a
4 part of this study, but I will present the
5 data related to the donor site in just a
6 moment.

7 The majority of adverse events
8 were grade I or grade II. The rates of grade
9 III or IV adverse events in both groups was
10 low. None of the grade III or IV adverse
11 events was related to InFuse bone graft.

12 There was one death among the 312
13 participants in these studies. The patient
14 was a 43-year-old woman at my site in San
15 Antonio who underwent an extraction socket
16 augmentation procedure. The operation was
17 uneventful and only had some expected mild
18 facial swelling in the immediate postoperative
19 period and no significant other adverse
20 events. The patient died three years
21 postoperatively. The cause of death was
22 judged not related to the study treatment by

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the investigator, which was me.

2 Antibody titers were measured
3 preoperatively and postoperatively in 184
4 patients who received rhBMP-2/ACS in 91
5 autogenous bone graft patients to monitor for
6 immune reactions to the components of InFuse.

7 There was no incidence of positive
8 anti-rhBMP-2 antibodies in the autogenous bone
9 graft and a 2.2 percent incidence in the
10 InFuse group. The titers, however, were low
11 and transient.

12 There was a 20 percent incidence
13 of antibodies to bovine collagen in the InFuse
14 patients. But, interestingly, the autogenous
15 bone graft group had an even higher incidence
16 at 31 percent, presumably due to exposure from
17 other bovine sources. Some titers continued
18 into the follow-up period.

19 Antibodies to human type I
20 collagen were not detected in either group of
21 patients. The presence of antibodies to
22 rhBMP-2 or bovine collagen was not associated

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 with immune-mediated adverse events, such as
2 allergic reactions.

3 As previously mentioned, there
4 were 91 patients who received autogenous bone
5 graft. And that bone graft was predominantly
6 harvested from three different areas: the
7 Iliac crest, the tibial plateau, and an
8 intra-oral bone site.

9 The pain and morbidity associated
10 with the harvest site is shown in this table.

11 Significant pain is experienced in a high
12 percentage of patients with all three bone
13 harvesting techniques.

14 Iliac crest harvesting is
15 associated with significant donor site pain
16 and gait disturbance out to ten days. Sensory
17 loss was present in 11.1 percent of these
18 patients, even at 6 months.

19 Even in the tibial plateau site,
20 there was significant pain and associated gait
21 disturbance, with pain present in 3.1 percent
22 of the patients out to 6 months.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 At 2 days, 27.6 percent of
2 patients experienced local sensory loss in the
3 perioral region, the gingiva, the lip, or in
4 the teeth, which can become nonviable due to
5 the harvesting of the intraoral bone. Sensory
6 loss was still present in 17 percent of these
7 patients at 6 months.

8 In summary, there is a large
9 percentage of patients that have problems with
10 these donor sites. The use of InFuse will
11 eliminate the morbidity associated with
12 autogenous bone graft harvesting.

13 This graph shows the time line of
14 the harvest site adverse events. The majority
15 of events occurred in the first 20 days after
16 the harvest procedure and dropped off
17 precipitously by 60 days with some morbidity
18 extending beyond 180 days. This was not
19 relevant to the InFuse group as these patients
20 did not undergo a harvest procedure.

21 In summary, InFuse bone graft has
22 an established safety profile through two

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 previous PMA approvals in more than 1,000
2 patients enrolled in rigorous FDA prospective
3 randomized human clinical trials.

4 In addition, the clinical data
5 from three prospective randomized IDE trials
6 that specifically evaluated its use in oral
7 and maxillofacial applications demonstrated
8 that there were significantly fewer adverse
9 events with InFuse bone graft than with
10 autogenous bone graft.

11 The use of InFuse as an
12 alternative to autogenous bone graft
13 eliminates the significant morbidity
14 associated with autogenous bone harvesting.
15 We believe that there is reasonable assurance
16 that InFuse is safe for these indications for
17 use.

18 I would now like to turn the
19 podium back to Ed Chin from Medtronic to
20 conclude our presentation.

21 DR. CHIN: Thank you, Dr. Cochran.

22 Members of the panel, based on the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 information presented to you today and the
2 information submitted in the PMA application,
3 we have demonstrated that there is a
4 reasonable assurance that the safety and
5 effectiveness for the use of InFuse bone graft
6 as an alternative to autogenous bone graft for
7 sinus augmentations and localized alveolar
8 ridge augmentations for defects associated
9 with extraction sockets.

10 Both indications are supported by
11 clinical data from prospective randomized
12 controlled FDA-approved clinical trials. We
13 believe our studies demonstrate that these
14 patients share a common clinical problem; that
15 is, the need to grow sufficient bone to
16 support the dental implant borne restoration,
17 bone that will respond to physiologic loading
18 over time.

19 Our studies of both the sinus
20 augmentation and extraction socket indications
21 demonstrated that InFuse induces significant
22 amounts of bone. The bone is histologically

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 normal bone. The bone was sufficient to allow
2 the dental implant and subsequent functional
3 loading over time. We believe that InFuse
4 produced nearly identical results.

5 When reviewing these studies and
6 analyzing the numbers, it is sometimes easy to
7 forget the significant clinical benefits these
8 patients derived from this technology.

9 This patient lost teeth and
10 subsequently the supporting bone.
11 Reconstruction using InFuse provided the bony
12 support for dental implantation and eventual
13 prosthetic restoration. This was accomplished
14 by eliminating the autogenous bone graft
15 harvest procedure and the associated pain and
16 morbidity with the harvest.

17 These non-restorable teeth were
18 extracted and the sockets were grafted with
19 InFuse, as shown in this slide. This
20 aesthetically pleasing and functional result
21 is only possible when adequate bone is grown
22 in the extraction socket, again eliminating

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the bone harvesting, the pain, and morbidity.

2 InFuse was effective in inducing
3 new bone in sinus augmentation and extraction
4 socket with buccal wall defects augmentation
5 procedures.

6 Our clinical studies demonstrated
7 that InFuse induced normal bone where no bone
8 existed before and this bone was
9 histologically and physiologically normal.
10 This bone responded to functional loading
11 stresses and supported dental implants under
12 physiologic loading conditions out to three
13 years.

14 The clinical data demonstrate that
15 InFuse is clinically effective in a
16 significant portion of the patient population
17 and has been shown to be an effective
18 alternative to bone graft.

19 We believe that we have
20 established that InFuse is safe. The
21 multitude of preclinical in vitro and in vivo
22 studies and extensive human clinical trials

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 performed on the product attest to the
2 established safety profile of InFuse.

3 Clearly the risks associated with
4 the use of InFuse are outweighed by the
5 benefits of the device, particularly when the
6 use of the device avoids the significant
7 morbidity associated with bone harvesting
8 procedures and/or the general anesthesia.

9 As clearly demonstrated in these
10 presentations and the information submitted in
11 the PMA application, a reasonable assurance of
12 safety and effectiveness of InFuse has been
13 provided. Functional animal model testing,
14 clinical data from two previously approved
15 PMAs, two large-scale IDE studies demonstrate
16 InFuse safety stimulates the formation of
17 bone.

18 The data are consistent. The data
19 are compelling. They are convincing. InFuse
20 can safely grow normal bone where none existed
21 before and is an effective alternative to a
22 bone graft. These data provide a reasonable

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 assurance that the device is safe and
2 effective for its intended use, the main
3 criterion for PMA approval.

4 We ask that you as members of the
5 panel acknowledge the significance and
6 validity of the information and make this
7 breakthrough technology available to surgeons
8 and their patients by recommending approval of
9 this PMA application.

10 This concludes Medtronic's
11 presentations. And we are available for
12 further questions.

13 CHAIRMAN BURTON: Thank you, Dr.
14 Chin, Dr. Cochran, and Dr. Marx.

15 At this time I would like to ask
16 the panel if there are any points of
17 clarification that they would like to have
18 from the three presentations that they would
19 like to have or we can also call them back
20 during our discussions later this morning and
21 in the afternoon, but I would certainly like
22 to entertain any questions at this time while

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the information is fresh. Yes, Dr. O'Brien?

2 MEMBER O'BRIEN: You may have
3 covered this already in your presentation, but
4 do you use antibiotics as part of the InFuse
5 procedure? Is that necessary for the
6 procedure or do you have any adverse effects
7 that you would like to avoid with antibiotics
8 or are they helpful in preventing adverse
9 effects?

10 DR. MARX: As per the protocol,
11 the individual site investigators were allowed
12 to use their standard antibiotic regimens in
13 both the bone graft groups and in the
14 treatment groups. No adverse reactions were
15 reported to the antibiotics related to the
16 study.

17 CHAIRMAN BURTON: Dr. Zuniga?

18 MEMBER ZUNIGA: I will address
19 this to Dr. Marx. In your presentation, you
20 introduced the tetracycline staining
21 protocols. Can you clarify, was that done in
22 the dosing study or the pivotal study? And

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 was there any quantitative analysis of that
2 tetracycline staining?

3 DR. MARX: The tetracycline
4 staining was done in both the dosing study and
5 the pivotal study, more in the pivotal study
6 because there were greater numbers. They were
7 given to a subset of the individuals.

8 The quantitation of that was not
9 recorded. There was mainly a qualitative
10 study with that only.

11 CHAIRMAN BURTON: Yes, Dr. Amar?

12 MEMBER AMAR: Was there any
13 attempt to look at a demineralized core biopsy
14 in terms of looking at the mineral content of
15 a recombinant human BMP-2 graft site, as
16 opposed to either original sites or site
17 grafted with DFDBA? I guess my question is,
18 was there any attempt to look at mineral
19 content and how would that progress over time?

20 DR. MARX: That was not part of
21 the original protocol to look at mineral
22 content. So that was accomplished in a few

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 patients, but it was not part of the original
2 protocol. Therefore, that data was not
3 presented.

4 MEMBER AMAR: Is there any data to
5 support that mineral content remains over
6 there? Because from what I saw, it is the
7 only demineralized section. Am I correct?

8 DR. MARX: Dr. Cook is behind me,
9 who was the histologic investigator. I think
10 he can address that question better than I.

11 DR. COOK: Yes. I'm Stephen Cook.
12 I'm a bioengineer. I'm professor of
13 orthopedic surgery at Tulane University. I am
14 a consultant for Medtronic, who are paying my
15 expenses to attend this meeting. I acted as
16 an independent histologic reviewer of the
17 sections and reports from all of these
18 studies. And that was my role in the project.

19 If I can go back maybe to your
20 first question, which was related to the
21 tetracycline labeling protocol. There was
22 quantification done on acquisition rates as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 part of the quantification of the sections in
2 the sinus pivotal study.

3 So there is data from the
4 tetracycline labels that was used more than
5 just looking at the change over time of the
6 mineralization, but also there was actually
7 some quantitative data that was performed.
8 And I believe that was in the packet of
9 information.

10 The second question was on
11 mineralization of the bone. And there were
12 bone density-type measurements that were shown
13 in the presentations by Dr. Marx and Dr.
14 Cochran based on CT examination of bone
15 density. And what that showed was that
16 earlier, at the four-month time period, when
17 the CTs were formed, mineralization was
18 slightly lower in the rhBMP-2 group.

19 But you have to remember that in
20 the autogenous bone group at four months and
21 indeed in the histologic sections, from the
22 biopsies that were generally 6 to 12 months

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 after placement of the graft, there was still
2 a significant amount of graft present in the
3 sections themselves.

4 So in the CT examination of bone
5 density, you're getting a false sense of
6 mineralization of new bone because you're
7 picking up the information, the mineralization
8 from the residual bone graft.

9 So although there were differences
10 at four months, as you got into later time
11 periods, as both the bone matured in both
12 groups, bone graft was reincorporated, they
13 became equivalent.

14 CHAIRMAN BURTON: Dr. Janosky?

15 MEMBER JANOSKY: I will direct my
16 questions to Dr. Cochran. You presented some
17 data about extraction socket studies. Am I
18 correct that the n is 21 that's in those
19 studies?

20 DR. COCHRAN: For the 1.5
21 milligram per ml concentration. It was an
22 80-patient study --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 MEMBER JANOSKY: Exactly.

2 DR. COCHRAN: -- with multiple
3 arms, yes.

4 MEMBER JANOSKY: Okay. For the
5 patients the n equals 21, what could you tell
6 us about either the provider data and/or the
7 patient data? I'm very interested in the mix
8 of providers and who provided that, training
9 and level, as well as patients.

10 You didn't present any data today
11 that shows us the heterogeneity of either the
12 patients or the providers. So what could you
13 tell us about that n equals 21, please?

14 DR. COCHRAN: I think that data
15 was analyzed between the different sites. I
16 think that's really what you're asking about.

17 Of the eight different sites that were
18 involved, what was the statistical variation
19 between the eight different sites? I have to
20 refer that to the statistical evaluator to do
21 that.

22 From a clinical investigator,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 ourselves, we did do training as a group of
2 investigators. So we did get together on
3 multiple occasions. And we did standardize
4 the procedure as much as we possibly could.
5 So the standardization from an investigator
6 training point of view was done in group
7 meetings.

8 I'll have to get the statistician
9 to give you the variability between the sites.

10 MEMBER JANOSKY: I'm also
11 interested in the n for each of those. So if
12 you have n equals 21, how many providers did
13 that represent?

14 DR. COCHRAN: Yes. I don't know
15 that information. Do you? Yes. We'll have
16 to get that. We'll look it up and get it back
17 to you.

18 CHAIRMAN BURTON: Dr. Li?

19 MEMBER LI: My question is for Dr.
20 Marx. The results of your pivotal study
21 showed that the success rate for the bone
22 graft group was fairly stable between the 6

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 months and 24 months, which was 90.8
2 throughout that period. But for the
3 rhBMP-2/ACS group, it declined slightly but
4 fairly consistently, which was 79 percent down
5 to 76 percent. It's slight, but each period
6 was a little bit lower.

7 Do you have any data beyond 24
8 months which indicates any trend that
9 continued or --

10 DR. MARX: The study concluded at
11 24 months. We don't have data beyond 24
12 months. The patients who declined were a
13 mixture. Many of them were dropouts whom we
14 couldn't get back for follow-up. And so
15 although the success rate declined, they were
16 successful up until the point we lost them to
17 follow-up. That explained a number of them.

18 But beyond 24 months, the study
19 was extinguished. So we don't have regular
20 follow-up on those. Many of them were
21 followed up outside the study, but that could
22 not be recorded.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 MEMBER LI: Okay. Thank you. And
2 I have a second question for Dr. Cochran. For
3 the extraction socket study, the completion
4 rate was 37.5 percent. You started with 80
5 subjects, right, total. Then the majority of
6 them dropped out or did not complete.

7 Do you have any thoughts on that,
8 the possible impact on the results?

9 DR. COCHRAN: Unfortunately, I can
10 probably tell you why that happened. The
11 problem was we provided the treatment up
12 front. We gave these patients new teeth, in
13 these cases really nice implant restorations.
14 And they were very satisfied with that
15 restoration. So they didn't like coming back
16 for the follow-up exams because they were
17 pretty satisfied. And that's a problem in
18 these studies, especially when you're
19 stretching it out pretty far.

20 And we saw these patients a lot.
21 I mean, we wanted to make sure that there were
22 no adverse events. And so Wyeth or Medtron --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 at that point it was Wyeth was very on top of
2 the study to make sure that we were following
3 these patients as much as possible.

4 But a lot of them, we wrote
5 letters. We sent certified letters to the
6 people. We would call them. And after a
7 while, they were pretty happy. And they,
8 unfortunately, didn't come back too often.

9 MEMBER LI: Thank you.

10 CHAIRMAN BURTON: Dr. Zuniga?

11 MEMBER ZUNIGA: I need to make a
12 follow-up question to Dr. Cook, I believe, who
13 earlier mentioned something about the
14 mineralization in tetracycline studies.

15 I thought I heard that both of
16 these analyses were done at the dosing and the
17 pivotal study. And, if so, was there a
18 quantitative difference between the two
19 dosages? And were there standard error or
20 standard deviation variances between the
21 autogenous bone graft in either of the BMP --

22 DR. COOK: In the earlier studies,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 which the dosing was with much smaller
2 numbers, they were qualitative assessments
3 made, rather than quantitative, in the pivotal
4 study when they n's were approximately 90 in
5 each group is where the quantification
6 actually took place and statistical
7 evaluations were performed. Qualitatively
8 assessments were done on a zero to three type
9 of scale in the pilot study as well as in the
10 dosing study.

11 MEMBER ZUNIGA: Was there much
12 variance, then, in the 1.5 milligrams per ml
13 in the autogenous bone graft comparators in
14 this scale?

15 DR. COOK: There's a vast
16 difference in the way the bone is formed in
17 the rhBMP group. It's a de novo bone
18 formation occurring very early. As you saw in
19 the slide that was presented, I believe, in
20 Dr. Marx's presentation, the ten-day stain,
21 the ten-day label, there was extensive
22 networking of new bone formed at that early

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 time period; whereas, in the autogenous group,
2 it was more of a connection of the bone
3 fragments to each other so that it formed a
4 network connecting the fragments themselves,
5 rather than a new network of de novo bone
6 formation; again, forming bone where there was
7 none present versus connecting pieces of
8 viable fragments.

9 CHAIRMAN BURTON: Dr. Patters?

10 DR. COCHRAN: Hang on one second,
11 Mark, if you don't mind. We found the data
12 for the 21 patients that you asked about a
13 little bit earlier.

14 Two of the sites had five patients
15 in that group each. Two sites had three
16 patients. Two sites had two patients. One
17 site had one patient. And one site didn't
18 have any patients in that group.

19 MEMBER JANOSKY: So am I correct
20 in there were 6 diverse sites representing 21
21 patients?

22 DR. COCHRAN: Seven, I believe.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 MEMBER JANOSKY: Seven. And the
2 largest number of patients treated at any one
3 site was?

4 DR. COCHRAN: Five.

5 MEMBER JANOSKY: Was five.

6 DR. COCHRAN: At two different
7 sites.

8 MEMBER JANOSKY: Okay. Thank you.

9 MEMBER PATTERS: Question for Dr.
10 Chin and Dr. Cochran. Dr. Chin, the proposed
11 indications that you have revolve around that
12 InFuse is an alternative to autogenous bone
13 grafts.

14 DR. CHIN: Yes.

15 MEMBER PATTERS: And you provide
16 direct data comparing autogenous bone grafts
17 and InFuse sinus augmentation with the pivotal
18 study.

19 DR. CHIN: That is correct.

20 MEMBER PATTERS: But that data
21 seems to be lacking with regard to extraction
22 sockets. And that's more implied, rather than

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 shown with the pivotal study. So what is the
2 justification for lumping those two together
3 under the indication of an alternative to
4 autogenous bone graft?

5 DR. CHIN: Okay. There are some
6 statistical implications there. So I would
7 maybe get some assistance there. But the data
8 for the sinus and the extraction socket are
9 looking at the 1.5-milligram concentration,
10 which is the commercial version of that
11 product today. So the analyses are done
12 comparing the 1.5 concentrations for that
13 efficacy.

14 DR. COCHRAN: I think I can
15 address that as well. In the sinus studies,
16 what we wanted to do was to take the standard
17 of care that existed at that time and compare
18 the treatment of InFuse bone graft to that
19 standard.

20 So at the times that we were doing
21 these studies, the standard there was
22 autogenous bone graft procedures mixed

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 occasionally with DFDBA as an extender or some
2 other type of material. So that was the
3 standard that we were trying to go against to
4 see what we could do as far as the standard of
5 care.

6 In the extraction socket defects,
7 the standard of care doesn't include bone
8 graft procedures, particularly at that time.
9 Most of the time when teeth are extracted,
10 there is nothing done in those cases. And so
11 that's why we did the comparison to the unfill
12 treatment because, really, that is the
13 standard of care.

14 We thought it was important to
15 have also as a control the carrier alone for
16 the BMP-2. So we had the unfilled cohort.
17 And then we had the collagen treatment alone
18 because we thought that was the most
19 scientific rigorous way to do that in the
20 carrier plus the protein. And, really, the
21 data showed exactly what we had hoped it did,
22 that the unfilled didn't form bone and that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the collagen formed a little bit but not
2 significantly. And it was very clear that the
3 BMP presence made a significant difference in
4 the outcome of that trial.

5 MEMBER PATTERS: Thank you, Dr.
6 Cochran.

7 I think I appreciate that. My
8 concern is really the indications as stated
9 here. I think that your data would clearly
10 show that InFuse has advantages over the
11 collagen carrier, but my question is, is there
12 any direct data to show that it has advantages
13 over autogenous bone being that's the
14 indication that the sponsor is looking for?

15 DR. COCHRAN: Yes. I would again
16 say that we didn't do that in the trial. So
17 we don't have that. But we can certainly
18 speculate on that. And I would speculate that
19 if you look at the histological specimen and
20 you look at the CAT scan data, you see the
21 growth of new bone in those areas.

22 People that do use some type of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 grafting material in extraction sockets, most
2 of the data indicates that in those cases,
3 these other materials, these osteoconductive
4 or really alloplastic materials, actually get
5 in the way of bone formation if you take
6 biopsies of that. And so that material will
7 stay around in the extraction defects. And,
8 really, it's not an optimal treatment for
9 placing implants.

10 The beauty of using a biological
11 growth factor like this is that you induce de
12 novo bone formation. So you don't have those
13 residual particles in that.

14 I thought the histology was really
15 nice with the fluorochrome stablin that showed
16 the particles of the autogenous bone stayed in
17 that area. And we really don't want to put
18 our implants in that type of bone.

19 Thanks.

20 CHAIRMAN BURTON: Dr. Diamond?

21 MEMBER DIAMOND: Thank you.

22 Dr. Cochran, don't sit down yet.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701