

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16

Appendix A

Standard Operating Procedures for the LLNA: DA Test Method

A1 Standard Operating Procedures/Protocol for the LLNA: DA Test MethodA-3
**A2 Results in the LLNA: DA Test Method for 1% Sodium Lauryl Sulfate (SLS)
Pretreatment versus without 1% SLS Pretreatment.....A-21**

17
18
19
20
21
22
23
24
25
26
27
28
29
30

[This Page Intentionally Left Blank]

31

32

33

34

35

36

37

38

39

40

41

Appendix A1

42

Standard Operating Procedures/Protocol for the LLNA: DA Test Method

43

44

44
45
46
47
48
49
50
51
52
53
54
55
56
57
58

[This Page Intentionally Left Blank]

59 These are the standard operating procedures performed during the two-phased interlaboratory
 60 test method validation study (Omori et al. 2008) for the murine local lymph node assay
 61 (LLNA) modified by Daicel Chemical Industries, Ltd., based on adenosine triphosphate
 62 content (ATP; referred to hereafter as the “LLNA: DA”) as confirmed by the LLNA: DA
 63 Validation Committee and provided by the study director.¹ These procedures are intended for
 64 tests conducted to evaluate a single test substance. Although the standard operating
 65 procedures detailed herein are specific for the interlaboratory test method validation study,
 66 the substances tested in the intralaboratory validation study followed a technically similar
 67 LLNA: DA test method protocol (Idehara et al. 2008; Idehara unpublished data).

68 **1.0 Preparation of Equipment and Materials**

69 Prepare the experimental equipment, materials, and reagents given in **Table A-1**.
 70 Luminometer tubes, 15 mL test tubes, 50 mL test tubes, petri dishes, and slide glass should
 71 be disposable. The underlined items will be provided by the LLNA: DA Validation
 72 Committee but in some cases, a luminometer will be furnished by the test facilities. All other
 73 materials will be provided by the test facilities.

74 **Table A-1 List of Required Equipment, Materials and Reagents**

Name of Equipment, Material, or Reagent	Manufacturer	Comment (Trade Name, Model Number, etc.)
<u>Luminometer</u>	Kikkoman Corporation, Japan	LUMITESTER C-100 Detection Range: 4×10^{-12} – 1×10^{-6} M Upper Limit: 1,000,000 RLU
<u>Luminometer tubes</u>	Kikkoman Corporation, Japan	Polypropylene, sterilized
<u>15 mL test tubes</u>	IWAKI brand	Polypropylene, sterilized
<u>50 mL test tubes</u>	IWAKI brand	Polypropylene, sterilized
<u>Petri dish</u>	Corning Incorporated	Cell culture dish, sterilized
<u>Cell scraper</u>	Costar brand	Disposable cell scraper, sterilized
<u>Slide glass</u>	Matsunami	Micro slide glass
<u>Vortex mixer</u>		
Analytical balance		For body weight measurements (readability of at least 0.1 g)

¹ 2/6/2006: Confirmed by LLNA: DA Validation Committee; 2/17/2006: Revised by Takashi Omori; 2/19/2006: Revised by Takashi Omori; 3/27/2006: Revised by Takashi Omori; 4/2/2006: Revised by Takashi Omori; 12/2/2006: Revised by Takashi Omori.

Name of Equipment, Material, or Reagent	Manufacturer	Comment (Trade Name, Model Number, etc.)
Analytical balance		For lymph node weight measurements (readability of at least 0.1 mg)
<u>Brush</u>	Ikkyuen	Osho
<u>Phosphate buffered saline</u>	Invitrogen Gibco™	pH 7.2, sterilized
<u>Luciferin-luciferase reagent</u>	Kikkoman Corporation, Japan	CheckLite™ 250 Plus ¹
Cages		Capable of housing four mice, with feed and water dispensers
Micropipette		For applying test solutions (25 µL), handling phosphate buffered saline (1000 µL), tissue suspension (20 µL), cell suspension (100 µL), and dissolved Luciferin-luciferase solution (100 µL)
Micropipette tips		Sterilized
Dissecting instruments		Large and small tweezers, scissors, surgical holder, injection needle and holder
Timer		With second display
General laboratory materials		Cotton, antiseptic solution, paper towel, clean sheet, test tube rack, microtube rack

75 Abbreviations: etc. = et cetera; g = grams; M = molar; mg = milligrams; µL = microliter; mL = milliliter; RLU
76 = relative luminescence units.

77 ¹For the substances tested in the intralaboratory validation study by Daicel Chemical Industries, Ltd. (Idehara et
78 al. 2008), the ATP content for potassium dichromate was measured by the CheckLite™ 250 Plus Kit
79 (Kikkoman Corporation, Japan) but that for all other substances was determined using the ViaLight® HS Kit
80 (Lonza Rockland, Inc., USA).

81 2.0 Preparations Prior to Delivery of Animals

82 The animals to be used in the tests are young adult female mice (nulliparous and non-
83 pregnant) of the CBA/JNCrlj strain, aged between eight to twelve weeks prior to application
84 of test and control substances. The animals will be provided by the LLNA: DA Validation
85 Committee. Preparations should be made according to the standards of the test facilities to
86 begin acclimatizing the animals once they have arrived on the previously agreed upon date of
87 delivery.

88 Six cages capable of holding four animals each should be prepared prior to the end of
 89 acclimatization.² The cages should be labeled as listed in **Table A-2**. The symbol “X”
 90 represents the code of the test substance to be provided. Mark the label using the letter
 91 indicated on the datasheets provided prior to the test. The animal test group numbers are also
 92 indicated on the datasheets. The numbers should be confirmed and the cages labeled with
 93 care. This test will be performed two or three times, so it is important to include the test
 94 number on the labels.

95 **Table A-2 Preparation of Test Group Cages**

Test Group Number	Label
Group 1	Acetone: Olive Oil (4:1)
Group 2	Positive Control
Group 3	Vehicle
Group 4	Test Substance “X” – Low Concentration
Group 5	Test Substance “X” – Medium Concentration
Group 6	Test Substance “X” – High Concentration

96 “X” represents the code of the test substance provided by the study management team.
 97

98 **3.0 Delivery, Acclimatization and Animal Assignment**

99 On the date of delivery, 25 animals will arrive and acclimatization should begin immediately.
 100 Acclimatization should be performed according to the standards of the test facilities. The
 101 animals should be acclimatized for at least five days, but no more than 16 days.

102 After acclimatization healthy animals with no observable skin lesions or other abnormalities
 103 should be randomly assigned to six groups of four animals each using randomly generated
 104 numbers. After assigning the animals to groups, four animals each should be placed in the six
 105 cages prepared as described in **Section 2.0**. Any animals remaining after the assignment of
 106 24 should be omitted from the test. Should there be fewer than 24 animals with no observed
 107 abnormalities, three animals should be assigned to each group beginning with the test group
 108 with the highest number until all of the animals are assigned.

² For the substances tested in the intralaboratory validation study by Daicel Chemical Industries, Ltd. (Idehara et al. 2008; Idehara unpublished data), at least three animals per dose group were used (i.e., in most cases, 4 animals per control group and three animals per test substance group).

109 From the delivery of the animals to the end of the test procedures the temperature of the
110 animal housing facility should be maintained at 22°C ($\pm 3^\circ\text{C}$) with a relative humidity of 30-
111 70%. The animals should be housed with a light: dark cycle of 12 hours light: 12 hours dark
112 and should be given food and water *ad libitum*. Any deviations from the standard housing
113 and feeding procedures should be recorded.

114 **4.0 Confirmation of Test Materials**

115 When the test materials sent by the LLNA: DA Validation Committee arrive, confirm that
116 the inventory document matches the contents.

117 The labels for each of the treatments (acetone: olive oil [4:1], positive control, vehicle, and
118 low, medium and high concentrations of test substances) include a test substance code and a
119 group number. After confirming that these codes match the datasheet, arrange the treatments
120 in a test tube rack according to group number. Sodium lauryl sulfate (SLS) solution will
121 arrive in one tube. Apportion 3 mL of SLS solution to each of the accompanying empty test
122 tubes, mark each tube with the group number, and arrange the tubes in order in the test tube
123 rack.

124 The treatments should be refrigerated immediately and only removed when beginning the
125 test. Refrigeration of the solutions used in these procedures should be between 0-10°C, and
126 preferably between 2-8°C, except when instructed differently. Should there be specific
127 instructions as to the handling of the solutions, the instructions will be included with the
128 materials shipment and they should be followed. For instance:

- 129 • SLS (CASRN: 151-21-3) is a 1% solution and should be kept at room
130 temperature
- 131 • Acetone: olive oil is 4:1 volume to volume ratio
- 132 • Positive control is a 25% acetone: olive oil (4:1) solution of hexyl cinnamic
133 aldehyde (CASRN: 101-86-0)

134 **5.0 Procedures on Test Days 1, 2, 3 and 7**

135 **5.1 Day 1**

136 Mark the animals on the tail with their test group number and a number from one to four.
137 Weigh the animals and record their weight to the nearest 0.1 g on the test forms.
138 Remove the test materials from the refrigerator. Should the materials arrive with instructions
139 to heat or sonicate the treatments prior to application, perform these procedures as instructed.

140 *5.1.1 Pre-treatment with 1% SLS Solution*

141 Beginning with Group 1 and proceeding in order to Group 6, the SLS solution should be
142 applied with a brush to the dorsum of both ears of the mice. The number of the SLS solution
143 used should match the test group number. The brush should be dipped in the SLS solution
144 and applied to the dorsum of one ear using a petting motion, covering the entire dorsum with
145 four to five strokes. Dip the brush again in the SLS solution and apply the solution to the
146 dorsum of the other ear in the same manner.

147 Record the time when beginning to apply SLS solution to Group 1 and when completing
148 application to Group 6. The application procedure should be performed continuously without
149 delay for Groups 1 through 6.

150 Six brushes should be prepared and numbered, using only one brush for each test group.
151 When performing the same application procedure on Days 2, 3, and 7 there is the possibility
152 of brush contamination due to residual solution on the mouse auricula. It is important to
153 switch brushes after finishing application for one group and check the number of the next
154 brush before proceeding to the next group. After use, the brushes should be washed
155 thoroughly and made available for the next day.

156 *5.1.2 Test Substance Application*

157 One hour after starting the SLS solution application, the numbered treatments should be
158 applied to the auriculae of the mice, beginning with Group 1 and ending with Group 6. Using
159 a micropipette or similar device, 25 µL of the test solution should be dripped slowly on the
160 dorsum of one of the mouse's ears, covering the dorsum entirely. Again take up 25µL of
161 treatment solution and apply it in the same manner to the dorsum of the mouse's other ear.

162 When applying the treatments, micropipette tips should be changed for each test group. After
163 completing application for one test group, remove the tip and spray the end of the
164 micropipette with an alcohol mist and wipe to avoid contamination.

165 Record the time when beginning to apply the test solution to Group 1 and when completing
166 application to Group 6. The application procedure should be performed continuously without
167 delay for Groups 1 through 6.

168 Immediately after completing application the test materials should be refrigerated.

169 5.1.3 *General Information on the 1% SLS Pre-treatment and Test Substance Application*

170 The objective of the application procedure is to first apply SLS solution to the entirety of the
171 dorsum of the ear and then to apply a prescribed amount of test solution to the same area.

172 Using ether anesthesia ensures ease and accuracy of the procedure. However, special care
173 should be taken to avoid taking the life of the animals in the course of anesthesia. If one
174 technician immobilizes the animal and extends the ear with tweezers while the other
175 technician applies the solution, the procedure can be performed with accuracy without using
176 anesthesia. If this approach is used six pairs of tweezers should be prepared, one for each
177 group, to avoid contamination. Alternatively, the tweezers should be wiped with an alcohol
178 swab after application is completed for each test group.

179 5.2 **Days 2 and 3**

180 Apply SLS solution and treatments using the same procedures as for Day 1.

181 When performing the application procedures the animals should be observed carefully for
182 necrosis, hardening, hyperplasia or erythema of the auricula, as well as piloerection, or a
183 decrease in locomotor activity. Any such abnormalities observed should be recorded on the
184 test forms.

185 5.3 **Day 7**

186 On Day 7 the same procedures should be performed as on Days 2 and 3.

187 Excision of the auricular lymph nodes will be performed from 24 to 30 hours after the start of
188 application on Day 7. It is therefore recommended that application procedures on Day 7
189 begin in the morning or early afternoon.

190 6.0 **Procedure on Test Day 8 (Excision of Auricular Lymph Nodes and 191 ATP Assay)**

192 6.1 **Laboratory Preparation**

193 Forty-eight 15 mL test tubes should each be filled with 1.98 mL of phosphate buffered saline
194 (PBS). The dispensing of PBS should be conducted under aseptic manipulation. Dispense a
195 minimum of 24 mL of PBS in a 50 mL test tube. Pipetting should be under aseptic
196 manipulation.

197 Dissolve the luciferin-luciferase reagent according to the ATP assay kit instructions (at least
198 4.8 mL are required). The ATP assay kit provided, CheckLite™ 250 Plus, includes five
199 bottles each of Luciferin-luciferase reagent, solvent water, and ATP releasing agent. Using
200 one bottle of each type, create a solution according to the instructions (approximately 5.5
201 mL). Shield the assay solutions from light using aluminum foil and refrigerate until the time
202 of use. Immediately before using, return to room temperature and remove the foil prior to
203 use. Dispense 0.1 mL of the ATP releasing agent included in the ATP assay kit to each of the
204 48 luminometer tubes. ATP assay kit reagents should be dispensed using sterilized pipette
205 tips under aseptic manipulation to avoid contamination with ATP and microorganisms.

206 **6.2 Body Weight Measurement**

207 Weigh the mice and record their body weights to the nearest 0.1 g on the test forms.

208 **6.3 Auricular Lymph Node Excision and Weight Measurement**

209 Perform procedures in **Sections 6.3, 6.4** and **6.5** within 24 to 30 hours after the start of
210 treatment application on Day 7. The necessary materials for procedures in **Sections 6.3, 6.4**
211 and **6.5** are given in **Annex I**.

212 Immediately after sacrificing the mice with ether anesthesia excise completely all auricular
213 lymph nodes for each ear (there can be one or two auricular lymph nodes) as illustrated in
214 **Figure A-1**. Place the excised lymph nodes for one animal in a disposable petri dish and
215 immediately measure the wet weight to the nearest 0.1 mg with an analytical balance.

216 **6.4 Preparation of Cell Suspension**

217 The lymph nodes from one animal should be sandwiched between two pieces of slide glass
218 and light pressure should be applied to crush the nodes (**Figure A-2**). After confirming that
219 the tissue has spread out thinly pull the two slides apart. Suspend the tissue on both pieces of
220 slide glass in 1 mL of PBS. As illustrated in **Figure A-3**, each piece of slide glass should be
221 held at an angle over the petri dish and rinsed with PBS while the tissue is scraped off of the



222 glass with repeated movements of a cell scraper. One mL of PBS should be used for rinsing
223 both slides.

224 The tissue suspension in the petri dish should be homogenized lightly with the cell scraper,
225 and 20 μ L of the suspension should be taken up with a micropipette, taking care not to take
226 up the membrane that is visible to the eye. The pipetted suspension should be added to 1.98
227 mL of PBS and homogenized well. This will be cell suspension No. 1. Again take up 20 μ L
228 of the suspension in the petri dish, add to 1.98 mL of PBS, and homogenize well. This will be
229 cell suspension No. 2.

230 These procedures should be performed while wearing gloves and a mask and micropipette
231 tips should be sterile. Detailed step-by-step procedures are given in **Annex II**.

232

232 **Figure A-1 Auricular lymph nodes³**




233

234

³ Taken from ICCVAM IWG LLNA Protocol (ICCVAM 2001)
A-13

234 **Figure A-2 Preparation of cell suspension**


235 Lymph nodes from each animal are sandwiched between two pieces of slide glass and light
236 pressure is applied to crush the nodes.



237

238 **Figure A-3 Preparation of cell suspension**

239 Rinse with PBS while scraping the tissue off of the glass with a cell scraper. Repeat the
240 scraping motion, scooping up liquid from the petri dish as need. Use 1 mL of PBS for the
241 nodes of each animal.



242

243

243 **6.5 ATP Assay**

244 Prepare 48 luminometer tubes in advance by dispensing 0.1 mL of the ATP releasing reagent
245 provided to each tube. Add 0.1 mL of each homogenized cell suspension to the luminometer
246 tubes and homogenize. After allowing the solution in the tube to stand for approximately 20
247 seconds, add 0.1 mL of the Luciferin-luciferase solution, promptly homogenize and place in
248 the luminometer. The amount of bioluminescence (RLU; relative luminescence units)
249 measured over 10 seconds will be displayed. Record this measurement on the test forms.

250 The amount of bioluminescence begins to decrease immediately after adding the Luciferin-
251 luciferase solution. It is therefore important that the series of procedures from the addition of
252 Luciferin-luciferase solution to switching on the luminometer are performed as quickly as
253 possible, ideally with the same rhythm.

254 These procedures should be performed while wearing gloves and a mask and micropipette
255 tips should be sterile. The detailed procedures are given in **Annex III**.

256 **7.0 Points of Caution on Procedures from Excision to ATP Assay**

257 The ATP content of the lymph node decreases over time after the sacrifice of the animal. It is
258 therefore desirable that the time elapsed between sacrifice of the animal and ATP assay is
259 uniform for each animal. The series of procedures from excision to ATP assay must be
260 performed rapidly and without delay.

261 If these procedures are performed by one technician, the animals should be sacrificed one at a
262 time. If there are multiple technicians, it is possible to divide tasks and sacrifice the animals
263 one group at a time. If two technicians perform the procedures, steps in **Section 6.3** should be
264 performed by one individual and steps in **Sections 6.4** and **6.5** should be performed by the
265 other. If three technicians perform the procedures, steps in **Sections 6.3**, **6.4** and **6.5** can each
266 be handled by one individual. If multiple technicians are involved, it is important that the
267 timing of excision is carefully planned so that there are no delays in subsequent steps.

268 **8.0 Data Entry**

269 Input the body weights on Day 1 and Day 8, the lymph node weight, and the amount of ATP
270 bioluminescence into the designated Excel file.

271 **Annex I: Equipment and Reagents Used for the Experimental Procedures**
272 **in Sections 6.3, 6.4, and 6.5**

273 For the equipment and reagents underlined below, the items provided by the LLNA: DA
274 Validation Committee should be used. In the event the test facility provides a luminometer, it
275 can be used.

276 **6.3 Auricular Lymph Node Excision and Weight Measurement**

277 Dissecting instruments set (Tweezers, scissors, surgical holder, injection needle and holder)

278 Antiseptic solution

279 Cotton

280 Petri dish (24)

281 Analytical balance (readability of at least 0.1 mg)

282 **6.4 Preparation of Cell Suspension**

283 15 mL test tubes with 1.98 mL phosphate buffered saline (PBS) (48)

284 50 mL test tubes with at least 24 mL PBS (1)

285 Slide glass (48)

286 Tweezers (1)

287 Micropipette 1000 μ L (1) (Volume to be measured: 1 mL)

288 Micropipette 100 μ L (1) (Volume to be measured: 20 μ L)

289 Cell scraper (1)

290 Sterilized pipette tips for 1000 μ L micropipette (24) and for 100 μ L micropipette (24)

291 Vortex mixer (1)

292 Paper towels

293 Clean sheet

294 Test tube rack

295 **6.5 ATP Assay**

296 Luminometer tubes with 0.1 mL ATP releasing agent (48)

297 15 mL test tube with dissolved luciferin-luciferase solution (1)

298 Micropipette - 100 μ L or 200 μ L (2) (Volume to be measured: 0.1 mL)

299 Sterilized micropipette tips (96)

300 Timer (with second display) (1)

301 Luminometer (1)

302 Vortex mixer (can use same mixer listed under **Section 6.4** Preparation of Cell Suspension)

303 Test tube rack

304 Luminometer tube rack (microtube rack)

305 **Annex II: Preparation of Cell Suspension for the Experimental Procedures**
306 **in Section 6.4**

- 307 1. Cover the laboratory bench with a clean sheet and place one piece of slide glass on the
308 sheet.
- 309 2. After measuring the lymph node weights, use tweezers to move the lymph nodes from
310 one animal from the petri dish to the center of the slide glass.
- 311 3. Place another piece of slide glass on top.
- 312 4. Pick up the two sandwiched pieces of slide glass. Squeeze the two pieces in the center to
313 crush the lymph nodes. (Apply only light pressure. Too much pressure can break the
314 cells.)
- 315 5. Confirm that the tissue has spread out thinly between the two slides and place the
316 sandwiched slides on the clean sheet.
- 317 6. Fasten a tip on the 1000 μ L micropipette and draw 1 mL phosphate buffered saline (PBS)
318 from the 50 mL tube.
- 319 7. Remove the upper slide glass from the sandwiched slides and place it on the clean sheet
320 with the side that was in contact with the lymph node tissue facing up. The other slide
321 glass should be held at an angle in the petri dish, the side with lymph node tissue affixed
322 facing forward, and washed with 1 mL PBS.
- 323 8. Dispose of the 1000 μ L micropipette tip.
- 324 9. Scrape the tissue off of the glass with a cell scraper, scooping up PBS from the petri dish
325 and repeating the scraping motion. Confirm that there is no tissue, or only trace amounts
326 of tissue, left on the slide before disposing of the slide glass.
- 327 10. Pick up the slide glass laid aside at step 7; scrape the tissue off in the same manner and
328 dispose of the slide glass. Note that it becomes difficult to scrape the tissue off of the
329 slide glass once it has dried. Perform steps 4 through 10 without delay. The scraping
330 should be performed while keeping the area of the slide glass to which the lymph node
331 tissue is affixed sufficiently wet with PBS from the petri dish.

- 332 11. The tissue suspension in the petri dish should be homogenized lightly with the cell
333 scraper. If large pieces of tissue are observed, stir with the cell scraper to break up the
334 pieces and obtain a uniform solution.
- 335 12. Wipe the cell scraper with a paper towel. (The cell scraper will be used for the next
336 animal.)
- 337 13. Fasten a tip to the 100 μ L micropipette, tilt the petri dish at an angle and mix the
338 suspension by pipetting in and out several times. Take up 20 μ L of the suspension with
339 the pipette, taking care not to take up any membrane that is visible to the eye.
- 340 14. Add the 20 μ L of suspension to a 15 mL test tube containing 1.98 mL PBS. Pipette the
341 solution and proceed to homogenize with the vortex mixer. (cell suspension No. 1)
- 342 15. Repeat steps 13 and 14 to prepare cell suspension No. 2.
- 343 16. Dispose of the 100 μ L micropipette tip.
- 344

345 Annex III: ATP Assay for the Experimental Procedures in Section 6.5

- 346 1. Fasten a tip on the 100 μ L (or 200 μ L) micropipette and draw 0.1 mL of vortex-
347 homogenized cell suspension No. 1.
- 348 2. To the luminometer tube filled with 0.1 mL ATP releasing reagent, add 0.1 mL of cell
349 suspension No. 1, making sure to note the time with a timer. Dispose of the tip.
- 350 3. Homogenize with the vortex mixer and place in the luminometer tube rack.
- 351 4. Fasten a tip on a separate 100 μ L (or 200 μ L) micropipette and draw 0.1 mL of solution
352 from the 15 mL tube containing dissolved Luciferin-luciferase reagent.
- 353 5. Take the luminometer tube from the rack and add 0.1 mL of Luciferin-luciferase solution
354 to the luminometer tube 20 seconds after the time noted in step 2.
- 355 6. Promptly homogenize in the vortex mixer, place in the luminometer and turn on the
356 switch. The amount of bioluminescence begins to decrease immediately after adding the
357 Luciferin-luciferase solution. Step 6 should be performed as quickly as possible, ideally
358 with the same rhythm.
- 359 7. Dispose of the tip.
- 360 8. After 10 seconds the amount of bioluminescence (RLU; relative luminescence units) will
361 be displayed. Record this measurement on the test forms.
- 362 9. Repeat steps 1 through 8 for cell suspension No. 2, measure the bioluminescence and
363 record.

364

365

366

367

368

369

370

371

372

373

374

375

Appendix A2

376

Results in the LLNA: DA Test Method for 1% Sodium Lauryl Sulfate (SLS)

377

Pretreatment versus without 1% SLS Pretreatment

378

379

380

381

382

383

384

385

386

387

388

389

390

391

[This Page Intentionally Left Blank]

392

392 **Appendix A2 Summary of Results in the LLNA: DA Test Method with 1% SLS**
 393 **Pretreatment versus without 1% SLS Pretreatment**

Substance Name	Vehicle	Concentration (%)	SI ¹ (+ SLS)	SI ¹ (- SLS)	Calculated EC3 ² (%) (+ SLS)	Calculated EC3 ² (%) (- SLS)
2, 4-Dinitrochlorobenzene	AOO	0.03	2.10	1.88	0.05	0.06
		0.10	5.02	4.46		
		0.30	9.74	14.61		
Potassium dichromate	DMSO	0.1	2.61	2.54	0.15	0.22
		0.3	4.24	3.34		
		1.0	5.51	5.66		
Isoeugenol	AOO	1.0	2.05	1.32	2.46	4.24
		2.5	3.02	2.21		
		5.0	2.85	3.35		
Citral	AOO	5	1.93	1.88	7.4	10.4
		10	4.15	2.91		
		25	6.97	5.90		
Hexyl cinnamic aldehyde	AOO	5	1.51	0.99	7.5	8.8
		10	4.52	3.64		
		25	4.84	3.79		
Cinnamic alcohol	AOO	10	2.46	2.44	14.1	18.5
		25	4.40	3.43		
		50	6.36	4.01		
Hydroxycitronellal	AOO	10	1.98	1.49	15.8	19.8
		25	4.61	3.81		
		50	6.59	6.74		
Imidazolidinyl urea	DMF	10	2.36	2.54	20.3	33.0
		25	3.29	2.38		
		50	6.02	4.31		
Methyl methacrylate	AOO	25	0.73	1.11	NA	NA
		50	0.68	0.92		

Substance Name	Vehicle	Concentration (%)	SI ¹ (+ SLS)	SI ¹ (- SLS)	Calculated EC3 ² (%) (+ SLS)	Calculated EC3 ² (%) (- SLS)
		100	1.31	1.83		
Nickel (II) chloride	DMSO	2.5	1.53	0.98	NA	NA
		5.0	1.57	1.16		
		10.0	2.24	1.87		
Methyl salicylate	AOO	5	0.89	0.83	NA	NA
		10	1.59	1.32		
		25	1.69	2.34		
Salicylic acid	AOO	5	1.21	1.13	NA	NA
		10	2.05	1.29		
		25	2.48	2.44		
Sulfanilamide	DMF	10	1.08	0.92	NA	NA
		25	1.03	0.90		
		50	0.94	0.84		

394 Abbreviations: AOO = acetone: olive oil (4:1); DMF = *N,N*-dimethylformamide; DMSO = dimethyl sulfoxide;
395 EC3 = estimated concentration required to produce a stimulation index of three; NA = not applicable; RLU =
396 relative luminescence units; SI = stimulation index; SLS = sodium lauryl sulfate.

397 ¹SI determined from mean ATP content (RLU).

398 ²EC3 value was calculated based on interpolation or extrapolation formulas discussed in Gerberick et al. 2004.

399 + SLS = with pretreatment of 1% SLS prior to test substance application

400 - SLS = without pretreatment of 1% SLS prior to test substance application

401

402