

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15

Appendix E
Accuracy Analyses Using Additional Approaches for Combining
Multiple Test Results

15
16
17
18
19
20
21
22
23
24
25
26
27
28
29

[This Page Intentionally Left Blank]

30 **1.0 LLNA: BrdU-ELISA Accuracy Analysis Using Alternative**
31 **Decision Criteria and Alternate Methods for Combining Data for**
32 **Substances Tested Multiple Times**

33 This appendix shows performance analyses for the LLNA: BrdU-ELISA using single
34 alternative decision criteria and two different approaches for combining test results for the 14
35 substances with multiple LLNA: BrdU-ELISA tests:

- 36 1. The positive/negative outcome for each substance for each criterion was
37 determined by the outcome of the test with the highest maximum SI of the
38 multiple tests.
- 39 2. The positive/negative outcome for each substance for each criterion was
40 determined by the outcome of the test with the lowest maximum SI of the
41 multiple tests.

42 **Section 6.4** provides the results for the analysis when the most prevalent outcome was used
43 as the result for each substance tested multiple times (for each criterion).

44 **1.1 Results of LLNA: BrdU-ELISA Accuracy Analysis Using Single Alternative**
45 **Decision Criteria and the Highest Maximum SI for the Outcome of Multiple**
46 **Tests**

47 When combining multiple test results for a single substance using the outcome of the test
48 with the highest maximum SI, the decision criterion of $SI \geq 2.0$ to identify sensitizers yielded
49 an accuracy of 77% (24/31), a sensitivity of 77% (17/22), a specificity of 78% (7/9), a false
50 positive rate of 22% (2/9), and a false negative rate of 23% (5/22) (**Table E-1**). $SI \geq 2.0$ was
51 the decision criterion used by the JSAAE interlaboratory validation study of the LLNA:
52 BrdU-ELISA. The performance for the additional decision criteria described in **Section 6.4**,
53 are also shown in **Table E-1**. Over the range of SI cutoffs evaluated, increasing the SI cutoff
54 decreased accuracy (84% at $SI \geq 1.3$ to 68% at $SI \geq 5.0$), decreased sensitivity (100% at
55 $SI \geq 1.3$ to 54% at $SI \geq 5.0$), increased specificity (44% at $SI \geq 1.3$ to 100% at $SI \geq 5.0$),
56 decreased the false positive rate (56% at $SI \geq 1.3$ to 0% at $SI \geq 5.0$), and increased the false
57 negative rate (0% at $SI \geq 1.3$ to 46% at $SI \geq 5.0$) (**Figure E-1** and **Table E-1**). Use of
58 ANOVA and summary statistics (i.e., mean absorbance values of treated groups $\geq 95\%$

59 confidence interval of the control group, or ≥ 2 or ≥ 3 SD from the control group mean),
60 yielded accuracy values of 74 to 81%, with sensitivity values of 86 to 100%, and false
61 negative rates of 0 to 14%. The specificity for these criteria ranged from 33 to 56% and the
62 false positive rates were 44 to 67%.

63 The highest accuracy and lowest false negative rate, for the approach using the highest
64 maximum SI for the substances with more than one test, was achieved using an $SI \geq 1.3$. The
65 accuracy for $SI \geq 1.3$ was 84% (26/31), with sensitivity of 100% (22/22), specificity of 44%
66 (5/9), a false positive rate of 56% (5/9), and false negative rate of 0% (0/22). However, using
67 an $SI \geq 1.3$ incorrectly classified lactic acid, isopropanol, and methyl salicylate, three of the
68 ICCVAM performance standards reference substances, as sensitizers. Use of mean
69 absorbance values of treated groups $\geq 95\%$ confidence interval of the control group to
70 identify sensitizers also produced the lowest false negative rate 0% (0/9), with a slightly
71 lower accuracy of 81% (25/31), and a higher false positive rate of 67% (6/9). This criterion
72 also incorrectly classified lactic acid and isopropanol as sensitizers. The lowest false positive
73 rates (0% [0/9]) were produced by SI cutoffs of 3.0 to 5.0, however the false negative rates at
74 those cutoffs were 27% (6/22) to 46% (10/22).

75 **Table E-1 Performance of the LLNA: BrdU-ELISA Compared with the Traditional LLNA Using Alternative Decision**
 76 **Criteria to Identify Sensitizers and the Highest Maximum SI for Substances with Multiple Tests**

Alternate Criterion	N ¹	Accuracy		Sensitivity		Specificity		False Positive Rate		False Negative Rate		Positive Predictivity		Negative Predictivity	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
Statistics ³	31	74	23/31	86	19/22	44	4/9	56	5/9	14	3/22	79	19/24	57	4/7
≥ 95% CI ⁴	31	81	25/31	100	22/22	33	3/9	67	6/9	0	0/22	79	22/28	100	3/3
≥ 2 SD ⁵	31	77	24/31	96	21/22	33	3/9	67	6/9	5	1/22	78	21/27	75	3/4
≥ 3 SD ⁶	31	81	25/31	91	20/22	56	5/9	44	4/9	9	2/22	83	20/24	71	5/7
SI ≥ 5.0	31	68	21/31	54	12/22	100	9/9	0	0/9	46	10/22	100	12/12	47	9/19
SI ≥ 4.5	31	71	22/31	59	13/22	100	9/9	0	0/9	41	9/22	100	13/13	50	9/18
SI ≥ 4.0	31	77	24/31	68	15/22	100	9/9	0	0/9	32	7/22	100	15/15	56	9/16
SI ≥ 3.5	31	77	24/31	68	15/22	100	9/9	0	0/9	32	7/22	100	15/15	56	9/16
SI ≥ 3.0	31	81	25/31	73	16/22	100	9/9	0	0/9	27	6/22	100	16/16	60	9/15
SI ≥ 2.5	31	81	25/31	77	17/22	89	8/9	11	1/9	23	5/22	94	17/18	62	8/13
SI ≥ 2.0	31	77	24/31	77	17/22	78	7/9	22	2/9	23	5/22	90	17/19	58	7/12
SI ≥ 1.5	31	81	25/31	91	20/22	56	5/9	44	4/9	9	2/22	83	20/24	71	5/7
SI ≥ 1.3	31	84	26/31	100	22/22	44	4/9	56	5/9	0	0/22	82	22/27	100	4/4

77 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA)
 78 detection of bromodeoxyuridine (BrdU); CI = confidence interval; No. = number; SD = standard deviation; SI = stimulation index

79 ¹ N = Number of substances included in this analysis.

80 ² The proportion on which the percentage calculation is based.


81 ³ Analysis of variance for difference of group means when substances were tested at multiple doses or *t*-test when substances were tested at one dose. The absorbance data
 82 were log-transformed prior to analysis of variance. Significance at *p* < 0.05 was further tested by Dunnett's test.

83 ⁴ The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.

84 ⁵ The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.

85 ⁶ The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group.

86 **Figure E-1 Performance of the LLNA: BrdU-ELISA with SI Compared to the**
87 **Traditional LLNA Using the Highest Maximum SI for Substances with**
88 **Multiple Tests**



89
90 As compared to traditional LLNA results, the lines show the change in performance characteristics
91 for the LLNA: BrdU-ELISA with the SI cutoff used to identify sensitizers. This analysis used LLNA:
92 BrdU-ELISA and traditional LLNA results for 31 substances (22 sensitizers and nine nonsensitizers
93 based on traditional LLNA results). For the 14 substances with multiple test results, the results for
94 each substance were combined by using the outcome for the test with the highest maximum SI value.
95 The solid line shows accuracy, the dashed line shows the false positive rate, and the dotted line shows
96 the false negative rate.

97
98

99 **1.2 Results of LLNA: BrdU-ELISA Accuracy Analysis Using Alternative**
100 **Decision Criteria and Lowest Maximum SI for the Outcome of Multiple Tests**

101 When combining multiple test results for a single substance using the outcome of the test
102 with the lowest maximum SI, the decision criterion of $SI \geq 2.0$ to identify sensitizers for
103 these 31 substances yielded an accuracy of 81% (25/31), a sensitivity of 73% (16/22), a
104 specificity of 100% (9/9), a false positive rate of 0% (0/9), and a false negative rate of 27%
105 (6/22) (**Table E-2**). $SI \geq 2.0$ was the decision criterion used by the JSAAE interlaboratory
106 validation study of the LLNA: BrdU-ELISA. The performance for the additional decision
107 criteria described in **Section 6.4**, are shown in **Table E-2**.

108 Over the range of SI cutoffs evaluated, increasing the SI cutoff decreased accuracy (90% at
109 $SI \geq 1.3$ to 45% at $SI \geq 5.0$), decreased sensitivity (100% at $SI \geq 1.3$ to 23% at $SI \geq 5.0$),
110 increased specificity (67% at $SI \geq 1.3$ to 100% at $SI \geq 5.0$), decreased the false positive rate
111 (33% at $SI \geq 1.3$ to 0% at $SI \geq 5.0$), and increased the false negative rate (0% at $SI \geq 1.3$ to
112 77% at $SI \geq 5.0$) (**Figure E-2** and **Table E-2**). Use of ANOVA and summary statistics (i.e.,
113 mean absorbance values of treated groups $\geq 95\%$ confidence interval of the control group, or
114 ≥ 2 or 3 SD from the control group mean), yielded accuracy of 84 to 90%, with sensitivity
115 values of 86 to 100%, and false negative rates of 0 to 14%. The specificity for these criteria
116 ranged from 67 to 89% and the false positive rates were 11 to 33%.

117 The highest accuracy and lowest false negative rate, for the approach using the lowest
118 maximum SI for the substances with more than one test, was achieved using an $SI \geq 1.3$ and
119 mean absorbance values of treated groups $\geq 95\%$ confidence interval of the control group.
120 Both criteria yielded an accuracy of 90% (26/31) and a false negative rate of 0% (0/22). Both
121 criteria also yielded sensitivity = 100% (22/22), specificity = 67% (6/9), and false positive
122 rate = 33% (3/9). However, these criteria incorrectly classified lactic acid, isopropanol, and
123 methyl salicylate, three of the ICCVAM performance standards reference substances, as
124 sensitizers. The lowest false positive rate (0% [0/9]) was produced by SI cutoffs of 2.0 to 5.0,
125 however the false negative rates at those cutoffs were 27% to 77%. Of those cutoffs, $SI \geq 2.0$
126 produced the highest accuracy, 81% (25/31), and the lowest false negative rate, 27% (6/22).

127 **Table E-2 Performance of the LLNA: BrdU-ELISA Compared with the Traditional LLNA Using Alternative Decision**
 128 **Criteria to Identify Sensitizers and the Lowest Maximum SI for Substances with Multiple Tests**

Alternate Criterion	N ¹	Accuracy		Sensitivity		Specificity		False Positive Rate		False Negative Rate		Positive Predictivity		Negative Predictivity	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
Statistics ³	31	84	26/31	86	19/22	78	7/9	22	2/9	14	3/22	91	19/21	70	7/10
≥ 95% CI ⁴	31	90	28/31	100	22/22	67	6/9	33	3/9	0	0/22	88	22/25	100	6/6
≥ 2 SD ⁵	31	87	27/31	96	21/22	67	6/9	33	3/9	5	1/22	88	21/24	86	6/7
≥ 3 SD ⁶	31	87	27/31	86	19/22	89	8/9	11	1/9	14	3/22	95	19/20	73	8/11
SI ≥ 5.0	31	45	14/31	23	5/22	100	9/9	0	0/9	77	17/22	100	5/5	35	9/26
SI ≥ 4.5	31	45	14/31	23	5/22	100	9/9	0	0/9	77	17/22	100	5/5	35	9/26
SI ≥ 4.0	31	55	17/31	36	8/22	100	9/9	0	0/9	64	14/22	100	8/8	39	9/23
SI ≥ 3.5	31	55	17/31	36	8/22	100	9/9	0	0/9	64	14/22	100	8/8	39	9/23
SI ≥ 3.0	31	64	20/31	50	11/22	100	9/9	0	0/9	50	11/22	100	11/11	45	9/20
SI ≥ 2.5	31	71	22/31	59	13/22	100	9/9	0	0/9	41	9/22	100	13/13	50	9/18
SI ≥ 2.0	31	81	25/31	73	16/22	100	9/9	0	0/9	27	6/22	100	16/16	60	9/15
SI ≥ 1.5	31	87	27/31	91	20/22	78	7/9	22	2/9	9	2/22	91	20/22	78	7/9
SI ≥ 1.3	31	90	28/31	100	22/22	67	6/9	33	3/9	0	0/22	88	22/25	100	6/6

129 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA)
 130 detection of bromodeoxyuridine (BrdU); CI = confidence interval; No. = number; SD = standard deviation; SI = stimulation index

131 ¹N = Number of substances included in this analysis.

132 ²The proportion on which the percentage calculation is based.


133 ³Analysis of variance for difference of group means when substances were tested at multiple doses or t-test when substances were tested at one dose. The absorbance data
 134 were log-transformed prior to analysis of variance. Significance at p < 0.05 was further tested by Dunnett's test.

135 ⁴The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.

136 ⁵The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.

137 ⁶The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group.

138 **Figure E-2 Performance of the LLNA: BrdU-ELISA with SI Compared to the**
 139 **Traditional LLNA Using the Lowest Maximum SI for Substances with**
 140 **Multiple Tests**



141
 142 As compared to traditional LLNA results, the lines show the change in performance characteristics for the
 143 LLNA: BrdU-ELISA with the SI cutoff used to identify sensitizers. This analysis used LLNA: BrdU-
 144 ELISA and traditional LLNA results for 31 substances (22 sensitizers and nine nonsensitizers based on
 145 traditional LLNA results). For the 14 substances with multiple test results, the results for each substance
 146 were combined by using the outcome for the test with the lowest maximum SI value. The solid line shows
 147 accuracy, the dashed line shows the false positive rate, and the dotted line shows the false negative rate.
 148

148 **2.0 Discordant Results For Accuracy Analysis of Alternative Decision** 149 **Criteria**

150 Using the decision criteria of $SI \geq 2.0$ to identify sensitizers and the most prevalent outcome for
151 the substances with multiple tests for the analysis of alternative decision criteria, the five
152 discordant substances (when compared to the traditional LLNA) were aniline, cyclamen
153 aldehyde, hydroxycitronellal, 2-mercaptobenzothiazole, and linalool (**Table 6-4**). As indicated in
154 **Section 6.4.1**, all five substances were misclassified as nonsensitizers when compared to the
155 traditional LLNA, which classified them as sensitizers.

156 **2.1 Discordant Results Using Alternative Decision Criteria and Highest Maximum SI** 157 **Outcome for Multiple Tests**

158 Using the decision criteria of $SI \geq 2.0$ to identify sensitizers and the test with the highest
159 maximum SI as the result for substances with multiple tests, yielded two additional discordant
160 substances: isopropanol and lactic acid, which were misclassified as sensitizers.

161 **Table E-3** shows how the number and identity of discordant substances changes with the
162 alternate decision criteria when using the test with the highest maximum SI as the result for
163 substances with multiple tests. Using an SI cutoff less than 2.0, $SI \geq 1.5$, to identify sensitizers
164 yielded six discordant substances. Two substances, hydroxycitronellal and linalool, were
165 misclassified as nonsensitizers, and four substances, hexane, isopropanol, lactic acid and
166 propylene glycol, were misclassified as sensitizers. Using an even lower SI to identify
167 sensitizers, $SI \geq 1.3$, yielded five discordant substances that were all misclassified as sensitizers
168 (hexane, isopropanol, lactic acid, methyl salicylate, and propylene glycol). Increasing the SI
169 cutoff to values greater than 2.0, increased the number of sensitizers that were misclassified as
170 nonsensitizers. At $SI \geq 2.0$, five sensitizers were misclassified as nonsensitizers while at $SI \geq 5.0$,
171 10 sensitizers were classified as nonsensitizers (**Table E-3**). At $SI \geq 2.0$, two nonsensitizers were
172 misclassified as sensitizers while, at $SI \geq 5.0$, no nonsensitizers were classified as sensitizers.

173 Use of a statistical test (i.e., ANOVA or *t*-test) or summary statistics (i.e., $\geq 95\%$ CI, ≥ 2 SD, or ≥ 3
174 SD) tended to misclassify more nonsensitizers than sensitizers. Using ANOVA or a *t*-test to
175 identify sensitizers misclassified three sensitizers (linalool, 2-mercaptobenzothiazole, and

176 aniline) as nonsensitizers and four nonsensitizers (glycerol, hexane, isopropanol, and lactic acid)
177 as sensitizers. Using treatment group absorbance $\geq 95\%$ CI or ≥ 2 SD of vehicle control mean
178 misclassified six nonsensitizers as sensitizers (glycerol, hexane, isopropanol, lactic acid, methyl
179 salicylate, and propylene glycol). Treatment group absorbance ≥ 2 SD of vehicle control mean
180 also misclassified one weak sensitizer as a nonsensitizer (linalool). Using treatment group
181 absorbance ≥ 3 SD of vehicle control mean misclassified three nonsensitizers as sensitizers
182 (hexane, isopropanol, and lactic acid) and three weak sensitizers as nonsensitizers
183 (hydroxycitronellal, linalool, and aniline).

184 Four ICCVAM performance standards reference substances were discordant for the analysis of
185 alternate decision criteria using the test with the highest maximum SI as the result for substances
186 with multiple tests (**Table E-3**). One sensitizer, 2-mercaptobenzothiazole, was misclassified by
187 some criteria as a nonsensitizer, and three nonsensitizers, isopropanol, lactic acid, and methyl
188 salicylate, were misclassified as sensitizers by some criteria. The criteria that yielded the correct
189 results for 2-mercaptobenzothiazole included summary statistics (i.e., $\geq 95\%$ CI, ≥ 2 SD, or
190 ≥ 3 SD), $SI \geq 1.5$, and $SI \geq 1.3$. The criteria that yielded the correct results for isopropanol
191 included $SI \geq 2.5$ to 5.0. The criteria that yielded the correct results for lactic acid included $SI \geq$
192 3.0 to 5.0. All criteria yielded the correct results for methyl salicylate except for treatment group
193 absorbance $\geq 95\%$ CI or ≥ 2 SD of vehicle control mean, and $SI \geq 1.3$.

194 **Table E-3 Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision Criteria Compared to the**
 195 **Traditional LLNA and the Highest Maximum SI for Substances with Multiple Tests**

Discordant Substance ¹	Alternate Decision Criterion ²												
	Statistics ³	≥ 95% CI ⁴	≥ 2 SD ⁵	≥ 3 SD ⁶	SI ≥ 5.0	SI ≥ 4.5	SI ≥ 4.0	SI ≥ 3.5	SI ≥ 3.0	SI ≥ 2.5	SI ≥ 2.0	SI ≥ 1.5	SI ≥ 1.3
2-Mercaptobenzothiazole (1.7%)	-				-	-	-	-	-	-	-		
Cinnamic aldehyde (2.4%)					-	-							
3-Aminophenol (3.2%)					-	-	-	-					
Nickel sulfate (4.8%)					-								
4-Chloroaniline (6.5%)					-	-	-	-	-				
Cyclamen aldehyde (22.3%)					-	-	-	-	-	-	-		
Hydroxycitronellal (24%)				-	-	-	-	-	-	-	-	-	
Linalool (30%)	-		-	-	-	-	-	-	-	-	-	-	
Isopropyl myristate (44%)					-	-							
Aniline (63%)	-			-	-	-	-	-	-	-	-		
Glycerol (-)	+	+	+										
Hexane (-)	+	+	+	+								+	+
Isopropanol (-)	+	+	+	+							+	+	+
Lactic acid (-)	+	+	+	+						+	+	+	+
Methyl salicylate (-)		+	+										+
Propylene glycol (-)		+	+									+	+

196 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA)
 197 detection of bromodeoxyuridine (BrdU); CI = confidence interval; SD = standard deviation; SI = stimulation index.

198 ¹Compared to the traditional LLNA. Traditional LLNA result in parentheses: “-” for nonsensitizers and EC3 (%) for sensitizers.

199 ²LLNA: BrdU result shown: “+” if the decision criterion was met and “-” if the decision criterion was not met.

200 ³Analysis of variance for difference of group means when substances were tested at multiple doses or *t*-test when substances were tested at one dose. The absorbance data
 201 were log-transformed prior to analysis of variance. Significance at *p* < 0.05 was further tested by Dunnett’s test.

202 ⁴The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.

203 ⁵The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.

204 ⁶The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group.

205 **2.2 Discordant Results Using Alternative Decision Criteria and Lowest Maximum SI**
206 **Outcome for Multiple Tests**

207 Using the decision criteria of $SI \geq 2.0$ to identify sensitizers and the most prevalent outcome for
208 the substances with multiple tests for the analysis of alternative decision criteria yielded five
209 discordant substances (when compared to the traditional LLNA). Aniline, cyclamen aldehyde,
210 hydroxycitronellal, 2-mercaptobenzothiazole, and linalool were misclassified as nonsensitizers
211 (**Table 6-4**), while the traditional LLNA classified them as sensitizers. Using the test with the
212 lowest maximum SI as the result for substances with multiple tests yielded six discordant
213 substances at $SI \geq 2.0$. One additional sensitizer, formaldehyde, was misclassified as a
214 nonsensitizer (**Table E-4**).

215 **Table E-4** shows how the number and identity of discordant substances changes with the
216 alternate decision criteria when using the test with the lowest maximum SI as the result for
217 substances with multiple tests. Using an SI cutoff less than 2.0, $SI \geq 1.5$, to identify sensitizers
218 yielded four discordant substances. Two substances, hydroxycitronellal and linalool, were
219 misclassified as nonsensitizers, and two substances, hexane and lactic acid, were misclassified as
220 sensitizers. Using an even lower SI to identify sensitizers, $SI \geq 1.3$, yielded three discordant
221 substances that were all misclassified as sensitizers (hexane, lactic acid, and methyl salicylate).
222 Increasing the SI cutoff to values greater than 2.0, increased the number of sensitizers that were
223 misclassified as nonsensitizers. At $SI \geq 2.0$, six sensitizers were misclassified as nonsensitizers
224 while at $SI \geq 4.5$ and $SI \geq 5.0$, 17 sensitizers were classified as nonsensitizers (**Table E-4**). From
225 $SI \geq 2.0$ to $SI \geq 5.0$ no nonsensitizers were misclassified as sensitizers.

226 Using the test with the lowest maximum SI as the result for substances with multiple tests caused
227 even potent sensitizers to be misclassified as nonsensitizers at the higher SI cutoffs. At $SI \geq 4.5$
228 and $SI \geq 5.0$, 2,4-dinitrochlorobenzene, glutaraldehyde, and formaldehyde were classified as
229 nonsensitizers. Glutaraldehyde was classified as a nonsensitizer at SI cutoffs as low as 2.5 and
230 formaldehyde was classified as a nonsensitizer at SI cutoffs as low as 2.0.

231 **Table E-4** **Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision Criteria Compared to the**
 232 **Traditional LLNA and the Lowest Maximum SI for Substances with Multiple Tests**

Discordant Substance ¹	Alternate Decision Criterion ²												
	Statistics ³	≥ 95% CI ⁴	≥ 2 SD ⁵	≥ 3 SD ⁶	SI ≥ 5.0	SI ≥ 4.5	SI ≥ 4.0	SI ≥ 3.5	SI ≥ 3.0	SI ≥ 2.5	SI ≥ 2.0	SI ≥ 1.5	SI ≥ 1.3
2,4-Dinitrochlorobenzene (0.049%)					-	-							
Glutaraldehyde (0.14%)					-	-	-	-	-	-			
Formaldehyde (0.53%)					-	-	-	-	-	-	-		
<i>trans</i> -Cinnamic aldehyde (1.4%)					-	-	-	-					
Isoeugenol (1.5%)					-	-	-	-	-	-			
2-Mercaptobenzothiazole (1.7%)	-				-	-	-	-	-	-	-		
Cinnamic aldehyde (2.4%)					-	-							
3-Aminophenol (3.2%)					-	-	-	-					
Nickel sulfate (4.8%)					-	-	-	-	-				
4-Chloroaniline (6.5%)					-	-	-	-	-				
Hexyl cinnamic aldehyde (9.7%)					-	-	-	-	-	-			
Eugenol (10.1%)					-	-	-	-					
Cyclamen aldehyde (22.3%)					-	-	-	-	-	-	-		
Hydroxycitronellal (24%)				-	-	-	-	-	-	-	-	-	
Linalool (30%)	-		-	-	-	-	-	-	-	-	-	-	
Isopropyl myristate (44%)					-	-							
Aniline (63%)	-			-	-	-	-	-	-	-	-		
Glycerol (-)	+	+	+										
Hexane (-)	+	+										+	+
Lactic acid (-)		+	+									+	+
Methyl salicylate (-)													+

233 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay (LLNA) with enzyme-linked immunosorbent assay (ELISA)
 234 detection of bromodeoxyuridine (BrdU); CI = confidence interval; SD = standard deviation; SI = stimulation index.

235 ¹Compared to the traditional LLNA. Traditional LLNA result in parentheses: “-” for nonsensitizers and EC3 (%) for sensitizers.

- 236 ²LLNA: BrdU result shown: “+” if the decision criterion was met and “-” if the decision criterion was not met.
- 237 ³Analysis of variance for difference of group means when substances were tested at multiple doses or t-test when substances were tested at one dose. The absorbance data
- 238 were log-transformed prior to analysis of variance. Significance at $p < 0.05$ was further tested by Dunnett’s test.
- 239 ⁴The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.
- 240 ⁵The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.
- 241 ⁶The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group.
- 242

243 Use of a statistical test (i.e., ANOVA or *t*-test) or summary statistics (i.e., $\geq 95\%$ CI, ≥ 2
244 SD, or $3 \geq \text{SD}$) more often misclassified nonsensitizers than sensitizers (**Table E-4**).
245 Using ANOVA or *t*-tests to identify sensitizers misclassified three sensitizers (2-
246 mercaptobenzothiazole, linalool, and aniline) as nonsensitizers and two nonsensitizers
247 (glycerol and hexane) as sensitizers. Using treatment group absorbance $\geq 95\%$ CI or ≥ 2
248 SD of vehicle control mean misclassified glycerol and lactic acid as sensitizers. Using
249 treatment group absorbance $\geq 95\%$ CI also misclassified hexane as a sensitizer.
250 Additionally, treatment group absorbance ≥ 2 SD of vehicle control mean misclassified
251 one weak sensitizer (linalool) as a nonsensitizer. Using treatment group absorbance ≥ 3
252 SD of vehicle control mean, however, misclassified three weak sensitizers as
253 nonsensitizers (hydroxycitronellal, linalool, and aniline).

254 Seven ICCVAM performance standards reference substances were discordant for the
255 analysis of alternate decision criteria using the test with the lowest maximum SI as the
256 result for substances with multiple tests (**Table E-4**). One strong sensitizer, 2,4-
257 dinitrochlorobenzene, was misclassified by some criteria as a nonsensitizer. Four
258 additional sensitizers, isoeugenol, 2-mercaptobenzothiazole, hexyl cinnamic aldehyde,
259 and eugenol, were also misclassified as nonsensitizers by some criteria. The criteria that
260 yielded the correct results for 2,4-dinitrochlorobenzene were all but the $\text{SI} \geq 4.5$ to 5.0
261 criteria. The criteria that yielded the correct results for 2-mercaptobenzothiazole included
262 summary statistics (i.e., $\geq 95\%$ CI, ≥ 2 SD, or ≥ 3 SD), $\text{SI} \geq 1.5$ and $\text{SI} \geq 1.3$. The criteria
263 that yielded the correct results for isoeugenol and hexyl cinnamic aldehyde were $\text{SI} \geq 1.5$
264 to 2.0 . The criteria that yielded the correct results for eugenol were $\text{SI} \geq 1.3$ to 3.0 .

265 Two nonsensitizers, lactic acid and methyl salicylate, from the list of ICCVAM
266 performance standards reference substances, were misclassified as sensitizers by some
267 criteria. The criteria that yielded the correct results for lactic acid were all except for
268 treatment group absorbance $\geq 95\%$ CI or ≥ 2 SD of vehicle control mean, $\text{SI} \geq 1.5$, and SI
269 ≥ 1.3 . The criteria that yielded the correct results for methyl salicylate were all except for
270 $\text{SI} \geq 1.3$.