

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15

**Non-radioactive Murine Local Lymph Node Assay: BrdU-ELISA
Test Method Protocol
(LLNA: BrdU-ELISA)
Revised Draft Background Review Document**

March 2009

15
16
17
18
19
20
21
22
23
24
25
26
27
28
29

[This Page Intentionally Left Blank]

30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53

Table of Contents

Page Number

List of Tables iv

List of Figures vi

List of Abbreviations and Acronyms vii

Interagency Coordinating Committee on the Validation Of Alternative Methods
(ICCVAM) Designated Agency Representatives ix

Acknowledgements x

Preface xiii

Executive Summary xv

1.0 Introduction 1

 1.1 Public Heath Perspective 1

 1.2 Historical Background for the Murine Local Lymph Node Assay 1

 1.3 The LLNA: BrdU-ELISA 3

2.0 LLNA: BrdU-ELISA Test Method Protocol 4

 2.1 Decision Criteria 4

3.0 LLNA: BrdU-ELISA Validation Database 6

4.0 Reference Data 9

5.0 Test Method Data and Results 10

6.0 Test Method Accuracy 11

 6.1 LLNA: BrdU-ELISA Database Used for the Accuracy Analysis 11

 6.2 Accuracy Analysis Using the $SI \geq 2.0$ Decision Criterion 13

 6.3 Accuracy Analysis ($SI \geq 2.0$) Based on the ICCVAM Performance
Standards Reference Substances 15

54 6.4 Discordant Results for Accuracy Analysis Using the $SI \geq 2.0$ Decision
 55 Criterion 18

56 6.5 LLNA: BrdU-ELISA Accuracy Analysis Using One Alternative
 57 Decision Criterion..... 22

58 6.6 Discordant Results for Accuracy Analysis of One Alternative Decision
 59 Criterion 30

60 6.7 LLNA: BrdU-ELISA Accuracy Analysis Using Multiple Alternative
 61 Decision Criteria..... 37

62 6.8 Discordant Results for Accuracy Analysis of Multiple Alternative
 63 Decision Criteria..... 38

64 7.0 Test Method Reliability 42

65 7.1 Intralaboratory Reproducibility 42

66 7.2 Interlaboratory Reproducibility 48

67 7.3 Reproducibility for the LLNA: BrdU-ELISA Using Multiple
 68 Alternative Decision Criteria..... 53

69 8.0 Data Quality 57

70 9.0 Other Scientific Reports and Reviews..... 58

71 10.0 Animal Welfare Considerations..... 59

72 10.1 Rationale for the Need to Use Animals 59

73 10.2 Basis for Determining the Number of Animals Used..... 59

74 10.3 Reduction Considerations 59

75 11.0 Practical Considerations 60

76 11.1 Transferability of the LLNA: BrdU-ELISA..... 60

77 11.2 Facilities and Major Fixed Equipment Required to Conduct the
 78 LLNA: BrdU-ELISA 60

79 11.3 LLNA: BrdU-ELISA Training Considerations 60

80 12.0 References..... 61

81 Appendix A LLNA: BrdU-ELISA Protocol..... A-1

82 Appendix B Physico-Chemical Properties of Substances Tested Using the

83 LLNA: BrdU-ELISA B-1

84 Appendix C Comparative LLNA: BrdU-ELISA, Traditional LLNA,

85 Guinea Pig, and Human Skin Sensitization Data..... C-1

86 Appendix D Individual Animal Data for the LLNA: BrdU-ELISA D-1

87 Appendix E Accuracy Analyses Using Additional Approaches for Combining

88 Multiple Test Results E-1

89 Appendix F Reproducibility Analyses for LLNA: BrdU-ELISA with Decision

90 Criterion of $SI \geq 1.5$ F-1

91

List of Tables

	Page Number
92	
93	
94	Table 3-1 Traditional LLNA EC3 Values and Chemical Classification of Substances
95	Tested in the LLNA: BrdU-ELISA 7
96	Table 6-1 Performance of the LLNA: BrdU-ELISA in Predicting Skin Sensitizing
97	Potential Using Decision Criteria of $SI \geq 2.0$ to Identify Sensitizers 14
98	Table 6-2 Performance of the LLNA: BrdU-ELISA ($SI \geq 2.0$) Using the
99	ICCVAM Performance Standards Reference Substances..... 16
100	Table 6-3 Characteristics of the Substances Tested in the LLNA: BrdU-ELISA vs.
101	the ICCVAM Performance Standards Reference Substances 18
102	Table 6-4 Discordant Results for LLNA: BrdU-ELISA (Using $SI \geq 2.0$ for Sensitizers)
103	Compared to Traditional LLNA and Guinea Pig Reference Data..... 20
104	Table 6-5 Discordant Results for LLNA: BrdU-ELISA ($SI \geq 2.0$) when Compared
105	to Traditional LLNA and Human Outcome Data..... 22
106	Table 6-6 Performance of the LLNA: BrdU-ELISA in Predicting Skin Sensitizing
107	Potential Using Alternative Decision Criteria to Identify Sensitizers
108	and the Most Prevalent Outcome for Substances with Multiple Tests 25
109	Table 6-7 Comparison of Performance for Decision Criteria of $SI \geq 1.5$ (Bold), $> 2SD$
110	(Bold Italics), and $SI \geq 2.0$ for Predicting Skin Sensitizing
111	Potential with LLNA: BrdU-ELISA..... 29
112	Table 6-8 Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision
113	Criteria Compared to the Traditional LLNA and the Most Prevalent
114	Outcome for Substances with Multiple Tests..... 31
115	Table 6-9 Discordant Results for LLNA: BrdU-ELISA (Using $SI \geq 1.5$ for
116	Sensitizers) Compared to the Traditional LLNA and Guinea Pig
117	Reference Data 34
118	Table 6-10 Discordant Results for LLNA: BrdU-ELISA ($SI \geq 1.5$) When
119	Compared to Traditional LLNA and Human Outcome Data 35

120

List of Tables (Continued)

121

Page Number

122

Table 6-11 Discordant Results for LLNA: BrdU-ELISA When Multiple Decision

123

Criteria Are Used 39

124

Table 7-1 Intralaboratory Reproducibility for the LLNA: BrdU-ELISA Outcome

125

of Substances Tested Multiple Times 44

126

Table 7-2 Intralaboratory Reproducibility for the SI of Tested Substances in the

127

LLNA: BrdU-ELISA - Coefficient of Variation 46

128

Table 7-3 Intralaboratory Reproducibility for the EC2 of Tested Substances in the

129

LLNA: BrdU-ELISA - Coefficient of Variation 47

130

Table 7-4 Intralaboratory Reproducibility for the EC3 of Tested Substances in the

131

Traditional LLNA 48

132

Table 7-5 Substances and Test Allocation for the Phase II Interlaboratory

133

Validation Study of the LLNA: BrdU-ELISA..... 49

134

Table 7-6 Qualitative Results for the Phase II Interlaboratory Validation Study

135

of the LLNA: BrdU-ELISA..... 50

136

Table 7-7 EC2 Values for the Phase II Interlaboratory Validation Study

137

of the LLNA: BrdU-ELISA..... 52

138

Table 7-8 Interlaboratory Reproducibility of the EC3 for Substances Tested in the

139

Traditional LLNA 53

140

Table 7-9 Frequency of Maximum SI for LLNA: BrdU-ELISA Tests

141

by Category and Traditional LLNA Outcome..... 55

142

Table 7-10 Concordance of LLNA: BrdU-ELISA Tests for Substances

143

with Multiple Tests by Maximum SI Category 56

144

145

		List of Figures	
145			
146			Page Number
147	Figure 6-1	Performance of the LLNA: BrdU-ELISA with SI Compared to the	
148		Traditional LLNA Using the Most Prevalent Outcome for Substances	
149		with Multiple Tests	27
150			

151

List of Abbreviations and Acronyms

152	ACD	Allergic contact dermatitis
153	ANOVA	Analysis of variance
154	AOO	Acetone: olive oil
155	BRD	Background review document
156	BrdU	Bromodeoxyuridine
157	CI	Confidence interval
158	CASRN	Chemical Abstracts Service Registry Number
159	Conc.	Concentration tested
160	CPSC	U.S. Consumer Product Safety Commission
161	CV	Coefficient of variation
162	DMF	<i>N,N</i> -dimethylformamide
163	DMSO	Dimethyl sulfoxide
164	DNA	Deoxyribonucleic acid
165	EC1.5	Estimated concentration needed to produce a stimulation index of 1.5
166		
167	EC2	Estimated concentration needed to produce a stimulation index of two
168		
169	EC3	Estimated concentration needed to produce a stimulation index of three
170		
171	ECt	Estimated concentration needed to produce a stimulation index equaling or greater than a specified threshold
172		
173	ELISA	Enzyme-linked immunosorbent assay
174	EPA	U.S. Environmental Protection Agency
175	GPMT	Guinea pig maximization test
176	HCA	Hexyl cinnamic aldehyde
177	HMT	Human Maximization Test
178	HPTA	Human patch test allergen
179	ICCVAM	Interagency Coordinating Committee on the Validation of Alternative Methods
180		
181	ISO	International Organization for Standardization
182	IWG	Immunotoxicity Working Group
183	JSAAE	Japanese Society for Alternatives to Animal Experiments
184	K _{ow}	Octanol-water partition coefficient
185	LLNA	Local Lymph Node Assay
186	LLNA: BrdU-ELISA	LLNA with enzyme-linked immunosorbent assay detection of bromodeoxyuridine
187		
188	MEK	Methyl ethyl ketone
189	MeSH	Medical Subject Headings
190	Min	Minimal
191	Mod	Moderate
192	MW	Molecular weight
193	NA	Not available
194	NC	Not calculated
195	NK	Not known

196	NICEATM	National Toxicology Program Interagency Center for the
197		Evaluation of Alternative Toxicological Methods
198	NT	Not tested
199	NTP	National Toxicology Program
200	OECD	Organisation for Economic Co-operation and Development
201	Res	Result
202	SD	Standard Deviation
203	SI	Stimulation Index
204	TG	Test Guideline
205	U.S.	United States
206	Unk	Unknown
207	Veh.	Vehicle
208	vs.	Versus
209	w/v	Weight to volume ratio

210
211

Interagency Coordinating Committee on the Validation of Alternative Methods: Agency Representatives

Agency for Toxic Substances And Disease Registry

- Moiz Mumtaz, Ph.D.

Consumer Product Safety Commission

- Marilyn L. Wind, Ph.D. (Chair)
- ◇ Kristina Hatlelid, Ph.D.
- Joanna Matheson, Ph.D.

Department of Agriculture

- Jodie Kulpa-Eddy, D.V.M. (Vice-Chair)
- ◇ Elizabeth Goldentyer, D.V.M.

Department of Defense

- Robert E. Foster, Ph.D.
- ◇ Patty Decot
- Peter J. Schultheiss, D.V.M., D.A.C.L.A.M.
- Harry Salem., Ph.D.

Department of Energy

- Michael Kuperberg, Ph.D.
- ◇ Marvin Stodolsky, Ph.D.

Department of the Interior

- Barnett A. Rattner, Ph.D.
- ◇ Sarah Gerould, Ph.D.

Department of Transportation

- George Cushmac, Ph.D.
- ◇ Steve Hwang, Ph.D.

Environmental Protection Agency

Office of Science Coordination and Policy

- Karen Hamernik, Ph.D.

Office of Research and Development

- ◇ Julian Preston, Ph.D.
- Suzanne McMaster, Ph.D.

Office of Pesticide Programs

- Amy Rispin, Ph.D.
- Deborah McCall

OECD Test Guidelines Program

- Jerry Smrcek, Ph.D.

- Principal agency representative
- ◇ Alternate principal agency representative

Food and Drug Administration

Office of Science

- Suzanne Fitzpatrick, Ph.D., D.A.B.T.

Center for Drug Evaluation and Research

- ◇ Abigail C. Jacobs, Ph.D.
- Paul C. Brown, Ph.D.

Center for Devices and Radiological Health

- Melvin E. Stratmeyer, Ph.D.
- Vasant G. Malshet, Ph.D., D.A.B.T.

Center for Biologics Evaluation and Research

- Richard McFarland, Ph.D., M.D.
- Ying Huang, Ph.D.

Center for Food Safety and Nutrition

- David G. Hattan, Ph.D.
- Robert L. Bronaugh, Ph.D.

Center for Veterinary Medicine

- Devaraya Jagannath, Ph.D.
- M. Cecilia Aguila, D.V.M.

National Center for Toxicological Research

- William T. Allaben, Ph.D.
- Paul Howard, Ph.D.

Office of Regulatory Affairs

- Lawrence D'Hoostelaere, Ph.D.

National Cancer Institute

- Alan Poland, M.D.
- ◇ T. Kevin Howcroft, Ph.D.

National Institute of Environmental Health Sciences

- William S. Stokes, D.V.M., D.A.C.L.A.M
- ◇ Raymond R. Tice, Ph.D.
- Rajendra S. Chhabra, Ph.D., D.A.B.T.
- Jerrold J. Heindel, Ph.D.

National Institute for Occupational Safety and Health

- Paul Nicolaysen, V.M.D.
- ◇ K. Murali Rao, M.D., Ph.D.

National Institutes of Health

- Margaret D. Snyder, Ph.D.

National Library of Medicine

- Pertti (Bert) Hakkinen, Ph.D.
- ◇ Jeanne Goshorn, M.S.

Occupational Safety and Health Administration

- Surender Ahir, Ph.D.

297

Acknowledgements

298

299

300

**Interagency Coordinating Committee on the Validation of
Alternative Methods (ICCVAM) Immunotoxicity Working Group (IWG)****301 U.S. Consumer Product Safety
302 Commission**

303 Joanna Matheson, Ph.D. (IWG Co-Chair)

304 Marilyn Wind, Ph.D.

305 U.S. Environmental Protection Agency306 *Office of Pesticide Programs*

307 Masih Hashim, D.V.M., Ph.D.

308 Marianne Lewis

309 Deborah McCall

310 Timothy McMahon, Ph.D.

311 Amy Rispin, Ph.D.

312 *Office of Prevention, Pesticides, and Toxic
313 Substances*

314 Ronald Ward, Ph.D.

315 *Office of Research and Development*

316 Marsha Ward, Ph.D.

317 *Office of Science Coordination and Policy*

318 Karen Hamernik, Ph.D.

319 U.S. Food and Drug Administration320 *Center for Devices and Radiological Health*

321 Daniel Lyle, Ph.D.

322 Vasant G. Malshet, Ph.D., D.A.B.T.

323 Jeffrey Toy, Ph.D.

324 *Center for Drug Evaluation and Research*

325 Paul Brown, Ph.D.

326 Abigail Jacobs, Ph.D. (IWG Co-Chair)

327 Jiaqin Yao, Ph.D.

328 *Center for Veterinary Medicine*

329 Ruth Barratt, Ph.D., D.V.M.

330

**330 National Institute of Environmental
331 Health Sciences**

332 Dori Germolec, Ph.D.

333 William S. Stokes, D.V.M., D.A.C.L.A.M.

334 Raymond R. Tice, Ph.D.

**335 National Institute for Occupational
336 Safety and Health**

337 B. Jean Meade, D.V.M., Ph.D.

338 National Library of Medicine

339 Pertti (Bert) Hakkinen, Ph.D.

**340 European Centre for the Validation of
341 Alternative Methods — Liaison**

342 Silvia Casati, Ph.D.

**343 Japanese Center for the Validation of
344 Alternative Methods — Liaison**

345 Hajime Kojima, Ph.D.

346 **National Toxicology Program Interagency Center For The Evaluation Of**
347 **Alternative Toxicological Methods**
348

349 **National Institute of Environmental Health Sciences**

350 William Stokes, D.V.M., D.A.C.L.A.M.
351 Director; Project Officer

352 Raymond Tice, Ph.D.

353 Deborah McCarley
354 Special Assistant; Assistant Project Officer

355

356 **NICEATM Support Contract Staff (Integrated Laboratory Systems, Inc.)**

347 David Allen, Ph.D.

348 Thomas Burns, M.S.

349 Linda Litchfield

350 Gregory Moyer, M.B.A.

351 Michael Paris

352 Eleni Salicru, Ph.D.

353 Catherine Sprankle

354 Frank Stack

355 Judy Strickland, Ph.D., D.A.B.T.

356 **Statistical Consultant for ILS, Inc.**

357 Joseph Haseman, Ph.D.

358

359

Other Acknowledgements

360 *ICCVAM and NICEATM gratefully acknowledge Masahiro Takeyoshi, Ph.D., of the*
361 *Chemicals Evaluation and Research Institute (CERI-Japan), Saitama, Japan, and Hajime*
362 *Kojima, Ph.D., Japanese Center for the Validation of Alternative Methods, Tokyo, Japan, for*
363 *submitting data to NICEATM used for the evaluation of the LLNA: BrdU-ELISA.*
364

365

Preface

366 In 1999, the U.S. Interagency Coordinating Committee on the Validation of Alternative
367 Methods (ICCVAM) recommended the murine (mouse) local lymph node assay (LLNA) as a
368 valid test method to assess the skin sensitization potential of most types of substances
369 (ICCVAM 1999). ICCVAM concluded that the LLNA (referred to herein as the “traditional
370 LLNA”) provided several advantages compared to the guinea pig method, including
371 elimination of potential pain and distress, use of fewer animals, less time required to perform,
372 and availability of dose-response information. United States and international regulatory
373 authorities subsequently accepted the traditional LLNA as an alternative test method for
374 allergic contact dermatitis testing. It is now commonly used around the world.

375 One disadvantage of the traditional LLNA is that it requires injection of a radioactive marker
376 to measure cell proliferation in lymph nodes. To avoid the use of radioactive markers,
377 scientists have recently developed several non-radioactive versions of the LLNA. In 2007,
378 the U.S. Consumer Product Safety Commission (CPSC) asked ICCVAM and the National
379 Toxicology Program Interagency Center for the Evaluation of Alternative Methods
380 (NICEATM) to evaluate the scientific validity of these non-radioactive versions. ICCVAM
381 assigned the nomination a high priority, and established the ICCVAM Immunotoxicity
382 Working Group (IWG) to work with NICEATM to review the current literature and evaluate
383 available data to assess the validity of three such test methods. A comprehensive draft
384 background review document (BRD) provided the information, data, and analyses supporting
385 the validation status of each of the non-radioactive test methods. ICCVAM also developed
386 draft test method recommendations for each test method regarding its usefulness and
387 limitations, test method protocol, performance standards, and future studies.

388 NICEATM and ICCVAM provided the draft BRDs and draft recommendations to an
389 international independent scientific peer review panel for their consideration at a public
390 meeting on March 4-6, 2008. A report of the Panel meeting was subsequently published on
391 the NICEATM-ICCVAM website¹. Both the Panel and ICCVAM concluded that more
392 information was needed before a recommendation on the usefulness and limitations of each
393 of the three test methods could be made. The Panel recommended that NICEATM obtain

¹ http://iccvam.niehs.nih.gov/methods/immunotox/llna_PeerPanel08.htm

394 additional existing data that were not available to the Panel and reanalyze the performance of
395 each non-radioactive LLNA method. NICEATM subsequently obtained additional data and
396 prepared updated BRDs. ICCVAM also prepared revised draft test method recommendations
397 based on the revised BRDs. This revised draft BRD addresses the validation database for the
398 LLNA: BrdU-ELISA.

399 The Panel will meet to consider the revised BRDs and to evaluate the extent to which the
400 available information supports the revised ICCVAM draft test method recommendations.
401 ICCVAM will consider the conclusions and recommendations of the Panel, along with
402 comments received from the public and the Scientific Advisory Committee for Alternative
403 Toxicological Methods, and then finalize the BRDs and test method recommendations. These
404 will then be forwarded to Federal agencies for their consideration and acceptance decisions
405 where appropriate.

406 We gratefully acknowledge the organizations and scientists who provided data and
407 information for this document. We also acknowledge the efforts of those individuals
408 contributing to the preparation of this BRD, including the following staff from the
409 NICEATM Support Contractor, Integrated Laboratory Systems, Inc.: David Allen, Ph.D.,
410 Thomas Burns, M.S., Gregory Moyer, M.B.A., Michael Paris, Eleni Salicru, Ph.D., Catherine
411 Sprankle, Frank Stack, and Judy Strickland, Ph.D. We also thank the members of the
412 ICCVAM IWG, chaired by Abigail Jacobs, Ph.D. (U.S. Food and Drug Administration) and
413 Joanna Matheson, Ph.D. (CPSC), and ICCVAM representatives who subsequently reviewed
414 and provided comments throughout the process leading to this final draft version.

415 Marilyn Wind, Ph.D.
416 Deputy Associate Executive Director
417 Directorate for Health Sciences
418 U.S. Consumer Product Safety Commission
419 Chair, ICCVAM

420
421 RADM William S. Stokes, D.V.M., D.A.C.L.A.M.
422 Assistant Surgeon General, U.S. Public Health Service
423 Director, NICEATM
424 Executive Director, ICCVAM

425
426 March 2009

427 **Executive Summary**

428 ***Background***

429 In 1999, the Interagency Coordinating Committee on the Validation of Alternative Methods
430 (ICCVAM) recommended to U.S. Federal agencies that the murine local lymph node assay
431 (LLNA) is a valid substitute for currently accepted guinea pig test methods to assess the
432 allergic contact dermatitis (ACD) potential of many, but not all, types of substances. ACD is
433 an allergic skin reaction characterized by redness, swelling, and itching that can result from
434 contact with a sensitizing chemical or product. The recommendation was based on a
435 comprehensive evaluation that included an independent scientific peer review panel (Panel)
436 assessment of the validation status of the LLNA. The Panel report and the ICCVAM
437 recommendations (ICCVAM 1999) are available at the National Toxicology Program
438 Interagency Center for the Evaluation of Alternative Toxicological Methods (NICEATM)-
439 ICCVAM website (http://iccvam.niehs.nih.gov/docs/immunotox_docs/llna/llnarep.pdf). The
440 LLNA was subsequently incorporated into national and international test guidelines for the
441 assessment of skin sensitization (Organisation for Economic Co-operation and Development
442 [OECD] Test Guideline 429 [OECD 2002]; International Organization for Standardization
443 [ISO] 10993-10: Tests for Irritation and Sensitization [ISO 2002]; U.S. Environmental
444 Protection Agency [EPA] Health Effect Testing Guidelines on Skin Sensitization [EPA
445 2003]).

446 In 2007, the U.S. Consumer Product Safety Commission (CPSC) formally nominated several
447 activities related to the LLNA for evaluation by ICCVAM and NICEATM (Available at
448 http://iccvam.niehs.nih.gov/methods/immunotox/llnadoocs/CPSC_LLNA_nom.pdf). One of
449 the nominated activities was an assessment of the validation status of non-radioactive
450 alternatives to the current version of the LLNA ([ICCVAM 1999; Dean et al. 2001] referred
451 to hereafter as the “traditional LLNA”), which uses radioactivity to detect sensitizers. The
452 information described in the original and this revised background review document (BRD)
453 was compiled by ICCVAM and NICEATM in response to this nomination. The BRD
454 provides a comprehensive review of available data and information regarding the usefulness
455 and limitations of one of these methods, the LLNA with detection of bromodeoxyuridine

456 (BrdU) incorporation by an enzyme-linked immunosorbent assay (ELISA) (referred to
457 hereafter as the “LLNA: BrdU-ELISA”).

458 ***Revisions to the LLNA: BrdU-ELISA Evaluation***

459 NICEATM and ICCVAM convened an independent scientific peer review panel meeting on
460 March 4-6, 2008. The Panel peer reviewed the draft BRD and commented on the extent that
461 it supported the draft ICCVAM test method recommendations on the usefulness and
462 limitations of the LLNA: BrdU-ELISA. Both ICCVAM and the Panel concluded that more
463 information was needed before a recommendation on the usefulness and limitations of the
464 LLNA: BrdU-ELISA could be made². The Panel indicated that the following information
465 was needed: a detailed protocol, individual animal data, and an evaluation of interlaboratory
466 reproducibility. The Panel recommended that additional data be obtained by NICEATM and
467 that a reanalysis of the performance of the LLNA: BrdU-ELISA be conducted. In response to
468 this recommendation, NICEATM obtained additional LLNA: BrdU-ELISA data from the test
469 sponsor, which were used to update the evaluation. These data include:

- 470 • LLNA: BrdU-ELISA data for six substances not previously provided to
471 NICEATM. (Note: The number of substances evaluated effectively increased
472 by seven with the location of reference data for one substance for which
473 LLNA: BrdU-ELISA data had been previously submitted). These data were
474 used in a reanalysis of test method accuracy, which is detailed in **Section 6.0**
475 of this BRD.
- 476 • Individual animal data for the LLNA: BrdU-ELISA studies included in the
477 interlaboratory validation study of 10 substances. These data were used in
478 additional quantitative analyses of test method reproducibility, which are
479 detailed in **Section 7.0** of this BRD.

480 ***Test Method Protocol***

481 The protocol in this draft BRD has been revised from the January 2008 draft BRD to include
482 the decision criterion of $SI \geq 2.0$, rather than $SI \geq 3.0$, to identify substances as sensitizers.
483 The LLNA: BrdU-ELISA was originally developed by Takeyoshi et al. (2001). While the

² http://iccvam.niehs.nih.gov/methods/immunotox/llna_PeerPanel08.htm

484 traditional LLNA assesses cellular proliferation by measuring the incorporation of
485 radioactivity into the deoxyribonucleic acid (DNA) of dividing lymph node cells, the LLNA:
486 BrdU-ELISA assesses the same endpoint by measuring the incorporation of the thymidine
487 analog BrdU using an ELISA. A stimulation index (SI), the ratio of the mean BrdU
488 incorporation into the lymph nodes of mice in the test substance group to the mean BrdU
489 incorporation into the lymph nodes of mice in the vehicle control group is used to identify a
490 substance as a sensitizer. Other than the procedure for measuring lymph node cell
491 proliferation, the protocol for the LLNA: BrdU-ELISA is similar to that of the traditional
492 LLNA (Dean et al. 2001; ICCVAM 1999).

493 ***Validation Database***

494 The validation database in this draft BRD has been revised from the January 2008 draft BRD
495 to include seven additional substances (six substances for which LLNA: BrdU-ELISA data
496 were not previously obtained and on previously included substance for which traditional
497 LLNA data were recently obtained). The accuracy and reliability of the LLNA: BrdU-ELISA
498 were assessed using the individual animal data for 31 substances from six published studies
499 (Takeyoshi et al. 2003; 2004a; 2004b; 2005; 2006; 2007a), one platform presentation
500 (Takeyoshi 2007b), and one poster presentation (Kojima et al. 2008). The reference test data
501 for these substances were obtained from the traditional LLNA, guinea pig (GP) skin
502 sensitization tests, and/or human skin sensitization tests or clinical information. Of the 31
503 substances with traditional LLNA data, 22 were classified by the traditional LLNA as skin
504 sensitizers and nine were classified as nonsensitizers.

505 ***Test Method Accuracy***

506 The accuracy evaluation in this draft BRD has been revised from the January 2008 draft
507 BRD to include the results for seven additional substances. Other revisions included the
508 evaluation of multiple decision criteria, including the $SI \geq 2.0$ recommended in the test
509 method protocol, and the evaluation of two different criteria used to classify sensitizers and
510 nonsensitizers. Based on the evaluation of multiple decision criteria, the optimal performance
511 was achieved using $SI \geq 2.0$ to classify sensitizers and $SI < 1.3$ to classify nonsensitizers.
512 When these two criteria are used, false positive results (0/9) and false negative results (0/22)
513 are eliminated compared with the traditional LLNA. However, using these criteria, 11

514 substances have an $SI \geq 1.3$ to < 2.0 , 6/11 substances were sensitizers and 5/11 substances
515 were nonsensitizers when tested in the traditional LLNA. Other available information, such
516 as peptide reactivity, could be used to interpret LLNA: BrdU-ELISA results when $1.3 \leq SI <$
517 2.0 . Sixty-seven percent (4/6) of the sensitizers in this range had peptide reactivity data and
518 all four had low to moderate peptide reactivity. All (5/5) of the nonsensitizers had minimal
519 peptide reactivity.

520 When a single decision criterion of $SI \geq 2.0$ was used to classify sensitizers vs.
521 nonsensitizers, compared to the traditional LLNA, accuracy was 84% (26/31), with a false
522 positive rate of 0% (0/9), and the false negative rate of 23% (5/22). Among the false negative
523 substances, no unique characteristics were identified that could be used as rationale for
524 excluding any particular types of substances from testing in the LLNA: BrdU-ELISA.

525 ***Test Method Reliability – Intralaboratory Reproducibility***

526 The intralaboratory reproducibility evaluation in this draft BRD has been revised from the
527 January 2008 draft BRD to include the results for a number of additional tests for which SI
528 values were newly available. Intralaboratory reproducibility was assessed using a
529 concordance analysis of sensitizer/nonsensitizer results, and a coefficient of variation (CV)
530 analysis of SI values and EC₂ values (estimated concentration needed to produce an SI of 2).
531 The qualitative analysis shows that multiple tests of eight substances (six sensitizers and two
532 nonsensitizers) yielded 100% concordance for sensitizer/nonsensitizer outcomes for seven of
533 the eight substances. In the quantitative analyses, the CVs for the SI values of nine
534 substance/concentration combinations that were tested up to five times each ranged from 1%
535 to 79%. The CVs for the EC₂ values of three substances that were tested up to five times at
536 multiple doses ranged from 16% to 73%.

537 ***Test Method Reliability – Interlaboratory Reproducibility***

538 The interlaboratory reproducibility evaluation is a new addition to this draft BRD because
539 interlaboratory data were not available for evaluation in the January 2008 draft BRD. This
540 draft BRD also includes a reproducibility analysis using separate SI criteria to identify
541 sensitizers and nonsensitizers. When using $SI \geq 2.0$ to classify sensitizers, the qualitative
542 interlaboratory reproducibility analysis of 10 substances (seven sensitizers and three
543 nonsensitizers), that were tested in up to seven laboratories indicated 100% agreement (3/3,

544 6/6, or 7/7) among the laboratories for seven substances (six sensitizers and one
545 nonsensitizer). There was 67% (2/3 or 4/6) agreement among the tests for the remaining
546 sensitizer and two nonsensitizers. Interlaboratory CV values for the EC2 values of the seven
547 sensitizers ranged from 20 to 101%.

548 When using $SI \geq 2.0$ to classify sensitizers and $SI < 1.3$ to classify nonsensitizers, the
549 concordance analysis for the 14 substances with multiple tests indicated that the SI results for
550 89% (8/9) of the sensitizers were 100% concordant (i.e., all yielded $SI \geq 2.0$). The SI results
551 for 40% (2/5) of the nonsensitizers were 100% (i.e., all yielded $1.3 \leq SI < 2.0$). The
552 concordance of the other three nonsensitizers was 50% (1/2) to 57% (4/7) for $SI < 1.3$ and
553 29% (2/7) to 33% (1/3) for $SI \geq 2.0$.

554 ***Animal Welfare Considerations***

555 The animal welfare considerations in this draft BRD have not changed from the January 2008
556 draft BRD. The LLNA: BrdU-ELISA will use the same number of animals when compared
557 to the updated ICCVAM-recommended LLNA protocol (Appendix A of ICCVAM 2009).
558 However, since use of the traditional LLNA is restricted in some institutions because it
559 involves radioactivity, availability and use of the non-radioactive LLNA: BrdU-ELISA may
560 lead to further reduction in use of the GP tests, which would provide for reduced animal use
561 and increased refinement due to the avoidance of pain and distress in the LLNA procedure.

562 ***Test Method Transferability***

563 The test method transferability considerations in this draft BRD have not changed from the
564 January 2008 draft BRD. The transferability of the LLNA: BrdU-ELISA is expected to be
565 similar to the traditional LLNA. Compared to the traditional LLNA, the LLNA: BrdU-
566 ELISA will not require facilities, equipment, and licensing permits for handling radioactive
567 materials. The level of training and expertise needed to conduct the LLNA: BrdU-ELISA
568 should be similar to the traditional LLNA except that the understanding and use of ELISA is
569 required.

570 ***ICCVAM Revised Draft Recommendations***

571 ICCVAM developed revised draft recommendations for the LLNA: BrdU-ELISA based on
572 the new data and analyses. Recommendations are provided for test method usefulness and

573 limitations, test method protocol, and future studies to further characterize its usefulness and
574 limitations. These are provided in a separate document, *Draft ICCVAM Test Method*
575 *Recommendations, Non-radioactive Murine Local Lymph Node Assay: BrdU-ELISA Test*
576 *Method Protocol (LLNA: BrdU-ELISA)*.

577 **1.0 Introduction**

578 **1.1 Public Health Perspective**

579 Allergic contact dermatitis (ACD) is a frequent occupational health problem. According to
580 the U.S. Department of Labor Bureau of Labor Statistics, in 2005, 980 cases of ACD
581 involved days away from work³.

582 ACD develops in two phases, induction and elicitation. The induction phase occurs when a
583 susceptible individual is exposed topically to a skin-sensitizing substance. Induction depends
584 on the substance passing through the epidermis, where it forms a hapten complex with
585 dermal proteins. The Langerhans cells, the resident antigen-presenting cells in the skin,
586 process the hapten complex. The processed hapten complex then migrates to the draining
587 lymph nodes. Antigen presentation to T-lymphocytes follows, which leads to the clonal
588 expansion of these cells. At this point, the individual is sensitized to the substance (Basketter
589 et al. 2003; Jowsey et al. 2006). Studies have shown that the magnitude of lymphocyte
590 proliferation correlates with the extent to which sensitization develops (Kimber and Dearman
591 1991, 1996).

592 The elicitation phase occurs when the individual is again topically exposed to the same
593 substance. As in the induction phase, the substance penetrates the epidermis, is processed by
594 the Langerhans cells, and presented to circulating T-lymphocytes. The T-lymphocytes are
595 then activated, which causes release of cytokines and other inflammatory mediators. This
596 release produces a rapid dermal immune response that can lead to ACD (ICCVAM 1999;
597 Basketter et al. 2003; Jowsey et al. 2006).

598 **1.2 Historical Background for the Murine Local Lymph Node Assay (LLNA)**

599 In 1999, the Interagency Coordinating Committee on the Validation of Alternative Methods
600 (ICCVAM) recommended that the LLNA is a valid substitute for currently accepted guinea
601 pig (GP) test methods to assess the ACD potential of many, but not all, types of substances.
602 The recommendation was based on a comprehensive evaluation that included an independent
603 scientific peer review panel (Panel) assessment of the validation status of the LLNA. The
604 Panel report and the ICCVAM recommendations (ICCVAM 1999) are available at the

³ Available at <http://www.bls.gov/>

605 National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
606 Toxicological Methods (NICEATM)-ICCVAM website
607 (http://iccvam.niehs.nih.gov/docs/immunotox_docs/llna/llnarep.pdf).

608 ICCVAM forwarded recommendations to U.S. Federal agencies that the LLNA should be
609 considered for regulatory acceptance or other non-regulatory applications for assessing the
610 ACD potential of substances, while recognizing that some testing situations would still
611 require the use of traditional GP test methods (ICCVAM 1999, Sailstad et al. 2001). The
612 LLNA was subsequently incorporated into national and international test guidelines for the
613 assessment of skin sensitization (Organisation for Economic Co-operation and Development
614 [OECD] Test Guideline 429 [OECD 2002]; International Standards Organization [ISO]
615 10993-10: Tests for Irritation and Sensitization [ISO 2002]; U.S. Environmental Protection
616 Agency [EPA] Health Effect Testing Guidelines on Skin Sensitization [EPA 2003]).

617 On January 10, 2007, the U.S. Consumer Product Safety Commission (CPSC) formally
618 nominated several activities related to the LLNA for evaluation by ICCVAM and NICEATM
619 (Available at
620 http://iccvam.niehs.nih.gov/methods/immunotox/llndocs/CPSC_LLNA_nom.pdf). One of
621 the nominated activities was an assessment of the validation status of non-radioactive
622 alternatives to the current version of the LLNA ([ICCVAM 1999, Dean et al. 2001] referred
623 to hereafter as the “traditional LLNA”), which uses radioactivity to detect sensitizers. The
624 information described in this background review document (BRD) was compiled by
625 ICCVAM and NICEATM in response to this nomination. The BRD provides a
626 comprehensive review of available data and information regarding the usefulness and
627 limitations of one of these methods, the LLNA with detection of bromodeoxyuridine (BrdU)
628 incorporation by enzyme-linked immunosorbent assay (ELISA) (referred to hereafter as the
629 “LLNA: BrdU-ELISA”). ICCVAM and its IWG evaluated this method in a draft background
630 review document (BRD) and developed draft test method recommendations based on this
631 evaluation. An independent peer review panel (Panel) reviewed the BRD in March 2008 to
632 evaluate the extent to which the information contained in the BRD supported the draft
633 recommendations. The Panel concluded that additional information was needed to evaluate
634 the method, including a detailed protocol, quantitative data for the method, and an evaluation

635 of interlaboratory reproducibility. After receiving the additional information, this revised
636 draft BRD was compiled for review by the Panel.

637 ICCVAM will consider the conclusions and recommendations of the Panel, along with
638 comments received from the public and the Scientific Advisory Committee for Alternative
639 Toxicological Methods, when developing the final BRD and final recommendations on the
640 usefulness and limitations of each non-radioactive alternative LLNA test method that is being
641 considered.

642 **1.3 The LLNA: BrdU-ELISA**

643 The LLNA: BrdU-ELISA was developed by Takeyoshi et al. (2001) as a non-radioactive
644 alternative to the traditional LLNA. While the traditional LLNA assesses cellular
645 proliferation by measuring the incorporation of radioactivity into the deoxyribonucleic acid
646 (DNA) of dividing lymph node cells, the LLNA: BrdU-ELISA assesses the same endpoint by
647 measuring the incorporation of the thymidine analog BrdU, which is detected and quantified
648 with an ELISA, which is available as a kit commercially from several sources.

649 This document provides:

- 650 • A comprehensive summary of the LLNA: BrdU-ELISA test method protocol
- 651 • The substances used in the validation of the test method and the test results
- 652 • The performance characteristics (accuracy and reliability) of the test method
- 653 • Animal welfare considerations
- 654 • Other considerations relevant to the usefulness and limitations of this test
655 method (e.g., transferability, cost of the test method).

656

656 **2.0 LLNA: BrdU-ELISA Test Method Protocol**

657 The protocol in this draft BRD has been revised from the January 2008 draft BRD to use the
658 decision criterion of $SI \geq 2.0$, rather than $SI \geq 3.0$, to identify substances as sensitizers. The
659 LLNA: BrdU-ELISA protocol (see **Appendix A**) is similar to the ICCVAM-recommended
660 protocol for the traditional LLNA (see Appendix A of ICCVAM [2009]), except for the
661 method used to assess lymphocyte proliferation. In both the LLNA: BrdU-ELISA and the
662 traditional LLNA, the test substance is administered on three consecutive days. In the
663 traditional LLNA, ^3H - thymidine or ^{125}I -iododeoxyuridine (in phosphate buffered saline; 250
664 $\mu\text{L}/\text{mouse}$) is administered via the tail vein two days after the final application of the test
665 substance. In the LLNA: BrdU-ELISA, 5 mg BrdU in a volume of 0.5 mL physiological
666 saline (concentration of 10 mg/mL) is administered via intraperitoneal injection two days
667 after the final application of the test substance. Takeyoshi et al. (2001) reported that one
668 injection of 5 mg BrdU was selected over two injections to minimize the incorporation of
669 BrdU in the control group. Injection of BrdU two days after topical treatment with test
670 substance yielded efficient incorporation of BrdU in comparison to injection one day or three
671 days after topical treatment with a test substance (Takeyoshi et al. 2001). On the day
672 following BrdU injection, lymph nodes are excised and a single cell suspension is prepared
673 from the lymph nodes of each animal. A standard aliquot of the cell suspension is added in
674 triplicate to the wells of a flat-bottom 96-well microplate and centrifuged. Supernatants are
675 then removed. FixDenat solution (Roche Applied Science), which fixes the cells and
676 denatures the DNA in one step, is added to each well, and the plate is incubated at room
677 temperature. The FixDenat solution is removed and the diluted anti-BrdU antibody solution
678 is added to each well. After each well is washed with phosphate-buffered saline, an aliquot of
679 substrate solution containing tetramethylbenzidine is added. After incubation at room
680 temperature, the absorbance is measured using a microplate reader.

681 **2.1. Decision Criteria**

682 Like the traditional LLNA, a stimulation index (SI) is used in the LLNA: BrdU-ELISA to
683 distinguish skin sensitizers from nonsensitizers. The SI is the ratio of the mean absorbance of
684 the incorporated BrdU in a lymph node suspension from individual mice in the test substance

685 group to the mean absorbance of the incorporated BrdU in a lymph node suspension from
686 individual mice in the vehicle control group as indicated by the formula below:

$$687 \quad SI = \frac{\text{Mean absorbance of the treatment group lymph nodes}}{\text{Mean absorbance of the vehicle control group lymph nodes}}$$

688 Consistent with the traditional LLNA, an $SI \geq 3.0$ was initially used as the threshold for
689 labeling a substance as a sensitizer. Takeyoshi et al. (2007b) evaluated the use of other
690 decision criteria such as specific differences in BrdU incorporation between treated and
691 control groups (i.e., greater than the 95% confidence interval [CI] of the control group,
692 greater than the two or three standard deviations [SD] from the control group mean, and
693 statistically significant differences by analysis of variance [ANOVA]) and other SI values to
694 distinguish sensitizers from nonsensitizers and found that lower cutoff values for the SI
695 improved accuracy when compared with the results of the traditional LLNA.

696 A multi-laboratory validation study of the LLNA: BrdU-ELISA organized by the Japanese
697 Society for Alternatives to Animal Experiments (JSAAE) used $SI \geq 2$ to classify sensitizers
698 (Kojima et al. 2008). The $SI \geq 2$ criterion was selected for the interlaboratory validation
699 study because prior studies (Takeyoshi et al. 2003; 2004a; 2004b; 2005; 2006; 2007a; 2007b)
700 indicated that the $SI \geq 3$ criterion was inadequate for reliably distinguishing sensitizers from
701 nonsensitizers (Kojima H, personal communication).

702

702 3.0 LLNA: BrdU-ELISA Validation Database

703 The validation database in this draft BRD has been revised from the January 2008 draft BRD
704 to include seven additional substances. To evaluate the validity of the LLNA: BrdU-ELISA,
705 data were available for 35 substances. Twenty-seven substances were tested in one laboratory
706 (Takeyoshi et al. 2003; 2004a; 2004b; 2005; 2006; 2007a; 2007b; unpublished data) and four
707 additional substances (along with six of the same substances tested by Takeyoshi et al.) were
708 tested in the multi-laboratory validation study coordinated by JSAAE (**Table 3-1**). Most of
709 these substances (31/35) had been previously tested in the traditional LLNA. No traditional
710 LLNA data were available for four substances, which include two dimers of eugenol
711 (dihydroxyl-3,3'-dimethoxy-5,5'-diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-2,3'-
712 dimethoxyphenyl ether) and two dimers of isoeugenol (4-[1-Hydroxy-2-(2-methoxy-4-
713 propenyl-phenyloxy)-propyl]-2-methoxy-phenol and 2-methoxy-4-(7-methoxy-3-methyl-5-
714 propenyl-2,3-dihydro-benzofuran-2yl)-phenol) (Takeyoshi et al. 2004a; 2007a). Of the 31
715 substances with traditional LLNA data, 22 were classified by the traditional LLNA as skin
716 sensitizers and nine were classified as nonsensitizers. The traditional LLNA EC3 values (i.e.,
717 estimated concentration needed to produce an SI = 3) for the 22 sensitizers ranged from
718 0.01% to 47.5% (**Table 3-1**).

719 **Appendix B** provides information on the physicochemical properties (e.g., physical form
720 tested), Chemical Abstracts Service Registry Number, and chemical class for each substance
721 tested. When available, chemical classes for each substance were retrieved from the National
722 Library of Medicine's ChemID Plus database. If chemical classes were unavailable, they
723 were assigned to each test substance using a standard classification scheme based on the
724 National Library of Medicine Medical Subject Headings classification system (available at
725 <http://www.nlm.nih.gov/mesh/meshhome.html>). A substance could be assigned to more than
726 one chemical class; however, no substance was assigned to more than three classes.
727 Chemical class information is presented only to provide an indication of the variety of
728 structural elements that are present in the structures that were evaluated in this analysis.
729 Classification of substances into chemical classes is not intended to indicate the impact of
730 structure on biological activity with respect to sensitization potential. **Table 3-1** shows that
731 18 chemical classes are represented by the substances tested in the LLNA: BrdU-ELISA.

732 Five substances are classified in more than one chemical class. The classes with the highest
 733 number of substances are carboxylic acids (12 substances) and aldehydes (six substances).

734 **Table 3-1 Traditional LLNA EC3 Values and Chemical Classification of Substances**
 735 **Tested in the LLNA: BrdU-ELISA**

Substance Name	Chemical Class ¹	Traditional LLNA EC3 (%) ²	N ³
p-Benzoquinone	Quinones	0.01	1
2,4-Dinitrochlorobenzene *	Hydrocarbon, Halogenated; Nitro Compounds; Hydrocarbons, Cyclic	0.049	15
Diphenylcyclopropenone	Hydrocarbons, Cyclic	0.05	1
Glutaraldehyde	Aldehydes	0.083 ⁴	3
4-Phenylenediamine *	Amines	0.11	6
Formaldehyde	Aldehydes	0.50 ⁴	4
<i>trans</i> -Cinnamaldehyde	Aldehydes	1.4	1
Isoeugenol *	Carboxylic Acids	1.5	47
2-Mercaptobenzothiazole *	Heterocyclic Compounds	1.7 ⁵	1
Cinnamic aldehyde	Aldehydes	1.9	6
3-Aminophenol	Amines; Phenols	3.2	1
Trimellitic anhydride	Anhydrides; Carboxylic Acids	4.7	2
Nickel sulfate	Inorganic Chemicals, Metals Inorganic Chemicals, Elements	4.8 ⁶	1
4-Chloroaniline	Amines	6.5	1
Citral *	Hydrocarbons, Other	9.2	6
Hexyl cinnamic aldehyde *	Aldehydes	9.7	21
Eugenol *	Carboxylic Acids	10.1	11
Cyclamen aldehyde	Aldehydes	22.3	1
Hydroxycitronellal	Hydrocarbons, Other	24.0	6
Linalool	Hydrocarbons, Other	30.0	1
Isopropyl myristate	Lipids	44.0	1
Aniline	Amines	47.5	3
2-Hydroxypropyl methacrylate	Carboxylic Acids	NA	1
Diethyl phthalate	Carboxylic Acids	NA	1
Dimethyl isophthalate	Carboxylic Acids	NA	1
Glycerol	Alcohols; Carbohydrates	NA ⁵	2
Hexane	Hydrocarbons, Acyclic	NA	1
Isopropanol *	Alcohols	NA	1
Lactic acid *	Carboxylic Acids	NA ⁶	1
Methyl salicylate *	Carboxylic Acids	NA	9
Propylene glycol	Alcohols	NA ⁷	1
2,2'-Dihydroxyl-3,3'-dimethoxy-5,5'-diallyl-biphenyl	Carboxylic Acids	NK	0

Substance Name	Chemical Class ¹	Traditional LLNA EC3 (%) ²	N ³
2-Methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-dihydro-benzofuran-2yl)-phenol	Carboxylic Acids	NK	0
4,5'-Diallyl-2'-hydroxy-2,3'-dimethoxyphenyl ether	Carboxylic Acids	NK	0
4-[1-Hydroxy-2-(2-methoxy-4-propenyl-phenoxy)-propyl]-2-methoxy-phenol (Synonym: □-O-4-Dilignol)	Carboxylic Acids	NK	0

736 Abbreviations: LLNA: BrdU-ELISA= Local lymph node assay with enzyme-linked immunosorbent assay
737 detection of bromodeoxyuridine; EC3 = Estimated concentration needed to produce a stimulation index (SI) =
738 3; NA = Not applicable since maximum SI < 3.0; NK = Not known (information not found).

739 *Reference substance from ICCVAM (2009).

740 ¹Chemical classifications based on the Medical Subject Headings classification for chemicals and drugs,
741 developed by the National Library of Medicine (<http://www.nlm.nih.gov/mesh/meshhome.html>).

742 ²Mean EC3 values from the NICEATM database of traditional LLNA studies. Vehicle for testing both
743 sensitizers and nonsensitizers was acetone: olive oil (4:1) unless otherwise noted.

744 ³Number of traditional LLNA studies from which the data were obtained.

745 ⁴Vehicle = Acetone.

746 ⁵Vehicle = *N,N*-Dimethylformamide.

747 ⁶Vehicle = Dimethyl sulfoxide.

748 ⁶Vehicle = Distilled water.

749

749 4.0 Reference Data

750 Twenty-six of the 31 substances previously tested in the traditional LLNA were considered
751 in the original evaluation of the LLNA by ICCVAM (ICCVAM 1999). The traditional LLNA
752 reference data used for the accuracy evaluation described in **Section 6.0** were obtained from
753 ICCVAM (1999) for twenty-four of these substances (**Appendix C**). The traditional LLNA
754 data for the two remaining substances included in the original LLNA evaluation (ICCVAM
755 1999), aniline and nickel sulfate, were obtained from more recent sources, Gerberick et al.
756 (2005) and Ryan et al. (2002), respectively. The traditional LLNA results in ICCVAM
757 (1999) for these two substances were negative, but the subsequent tests at higher
758 concentrations produced positive results. The traditional LLNA data for the remaining five
759 substances that were not considered in the original ICCVAM evaluation (ICCVAM 1999),
760 *trans*-cinnamaldehyde, cyclamen aldehyde, glutaraldehyde, isopropyl myristate, and linalool,
761 were obtained from Gerberick et al. (2005), Basketter et al. (2005), Hilton et al. (1998), Ryan
762 et al. (2000), and Gerberick et al. (2005), respectively.

763 The reference data for the GP tests (guinea pig maximization test [GPMT] or Buehler test)
764 and human tests (human maximization test, human patch test allergen, or other human data)
765 were obtained from Marzulli and Maibach (1974), Opdyke (1976), Bjorkner (1984), Gad et
766 al. (1986), Klecak et al. (1997), ICCVAM (1999), Basketter et al. (1999, 2005), Kwon et al.
767 (2003), Takeyoshi et al. (2004a), and Takeyoshi et al. (2007a). Although there were no
768 traditional LLNA data available for the eugenol dimers (dihydroxyl-3,3'-dimethoxy-5,5'-
769 diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-2,3'-dimethoxyphenyl ether) or the isoeugenol
770 dimers (4-[1-Hydroxy-2-(2-methoxy-4-propenyl-phenoxy)-propyl]-2-methoxy-phenol and
771 2-Methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-dihydro-benzofuran-2yl)-phenol),
772 Takeyoshi et al. (2004a and 2007a, respectively) provided results from the GPMT for these
773 compounds.

774 An independent quality assurance contractor for the NTP audited the traditional LLNA data
775 provided in ICCVAM (1999). Audit procedures and findings are presented in the quality
776 assurance report on file at the National Institute of Environmental Health Sciences. The audit
777 supports the conclusion that the transcribed test data in the submission were accurate,
778 consistent, and complete as compared to the original study records.

779 **5.0 Test Method Data and Results**

780 The test method data in this draft BRD has been revised from the January 2008 draft BRD to
781 include the individual animal results for all of the LLNA: BrdU-ELISA results evaluated in
782 this BRD. The LLNA: BrdU-ELISA data evaluated in this technical summary were obtained
783 from individual animal data that were submitted to NICEATM. These data supported six
784 published studies (Takeyoshi et al. 2003; 2004a; 2004b; 2005; 2006; 2007a), one platform
785 presentation (Takeyoshi et al. 2007b), and one poster presentation (Kojima et al. 2008). Dr.
786 Takeyoshi also submitted unpublished data to NICEATM in January 2009. All test results
787 were obtained using the protocol in **Appendix A**. The substances tested by Takeyoshi et al.
788 were not coded to prevent the possibility of bias in the interpretation of test results. The
789 interlaboratory validation study reported by Kojima et al. (2008), however, used coded test
790 substances to mask the identity of the test substances from the testing laboratories. **Appendix**
791 **C** contains summary data for the LLNA: BrdU-ELISA and comparative reference data for
792 the 35 substances tested in these studies and **Appendix D** contains the individual animal data
793 for the LLNA: BrdU-ELISA.

794

794 **6.0 Test Method Accuracy**

795 The accuracy evaluation in this draft BRD has been revised from the January 2008 draft
796 BRD to include the results for seven additional substances. Other revisions included the
797 evaluation of multiple decision criteria, including the $SI \geq 2.0$ recommended in the test
798 method protocol, and the evaluation of two different criteria used simultaneously to classify
799 sensitizers and nonsensitizers.

800 A critical component of a formal evaluation of the validation status of a test method is an
801 assessment of the accuracy of the proposed tested method when compared to the current
802 reference test method (ICCVAM 2003). Additional comparisons should also be made against
803 available human data, including experience from testing or accidental exposures. This aspect
804 of assay performance is typically evaluated by calculating:

- 805 • Accuracy (concordance): the proportion of correct outcomes (positive and
806 negative) of a test method
- 807 • Sensitivity: the proportion of all positive substances that are classified as
808 positive
- 809 • Specificity: the proportion of all negative substances that are classified as
810 negative
- 811 • False positive rate: the proportion of all negative substances that are
812 incorrectly identified as positive
- 813 • False negative rate: the proportion of all positive substances that are
814 incorrectly identified as negative.

815 **6.1 LLNA: BrdU-ELISA Database Used for the Accuracy Analysis**

816 Thirty-one of the 35 substances listed in **Table 3-1** had sufficient LLNA: BrdU-ELISA and
817 traditional LLNA data to conduct an accuracy analysis. The eugenol dimers (dihydroxyl-3,3'-
818 dimethoxy-5,5'-diallyl-biphenyl and 4,5'-diallyl-2'-hydroxy-2,3'-dimethoxyphenyl ether), and
819 the isoeugenol dimers (4-[1-Hydroxy-2-(2-methoxy-4-propenyl-phenoxy)-propyl]-2-
820 methoxy-phenol and 2-methoxy-4-(7-methoxy-3-methyl-5-propenyl-2,3-dihydro-

821 benzofuran-2yl)-phenol) were excluded from the accuracy analyses because traditional
822 LLNA data for these substances were not identified.

823 Of the 31 substances tested with both LLNA: BrdU-ELISA and the traditional LLNA, 24 had
824 had GP data for a comparison of the performance of the LLNA: BrdU-ELISA vs. GP data
825 with that of the traditional LLNA vs. GP data. No GP data were found for *trans*-
826 cinnamaldehyde, cyclamen aldehyde, diphenylcyclopropenone, hexane, isopropyl myristate,
827 or linalool. Additionally, 3-aminophenol was excluded from the accuracy analyses for the
828 dataset with LLNA: BrdU-ELISA, traditional LLNA, and GP data since the available GP
829 data were generated with a nonstandard GPMT protocol⁴.

830 Of the 31 substances tested with both LLNA: BrdU-ELISA and the traditional LLNA, 29 had
831 human data for a comparison of the performance of the LLNA: BrdU-ELISA vs. human data
832 with that of the traditional LLNA vs. human data. No human data for *trans*-cinnamaldehyde
833 or trimellitic anhydride were located. The complete set of comparative data for each
834 substance is located in **Appendix C**.

835 Multiple tests were available for 14 substances tested with the LLNA: BrdU-ELISA. For the
836 accuracy analyses, results for multiply tested substances were combined so that each
837 substance was represented by one result for the accuracy analysis. In this case, the single
838 result used for each substance represented the outcome that was most prevalent. For example,
839 at $SI \geq 2.0$, isopropanol was a nonsensitizer because five of the seven tests for isopropanol
840 were negative.

841 Discordant test results were noted for three of the substances with multiple test results:
842 formaldehyde, isopropanol, and lactic acid. For all three substances, the solvents used for
843 each test were the same. One of the three laboratories in the interlaboratory validation study
844 reported an SI of 1.97 for formaldehyde; while the others produced $SI > 2$ (Kojima et al.
845 2008). Two of the seven tests of isopropanol yielded $SI \geq 2$ ($SI = 2.0$ and $SI = 2.2$), while the
846 others yielded negative results. These discordant tests were obtained by two of the six
847 laboratories in the interlaboratory validation study. The seventh test of isopropanol yielded SI
848 < 2 (Takeyoshi et al. 2007b). One of the three tests for lactic acid from the interlaboratory

⁴ The nonstandard GP protocol did not include the 48-hour topical patch induction that should follow induction by intradermal injections and it replaced the 24-hour skin patch challenge (usually two weeks after topical induction) with a 6-hour skin patch challenge (Basketter D, personal communication).

849 validation study produced $SI \geq 2$ (i.e., $SI = 2.5$), while the others yielded $SI > 2$ (Kojima et al.
850 2008).

851 **6.2 Accuracy Analysis Using the $SI \geq 2.0$ Decision Criterion**

852 The performance characteristics of the LLNA: BrdU-ELISA were first evaluated using the
853 criterion of $SI \geq 2.0$ to identify sensitizers, which was the threshold for a positive response
854 used in the interlaboratory validation study (the complete protocol used in the validation
855 study is included in **Appendix A**).

856 *6.2.1 Accuracy vs. the Traditional LLNA*

857 When compared to the traditional LLNA and using a decision criteria of $SI \geq 2.0$ to identify
858 sensitizers, the LLNA: BrdU-ELISA had an accuracy of 84% (26/31), a sensitivity of 77%
859 (17/22), a specificity of 100% (9/9), a false positive rate of 0% (0/9), and a false negative rate
860 of 23% (5/22) (**Table 6-1**).

861 *6.2.2 Accuracy vs. Guinea Pig Data*

862 When the accuracy of the LLNA: BrdU-ELISA ($SI \geq 2.0$) and the traditional LLNA were
863 compared based on their performance relative to the GP test, the LLNA: BrdU-ELISA had a
864 lower accuracy (88% [21/24] vs. 100% [24/24]) and sensitivity (81% [13/16] vs. 100%
865 [16/16]), and higher false negative rate (19% [3/16] vs. 0% [0/16]; **Table 6-1**). The
866 specificity (100% [8/8]) and the false positive rate (0% [0/8]) for the LLNA: BrdU-ELISA
867 and the traditional LLNA were the same when they were compared with GP data.

868 *6.2.3 Accuracy vs. Human Data*

869 When the accuracy of the LLNA: BrdU-ELISA ($SI \geq 2.0$) and the traditional LLNA were
870 compared based on their performance relative to the available human data, the LLNA: BrdU-
871 ELISA had a lower accuracy (72% [21/29] vs. 76% [22/29]) and sensitivity (67% [14/21] vs.
872 81% [17/21]) and a higher false negative rate (33% [7/21] vs. 19% [4/21]) than the traditional
873 LLNA (**Table 6-1**). The specificity for the LLNA: BrdU-ELISA was higher (88% [7/8] vs.
874 63% [5/8]) and the false positive rate was lower (12% [1/8] vs. 38% [3/8]) for the LLNA:
875 BrdU-ELISA than that for the traditional LLNA.

876

877 **Table 6-1 Performance of the LLNA: BrdU-ELISA in Predicting Skin Sensitizing Potential Using Decision Criteria of**
 878 **SI ≥ 2.0 to Identify Sensitizers**

Comparison	n ¹	Accuracy		Sensitivity		Specificity		False Positive Rate		False Negative Rate		Positive Predictivity		Negative Predictivity	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
BrdU-ELISA vs. Traditional LLNA	31	84	26/31	77	17/22	100	9/9	0	0/9	23	5/22	100	17/17	64	9/14
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and GP Data</i>															
BrdU-ELISA vs. Traditional LLNA	24	88	21/24	81	13/16	100	8/8	0	0/8	19	3/16	100	13/13	73	8/11
LLNA: BrdU-ELISA vs. GP³	24	88	21/24	81	13/16	100	8/8	0	0/8	19	3/16	100	13/13	73	8/11
Traditional LLNA vs. GP³	24	100	24/24	100	16/16	100	8/8	0	0/8	0	0/16	100	16/16	100	8/8
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and Human Data</i>															
BrdU-ELISA vs. Traditional LLNA	29	83	24/29	75	15/20	100	9/9	0	0/9	25	5/20	100	15/15	64	9/14
LLNA: BrdU-ELISA vs. Human⁴	29	72	21/29	67	14/21	88	7/8	12	1/8	33	7/21	93	14/15	50	7/14
Traditional LLNA vs. Human⁴	29	76	22/29	81	17/21	63	5/8	38	3/8	19	4/21	85	17/20	56	5/9

879 Abbreviations: LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; GP =
 880 guinea pig skin sensitization outcomes; ICCVAM = Interagency Coordinating Committee on the Validation of Alternative Methods; LLNA = murine local
 881 lymph node assay; No. = number.
 882 ¹n = Number of substances included in this analysis.
 883 ²The data on which the percentage calculation is based.
 884 ³GP refers to outcomes obtained by studies conducted using either the Guinea Pig Maximization Test or the Buehler Test.
 885 ⁴Human refers to outcomes obtained by studies conducting using the human maximization test, inclusion of the test substance in a human patch test allergen
 886 kit, and/or published clinical case studies/reports.

887 **6.3 Accuracy Analysis (SI \geq 2.0) Based on the ICCVAM Performance Standards**
888 **Reference Substances**

889 ICCVAM has developed recommended test method performance standards for the traditional
890 LLNA (ICCVAM 2009)⁵, which are proposed to evaluate the performance of modified
891 LLNA test methods that are mechanistically and functionally similar to the traditional
892 LLNA. Because the validation studies for the LLNA: BrdU-ELISA test method were
893 completed prior to the development of LLNA performance standards, the LLNA: BrdU-
894 ELISA is not being evaluated using the ICCVAM-recommended LLNA performance
895 standards. Thus, evaluations of the LLNA: BrdU-ELISA test substances to the ICCVAM-
896 recommended LLNA performance standards test substances are shown to provide a general
897 comparison to a set list of reference substances (18 required reference substances and four
898 optional reference substances) that represent a diverse substance group. As shown in **Table**
899 **6-2**, 10 of the 18 minimum reference substances included in the ICCVAM LLNA
900 Performance Standards have been tested in the LLNA: BrdU-ELISA. Nine of the ten
901 substances yielded the same sensitizer/nonsensitizer outcome in the LLNA: BrdU-ELISA as
902 in the traditional LLNA.

903 **Table 6-3** provides the range and characteristics for 31 substances tested in the LLNA:
904 BrdU-ELISA based on traditional LLNA data. These substances are compared to the range of
905 18 required reference substances included on the ICCVAM-recommended LLNA
906 performance standards reference substances list (ICCVAM 2009). The table indicates that
907 although not all of the 18 required reference substances from the ICCVAM-recommended
908 performance standards reference substances have been tested, the range of the substances
909 tested in the LLNA: BrdU-ELISA is similar to that included in the performance standards
910 list. In general, there are a proportionally increased number of substances tested in the
911 LLNA: BrdU-ELISA in each of the categories included in the table.

912

⁵ Available at http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm.

913 **Table 6-2 Performance of the LLNA: BrdU-ELISA (SI ≥ 2.0) Using the ICCVAM Performance Standards Reference**
 914 **Substances¹**

Substance	Recommended Performance Standards ¹				LLNA: BrdU-ELISA ²			
	Vehicle	Result	EC3 (%) ¹	N ³	Vehicle	Result	EC2 (%)	N ³
5-Chloro-2-methyl-4-isothiazolin-3-one	DMF	+	0.009	1	NT	NT	NT	NT
2, 4-Dinitrochlorobenzene	AOO	+	0.049	15	AOO	+	0.044	8
4-Phenylenediamine	AOO	+	0.11	6	AOO	+	NC	2
Methyl methacrylate	DMF	+	90	1	NT	NT	NT	NT
Isoeugenol	AOO	+	1.5	47	AOO	+	7.6	2
2-Mercaptobenzothiazole	DMF	+	1.7	1	DMF	-	NA (-)	1
Cobalt chloride	DMSO	+	0.6	2	NT	NT	NT	NT
Citral	AOO	+	9.2	6	AOO	+	NC	1
Hexyl cinnamic aldehyde	AOO	+	9.7	21	AOO	+	17.4	11
Eugenol	AOO	+	10.1	11	AOO	+	9.8	8
Phenyl benzoate	AOO	+	13.6	3	NT	NT	NT	NT
Cinnamic alcohol	AOO	+	21	1	NT	NT	NT	NT
Imidazolidinyl urea	DMF	+	24	1	NT	NT	NT	NT
Chlorobenzene	AOO	-	NA	1	NT	NT	NT	NT
Isopropanol	AOO	-	NA	1	AOO	-	NA (-) ⁴	7
Lactic acid	DMSO	-	NA	1	DMSO	+	NA (-) ⁵	3
Methyl salicylate	AOO	-	NA	9	AOO	NT	NA (-)	3
Salicylic acid	AOO	-	NA	1	NT	NT	NT	NT
Ethylene glycol dimethacrylate	MEK	False +	28 (FP)	1	NT	NT	NT	NT
Sodium lauryl sulfate	DMF	False +	8.1 (FP)	5	NT	NT	NT	NT
Nickel chloride	DMSO	False -	NA (FN)	2	NT	NT	NT	NT
Xylene	AOO	False -	95.8 (FP)	1	NT	NT	NT	NT

915 Bolded italics text highlights discordant LLNA: BrdU-ELISA vs. traditional LLNA test results.
 916 Abbreviations: AOO = acetone: olive oil (4: 1); LLNA: BrdU-ELISA= murine local lymph node assay with enzyme-linked immunosorbent assay detection
 917 of bromodeoxyuridine; DMF = *N,N*-dimethylformamide; DMSO = dimethyl sulfoxide; EC3 = estimated concentration needed to produce a stimulation index
 918 of 3; EC2 = estimated concentration needed to produce a stimulation index of 2; FN = false negative in traditional LLNA when compared to guinea pig
 919 and/or human results; FP = false positive in traditional LLNA when compared to guinea pig and/or human results; LLNA = murine local lymph node assay;
 920 MEK = methyl ethyl ketone; NA = not applicable; NC = not calculated; only one concentration tested; NT = not tested; SI = Stimulation index.
 921 + = Sensitizer.
 922 - = Nonsensitizer.

923 ¹From *Recommended Performance Standards: Murine Local Lymph Node Assay* (ICCVAM 2009; available:
924 http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm.
925 ²Calculated from data supporting Takeyoshi et al. (2003, 2004b, 2005, 2006, 2007a, 2007b, and unpublished) and Kojima et al (2008). Substances for which
926 EC2 values were not available include the outcome of the LLNA: BrdU-ELISA test (+ = sensitizer; - = nonsensitizer) in parentheses.
927 ³Number of values used to derive the mean EC3 or EC2.
928 ⁴Based on the most prevalent outcome (i.e., 5/7 tests yielded SI < 2).
929 ⁵Based on the most prevalent outcome (i.e., 2/3 tests yielded SI < 2).
930

931 **Table 6-3 Characteristics of the Substances Tested in the LLNA: BrdU-ELISA vs.**
 932 **the ICCVAM Performance Standards Reference Substances¹**

EC3 Range (%)	No. Chems	Solid/Liquid	Actual EC3 Range (%)	Maximum SI Range	Human Data	Peptide Reactivity (Hi/Mod/Min/Lo/Unk) ³
<0.1	4	3/1	0.01 - 0.083	18.0 –59.0	4	4/0/0/0/0
	2	1/1	0.009 - 0.05	22.6 - 52.3	2	2/0/0/0/0
≥ 0.1 to <1	2	1/1	0.11 - 0.50	4.0 – 26.4	2	0/1/0/0/1
	2	2/0	0.11 - 0.6	6.7 - 75.3	2	0/0/0/0/2
≥ 1 to <10	10	4/6	1.4 - 9.7	3.1 – 31.0	8	2/0/1/1/6
	4	1/3	1.5 - 9.7	8.6 - 29.5	4	1/0/1/0/2
≥ 10 to <100	5	0/5	10.1 – 47.5	3.4 - 17.0	5	0/0/1/2/2
	5	3/2	10.1 - 90	5.5 - 70.3	5	0/1/0/0/4
Negative	10	2/8	NC	1.0 – 2.9	10	0/0/7/1/2
	5	1/4	NC	0.9 - 2.8	3	0/0/2/0/3
Overall	31	10/21	0.01 – 47.5	0.9 - 28.6	29	6/1/9/4/11
	18	10/8	0.009 - 24	0.9 - 75.3	16	3/1/3/0/11

933 Bolded text represents characteristics of the LLNA: BrdU-ELISA database.

934 Abbreviations: Chems = chemicals; EC3 = Estimated concentration needed to produce SI = 3; LLNA: BrdU-
 935 ELISA= murine local lymph node assay with enzyme-linked immunosorbent assay detection of
 936 bromodeoxyuridine; NC = Not calculated because maximum SI < 3; No. = number; Lo = low; Min = minimal;
 937 Mod = moderate; SI = stimulation index; Unk = unknown.

938 ¹From *Recommended Performance Standards: Murine Local Lymph Node Assay* (ICCVAM 2009; available:
 939 http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm. Includes the 18 "required" substances for
 940 testing.

941 ²Data obtained from: Gerberick et al. (2007)

942 **6.4 Discordant Results for Accuracy Analysis Using the SI ≥ 2.0 Decision** 943 **Criterion**

944 6.4.1 *Discordance Between the LLNA: BrdU-ELISA and the Traditional LLNA*

945 When the outcomes for the 31 substances tested in the LLNA: BrdU-ELISA (using SI ≥ 2.0)
 946 and the traditional LLNA were compared, the classifications for five substances were
 947 different. The LLNA: BrdU-ELISA classified aniline, cyclamen aldehyde,
 948 hydroxycitronellal, 2-mercaptobenzothiazole, and linalool as nonsensitizers while the
 949 traditional LLNA classified them as sensitizers (i.e., false negative outcome) (**Table 6-4**).
 950 The substances were tested in the same vehicle in both the LLNA: BrdU-ELISA and the

951 traditional LLNA tests. The only commonality noted among these four substances was their
952 molecular weights (MW), which range from 93 to 172 g/mole. No commonalities in
953 chemical class, physical form, peptide reactivity (see **Appendix B** for physical/chemical
954 information), or potential for skin irritation were noted among these substances.

- 955 • Aniline (MW = 93.13 g/mole) is an amine, cyclamen aldehyde is a carboxylic
956 acid (MW = 190.28 g/mole), 2-mercaptobenzothiazole (MW = 167.26 g/mole)
957 is a heterocyclic compound, and hydroxycitronellal (MW = 172.26 g/mole)
958 and linalool (MW = 154.25 g/mole) are hydrocarbons.
- 959 • Aniline, cyclamen aldehyde, hydroxycitronellal, and linalool are liquids, while
960 2-mercaptobenzothiazole is a solid.
- 961 • Of the three substances for which peptide reactivity information was
962 available, hydroxycitronellal and cyclamen aldehyde had low peptide
963 reactivity and 2-mercaptobenzothiazole had high peptide reactivity.
- 964 • Cyclamen aldehyde, 2-mercaptobenzothiazole and linalool are skin irritants at
965 the concentrations tested in the LLNA: BrdU-ELISA, and linalool is also a
966 skin irritant at the concentrations tested the traditional LLNA.
- 967 • None of the five discordant substances generated strongly positive result in
968 the traditional LLNA (EC3 = 1.7% to 47.5%).

969

969 **Table 6-4** **Discordant Results for LLNA: BrdU-ELISA (Using SI \geq 2.0 for**
 970 **Sensitizers) Compared to Traditional LLNA and Guinea Pig Reference**
 971 **Data¹**

Substance Name ²	Vehicle ³	LLNA: BrdU-ELISA ⁴	Traditional LLNA ⁴	Guinea Pig Studies	Skin Irritant?
Aniline (47.5%)	AOO	- (1.5, 50%)	+ (3.6, 100%) ⁵	+	Negative at \leq 100%
Hydroxycitronellal (24.0%)	AOO	- (1.3, 100%)	+ (8.5, 100%)	+	Negative at \leq 100%
Cyclamen aldehyde (22.3%)	AOO	- (1.97, 100%)	+ (5.2, 50%)	NA	Irritant at 100%
Linalool (30.0%)	AOO	- (1.45, 100%) ⁵	+ (8.3, 100%)	NA	Mild irritant at 10%
2-Mercaptobenzothiazole (1.7%)	DMF	- (1.62, 50%) ⁶	+ (8.6, 10%)	+	Negative at \leq 10%

972 Abbreviations: LLNA: BrdU-ELISA= murine local lymph node assay with enzyme-linked immunosorbent
 973 assay detection of bromodeoxyuridine; GP = outcomes of guinea pig skin sensitization tests; LLNA = murine
 974 local lymph node assay; NA = not available; SI = stimulation index.

975 + = Sensitizer.

976 - = Nonsensitizer.

977 ¹Data sources provided in **Appendix C-1**.

978 ²Numbers in parentheses are the EC3 values for the traditional LLNA (from **Table 3-1**).

979 ³Vehicles apply to tests for both the LLNA: BrdU-ELISA and the traditional LLNA.

980 ⁴Numbers in parentheses are highest SI values and maximum concentrations tested.

981 ⁵Highest SI occurred at concentration of 50%.

982 ⁶Highest SI occurred at concentration of 12.5%.

983 6.4.2 *Discordance Among the LLNA: BrdU-ELISA, the Traditional LLNA, and/or the* 984 *Guinea Pig Test*

985 For the 24 substances with LLNA: BrdU-ELISA, traditional LLNA, and GP test results, the
 986 results for aniline, hydroxycitronellal, and 2-mercaptobenzothiazole were also discordant
 987 with the GP test results (**Table 6-4**). The LLNA: BrdU-ELISA results for aniline,
 988 hydroxycitronellal, and 2-mercaptobenzothiazole were negative, while the traditional LLNA
 989 and GP results were positive. No guinea pig results were available for linalool or cyclamen
 990 aldehyde, which were negative in the LLNA: BrdU-ELISA and positive in the traditional
 991 LLNA. As noted in **Section 6.3.1**, there were no commonalities associated with these
 992 discordant substances.

993

993 6.4.3 *Discordance Among the LLNA: BrdU-ELISA, the Traditional LLNA, and/or the*
994 *Human Outcome*

995 When analyses were restricted to the 29 substances with LLNA: BrdU-ELISA, traditional
996 LLNA, and human outcomes, the LLNA: BrdU-ELISA misclassified eight substances. Both
997 the LLNA: BrdU-ELISA and the traditional LLNA misclassified four human sensitizers
998 (diethyl phthalate, 2-hydroxypropylmethacrylate, isopropanol, and propylene glycol) as
999 nonsensitizers (**Table 6-5**). The LLNA: BrdU-ELISA also misclassified three other
1000 sensitizers as nonsensitizers that were correctly classified by the traditional LLNA (aniline,
1001 2-mercaptobenzothiazole, and hydroxycitronellal).

1002 The eighth misclassified substance was isopropyl myristate, which was misclassified as a
1003 sensitizer by the LLNA: BrdU-ELISA and the traditional LLNA. Isopropyl myristate
1004 exhibited a weak response in the traditional LLNA (EC3 = 44%). It was tested in both
1005 methods at concentrations that were not irritating to skin (**Table 6-5**). Isopropyl myristate
1006 (MW = 270.46 g/mole) is a liquid lipid that exhibits low peptide reactivity.

1007 6.4.4 *Discordance Between the LLNA: BrdU-ELISA and the Traditional LLNA When*
1008 *Testing the LLNA Performance Standards Substances*

1009 There was one discordant substance (2-mercaptobenzothiazole) noted among the 10
1010 performance standards reference substances that were tested in LLNA: BrdU-ELISA. The
1011 LLNA: BrdU-ELISA classified this substance as a nonsensitizer, while the traditional LLNA,
1012 GP, and human tests classified it as a sensitizer. The EC3 value for the traditional LLNA,
1013 1.7%, was derived from a test of 1, 3, and 10% 2-mercaptobenzothiazole in *N,N*-
1014 dimethylformamide (Gerberick et al 2005). The maximum SI was 8.6 at 10%. The LLNA:
1015 BrdU-ELISA test used the same vehicle and tested concentrations of 12.5%, 25%, 50% 2-
1016 mercaptobenzothiazole which yielded SI values of 1.6, 1.4 and 1.5, respectively.

1017

1017 **Table 6-5 Discordant Results for LLNA: BrdU-ELISA (SI \geq 2.0) When Compared**
 1018 **to Traditional LLNA and Human Outcome Data¹**

Substance Name ²	Vehicle ³	LLNA: BrdU-ELISA ⁴	Traditional LLNA ⁴	Human Outcome ⁵	Skin Irritant?
Diethyl phthalate	AOO	- (0.9, 50%)	- (1.5, 100%)	+ (HPTA)	Negative at \leq 100%
2-Hydroxypropylmethacrylate	AOO	- (1.1, 50%)	- (1.3, 50%)	+ (case study, 0.1%)	Negative at \leq 10%
Isopropanol	AOO	- (2.2, 50%) ⁶	- (1.7, 50%) ⁷	+ (case study, 0.001%)	Negative at \leq 100%
Propylene glycol	AOO	- (1.6, 50%)	- (1.6, 100%) ⁸	+ (HPTA)	Negative at \leq 25%
Aniline (47.5%)	AOO	- (1.5, 50%)	+ (3.6, 100%) ⁹	+ (7/25, 20%)	Negative at \leq 100%
2-Mercaptobenzothiazole (1.7%)	DMF	- (1.62, 50%) ¹⁰	+ (8.6, 10%)	+ (5/24, 10%)	Negative at \leq 10%
Hydroxycitronellal (24.0%)	AOO	- (1.3, 100%)	+ (8.5, 100%)	+ (14/73, 20%)	Negative at \leq 100%
Isopropyl myristate (44%)	AOO	+ (4.2, 50%)	+ (3.4, 100%)	- (0/25, 20%)	Negative at \leq 100%

1019 Abbreviations: LLNA: BrdU-ELISA= murine local lymph node assay with enzyme-linked immunosorbent
 1020 assay detection of bromodeoxyuridine; HPTA = human patch test allergen; LLNA = murine local lymph node
 1021 assay.

1022 + = Sensitizer.

1023 - = Nonsensitizer.

1024 ¹Data sources listed in **Appendix C-1**.

1025 ²Numbers in parentheses are EC3 values for the traditional LLNA (from **Table 3-1**).

1026 ³Vehicles apply to tests for both the LLNA: BrdU-ELISA and the traditional LLNA unless otherwise noted.

1027 ⁴Numbers in parentheses are highest SI values and maximum concentrations tested.

1028 ⁵Information in parentheses indicates the basis for the human outcome. Numbers indicate the incidence of
 1029 positive human response and concentration.

1030 ⁶Negative based on most prevalent call. Highest SI of any test is shown. Highest SIs for most tests occurred at <
 1031 50%.

1032 ⁷Highest SI occurred at 10%.

1033 ⁸Vehicle for the traditional LLNA was distilled water.

1034 ⁹Highest SI occurred at 50%.

1035 ¹⁰Highest SI occurred at 12.5%.

1036 **6.5 LLNA: BrdU-ELISA Accuracy Analysis Using One Alternative Decision** 1037 **Criterion**

1038 In addition to the accuracy analysis using SI \geq 2.0 to classify substances as sensitizers, other
 1039 decision criteria were evaluated for test method performance with the traditional LLNA
 1040 serving as the reference test. The performance characteristics for 13 different decision criteria
 1041 for determining whether the skin sensitization potential for the substances were positive or

1042 negative are reported in this section. The substances evaluated were the 31 substances with
1043 both LLNA: BrdU-ELISA and traditional LLNA data discussed in **Section 6.1**. The decision
1044 criteria included:

- 1045 1. SI values ≥ 1.3 , ≥ 1.5 , ≥ 2.0 , ≥ 2.5 , ≥ 3.0 , ≥ 3.5 , ≥ 4.0 , ≥ 4.5 , or ≥ 5.0
- 1046 2. Statistically significant difference between any treatment group and the
1047 vehicle control group. Absorbance values of treated groups were compared
1048 with vehicle control group using ANOVA with a post-hoc Dunnett's test,
1049 when multiple treatment groups were tested, or Student's t-test when only
1050 there was only one treatment group.
- 1051 3. Mean absorbance values of treated groups $\geq 95\%$ CI of the control group
- 1052 4. Mean absorbance values of treated groups ≥ 2 SD or ≥ 3 SD from the control
1053 group mean

1054 Multiple tests were available for 14 substances tested with the LLNA: BrdU-ELISA. The
1055 results for each of these substances were combined so that each substance was represented by
1056 one sensitizer or nonsensitizer result for each criterion evaluated for the accuracy analysis.
1057 The results were combined in three ways and a separate accuracy analysis was performed for
1058 each approach.

- 1059 1. The sensitizer/nonsensitizer outcome for each substance was most prevalent
1060 outcome for each criterion. For example, for the criterion for a statistical
1061 difference between control and treatment groups, two of the three lactic acid
1062 tests produced sensitizer results. Thus, the single outcome for lactic acid for
1063 the accuracy analysis was a sensitizer result. If the number of positive and
1064 negative outcomes were equal, the most conservative (i.e., positive) result was
1065 used for the accuracy analyses.
- 1066 2. The positive/negative outcome for each substance for each criterion was
1067 determined by the outcome of the test with the highest maximum SI of the
1068 multiple tests.

1069 3. The positive/negative outcome for each substance was determined by the
1070 outcome of the test with the lowest maximum SI of the multiple tests.

1071 The analysis presented here is based on using the most prevalent outcome for substances with
1072 multiple tests, as this is representative of the most likely outcome for a given chemical. The
1073 analyses using the highest maximum SI and the lowest maximum SI are detailed in
1074 **Appendix E.**

1075 As shown in **Section 6.1**, using the most prevalent outcome and the decision criterion of $SI \geq$
1076 2.0 resulted in an accuracy of 84% (26/31), a sensitivity of 77% (17/22), a specificity of
1077 100% (9/9), a false positive rate of 0% (0/9), and a false negative rate of 23% (5/22) (**Tables**
1078 **6-1** and **6-6**). Although using $SI \geq 2.5$ produced the same results as $SI \geq 2.0$, using higher SI
1079 values (i.e., $SI \geq 3.0$ to $SI \geq 5.0$) resulted in reduced overall accuracy, higher false negative
1080 rates, and lower false positive rates as compared to $SI \geq 2.0$ (**Figure 6-1** and **Table 6-6**).
1081 Using a lower SI value ($SI \geq 1.5$) as the decision criterion produced the same accuracy as SI
1082 ≥ 2.0 (84% [26/31]), but the false positive rate increased to 33% (3/9), and the false negative
1083 rate decreased to 9% (2/22). $SI \geq 1.3$ is shown for comparison as it was previously
1084 recommended by ICCVAM, but was considered to be inadequate by the March 2008 Peer
1085 Review Panel (ICCVAM 2008). Use of ANOVA and summary statistics (i.e., mean
1086 absorbance values of treated groups \geq 95% confidence interval of the control group, or ≥ 2 or
1087 3 SD from the control group mean), yielded accuracy values of 81 to 87%, with false
1088 negative rates of 0 to 14%, and false positive rates were 11 to 56%.

1089 **Table 6-6 Performance of the LLNA: BrdU-ELISA in Predicting Skin Sensitizing Potential Using Alternative Decision**
 1090 **Criteria to Identify Sensitizers and the Most Prevalent Outcome for Substances with Multiple Tests**

Alternate Criterion	N ¹	Accuracy		Sensitivity		Specificity		False Positive Rate		False Negative Rate		Positive Predictivity		Negative Predictivity	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
Statistics ³	31	81	25/31	86	19/22	67	6/9	33	3/9	14	3/22	86	19/22	67	6/9
≥ 95% CI ⁴	31	84	26/31	100	22/22	44	4/9	56	5/9	0	0/22	82	22/27	100	4/4
≥ 2 SD ⁵	31	84	26/31	96	21/22	56	5/9	44	4/9	5	1/22	84	21/25	83	5/6
≥ 3 SD ⁶	31	87	27/31	86	19/22	89	8/9	11	1/9	14	3/22	95	19/20	73	8/11
SI ≥ 5.0	31	58	18/31	41	9/22	100	9/9	0	0/9	59	13/22	100	9/9	41	9/22
SI ≥ 4.5	31	58	18/31	41	9/22	100	9/9	0	0/9	59	13/22	100	9/9	41	9/22
SI ≥ 4.0	31	71	22/31	59	13/22	100	9/9	0	0/9	41	9/22	100	13/13	50	9/18
SI ≥ 3.5	31	74	23/31	64	14/22	100	9/9	0	0/9	36	8/22	100	14/14	53	9/17
SI ≥ 3.0	31	77	24/31	68	15/22	100	9/9	0	0/9	32	7/22	100	15/15	56	9/16
SI ≥ 2.5	31	84	26/31	77	17/22	100	9/9	0	0/9	23	5/22	100	17/17	64	9/14
SI ≥ 2.0	31	84	26/31	77	17/22	100	9/9	0	0/9	23	5/22	100	17/17	64	9/14
SI ≥ 1.5	31	84	26/31	91	20/22	67	6/9	33	3/9	9	2/22	87	20/23	75	6/8
SI ≥ 1.3	31	87	27/31	100	22/22	56	5/9	44	4/9	0	0/22	85	22/26	100	5/5

1091 Abbreviations: LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine (BrdU); CI = confidence
 1092 interval; No. = number; SD = standard deviation; SI = stimulation index

1093 ¹ N = Number of substances included in this analysis.

1094 ² The proportion on which the percentage calculation is based.


1095 ³ Analysis of variance for difference of group means when substances were tested at multiple doses or *t*-test when substances were tested at one dose. The absorbance data were
 1096 log-transformed prior to analysis of variance. Significance at *p* < 0.05 was further tested by Dunnett's test.

1097 ⁴ The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.

1098 ⁵ The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.

1099 ⁶The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group

1100 **Figure 6-1 Performance of the LLNA: BrdU-ELISA with SI Compared to the**
 1101 **Traditional LLNA Using the Most Prevalent Outcome for Substances**
 1102 **with Multiple Tests**



1103
 1104 As compared to traditional LLNA results, the lines show the change in performance characteristics
 1105 for the LLNA: BrdU-ELISA with the SI cutoff used to identify sensitizers. This analysis used LLNA:
 1106 BrdU-ELISA and traditional LLNA results for 29 substances (20 sensitizers and nine nonsensitizers).
 1107 For the 14 substances with multiple test results, the results for each substance were combined using
 1108 the most prevalent outcome. The solid line shows accuracy, the dashed green line shows the false
 1109 positive rate, and the dotted red line shows the false negative rate.

1110
 1111 The decision criteria of $SI \geq 1.5$ and mean absorbance of treated group ≥ 2 SD from the
 1112 vehicle control group are compared with $SI \geq 2.0$ for accuracy of the LLNA: BrdU-ELISA
 1113 against GP and human data in **Table 6-7**. When GP test results were used as the reference
 1114 data, $SI \geq 1.5$ and mean absorbance of treated group ≥ 2 SD from the vehicle control group
 1115 had the same accuracy (88% [21/24]), lower false negative rates (6% [1/16] for $SI \geq 1.5$ vs.
 1116 0% [0/16] for mean absorbance of treated group ≥ 2 SD from the vehicle control group vs.

1117 19% [3/16] for $SI \geq 2.0$), and increased false positive rate (25% [2/8] for $SI \geq 1.5$ vs. 38%
1118 [3/8] for mean absorbance of treated group ≥ 2 SD from the vehicle control group vs. 0%
1119 [0/8] for $SI \geq 2.0$) when compared with $SI \geq 2.0$. Using mean absorbance of treated group ≥ 2
1120 SD from the vehicle control group had the most impact on the false negative rate. It
1121 decreased the number of false negatives from three (for $SI \geq 2.0$) to zero, but the number of
1122 false positives increased from zero (for $SI \geq 2.0$) to three.

1123 When results were compared to human data, $SI \geq 1.5$ and mean absorbance of treated group
1124 ≥ 2 SD from the vehicle control group produced the same the accuracy (72% [21/29]),
1125 decreased the false negative rate (19% [4/21] for $SI \geq 1.5$ vs. 14% [3/21] for mean
1126 absorbance of treated group ≥ 2 SD from the vehicle control group vs. 33% [7/21] for $SI \geq$
1127 2.0), and increased the false positive rate (50% [4/8] for $SI \geq 1.5$ vs. 63% [5/8] for mean
1128 absorbance of treated group ≥ 2 SD from the vehicle control group vs. 12% [1/8] for $SI \geq$
1129 2.0) compared with $SI \geq 2.0$. Using mean absorbance of treated group ≥ 2 SD from the
1130 vehicle control group had the most impact on the false negative rate. It decreased the number
1131 of false negatives from seven (for $SI \geq 2.0$) to three, but the number of false positives
1132 increased from one (for $SI \geq 2.0$) to five.

1133

1134 **Table 6-7 Comparison of Performance for Decision Criteria of SI ≥ 1.5 (Bold), > 2 SD (Bold Italics), and SI ≥ 2.0 for**
 1135 **Predicting Skin Sensitizing Potential with LLNA: BrdU-ELISA**

Comparison	n ¹	Accuracy		Sensitivity		Specificity		False Positive Rate		False Negative Rate		Positive Predictivity		Negative Predictivity	
		%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²	%	No. ²
BrdU-ELISA vs. Traditional LLNA	31	84	26/31	91	20/22	67	6/9	33	3/9	9	2/22	87	20/23	75	6/8
		<i>84</i>	<i>26/31</i>	<i>96</i>	<i>21/22</i>	<i>56</i>	<i>5/9</i>	<i>44</i>	<i>4/9</i>	<i>5</i>	<i>1/22</i>	<i>84</i>	<i>21/25</i>	<i>83</i>	<i>5/6</i>
		84	26/31	77	17/22	100	9/9	0	0/9	23	5/22	100	17/17	64	9/14
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and GP Data</i>															
BrdU-ELISA vs. Traditional LLNA	24	88	21/24	94	15/16	75	6/8	25	2/8	6	1/16	88	15/17	86	6/7
		<i>88</i>	<i>21/24</i>	<i>100</i>	<i>16/16</i>	<i>63</i>	<i>5/8</i>	<i>38</i>	<i>3/8</i>	<i>0</i>	<i>0/16</i>	<i>84</i>	<i>16/19</i>	<i>100</i>	<i>5/5</i>
		88	21/24	81	13/16	100	8/8	0	0/8	19	3/16	100	13/13	73	8/11
LLNA: BrdU-ELISA vs. GP ³	24	88	21/24	94	15/16	75	6/8	25	2/8	6	1/16	88	15/17	86	6/7
		<i>88</i>	<i>21/24</i>	<i>100</i>	<i>16/16</i>	<i>63</i>	<i>5/8</i>	<i>38</i>	<i>3/8</i>	<i>0</i>	<i>0/16</i>	<i>84</i>	<i>16/19</i>	<i>100</i>	<i>5/5</i>
		88	21/24	81	13/16	100	8/8	0	0/8	19	3/16	100	13/13	73	8/11
Traditional LLNA vs. GP ³	24	100	24/24	100	16/16	100	8/8	0	0/8	0	0/16	100	16/16	100	8/8
<i>Substances with LLNA: BrdU-ELISA, Traditional LLNA, and Human Data</i>															
BrdU-ELISA vs. Traditional LLNA	29	83	24/29	90	18/20	67	6/9	33	3/9	10	2/20	86	18/21	75	6/8
		<i>83</i>	<i>24/29</i>	<i>95</i>	<i>19/20</i>	<i>56</i>	<i>5/9</i>	<i>44</i>	<i>4/9</i>	<i>5</i>	<i>1/20</i>	<i>83</i>	<i>19/23</i>	<i>83</i>	<i>5/6</i>
		83	24/29	75	15/20	100	9/9	0	0/9	25	5/20	100	15/15	64	9/14
LLNA: BrdU-ELISA vs. Human ⁴	29	72	21/29	81	17/21	50	4/8	50	4/8	19	4/21	81	17/21	50	4/8
		<i>72</i>	<i>21/29</i>	<i>86</i>	<i>18/21</i>	<i>38</i>	<i>3/8</i>	<i>63</i>	<i>5/8</i>	<i>14</i>	<i>3/21</i>	<i>78</i>	<i>18/23</i>	<i>50</i>	<i>3/6</i>
		72	21/29	67	14/21	88	7/8	12	1/8	33	7/21	93	14/15	50	7/14
Traditional LLNA vs. Human ⁴	29	76	22/29	81	17/21	63	5/8	38	3/8	19	4/21	85	17/20	56	5/9

1136 Abbreviations: LLNA: BrdU-ELISA= murine local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; GP = guinea pig skin
 1137 sensitization outcomes; LLNA = murine local lymph node assay; No. = number.

1138 ¹n = Number of substances included in this analysis.

1139 ²The data on which the percentage calculation is based.

1140 ³GP refers to outcomes obtained by studies conducted using either the guinea pig maximization test or the Buehler test.

1141 ⁴Human refers to outcomes obtained by studies conducting using the human maximization test, inclusion of the test substance in a human patch test allergen kit, and/or
 1142 published clinical case studies/reports.

1143 **6.6 Discordant Results for Accuracy Analysis Using One Alternative Decision**
1144 **Criterion**

1145 This section discusses the discordant results obtained for the analyses using the alternative
1146 decision criteria shown in **Tables 6-6** and **6-7** to provide a comparison to the discordant
1147 substances identified using the decision criteria of $SI \geq 2.0$ to identify sensitizers. Discordant
1148 results are first discussed using the traditional LLNA as the reference test (**Section 6.5.1**) and
1149 then discordant results for $SI \geq 1.5$, the optimized single criterion, are discussed using the
1150 traditional LLNA, GP, and human outcomes as references (**Section 6.5.2**).

1151 6.6.1 *Discordant Results Using Alternative Decision Criteria Compared with the*
1152 *Traditional LLNA*

1153 Using decision criteria of $SI \geq 2.0$ and the most prevalent outcome for the substances with
1154 multiple tests, the five discordant substances, when compared to the traditional LLNA, were
1155 aniline, cyclamen aldehyde, hydroxycitronellal, linalool, and 2-mercaptobenzothiazole
1156 (**Table 6-4**). As indicated in **Section 6.3**, all five substances were false negatives when
1157 compared to the traditional LLNA.

1158 **Table 6-8** shows how the number and identity of discordant substances changes with the
1159 alternate decision criteria when using the most prevalent outcome for the substances with
1160 multiple tests. Use of a statistical test (i.e., ANOVA or *t*-test; “Statistics” in **Table 6-6**) or
1161 summary statistics (i.e., $\geq 95\%$ CI, or ≥ 2 or 3 SD in **Table 6-6**) did not result in
1162 substantively improved performance relative to using $SI \geq 1.5$.

1163 **Table 6-8** **Discordant Results for LLNA: BrdU-ELISA Using Alternative Decision Criteria Compared to the**
 1164 **Traditional LLNA and the Most Prevalent Outcome for Substances with Multiple Tests**

Discordant Substance ¹	Alternate Decision Criterion ²												
	Statistics ³	\geq 95% CI ⁴	\geq 2 SD ⁵	\geq 3 SD ⁶	SI \geq 5.0	SI \geq 4.5	SI \geq 4.0	SI \geq 3.5	SI \geq 3.0	SI \geq 2.5	SI \geq 2.0	SI \geq 1.5	SI \geq 1.3
Formaldehyde (0.53%)					-	-							
<i>trans</i> -Cinnamic aldehyde (1.4%)					-	-							
2-Mercaptobenzothiazole (1.7%)	-				-	-	-	-	-	-	-		
Cinnamic aldehyde (2.4%)					-	-							
3-Aminophenol (3.2%)					-	-	-	-					
Nickel sulfate (4.8%)					-	-	-	-	-				
4-Chloroaniline (6.5%)					-	-	-	-	-				
Hexyl cinnamic aldehyde (9.7%)					-	-	-						
Cyclamen aldehyde (22.3%)					-	-	-	-	-	-	-		
Hydroxycitronellal (24%)				-	-	-	-	-	-	-	-	-	
Linalool (30%)			-	-	-	-	-	-	-	-	-	-	
Isopropyl myristate (44%)					-	-							
Aniline (63%)	-			-	-	-	-	-	-	-	-		
Glycerol (-)	+	+	+										
Hexane (-)	+	+	+	+								+	+
Lactic acid (-)	+	+	+									+	+
Methyl salicylate (-)		+											+
Propylene glycol (-)		+	+									+	+

1165 Abbreviations: LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; CI = confidence interval;
1166 SD = standard deviation; SI = stimulation index.
1167 ¹LLNA: BrdU-ELISA outcomes are indicated by “+” for sensitizer results and “-” for nonsensitizer results.
1168 ²Compared to the traditional LLNA. Traditional LLNA result in parentheses: “-” for nonsensitizers and EC3 (%) for sensitizers.
1169 ³Analysis of variance for difference of group means when substances were tested at multiple doses or t-test when substances were tested at one dose. The absorbance data
1170 were log-transformed prior to analysis of variance. Significance at $p < 0.05$ was further tested by Dunnett’s test.
1171 ⁴The mean absorbance of at least one treatment group was outside the 95% confidence interval for the mean absorbance of the vehicle control group.
1172 ⁵The mean absorbance of at least one treatment group was greater than 3 SD from the mean absorbance of the vehicle control group.
1173 ⁶The mean absorbance of at least one treatment group was greater than 2 SD from the mean absorbance of the vehicle control group.
1174

1175 Four ICCVAM performance standards reference substances were discordant for the analysis
1176 of alternate decision criteria using the most prevalent outcome for substances with multiple
1177 tests (**Table 6-6**). Two sensitizers (2-mercaptobenzothiazole and hexyl cinnamic aldehyde)
1178 were misclassified by some criteria as nonsensitizers, and two nonsensitizers (lactic acid and
1179 methyl salicylate) were misclassified as sensitizers by some criteria. The criteria that yielded
1180 the correct results for 2-mercaptobenzothiazole included summary statistics (i.e., $\geq 95\%$ CI, \geq
1181 2 SD, or $3 \geq$ SD) and $SI \geq 1.5$. The criteria that yielded the correct results for hexyl cinnamic
1182 aldehyde included statistical tests (i.e., ANOVA or *t*-test), summary statistics (i.e., $\geq 95\%$ CI,
1183 ≥ 2 SD, or $3 \geq$ SD), and $SI \geq 1.5$ to 3.5. The criteria that yielded the correct results for lactic
1184 acid included treatment group mean ≥ 3 SD from the vehicle control, and $SI \geq 2.0$ to 5.0. All
1185 criteria yielded the correct results for methyl salicylate except for treatment group absorbance
1186 $\geq 95\%$ CI of vehicle control mean.

1187 6.6.2 *Discordant Results for Accuracy Analysis of the $SI \geq 1.5$ Decision Criterion*

1188 When the outcomes for the 31 substances tested in the LLNA: BrdU-ELISA (using $SI \geq 1.5$)
1189 and the traditional LLNA were compared, the classifications for five substances were
1190 different. For the three substances with GP data, the GP tests and traditional LLNA yielded
1191 the same sensitizer/nonsensitizer results (**Table 6-9**). Two substances were misclassified in
1192 the LLNA: BrdU-ELISA as nonsensitizers (hydroxycitronellal and linalool) and three were
1193 misclassified as sensitizers (hexane, lactic acid and propylene glycol). Chemical class,
1194 physical form, MW, peptide reactivity (see **Appendix B** for physical/chemical properties),
1195 traditional LLNA EC3 range, or potential for skin irritation were examined to identify
1196 commonalities among the discordant substances. For the two substances misclassified as
1197 nonsensitizers:

- 1198 • Hydroxycitronellal (MW = 172.26 g/mole) and linalool (MW = 154.25
1199 g/mole) are hydrocarbons in a liquid form with similar molecular weights.
- 1200 • Hydroxycitronellal exhibits low peptide reactivity; peptide reactivity
1201 information is not available for linalool.
- 1202 • Hydroxycitronellal and linalool were not strongly positive in the traditional
1203 LLNA (EC3 = 24% and 30% with maximum $SI = 8.5$ and 8.3 , respectively, at
1204 100%).

1205 • Linalool is a skin irritant at the concentrations tested in the LLNA: BrdU-
1206 ELISA and traditional LLNA, but hydroxycitronellal was not.

1207 For the three substances misclassified as sensitizers in the LLNA: BrdU-ELISA (hexane,
1208 lactic acid and propylene glycol), although they represented different chemical classes
1209 (acyclic hydrocarbons, carboxylic acids, and alcohols, respectively) all three:

- 1210 • Are liquids
- 1211 • Have minimal peptide reactivity
- 1212 • Have molecular weights below 100 g/mole
- 1213 • Were tested at concentrations that are irritating to skin.

1214 **Table 6-9 Discordant Results for LLNA: BrdU-ELISA (Using SI \geq 1.5 for**
1215 **Sensitizers) Compared to Traditional LLNA and Guinea Pig Reference**
1216 **Data¹**

Substance Name ²	Vehicle ³	LLNA: BrdU-ELISA ⁴	Traditional LLNA ⁴	Guinea Pig Studies	Skin Irritant?
Hydroxycitronellal (24.0%)	AOO	- (1.3, 100%)	+ (8.5, 100%)	+	Negative at \leq 100%
Linalool (30.0%)	AOO	- (1.45, 100%) ⁵	+ (8.3, 100%)	NA	Mild irritant at 10%
Hexane	AOO	+ (1.8, 100%) ⁶	- (2.2, 100%)	NA	Irritant at 100%
Lactic acid	DMSO	+ (2.5, 50%)	- (2.2, 25%)	-	Slightly irritating at 10%
Propylene glycol	AOO	+ (1.6, 50%)	- (1.6, 100%)	-	Negative at \leq 25%

1217 Abbreviations: AOO = acetone: olive oil (4: 1); LLNA: BrdU-ELISA= murine local lymph node assay with
1218 enzyme-linked immunosorbent assay detection of bromodeoxyuridine; DMSO = dimethyl sulfoxide;; LLNA =
1219 murine local lymph node assay; NA = not available.

1220 + = Sensitizer.

1221 - = Nonsensitizer.

1222 ¹Data sources provided in **Appendix C-1**.

1223 ²Numbers in parentheses are the EC3 values (estimated concentration needed to produce a stimulation index
1224 [SI] of 3) for the traditional LLNA (from **Table 3-1**).

1225 ³Vehicles apply to tests for both the LLNA: BrdU-ELISA and the traditional LLNA.

1226 ⁴Numbers in parentheses are highest SI values and maximum concentrations tested.

1227 ⁵Highest SI occurred at concentration of 50%.

1228 ⁶An additional test yielded an SI of 1.9 at 50%.

1229 When the outcomes for the 29 substances with LLNA: BrdU-ELISA (using SI \geq 1.5) and
1230 human outcome data were compared, the classifications for eight substances were different

1231 (Table 6-10). The LLNA: BrdU-ELISA results for three of these substances
 1232 (hydroxycitronellal, hexane, and lactic acid) were discordant with the traditional LLNA. The
 1233 LLNA: BrdU-ELISA classified four human sensitizers as nonsensitizers (diethyl phthalate,
 1234 2-hydroxypropylmethacrylate, isopropanol, and hydroxycitronellal) and four human
 1235 nonsensitizers as sensitizers (cyclamen aldehyde, isopropyl myristate, hexane, and lactic
 1236 acid).

1237 **Table 6-10 Discordant Results for LLNA: BrdU-ELISA (SI \geq 1.5) When Compared**
 1238 **to Traditional LLNA and Human Outcome Data¹**

Substance Name ²	Vehicle ³	LLNA: BrdU-ELISA ⁴	Traditional LLNA ⁴	Human Outcome ⁵	Skin Irritant?
Diethyl phthalate	AOO	- (0.9, 50%)	- (1.5, 100%)	+ (HPTA)	Negative at \leq 100%
2-Hydroxypropylmethacrylate	AOO	- (1.1, 50%)	- (1.3, 50%)	+ (case study, 0.1%)	Negative at \leq 10%
Isopropanol	AOO	- (2.2, 50%) ⁶	- (1.7, 50%) ⁷	+ (case study, 0.001%)	Negative at \leq 100%
Hydroxycitronellal (24.0%)	AOO	- (1.3, 100%)	+ (8.5, 100%)	+ (14/73, 20%)	Negative at \leq 100%
Cyclamen aldehyde (22.3%)	AOO	+ (1.97, 100%)	+ (5.2, 50%)	- (0/64, 4%)	Irritant at 100%
Isopropyl myristate (44%)	AOO	+ (4.2, 50%)	+ (3.4, 100%)	- (0/25, 20%)	Negative at \leq 100%
Hexane	AOO	+ (1.8, 100%) ⁸	- (2.2, 100%)	- (0/25, 100%)	Irritant at 100%
Lactic acid	DMSO	+ (2.5, 50%)	- (2.2, 25%)	- (no data)	Slightly irritating at 10%

1239 Abbreviations: AOO = acetone: olive oil (4: 1); DMSO = dimethyl sulfoxide; LLNA: BrdU-ELISA = Murine
 1240 local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; HPTA =
 1241 human patch test allergen; LLNA = murine local lymph node assay.

1242 + = Sensitizer.

1243 - = Nonsensitizer.

1244 ¹Data sources provided in **Appendix C-1**.

1245 ²Numbers in parentheses are EC3 values (estimated concentration needed to produce a stimulation index [SI] of
 1246 3) for substances that are sensitizers in the traditional LLNA; from **Table 3-1**.

1247 ³Vehicles apply to tests for both the LLNA: BrdU-ELISA and the traditional LLNA unless otherwise noted.

1248 ⁴Numbers in parentheses are highest SI values and maximum concentrations tested.

1249 ⁵Information in parentheses indicates the basis for the human outcome. Numbers indicate the incidence of
 1250 positive human response and concentration.

1251 ⁶Negative based on most prevalent call. Highest SI of any test is shown. Highest SIs for most tests occurred at <
 1252 50%.

1253 ⁷Highest SI occurred at 10%.

1254 ⁸An additional test yielded SI=1.9 at 50%.

1255 Few commonalities in chemical class, physical form, molecular weight, peptide reactivity,
1256 traditional LLNA EC3 range, or potential for skin irritation were noted among the discordant
1257 substances. The four human sensitizers that were misclassified as nonsensitizers:

- 1258 • Represented three different chemical classes: carboxylic acids (diethyl
1259 phthalate and 2-hydroxypropylmethacrylate), alcohols (isopropanol), and
1260 hydrocarbons (hydroxycitronellal) (**Tables 3-1 and 6-10**).
- 1261 • Three substances were liquids (diethyl phthalate, isopropanol, and
1262 hydroxycitronellal) and one was a solid (2-hydroxypropylmethacrylate).
- 1263 • Molecular weights ranged from 60.10 (isopropanol) to 222.24 g/mole (diethyl
1264 phthalate).
- 1265 • All four substances exhibited low peptide reactivity.
- 1266 • Three were classified as nonsensitizers by the traditional LLNA and one was
1267 classified as a sensitizer (hydroxycitronellal with EC3 = 24.0%).
- 1268 • Although 2-hydroxypropylmethacrylate is a skin irritant at the concentrations
1269 tested in the LLNA: BrdU-ELISA, the other three substances were not
1270 irritating to skin at the concentrations tested (**Table 6-10**).

1271 There were few commonalities in chemical class, physical form, molecular weight, peptide
1272 reactivity, EC3 range (based on the traditional LLNA), or potential for skin irritation noted
1273 among the four human nonsensitizers that were misclassified as sensitizers by the LLNA:
1274 BrdU-ELISA.

- 1275 • The four substances represented three different chemical classes: carboxylic
1276 acids (cyclamen aldehyde and lactic acid), lipids (isopropyl myristate), and
1277 acyclic hydrocarbons (hexane) (**Tables 3-1 and 6-10**).
- 1278 • While all four substances are liquids, with minimal to low peptide reactivity,
1279 molecular weights ranged from 86.15 g/mole for hexane to 270.46 g/mole for
1280 isopropyl myristate.
- 1281 • Cyclamen aldehyde and isopropyl myristate were also classified as sensitizers
1282 by the traditional LLNA (EC3 values were 22.3% and 44%, respectively), but

1283 hexane and lactic acid were classified as nonsensitizers by the traditional
1284 LLNA.

1285 • Two of the substances (cyclamen aldehyde and lactic acid) misclassified as
1286 sensitizers were tested at concentrations that are irritating to skin, but two
1287 were not (isopropyl myristate and hexane) (**Table 6-10**).

1288 **6.7 LLNA: BrdU-ELISA Accuracy Analysis Using Multiple Alternative Decision** 1289 **Criteria**

1290 As detailed in **Section 6.5**, the accuracy of the LLNA: BrdU-ELISA when using a number of
1291 alternative decision criteria was evaluated using the traditional LLNA as the reference test.
1292 The lowest decision criterion with a 0% (0/9) false positive rate was $SI \geq 2.0$, which was
1293 used by the JSAAE interlaboratory validation study. The accuracy at $SI \geq 2.0$ was 84%
1294 (26/31) and the false negative rate was 23% (5/22) (**Table 6-6**). Higher SI values also
1295 produced false positive rates of 0% (0/9), but the false negative rate increased as the SI
1296 increased. The lowest false negative rate was produced at $SI \geq 1.3$ (0% [0/22]), but the false
1297 positive rate at $SI \geq 1.3$ was 44% (4/9).

1298 The 0% false positive rate using $SI \geq 2.0$ and the 0% false negative rate using $SI \geq 1.3$
1299 prompted an evaluation of using two decision criteria for LLNA: BrdU-ELISA results: one
1300 criterion to classify substances as sensitizers (i.e., $SI \geq 2.0$) and one criterion to classify
1301 substances as nonsensitizers (i.e., $SI < 1.3$). The $SI \geq 1.3$ criterion, when used to classify
1302 sensitizers, resulted in no false negative results with respect to the traditional LLNA results.
1303 However, using $SI \leq 1.3$ to classify substances as nonsensitizers resulted in one false
1304 negative result (4% [1/22]), which was for hydroxycitronellal (at 100% the LLNA: BrdU-
1305 ELISA $SI = 1.30$, while the traditional LLNA $SI = 8.5$). Thus, $SI < 1.3$ is proposed to classify
1306 substances as nonsensitizers because this criterion results in no false negative results.

1307 It should be noted that this analyses was based on the same strategy for combining results
1308 from multiply tested substances described in **Section 6.5** (i.e., the sensitizer/nonsensitizer
1309 outcome for each substance was most prevalent outcome). **Section 7.3** details the
1310 reproducibility of multiply tested substances and indicates that, while there were isolated
1311 instances of false positive results for nonsensitizers (i.e., $SI \geq 2.0$), there were no false
1312 negatives. Among the 78 tests that produced a maximum $SI \geq 2.0$, 4% (3/78) were

1313 nonsensitizers (i.e., produced a false positive result). These results were obtained for
1314 isopropanol and lactic acid, which produced SI = 2.0 and 2.2 in two different tests in the
1315 LLNA: BrdU-ELISA and one test of lactic acid, which produced an SI = 2.5. Isopropanol
1316 and lactic acid are classified as nonsensitizers based on maximum SI values of 1.7 and 2.2,
1317 respectively in the traditional LLNA. See **Section 7.3** from more details regarding these
1318 results.

1319 **6.8 Discordant Results for Accuracy Analysis Using Multiple Alternative** 1320 **Decision Criteria**

1321 While optimum false positive and false negative rates can be achieved using these two
1322 different decision criteria, a range of SI values ($1.3 \leq SI < 2.0$) now exists for which the
1323 correct classification is not definitive (i.e., there is a chance for false positives or false
1324 negatives for substances in this range). Chemical class, physical form, MW, peptide
1325 reactivity (see **Appendix B** for physical/chemical properties), traditional LLNA EC3 range,
1326 or potential for skin irritation were examined to identify commonalities among the substances
1327 that produced SI values of 1.3 to < 2.0 in an attempt to identify common characteristics
1328 among these substances that could be used to correctly classify such substances.

1329 Eleven substances produced SI values from 1.3 to < 2.0 (see **Table 6-11**). Five of the 11
1330 substances are nonsensitizers and six are sensitizers based on traditional LLNA results. The
1331 five substances classified by the traditional LLNA as nonsensitizers (methyl salicylate,
1332 isopropanol, propylene glycol, hexane and lactic acid), represented four chemical classes
1333 (carboxylic acids, alcohols, phenols and acyclic hydrocarbons).

- 1334 • Two substances are classified as carboxylic acids (methyl salicylate [also a
1335 phenol] and lactic acid) and two were classified as alcohols (isopropanol, and
1336 propylene glycol).
- 1337 • Hexane is an acyclic hydrocarbon.

1338

1338
1339**Table 6-11 Discordant Results for LLNA: BrdU-ELISA When Multiple Decision Criteria Are Used ¹**

Substance Name ²	Vehicle ³	LLNA: BrdU-ELISA ⁴	Traditional LLNA ⁵	Skin Irritant?
Hexane	AOO	1.76, 100% ⁷ 1.9, 50% (2/2 tests)	- (2.2, 100%)	Irritant at 100%
Isopropanol	AOO	1.6, 50% (1/7 tests)	- (1.7, 50%) ⁹	No, up to 100%
Lactic acid	DMSO	1.8, 50% 1.9, 50% (2/2 tests)	- (2.2, 25%)	Slightly irritating at 10%
Methyl salicylate	AOO	1.4, 50% (3/3 tests at SI = 1.4)	- (2.9, 20%)	Irritant at 10%
Propylene glycol	AOO	1.6, 50% (1/2 tests)	- (1.6, 100%)	No, up to 25%
Aniline (47.5%)	AOO	1.5, 50%	+ (3.6, 100%) ⁷	No, up to 100%
Hydroxycitronellal (24.0%)	AOO	1.30, 100%	+ 8.5, 100%)	No, up to 50%
Linalool (30.0%)	AOO	1.45, 100% ⁷	+ (8.3, 100%)	Mild irritant at 100%
2-Mercaptobenzo-thiazole (1.7%)	DMF	1.6, 50% ¹⁰	+ (8.6, 10%)	No, up to 10%
Cyclamen aldehyde (22.3%)	AOO	1.97, 100%	+ 5.2, 50%	Irritant at 100%
Formaldehyde (0.50%)	ACE	1.97, 10% (1/3 tests)	+ (4.0, 1.8%)	No, up to 2%

1340 Abbreviations: ACE = acetone; AOO = acetone: olive oil (4: 1); DMF = *N,N*-Dimethylformamide; DMSO =
1341 dimethyl sulfoxide; LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent
1342 assay detection of bromodeoxyuridine; HPTA = human patch test allergen; LLNA = murine local lymph node
1343 assay; NA = not available; + = Sensitizer; - = Nonsensitizer.

1344 ¹Data sources provided in **Appendix C-1**.

1345 ²Numbers in parentheses are EC3 values (estimated concentration needed to produce a stimulation index [SI] of
1346 3) for substances that are sensitizers in the traditional LLNA; from **Table 3-1**.

1347 ³Vehicles apply to tests for both the LLNA: BrdU-ELISA and the traditional LLNA unless otherwise noted.

1348 ⁴Numbers are highest SI values achieved and maximum concentration tested.

1349 ⁵Information in parentheses indicates the basis for the human outcome. Numbers indicate the incidence of
1350 positive human response and concentration tested.

1351 ⁷Highest SI occurred at 50%.

1352 ⁸The solvent for the traditional LLNA was *N,N*-dimethylformamide.

1353 ⁹Highest SI occurred at 10%.

1354 ¹⁰Highest SI occurred at 12.5%.

1355

1356 Other characteristics of the nonsensitizers (by the traditional LLNA) include:

- 1357
- All of the five substances are liquids and have minimal peptide reactivity.
- 1358
- Four substances have MW < 100 g/mole. The other substances, methyl
- 1359 salicylate, has a MW of 152.15 g/mole, respectively.
- 1360
- Four of the five substances were tested at irritating concentrations in both the
- 1361 traditional LLNA and in the LLNA: BrdU-ELISA: methyl salicylate,
- 1362 propylene glycol, hexane and lactic acid. Isopropanol was tested at
- 1363 concentrations nonirritating to skin.
- 1364
- Two of the five substances yielded SI < 2 in the traditional LLNA:
- 1365 isopropanol and propylene glycol. The other three substances yielded SI
- 1366 values between 2 and 3 (exclusive): hexane, lactic acid and methyl salicylate.

1367 The six substances classified by the traditional LLNA as sensitizers (aniline,
1368 hydroxycitronellal, linalool, 2-mercaptobenzothiazole, formaldehyde, and cyclamen
1369 aldehyde) represent five chemical classes. Aniline is an amine, hydroxycitronellal and
1370 linalool are hydrocarbons (other), 2-mercaptobenzothiazole is a heterocyclic compound,
1371 formaldehyde is an aldehyde, and cyclamen aldehyde is a carboxylic acid. Other
1372 characteristics of the discordant substances that are classified as sensitizers by the traditional
1373 LLNA include:

- 1374
- Five are liquids and one is a solid (2-mercaptobenzothiazole).
- 1375
- Three substances have MW between 150 and 200 g/mole. Formaldehyde and
- 1376 aniline both have MW less than 100 g/mole (MW =30 g/mole and 93.13
- 1377 g/mole, respectively).
- 1378
- Hydroxycitronellal and cyclamen aldehyde exhibit low peptide reactivity,
- 1379 formaldehyde exhibits moderate peptide reactivity, 2-mercaptobenzothiazole
- 1380 exhibits high peptide reactivity, and peptide reactivity information is not
- 1381 available for the other two substances.

- 1382
- 1383
- 1384
- 1385
- 1386
- 1387
- 1388
- 1389
- 1390
- 1391
- Aniline, linalool, hydroxycitronellal, and cyclamen aldehyde were not strongly positive in the traditional LLNA (EC3 = 47.5%, 30%, 24%, and 22.3%) respectively, with maximum SI = 3.6, 8.3, 8.5, and 5.2, respectively, when tested at concentrations up to 100%.
 - Hydroxycitronellal, linalool, and 2-mercaptobenzothiazole were tested at concentrations in the LLNA: BrdU-ELISA and traditional LLNA that were irritating to skin, but aniline was not. Formaldehyde and cyclamen aldehyde were tested at concentrations in the LLNA: BrdU-ELISA that were irritating to skin, but was not tested at irritating concentrations in the traditional LLNA.

1391 **7.0 Test Method Reliability**

1392 An assessment of test method reliability (intra- and inter-laboratory reproducibility) is an
1393 essential element of any evaluation of the performance of an alternative test method
1394 (ICCVAM 2003). Intralaboratory reproducibility refers to the extent to which qualified
1395 personnel within the same laboratory can replicate results using a specific test protocol at
1396 different times. Interlaboratory reproducibility refers to the extent to which different
1397 laboratories can replicate results using the same protocol and test substances, and indicates
1398 the extent to which a test method can be transferred successfully among laboratories.

1399 The reproducibility evaluation in this draft BRD has been revised from the January 2008
1400 draft BRD to include the results for a number of additional intralaboratory tests for which SI
1401 values were newly available. The interlaboratory reproducibility evaluation is a new addition
1402 to this draft BRD because interlaboratory data were not available for evaluation in the
1403 January 2008 draft BRD. This draft BRD also includes a reproducibility analysis using
1404 separate SI criteria to identify sensitizers and nonsensitizers.

1405 The available LLNA: BrdU-ELISA data were amenable to both intralaboratory and
1406 interlaboratory reproducibility analyses. This section provides an assessment of
1407 reproducibility for the decision criterion of $SI \geq 2.0$ to identify sensitizers. The evaluation of
1408 single decision criteria in **Section 6.6** showed that this was the lowest SI value that produced
1409 a false positive rate of 0% (0/9) when the traditional LLNA was the reference test (**Table 6-**
1410 **6**). The evaluation of multiple decision criteria in **Section 6.7** evaluated $SI \geq 2.0$ as the
1411 decision criterion for classifying substances as sensitizers when used with a decision criterion
1412 of <1.3 to identify nonsensitizers. In addition, the decision criterion of $SI \geq 2.0$ to identify
1413 sensitizers was used in the JSAAE interlaboratory validation study. **Appendix F** describes
1414 the evaluation of reproducibility for the decision criterion of $SI \geq 1.5$ to identify sensitizers,
1415 which was evaluated in **Section 6.6**.

1416 **7.1 Intralaboratory Reproducibility**

1417 The test results for the LLNA: BrdU-ELISA were amenable to intralaboratory reproducibility
1418 analyses for three endpoints: sensitizer or nonsensitizer classification, SI values, and EC2
1419 values. Analyses of intralaboratory reproducibility were performed using a concordance
1420 analysis for the qualitative results (sensitizer vs. nonsensitizer) (**Section 7.1.1**) and a

1421 coefficient of variation (CV) analysis for the quantitative results (SI values and EC3 values)
1422 (**Sections 7.1.2** and **7.1.3**, respectively).

1423 7.1.1 *Intralaboratory Reproducibility – Qualitative Results*

1424 The dataset available for an intralaboratory concordance analysis of the qualitative test
1425 results for the LLNA: BrdU-ELISA included eight substances that were tested multiple times
1426 and classified as sensitizers or nonsensitizers. Hexyl cinnamic aldehyde was tested six times,
1427 eugenol was tested five times, and isoeugenol was tested three times, and 2,4-
1428 dinitrochlorobenzene, glutaraldehyde, hexane, 4-phenylenediamine, and propylene glycol
1429 were each tested twice (Takeyoshi et al. 2003, 2004a, 2005, 2006, 2007a; unpublished data)
1430 (**Table 7-1**). All substances were sensitizers in the traditional LLNA except for propylene
1431 glycol and hexane. The multiple test results for 8/8 substances were 100% concordant when
1432 $SI \geq 2.0$ was used to classify substances as sensitizers.

1433 By comparison, the qualitative intralaboratory concordance analysis for the traditional LLNA
1434 (ICCVAM 1999) was based on a dataset of six substances that included six results each for
1435 benzocaine and hexyl cinnamic aldehyde, five results for eugenol, four results each for
1436 isoeugenol and methyl salicylate, and three results for 2,4-dinitrochlorobenzene.
1437 Intralaboratory results for each substance were 100% concordant with the exception of
1438 benzocaine. One of the six benzocaine (5/6 or 83% concordance) results for the traditional
1439 LLNA was reported as equivocal because SI increased with dose, but did not reach the
1440 criterion of $SI \geq 3.0$. Thus, the proportion of substances for which intralaboratory
1441 concordance of qualitative results was 100% was similar for LLNA: BrdU-ELISA (7/8) and
1442 the traditional LLNA (5/6).

1443

1443 **Table 7-1 Intralaboratory Reproducibility for the LLNA: BrdU-ELISA Outcome of**
 1444 **Substances Tested Multiple Times**

Substance	Highest Concentration Tested (%)	Highest SI	Outcome ¹	Takeyoshi et al. Reference
2,4-Dinitro-chlorobenzene	2	17.9	+	2005
	2	6.8	+	2006, 2007b
Eugenol	30	3.3	+	2004a
	30	3.8	+	2007a
	50	12.3	+	2005
	50	3.1	+	2006
	50	17.7	+	2007b
Glutaraldehyde	2	14.6	+	2005, 2007b
	10	15.5	+	2005, 2007b
Hexane	50	1.9	-	2005
	100	1.8	-	unpublished data
Hexyl cinnamic aldehyde	25	2.4	+	2003
	50	3.6	+	2003
	50	5.9	+	2005
	50	3.6	+	2006
	50	2.7	+	2006
	50	3.0	+	2007b
Isoeugenol	10	8.4	+	2005
	10	2.4	+	2006, 2007b
	30	6.7	+	
4-Phenylenediamine	2	11.7	+	2005, 2007b
	10	14.7	+	2005, 2007b
Propylene glycol	50	1.6	-	2005
	50	0.9	-	2006, 2007b

1445 Abbreviations: LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent
 1446 assay detection of bromodeoxyuridine; SI = stimulation index.

1447 ¹(+) = Sensitizer; (-) = nonsensitizer

1448

1449 7.1.2 *Intralaboratory Reproducibility – SI*

1450 There were six substances that were tested multiple times by Takeyoshi et al. (2003, 2004a,
 1451 2005, 2006, 2007a, 2007b, unpublished data). Because two substances had multiple tests for
 1452 more than one concentration, there were nine substance/concentration combinations that
 1453 were tested two to five times in separate experiments. The multiple SI values for each
 1454 substance/concentration were used to calculate a CV for the assessment of intralaboratory
 1455 variability. As shown by **Table 7-2**, the CVs ranged from 1% (25% hexyl cinnamic
 1456 aldehyde) to 79% (10% isoeugenol). The intralaboratory reproducibility of the traditional
 1457 LLNA was not assessed by CV analysis of SI values (ICCVAM 1999).

1458 7.1.3 *Intralaboratory Reproducibility – EC2*

1459 CV values were also calculated for the EC2 values for the three sensitizers that were tested
 1460 more than once using multiple doses by Takeyoshi et al. (2003; 2004a, 2005, 2006, 2007a,
 1461 2007b). The individual animal data for eugenol, hexyl cinnamic aldehyde, and isoeugenol,
 1462 were used to calculate EC2 values for the LLNA: BrdU-ELISA. The methods for calculating
 1463 EC2 values for each sensitizer were modified from those used by Ryan et al. (2007) to
 1464 calculate EC3 values. Linear interpolation was used to calculate EC2 values for each test
 1465 with SI values higher or lower than two and extrapolation was used to calculate EC2 values
 1466 for tests with no SI values below two. The equation for linear interpolation was:

$$1467 \quad EC2 = c + \left[\frac{(2-d)}{(b-d)} \right] \times (a-c)$$

1468 The linear interpolation equation uses the points immediately above and below SI = 2, with
 1469 the (dose, SI) coordinates of (a, b) immediately above SI = 2 and (c, d) immediately below SI
 1470 = 2. The equation for extrapolation was:

$$1471 \quad EC2_{ex} = 2 \left\{ \log_2(c) + \frac{(2-d)}{(b-d)} \times [\log_2(a) - \log_2(c)] \right\}$$

1472

1472 **Table 7-2 Intralaboratory Reproducibility for the SI of Tested Substances in**
 1473 **LLNA: BrdU-ELISA - Coefficient of Variation**

Substance	Concentration Tested (%)	SI	Mean	SD	CV (%)	Takeyoshi et al. Reference
2,4-Dinitrochlorobenzene	2	17.9	12.4	7.8	64	2005
	2	6.8				2006, 2007b
Eugenol	30	3.3	3.6	0.4	10	2004a
	30	3.8				2007a
Eugenol	50	12.3	11.0	7.4	67	2005
	50	3.1				2006
	50	17.7				2007b
Hexane	50	1.9	1.6	0.4	22	2005
	50	1.4				unpublished
Hexyl cinnamic aldehyde	12.5	1.87	1.73	0.21	12	2003
	12.5	1.58				2003
Hexyl cinnamic aldehyde	25	2.42	2.4	0.01	1	2003
	25	2.40				2003
Hexyl cinnamic aldehyde	50	3.6	3.8	1.3	34	2003
	50	5.9				2005
	50	3.6				2006
	50	2.7				2006
	50	3.0				2007b
Isoeugenol	10	8.4	5.4	4.2	79	2005
	10	2.4				2006, 2007b
Propylene glycol	50	1.6	1.1	0.6	55	2005
	50	0.7				2006, 2007b

1474 Abbreviations: CV = coefficient of variation; LLNA: BrdU-ELISA = murine local lymph node assay with
 1475 enzyme-linked immunosorbent assay detection of bromodeoxyuridine; SD = standard deviation, SI =
 1476 stimulation index.
 1477

1477 The extrapolation equation uses the two points immediately above SI = 2, with the
 1478 coordinates of (a, b) for the point closest to SI = 2, and (c, d) for the higher point. As shown
 1479 in **Table 7-3**, there were five EC2 values for hexyl cinnamic aldehyde, four EC2 values for
 1480 eugenol, and two EC2 values for isoeugenol. The CV values were 73% for eugenol, 25% for
 1481 hexyl cinnamic aldehyde, and 16% for isoeugenol. The ICCVAM LLNA *Performance*
 1482 *Standards* criteria for demonstrating adequate intralaboratory reproducibility is based on
 1483 results from at least four independent tests of hexyl cinnamic aldehyde (ICCVAM 2009).
 1484 Intralaboratory reproducibility is considered adequate when each test yields an EC_t value
 1485 (i.e., the estimated concentration needed to produce an SI of a specific threshold value; in this
 1486 case, SI = 1.5) within 5% to 20% (ICCVAM 2009). Two of the five EC2 values for hexyl
 1487 cinnamic aldehyde were within the acceptable range for intralaboratory reproducibility.

1488 **Table 7-3 Intralaboratory Reproducibility for the EC2 of Tested Substances in**
 1489 **LLNA: BrdU-ELISA - Coefficient of Variation**

Substance	EC2	Mean	SD	CV (%)	Takeyoshi et al. Reference
Eugenol	11.2	12.6	9.2	73	2004a
	23.6				2006
	1.2				2007b
	14.6				2007a
Hexyl cinnamic aldehyde	15.2	22.6	5.7	25	2003
	18.8				2003
	29.9				2006
	25.5				2006
	23.4				2007b
Isoeugenol	8.4	7.6	1.2	16	2006; 2007b
	6.7				2007a

1490 Abbreviations: CV = coefficient of variation; EC2 = estimated concentration needed to produce a stimulation
 1491 index of two; LLNA: BrdU-ELISA = murine local lymph node assay with enzyme-linked immunosorbent assay
 1492 detection of bromodeoxyuridine ; SD = standard deviation.

1493 The intralaboratory reproducibility of the traditional LLNA was assessed by CV analysis of
 1494 EC3 values using a larger dataset (ICCVAM 1999) than that available for the LLNA: BrdU-
 1495 ELISA analysis. Two EC3 values were reported by each of five laboratories for
 1496 2, 4-dinitrochlorobenzene, five EC3 values were reported by one laboratory for isoeugenol,
 1497 six EC3 values were reported for hexyl cinnamic aldehyde by two laboratories, and five EC3
 1498 values were reported for eugenol by one laboratory (**Table 7-4**).

1499 **Table 7-4 Intralaboratory Reproducibility for the EC3 of Tested Substances in the**
 1500 **Traditional LLNA¹**

Substance	Number of Laboratories	Number of Tests per Laboratory	CV (%)
2, 4-Dinitrochlorobenzene	5	2	13-47
Isoeugenol	1	5	26
Hexyl cinnamic aldehyde	2	6	19-27
Eugenol	1	5	18

1501 Abbreviations: CV = coefficient of variation; EC3 = estimated concentration needed to produce a stimulation
 1502 index of three; LLNA = murine local lymph node assay;
 1503 ¹From ICCVAM (1999).

1504 For two of three substances, the intralaboratory CV values for the EC2 values from LLNA:
 1505 BrdU-ELISA tests were higher than EC3 values for the same substances from the traditional
 1506 LLNA reported in ICCVAM (1999). The intralaboratory EC2 CV from the LLNA: BrdU-
 1507 ELISA tests of eugenol was higher than that reported by ICCVAM (1999) (73% vs. 18%).
 1508 The intralaboratory EC2 CV from the LLNA: BrdU-ELISA tests of isoeugenol was greater
 1509 than that from ICCVAM (1999) (16% vs. 26%). However, the intralaboratory EC2 CV from
 1510 the LLNA: BrdU-ELISA tests of hexyl cinnamic aldehyde was within the range reported by
 1511 ICCVAM (1999) (25% vs. 19 to 27%).

1512 7.2 Interlaboratory Reproducibility

1513 The interlaboratory reproducibility of the LLNA: BrdU-ELISA was assessed using the
 1514 individual animal data from the multi-laboratory validation study organized by the JSAAE
 1515 (Kojima et al. 2008). Phase I of the study evaluated the reliability and transferability of the
 1516 test method protocol by testing 12 substances in three to nine laboratories. With the
 1517 exception of the positive control data, neither the summary results nor the individual animal
 1518 data from Phase I of the validation study have been released. Phase II of the study tested 10
 1519 substances in three to seven laboratories as shown in **Table 7-5**. All the laboratories that
 1520 participated in the validation study used the same experimental protocol (**Appendix A**) and
 1521 participated in a one-day seminar that explained the protocol and execution of the test
 1522 method. The same commercial ELISA kit, test materials, and the same doses of the test
 1523 substances were used in all of the laboratories. The Validation Management Team

1524 determined the doses and vehicles for testing and coded the identity of the test substances
 1525 prior to distribution to the test laboratories. Seven substances were sensitizers and three
 1526 substances were nonsensitizers according to the traditional LLNA. Six substances were
 1527 ICCVAM *Recommended Performance Standards* reference substances: 2,4-
 1528 dinitrochlorobenzene, eugenol, hexyl cinnamic aldehyde, lactic acid, isopropanol, and methyl
 1529 salicylate (ICCVAM 2009).

1530 **Table 7-5 Substances and Test Allocation for the Phase II Interlaboratory**
 1531 **Validation Study of the LLNA: BrdU-ELISA**
 1532

Substance ¹	Vehicle	Concentrations Tested			Laboratory ²						
					1	2	3	4	5	6	7
Nickel sulfate (+)	DMSO	1%	3%	10%			X	X			X
Isopropanol (-)	AOO	10%	25%	50%	X	X	X	X	X	X	X
Eugenol (+)	AOO	10%	25%	50%		X				X	X
Cinnamic aldehyde (+)	AOO	1%	3%	10%		X		X	X		
2,4-Dinitrochlorobenzene (+)	AOO	0.1%	0.3%	1%	X	X	X	X	X	X	X
Glutaraldehyde (+)	ACE	0.1%	0.3%	1%	X				X	X	
Methyl salicylate (-)	AOO	10%	25%	50%	X	X	X				
Hexyl cinnamic aldehyde (+)	AOO	10%	25%	50%	X	X	X	X	X	X	X
Lactic acid (-)	DMSO	10%	25%	50%			X	X			X
Formaldehyde (+)	ACE	1%	3%	10%	X				X	X	

1533 Abbreviations: ACE = acetone; AOO = acetone: olive oil (4: 1); DMSO = dimethyl sulfoxide; LLNA: BrdU-ELISA =
 1534 murine local lymph node assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine

1535 ¹(+) indicates sensitizers and (-) indicates nonsensitizers according to traditional LLNA tests.

1536 ²X indicates that a substance was tested in a particular laboratory. 1 = Daicel Chemical Industries Ltd.; 2 = Food and Drug
 1537 Safety Center; 3 = Otsuka Pharmaceutical Co. Ltd.; 4 = Taisho Pharmaceutical Co. Ltd.; 5 = Fuji Film Co. Ltd.; 6 =
 1538 Biosafety Research Center, Foods, Drugs and Pesticides; 7 = National Institute of Health Sciences.

1539 The LLNA: BrdU-ELISA test results from the JSAAE validation study were used for
 1540 interlaboratory reproducibility analyses for three endpoints: sensitizer or nonsensitizer
 1541 classification and EC2 values. Analyses of interlaboratory reproducibility were performed
 1542 using a concordance analysis for the qualitative results (sensitizer vs. nonsensitizer) (**Section**
 1543 **7.2.1**) and a CV analysis for the quantitative results (EC2 values) (**Sections 7.2.2 and 7.2.3,**
 1544 respectively).

1545

1545 7.2.1 *Interlaboratory Reproducibility – Qualitative Results*

1546 The available quantitative absorbance data for interlaboratory reproducibility analysis were
 1547 used to calculate SI values for each substance and dose tested. Substances with $SI \geq 2.0$ at
 1548 any dose were classified as sensitizers. The qualitative (sensitizer/nonsensitizer)
 1549 interlaboratory concordance analysis for the 10 substances tested during Phase II of the
 1550 JSAAE interlaboratory validation study is shown in **Table 7-6**. The qualitative comparison
 1551 of LLNA: BrdU-ELISA results (i.e., positive vs. negative) for 10 substances tested among up
 1552 to 7 laboratories were consistent. The concordance results show that interlaboratory
 1553 concordance was 100% (3/3, 6/6, or 7/7) for seven substances. There were three discordant
 1554 substances (formaldehyde, isopropanol, and lactic acid) for which interlaboratory
 1555 concordance was 67% (2/3 or 4/6). One of the three laboratories reported an SI of 1.97 for
 1556 formaldehyde; while the others produced $SI > 2$. Two of the six tests of isopropanol yielded
 1557 $SI \geq 2.0$ ($SI = 2.0$ and $SI = 2.2$); while the others yielded negative results. One of the three
 1558 tests for lactic acid produced $SI \geq 2.0$ (i.e., $SI = 2.5$), while the others yielded $SI < 2.0$. The
 1559 Validation Management Team considered the interlaboratory reproducibility to be acceptable
 1560 (Kojima et al. 2008). Because the evaluation of interlaboratory reproducibility for the
 1561 traditional LLNA did not include an evaluation of qualitative results (ICCVAM 1999), there
 1562 were no traditional concordance data for comparison with the BrdU-ELISA concordance.

1563 **Table 7-6 Qualitative Results for the Phase II Interlaboratory Validation Study on**
 1564 **the LLNA: BrdU-ELISA¹**

Substance	Laboratory							Concordance
	1	2	3	4	5	6	7	
2,4-Dinitrochlorobenzene	+	+	+	+	+	+	+	7/7
Glutaraldehyde	+				+	+		3/3
Nickel sulfate			+	+			+	3/3
<i>trans</i> -Cinnamic aldehyde		+		+	+			3/3
Formaldehyde	+				+	- ⁴		2/3
Eugenol		+				+	+	3/3
Hexyl cinnamic aldehyde	+	- ³	+	+	+	+ ⁵	+	6/6
Isopropanol	+ ²	- ³	-	-	-	+ ^{2,6}	-	4/6
Lactic acid			-	-			+	2/3
Methyl salicylate	-	-	-					3/3

1565 Abbreviation: LLNA: BrdU-ELISA = Murine local lymph node assay with enzyme-linked immunosorbent assay detection
 1566 of bromodeoxyuridine.

1567 ¹(+) indicates sensitizer result; (-) indicates nonsensitizer result.

1568 ²Stimulation index [SI] ≥ 2 at lowest dose tested, but < 2 at the higher doses. The Validation Management Team considered

1569 these to be nonsensitizer results (Kojima et al. 2008).

1570 ³Test failed because concurrent positive control failed (i.e., SI < 2). Result not included in the concordance analysis.

1571 ⁴Maximum SI = 1.97.

1572 ⁵Three mice tested at highest dose.

1573 ⁶Three mice per dose group.

1574 7.2.2 *Interlaboratory Reproducibility – EC2 Values*

1575 The SI values from the interlaboratory validation study were used to calculate EC2 values for
1576 each sensitizer according to the methods reported in Section 7.1.3. The EC2 values from each
1577 laboratory were then used to calculate CV values for each substance. The resulting values are
1578 shown in **Table 7-7**. CV values ranged from 20% (formaldehyde) to 101% (glutaraldehyde).
1579 The mean CV was 58%.

1580 The ICCVAM LLNA performance standards indicates that interlaboratory reproducibility
1581 should be evaluated with at least two sensitizing chemicals with well-characterized activity in
1582 the traditional LLNA (ICCVAM 2009). Acceptable reproducibility is attained when each
1583 laboratory obtains EC_t values within 0.025% to 0.1% for 2,4-dinitrochlorobenzene and
1584 within 5% to 20% for hexyl cinnamic aldehyde (ICCVAM 2009). EC2 values from two
1585 laboratories were outside these ranges for both substances. Laboratory 2 and Laboratory 5
1586 reported EC2 values that were lower than the specified acceptance range for 2,4-
1587 dinitrochlorobenzene (0.019% and 0.0025%, respectively). For hexyl cinnamic aldehyde,
1588 Laboratory 3 obtained an EC2 value of 24.0%, which was higher than the acceptance range
1589 and Laboratory 5 obtained an EC2 value of 4.07%, which was lower than the acceptance
1590 range.

1591

1592 **Table 7-7 EC2 Values from the Phase II Interlaboratory Validation Study on the LLNA: BrdU-ELISA¹**

Substance	Laboratory							Mean	% CV
	1	2	3	4	5	6	7		
2,4-Dinitro-chlorobenzene	0.084 (4.3 @ 1%)	0.019 (8.37 @ 1%)	0.029 (5.99 @ 0.3%)	0.030 (5.50 @ 1%)	0.0025 (18.80 @ 0.3%)	0.025 (4.83 @ 0.3%)	0.053 (12.18 @ 1%)	0.035	76
Hexyl cinnamic aldehyde	16.2 (3.4 @ 50%)	- ¹ (1.83 @ 50%)	24.0 (2.87 @ 50%)	9.36 (3.34 @ 50%)	4.07 (13.5 @ 50%)	13.0 ² (3.27 @ 50%)	14.2 (3.84 @ 50%)	13.5	50
Glutaraldehyde	0.18	NT	NT	NT	0.034	0.51	NT	0.24	101
Nickel sulfate	NT	NT	3.85	0.95	NT	NT	1.31	2.0	78
<i>trans</i> -Cinnamic aldehyde	NT	2.59	NT	1.63	2.79	NT	NT	2.3	27
Formaldehyde	0.41	NT	NT	NT	0.31	- ³	NT	0.36	20
Eugenol	NT	19.1	NT	NT	NT	16.4	5.06	13.5	55

1593 Note: Bolded font indicates substances recommended for assessing interlaboratory reproducibility in *Recommended Performance Standards* (ICCVAM 2009). Bolded EC2
1594 values are outside of the acceptable range from the ICCVAM LLNA performance standards: 5 - 20% for hexyl cinnamic aldehyde and 0.025 - 0.1% for 2,4-
1595 dinitrochlorobenzene. Values in parentheses are the highest SI values achieved.

1596 Abbreviations: CV = coefficient of variation; EC2 = estimated concentration needed to produce a stimulation index of two; LLNA: BrdU-ELISA = murine local lymph node
1597 assay with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; NT = not tested; SI = stimulation index.

1598 ¹Test failed because associated positive control failed (i.e., SI < 2; vehicle control absorbance was unusually high). Result not included in the mean EC2 and CV.

1599 ²Three mice tested at highest dose.

1600 ³Maximum SI = 1.97.

1601

1602 The interlaboratory CV values for the LLNA: BrdU-ELISA EC2 values were higher than that
 1603 for the traditional LLNA EC3 values. The analysis of interlaboratory variation of EC3 values
 1604 for the traditional LLNA reported CV values of 7 to 84% for five substances tested in five
 1605 laboratories (**Table 7-8**; ICCVAM 1999). Three of the same substances were evaluated in the
 1606 traditional LLNA and the LLNA: BrdU-ELISA. All interlaboratory CV values for LLNA:
 1607 BrdU-ELISA were greater than that for the traditional LLNA. The CV of 76% for 2,4-
 1608 dinitrochlorobenzene was greater than the two CV values of 37% and 27%, calculated from
 1609 five values each, reported by ICCVAM (1999). The CV of 50% for hexyl cinnamic aldehyde
 1610 tested in the LLNA: BrdU-ELISA was greater than the 7% reported by ICCVAM (1999).
 1611 The CV of 55% for eugenol tested in the LLNA: BrdU-ELISA was greater than the 42%
 1612 reported by ICCVAM (1999).

1613 **Table 7-8 Interlaboratory Reproducibility of the EC3 for Substances Tested in the**
 1614 **Traditional LLNA¹**

Substance	Laboratory					CV (%)
	1	2	3	4	5	
2, 4-Dinitrochlorobenzene	0.3	0.5	0.6	0.9	0.6	37
	0.5	0.6	0.4	0.6	0.3	27
Hexyl cinnamic aldehyde	7.9	7.6	8.4	7.0	8.1	7
Isoeugenol	1.3	3.3	1.8	3.1	1.6	41
Eugenol	5.8	14.5	8.9	13.8	6.0	42
Sodium lauryl sulfate	13.4	4.4	1.5	17.1	4.0	84

1615 Abbreviations: CV = coefficient of variation; EC3 = estimated concentration needed to produce a stimulation
 1616 index of three; LLNA = murine local lymph node assay.

1617 ¹From ICCVAM (1999).

1618 **7.3 Reproducibility for the LLNA: BrdU-ELISA Using Multiple Alternative** 1619 **Decision Criteria**

1620 **Section 6.7** details the accuracy analysis for the LLNA: BrdU-ELISA (using the most
 1621 prevalent outcome for substances with multiple tests) when using two decision criteria for
 1622 LLNA: BrdU-ELISA results: one criterion to classify substances as sensitizers ($SI \geq 2.0$) and
 1623 one criterion to classify substances as nonsensitizers (i.e., $SI < 1.3$). $SI \geq 2.0$ was evaluated
 1624 for classifying sensitizers because it resulted in no false positives with respect to the

1625 traditional LLNA. $SI < 1.3$ was evaluated for classifying substances as nonsensitizers
1626 because it resulted in no false negatives. This section evaluates reproducibility of the
1627 concordance with the traditional LLNA results by examining the frequency with which SI
1628 values in the validation database of 31 substances occurred in one of three SI categories. The
1629 three SI categories were:

- 1630 • $SI < 1.3$ for classifying nonsensitizers
- 1631 • $SI \geq 1.3$ and < 2.0 , the range of uncertainty with respect to classification by the
1632 traditional LLNA
- 1633 • $SI \geq 2.0$ to classify substances as sensitizers

1634 The validation database for the LLNA: BrdU-ELISA consists of 102 tests of 31 substances.
1635 The maximum SI achieved by each test and the traditional LLNA outcome (sensitizer vs.
1636 nonsensitizer) were used to determine the frequency of the maximum SI. **Table 7-9** shows
1637 the proportion of sensitizers and nonsensitizers, according to the traditional LLNA for each
1638 SI category. All of the tests (9/9 [100%]) that yielded $SI < 1.3$ were for substances that were
1639 classified as nonsensitizers by the traditional LLNA. Forty percent (6/15) of the tests that
1640 yielded SI values of $1.3 \leq SI < 2.0$ were for substances that were classified as sensitizers by
1641 the traditional LLNA. Three tests produced SI values at either end of this range (i.e., $SI = 1.3$
1642 or $SI = 2.0$). Hydroxycitronellal produced $SI = 1.3$ and the cyclamen aldehyde test and one
1643 formaldehyde test produced $SI = 1.97$. The remainder of the tests in this category, 60%
1644 (9/15), were classified as nonsensitizers by the traditional LLNA. Ninety-six percent (75/78)
1645 of the tests that yielded $SI \geq 2.0$ were for substances that were classified as sensitizers by the
1646 traditional LLNA and only 4% (3/78) were classified as nonsensitizers. The three
1647 nonsensitizer tests were two tests of isopropanol, which yielded $SI = 2.0$ and 2.2 in the
1648 LLNA: BrdU-ELISA, and one test of lactic acid, which produced an $SI = 2.5$.

1649

1649 **Table 7-9 Frequency of Maximum SI for LLNA: BrdU-ELISA Tests by Category**
 1650 **and Traditional LLNA Outcome**

Classification Based on Traditional LLNA	Classification Concordance with Traditional LLNA ¹			Total
	Maximum SI < 1.3	1.3 ≤ Maximum SI < 2.0	Maximum SI ≥ 2.0	
Sensitizer	0 (0%)	6 (40%)	75 (96%)	86
Nonsensitizer	9 (100%)	9 (60%)	3 (4%)	16
Total	9	15	78	102

1651 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay
 1652 with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; SI = stimulation index.

1653 ¹Numbers shown reflect number of tests. Includes all tests of substances that were tested multiple times.

1654 Percentage in parentheses reflects percentage of the total number of tests for each SI category.

1655 The 102 tests evaluated in **Table 7-9** include multiple tests for 14 substances. For the 14
 1656 substances, two to 31 tests were available. **Table 7-10** shows the proportion of the tests for
 1657 each substance that produced SI values in each category. For the nine sensitizers with
 1658 multiple test results, there were no tests that produced SI < 1.3 and only one test that
 1659 produced an SI of 1.3 to <2. This was a formaldehyde test that produced SI = 1.97. For the
 1660 five nonsensitizers with multiple test results, however, SI values occurred in all three SI
 1661 categories. The results for isopropanol were particularly variable: 57% (4/7) produced SI <
 1662 1.3 (two tests with SI= 0.9 and two tests with SI = 1.0), 14% (1/7) produced 1.3 ≤ SI < 2 (SI
 1663 = 1.6), and 29% (2/7) produced SI ≥ 2 (SI = 2.0 and 2.2). Lactic acid tests produced SI values
 1664 in two categories: 67% (2/3) of the tests had 1.3 ≤ SI < 2 (SI = 1.8 and 1.9), and 33% (1/3) of
 1665 the tests had SI ≥ 2 (SI = 2.5). Propylene glycol tests produced SI values in two categories:
 1666 50% (1/2) of the tests had SI < 1.3 (0.9) and one test produced 1.3 ≤ SI < 2 (SI = 1.9). The
 1667 multiple test results for hexane and methyl salicylate were 100% concordant. The two hexane
 1668 tests produced SI values in the 1.3 ≤ SI < 2 category (SI = 1.76 and 1.9) and the three methyl
 1669 salicylate tests also produced SI values in that category (all three SI = 1.4).

1670

1670 **Table 7-10 Concordance of LLNA: BrdU-ELISA Tests for Substances with Multiple**
 1671 **Tests by Maximum SI Category**

Substance	Concordance Among Multiple Tests ¹			Total
	Maximum SI < 1.3	1.3 ≤ Maximum SI < 2.0	Maximum SI ≥ 2.0	
Sensitizers²				
2,4-Dinitrochlorobenzene	0 (0%)	0 (0%)	9 (100%)	9
Eugenol	0 (0%)	0 (0%)	8 (100%)	8
Formaldehyde	0 (0%)	1 (33%)	2 (67%)	3
Glutaraldehyde	0 (0%)	0 (0%)	5 (100%)	5
Hexyl cinnamic aldehyde	0 (0%)	0 (0%)	31 (100%)	31
Isoeugenol	0 (0%)	0 (0%)	3 (100%)	3
Nickel sulfate	0 (0%)	0 (0%)	3 (100%)	3
1,4-Phenylenediamine	0 (0%)	0 (0%)	2 (100%)	2
<i>trans</i> -Cinnamaldehyde	0 (0%)	0 (0%)	4 (100%)	4
Nonsensitizers²				
Hexane	0 (0%)	2 (100%)	0 (%)	2
Isopropanol	4 (57%)	1 (14%)	2 (29%)	7
Lactic acid	0 (0%)	2 (67%)	1 (33%)	3
Methyl salicylate	0 (0%)	3 (100%)	0 (0%)	3
Propylene glycol	1 (50%)	1 (50%)	0 (0%)	2

1672 Abbreviations: LLNA = murine local lymph node assay; LLNA: BrdU-ELISA = murine local lymph node assay
 1673 with enzyme-linked immunosorbent assay detection of bromodeoxyuridine; SI = stimulation index.

1674 ¹Numbers shown reflect number of tests. Percentage in parentheses reflects percentage of the total number of
 1675 tests for each substance.

1676 ²According to traditional LLNA results.

1677

1677 **8.0 Data Quality**

1678 The data quality section in this draft BRD has been revised from the January 2008 draft BRD
1679 only to include data quality information about the interlaboratory validation study organized
1680 by the JSAAE.

1681 The data submitted by Dr. Takeyoshi were generated at the Hita Laboratory of the Chemicals
1682 Evaluation and Research Institute, Japan. Although the Hita Laboratory is a Good Laboratory
1683 Practice (GLP)-conforming facility, the studies on the LLNA: BrdU-ELISA did not conform
1684 fully with GLP guidelines since they were not intended for regulatory purposes. However, all
1685 systems employed for these studies (i.e., test facilities, study staff, reagents, and the other
1686 study elements) were reportedly the same as those employed in the fully GLP-compliant
1687 studies conducted in the laboratory. Although multiple staff members checked the reported
1688 data for consistency with the raw data, no audit report is available (Takeyoshi M, personal
1689 communication). The raw data are also not available for audit.

1690 The data from the interlaboratory validation study (Kojima et al. 2008) were generated in
1691 GLP laboratories, but the LLNA: BrdU-ELISA studies were not fully GLP-compliant. The
1692 data from each laboratory were reviewed by the chief of the Validation Management Team
1693 and the biostatistician.

1694

1694 **9.0. Other Scientific Reports and Reviews**

1695 This section has been revised from the January 2008 draft BRD only to include information
1696 about the interlaboratory validation study of the LLNA: BrdU-ELISA that was organized by
1697 the JSAAE. The Validation Management Team for the multi-laboratory validation study
1698 concluded that the LLNA: BrdU-ELISA, using the $SI \geq 2$ criterion to identify sensitizers, had
1699 sufficient relevance compared with the traditional LLNA and acceptable interlaboratory
1700 reproducibility (Kojima et al. 2008). The validation study has been peer reviewed in Japan.
1701 The peer review report is expected to be completed by the end of February 2009 (Kojima H,
1702 personal communication).

1703 A set of studies were conducted by Yamano et al. using a similar LLNA: BrdU-ELISA based
1704 method (Yamano et al. 2003, 2004, 2005, 2006, 2007). The test method protocol (e.g.,
1705 application of test substance to ear of mouse) was similar to what was described in the
1706 Takeyoshi et al. studies discussed above. Compared to the method Takeyoshi et al., which
1707 administered 5 mg BrdU/mouse, the concentration of BrdU administered (via intraperitoneal
1708 injection) was 150 mg/kg/15 mL saline, which would be approximately 3 mg BrdU/mouse
1709 (based on a 20 g mouse). The studies discussed the use of a BrdU-ELISA based method to
1710 assess the skin sensitization potential of a variety of substances including metal salts of
1711 naphthenic acid, methylated phenols, industrial biocides, and preservatives.

1712 The outcomes of these studies were not included in this evaluation since comparative
1713 traditional LLNA data were not available for the substances evaluated. Therefore, a
1714 comparison of the accuracy of the LLNA: BrdU-ELISA versus the traditional LLNA, when
1715 outcomes were compared to guinea pig or human results, could not be conducted.

1716

1716 **10.0 Animal Welfare Considerations**

1717 This section of the draft BRD has not changed from the January 2008 draft BRD. The
1718 LLNA: BrdU-ELISA will require the use of the same number of animals when compared to
1719 the updated ICCVAM LLNA protocol (ICCVAM 2009). However, since the traditional
1720 LLNA uses radioactivity and as such its use is restricted in some institutions, broader use of
1721 the non-radioactive LLNA: BrdU-ELISA protocol in place of the GP test could further
1722 reduce the number of guinea pigs that are still being used to assess skin sensitization.

1723 **10.1 Rationale for the Need to Use Animals**

1724 The rationale for the use of animals in the LLNA: BrdU-ELISA is the same as the rationale
1725 for the traditional LLNA; there are no valid and accepted non-animal ways to determine the
1726 ACD potential of substances and products, except for situations where human studies could
1727 be conducted ethically and where such studies would meet regulatory safety assessment
1728 requirements. The most detailed information about the induction and regulation of
1729 immunological responses are available for mice (ICCVAM 1999).

1730 **10.2 Basis for Determining the Number of Animals Used**

1731 The number of animals used for the experimental, vehicle, and positive control groups is
1732 based on the number of animals specified in the updated ICCVAM LLNA protocol
1733 (Appendix A of ICCVAM 2009).

1734 **10.3 Reduction Considerations**

1735 A further reduction of 40% (12 vs. 20) could be achieved by using a limit dose version of the
1736 LLNA: BrdU-ELISA in cases where dose response information is not needed for hazard
1737 identification purposes. In such an approach, only the highest soluble dose of the test article
1738 that does not produce skin irritation or systemic toxicity would be administered, and the two
1739 lower dose groups would not be used. Additional reductions could be achieved by testing
1740 more substances concurrently, so that the same vehicle and positive control group could be
1741 used for multiple substances, thus further reducing the number of animals for each additional
1742 substance by eight animals, or 40% (12 vs. 20).

1743

1743 **11.0 Practical Considerations**

1744 This section of the draft BRD has not changed from the January 2008 draft BRD. Several
1745 issues are taken into account when assessing the practicality of using an alternative to an
1746 existing test method. In addition to performance evaluations, assessments of the laboratory
1747 equipment and supplies needed to conduct the alternative test method, level of personnel
1748 training, labor costs, and the time required to complete the test method relative to the existing
1749 test method are necessary. The time, personnel cost, and effort required to conduct the
1750 proposed test method(s) must be considered to be reasonable when compared to the existing
1751 test method it is intended to replace.

1752 **11.1 Transferability of the LLNA: BrdU-ELISA**

1753 Test method transferability addresses the ability of a method to be accurately and reliably
1754 performed by multiple laboratories (ICCVAM 2003), including those experienced in the
1755 particular type of procedure as well as laboratories with less or no experience in the
1756 particular procedure. It would be expected that the transferability of the LLNA: BrdU-ELISA
1757 would similar to the traditional LLNA, since the protocols of the two methods (except for the
1758 detection of lymphocyte proliferation) are similar.

1759 **11.2 Facilities and Major Fixed Equipment Required to Conduct the LLNA:** 1760 **BrdU-ELISA**

1761 Compared to the traditional LLNA, the LLNA: BrdU-ELISA will not require facilities,
1762 equipment, and licensing permits for handling radioactive materials. The remaining facilities
1763 (e.g., animal care facilities) are the same between the two methods.

1764 **11.3 LLNA: BrdU-ELISA Training Considerations**

1765 The level of training and expertise needed to conduct the LLNA: BrdU-ELISA should be
1766 similar to the traditional LLNA. Additionally, individuals will need to understand and know
1767 how to perform ELISAs.

1768 **12.0 References**

- 1769 Ashby J, Basketter DA, Patona D, Kimber I. 1995. Structure activity relationships in skin
1770 sensitization using the murine local lymph node assay. *Toxicology* 103:177-194.
- 1771 Basketter DA, Blaikie L, Dearman RJ, Kimber I, Ryan CA, Gerberick GF, et al. 2000. Use of
1772 the local lymph node assay for the estimation of relative contact allergenic potency. *Contact*
1773 *Dermatitis* 42:344-348.
- 1774 Basketter DA, Cadby P. 2004. Reproducible prediction of contact allergenic potency using
1775 the local lymph node assay. *Contact Dermatitis* 50:15-17.
- 1776 Basketter DA, Clapp C, Jefferies D, Safford B, Ryan CA, Gerberick F, et al. 2005. Predictive
1777 identification of human skin sensitization thresholds. *Contact Dermatitis* 53:260-267.
- 1778 Basketter DA, Gerberick GF, Kimber I. 1998. Strategies for identifying false positive
1779 responses in predictive skin sensitization tests. *Food Chem Toxicol* 36:327-333.
- 1780 Basketter DA, Gerberick GF, Kimber I, Loveless SE. 1996. The local lymph node assay: A
1781 viable alternative to currently accepted skin sensitization tests. *Food Chem Toxicol* 34:985-
1782 997.
- 1783 Basketter DA, Kan-King-Yu D, Dierkes P, Jowsey IR. 2007. Does irritation potency
1784 contribute to the skin sensitization potency of contact allergens? *Cutan Ocul Toxicol*
1785 26(4):279-286.
- 1786 Basketter DA, Kimber I. 2001. Predictive testing in contact allergy: Facts and future. *Allergy*
1787 56: 937-943.
- 1788 Basketter DA, Kimber I. 2006. Predictive tests for irritants and allergens and their use in
1789 quantitative risk assessment. In “Contact Dermatitis”, 4th Edition. (Frosch P J, Menné T,
1790 Lepoittevin J-P, eds.) Springer Verlag, Heidelberg, pp 179 – 188.
- 1791 Basketter DA, Lea LJ, Cooper KJ, Ryan CA, Gerberick GF, Dearman RJ, et al. 1999a.
1792 Identification of metal allergens in the local lymph node assay. *Am J Contact Dermat*
1793 10:207-212.

- 1794 Basketter DA, Lea LJ, Cooper K, Stocks J, Dickens A, Pate I, et al. 1999b. Threshold for
1795 classification as a skin sensitizer in the local lymph node assay: A statistical evaluation. *Food*
1796 *Chem Toxicol* 37:1167-1174.
- 1797 Basketter DA, Lea LJ, Dickens A, Briggs D, Pate I, Dearman RJ, et al. 1999c. A comparison
1798 of statistical approaches to the derivation of EC3 values from local lymph node assay dose
1799 responses. *J Appl Toxicol* 19:261-266.
- 1800 Basketter DA, Scholes EW, Kimber I, Botham PA, Hilton J, Miller K, et al. 1991.
1801 Interlaboratory Evaluation of the Local Lymph Node Assay with 25 Chemicals and
1802 Comparison with Guinea Pig Test Data. *Toxicology Methods* 1:30-43.
- 1803 Basketter DA, Scholes EW. 1992. Comparison of the Local Lymph Node Assay with the
1804 Guinea-Pig Maximization Test for the Detection of a Range of Contact Allergens. *Food*
1805 *Chem Toxicol* 30:65-69.
- 1806 Basketter DA, Scholes EW, Cumberbatch M, Evans CD, Kimber I. 1992. Sulphanilic acid:
1807 Divergent results in the guinea pig maximization test and the local lymph node assay.
1808 *Contact Dermatitis* 27:209-213.
- 1809 Basketter DA, Smith Pease CK, Patlewicz GY. 2003. Contact allergy: the local lymph node
1810 assay for the prediction of hazard and risk. *Clin Exp Dermatol* 28:218-221.
- 1811 Betts CJ, Dearman RJ, Heylings JR, Kimber I, Basketter DA. 2006. Skin sensitization
1812 potency of methyl methacrylate in the local lymph node assay: Comparisons with guinea-pig
1813 data and human experience. *Contact Dermatitis* 55:140-147.
- 1814 Bjorkner B. 1984. Contact allergy to 2-hydroxypropyl methacrylate (2-HPMA) in an
1815 ultraviolet curable ink. *Acta Derm Venereol* 64: 264-267.
- 1816 Buehler EV. 1965. Delayed contact hypersensitivity in the guinea pig. *Arch Dermatol*
1817 91:171-177.
- 1818 Cosmetic Ingredient Review Expert Panel. 1998. Final report on the safety assessment of
1819 glycolic acid, ammonium, calcium, potassium, and sodium glycolates, methyl, ethyl, propyl,
1820 and butyl glycolates, and lactic acid, ammonium, calcium, potassium, sodium, and tea-
1821 lactates, methyl, ethyl, isopropyl, and butyl lactates, and lauryl, myristyl, and cetyl lactates.
1822 *Int J Toxicol* 17(Suppl 1):1-203.

- 1823 Dean JH, Twerdok LE, Tice RR, Dailstad DM, Hattan DG, Stokes WS. 2001. ICCVAM
1824 Evaluation of the Murine Local Lymph Node Assay II. Conclusions and Recommendations
1825 of an Independent Scientific Peer Review Panel. *Regul Toxicol Pharmacol* 34:258-273.
- 1826 Dearman RJ, Betts CJ, Farr C, McLaughlin J, Berdasco N, Wiench K, et al. 2007.
1827 Comparative analysis of skin sensitization potency of acrylates (methyl acrylate, ethyl
1828 acrylate, butyl acrylate, and ethylhexyl acrylate) using the local lymph node assay. *Contact*
1829 *Dermatitis* 57:242-247.
- 1830 Dearman DJ, Kimber I. 2001. Development of the Local Lymph Node Assay for risk
1831 assessment of chemical and formulations. HSE contract research report 399. Syngenta
1832 Central Toxicology Laboratory, UK.
- 1833 ECETOC. 1995. Skin Irritation and Corrosion: Reference Chemicals Data Bank. Technical
1834 Report Number 66. Brussels: European Centre for Ecotoxicology and Toxicology of
1835 Chemicals.
- 1836 EPA. 2003. Health Effects Test Guidelines: OPPTS 870.2600 - Skin Sensitization. EPA 712–
1837 C–98–197. Washington, DC: U.S. Environmental Protection Agency. Available:
1838 [http://www.epa.gov/opptsfrs/publications/OPPTS_Harmonized/870_Health_Effects_Test_G](http://www.epa.gov/opptsfrs/publications/OPPTS_Harmonized/870_Health_Effects_Test_Guidelines/Series/870-2600.pdf)
1839 [uidelines/Series/870-2600.pdf](http://www.epa.gov/opptsfrs/publications/OPPTS_Harmonized/870_Health_Effects_Test_Guidelines/Series/870-2600.pdf)
- 1840 Gad SC, Dunn BJ, Dobbs DW, Reilly C, Walsh RD. 1986. Development and validation of an
1841 alternative dermal sensitization test: the mouse ear swelling test (MEST). *Toxicol Appl*
1842 *Pharmacol* 84:93-114.
- 1843 Gerberick GF, Cruse LW, Ryan CA, Hulette BC, Chaney JG, Skinner RA, et al. 2002. Use of
1844 a B cell marker (B220) to discriminate between allergens and irritants in the local lymph
1845 node assay. *Toxicol Sci* 68(2):420-428.
- 1846 Gerberick GF, House RV, Fletcher ER, Ryan CA. 1992. Examination of the local lymph
1847 node assay for use in contact sensitization risk assessment. *Fundam Appl Toxicol* 19:438-
1848 445.
- 1849 Gerberick GF, Ryan CA, Kern PS, Dearman RJ, Kimber I, Patlewicz GY, et al. 2004. A
1850 chemical dataset for evaluation of alternative approaches to skin-sensitization testing.
1851 *Contact Dermatitis* 50:274-288.

- 1852 Gerberick GF, Ryan CA, Kern PS, Schlatter H, Dearman RJ, Kimber I, et al. 2005.
1853 Compilation of historical local lymph node data for evaluation of skin sensitization
1854 alternative methods. *Dermatitis* 16:157-202.
- 1855 Gerberick GF, Ryan CA, Dearman RJ, Kimber I. 2007. Local lymph node assay (LLNA) for
1856 detection of sensitization capacity of chemicals. *Methods* 41:54-60.
- 1857 Gerberick GF, Vassallo JD, Foertsch LM, Price BB, Chaney JG, Lepoittevin JP. 2007.
1858 Quantification of chemical peptide reactivity for screening contact allergens: a classification
1859 tree model approach. *Toxicol Sci* 97:417-27.
- 1860 Hilton J, Dearman RJ, Harvey P, Evans P, Basketter DA, Kimber I. 1998. Estimation of
1861 relative skin sensitizing potency using the local lymph node assay: A comparison of
1862 formaldehyde with glutaraldehyde. *Am J Contact Dermat* 9:29-33.
- 1863 Hutchings CV, Shum KW, Gawkrödger DJ. 2001. Occupational contact dermatitis has an
1864 appreciable impact on quality of life. *Contact Dermatitis* 45:17-20.
- 1865 ICCVAM. 1997. Validation and Regulatory Acceptance of Toxicological Test Methods: A
1866 Report of the ad hoc Interagency Coordinating Committee on the Validation of Alternative
1867 Methods. NIH Publication No.: 97-3981. Research Triangle Park: National Institute of
1868 Environmental Health Sciences.
- 1869 ICCVAM 1999. The murine local lymph node assay: A test method for assessing the allergic
1870 contact dermatitis potential of chemical/compounds. NIH Publication No. 99-4494. Research
1871 Triangle Park, NC: National Institute of Environmental Health Sciences.
- 1872 ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised,
1873 and Alternative Test Methods. NIH Publication No: 03-4508. Research Triangle Park:
1874 National Institute of Environmental Health Sciences.
- 1875 ICCVAM. 2009. Recommended Performance Standards: Murine Local Lymph Node Assay.
1876 NIH Publication Number 09-7357. Research Triangle Park, NC: National Institute of
1877 Environmental Health Sciences. Available at
1878 http://iccvam.niehs.nih.gov/methods/immunotox/llna_PerfStds.htm.
- 1879 ICCVAM. 2008. Independent Scientific Peer Review Panel Report. Validation Status of New

- 1880 Versions and Applications of the Murine Local Lymph Node Assay: A Test Method for
1881 Assessing the Allergic Contact Dermatitis Potential of Chemicals and Products. Research
1882 Triangle Park, NC: National Institute of Environmental Health Sciences. Available at:
1883 http://iccvam.niehs.nih.gov/docs/immunotox_docs/LLNAPRPRpt2008.pdf.
- 1884 ISO. 2002. Biological evaluation of medical devices -- Part 10: Tests for irritation and
1885 delayed-type hypersensitivity. Available for purchase at: <http://www.iso.org/iso/home.htm>.
- 1886 Jowsey IR, Basketter DA, Westmoreland C, Kimber I. 2006. A future approach to measuring
1887 relative skin sensitising potency: a proposal. *J Appl Toxicol* 26:341-350.
- 1888 Kimber I, Basketter DA, Butler M, Gamer A, Garrigue J-L, Gerberick GF, Newsome C,
1889 Steiling W, Vohr H-W. 2003. Classification of contact allergens according to potency:
1890 proposals. *Food Chem Toxicol* 41:1799-1809.
- 1891 Kimber I, Dearman RJ. 1991. Investigation of lymph node cell proliferation as a possible
1892 immunological correlate of contact sensitizing potential. *Food Chem Toxicol* 29: 125-129.
- 1893 Kimber I, Dearman RJ. 1996. Contact hypersensitivity: immunological mechanisms. In
1894 *Toxicology of Contact Hypersensitivity*. (Kimber I, Maurer T, eds). London:Taylor and
1895 Francis, 4-25.
- 1896 Kimber I, Dearman RJ, Scholes EW, Basketter DA. 1994. The local lymph node assay:
1897 developments and applications. *Toxicology* 93:13-31.
- 1898 Kimber I, Hilton J, Botham PA, Basketter DA, Scholes EW, Miller K, et al. 1991. The
1899 murine local lymph node assay: results of an inter-laboratory trial. *Toxicol Lett* 55:203-213.
- 1900 Kimber I, Hilton J, Dearman RJ, Gerberick GF, Ryan CA, Basketter DA, et al. 1995. An
1901 international evaluation of the murine local lymph node assay and comparison of modified
1902 procedures. *Toxicology* 103:63-73.
- 1903 Kimber I, Hilton J, Dearman RJ, Gerberick GF, Ryan CA, Basketter DA, et al. 1998.
1904 Assessment of the skin sensitization potential of topical medicaments using the local lymph
1905 node assay: An interlaboratory evaluation. *J Toxicol Environ Health - Part A* 53:563-579.

- 1906 Kimber I, Hilton J, Weisenberger C. 1989. The murine local lymph node assay for
1907 identification of contact allergens: A preliminary evaluation of in situ measurement of
1908 lymphocyte proliferation. *Contact Dermatitis* 21:215-220.
- 1909 Klecak G, Geleicke H, Frey JR. 1977. Screening of fragrance materials for allergenicity in
1910 the guinea pig. I. Comparison of four testing methods. *J Soc Cosmet Chem* 28:53-64.
- 1911 Kligman AM. 1966. The Identification of Contact Allergens by Human Assay. *J Invest*
1912 *Dermatol* 47:369-409.
- 1913 Kojima H, Takeyoshi M, Omori T, Sozu, T, Arima K, Idehara K, Ikarashi Y, Kanazawa Y,
1914 Maki E, Nakagiri N, Tanaka M, Yuasa A, Yoshimura I. 2008. Inter-laboratory validation
1915 study on LLNA-BrdU. *Toxicologist* 102:216.
- 1916 Kwon JA, Lee MS, Kim MI, Park YM, Kim HO, Kim CW. 2003. Allergic contact dermatitis
1917 from dodecyl diaminoethylglycine and isopropyl alcohol in a commercial disinfectant swab.
1918 *Contact Dermatitis* 48:339-340.
- 1919 Loveless SE, Ladics GS, Gerberick GF, Ryan CA, Basketter DA, Scholes EW, et al. 1996.
1920 Further evaluation of the local lymph node assay in the final phase of an international
1921 collaborative trial. *Toxicology* 108:141-152.
- 1922 Magnusson B, Kligman AM. 1970. *Allergic Contact Dermatitis in the Guinea Pig*.
1923 Springfield, IL: Charles C. Thomas.
- 1924 Manetz TS, Meade BJ. 1999. Development of a flow cytometry assay for the identification
1925 and differentiation of chemicals with the potential to elicit irritation, IgE- mediated, or T cell-
1926 mediated hypersensitivity responses. *Toxicol Sci* 48:206-217.
- 1927 Marzulli FN, Maibach HI. 1974. The use of graded concentrations in studying skin
1928 sensitizers: experimental contact sensitization in man. *Food Cosmet Toxicol* 12:219-227.
- 1929 Marzulli FN, Maibach HI. 1980. Contact Allergy: Predictive Testing of Fragrance
1930 Ingredients in Humans by Draize and Maximization Methods. *J Environ Pathol Toxicol*
1931 3:235-245.
- 1932 Meylan WM, Howard PH. 1995. Atom/fragment contribution method for estimating octanol-
1933 water partition coefficients. *J Pharm Sci.* Jan;84(1):83-92.

- 1934 Moriguchi I, Hirono S, Nakagome I, Hirano H. 1994. Comparison of Reliability of log *P*
1935 Values for Drugs Calculated by Several Methods. *Chem Pharm Bull* 42(3):976-978.
- 1936 OECD. 2002. Guideline For Testing of Chemicals - Test Guideline 429: Skin Sensitisation:
1937 Local Lymph Node Assay. Paris:OECD. Available:
1938 [http://puck.sourceoecd.org/vl=4134815/cl=32/nw=1/rpsv/cw/vhosts/oecdjournals/1607310x/
1939 v1n4/contp1-1.htm](http://puck.sourceoecd.org/vl=4134815/cl=32/nw=1/rpsv/cw/vhosts/oecdjournals/1607310x/v1n4/contp1-1.htm)
- 1940 OECD. 2002. Test guideline 429. Skin Sensitisation: Local Lymph Node Assay, adopted
1941 April 24, 2002. In: OECD Guidelines for Testing of Chemicals. Paris:OECD.
- 1942 Opdyke DLJ. 1976. Fragrance raw materials monographs. Isopropyl myristate. *Food Cosmet*
1943 *Toxicol* 14:323-325.
- 1944 Robinson MK, Nusair TL, Fletcher ER, Ritz HL. 1990. A review of the Buehler guinea pig
1945 skin sensitization test and its use in a risk assessment process for human skin sensitization.
1946 *Toxicology* 61:91-107.
- 1947 Ryan CA, Chaney JG, Gerberick GF, Kern PS, Dearman RJ, Kimber I, Basketter DA. 2007.
1948 Extrapolating local lymph node assay EC3 values to estimate relative sensitizing potency.
1949 *Cutan Ocul Toxicol* 26:135-145.
- 1950 Ryan CA, Cruse LW, Skinner RA, Dearman RJ, Kimber I, Gerberick GF. 2002. Examination
1951 of a vehicle for use with water soluble materials in the murine local lymph node assay. *Food*
1952 *Chem Toxicol* 40:1719-1725.
- 1953 Ryan CA, Gerberick GF, Cruse LW, Basketter DA, Lea L, Blaikie L, et al. 2000. Activity of
1954 human contact allergens in the murine local lymph node assay. *Contact Dermatitis* 43:95-
1955 102.
- 1956 Sailstad DM, Hattan D, Hill RN, Stokes WS. 2001. ICCVAM Evaluation of the Murine
1957 Local Lymph Node Assay (LLNA) I: The ICCVAM Review Process. *Regul Toxicol*
1958 *Pharmacol* 34:249-257.
- 1959 Schneider K, Akkan Z. 2004. Quantitative relationship between the local lymph node assay
1960 and human skin sensitization assays. *Regul Toxicol Pharmacol* 39:245-255.

- 1961 Scholes EW, Basketter DA, Sarll AE, Kimber I, Evans CD, Miller K, et al. 1992. The local
1962 lymph node assay: Results of a final inter-laboratory validation under field conditions. *J Appl*
1963 *Toxicol* 12:217-222.
- 1964 Skoet R, Zachariae R, Agner T. 2003. Contact dermatitis and quality of life: a structured
1965 review of the literature. *Br J Dermatol* 149:452-456.
- 1966 Takeyoshi M, Iida K, Shiraishi K, Hoshuyama S. 2005. Novel approach for classifying
1967 chemicals according to skin sensitizing potency by non-radioisotopic modification of the
1968 local lymph node assay. *J Appl Toxicol* 25:129-34.
- 1969 Takeyoshi M, Iida K, Suzuki K, Yamazaki S. 2007a. Skin sensitization potency of
1970 isoeugenol and its dimers evaluated by a non-radioisotopic modification of the local lymph
1971 node assay and guinea pig maximization test. *J Appl Toxicol* DOI: 10.1002/jat.1305 [Online
1972 10 October 2007].
- 1973 Takeyoshi M, Minobe Y, Takatsuki M. 2007b. Promising endpoint and assay performance of
1974 non-radioisotopic Local Lymph Node Assay based on BrdU incorporation. Presentation at
1975 the 6th World Congress on Alternatives and Animal Use in the Life Sciences. Tokyo, Japan.
- 1976 Takeyoshi M, Noda S, Yamasaki K. 2004a. Differences in responsiveness of mouse strain
1977 against p-benzoquinone as assessed by non-radioisotopic murine local lymph node assay.
1978 *Exp Anim* 53:171-3.
- 1979 Takeyoshi M, Noda S, Yamazaki S, Kakishima H, Yamasaki K, Kimber I. 2004b.
1980 Assessment of the skin sensitization potency of eugenol and its dimers using a non-
1981 radioisotopic modification of the local lymph node assay. *J Appl Toxicol* 24:77-81.
- 1982 Takeyoshi M, Noda S, Yamasaki K, Kimber I. 2006. Advantage of using CBA/N strain mice
1983 in a non-radioisotopic modification of the local lymph node assay. *J Appl Toxicol* 26:5-9.
- 1984 Takeyoshi M, Sawaki M, Yamasaki K, Kimber I. 2003. Assessment of statistic analysis in
1985 non-radioisotopic local lymph node assay (non-RI-LLNA) with alpha-hexylcinnamic
1986 aldehyde as an example. *Toxicology* 191:259-63.
- 1987 Takeyoshi M, Yamasaki K, Yakabe Y, Takatsuki M, and Kimber I. 2001. Development of
1988 non-radio isotopic endpoint of murine local lymph node assay based on 5-bromo-2'-

- 1989 deoxyuridine (BrdU) incorporation. *Toxicology Letters* 119:203-208.
- 1990 Van Der Walle HB, Klecak G, Geleick H, Bensink T. 1982. Sensitizing potential of 14 mono
1991 (meth) acrylates in the guinea pig. *Contact Dermatitis* 8:223-235.
- 1992 Van Och FMM, Vandebriel RJ, Prinsen MK, De Jong WH, Slob W, Van Loveren H. 2001.
1993 Comparison of dose-responses of contact allergens using the guinea pig maximization test
1994 and the local lymph node assay. *Toxicology* 167:207-215.
- 1995 Wahlberg JE, Boman A. 1985. Guinea pig maximization test. *Curr Probl Dermatol* 14:59-
1996 106.
- 1997 Yamano T, Ichihara M, Shimizu M, Noda T, Tsujimoto Y. 2007. Immunomodulatory effects
1998 of mono-, di-, and trimethylphenols in mice. *Toxicology* 232:132-7.
- 1999 Yamano T, Shimizu M, Noda T. 2003. Allergenicity evaluation of p-chloro-m-cresol and p-
2000 chloro-m-xyleneol by non-radioactive murine local lymph-node assay and multiple-dose
2001 guinea pig maximization test. *Toxicology* 190:259-66.
- 2002 Yamano T, Shimizu M, Noda T. 2004. Allergenicity evaluation of Bioban CS-1135 in
2003 experimental animals. *Contact Dermatitis* 50:339-43.
- 2004 Yamano T, Shimizu M, Noda T. 2005. Quantitative comparison of the results obtained by the
2005 multiple-dose guinea pig maximization test and the non-radioactive murine local lymph-node
2006 assay for various biocides. *Toxicology* 211:165-75.
- 2007 Yamano T, Shimizu M, Noda T. 2006. Allergenicity and cross-reactivity of naphthenic acid
2008 and its metallic salts in experimental animals. *Contact Dermatitis* 54:25-8.