

1 **NOTE: The European Centre for the Validation of Alternative Methods (ECVAM)**
2 **has prepared a comprehensive background review document (BRD) for the**
3 **Cytosensor (CM) test method and has agreed to make it publicly available for**
4 **review and comment. However, they have not been able to obtain agreement for**
5 **release of this document from companies that produced unpublished data that are**
6 **included in the BRD. Therefore, the following draft proposed recommendations are**
7 **based on data obtained from three peer reviewed publications (Balls et al. 1995,**
8 **Gettings et al. 1996, Brantom et al. 1997). Subsequent to these peer review activities,**
9 **ICCVAM will work with ECVAM (which is conducting a peer review of CM by the**
10 **ECVAM Scientific Advisory Committee), the Japanese Center for the Validation of**
11 **Alternative Methods, and Health Canada to develop final test recommendations for**
12 **CM.**

13 **Draft Proposed ICCVAM Recommendations for the Cytosensor Test**

14 **Method: Uses and Limitations**

15 **Use of the Cytosensor Test Method (INVITTOX Protocol Number 102) to Identify** 16 **Ocular Corrosives and Severe Irritants**

17 The database of 53 water-soluble surfactants tested using INVITTOX Protocol 102
18 includes 21 surfactant chemicals and 32 surfactant-containing formulations tested across
19 seven different laboratories. Most of the 32 formulations, which are limited to cosmetic
20 and personal care products, contain one or more surfactants at a final concentration of
21 greater than five percent. There were no pesticide formulations included in the validation
22 database. Using INVITTOX 102 to identify ocular corrosives and severe irritants among
23 these surfactant-containing substances, the false positive rate ranged from 3-10% (1/29 to
24 3/29) when compared to *in vivo* results. The three false positives when using the EPA
25 classification system are classified as Category II (n=2) or III (n=1) based on *in vivo* data.
26 The one false positive when using the GHS and EU classification systems are Not
27 Labeled based on *in vivo* data. The false negative rate ranged from 9-22% (2/23 to 5/23)
28 when compared to *in vivo* results. In each case, these substances were classified as mild
29 or moderate irritants *in vitro* based on the EPA, EU, and GHS classification systems (i.e.,
30 Category II/III, R36, or Category 2A/2B, respectively).

31 The nonsurfactant substances database for INVITTOX 102 consisted of 29 water-soluble
32 nonsurfactant chemicals (n=27), which included a range of chemical classes (e.g., acids,
33 alcohols, alkalis, and ketones), and nonsurfactant formulations (n=2) tested in seven
34 laboratories. Using INVITTOX 102 to identify ocular corrosives and severe irritants
35 among these nonsurfactant substances, the false positive rate ranged from 0-6% (0/18 to
36 1/18) when compared to *in vivo* results. The one false positive when using the EPA or EU
37 classification systems was Category III and R36, respectively based on *in vivo* data.
38 There were no false positives when using the GHS classification system. The false
39 negative rate ranged from 43-55% (3/7 to 6/11) when compared to *in vivo* results. Three
40 substances were false negatives when using the EPA classification system and were
41 classified *in vitro* as either Category II/III (n=2) or IV (n=1). Five substances were false
42 negatives when using the GHS classification system and were classified *in vitro* as either

43 Category 2A/2B (n=4) or Not Labeled (n=1). Six substances were false negatives when
44 using the EU classification system and were classified *in vitro* as either R36 (n=5) or Not
45 Labeled (n=1).

46 Based on these data and test method performance, ICCVAM proposes that the
47 Cytosensor test method can be used as a screening test to identify water-soluble
48 substances as ocular corrosives and severe irritants (i.e., EPA Category I, GHS Category
49 1, EU R41) in a tiered-testing strategy, as part of a weight-of-evidence approach. A
50 substance that tests negative with Cytosensor would need to be tested in another test
51 method that is capable of identifying possible *in vitro* false negative severe irritants and
52 ocular corrosives and to distinguish between moderate and mild ocular irritants.
53 Currently, the *in vivo* rabbit eye test is the only test method capable of making such a
54 distinction.

55 **Use of the Cytosensor Test Method (INVITTOX Protocol Number 102) to Identify** 56 **Substances Not Labeled as Irritants**

57 Using INVITTOX 102 to identify substances not labeled as irritants among the database
58 of 53 water soluble surfactants and surfactant-containing formulations, the false negative
59 rate ranged from 0-2% (0/27 to 1/46) when compared to *in vivo* results. The one false
60 negative, which occurred only for the EPA classification system, was classified as
61 Category III based on *in vivo* data. For this substance, six rabbits were included in the *in*
62 *vivo* test. One rabbit had no observable effects, three rabbits had conjunctival redness
63 (score = 1) that cleared after one (n=1) or two days (n=2), and two rabbits had corneal
64 opacity (score = 1) that cleared after one day. The false positive rate ranged from 50-69%
65 (3/6 to 18/26) when compared to *in vivo* results. Three substances were false positives
66 when using the EPA classification system and were classified *in vitro* as Category II/III.
67 Seventeen substances were false positives when using the GHS classification system and
68 were classified *in vitro* as Category 2A/2B (n=16) or Category 1 (n=1). Eighteen
69 substances were false positives when using the EU classification system and classified *in*
70 *vitro* as R36 (n=17) or R41 (n=1).

71 Using INVITTOX 102 to identify substances not labeled as irritants among the database
72 of 29 nonsurfactant substances, the false negative rate ranged from 24-38% (5/21 to 8/21)

73 when compared to *in vivo* results. The false positive rate ranged from 25-40% (1/4 to 2/5)
74 when compared to *in vivo* results.

75 Based on these data, ICCVAM proposes that the Cytosensor test method can be used as a
76 screening test to identify water-soluble surfactant chemicals and certain types of
77 surfactant-containing formulations (e.g., cosmetics and personal care product
78 formulations, but not pesticide formulations) as substances not labeled as irritants (i.e.,
79 EPA Category IV, GHS Category NL, EU Category NL) in a tiered-testing strategy, as
80 part of a weight-of-evidence approach. However, based on the false positive rate, a
81 substance that tests positive with Cytosensor would need to be tested in another test
82 method that is capable of correctly identifying possible *in vitro* false positives. Positives
83 would also need to be additionally tested with methods that can correctly identify severe,
84 moderate, and mild ocular irritants.

85 Because of the high false negative rate for Cytosensor when testing water-soluble
86 nonsurfactant substances and formulations, Cytosensor is not recommended as a
87 screening test to identify substances not labeled as irritants among these types of
88 substances.

89 **Use of the Cytosensor Test Method (INVITTOX Protocol Number 102) to Identify**
90 **Either Ocular Corrosives/Severe Irritants or Substances Not Labeled as Irritants**

91 Given that the Cytosensor test method (INVITTOX Protocol Number 102) is proposed
92 for use as a screening test to identify both ocular corrosive/severe irritants and
93 nonirritants, specifically for water-soluble surfactant chemicals and specific types of
94 surfactant-containing formulations (e.g., cosmetics and personal care product
95 formulations, but not pesticide formulations), users may want to consider using
96 Cytosensor prior to another *in vitro* ocular test method for testing these types of
97 substances. However, water-soluble surfactant chemicals and surfactant formulations that
98 are not identified as ocular corrosive/severe irritants or as nonirritants with Cytosensor
99 would need to be tested in another test method(s) capable of correctly classifying
100 substances into each of the four hazard classification categories for EPA or GHS.
101 Currently, the only test method accepted for these purposes is the *in vivo* Draize test.
102 Because of the high false positive rate (> 50%) for the non-irritant decision criteria, users

103 may not want to use Cytosensor if the intended use is to start with identifying
104 nonirritants.

105 **References**

106 Balls M, Botham PA, Bruner LH, Spielmann H. 1995. The EC/HO international
107 validation study on alternatives to the Draize eye irritation test. *Toxicol In Vitro* 9:871-
108 929.

109 Brantom, PG, Bruner, LH, Chamberlain M, de Silva O, Dupuis J, Earl LK, Lovell DP,
110 Pape WJW, Uttley M, Bagley DM, Baker FW, Bracher, M, Courtellemont P, Declercq L,
111 Freeman S, Steiling W, Walker AP, Carr GJ, Dami N, Thomas G, Harbell J, Jones PA,
112 Pfannenbecker U, Southee JA, Tchong M, Argembeaux H, Castelli D, Clothier R, Esdaile
113 DJ, Itigaki H, Jung K, Kasai Y, Kojima H, Kristen U, Larnicol M, Lewis RW, Marenus
114 K, Moreno O, Peterson A, Rasmussen ES, Robles C, Stern M. 1997. Summary report of
115 the COLIPA international validation study on alternatives to the Draize rabbit eye
116 irritation test. *Toxicol In Vitro* 11:141-179.

117 Gettings S, Lordo R, Hintze K, Bagley D, Casterton P, Chudkowski M, Curren RD,
118 Demetrulias JL, Dipasquale LC, Earl LK, Feder PI, Galli CL, Glaza SM, Gordon VC,
119 Janus J, Kurtz PJ, Marenus KD, Moral J, Pape WJW, Renskers KJ, Rheins LA, Roddy
120 MT, Rozen MG, Tedeschi JP, Zyrocki J. 1996. The CTFA evaluation of alternatives
121 program: An evaluation of in vitro alternatives to the Draize primary eye irritation test.
122 (Phase III) surfactant-based formulations. *Food Chem Toxic* 34:79-117.