

INSTRUMENT

HB-1


BEAM LINE HIGH FLUX ISOTOPE REACTOR

Fact Sheet



POLARIZED TRIPLE-AXIS SPECTROMETER

The HB-1 Polarized Triple-Axis Spectrometer is designed primarily for the study of excitations in crystalline solids at intermediate energies. Thanks to the vertical beam focusing and the very high time-averaged flux at HFIR, its geometry is optimal for investigating small samples and weak scattering in specific areas of energy-momentum space. The sample goniometers and a full software implementation of the three-dimensional sample orientation matrix allow measurements outside the traditional single-scattering plane. The unique capability of HB-1 is the polarized configuration for studies of excitations, phase transitions, structures, and density distributions in magnetic materials. Use of a beryllium low-pass filter makes the instrument suitable for hypothermal neutron measurements with incident or final energy below 5 meV.



SPECIFICATIONS

Beam spectrum	Thermal
Monochromators	Unpolarized PG(002) Polarized Heusler (1 1 1)
Analyzers	Unpolarized PG(002), Be(101), Be(002) Polarized Heusler (1 1 1)
Monochromator takeoff angle	$2\theta_M = 18$ to 75°
Sample angles	360°
Scattering angle	-90 to 140°
Analyzer angles	-40 to 140°
Collimations (FWHM)	Premonochromator: 15', 30', 48' Monochromator-sample: 20', 40', 60', 80' Sample-analyzer: 20', 40', 60', 80' Analyzer-detector: 70', 120', 240'

Status: Operational

APPLICATIONS

The following are some of the scientific applications for which the Triple-Axis Spectrometer is particularly well suited.

- Spin waves in ordered magnetic materials
- Exotic excitations in low-dimensional, molecular, itinerate, and other "quantum" magnets
- Spin and lattice excitations in high- T_c superconductivity, colossal magnetoresistance materials, and multiferroic systems
- Spin density distributions in magnetic compounds
- Phonon dispersion curves in alloys and phonon-driven phase transitions

FOR MORE INFORMATION, CONTACT

Instrument Scientist: Andrey Zheludev, zheludevai@ornl.gov, 865.241.0098

http://neutrons.ornl.gov/hfir_instrument_systems/HB-1.shtml



February 2009