

1
2
3
4
5
6
7
8
9
10
11
12
13

DRAFT ICCVAM TEST METHOD RECOMMENDATIONS

***In Vitro* Pyrogenicity Test Methods**

December 1, 2006

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

[This Page Intentionally Left Blank]

28

TABLE OF CONTENTS

29

Page Number

30	1.0	ICCVAM Draft Recommendations for <i>In Vitro</i> Pyrogenicity Test Methods	7
31	1.1	Draft Recommended Test Method Uses	7
32	1.2	Draft Recommended Future Studies	9
33			
34	Appendix A	Draft Performance Standards for <i>In Vitro</i>	
35		Pyrogenicity Test Methods	10
36	1.0	Purpose and Background of Performance Standards	12
37	1.1	Introduction	12
38	1.2	Elements of ICCVAM Performance Standards.....	13
39	1.3	ICCVAM Process for the Development of Performance Standards.....	13
40	1.4	ICCVAM Development of Performance Standards for <i>In Vitro</i>	
41		Pyrogenicity Test Methods.....	14
42	1.4.1	Current Regulatory Testing Requirements for Pyrogenicity	15
43	1.4.2	Test Methods for Assessing Pyrogenicity	15
44	1.4.3	Intended Regulatory Uses for <i>In Vitro</i> Pyrogenicity Test	
45		Methods	16
46	1.4.4	Similarities and Differences in the Endpoints of <i>In Vitro</i>	
47		Pyrogenicity Test Methods and Currently Recognized	
48		Pyrogenicity Test Methods	16
49	2.0	<i>In Vitro</i> Pyrogenicity Test Methods	17
50	2.1	Background.....	17
51	2.2	Principles of <i>In Vitro</i> Pyrogenicity Test Methods	17
52	2.3	Essential Test Method Components for <i>In Vitro</i> Pyrogenicity	
53		Test Methods	18
54	2.3.1	<i>In Vitro</i> Cell Culture Conditions	19
55	2.3.2	Application of the Test Substances	20
56	2.3.2.1	Test Substance Preparation.....	20
57	2.3.2.2	Test Substance Application and Sample Collection	20

58	2.3.3	Control Substances	21
59	2.3.3.1	Negative Control.....	21
60	2.3.3.2	Positive Control	21
61	2.3.3.3	Benchmark Controls.....	21
62	2.3.4	Cytokine Measurements	22
63	2.3.5	Interpretation of Results.....	22
64	2.3.6	Test Report.....	22
65	2.4	Reference Substances for <i>In vitro</i> Pyrogenicity Test Methods	24
66	2.5	Accuracy and Reliability	25
67	2.5.1	Accuracy	25
68	2.5.2	Reliability.....	26
69	2.5.2.1	Intralaboratory Repeatability	27
70	2.5.2.2	Intralaboratory Reproducibility	27
71	2.5.2.3	Interlaboratory Reproducibility	29
72	3.0	References	31
73	Appendix B	Proposed <i>In Vitro</i> Pyrogenicity Test Method Protocols.....	B-1
74	B1	Cryopreserved Whole Blood (Cryo WB)/Interleukin-1 (IL-1 β) <i>In Vitro</i>	
75		Pyrogen Test.....	B1-1
76	B2	Mono Mac 6 (MM6)/Interleukin-6 (IL-6) <i>In Vitro</i> Pyrogen Test	B2-1
77	B3	Peripheral Blood Mononuclear Cell (PBMC)/Interleukin-6 (IL-6)	
78		<i>In Vitro</i> Pyrogen Test.....	B3-1
79	B4	Whole Blood (WB)/Interleukin-1 β (IL-1 β) <i>In Vitro</i> Pyrogen Test.....	B4-1
80	B5	Whole Blood (WB)/Interleukin -6 (IL-6) <i>In Vitro</i> Pyrogen Test	B5-1
81			

82

83

LIST OF TABLES

84

Page Number

85 Table 1 Performance Statistics for *In Vitro* Pyrogenicity Test Methods..... 8

86 Table 2 Test Substances (Parenteral Drugs) Used in the Validation Studies

87 for Determining Test Method Performance 8

88 Table 1-1 Summary of U.S. and European Legislation and Statutory Protocol

89 Requirements for Pyrogenicity Testing..... 16

90 Table 2-1 Recommended Reference Substances for the *In Vitro* Pyrogenicity

91 Test Methods 26

92 Table 2-2 Performance Statistics for *In Vitro* Pyrogenicity Test Methods..... 2793 Table 2-3 Intralaboratory Reproducibility of *In Vitro* Pyrogenicity

94 Test Methods 29

95 Table 2-4 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity

96 Test Methods 30

97 Table 2-5 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity

98 Test Methods 31

99

100

101

101

102

103

104

105

106

107

108

109

[This Page Intentionally Left Blank]

110

111

112 **1.0 ICCVAM DRAFT RECOMMENDATIONS FOR *IN VITRO***
113 **PYROGENICITY TEST METHODS**

114

115 **1.1 Draft Recommended Test Method Uses**

116

117 ICCVAM has evaluated the validation status of the following *in vitro* test methods
118 proposed as replacements for the *in vivo* rabbit pyrogen test (RPT) only (i.e., not for the
119 bacterial endotoxin test [BET]):

- 120 • cryo WB/IL-1 (The Human Whole Blood/IL-1 *In Vitro* Pyrogen Test:
121 Application of cryopreserved human whole blood)
- 122 • MM6/IL6 (An Alternative *In Vitro* Pyrogen Test Using the Human
123 Monocytoid Cell Line MONO MAC-6 [MM6])
- 124 • PBMC/IL-6 (The Human Peripheral Blood Mononuclear Cell [PBMC]/IL-
125 6 *In Vitro* Pyrogen Test)
- 126 • WB/IL-1 (The Human Whole Blood [WB]/IL-1 *In Vitro* Pyrogen Test)
- 127 • WB/IL-6 (The Human Whole Blood/IL-6 *In Vitro* Pyrogen Test)

128

129 There is sufficient information (see **Table 1**), based on validation studies with a limited
130 number of pharmaceuticals (see **Table 2**), to substantiate the use of these test methods
131 (PBMC/IL-6, cryo WB/IL-1 [96 well plate method], WB/IL-6, and MM6/IL-6) for the
132 detection of pyrogenicity mediated by Gram-negative endotoxin in materials that are
133 currently tested in the RPT, subject to product-specific validation to demonstrate
134 equivalency^{1,2}. While the scientific basis of these test methods suggests that they have the
135 capability to detect pyrogenicity produced by a wider range of pyrogens (i.e., those
136 mediated by non-endotoxin sources), there is insufficient data to support this broader
137 application.

138

¹ Equivalent methods can be regulated under 21 CFR 610.9 as alternatives to the currently accepted test method(s).

² There are substances other than endotoxin that may induce the cellular release of IL-1 β and/or IL-6. For this reason, users of these test methods should be aware of the potential for a false positive result, suggesting that endotoxin is present, which actually is due to the presence of another pyrogenic material.

139

140 **Table 1 Performance Statistics for *In Vitro* Pyrogenicity Test Methods¹**

Test Method	Concordance ²	Sensitivity	Specificity	False Negative Rate	False Positive Rate
PBMC/IL-6	93.3% (140/150)	92.2% (83/90)	95.0% (57/60)	7.8% (7/90)	5.0% (3/60)
cryo WB/IL-1	91.7% (110/120)	97.4% (75/77)	81.4% (35/43)	2.6% (2/77)	18.6% (8/43)
WB/IL-6	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
MM6/IL-6	93.2% (138/148)	95.5% (85/89)	89.8% (53/59)	4.5% (4/89)	10.2% (6/59)
WB/IL-1 (plate method)	92.0% (129/139)	98.8% (83/84)	83.6% (46/55)	1.2% (1/84)	16.4% (9/55)

141

¹Based on combined results of 10 different parenteral drugs tested in each of three different laboratories; samples of each drug were tested with or without being spiked with a Gram-negative endotoxin standard (0, 0.25, 0.5, 0.5, and 1.0 EU/mL).

142

²Percentage (Number of correct runs/total number of runs)

143

144

145

146 **Table 2 Test Substances (Parenteral Drugs) Used in the Validation Studies for**
 147 **Determining Test Method Performance¹**

Test Substance ²	Source	Active Ingredient	Indication	MVD (-fold)
Beloc®	Astra Zeneca	Metoprolol tartrate	Heart dysfunction	140
Binotal®	Aventis	Ampicillin	Antibiotic	140
Ethanol 13% (w/w)	B. Braun	Ethanol	Diluent	35
Fenistil®	Novartis	Dimetindenmaleat	Antiallergic	175
Glucose 5% (w/v)	Eifel	Glucose	Nutrition	70
MCP®	Hexal	Metoclopramid	Antiemetic	350
Orasthin®	Aventis	Oxytocin	Initiation of Delivery	700
Sostril®	GSK	Ranitidine	Antiacidic	140
Drug A - 0.9% NaCl	-	0.9% NaCl	-	35
Drug B - 0.9% NaCl	-	0.9% NaCl	-	70

148

¹Each substance was tested in all five *in vitro* pyrogenicity test methods.

149

²Each test substance was spiked with 0, 0.25, 0.5, 0.5, 1.0 EU/mL of endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]). Each sample contained the appropriate spike concentration when tested at its Maximum Valid Dilution (MVD).

150

151

152

153 Users should be aware that the performance characteristics for these *in vitro* pyrogenicity
 154 test methods could be revised as additional data become available. Therefore, test method
 155 users should consult the ICCVAM/NICEATM website (<http://iccvam.niehs.nih.gov/>) and
 156 other sources to obtain the most current information relevant to the current performance
 157 and usefulness and limitations of these test methods.

158

159

159 **1.2 Draft Recommended Future Studies**

160

161 To further the use of these five test methods (cryo WB/IL-1, MM6/IL-6, PBMC/IL-6,
162 WB/IL-1 [plate method], and WB/IL-6) as potential replacements for the RPT for
163 detecting non-endotoxin pyrogens, additional studies that include a broader range of
164 pyrogenic materials are recommended. For a direct comparison between the *in vitro*
165 pyrogen test(s) and the RPT, such studies should include parallel RPT testing³.

166

167 The hazards associated with human blood products should be carefully considered and all
168 technical staff must be trained to observe all necessary safety precautions.

169

170 **Appendix A** provides Draft Performance Standards for *In Vitro* Pyrogenicity Test
171 Methods that are based on ICCVAM guidelines (ICCVAM 2003⁴). **Appendix B** provides
172 five draft proposed *in vitro* pyrogenicity test method protocols that are based on those
173 used in the ECVAM validation study. **Appendix B1** is the *Proposed Test Method*
174 *Protocol for the Human Whole Blood/IL-1 In Vitro Pyrogen Test: Application of*
175 *Cryopreserved Human Whole Blood*. **Appendix B2** is the *Proposed Test Method*
176 *Protocol for the In Vitro Pyrogen Test Using the Human Monocytoid Cell Line MONO*
177 *MAC-6 (MM6)*. **Appendix B3** is the *Proposed Test Method Protocol for the Human*
178 *PBMC/IL-6 In Vitro Pyrogen Test*. **Appendix B4** is the *Proposed Test Method Protocol*
179 *for the Human Whole Blood/IL-1 In Vitro Pyrogen Test*. **Appendix B5** is the *Proposed*
180 *Test Method Protocol for the Human Whole Blood/IL-6 In Vitro Pyrogen Test*.

³ In order to demonstrate the utility of these test methods for the detection of non-endotoxin pyrogens, either an international standard is needed (as is available for endotoxin [i.e., WHO-LPS 94/580 *E. coli* O113:H10:K-]) or, when a positive non-endotoxin-mediated RPT result is encountered, this same sample should be subsequently tested *in vitro*.

⁴ ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised, and Alternative Test Methods. NIH Publication No. 03-4508. National Institute for Environmental Health Sciences, Research Triangle Park, NC. Available: <http://iccvam.niehs.nih.gov/>. [accessed 2 June 2005].

181

182

183

184

185

186

187

188

189

190

191

192

193

194

195

[This Page Intentionally Left Blank]

195

196

197

198

199

200

201

Appendix A

202

203

204

Draft Performance Standards for *In Vitro* Pyrogenicity Test Methods

205

205

206

207

208

209

210

211

212

213

214

215

216

217

218

219

[This Page Intentionally Left Blank]

220 1.0 PURPOSE AND BACKGROUND OF PERFORMANCE STANDARDS

221

222 1.1 Introduction

223

224 Prior to the acceptance of a new test method for regulatory testing applications, validation
225 studies are conducted to assess its reliability (i.e., the extent of intra- and inter-laboratory
226 reproducibility) and its relevance (i.e., the ability of the test method to correctly predict or
227 measure the biological effect of interest) (OECD 1996, 2002; ICCVAM 1997, 2003). The
228 purpose of performance standards is to communicate the basis by which new proprietary
229 (i.e., copyrighted, trademarked, registered) and nonproprietary test methods have been
230 determined to have sufficient relevance and reliability for specific testing purposes. These
231 performance standards, based on test methods accepted by regulatory agencies, can be
232 used to evaluate the reliability and relevance of other test methods that are based on
233 similar scientific principles and measure or predict the same biological or toxic effect.
234 Five *in vitro* pyrogenicity test methods, cryo WB/IL-1 (The Human Whole Blood/IL-1 *In*
235 *Vitro* Pyrogen Test: Application of cryopreserved human whole blood), MM6/IL6 (An
236 Alternative *In Vitro* Pyrogen Test Using the Human Monocytoid Cell Line MONO
237 MAC-6 [MM6]), PBMC/IL-6 (The Human PBMC/IL-6 *In Vitro* Pyrogen Test, WB/IL-1
238 (The Human Whole Blood/IL-1 *In Vitro* Pyrogen Test), and WB/IL-6 (The Human
239 Whole Blood/IL-6 *In Vitro* Pyrogen Test) were included in a validation study to evaluate
240 the correlation between *in vitro* cytokine release and the *in vivo* rabbit fever response and
241 the feasibility of using *in vitro* cytokine assays to predict a pyrogenic response.

242

243 This section describes the three elements of performance standards identified by
244 ICCVAM (2003) and the ICCVAM process used to develop performance standards
245 during a test method evaluation. These test method performance standards are proposed
246 as standards that can be used to evaluate future *in vitro* pyrogenicity test methods. If other
247 *in vitro* pyrogenicity test methods are adequately validated and demonstrate significantly
248 improved performance, then the test method performance standards may be revised
249 accordingly.

250

251 1.2 Elements of ICCVAM Performance Standards

252

253 Performance standards are standards based on a validated test method that provide a basis
254 for evaluating the comparability of a proposed test method that is mechanistically and
255 functionally similar (ICCVAM 2003). The three elements of performance standards are:

- 256 • **Essential test method components:** These consist of essential structural,
257 functional, and procedural elements of a validated test method that should
258 be included in the protocol of a proposed, mechanistically and functionally
259 similar test method. Essential test method components include unique
260 characteristics of the test method, critical procedural details, and quality
261 control measures.
- 262 • **A minimum list of reference substances:** Reference substances are used to
263 assess the accuracy and reliability of a proposed, mechanistically and
264 functionally similar test method. These substances are a representative
265 subset of those used to demonstrate the reliability and the accuracy of the
266 validated test method, and are the minimum number that should be used to
267 evaluate the performance of a proposed, mechanistically and functionally
268 similar test method.
- 269 • **Accuracy and reliability values:** These are the accuracy and reliability
270 characteristics that the proposed test method should be comparable to or
271 exceed when evaluated using the minimum list of reference chemicals.

272

273 1.3 ICCVAM Process for the Development of Performance Standards

274

275 The process followed by ICCVAM for developing performance standards for new test
276 methods is as follows:

- 277 • NICEATM and the ICCVAM Pyrogenicity Working Group (PWG) develop
278 proposed performance standards for consideration during the ICCVAM
279 evaluation process. If performance standards are proposed by a test method
280 sponsor, they will be considered by ICCVAM at this stage. Generally, the

- 281 proposed performance standards are based on the information and data
282 provided in the test method submission or on other available applicable data.
- 283 • The ICCVAM/NICEATM Peer Review Panel evaluates the proposed
284 performance standards for completeness and appropriateness during its
285 evaluation of the validation status of the proposed test method. The
286 proposed performance standards, as well as the test method submission, are
287 made available to the public for comment prior to and during the Peer
288 Review Panel meeting.
 - 289 • The PWG, with the assistance of NICEATM, prepares the final performance
290 standards for ICCVAM approval, taking into consideration the
291 recommendations of the Peer Review Panel and public comments.

292

293 Performance standards recommended by ICCVAM are incorporated into ICCVAM test
294 method evaluation reports, which are then provided to U.S. Federal agencies and made
295 available to the public. Regulatory authorities can then reference the performance
296 standards in the ICCVAM report when they communicate their acceptance of a new test
297 method. In addition, performance standards adopted by U.S. Federal regulatory
298 authorities can be provided in guidelines issued for new test methods. Availability of
299 ICCVAM test method evaluation reports are announced routinely in the *Federal Register*,
300 in NTP Newsletters, and by e-mail to ICCVAM/NICEATM listserve groups.

301

302 **1.4 ICCVAM Development of Performance Standards for *In Vitro*** 303 **Pyrogenicity Test Methods**

304

305 **1.4.1 Current Regulatory Testing Requirements for Pyrogenicity**

306 The major regulatory requirement for pyrogenicity testing is for end product release of
307 human and animal parenteral drugs, medical devices, and human biological products.
308 Results from pyrogenicity testing are used to limit to an acceptable level the risks of
309 febrile reaction in the patient to injection and/or implantation of the product of concern.
310 The current U.S. legislation requiring the use of pyrogenicity testing is stated in the
311 Federal Food, Drug, and Cosmetic Act (U.S.C., Title 21, Chapter 9). In addition, the U.S.

312 Pharmacopeia (USP) maintains sterility requirements for pharmaceuticals that include
 313 pyrogenicity testing. As detailed in **Table 1-1**, the U.S. Food and Drug Administration
 314 (FDA) is the principal U.S. regulatory agency that requires pyrogenicity testing, with
 315 different Centers within the FDA regulating the affected products (i.e., human and animal
 316 parenteral drugs, biological products, and medical devices). **Table 1-1** also shows the
 317 statutory protocol requirements used by each FDA Center, along with the comparable
 318 international standards presently required by European Union member nations.

319

320 **Table 1-1 Summary of U.S. and European Legislation and Statutory Protocol**
 321 **Requirements for Pyrogenicity Testing**

Agency	Regulated Products	Legislation	Statutory Protocol Requirements	Non-Governmental Standards
United States				
FDA-CBER	Biological products	- Federal Food, Drug, and Cosmetic Act (U.S.C. Title 21, Chapter 9)	- 21 CFR 610.13	- USP28 NF23<85> - USP28 NF23<151> - ISO 10993-11
FDA-CDER	Human parenteral pharmaceuticals			
FDA-CDRH	Medical devices			
FDA-CVM	Veterinary pharmaceuticals			
Europe				
EDQM	Human/veterinary parenteral pharmaceuticals, biological products, medical devices	- Council Regulation (EEC) 230/9/93	- EP5.0 2.6.8 - EP5.0 2.6.14	- ISO 10993-11
EMA		- Council Directive 93/39/EEC		
Regulatory Authorities for Individual EU Countries		- Council Directive 93/40/EEC		

322 Abbreviations: CBER = Center for Biologics Evaluation and Research; CDER = Center for Drug Evaluation and
 323 Research; CDRH = Center for Devices and Radiological Health; CFR = Code of Federal Regulations; CVM = Center
 324 for Veterinary Medicine; EDQM = European Directorate for the Quality of Medicines; EMA = European Medicines
 325 Agency; EP = European Pharmacopoeia; EU = European Union; FDA = U.S. Food and Drug Administration; US =
 326 U.S. Pharmacopeia

327

328 1.4.2 Test Methods for Assessing Pyrogenicity

329 The currently recognized test methods for evaluating pyrogenicity are the *in vivo* rabbit
 330 pyrogen test (RPT) and the bacterial endotoxin test (BET). The RPT (USP28
 331 NF23<151>, EP5.0 2.6.8) involves measuring the rise in body temperature evoked in
 332 rabbits by the intravenous injection of a test solution. The RPT is a sequential test, using
 333 the response of the first three rabbits tested to determine the need for additional testing.
 334 The BET (USP28 NF23<85>, EP5.0 2.6.14) is used to detect or quantify the presence of

335 gram-negative bacterial endotoxins using amoebocyte lysate from the horseshoe crab
336 (*Limulus polyphemus*). There are three different BET techniques: the gel-clot technique
337 (based on gel formation due to clotting of *Limulus* amoebocyte lysate [LAL] in the
338 presence of endotoxin); the turbidometric technique (based on the development of
339 turbidity after cleavage of an endogenous substrate); and the chromogenic technique
340 (based on color development resulting from cleavage of a synthetic peptide-chromagen
341 complex).

342

343 1.4.3 Intended Regulatory Uses for *In Vitro* Pyrogenicity Test Methods

344 The *in vitro* pyrogenicity test methods are not intended as replacements for the BET.
345 However, five of these methods (cryo WB/IL-1, MM6/IL-6, PBMC/IL-6, WB/IL-1 [plate
346 method], and WB/IL-6) may be considered for the detection of the presence of Gram-
347 negative endotoxin in materials that are currently tested in the RPT, subject to product-
348 specific validation to demonstrate equivalency.

349

350 1.4.4 Similarities and Differences in the Endpoints of *In Vitro* Pyrogenicity Test 351 Methods and Currently Recognized Pyrogenicity Test Methods

352 The endpoint measured in the *in vitro* pyrogenicity test methods is cytokine release,
353 either IL-1 β or IL-6, depending on the test method employed. The RPT involves
354 measuring the rise in body temperature evoked in rabbits by the intravenous injection of a
355 test solution. While there is not a direct association between the endpoints measured in
356 these assays, cytokine release is involved in the development of an inflammatory
357 response, which can result in an increase in temperature. Therefore, the *in vitro* release of
358 pro-inflammatory cytokines, such as IL-1 β and IL-6, is intended to predict the onset of
359 such a response. The cell types used for the various *in vitro* methods presumably are
360 those that would be directly associated with an inflammatory response. Both the *in vitro*
361 and *in vivo* tests provide quantitative endpoints.

362

363 There are different endpoints for the BET, depending on the technique used. The gel-clot
364 technique is based on the observation of gel formation due to clotting of LAL in the
365 presence of endotoxin. The turbidometric technique evaluates the development of turbidity

366 after cleavage of an endogenous substrate. Finally, the chromogenic technique measures
367 color development resulting from cleavage of a synthetic peptide-chromogen complex.
368 Clearly, there are no biological similarities between the endpoints measured in the *in*
369 *vitro* test methods and the various BET techniques. However, like the *in vitro* test
370 methods, the turbidometric and chromogenic techniques provide quantitative
371 measurements, while the gel-clot technique is qualitative.

372

373 **2.0 IN VITRO PYROGENICITY TEST METHODS**

374

375 **2.1 Background**

376

377 Pre-validation and validation studies have been completed to evaluate the ability of the
378 five *in vitro* pyrogenicity test methods to be used as alternatives to the RPT. This section
379 briefly describes the principles of *in vitro* pyrogenicity test methods followed by the
380 recommended performance standards that would be used to evaluate test methods that are
381 functionally and mechanistically similar to these methods. The performance standards
382 consist of 1) essential test method components, 2) reference substances, and 3) the
383 comparable accuracy and reliability that should be achieved.

384

385 **2.2 Principles of *In Vitro* Pyrogenicity Test Methods**

386

387 Although there are differences among the *in vitro* pyrogenicity test methods based
388 predominantly on the cell type used, there are some basic steps that are consistent across
389 all methods as follows:

- 390 • The test substance is applied to the specific human-derived cells used in
391 the *in vitro* test method (i.e., mixed with a suspension of cells).
- 392 • The test substance is incubated with the cells for a specified period of
393 time.
- 394 • The concentration of pro-inflammatory cytokines (e.g., IL-1 β , IL-6) is
395 quantified via a cytokine-specific enzyme-linked immunosorbent assay
396 (ELISA) by comparison to a standard curve.

- 397 • Using an endotoxin standard curve, the endotoxin content of the product
398 is calculated.
- 399 • A product “passes” (i.e., is considered negative for endotoxin) if the
400 endotoxin content is < 0.5 endotoxin units (EU)/mL.

401

402 While it is possible that comparable cell types derived from other species would also be
403 capable of detecting the presence of Gram-negative endotoxin, one of the strengths of
404 these test methods is that they are derived from human tissues, and thus avoid the
405 potential uncertainty associated with cross-species interpretation. Therefore, test method
406 developers are encouraged to focus on human cell-based systems.

407

408 Investigators using an *in vitro* pyrogenicity test method for detecting the presence of
409 Gram-negative endotoxin must be able to demonstrate that the assay is valid for its
410 intended use. This includes demonstrating that any modification to the existing validated
411 reference test method does not adversely affect its performance characteristics. *In vitro*
412 pyrogenicity test methods may be used to test pharmaceuticals, biological products, and
413 medical devices. Interference testing must be included to demonstrate that the properties
414 of the test substance do not impede the release and detection of proinflammatory
415 cytokines.

416

417 **2.3 Essential Test Method Components for *In Vitro* Pyrogenicity Test**

418 **Methods**

419

420 Essential test method components consist of essential structural, functional, and
421 procedural elements of a validated test method that should be included in the protocol of
422 a mechanistically and functionally similar proposed test method. These components
423 include unique characteristics of the test method, critical procedural details, and quality
424 control measures. Adherence to these components will help assure that a proposed test
425 method is based on the same concepts as the corresponding validated test method.

426

427 The following is a description of the essential test method components for *in vitro*
428 pyrogenicity test methods.

429

430 2.3.1 *In Vitro* Cell Culture Conditions

- 431 • A mammalian cell line, primary cells, or heparinized whole blood (as
432 described above, either would preferably be of human origin) is used.
433 Cryopreserved cells/whole blood may be used, where it has been
434 demonstrated that cryopreservation is not detrimental to the test method.
- 435 • Fresh whole blood may be stored at room temperature, but should be used
436 within four hours of collection. Blood donors should be in good health
437 (i.e., not suffering from bacterial or viral infections for at least one week
438 prior to donation), and not taking any medications known to influence
439 cytokine production (e.g., immunosuppressant or anti-inflammatory
440 drugs). As an additional measure, assay acceptance criteria allow for the
441 identification of low or high responders such that results with
442 Blood/PBMCs from a compromised donor are omitted.
- 443 • Where necessary, cells are propagated in sterile tissue culture flasks and
444 then subcultured to sterile 96-well plates for use in testing. Initial cell
445 seeding should be done at a density that allows rapid growth throughout
446 the exposure period. However, cell density should not reach confluency by
447 the end of the test exposure period.
- 448 • Appropriate cell culture growth conditions (e.g., 37°C ± 1°C, 90% ± 10%
449 humidity, 5.0% ± 1% CO₂ in ambient air) should be maintained
450 throughout the testing period. The cell cultures should be free of
451 contamination with bacteria, mycoplasma, or fungi.

452

453 Cell culture media should be prequalified by the testing laboratory via a standardized
454 protocol before initiating the test to guarantee that the media provide cells with
455 appropriate nutrients to meet the growth criteria needed for the test method.

456

457

457 2.3.2. Application of the Test Substances458 2.3.2.1 *Test Substance Preparation*

- 459 • All disposables (e.g., pipette tips, pipettes, culture ware, etc.) should be
460 labeled sterile, pyrogen free.
- 461 • Test substances (i.e., pharmaceuticals, biological products) should be
462 diluted to their respective maximum valid dilution (MVD) in sterile,
463 pyrogen-free 0.9% NaCl.
- 464 • Medical devices can be directly incubated with to the cells in suspension.
465 Alternatively, and if necessary, eluates/extracts from medical devices may
466 be prepared with a volume of pyrogen-free water appropriate to their use
467 and, where applicable, to the surface area that comes in contact with body
468 tissues or fluids.
- 469 • Each test should contain a range of concentrations of either the
470 international reference standard endotoxin (i.e., WHO LPS 94/580), or an
471 LPS standard that has been calibrated against this standard, with which to
472 generate a standard curve. *NOTE: In the ECVAM validation study,*
473 *concentrations of 0, 0.25, 0.5, and 1.0 endotoxin units (EU)/mL were used*
474 *to establish the decision criteria for a pyrogenic response based on the in*
475 *vivo rabbit threshold fever concentration (i.e., the concentration at which*
476 *an increase in temperature was recorded in 50% animals tested) (see*
477 **Section 2.3.5).**
- 478 • Test substances should be fully solubilized (i.e., no visual observation of
479 test substance in the dosing solution).

480

481 2.3.2.2 *Test Substance Application and Sample Collection*

- 482 • Whole blood samples may be dosed in either 96-well plates or
483 microcentrifuge tubes.
- 484 • The cells should be exposed for from 16 to 24 hours
- 485 • Each substance should be tested in a minimum of three replicates.

- 486 • At the end of the exposure period, supernatants may be collected either
487 directly from each well, or following centrifugation for microcentrifuge
488 tubes.

489

490 2.3.3 Control Substances

491 2.3.3.1 *Negative Control*

492 To ensure that the test system is functioning properly and that the specific test is valid,
493 the negative control (i.e., 0.9% NaCl) should not induce a significant increase in IL-1 β or
494 IL-6 release.

495

496 2.3.3.2 *Positive Control*

497 The purpose of a positive control chemical is to demonstrate that the cell system is
498 responding with adequate sensitivity to a pyrogenic substance for which the magnitude of
499 the pyrogenic response is well characterized. Each test should generate a response that is
500 comparable to the historical range generated by the laboratory. Therefore, the positive
501 control should be the international reference standard endotoxin (i.e., WHO-LPS 94/580
502 [*E. coli* 0113:h10:K-]), or an endotoxin standard that has been calibrated against this
503 standard. A laboratory should perform a minimum of 10 *in vitro* pyrogenicity tests using
504 the positive control over a number of days to develop a minimum historical database of
505 cytokine data. Typically, for biologically based test methods, suggested acceptable ranges
506 for the positive control response are within two to three standard deviations of the
507 historical mean response, but developers of proprietary test methods may establish tighter
508 ranges. The positive control chemical should be tested concurrently with (and
509 independent of) the test substance. Test substances spiked with known quantities of the
510 positive control should be used for interference testing.

511

512 2.3.3.3 *Benchmark Controls*

513 Benchmark controls may be useful to demonstrate that the test method is functioning
514 properly for detecting the pyrogenic potential of chemicals (e.g., parenterals or medical
515 device eluates) of a specific chemical class or a specific range of responses, or for

516 evaluating the relative pyrogenic potential of a test substance. Appropriate benchmark
517 controls should have the following properties:

- 518 • structural and functional similarity to the class of the substance being tested
- 519 • known physical/chemical characteristics
- 520 • supporting data on known effects in animal models
- 521 • known potency in the range of response

522

523 2.3.4 Cytokine Measurements

524 Only standardized, quantitative methods should be used to measure cytokine release (e.g.,
525 an enzyme immunoassay with a species-specific antibody for the relevant
526 proinflammatory cytokine). Each assay should contain a range of concentrations of the
527 relevant cytokine standard (e.g., IL-1 β , IL-6) in order to generate a standard curve for the
528 analytical assay. The protocol should be compatible with analytical laboratory equipment
529 (e.g., spectrophotometer) that allows a quick and precise measurement of the endpoint.
530 Colorimetric, fluorometric, or luminometric endpoints should have the optical density
531 (OD) measured at the appropriate wavelength, and OD values for blanks should be
532 subtracted from all measurements. Each supernatant should be assayed with a minimum
533 of three replicates.

534

535 2.3.5 Interpretation of Results

536 The endpoint values obtained for each test sample can be used to calculate the level of
537 cytokine release relative to the positive control samples (i.e., the endotoxin standard
538 curve). A sample is considered positive for a pyrogenic response if the level of cytokine
539 release is greater than or equal to that induced by the 0.5 EU/mL endotoxin standard, the
540 reported threshold fever concentration for the *in vivo* rabbit test (see **Section 2.3.2.1**).

541

542 2.3.6 Test Report

543 The test report should include the following information, if relevant to the conduct of the
544 study:

545 *Test Substances and Control Substances*

- 546 • Name of pharmaceutical, biological product, medical device eluate, etc.

- 547 • Purity and composition of the substance or preparation
- 548 • Physicochemical properties (e.g., physical state, water solubility) relevant to the
- 549 conduct of the study
- 550 • Treatment of the test/control substances prior to testing, if applicable (e.g.,
- 551 vortexing, sonication, warming; resuspension solvent)

552 *Justification of the In Vitro Test Method and Protocol Used*

553 *Test Method Integrity*

- 554 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
- 555 test method over time
- 556 • If the test method employs proprietary components, documentation on the
- 557 procedure used to ensure their integrity from “lot-to-lot” and over time
- 558 • The procedures that the user may employ to verify the integrity of the proprietary
- 559 components

560 *Criteria for an Acceptable Test*

- 561 • Acceptable concurrent positive control ranges based on historical data
- 562 • Acceptable negative control data

563 *Test Conditions*

- 564 • Cell system used
- 565 • Calibration information for the spectrophotometer used to read the ELISA
- 566 • Details of test procedure
- 567 • Description of any modifications of the test procedure
- 568 • Reference to historical data of the model
- 569 • Description of evaluation criteria used

570 *Results*

- 571 • Tabulation of data from individual test samples

572 *Description of Other Effects Observed*

573 *Discussion of the Results*

574 *Conclusion*

575 *A Quality Assurance Statement for Good Laboratory Practice (GLP)-Compliant Studies*

- 576 • This statement should indicate all inspections made during the study and the dates
577 any results were reported to the Study Director. This statement should also
578 confirm that the final report reflects the raw data

579 If GLP-compliant studies are performed, then additional reporting requirements provided
580 in the relevant guidelines (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003) should be
581 followed.

582

583 **2.4 Reference Substances for *In Vitro* Pyrogenicity Test Methods**

584

585 Reference substances are used to assess the accuracy and reliability of a proposed,
586 mechanistically and functionally similar test method and are a representative subset of
587 those used to demonstrate the reliability and the accuracy of the validated test method.

588 These substances are:

- 589 • representative of the range of responses that the validated test method is
590 capable of measuring or predicting
- 591 • have produced consistent results in the validated test method
- 592 • reflect the accuracy of the validated test method
- 593 • have well-defined chemical structures
- 594 • are readily available
- 595 • are not associated with excessive hazard or prohibitive disposal costs

596

597 To demonstrate technical proficiency with the validated test method, the user should
598 evaluate his/her ability to identify Gram-negative endotoxin that has been spiked into
599 each of the reference substances listed in **Table 2-1**. These eight substances are marketed
600 pharmaceuticals that were tested in the ECVAM *in vitro* pyrogenicity test methods
601 validation study. Only released clinical lots that have been labeled as having no
602 detectable pyrogens should be used as reference substances to be spiked with Gram-
603 negative endotoxin. As indicated in **Section 2.3.3**, the spike should be either the
604 international reference standard endotoxin (WHO-LPS 94/580 [*E. coli* 0113:h10:K-]) or
605 an endotoxin standard that has been calibrated against this standard. Each reference
606 substance should be tested clean (i.e., unspiked) and spiked with endotoxin (0.5 EU/mL).

607

608 **Table 2-1 Recommended Reference Substances for *In Vitro* Pyrogenicity Test**
 609 **Methods**

Test Substance ¹	Source	Active Ingredient	MVD (-fold)
Beloc®	Astra Zeneca	Metoprolol tartrate	140
Binotal®	Aventis	Ampicillin	140
Ethanol 13% (w/w)	B. Braun	Ethanol	35
Fenistil®	Novartis	Dimetindenmaleat	175
Glucose 5% (w/v)	Eifel	Glucose	70
MCP®	Hexal	Metoclopramid	350
Orasthin®	Aventis	Oxytocin	700
Sostril®	GSK	Ranitidine	140

610 ¹Each reference substance should be spiked with 0.5 EU/mL of endotoxin (WHO-LPS 94/580 [*E. coli*
 611 O113:H10:K-]). Each sample should contain the appropriate spike concentration when tested at its
 612 Maximum Valid Dilution (MVD). MVD = (endotoxin limit concentration)/(detection limit of the assay)
 613 Positive Control: 0.9% NaCl spiked with 0.5 EU/mL
 614 Negative Control: 0.9% NaCl
 615

616 2.5 Accuracy and Reliability

617

618 The third element of the performance standards is the determination of accuracy (also
 619 known as relevance) and reliability values. A proposed test method, functionally and
 620 mechanistically similar to the *in vitro* pyrogenicity test methods described above, will use
 621 selected reference substances to assess accuracy and reliability.

622

623 2.5.1 Accuracy

624 When evaluated using the minimum list of recommended reference substances (**Table 2-**
 625 **1**), the proposed test method should have performance characteristics that are comparable
 626 to the performance of the validated *in vitro* pyrogenicity test methods. Accuracy is
 627 defined as the closeness of agreement between a test method result and an accepted
 628 reference value (ICCVAM 2003). The substances tested in the ECVAM validation
 629 studies are included so that the performance of the proposed test method can be
 630 determined and compared to that of the validated reference test methods.

631

632 The accuracy of these assays to identify a pyrogenic concentration of Gram-negative
 633 endotoxin was evaluated. This accuracy evaluation characterizes the extent that
 634 additional test methods will be necessary to achieve accurate *in vitro* predictions of

635 contamination by Gram-negative endotoxin for labeling and lot release purposes. **Table**
 636 **2-2** shows that overall accuracy among the test methods is comparable (91.7% to 93.3%),
 637 with false negative rates ranging from 1.2% to 11.2%, and false positive rates ranging
 638 from 3.4% to 18.6%.

639

640

641 **Table 2-2 Performance Statistics for *In Vitro* Pyrogenicity Test Methods¹**

Test Method	Accuracy ²	Sensitivity	Specificity	False Negative Rate	False Positive Rate
PBMC/IL-6	93.3% (140/150)	92.2% (83/90)	95.0% (57/60)	7.8% (7/90)	5.0% (3/60)
cryo WB/IL-1	91.7% (110/120)	97.4% (75/77)	81.4% (35/43)	2.6% (2/77)	18.6% (8/43)
WB/IL-6	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
MM6/IL-6	93.2% (138/148)	95.5% (85/89)	89.8% (53/59)	4.5% (4/89)	10.2% (6/59)
WB/IL-1 (plate method)	92.8% (129/139)	98.8% (83/84)	83.6% (46/55)	1.2% (1/84)	16.4% (9/55)

642

643

644

¹Based on combined results of 10 different substances tested in three different laboratories

²Percentage (Number of correct runs/total number of runs)

645

2.5.2 Reliability

646

647

648

649

650

651

652

653

654

655

656

657

658

Test method reliability (intralaboratory repeatability and intra- and inter-laboratory reproducibility) is the degree to which a test method can be performed reproducibly within and among laboratories over time (ICCVAM 2003). Repeatability refers to the closeness of agreement between test results obtained within a single laboratory when the procedure is performed on the same substance under identical conditions within a given time period. Intralaboratory reproducibility refers to the determination of the extent to which qualified personnel within the same laboratory can replicate results using a specific test protocol at different times. Interlaboratory reproducibility refers to the determination of the extent to which different laboratories can replicate results using the same protocol and test chemicals, and indicates the extent to which a test method can be transferred successfully among laboratories.

659 The reliability of the proposed test method for the reference substances should be
660 comparable to or better than that of the validated *in vitro* pyrogenicity test methods. The
661 following sections provide these reference statistics.

662

663 2.5.2.1 *Intralaboratory Repeatability*

664 In the ECVAM validation study, intralaboratory repeatability was evaluated in each test
665 method by testing saline and various endotoxin spikes (0.06 to 0.5 EU/mL) in saline and
666 evaluating the closeness of agreement between optical density readings for cytokine
667 measurements for each concentration. Up to 20 replicates per concentration were tested
668 and results indicated that variability in measurements increased with endotoxin
669 concentration, but that the 0.5 EU/mL spike concentration (i.e., the threshold for
670 pyrogenicity for each test method) was clearly distinguishable from lower concentrations.

671

672 2.5.2.2 *Intralaboratory Reproducibility*

673 Intralaboratory reproducibility was evaluated with three marketed pharmaceuticals spiked
674 with various concentrations of endotoxin. Reproducibility was assessed from three
675 identical, independent runs conducted in each of the three testing laboratories (with the
676 exception of the cryo WB/IL-1 test method¹). From these results, agreement between
677 different runs was determined for each substance in three laboratories. As shown in
678 **Table 2-3**, the agreement across three runs in an individual lab ranged from 75% to
679 100%.

680

¹ The cryo WB/IL-1 test method BRD states that an assessment of intralaboratory reproducibility was performed in the WB IL-1 (fresh blood) test method, and it is assumed that variability is not affected by the change to cryopreserved blood assayed in 96-well plates.

581 **Table 2-3 Intralaboratory Reproducibility of *In Vitro* Pyrogenicity Test Methods**

Run Comparison ¹	WB/IL-1			Cryo WB/IL-1			WB/IL-6			PBMC/IL-6			MM6/IL-6		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
1 vs 2	92%	100%	100%	ND	ND	ND	75%	92%	100%	92%	100%	100%	100%	92%	100%
1 vs 3	83%	88%	92%	ND	ND	ND	100%	92%	100%	100%	100%	92%	100%	92%	92%
2 vs 3	92%	NA	92%	ND	ND	ND	75%	100%	100%	92%	100%	92%	100%	100%	92%
Mean	89%	-	94%	ND	ND	ND	83%	94%	100%	94%	100%	94%	100%	94%	94%
Agreement across 3 runs	83%	-	92%	ND	ND	ND	75%	92%	100%	92%	100%	94%	100%	92%	92%

582 NA: Not assessed due to lack of sufficient data. The sensitivity criteria were not met for 1/3 substances in run 2, and 1/3 substances in run 3; ND: Not done. The cryo WB/IL-1 test method BRD states
 583 that an assessment of intralaboratory reproducibility was performed in the WB IL-1 (fresh blood) test method, and it is assumed that variability is not affected by the change to cryopreserved blood
 584 assayed in 96-well plates.

585 ¹Comparison between 3 individual runs within each laboratory

586

587

688 2.5.2.3 *Interlaboratory Reproducibility*

689 Interlaboratory reproducibility was evaluated in two different studies. In both studies,
 690 each run from one laboratory was compared with all other runs of another laboratory. The
 691 proportion of equally classified samples provides a measure of reproducibility. In the first
 692 study, like the intralaboratory reproducibility evaluation, interlaboratory reproducibility
 693 was evaluated with the three marketed pharmaceuticals spiked with endotoxin, and tested
 694 three times in three different laboratories. As shown in **Table 2-4**, the agreement across
 695 three laboratories for each test method (where three runs per laboratory were conducted)
 696 ranged from 72% to 86%, depending on the test method considered. In comparison, the
 697 agreement across three laboratories for the WB/IL-1 (96-well plate method) and cryo
 698 WB/IL-1 test methods, for which only one run per laboratory was conducted, was 83%
 699 and 92%, respectively.

700

701 **Table 2-4 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test**702 **Methods**

Lab Comparison ¹	Agreement Between Laboratories ¹				
	WB/IL-1	Cryo WB/IL-1	WB/IL-6	PBMC/IL-6	MM6/IL-6
1 vs 2	92% (11/12)	92% (11/12)	72% (78/108) ²	81% (87/108)	97% (105/108)
1 vs 3	83% (10/12)	92% (11/12)	75% (81/108) ²	86% (93/108)	89% (96/108)
2 vs 3	92% (11/12)	92% (11/12)	97% (105/108) ²	89% (96/108)	86% (93/108)
Mean	89%	92%	81%	85%	90%
Agreement across 3 labs ³	83% (10/12)	92% (11/12)	72% (234/324) ²	78% (252/324)	86% (279/324)

703 ¹Data from three substances (see **Table 2-3**) spiked with endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]) at 0, 0, 0.5 and 1.0
 704 EU/mL tested three times in three different laboratories, with the exception of the WB/IL-1 (96-well plate method) and the cryo
 705 WB/IL-1 (only the preliminary run from each laboratory used for analysis)

706 ²Some of the runs did not meet the assay acceptance criteria and therefore were excluded from the analysis.

707 ³All possible combinations of runs among the 3 laboratories were compared (with the exception of cryo WB/IL-1, which was only
 708 tested once in each laboratory, resulting in only one possible combination per substance).
 709

710 In a second study, reproducibility was evaluated with the same ten substances used for
 711 evaluating accuracy. In this study, the ten substances were spiked with five
 712 concentrations of endotoxin and tested once in each of three laboratories. As indicated in
 713 **Table 2-5**, the agreement across three laboratories for each test method ranged from 79%
 714 to 88%, depending on the test method considered.

715

716

716 **Table 2-5 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test**
 717 **Methods**

Lab Comparison ¹	Agreement Between Laboratories ¹				
	WB/IL-1	Cryo WB/IL-1	WB/IL-6	PBMC/IL-6	MM6/IL-6
1 vs 2	88% (37/42)	84% (38/45)	85% (41/48)	84% (42/50)	90% (45/50)
1 vs 3	90% (35/39)	88% (21/24)	85% (41/48)	86% (43/50)	90% (43/48)
2 vs 3	92% (43/47)	100% (25/25)	88% (44/50)	90% (45/50)	83% (40/48)
Mean	90%	91%	86%	87%	88%
Agreement across 3 labs	85% (33/39)	88% (21/24)	79% (38/48)	80% (40/50)	81% (39/48)

718 ¹Data from 10 substances spiked with endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]) at 0, 0.25, 0.5, 0.5, and 1.0 EU/mL
 719 tested once in three different laboratories

720 **3.0 REFERENCES**

721

722 ICCVAM. 1997. Validation and Regulatory Acceptance of Toxicological Test Methods.
723 A Report of the *ad hoc* Interagency Coordinating Committee on the Validation of
724 Alternative Methods. NIH Publication No. 97-3981. Research Triangle Park, NC.
725 National Institute of Environmental Health Sciences. Available:
726 <http://iccvam.niehs.nih.gov/docs/guidelines/validate.pdf> [accessed 20 July 2006].

727

728

729 ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New,
730 Revised, and Alternative Test Methods. NIH Publication No. 03-4508. National Institute
731 for Environmental Health Sciences, Research Triangle Park, NC. Available:
732 <http://iccvam.niehs.nih.gov/>. [accessed 20 July 2006].

733

734 OECD. 1996. Final Report of the OECD Workshop on Harmonisation of Validation and
735 Acceptance Criteria for Alternative Toxicological Methods. (Solna Report) Paris:
736 Organisation for Economic Co-operation and Development. Available:
737 <http://www.oecd.org> [accessed 20 July 2006].

738

739 OECD. 2002. Report of the Stockholm Conference on Validation and Regulatory
740 Acceptance of New and Updated Methods in Hazard Assessment. Paris: Organisation for
741 Economic Co-operation and Development. Available: <http://www.oecd.org> [accessed 20
742 July 2006]

743

744

745

746

746
747
748
749
750
751
752
753
754
755
756
757
758
759
760
761
762
763
764
765

[This Page Intentionally Left Blank]

766
767
768
769
770
771
772
773
774
775
776
777
778
779
780
781
782
783
784
785

Appendix B

Individual BRDs Submitted by ECVAM on Five *In Vitro* Pyrogenicity Test Methods

B1 The Human Whole Blood (WB)/IL-1 *In Vitro* Pyrogen Test:
Application of Cryopreserved Human WB B1-1

**B2 An Alternative *In Vitro* Pyrogenicity Test Using the Monocytoid
Cell Line Mono Mac 6 (MM6)/IL-6 B2-1**

**B3 The Human Peripheral Blood Mononuclear Cell (PBMC)/IL-6
In Vitro Pyrogen Test..... B3-1**

B4 The Human WB/IL-1 *In Vitro* Pyrogen Test B4-1

B5 The Human WB/IL-6 *In Vitro* Pyrogen Test B5-1

785

786

787

788

789

790

791

792

793

794

795

796

797

798

799

[This Page Intentionally Left Blank]