

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY TEMPLATE**

A. 510(k) Number:

K043224

B. Purpose for Submission:

The APTIMA Combo 2 (AC2) assay is a nucleic acid amplification test (NAAT) intended for the qualitative detection and differentiation of ribosomal RNA from *Chlamydia trachomatis* (CT) and/or *Neisseria gonorrhoeae* (GC) in endocervical, male urethral and vaginal swab specimens and in female and male urine specimens. The assay originally received FDA clearance on May 21, 2001 (K032554). The current application is for the additional indication of testing specimens collected and processed with the Cytoc ThinPrep 2000 System. An ancillary kit called the “GEN-PROBE APTIMA Specimen Transfer Kit” is included in the current submission to facilitate specimen transport and processing.

C. Measurand:

Chlamydia trachomatis (CT) and/or *Neisseria gonorrhoeae* (GC) ribosomal RNA

D. Type of Test:

NAAT

E. Applicant:

Gen-Probe Incorporated

F. Proprietary and Established Names:

GEN-PROBE[®] APTIMA[®] Combo 2 Assay

G. Regulatory Information:

1. Regulation section:

CT: 21 CFR 866.3120

GC: 21 CFR 866.3390

2. Classification:

CT: Class I

GC: Class II

3. Product code:

CT: MKZ

GC: LSL

4. Panel:

Microbiology (83)

H. Intended Use:

1. Intended use(s):

AC2 Assay package insert:

The APTIMA Combo 2 Assay is a target amplification nucleic acid probe test that utilizes target capture for the in vitro qualitative detection and differentiation of ribosomal RNA (rRNA) from *Chlamydia trachomatis* (CT) and/or *Neisseria gonorrhoeae* (GC) in clinician-collected endocervical, vaginal, and male urethral swab specimens, patient-collected vaginal swab specimens*, and female and male urine specimens. The assay is also intended for use with testing of gynecological specimens collected in the PreservCyt Solution and processed with the Cytoc ThinPrep 2000 System. The assay may be used to test specimens from symptomatic and asymptomatic individuals to aid in the diagnosis of gonococcal and/or chlamydial urogenital disease.

*Patient-collected vaginal swab specimens are an option for screening women when a pelvic exam is not otherwise indicated. The vaginal swab specimen collection kit is not for home use.

Ancillary Kit package insert:

The GEN-PROBE® APTIMA® Specimen Transfer Kit is only for use with GEN-PROBE APTIMA assays for the detection of *Chlamydia trachomatis* and/or *Neisseria gonorrhoeae*. The GEN-PROBE APTIMA Specimen Transfer Kit allows for APTIMA Assay testing of gynecological specimens collected and processed by the Cytoc ThinPrep 200 Processor according to the instructions provided.

2. Indication(s) for use:

See intended use above

3. Special conditions for use statement(s):

This device is for prescription use only.

4. Special instrument requirements:

Gen-Probe Leader HC+ liminometer and Gen-Probe Target Capture System

I. Device Description:

The APTIMA Combo 2 (AC2) assay is a nucleic acid amplification test (NAAT).

See Test Principle below for more details.

J. Substantial Equivalence Information:

1. Predicate device name(s):

APTMA Combo 2 Assay

2. Predicate 510(k) number(s):

K003395

3. Comparison with predicate:

Same device as current submission, but with added PreservCyt (PC) specimen indication. Collection device and media are different, as are specimen handling and storage instructions.

K. Standard/Guidance Document Referenced (if applicable):

Not Applicable

L. Test Principle:

The GEN-PROBE APTIMA Combo 2 Assay combines the technologies of target capture, Transcription-Mediated Amplification (TMA), and Dual Kinetic Assay (DKA). During target capture, rRNA molecules are isolated from specimens by capture oligomers on magnetic microparticles. After target capture, the specimens are ready for TMA. The GEN-PROBE APTIMA Combo 2 Assay reaction replicates a specific region of the 23S rRNA from *C. trachomatis* and a specific region of the 16S rRNA from *N. gonorrhoeae* via DNA intermediates. Detection of the rRNA amplicons is achieved using single-stranded chemiluminescent DNA probes, which are labeled with different acridinium ester molecules. The labeled DNA probes combine with amplicon to form stable RNA:DNA hybrids and light emitted from the labeled RNA:DNA hybrids is reported as Relative Light Units (RLU). In DKA, differences in the kinetic profiles of the CT and GC probes allow for the differentiation of signal. The chemiluminescent detection reaction for CT signal has the “flasher” kinetic type. The chemiluminescent detection reaction for the GC signal has the “glower” kinetic type. Assay results are determined by a cut-off based on the total RLU and the kinetic curve type.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. *Precision/Reproducibility:*

Precision testing was performed at three sites to obtain measures of repeatability and reproducibility. Reproducibility was established with a 12-member panel generated by spiking PreservCyt Solution with 0 to 2000 fg/assay of *C. trachomatis* and 0 to 5,000 fg/assay of *N. gonorrhoeae* rRNA

and aliquotting 1mL into the APTIMA Specimen Transfer Kit collection tube. Two (2) operators at each of the three sites performed one run per day on each of three days, totaling three valid runs per operator. Testing was performed using one assay kit lot. The results of this precision study are summarized below.

Reproducibility when testing PreservCyt liquid Pap clinical specimens containing target organism has not been determined.

Concentration (fg/assay)															
				Intra-Run		Inter-Run		Inter-Site		Inter-Operator					
CT	GC	N	Agreement	Mean RLU (x 1,000)	SD (x 1,000)	CV (%)	SD (x 1,000)	CV (%)	SD (x 1,000)	CV (%)	SD (x 1,000)	CV (%)			
0	0	162	97.5%	9.7	31.6	N/A	3.4	N/A	6.4	N/A	4.7	N/A			
0	5,000	54	98.3%	1298	146	11.3	54.8	4.2	0.0	0.0	0.0	0.0			
2,000	0	54	100%	1140	54.1	4.7	79.8	7.0	101	8.9	2.4	0.2			
2,000	5,000	54	100%	2345	79.8	3.4	78.0	3.3	94.7	4.0	37.9	1.6			
0	250	54	100%	953	114	12.0	0.0	0.0	161	18.9	90.7	9.5			
5	0	54	100%	971	58.3	6.0	71.7	7.4	22.8	2.4	85.0	8.8			
1,000	2,500	54	100%	2294	114	5.0	88.9	3.9	153	6.7	0.0	0.0			
100	250	54	98.1%	1911	139	7.3	130	6.8	348	18.2	39.7	2.1			
5	5,000	54	100%	2138	113	5.3	130	6.1	98.8	4.6	166	7.8			
2,000	250	54	98.3%	2044	138	6.7	169	8.3	360	17.6	26.9	1.3			

RLU - Relative Light Units
SD = Standard Deviation
CV = Coefficient of Variation
N/A = Not applicable for negative panel members

b. *Linearity/assay reportable range:*

Not Applicable

c. *Traceability, Stability, Expected values (controls, calibrators, or methods):*

Not Applicable

d. *Detection limit:*

A study was performed that showed the AC2 Assay detected CT cells at least 4 fold below the analytical sensitivity claim (claim = 1 IFU/assay) for 3 replicates of each of 15 CT serovars tested in PC media. Likewise, 3 replicates of each of 20 GC clinical isolates were detected in PC media at 10 fold below the analytical sensitivity claim (claim = 50 cells/assay).

Table 5.5-10: Analytical Sensitivity for Detection of CT

Serovar	IFU/ Assay	Replicate 1		Replicate 2		Replicate 3	
		RLU	Result	RLU	Result	RLU	Result
A	10	1,006,000	CT+	1,019,000	CT+	964,000	CT+
	1.0	402,000	CT+	295,000	CT+	455,000	CT+
	0.1	53,000	CT*	82,000	CT*	41,000	CT*
	0.01	10,000	CT-	19,000	CT-	12,000	CT-
B	10	1,044,000	CT+	1,127,000	CT+	1,005,000	CT+
	1.0	607,000	CT+	725,000	CT+	724,000	CT+
	0.1	165,000	CT+	44,000	CT*	153,000	CT+
	0.01	68,000	CT*	27,000	CT-	18,000	CT-
Ba	10	1,139,000	CT+	1,201,000	CT+	1,275,000	CT+
	1.0	973,000	CT+	1,045,000	CT+	1,083,000	CT+
	0.1	699,000	CT+	826,000	CT+	783,000	CT+
	0.01	422,000	CT+	271,000	CT+	262,000	CT+
C	10	1,179,000	CT+	1,205,000	CT+	1,078,000	CT+
	1.0	740,000	CT+	716,000	CT+	576,000	CT+
	0.1	175,000	CT+	172,000	CT+	240,000	CT+
	0.01	7,000	CT-	17,000	CT-	40,000	CT*
D	10	989,000	CT+	1,005,000	CT+	993,000	CT+
	1.0	577,000	CT+	498,000	CT+	473,000	CT+
	0.1	48,000	CT*	56,000	CT*	91,000	CT*
	0.01	15,000	CT-	11,000	CT-	7,000	CT-
E	10	1,182,000	CT+	1,168,000	CT+	1,262,000	CT+
	1.0	1,134,000	CT+	1,150,000	CT+	1,204,000	CT+
	0.1	845,000	CT+	791,000	CT+	773,000	CT+
	0.01	442,000	CT+	403,000	CT+	268,000	CT+
F	10	1,167,000	CT+	1,179,000	CT+	1,169,000	CT+
	1.0	831,000	CT+	838,000	CT+	804,000	CT+
	0.1	241,000	CT+	6,000	CT-	19,000	CT-
	0.01	6,000	CT-	5,000	CT-	5,000	CT-
G	10	1,190,000	CT+	1,236,000	CT+	1,175,000	CT+
	1.0	1,018,000	CT+	1,030,000	CT+	979,000	CT+
	0.1	816,000	CT+	916,000	CT+	675,000	CT+
	0.01	71,000	CT*	412,000	CT+	454,000	CT+

* = Equivocal Result

Serovar	IFU/ Assay	Replicate 1		Replicate 2		Replicate 3	
		RLU	Result	RLU	Result	RLU	Result
H	10	1,136,000	CT+	1,120,000	CT+	1,111,000	CT+
	1.0	872,000	CT+	937,000	CT+	998,000	CT+
	0.1	795,000	CT+	733,000	CT+	657,000	CT+
	0.01	369,000	CT+	251,000	CT+	114,000	CT+
I	10	1,039,000	CT+	1,015,000	CT+	1,039,000	CT+
	1.0	1,025,000	CT+	992,000	CT+	896,000	CT+
	0.1	907,000	CT+	884,000	CT+	938,000	CT+
	0.01	899,000	CT+	924,000	CT+	928,000	CT+
J	10	1,159,000	CT+	1,221,000	CT+	1,297,000	CT+
	1.0	1,160,000	CT+	1,086,000	CT+	1,049,000	CT+
	0.1	802,000	CT+	803,000	CT+	788,000	CT+
	0.01	524,000	CT+	487,000	CT+	486,000	CT+
K	10	1,175,000	CT+	1,189,000	CT+	1,117,000	CT+
	1.0	886,000	CT+	844,000	CT+	909,000	CT+
	0.1	469,000	CT+	672,000	CT+	478,000	CT+
	0.01	5,000	CT-	5,000	CT-	7,000	CT-
L1	10	1,192,000	CT+	1,147,000	CT+	1,127,000	CT+
	1.0	947,000	CT+	935,000	CT+	942,000	CT+
	0.1	276,000	CT+	806,000	CT+	787,000	CT+
	0.01	504,000	CT+	334,000	CT+	294,000	CT+
L2	10	1,072,000	CT+	1,098,000	CT+	1,203,000	CT+
	1.0	979,000	CT+	839,000	CT+	901,000	CT+
	0.1	509,000	CT+	213,000	CT+	471,000	CT+
	0.01	11,000	CT-	5,000	CT-	230,000	CT+
L3	10	1,205,000	CT+	1,250,000	CT+	1,133,000	CT+
	1.0	999,000	CT+	967,000	CT+	951,000	CT+
	0.1	566,000	CT+	582,000	CT+	637,000	CT+
	0.01	196,000	CT+	9,000	CT-	111,000	CT+

* = Equivocal Result

Table 5.5-11: Additional Dilution Testing of CT Serovars A, B, D and F.

Serovar A					
IFU/Assay	Rep #	PreservCvt/STM		STM only	
		RLU	Results	RLU	Results
1.0	1	713,000	CT+	769,000	CT+
	2	632,000	CT+	675,000	CT+
	3	594,000	CT+	741,000	CT+
0.75	1	412,000	CT+	692,000	CT+
	2	287,000	CT+	674,000	CT+
	3	618,000	CT+	606,000	CT+
0.50	1	565,000	CT+	365,000	CT+
	2	344,000	CT+	109,000	CT+
	3	316,000	CT+	120,000	CT+
0.25	1	322,000	CT+	170,000	CT+
	2	272,000	CT+	279,000	CT+
	3	189,000	CT+	309,000	CT+
0.10	1	14,000	CT-	38,000	CT*
	2	43,000	CT*	44,000	CT*
	3	45,000	CT*	83,000	CT*

Serovar B					
IFU/Assay	Rep #	PreservCvt/STM		STM only	
		RLU	Results	RLU	Results
1.0	1	771,000	CT+	533,000	CT+
	2	754,000	CT+	572,000	CT+
	3	663,000	CT+	546,000	CT+
0.75	1	616,000	CT+	562,000	CT+
	2	576,000	CT+	539,000	CT+
	3	581,000	CT+	640,000	CT+
0.50	1	473,000	CT+	558,000	CT+
	2	454,000	CT+	620,000	CT+
	3	483,000	CT+	557,000	CT+
0.25	1	439,000	CT+	426,000	CT+
	2	514,000	CT+	360,000	CT+
	3	376,000	CT+	527,000	CT+
0.10	1	154,000	CT+	295,000	CT+
	2	216,000	CT+	175,000	CT+
	3	283,000	CT+	280,000	CT+

* = Equivocal Result

Table 5.5-11 (cont). Additional Dilution Testing of CT Serovars A, B, D and F.

Serovar D					
IFU/Assay	Rep #	PreservCvt/STM		STM only	
		RLU	Results	RLU	Results
1.0	1	873,000	CT+	897,000	CT+
	2	904,000	CT+	928,000	CT+
	3	701,000	CT+	764,000	CT+
0.75	1	609,000	CT+	731,000	CT+
	2	422,000	CT+	649,000	CT+
	3	859,000	CT+	542,000	CT+
0.50	1	803,000	CT+	388,000	CT+
	2	628,000	CT+	465,000	CT+
	3	714,000	CT+	531,000	CT+
0.25	1	429,000	CT+	147,000	CT+
	2	366,000	CT+	250,000	CT+
	3	348,000	CT+	122,000	CT+
0.10	1	28,000	CT-	153,000	CT+
	2	92,000	CT*	170,000	CT+
	3	122,000	CT+	69,000	CT*

Serovar F					
IFU/Assay	Rep #	PreservCvt/STM		STM only	
		RLU	Results	RLU	Results
1.0	1	741,000	CT+	552,000	CT+
	2	863,000	CT+	133,000	CT+
	3	855,000	CT+	168,000	CT+
0.75	1	734,000	CT+	341,000	CT+
	2	716,000	CT+	168,000	CT+
	3	485,000	CT+	26,000	CT-
0.50	1	676,000	CT-	1,049,000	CT+
	2	310,000	CT+	911,000	CT+
	3	408,000	CT+	919,000	CT+
0.25	1	631,000	CT+	439,000	CT+
	2	561,000	CT+	725,000	CT+
	3	656,000	CT+	526,000	CT+
0.10	1	6,000	CT-	5,000	CT-
	2	7,000	CT-	6,000	CT-
	3	6,000	CT-	5,000	CT-

* = Equivocal Result

Table 5.5-12: Analytical Sensitivity for Detection of GC

GP C/μ	GC Cells Per Assay	Replicate 1		Replicate 2		Replicate 3	
		RLU	Result	RLU	Result	RLU	Result
760	500	1,195,000	GC+	1,181,000	GC+	1,051,000	GC+
	50	707,000	GC+	839,000	GC+	659,000	GC+
	5	326,000	GC+	420,000	GC+	466,000	GC+
	0.5	50,000	GC-	65,000	GC±	70,000	GC±
764	500	1,203,000	GC+	1,169,000	GC+	1,198,000	GC+
	50	1,073,000	GC+	1,049,000	GC+	1,120,000	GC+
	5	665,000	GC+	679,000	GC+	633,000	GC+
	0.5	79,000	GC±	107,000	GC±	45,000	GC-
772	500	1,277,000	GC+	1,255,000	GC+	1,234,000	GC+
	50	912,000	GC+	631,000	GC+	1,085,000	GC+
	5	532,000	GC+	394,000	GC+	652,000	GC+
	0.5	100,000	GC±	135,000	GC±	66,000	GC±
783	500	1,075,000	GC+	1,064,000	GC+	1,101,000	GC+
	50	983,000	GC+	938,000	GC+	963,000	GC+
	5	693,000	GC+	719,000	GC+	786,000	GC+
	0.5	226,000	GC+	237,000	GC+	283,000	GC+
787	500	1,146,000	GC+	1,046,000	GC+	1,093,000	GC+
	50	959,000	GC+	871,000	GC+	875,000	GC+
	5	295,000	GC+	397,000	GC+	396,000	GC+
	0.5	95,000	GC±	47,000	GC-	58,000	GC-
789	500	1,186,000	GC+	1,181,000	GC+	1,118,000	GC+
	50	954,000	GC+	1,014,000	GC+	941,000	GC+
	5	209,000	GC+	551,000	GC+	628,000	GC+
	0.5	264,000	GC+	129,000	GC±	93,000	GC±
790	500	1,041,000	GC+	1,011,000	GC+	1,034,000	GC+
	50	791,000	GC+	793,000	GC+	814,000	GC+
	5	365,000	GC+	313,000	GC+	275,000	GC+
	0.5	25,000	GC-	18,000	GC-	16,000	GC-
793	500	1,259,000	GC+	1,246,000	GC+	1,239,000	GC+
	50	1,182,000	GC+	1,150,000	GC+	1,139,000	GC+
	5	830,000	GC+	792,000	GC+	915,000	GC+
	0.5	343,000	GC+	386,000	GC+	349,000	GC+
794	500	1,192,000	GC+	1,117,000	GC+	1,237,000	GC+
	50	1,166,000	GC+	1,132,000	GC+	1,157,000	GC+
	5	797,000	GC+	159,000	GC+	706,000	GC+
	0.5	220,000	GC+	205,000	GC+	200,000	GC+
795	500	1,148,000	GC+	1,116,000	GC+	1,113,000	GC+
	50	923,000	GC+	1,038,000	GC+	1,053,000	GC+
	5	702,000	GC+	547,000	GC+	714,000	GC+
	0.5	177,000	GC+	93,000	GC+	158,000	GC+

GP Cl#	GC Cells Per Assay	Replicate 1		Replicate 2		Replicate 3	
		RLU	Result	RLU	Result	RLU	Result
799	500	1,170,000	GC+	1,134,000	GC+	1,109,000	GC-
	50	976,000	GC+	950,000	GC+	992,000	GC-
	5	544,000	GC+	580,000	GC+	524,000	GC-
	0.5	62,000	GCe	54,000	GC-	103,000	GCe
800	500	1,231,000	GC+	1,221,000	GC+	1,213,000	GC-
	50	1,164,000	GC+	1,111,000	GC+	1,118,000	GC+
	5	903,000	GC+	931,000	GC+	900,000	GC-
	0.5	147,000	GCe	214,000	GC+	191,000	GC-
801	500	1,059,000	GC+	1,091,000	GC+	1,175,000	GC-
	50	1,002,000	GC+	1,025,000	GC+	974,000	GC+
	5	442,000	GC+	469,000	GC+	491,000	GC-
	0.5	62,000	GCe	69,000	GCe	40,000	GC-
802	500	1,156,000	GC+	1,187,000	GC+	1,168,000	GC-
	50	1,026,000	GC+	1,122,000	GC+	1,100,000	GC+
	5	709,000	GC+	622,000	GC+	323,000	GC-
	0.5	72,000	GCe	134,000	GCe	153,000	GC+
806	500	1,145,000	GC+	1,088,000	GC+	1,010,000	GC-
	50	781,000	GC+	851,000	GC+	822,000	GC-
	5	255,000	GC+	247,000	GC+	178,000	GC-
	0.5	61,000	GCe	31,000	GC-	39,000	GC-
813	500	1,112,000	GC+	1,134,000	GC-	1,163,000	GC-
	50	1,006,000	GC+	1,030,000	GC+	1,108,000	GC+
	5	669,000	GC+	688,000	GC+	608,000	GC+
	0.5	141,000	GCe	135,000	GCe	176,000	GC+
827	500	1,192,000	GC+	1,089,000	GC+	1,250,000	GC-
	50	1,168,000	GC+	1,192,000	GC+	1,133,000	GC+
	5	744,000	GC+	841,000	GC+	712,000	GC+
	0.5	178,000	GC+	245,000	GC+	202,000	GC-
3043	500	1,235,000	GC+	1,238,000	GC+	1,232,000	GC-
	50	1,090,000	GC+	1,036,000	GC+	1,125,000	GC+
	5	292,000	GC+	391,000	GC+	493,000	GC-
	0.5	78,000	GCe	99,000	GCe	73,000	GCe
3045	500	1,184,000	GC+	1,110,000	GC+	1,211,000	GC-
	50	1,116,000	GC+	1,077,000	GC+	978,000	GC+
	5	641,000	GC+	808,000	GC+	513,000	GC-
	0.5	110,000	GCe	111,000	GCe	151,000	GC-
3047	500	1,189,000	GC+	1,174,000	GC+	1,163,000	GC-
	50	1,071,000	GC+	1,148,000	GC+	1,182,000	GC+
	5	880,000	GC+	896,000	GC+	794,000	GC-
	0.5	70,000	GCe	207,000	GC+	177,000	GC-

e = Equivocal Result

e. Analytical specificity:

The Chlamydia and Neisseria species were used to evaluate the analytical specificity of the AC2 Assay. A total of 50 culture isolates were tested in the liquid Pap media. None of the 50 culture isolates produced a positive result in the AC2 Assay. See results below:

Table 5.5-13: Specificity of the AC2 Assay

PHYLOGENETIC CROSS-SECTION	GP No.	ATCC No.	Concentration Tested/Assay	Rep #	Results (RLU)
<i>Chlamydia psittaci</i>	1537	VR601	7.9 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Chlamydia psittaci</i>	768	VR629	1 x 10 ⁸ CELD ₉₉ /0.2 ml	1	4,000
				2	5,000
<i>Chlamydia pneumoniae</i>	1404	VR1360	4.0 x 10 ⁸ cells	1	5,000
				2	6,000
<i>Neisseria elongata</i>	CI1502	49377	1.2 x 10 ⁹ cells	1	6,000
				2	6,000
<i>Neisseria elongata</i>	CI1503	49378	1.2 x 10 ⁸ cells	1	5,000
				2	5,000
<i>Neisseria elongata</i>	CI1504	49379	1.8 x 10 ⁸ cells	1	9,000
				2	10,000
<i>Neisseria flava</i>	1538	14221	2.5 x 10 ⁸ cells	1	4,000
				2	4,000
<i>Neisseria mucosa</i>	190	19696	2.4 x 10 ⁸ cells	1	5,000
				2	5,000
<i>Neisseria mucosa</i>	791	25999	1.3 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Neisseria perflava</i>	1559	10555	1.0 x 10 ⁸ cells	1	6,000
				2	4,000
<i>Neisseria sicca</i>	272	9913	5.9 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Neisseria sicca</i>	762	29193	7.2 x 10 ⁷ cells	1	4,000
				2	4,000
<i>Neisseria subflava</i>	CI3113	NA	4.5 x 10 ⁷ cells	1	4,000
				2	4,000
<i>Neisseria subflava</i>	NH1	NA	2.2 x 10 ⁷ cells	1	4,000
				2	3,000
<i>Neisseria subflava</i>	NH5	NA	2.4 x 10 ⁷ cells	1	3,000
				2	4,000
<i>Neisseria subflava</i>	NH6	NA	2.7 x 10 ⁸ cells	1	4,000
				2	4,000
<i>Neisseria subflava</i>	NH7	NA	1.2 x 10 ⁸ cells	1	4,000
				2	5,000
<i>Neisseria subflava</i>	NH8	NA	4.7 x 10 ⁸ cells	1	6,000
				2	4,000
<i>Neisseria subflava</i>	NH11	NA	1.0 x 10 ⁸ cells	1	5,000
				2	5,000
<i>Neisseria subflava</i>	NH13	NA	3.0 x 10 ⁷ cells	1	6,000
				2	6,000

PHYLOGENETIC CROSS-SECTION	GP No.	ATCC No.	Concentration Tested/Assay	Rep #	Results (RLU)
<i>Neisseria subflava</i>	NH14	NA	4.5 x 10 ⁷ cells	1	6,000
				2	5,000
<i>Neisseria subflava</i>	NH15	NA	8.7 x 10 ⁷ cells	1	4,000
				2	4,000
<i>Neisseria subflava</i>	NH17	NA	1.1 x 10 ⁸ cells	1	3,000
				2	4,000
<i>Neisseria subflava</i>	NH18	NA	5.7 x 10 ⁷ cells	1	5,000
				2	5,000
<i>Neisseria subflava</i>	NH20	NA	2.1 x 10 ⁸ cells	1	6,000
				2	5,000
<i>Neisseria cinerea</i>	761	14685	1.6 x 10 ⁸ cells	1	8,000
				2	8,000
<i>Neisseria cinerea</i>	CI3051	NA	2.9 x 10 ⁸ cells	1	8,000
				2	8,000
<i>Neisseria cinerea</i>	CI4543	NA	1.3 x 10 ⁸ cells	1	9,000
				2	10,000
<i>Neisseria cinerea</i>	CI4546	NA	4.9 x 10 ⁸ cells	1	7,000
				2	7,000
<i>Neisseria dentrificans</i>	763	14686	3.0 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Neisseria lactamica</i>	760	23970	1.1 x 10 ⁸ cells	1	6,000
				2	6,000
<i>Neisseria lactamica</i>	CI3013	NA	9.5 x 10 ⁸ cells	1	5,000
				2	5,000
<i>Neisseria lactamica</i>	CI3018	NA	1.5 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Neisseria lactamica</i>	CI3021	NA	5.0 x 10 ⁸ cells	1	4,000
				2	14,000
<i>Neisseria lactamica</i>	CI3022	NA	9.4 x 10 ⁸ cells	1	5,000
				2	4,000
<i>Neisseria lactamica</i>	CI3049	NA	2.3 x 10 ¹⁰ cells	1	5,000
				2	4,000
<i>Neisseria lactamica</i>	CI3065	NA	1.5 x 10 ¹⁰ cells	1	4,000
				2	5,000
<i>Neisseria lactamica</i>	CI3067	NA	2.9 x 10 ¹⁰ cells	1	4,000
				2	4,000
<i>Neisseria lactamica</i>	CI834	NA	6.0 x 10 ⁸ cells	1	6,000
				2	5,000
<i>N. meningitidis</i> Serogroup A	755	13077	4.0 x 10 ¹⁰ cells	1	9,000
				2	8,000
<i>N. meningitidis</i> Serogroup B	756	13090	3.1 x 10 ⁸ cells	1	6,000
				2	6,000

PHYLOGENETIC CROSS-SECTION	GP No.	ATCC No.	Concentration Tested/Assay	Rep #	Results (RLU)
<i>N. meningitidis</i> Serogroup C	757	13102	5.0 x 10 ⁸ cells	1	8,000
				2	8,000
<i>N. meningitidis</i> Serogroup C	1388	13109	2.9 x 10 ¹¹ cells	1	5,000
				2	4,000
<i>N. meningitidis</i> Serogroup C	1389	13100	8.0 x 10 ¹⁰ cells	1	5,000
				2	4,000
<i>N. meningitidis</i> Serogroup C	1390	13112	2.7 x 10 ⁸ cells	1	10,000
				2	7,000
<i>N. meningitidis</i> Serogroup D	401	13113	1.7 x 10 ¹¹ cells	1	8,000
				2	7,000
<i>N. meningitidis</i> Serogroup Y	787	35561	3.0 x 10 ¹⁰ cells	1	8,000
				2	9,000
<i>N. meningitidis</i> Serogroup W135	1387	43744	1.6 x 10 ⁸ cells	1	10,000
				2	10,000
<i>Neisseria polysacchara</i>	1489	43768	3.0 x 10 ⁸ cells	1	12,000
				2	12,000

f. Assay cut-off:

Not Applicable

2. Comparison studies:

a. Method comparison with predicate device:

Not Applicable

3. Clinical studies:

A prospective multi-center clinical study was conducted to evaluate the use of the PreservCyt Solution (a component of the ThinPrep 2000 System) as an alternative medium for gynecological specimens for the detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae*. One thousand six hundred forty-seven (1,647) symptomatic and asymptomatic female subjects attending OB/GYN, family planning, public health, women's and STD clinics were evaluated in the clinical study. Of the 1,647 evaluable subjects, 1,288 were asymptomatic subjects and 359 were symptomatic subjects. Subjects were enrolled from sites with CT prevalence that ranged from 3.2 to 14.0% and GC prevalence that ranged from 0 to 5.0%. Two specimens were collected from each eligible subject: one PreservCyt Solution specimen and one endocervical swab. PreservCyt Solution specimens were processed in accordance with the ThinPrep 2000 Processor Operator's Manual and APTIMA Specimen Transfer Kit Package Insert. After processing the PreservCyt Solution specimen with the ThinPrep 2000 Processor, the specimen was transferred into the APTIMA Specimen Transfer Kit for testing with the APTIMA Combo 2 Assay. The PreservCyt liquid Pap specimens and endocervical swab specimens were tested with the APTIMA Combo 2 Assay.

Sensitivity and specificity for PreservCyt liquid Pap specimens were calculated by comparing results to a patient infected status algorithm. In the algorithm, the designation of a subject as being infected or non-infected with *C. trachomatis* or *N. gonorrhoeae* was based on endocervical swab specimen results from two commercially-available NAATs. For *C. trachomatis*, the reference NAATs included the APTIMA Combo 2 Assay and the APTIMA CT Assay. For *N. gonorrhoeae*, the reference NAATs included the APTIMA Combo 2 Assay and the APTIMA GC Assay. Positive results from both reference NAATs were required to establish an infected patient. A non-infected patient was established if the results from the two reference NAATs disagreed or were negative.

Sensitivity and specificity for *C. trachomatis* in PreservCyt liquid Pap specimens tested in the APTIMA Combo 2 Assay, by symptom status and overall, is presented in Table 5c. For *C. trachomatis*, overall sensitivity was 96.7% (87/90). In symptomatic and asymptomatic subjects, sensitivity was 96.7% (29/30) and 96.7% (58/60), respectively. Overall specificity for *C. trachomatis* PreservCyt liquid Pap specimens was 99.2% (1545/1557). In symptomatic and asymptomatic subjects, specificity was 98.5% (324/329) and 99.4% (1221/1228), respectively. Table 6c shows the APTIMA Combo 2 Assay sensitivity and specificity values for *C. trachomatis* in PreservCyt

liquid Pap specimens by clinical site and overall. For *C. trachomatis*, the sensitivity ranged from 92.9% to 100%. The specificity ranged from 97.7% to 100%. Sensitivity and specificity for *N. gonorrhoeae* in PreservCyt liquid Pap specimens tested in the APTIMA Combo 2 Assay, by symptom status and overall, is presented in Table 9c. For *N. gonorrhoeae*, overall sensitivity was 92.3% (12/13). In symptomatic and asymptomatic subjects, sensitivity was 100% (7/7) and 83.3% (5/6), respectively. Overall specificity for *N. gonorrhoeae* PreservCyt liquid Pap specimens was 99.8% (1630/1633). In symptomatic and asymptomatic subjects, specificity was 100% (352/352) and 99.8% (1278/1281), respectively. Table 10c shows the APTIMA Combo 2 Assay sensitivity and specificity values for *N. gonorrhoeae* in PreservCyt liquid Pap specimens by clinical site and overall. For *N. gonorrhoeae*, the sensitivity ranged from 80.0% to 100%. Specificity ranged from 99.0% to 100%

Table 7c. *C. trachomatis* PreservCyt Liquid Pap Specimen Analysis for Female Patient Infected Status

Patient Infected Status	Endocervical Swab		Symptom Status	
	APTIMA Combo 2 Assay	APTIMA CT Assay	Symptomatic	Asymptomatic
Infected	+	+	30	60
Non-Infected	-	+	4	12
Non-Infected	+	-	3	2
Non-Infected	-	-	322	1214
Total			359	1288

Table 5c. *C. trachomatis* Sensitivity and Specificity: APTIMA Combo 2 Assay PreservCyt Specimens vs. Patient Infected Status

Symptom Status	AC2/CT ThinPrep Result	AC2 Assay / ACT Assay				Sensitivity (95% C.I.)	Specificity (95% C.I.)
		+/+	+/-	-/+	-/-		
Asymptomatic	Positive	58	1	0	6	98.7% (88.5 - 99.8)	99.4% (98.8 - 99.8)
	Negative	2	1	12	1208		
	Total	60	2	12	1214		
Symptomatic	Positive	29	0	0	5	98.7% (82.8 - 99.9)	98.5% (96.5 - 99.5)
	Negative	1	3	4	317		
	Total	30	3	4	322		
All	Positive	87	1	0	11	98.7% (90.6 - 99.3)	99.2% (98.7 - 99.8)
	Negative	3	4	16	1525		
	Total	90	5	16	1536		

+/+ = Positive Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the ACT Assay
 +/- = Positive Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the ACT Assay
 -/+ = Negative Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the ACT Assay
 -/- = Negative Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the ACT Assay

Table 6c. *C. trachomatis* Performance by Clinical Site: APTIMA Combo 2 Assay PreservCyt Specimens vs. Patient Infected Status

Site	AC2/CT ThinPrep Result	+/+	+/-	-/+	-/-	Prev (%)	Sensitivity (95% C.I.)	Specificity (95% C.I.)	PPV (%)	NPV (%)
	Negative	0	0	1	83					
	Total	14	0	1	85					
2	Positive	4	0	0	0	3.2	100% (39.8 - 100)	100% (97.0 - 100)	100	100
	Negative	0	0	2	118					
	Total	4	0	2	118					
3	Positive	29	0	0	2	6.5	93.5% (78.6 - 99.2)	99.5% (98.4 - 99.9)	93.5	99.5
	Negative	2	0	2	440					
	Total	31	0	2	442					
4	Positive	8	1	0	4	2.8	100% (63.1 - 100)	98.2% (95.9 - 99.4)	61.5	100
	Negative	0	2	1	271					
	Total	8	3	1	275					
5	Positive	13	0	0	2	4.7	92.9% (66.1 - 99.8)	99.3% (97.5 - 99.9)	86.7	99.6
	Negative	1	1	4	276					
	Total	14	1	4	278					
6	Positive	19	0	0	1	5.2	100% (82.4 - 100)	99.7% (98.4 - 100)	95.0	100
	Negative	0	1	6	337					
	Total	19	1	6	338					
All	Positive	87	1	0	11	5.5	96.7% (90.6 - 99.3)	99.2% (98.7 - 99.6)	87.9	99.8
	Negative	3	4	16	1525					
	Total	90	5	16	1536					

+/+ = Positive Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the ACT Assay

+/- = Positive Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the ACT Assay

-/+ = Negative Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the ACT Assay

-/- = Negative Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the ACT Assay

Table 11c. *N. gonorrhoea* PreservCyt Liquid Pap Specimen Analysis for Female Patient Infected Status

Patient Infected Status	Endocervical Swab		Symptom Status	
	APTIMA Combo 2 Assay	APTIMA GC Assay	Symptomatic	Asymptomatic
Infected	+	+	7	6
Non-Infected	=	+	0	1
Non-Infected	-	+	0	5
Non-Infected	-	-	352	1276
Total			359	1288

The equal symbol ("=") represents an equivocal result.

Table 9c. *N. gonorrhoeae* Sensitivity and Specificity: APTIMA Combo 2 Assay PreservCyt Specimens vs. Patient Infected Status

Symptom Status	AC2/GC ThinPrep Result	+/+	+/-	-/+	-/-	Sensitivity (95% C.I.)	Specificity (95% C.I.)
Asymptomatic	Positive	5	0	0	3	83.3% (35.9 - 99.6)	99.8% (99.3 - 100)
	Negative	1	0	5	1273		
	Total	6	0	5	1276		
Symptomatic	Positive	7	0	0	0	100% (59.0 - 100)	100% (99.0 - 100)
	Negative	0	0	0	352		
	Total	7	0	0	352		
All	Positive	12	0	0	3	92.3% (64.0 - 99.8)	99.8% (99.5 - 100)
	Negative	1	0	5	1625		
	Total	13	0	5	1628		

+/+ = Positive Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the AGC Assay
 +/- = Positive Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the AGC Assay
 -/+ = Negative Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the AGC Assay
 -/- = Negative Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the AGC Assay

Table 10c. *N. gonorrhoeae* Performance by Clinical Site: APTIMA Combo 2 Assay PreservCyt Specimens vs. Patient Infected Status

Site	AC2/GC ThinPrep Result	+/+	+/-	-/+	-/-	Prev (%)	Sensitivity (95% C.I.)	Specificity (95% C.I.)	PPV (%)	NPV (%)
1	Positive	5	0	0	0	5.0	100% (47.8 - 100)	100% (96.2 - 100)	100	100
	Negative	0	0	0	95					
	Total	5	0	0	95					
2	Positive	1	0	0	0	0.8	100% (2.5 - 100)	100% (97.0 - 100)	100	100
	Negative	0	0	0	123					
	Total	1	0	0	123					
3	Positive	4	0	0	0	1.1	80.0% (28.4 - 99.5)	100% (99.2 - 100)	100	99.8
	Negative	1	0	0	470					
	Total	5	0	0	470					
4	Positive	1	0	0	0	0.3	100% (2.5 - 100)	100% (98.7 - 100)	100	100
	Negative	0	0	3	283					
	Total	1	0	3	283					
5	Positive	0	0	0	3	0.0	N/A	99.0% (97.1 - 99.8)	0.0	100
	Negative	0	0	0	294					
	Total	0	0	0	297					
6	Positive	1	0	0	0	0.3	100% (2.5 - 100)	100% (99.0 - 100)	100	100
	Negative	0	0	2	360					
	Total	1	0	2	360					
All	Positive	12	0	0	3	0.8	92.3% (64.0 - 99.8)	99.8% (99.5 - 100)	80.0	99.9
	Negative	1	0	5	1625					
	Total	13	0	5	1628					

+/+ = Positive Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the AGC Assay
 +/- = Positive Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the AGC Assay
 -/+ = Negative Endocervical Swab Specimen Result in the AC2 Assay / Positive Endocervical Swab Specimen Result in the AGC Assay
 -/- = Negative Endocervical Swab Specimen Result in the AC2 Assay / Negative Endocervical Swab Specimen Result in the AGC Assay

4. Clinical cut-off:

Not Applicable

5. Expected values/Reference range:

The prevalence of *C. trachomatis* and/or *N. gonorrhoeae* disease in patient populations depends on risk factors such as age, gender, the presence of symptoms,

the type of clinic, and the test method. A summary of the prevalence of *C. trachomatis* and *N. gonorrhoeae* as determined by the APTIMA Combo 2 Assay results on PreservCyt liquid Pap specimens is shown below by clinical site and overall.

Table 1c. Prevalence of *C. trachomatis* and/or *N. gonorrhoeae* Disease by Clinical Site: PreservCyt liquid Pap Specimen

Site	PreservCyt liquid Pap % Prevalence (# positive/# tested)		
	CT+/GC+	CT+/GC-	CT-/GC+
1	3.0 (3/100)	13.0 (13/100)	2.0 (2/100)
2	0 (0/124)	3.2 (4/124)	0.8 (1/124)
3	0.4 (2/475)	6.1 (29/475)	0.4 (2/475)
4	0.4 (1/287)	4.2 (12/287)	0 (0/287)
5	0 (0/297)	5.1 (15/297)	1.0 (3/297)
6	0 (0/384)	5.5 (20/384)	0.8 (2/384)
ALL	0.4 (6/1647)	5.6 (93/1647)	0.6 (10/1647)

Note: The CT and GC prevalence were calculated using the APTIMA Combo 2 Assay results of PreservCyt liquid Pap specimen.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

1. The submitted information in this premarket notification is complete and supports a substantial equivalence decision.