

# CHROMIUM, Hexavalent, in Settled Dust Samples

9101

Cr(VI)

MW: 52.00

CAS: 18540-29-9

RTECS: GB6262000

---

METHOD: 9101, Issue 1

EVALUATION: Not Applicable

Issue 1: 15 May 1996

---

**PURPOSE:** Estimation of soluble hexavalent chromium content of settled dust.

**LIMIT OF**

**DETECTION:** 1 µg Cr(VI) per sample

**FIELD**

**EQUIPMENT:**

1. Chromate (Diphenylcarbazide reagent) test kit (Chemetrics Chromate Kit, or equivalent)
2. Sulfuric acid, 20% w/v (included in test kit)
3. Extraction solution, 2% NaOH/3% Na<sub>2</sub>CO<sub>3</sub> in deionized water
4. Deionized water
5. Centrifuge tubes, 15-mL, graduated, clear plastic with screw-caps, disposable
6. Spatula, ~0.1 cm<sup>3</sup> capacity
7. pH paper

**PROCEDURE:**

1. Place 1 spatula full of dust (approximately 0.1 cm<sup>3</sup>; the size of a small pea) to be tested in a 15-mL clear plastic centrifuge tube. Add extraction solution up to the 2-mL mark. Cap the tube and shake vigorously.
2. Allow the tube to stand for 10 minutes, or longer, with occasional shaking.  
NOTE: Gently heating the tube in hot water will increase the sensitivity of the test.
3. Uncap the centrifuge tube and add deionized water to the 7-mL mark. Mix and allow the residue to settle.
4. Decant or pipet approximately 3 mL of the supernatant liquid into a second tube.  
NOTE: The sample may be filtered, if excessively turbid.
5. Add 9 drops of 20% sulfuric acid (3 drops /mL of decanted liquid), cap the tube, and invert to mix the contents.
6. Check the pH of the liquid with pH paper. If necessary, add 20% sulfuric acid dropwise to bring to pH <1.
7. Follow the instructions for color development.

NOTE: For more accurate determination of total hexavalent chromium in the dust, send a sample to the laboratory for analysis by Method 7600 or 7604.

**METHOD WRITTEN BY:** Mark Millson, MRSB/DPSE, and Peter Eller, QASA/DPSE