

# Detection of Ki-67 (TEC-3) in Formalin Fixed Paraffin-Embedded Mouse Tissue

## Reagents:

[1X Automation Buffer](#)

[3% Hydrogen Peroxide](#)

[Antibody Diluent](#)

[Citrate Buffer](#)

[DAB Chromagen](#)

[Hematoxylin](#)

## Antibody Information

Block: Normal Rabbit Serum

Jackson Immunoresearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Catalog #011-000-001

Avidin/Biotin Blocking Kit

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog #SP-2001

Negative control: Normal Rat Serum

Jackson Immunoresearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Primary antibody: Rat anti- mouse Ki-67 (TEC 3)

Dako Corporation

Carpinteria, CA 93013

[www.dakousa.com](http://www.dakousa.com)

1-800-235-5763

Catalog #M7249

Secondary antibody: Biotinylated rabbit anti-rat IgG  
Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666

Label antibody :Vector Elite Label  
Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666  
Catalog #PK-6100

### **Staining Procedure**

- Positive Control Tissue: Mouse GI tissue
- Stain Localization: Restricted to the nucleus of proliferating cells with a distinct, prominent staining of nucleoli.

Deparaffinize and hydrate slides through the following solutions.

Xylene	2 times	5 minutes
100% EtOH	2 times	3 minutes
95% EtOH	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1 Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.

2 Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

3. Unmasking Techniques using the decloaker.

Add 500 ml D/W to the pan of the decloaker.

Place full rack of slides in 200 ml of 1X citrate buffer and place in the decloaker.

Decloak for 5 minutes. Pressure \_\_\_\_\_

Depressurize for 10 minutes.

Remove pan top and cool for 10 min. Temp \_\_\_\_\_

Rinse in D/W, 2x for 3 min each

4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

5. Apply 5% Normal Rabbit Serum solution for 20 minutes at room temperature.  
lot# \_\_\_\_\_ Reconstituted Date \_\_\_\_\_

6. Apply Avidin/Biotin block

Lot# \_\_\_\_\_ Exp Date \_\_\_\_\_ New Kit: yes / no

Apply avidin block - 15 min at RT.

Quick rinse in 1X AB.

Apply biotin block - 15 min at RT.

No wash, wipe excess block

**DO NOT RINSE SECTIONS WITH BUFFER.**

7. Apply primary antibody (TEC-3 Ki-67) at 1:80 and incubate 1 hour.

Lot# \_\_\_\_\_ Exp Date \_\_\_\_\_

For negative control slides, normalize the protein concentration of normal rat serum to the protein concentration of the primary antibody and use this to make the 1:80 dilution. Apply normal rat serum to the slides and incubate for one hour.

Lot# \_\_\_\_\_ Reconstituted Date \_\_\_\_\_

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply the Rabbit anti-rat secondary antibody at a 1:300 dilution and incubate for 30 minutes.

Lot# \_\_\_\_\_ Reconstituted Date \_\_\_\_\_

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply the label antibody (Vector elite) and incubate for 30 minutes.

Exp. Date \_\_\_\_\_ New Kit: yes / no

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.

(Add 1 drop of DAB per ml of substrate)

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes no

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30 seconds.

16. Rinse in tap water until water is clear.

17. Place slides in 1X Automation buffer for 1 minute with gentle agitation to blue slides.

18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% Ethanol	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip.

last updated 10/12/04