

Detection of Pan-MCM in Formalin-Fixed, Paraffin Embedded Rat Tissue Using the Ventana

Staining Procedure

-Positive Control Tissue: GI tract

-Stain Localization: Nuclear

Cut 2 sequential 6 micron sections onto tissue adhesive slides (Daigger ? Vernon Hills, Illinois - Superfrost Plus Catalog #GX15978Z). Do not pre-heat treat the slides in histology. The slides should be allowed to air dry.

Deparaffinize and hydrate slides through the following solutions

Xylene	2 times	5 minutes
100% EtOH	2 times	3 minutes
95% EtOH	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
2. Perform Heat Induced Epitope Retrieval using Microwave Oven.
Place a full rack of slides in Tissue Tek[®] container containing 200 mls distilled water.

MWO for 5 minutes at 50% power
Cool for 1 minute (Add 50 mls distilled water to container)
MWO for 5 minutes at 50% power
Cool 20 minutes at room temperature
Rinse in distilled water 3 X 2 minutes each

3. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
4. Place slides on Ventana

Protocol # 76

No Amplification

A/B Block Master Lot # _____

DAB Detection Kit Master Lot # _____

Titration run ? 32 minute incubation

Primary antibody : Polyclonal Rabbit anti-Human Pan-Minichromosome

Maintenance Protein (MCM)

BDPharmingen

San Deigo, CA 92121

www.bdbiosciences.com

1-877-232-8995

Catalog #559541

Lot # _____ Date received _____

1:300

Upon completion of run, remove slides from machine and rinse in warm detergent water to remove Ventana coverslip from slides.

5. Rinse in tap water 3 minutes.
6. Counterstain with Modified Harris Hematoxylin for 45 seconds.
7. Rinse in tap water until water is clear.
8. Place slides in 1X Automation buffer for 1 minute with gentle agitation to blue slides.
9. Dehydrate through the following solutions.

1. 95% Ethanol	2. 1 change	3. 3 minutes
4. 100% Ethanol	5. 3 changes	6. 3 minutes
7. Xylene	8. 2 changes	9. 5 minutes

10. Coverslip.

last updated 4/15/2003