




CONTACT

Heinrich Nakotte
 505-664-0284
 505-665-2676 (fax)
 hnakotte@nmsu.edu

instruments


SCD is used to determine the crystal and magnetic structures of a wide variety of materials. Neutrons are scattered from the crystalline sample onto a ^3He -gas-filled area detector that is position sensitive and measures $25 \times 25 \text{ cm}^2$. The wavelengths of the neutrons are determined by their time of flight from the source to the detector. To cover all of the reciprocal space for a particular crystal, the sample is mounted on a goniometer. For most crystals, about 10-15 histograms will provide a full data set.

Crystal Structure of $\text{La}_{2-2x}\text{Sr}_{1+2x}\text{Mn}_2\text{O}_7$



SCD can be used to study the structure of organometallic molecules; crystal-structure changes at solid/solid phase transitions; magnetic spin structures; twinned or multiple crystals; texture in polycrystalline materials; and structural or magnetic modifications of materials under pressures up to 20,000 atm. The instrument measures a large fraction of reciprocal space at one time and therefore can be used for studies of unknown incommensurate structures and diffuse scattering.

Orbital Ordering in $\text{LaSr}_2\text{Mn}_2\text{O}_7$



Super-lattice reflections in $\text{LaSr}_2\text{Mn}_2\text{O}_7$ with propagation vector $q = (h+1/4, k-1/4, l)$ are observed between 210-100 K.

SCD Specifications	
Wavelength range	0.5–10 Å
Beam diameter at sample	1–5 mm
Time resolution	1%
Maximum lattice constant	20 Å
Detector	1 multiwire proportional counter (25 cm x 25 cm) at 90°
Detector resolution	2.5 mm
Moderator	Chilled water at 283 K
Sample environment	10–300 K
Sample size	0.5–10 cm ³
Experiment duration	2 to 4 days per octant of reciprocal space