

X8A SAFETY CHECKLIST			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
Rev: C	Effective: 5/24/2010	Page 1 of 1	Number: LS-SCL-0017
Prepared by: Philip Marino		Reviewed by: P. Marino	Approved by: B. Chmiel

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

- 1. Valid Padlock and Index, no padlocks open
- 2. Exclusion Zone (EZ #2) wire mesh in place and secure between Valve V2A and BS #2.
- 3. Bremsstrahlung shield (BS #2) in place and banded, and bellows wrapped in lead as per photo.
- 4. Bremsstrahlung shield (BS #3) in place and banded, and bellows wrapped in lead as per photo.
- 5. Leaded glass on 8" View port on Bench #3. View port wrapped in lead.
- 6. Bremsstrahlung shield (BS #4) in place and banded, and bellows wrapped in lead as per photo.
- 7. Leaded glass on 2-3/4" View port downstream of Bench #3. View port wrapped in lead.
- 8. Leaded glass on 8" View port on Bench # 4. View port wrapped in lead.
- 9. Leaded glass on 4-1/2" View port located upstream of monochromator. View port wrapped in lead.
- 10. Leaded glass on 2-3/4" View port on upstream side of monochromator. View port wrapped in lead.
- 11. Leaded glass on 4-1/2" View port downstream of Photon shutter. View port wrapped in lead.

Note: Bremsstrahlung shield (BS #1) and Exclusion Zone (EZ #1) located downstream of shield wall are not accessible for this checklist.

OPCO / Beamline Staff: _____ Date: _____