Brookhaven National Laboratory – Photon Sciences - National Synchrotron Light Source					
Beamline Hazard Analysis - Beamline					
LS-BHA-X17A	Revision: 01	Effective: 10/17/11	Page 1 of 1		
The following tacks are included in the beamline safety review. Any operations not listed here, or performed by others than					

The following tasks are included in the beamline safety review. Any operations not listed here, or performed by others than those listed, need additional review by ES&H personnel and/or beamline personnel.

Prepared By:	L. Ehm	Approved By:	Lori Stiegler
· · epai ca byi	2. 2000	, pp. 01 cu b j .	Lon Stiegler

* Approval signatures on file with master copy

Tasks / Equipment	Performed by:		Controls			
	Beamline Users Staff					
Inert compressed gas			 Compressed gas training – TQ-COMPGAS1 or NSLS User training Transport cylinders using a cylinder cart Properly secure cylinders Wear safety glasses with side shields while connecting lines Leak check gas lines 			
Pouring cryogens ≤ 5L			 Use goggles Long pants and sleeves - avoid clothing which can trap cryogens (pockets, cuffs) Cryogen or heavy gloves Enclosed shoes Never pour from above chest level 			
Using the LN2 Fill Station or pressurized transfer			 LS-PROC-LN2FILL-MANL training for fill station Use full face shield along with goggles or safety glasses Long pants and sleeves - avoid clothing which can trap cryogens (pockets, cuffs) Cryogen or heavy gloves Enclosed shoes 			
Cryogenic autofill system			Ensure hutch ventilation system is operating			
High temperature			Warning signsTraining			
Cutting/razor blades	Y	Y	 Be aware of hand positioning Use safety knife when possible Consider using cut resistant gloves for frequent cutting tasks 			