

Brookhaven National Laboratory/ LIGHT SOURCES DIRECTORATE			
<b>Subject:</b>	<b>VACUUM PROCEDURES FOR BEAMLINE X-14A</b>		
<b>Number:</b>	LS-OPS-0139	<b>Revision:</b>	B
		<b>Effective:</b>	03/15/2011
			<b>Page 1 of 4</b>
Prepared By: J. Bai	Reviewed By: J. Klug	Approved By: S. Ehrlich	Approved By: E. Hu

\*Approval signatures on file with master copy.

**The only official copy of this file is the one on-line in the NSLS website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS website.**

The following procedures must be followed when bleeding up different beamline sections and when returning these sections to operation (refer to Beam Line Layout Drawing):

## **I. FRONT-END (PROCEDURE TO BE PERFORMED BY NSLS VACUUM GROUP ONLY)**

### **A. Bleed-Up**

1. Notify the Coordinator (Beeper 5824).
2. Refer to Front-End Vacuum Procedures (SLS-07.19-13-1).

### **B. Return to Operation**

1. Notify the Coordinator (Beeper 5824).
2. Refer to Front-End Vacuum Procedures (SLS-07.19-13-1).

## **II. SECTION BETWEEN ISOLATION VALVE 1A AND Be WINDOW 1A**

### **A. Bleed-Up**

1. Notify the Coordinator (Beeper 5824).
2. Close and seal the Front-End High Vacuum Valve and Valve 1A.
3. Hook up turbo pump to this section and isolate turbo.
4. Coordinator places Yellow Tags on Valve 1A and the Front-End High Vacuum Valve.
5. Slowly bleed-up with boil-off N<sub>2</sub> while Coordinator monitors the Front-End pressure.

### **B. Return to Operation**

1. Bake and pump to  $< 2 \times 10^{-9}$  Torr.
2. Notify the Coordinator (Beeper 5824).
3. Prepare for RGA Scan.\*
4. Open Valve 1A provided pressure is  $< 2 \times 10^{-9}$  Torr downstream of the valve.
5. Perform RGA Scan.\*
6. If RGA Scan or pressure reading (if no RGA Scan required) is satisfactory, Coordinator removes Yellow Tags from Valve 1A and the Front-End High Vacuum Valve.
7. Remove any unprotected turbo pump from this section or valve off the turbo pump and place a Tag on the valve. \*\*

## **III. SECTION BETWEEN Be WINDOW 1A AND HUTCH KAPTON WINDOW, MONOCHROMATOR**

### **A. Bleed-Up**

1. Notify the Coordinator (Beeper 5824).
2. Close and seal Valve 1A.
3. Coordinator places Yellow Tag on Valve 1A.
4. Slowly bleed-up through Leak Valve #5 while Coordinator monitors pressure between Valve 1A and Be Window 1A.

### **B. Return to Operation**

1. Pump to  $< 5 \times 10^{-2}$  Torr.
2. Notify the Coordinator (Beeper 5824).
3. Open Valve 1A if pressure is  $< 2 \times 10^{-9}$  Torr downstream of the valve.
4. Coordinator removes Yellow Tag from Valve 1A.

Brookhaven National Laboratory/ LIGHT SOURCES DIRECTORATE			
<b>Subject:</b>	<b>VACUUM PROCEDURES FOR BEAMLINE X-14A</b>		
<b>Number:</b>	LS-OPS-0139	<b>Revision:</b>	B
		<b>Effective:</b>	03/15/2011
			Page 2 of 4

#### IV. SECTION BETWEEN Be WINDOW 1A AND VALVE 6A, MONOCHROMATOR

##### A. Bleed-Up

1. Notify the Coordinator (Beeper 5824).
2. Close and seal Valve 1A.
3. Coordinator places Yellow Tag on Valve 1A.
4. Close Valve 6A.
5. Slowly bleed-up through Leak Valve #5 while Coordinator monitors pressure between Valve 1A and Be Window 1A.

##### B. Return to Operation

1. Pump to  $< 5 \times 10^{-2}$  Torr.
2. Open Valve 6A.
2. Notify the Coordinator (Beeper 5824).
3. Open Valve 1A if pressure is  $< 2 \times 10^{-9}$  Torr downstream of the valve.
4. Coordinator removes Yellow Tag from Valve 1A.

#### V. SECTION BETWEEN VALVE 3A AND HUTCH KAPTON WINDOW, BEAMLINE

##### A. Bleed-Up

1. Notify the Coordinator (Beeper 5824).
2. Close and seal Valve 3A.
3. Coordinator places Yellow Tag on Valve 3A.
4. Slowly bleed-up through Leak Valve #5 while Coordinator monitors pressure in the Monochromator.

##### B. Return to Operation

1. Pump to  $< 5 \times 10^{-2}$  Torr.
2. Open Valve 6A.
2. Notify the Coordinator (Beeper 5824).
3. Open Valve 3A if pressure is  $< 5 \times 10^{-2}$  Torr downstream of the valve.
4. Coordinator removes Yellow Tag from Valve 3A.

#### VI. SECTION BETWEEN VALVE 3A AND VALVE 6A, BEAMLINE

##### A. Bleed-Up

1. Notify the Coordinator (Beeper 5824).
2. Close and seal Valve 3A.
3. Close and seal Valve 6A.
3. Coordinator places Yellow Tag on Valve 3A.
4. Slowly bleed-up through Leak Valve #5 while Coordinator monitors pressure in the Monochromator.

##### B. Return to Operation

1. Pump to  $< 5 \times 10^{-2}$  Torr.
2. Open Valve 6A.
2. Notify the Coordinator (Beeper 5824).
3. Open Valve 3A if pressure is  $< 5 \times 10^{-2}$  Torr downstream of the valve.
4. Coordinator removes Yellow Tag from Valve 3A.

#### VII. SECTION BETWEEN VALVE 6A AND HUTCH KAPTON WINDOW, BEAMLINE

##### A. Bleed-Up

1. Notify the Coordinator (Beeper 5824).
2. Close and seal Valve 3A.

<b>Subject:</b>	<b>VACUUM PROCEDURES FOR BEAMLINE X-14A</b>		
<b>Number:</b>	LS-OPS-0139	<b>Revision:</b>	B
		<b>Effective:</b>	03/15/2011
			<b>Page 3 of 4</b>

3. Close and seal Valve 6A.
3. Coordinator places Yellow Tag on Valve 6A.
4. Slowly bleed-up through Leak Valve #7 while Coordinator monitors pressure in the beamline upstream Of Valve 6A.

**B. Return to Operation**

1. Pump to  $< 5 \times 10^{-2}$  Torr.
2. Open Valve 6A.
2. Notify the Coordinator (Beeper 5824).
3. Open Valve 3A if pressure is  $< 5 \times 10^{-2}$  Torr downstream of the valve.
4. Coordinator removes Yellow Tag from Valve 6A.

**\* NSLS POLICY FOR RGA SCANS (24 HOUR NOTICE REQUIRED)**

An RGA scan is required before returning to operation if there is a major change of hardware in the vacuum system, i.e. changing of samples, mirrors, windows, monochromator crystals or gratings, manipulators, detectors, etc., **with the following two exceptions:**

1. After UHV sample chambers have been bled up for replacing components, an RGA scan will not be required if the chamber pressure is returned to  $< 2 \times 10^{-9}$  Torr and the Front End pressure remains  $< 2 \times 10^{-9}$  Torr when vacuum sections upstream of the chamber are opened into the Front End.
2. If any vacuum section upstream of the bled-up section remains at a pressure of  $< 9 \times 10^{-10}$  Torr as read using a hot-filament ion gauge, when the entire beamline is opened into the Front End, and the Front End pressure does not increase, no RGA is required.

**\*\* NSLS TURBO PUMP POLICY**

An unprotected turbo pump is one not separated from the Front End by a beamline valve which automatically closes in the event of a power loss or a pressure increase at the turbo pump. **No unprotected turbo pump can share a contiguous vacuum with the Front End.**

<b>Subject:</b>	<b>VACUUM PROCEDURES FOR BEAMLINE X-14A</b>		
<b>Number:</b>	LS-OPS-0139	<b>Revision:</b>	B
		<b>Effective:</b>	03/15/2011
			<b>Page 4 of 4</b>

<b>Document Review Frequency</b>
<b>3</b> Years

Review signatures on file  
with master copy of  
controlled document

<b>LIGHT SOURCES DIRECTORATE REVISION LOG</b>		
<b>Document Number:</b>	LS-OPS-0139	
<b>Subject:</b>	VACUUM PROCEDURES FOR BEAMLINE X-14A	
<b>Rev</b>	<b>Description</b>	<b>Date</b>
B	Initial release in to the Controlled Document System.	03/15/2011