

# **Protocol for Reagent Preparation for Use in the Neutralizing Antibody Assay for HIV-1 in A3R5 Cells (December 2011)**

## **I. Introduction**

The preparation and maintenance of key reagents used for the Neutralizing Antibody Assay for HIV-1 in A3R5 Cells is crucial for obtaining accurate and reproducible results. Reagents must be created and stored as per manufacturer's guidelines and must be used within pre-established expiration dates.

## **II. Definitions**

**GM: Growth Medium**

**RPMI: Roswell Park Memorial Institute**

**FBS: Fetal Bovine Serum**

**DEAE-Dextran: Diethylaminoethyl-Dextran**

**HEPES: N-2-Hydroxyethylpiperazine-N'-2-Ethanesulfonic Acid**

**DMSO: Dimethyl Sulfoxide**

## **III. Reagents and Materials**

Recommended vendors are listed. Unless otherwise specified, products of equal or better quality than the recommended ones can be used whenever necessary.

### **Complete Growth Medium for A3R5 Assay**

RPMI-1640, with 25 mM HEPES and L-glutamine

*Vendor:* Gibco BRL Life Technologies

Sterile, store at 4°C

Fetal Bovine Serum

*Vendor:* Hyclone

Heat-inactivated 56°C for 30 minutes, 500 ml bottle, sterile. Store at -20°C. Once thawed, store at 4°C.

Gentamicin

*Vendor:* Sigma

Sterile, store at 4°C

**DEAE-Dextran, hydrochloride, average Mol. Wt. 500,000**

*Vendor:* Sigma

**ViviRen Live Cell Substrate**

*Vendor:* Promega

Store at  $-20^{\circ}\text{C}$

**DMSO**

*Vendor:* Sigma

Sterile, store at room temperature

**IV. Instrumentation**

Recommended manufacturers are listed. Unless otherwise specified, equipment of equal or better quality than the recommended ones can be used whenever necessary.

**Pipettor**

*Manufacturer:* Biohit

*Manufacturer:* Rainin

**Scale**

*Manufacturer:* Mettler

**$4^{\circ}\text{C}$  Refrigerator**

*Manufacturer:* Sci-Cool

**$-20^{\circ}\text{C}$  Freezer**

*Manufacturer:* Sci-Cool

**Low Temperature Freezer (at least  $-70^{\circ}\text{C}$ )**

*Manufacturer:* Thermo LabSystems

**V. Protocol**

**1. Complete Growth Medium for A3R5 Assay**

1.1. Complete GM for use in the A3R5 Neutralization Assay consists of RPMI-1640 medium containing  $50\ \mu\text{g/ml}$  gentamicin,  $25\text{mM}$  HEPES, and 10% heat-inactivated FBS.

1.2. To make 500 ml of Complete GM, combine and mix in a sterile bottle:

447.5 ml RPMI-1640  
50.0 ml FBS  
2.5 ml gentamicin ( $10\text{mg/ml}$ )

1.3. Store the Complete GM at  $4^{\circ}\text{C}$  for up to 2 months (or to the earliest expiration date of any one of the constituent reagents, whichever comes first). Before use in the assay, warm medium to  $20^{\circ}$ - $37^{\circ}\text{C}$ .

**2. DEAE-Dextran**

- 2.1. To prepare a 7.5 mg/ml solution, dissolve 3.75 gm of DEAE-Dextran in 500 ml of sterile water.
- 2.2. Create 10 ml aliquots in 15 ml sterile conical polypropylene tubes.
- 2.3. Store aliquots at -80°C.
- 2.4. DEAE-Dextran from some manufacturers does not have a listed expiration date. Contact the manufacturer for the stability of each DEAE-Dextran lot.

### **3. ViviRen Live Cell Substrate**

- 3.1. Add 1 ml of DMSO to a vial containing 37 mg of the ViviRen powder. The concentration of the ViviRen substrate will be 60 mM after this step.
- 3.2. Aliquot the substrate in increments of 10 µl (i.e. 10, 30, 50) because 10 µl is needed for one 96-well plate.
- 3.3. Store the ViviRen Substrate at -20°C for up to 2.5 months (or to the earliest expiration date of any one of the constituent reagents, whichever comes first).
- 3.4. Thaw aliquots immediately prior to use in the assay at room temperature.